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This is a true and accurate
copy of the original

Leuzli
Regulatory Affairs Coordinator

July 5 2023

William Steinmetz
William Steinmetz
Notary Public, State of Ohio
My Commission Expires 04-20-2027
05 July 2023

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«The Cutter Reloads with staples for stapling
and cutting devices, variations:

1. The Echelon Contour Curved Cutter Reloads,
standard
2. The Echelon Contour Curved Cutter Reloads, thick»

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DEVICE DESCRIPTION

The Echelon Contour Curved Cutter with Reload is a multifire, single-patient-use device with a curved head that cuts and staples. The device will provide ligation of colorectal structures when permanent ligation is required.

The device delivers four staggered rows of titanium staples, with a knife between the second and third row of staples and creates a 40 mm curved transection. The device is designed with a feature which prevents closing if a used reload or no reload is in the instrument. Another feature is provided to prevent firing unless the closure trigger is latched in the closed position. A retaining pin holds tissue in place and can be positioned either manually or by squeezing the closure trigger.

The instrument is shipped with a reload. The Echelon Contour Curved Cutter may be reloaded five times for a maximum of six firings per instrument during a single procedure.

Each reload includes a knife blade with two staggered rows of staples on each side, an anvil, a yellow cutting washer, a retaining pin, and a staple retainer. Reloads are available in two sizes. Refer to the table below.

Gripping Surface Technology (GST), a proprietary ridged surface, precisely compresses tissue where needed. GST precisely compresses only where it is needed to prepare the tissue for stapling compressive forces on tissue and less tissue trauma compared to traditional staplers.

3D Stapling Technology, designed with offset closure of the staple legs, creates a wider staple footprint. This distributes compression more evenly throughout the staple line.

The device design ensures parallel reduction of the jaws, which guarantees the high-quality formation of staples of the same height along the entire staple line length.

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Product Codes Table

WARNING: The blue (standard) and green (thick) Echelon Contour Curved Cutter Reloads may only be used with The Echelon Contour Curved Cutter but can be used interchangeably.

Reload Code	Tissue Thickness	Reload Color	Number of staples	Closed Staple Height
GCR40B	The Echelon Contour Curved Cutter Reloads, standard	Blue	46	1.5 mm
GCR40G	The Echelon Contour Curved Cutter Reloads, thick	Green	46	2.0 mm

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INSTRUCTIONS FOR USE

Verify compatibility of all instruments and accessories prior
to using the instrument.

1. Remove the used reload from the instrument. Grasp the top
of the reload near the grasping feature and pull out in a slight
arc motion, unsnapping the reload from the jaws. Properly
discard the used reload.

2. Thoroughly clean the instrument by completely submerging
the curved head and reload housing in a sterile solution and
swish vigorously. After swishing, visually inspect the curved
head and reload housing and remove all residual staples
and/or foreign matter from instrument prior to reloading.

WARNING: Failure to properly clean the instrument prior to
reloading could damage the instrument and compromise cut
line and staple formation of subsequent firings.

3. Determine the proper reload. Tissue thickness should be
carefully evaluated prior to using the instrument and reload.

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Refer to the Product Codes Table for proper reload selection.

WARNING: Use of an inappropriate reload may result in improper staple formation.

4. Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.

5. Examine the new reload for the presence of a staple retainer (Illustration 1-3). If the retainer is not in place, discard the reload. Do not remove the staple retainer until the reload is inserted into the device. The staple retainer cannot be removed until the reload is correctly loaded into the device. The loaded device cannot be used until the staple retainer is removed.

6. Depress the release button to ensure the instrument is in the open position. Ensure that the retaining pin actuator has been retracted. Load the device with the appropriate reload by aligning the reload with the device. The tracks on each side of the reload should be used as guides to align the reload within the jaws of the instrument (Illustration 3). When the reload (staple retainer still in place) is properly aligned, push the reload into the instrument until it is fully seated. Check that the reload is held securely within the jaws. The user may notice audible and/or tactile feedback when the Echelon Contour Curved Cutter Reload is correctly seated (Illustration 4).

WARNING: Do not remove the yellow washer (Illustration 1-2) as it is critical to the structural integrity of the reload. Removing the washer may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line.

WARNING: The staple retainer ensures proper reload alignment and protects the staple leg points during shipping and transportation. If the reload is removed from the device without the staple retainer in place, the reload cannot be used and should be discarded. Failure to do so may result in malformed staples, incomplete cut line, bleeding, leakage from the staple line, and/or difficulty removing the device. In addition, instrument damage or failure may result.

7. Remove the staple retainer by pushing on the pad labeled "Push to Remove" (Illustration 5) and discard. The loaded device cannot be used until the staple retainer is removed.

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WARNING: Do not remove the yellow washer (Illustration 1-2) as it is critical to the structural integrity of the reload. Removing the washer may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line.


Caution: After removing the staple retainer, observe the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. If colored drivers are visible, the reload may not contain staples. Attempting to close the device when the colored drivers are visible will initiate the lockout feature and prevent closure.

8. Check that the reload is held securely within the jaws. To ensure the reload is fully seated, the user may push on the blue or green cartridge (Illustration 1-1) where indicated by the arrow in Illustration 4. The user may notice audible and/or tactile feedback when the reload is correctly seated.


Caution: Attempting to close the device with an improperly seated reload will initiate the lockout feature and prevent closure.

9. The instrument is now reloaded and ready to use.

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

WARNINGS AND PRECAUTIONS

- Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- Minimally invasive instruments may vary in diameter from manufacturer to manufacturer. When minimally invasive

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instruments and accessories from different manufacturers are employed together in a procedure, verify compatibility prior to initiation of the procedure.

- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments. Ensure that electrical insulation or grounding is not compromised. Do not immerse electrosurgical instruments in liquid unless the instruments are designed and labeled to be immersed.

- Pre-operative radiotherapy may result in changes to tissue characteristics. These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected staple. Careful consideration should be given to any pre-surgical treatment the patient may have undergone, which may require alterations to surgical technique or alternative surgical procedures.

- The blue (standard) and green (thick) Echelon Contour Curved Cutter Reloads may only be used with The Echelon Contour Curved Cutter but can be used interchangeably.

- Do not fire device without a loaded reload.

- Use of an inappropriate reload may result in improper staple formation.

- Attempting to load the device without a staple retainer could result in misalignment, creating malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device. In addition, instrument damage or failure may result.

- The staple retainer ensures proper reload alignment and protects the staple leg points during shipping and transportation. If the reload is removed from the device without the staple retainer in place, the reload cannot be used and should be discarded. Failure to do so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device. In addition, instrument damage or failure may result.

- Do not remove the yellow washer as it is critical to the structural integrity of the reload. Removing the washer may

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result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line. If the reload is removed from the device without the staple retainer in place, the reload cannot be used and should be discarded. Failure to do so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device. In addition, instrument damage or failure may result.

- After removing the staple retainer, examine the surface of the new reload. The reload must be replaced with another reload if any colored drivers are visible. (If colored drivers are visible, the reload may not contain staples.) Attempting to close the device when the colored drivers are visible will initiate the lockout feature and prevent closure.

- Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.

- Ensure that unintentional tissue is not cut.

- Care should be taken not to fire this instrument when the resulting staple line will be in tension.

- The Echelon Contour Curved Cutter is designed with a feature which prevents closure if a used reload, no reload, or an incorrectly inserted reload is in the instrument. The used or misloaded reload module should then be replaced with a new reload or, if no reload was present, a reload should be loaded in the instrument. After following the above steps, any device which does not close either partially or completely should not be used.

- If the reload is removed from the device without the staple retainer in place, the reload cannot be used and should be discarded. Failure to do so may result in malformed staples, incomplete cut line, bleeding, and leakage from the staple line and/or difficulty removing the device. In addition, instrument damage or failure may result.

- Make sure tissue to be stapled is properly positioned in the jaws before stapling. Bunching, stretching or uneven loading of tissue could result in leakage, lack of hemostasis, or disruption of staple line.

- Failure to properly clean the instrument prior to reloading

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could damage the instrument; compromise cut line and staple formation of subsequent firings.

- When using other technologies in the procedure, observe the precautions suggested by the original equipment manufacturer to avoid the hazards associated with their use.

- Do not use the Echelon Contour Curved Cutter or Reloads after the expiration date or if the package is damaged. This may render the equipment inoperable or non-sterile.

- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.

- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death.

- Reuse and improper reprocessing or resterilization of single-use devices may create a risk of contamination and/or cause infection or cross-infection, including, but not limited to, the transmission of infectious diseases. Contamination may lead to injury, illness, or death.

- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination. Disposal of all medical device instruments and accessories should be performed according to local requirements and regulations.

- Attempting to close the device with an improperly seated reload will initiate the lockout feature and prevent closure.

- Discard used items containing knives or blades in "sharps" containers.

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MRI Safety Information

Non-clinical testing has demonstrated the staples made of titanium (Ti3Al2.5V) alloy are MR Conditional. A patient with the implanted staples can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 3,000 gauss/cm (30 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the staples are expected to produce a maximum temperature rise of less than 2° C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device

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extends approximately 4 mm from the staple line when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

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TRANSPORT AND STORAGE

Storage and transportation conditions:

Designation	GCR40B GCR40G
Storage Temperature	Less than 25°C
Transportation Temperature	Less than 60°C
Relative humidity	Less than 60 %
Barometric pressure	No specific conditions

Keep Dry.

Store away from direct sunlight and heat.

OPERATING ENVIRONMENT CONDITIONS

Hospital, outpatient surgical center.

Operation conditions	GCR40B GCR40G
Temperature	20 °C to 24 °C
Relative Humidity	20 % to 60 %
Barometric pressure	No specific conditions

Do not use MD after the expiration date or if the package is damaged.

HOW SUPPLIED

The Echelon Contour Curved Cutter Reloads are supplied sterile for single-patient use. Discard after use. Sales unit contains 6 primary packaging in one box with instruction for use.

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THE POSSIBILITY AND WAYS OF INTEGRATING WITH OTHER MEDICAL DEVICES

The blue (standard) and green (thick) Echelon Contour Curved Cutter Reloads may only be used with The Echelon Contour Curved Cutter but can be used interchangeably. The Echelon Contour Curved Cutter may be reloaded five times for a maximum of six firings per instrument during a single procedure.

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DL-JNJRUMedicalQA@ITS.JNJ.com

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RUSSIA:

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Instruction for Use for the medical device:

«The Cutter Reloads with staples for stapling
and cutting devices, variation: The ECHELON
Flex 35 Cutter Reloads»

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INDICATIONS

The ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip (35mm, 4 Row) and ENDOPATH ECHELON™ Vascular White Reload are intended for transection and resection of tissue and vasculature. The instruments have application in multiple open or minimally invasive general, urologic, thoracic, and pediatric surgical procedures.

CONTRAINDICATIONS

- Do not use the instrument on the aorta.
- Do not use the instrument on ischemic or necrotic tissue.
- Do not use any linear cutter on major vessels without making provision for proximal and distal control.
- Tissue thickness should be carefully evaluated before firing any stapler. Refer to the Reload Product Code Chart for the tissue compression requirements (closed staple height). If tissue cannot comfortably compress to the closed staple height, or easily compresses to less than the closed staple height, the tissue is contraindicated as it may be too thick or too thin for the staple identified in the Reload Product Code Chart.
- This instrument is not intended for use when surgical stapling is contraindicated.
- Do not use the instrument where use of buttress material is indicated.
- Do not use the instrument if upon closing, the tissue extends beyond the cut line marking, which represents the minimum cut distance, such as in creating an anastomosis.
- Do not use the instrument on the parenchyma of the liver, pancreas, kidney, or spleen, or other tissues which cannot comfortably compress to the closed staple height.

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UNDESIRABLE SIDE EFFECTS/RESIDUAL RISKS

Undesirable side effects and risks associated with surgical stapling include the potential for bleeding, tissue injury, introduction of non-sterile surfaces or pathogen transfer, inflammatory or unintended tissue reaction, hypersensitivity, property or environmental damage, electrical shock, and foreign body magnetic resonance incompatibility. Also, unintended harm, extended surgery, or altered surgical approach may result from staple line failures, inability to cut, or damaged devices.

DEVICE DESCRIPTION

The Instrument Reloads are sterile, single patient use devices that simultaneously cut and staple tissue. There are four staggered rows of staples, two on either side of the cut line. The Instrument has a staple line that is approximately 35 mm long and a cut line that is approximately 30 mm long. The shaft can rotate freely in both directions and an articulation mechanism enables the distal portion of the shaft to pivot to facilitate lateral access to the operative site.

The instrument is packaged with a primary lithium battery pack that must be installed prior to use. There are specific requirements for disposing of the battery pack.

The instrument is packaged without a reload and must be loaded prior to use. A staple retaining cap on the reload protects the staple leg points during shipping, transportation, and installation of the reload. The instrument's lockout feature is designed to prevent a used or improperly installed reload from being refired or an instrument from being fired without a reload.

Caution: Do not load the instrument more than 12 times for a maximum of 12 firings per instrument.

Reload Product Code

WARNING: The ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip may only be used with the ENDOPATH ECHELON Vascular White Reload for Advanced Placement Tip, (35mm, 4 Row).

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Reload Code	Open Staple Height	Closed Staple Height	Staple Line Length	Reload Color	Number of staples	Staple Rows
VASECR35	2.5 mm	1.0 mm	35 mm	White	36	4

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INSTRUCTIONS FOR USE

Prior to Using Reload

1 Verify that the reload size matches the instrument size to be used (e.g. use an ENDOPATH ECHELON Vascular White Re-load for Advanced Placement Tip, (35mm, 4 Row) with an ECHELON FLEX Powered Vascular Stapler instrument..

2 Verify compatibility of all instruments and accessories prior to using the instrument.

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2 Prior to inserting the reload, ensure the instrument is in the open position (Illustration 2).

3 Examine the new reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

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4 Insert the new reload by sliding it against the bottom of the channel jaw until the reload alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now reloaded and ready for use (Illustration 3).

WARNING: Failure to remove the staple retaining cap prior to position-

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ing the instrument around the tissue to be stapled will result in improperly formed staples.

Caution: After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any blue-colored drivers are visible. (If blue-colored drivers are visible, the reload may not contain staples.)

5 Refer to the ETHICON ECHELON FLEX Powered Vascular Articulating Stapler package insert for instructions for using the instrument after it is loaded.

Reloading the Instrument

1 Using sterile technique, remove the reload from the package. To avoid damage, do not flip the reload into the sterile field.

2 Prior to reloading, ensure the instrument is in the open position (Illustration 2).

3 Push upward (toward the anvil) to unsnap the reload from the channel jaw. Discard the used reload (Illustration 4).

WARNING: Prior to reloading the instrument, hold the instrument in a vertical position, with anvil and channel jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the anvil and channel jaw to clear any unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and channel jaw.

4 Examine the new reload for the presence of a staple retaining cap. If the retaining cap is not in place, discard the reload.

Caution: Tissue thickness should be carefully evaluated prior to using the instrument. Refer to the Reload Product Code Chart.

WARNING: Loading the incorrect reload for the instrument size or model may result in transecting the tissue without stapling.

5 Insert the new reload by sliding it against the bottom of the channel jaw until the reload alignment tab stops in the reload alignment slot. Snap the reload securely in place. Remove the staple retaining cap and discard. The instrument is now reloaded and ready for use (Illustration 3).

WARNING: Failure to remove the staple retaining cap prior to positioning the instrument around the tissue to be stapled will result in improperly formed staples.

Caution: After removing the staple retaining cap, observe the surface of

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strument. Refer to the Reload Product Code Chart.

- Failure to remove the staple retaining cap prior to positioning the instrument around the tissue to be stapled will result in improperly formed staples.
- After removing the staple retaining cap, observe the surface of the new reload. The reload must be replaced with another reload if any blue-colored drivers are visible. (If blue-colored drivers are visible, the reload may not contain staples.)
- Prior to reloading the instrument, hold the instrument in a vertical position, with anvil and channel jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the anvil and channel jaw to clear any unused staples from the instrument. Do not use the instrument until it has been visually inspected to confirm there are no staples on the anvil and channel jaw.
- Dispose of all opened instruments whether used or unused. This device is packaged and sterilized for single use only.
- Reuse and improper reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the country in which the user and/or patient is established.

MRI Safety Information

Non-clinical testing has demonstrated the staples made of titanium (Ti3Al2.5V) alloy are MR Conditional. A patient with the implanted staples can be scanned safely under the following conditions:

- Static magnetic field of 1.5 T and 3.0 T.
- Maximum spatial field gradient of 8,900 gauss/cm (89 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the staples are expected to

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produce a maximum temperature rise of less than 2° C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 4 mm from the staple line when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

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It is used only in medical and preventive health care institutions.

Operation conditions	VASECR35
Temperature	10 °C to 40 °C
Relative Humidity	30% – 75%
Barometric pressure	No specific conditions

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HOW SUPPLIED

The reloads are supplied sterile for single-patient use. Discard after use. Sales unit contains 12 primary packaging in one box with instruction for use.

STERILIZATION METHOD

Sterilization is made by gamma irradiation.

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



The ECHELON FLEX 35 Stapler with Advanced Placement Tip may only be used with the ENDOPATH ECHELON Vascular White Reload for Advanced Placement Tip, (35mm, 4 Row).

Do not load the instrument more than 12 times for a maximum of 12 firings per instrument.

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**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN
RUSSIA:**

Johnson & Johnson LLC,
121614, Moscow, Krylatskaya str., 17/2.
Tel. +7 (495) 580-7777, Fax: +7 (495) 580-7878
DL-JNJRUMedicalQA@ITS.JNJ.com,

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