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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SHENZHEN NEW INDUSTRIES BIOMEDICAL ENGINEERING CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized

Wei Shanshan

Signature:

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(Date: Oct. 23, 2023)

证明书网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Starter reagent kit (Starter 1, Starter 2) for the automatic analyzers MAGLUMI[®], in versions

PURPOSE

The medical device is used as an aid for in vitro diagnostics and monitoring on human whole blood, serum, plasma or urine samples and is designed to trigger the chemiluminescence reaction when used on the chemiluminescence immunoassay analyzers MAGLUMI[®].

BRIEF REVIEW

Starter 1 provides an alkaline environment, and then starter 2 starts N-(4-aminobutyl)-N-ethylisoluminol (ABEL) to produce a light signal that will be detected by the analyzer.

DEVICE COMPOSITION

Starter reagent kit (Starter 1, Starter 2) for the automatic analyzers MAGLUMI[®], in versions:

Version 1

Starter reagent kit (Starter 1, Starter 2) for the automatic analyzers MAGLUMI[®], 230 ml, as part of:

1. Starter 1 – 230 ml
2. Starter 2 – 230 ml
3. Manufacturer's quality control certificate – 1 pc.
4. Instructions for use – 1 pc.

Version 2

Starter reagent kit (Starter 1, Starter 2) for the automatic analyzers MAGLUMI[®], 1500 ml, as part of:

1. Starter 1 – 1500 ml
2. Starter 2 – 1500 ml
3. Manufacturer's quality control certificate – 1 pc.
4. Instructions for use – 1 pc.

COMPONENTS

Catalog number	Description	Specification
130299004M	1.5% NaOH solution (Starter 1)	Starter 1: 1×230 ml Starter 2: 1×230 ml
130299027M	0.18% H ₂ O ₂ solution (Starter 2)	Starter 1: 1×1500 ml Starter 2: 1×1500 ml

Devices used in combination:

Material	No. RU/REF
Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI [®] with accessories	RZN 2021/15740
In vitro Reagent Kit for the automatic analyzers MAGLUMI [®]	Individually
Reaction Modules for the automatic analyzers MAGLUMI [®]	RZN 2021/14862
Tip for automatic analyzers MAGLUMI [®]	130208000201
Reaction Cups for the automatic analyzers MAGLUMI [®]	130105000101
Luminescence measurement control (Light Check) for the automatic analyzers MAGLUMI [®]	RZN 2021/14861
Wash Concentrate for the automatic analyzers MAGLUMI [®]	RZN 2021/14863
System Tubing Cleaning Solution for the automatic analyzers MAGLUMI [®]	RZN 2021/14881

For the specific accessories required and the specifications of accessories for each model, see in the Operating Manual of the analyzer MAGLUMI[®].

Order these materials from Shenzhen New Industries Biomedical Engineering Co., Ltd (SNIBE Co., Ltd.) or from our authorized representatives in the Russian Federation.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostics.

or professional use only.

Do not use the medical device after the specified expiration date (12 months from the date of manufacture).

Follow the usual precautions necessary for handling all laboratory reagents.

Whenever replacing the starting reagents, follow the procedure described in the Analyzer Operating Manual MAGLUMI®.

- Do not shake the vials to avoid the formation of bubbles.
- Always close the starter cover after loading the reagents.
- Never leave bottles of starter reagents open.
- Do not change the bottle caps of starting reagents 1 and 2.
- Avoid ingestion, contact with skin and eyes, if this happens, please drink or rinse with plenty of water, consult a doctor if necessary.
- Safety data sheets are available on request.
- All waste should be disposed of in accordance with local guidelines.

Note: If any serious incident has occurred with respect to the product, please report it to Shenzhen New Industries Biomedical Engineering Co., Ltd. (Snibe) or our authorized representative, as well as to the competent authority of the State in which you are located.

In accordance with Regulation EC 1272/2008 (CLP), hazardous reagents are classified and labelled as indicated below:

Reagents:	Starter 1
Classification:	Skin irritation, Category 2
Signal word:	Warning:
Symbols /pictograms:	 GHS07 Exclamation Mark
Hazard statements:	H315 Causes skin irritation.
Precautionary measures:	P264: Wash hands thoroughly after handling. P280: Wear protective gloves/protective clothing. P362+P364: Remove contaminated clothing and wash before reuse. P302+P352: IN CONTACT WITH SKIN: Rinse with plenty of soap and water.
Contains: (only substances prescribed in accordance with Article 18 of Regulation EC 1272/2008).	Sodium hydroxide; caustic soda (Starter 1)

APPLICATION PROCEDURE

Place the starter 1 and starter 2 vials in the analyzer compartment "1" and "2". To avoid confusion, the starting reagents are labeled. The starters must be positioned correctly.

For more information, see refer to the Analyzer User Manual MAGLUMI®.

HEIGHT AND SIZE CHARACTERISTICS

Parameter	Starter 1 (230 ml)	Starter 2 (230 ml)	Starter 1 (1500 ml)	Starter 2 (1500 ml)
Bottle height (with lid)	97.4±0.8 mm	97.4±0.8 mm	208.5±1.5 mm	208.5±1.5 mm
Bottle length	61.5±0.5 mm	61.5±0.5 mm	150±1.2 mm	150±1.2 mm
Bottle width	61.5±0.5 mm	61.5±0.5 mm	82±1.0 mm	82±1.0 mm
Weight of the filled reagent bottle	265.1±4.7 g	259.9±4.6 g	1702.1±30.6 g	1667.4±29.9 g
Reagent volume	230±4.6 ml	230±4.6 ml	1500±50 ml	1500±50 ml
Height of the lid	18.2±0.86 mm	18.2±0.86 mm	18±0.2 mm	18±0.2 mm
Diameter of the lid	32.3±0.2 mm	32.3±0.2 mm	41±0.3 mm	41±0.3 mm

ORGANOLEPTIC AND CHEMICAL-PHYSICAL CHARACTERISTICS

Indicator	Acceptance criteria
Appearance	The packaging is clean, the label is clear, precise and clear. Inside the bottle there is a colorless liquid without sediment and flakes.
pH	Starter 1: from 12 to 13 Starter 2: from 7 to 8

FUNCTIONAL CHARACTERISTICS OF THE DEVICE

Indicator	Acceptance criteria
Background RLU	The Background RLU is in the range of 100-1200
Precision	The coefficient of variation is less than 10%.

STORAGE AND STABILITY

- Store at a temperature of 2-30 ° C: until the expiration date indicated on the package (12 months from the date of manufacture).
- Open at a temperature of 2-30°C: 45 days.
- Stability on board: 4 weeks.
- Keep away from direct sunlight.

ENVIRONMENTAL REQUIREMENTS

During use, transportation and storage, the medical device does not have a negative impact on humans and the environment, subject to the conditions of use and precautions specified in the instructions.

TRANSPORTATION

Transportation is carried out by all types of transport in covered vehicles in accordance with the rules of cargo transportation applicable to each type of transport. Transport at a temperature of 2-30 °C.

PACKAGING

Special transport packaging.

POSAL

To avoid possible environmental pollution, expired products, waste of used reagents, as well as products with visible damage to the individual packaging must be disposed of in accordance with SanPiN 2.1.3684-21. Proper destruction and disposal will support the conservation of natural resources and contribute to the protection of human health and the environment.

PACKAGING WASTE DISPOSAL

The packaging of the products is designed to minimize environmental pollution while maintaining the integrity of the product during transportation and storage. Packaging waste should be disposed of at local packaging collection points, in specialized containers located in municipalities.

MAINTENANCE AND REPAIR

Not applicable to this medical device.

WARRANTY OBLIGATIONS


SNIBE Co., Ltd. guarantees that this product meets the stated specifications, provided that the product is used by qualified personnel in accordance with the instructions in this manual and the materials are used before the expiration date indicated on the label.

Shenzhen New Industries Biomedical Engineering Co.,Ltd
No. 23, Jinxiu East Road, Pingshan District, 518122 Shenzhen,
People's Republic of China
Phone: +86-755-21536601 Fax: +86-755-28292740

Authorized
Representative for the
Russian Federation:

LLC "CLS" , 125363, MOSCOW, NOVOPOSELKOVA STREET, BUILDING 6,
BUILDING 216.1, ROOM 58
Phone: +7-495-9469597

Shanghai International Holding Corp. GmbH (Europe)
Eifestrasse 80, 20537 Hamburg, Germany
Tel: + 49-40-2513175 Fax: + 49-40-255726

SYMBOLS

	Refer to the instructions for use		Manufacturer
	Temperature range (Store at 2-30°C)		Avoid exposure to sunlight
	Catalog number		Use up to
	Medical device for in vitro diagnostics		Batch code
	European Quality Certificate		Composition
	Registration certificate		The top of the device
	Warning	-	

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