17, Krylatskaya str., bld. 2, Moscow, 121614  
Tel.: (495) 580-77-77 Fax: (495) 580-78-78

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EN

### Device Description

The TECNIS Eyhance IOL (intraocular lens), model ICB00, is an ultraviolet-light absorbing posterior chamber IOL which provides distance vision and a dysphotopsia profile comparable to a standard aspheric monofocal IOL. The Model ICB00 lens has a modified aspheric anterior surface that is designed to extend depth of focus, which improves intermediate vision as compared to a standard aspheric monofocal IOL. Accommodation will not be restored. The TECNIS Eyhance IOL, model ICB00, is designed to be positioned in the lens capsule to replace the optical function of the natural crystalline lens. The IOL contains a squared posterior optic edge that provides a 360-degree barrier. The edge of the optic has a frosted design to reduce potential edge glare effects.

### Indications For Use

The TECNIS Eyhance IOL, model ICB00, is indicated for the visual correction of aphakia in adult patients in whom a cataractous lens has been removed by extracapsular cataract extraction. The lens extends the depth of focus, which improves vision for intermediate tasks, and provide similar distance vision as compared to a standard aspheric monofocal IOL. The lens is indicated for placement in the capsular bag only.

### Precautions

1. Autorefractors may not provide optimal postoperative refraction of patients with TECNIS Eyhance IOL, model ICB00, lenses. Manual refraction with maximum plus technique is strongly recommended.
2. Recent contact lens usage may affect the patient's refraction; therefore, for patients that wear contact lenses, surgeons should establish corneal stability without contact lenses prior to determining IOL power.
3. Do not resterilize the lens. Most sterilizers are not equipped to sterilize the soft acrylic material without producing undesirable side effects.
4. Do not soak or rinse the intraocular lens with any solution other than sterile balanced salt solution or sterile normal saline.
5. Do not store the lens in direct sunlight or at a temperature greater than 113°F (45°C).

ametropia is targeted. The A-constants listed on the outer label are presented as a guideline and serve as a starting point for implant power calculations. The physician should determine preoperatively the dioptric power of the lens to be implanted.

Physicians requiring additional information on lens power calculations may contact the local Johnson & Johnson Surgical Vision, Inc. representative. Lens power calculation methods are described in the following references:

- Haeigis W. "The Haeigis formula". In: Shammas HJ, ed. Intraocular Lens Power Calculations. Thorofare, NJ: Slack; 2004; 41-57.
- Hoffer K.J. "The Hoffer Q formula: a comparison of theoretic and regression formulas". J Cataract Refract Surg. 19; 700-712 (1993). Erratum in: J Cataract Refract Surg 1994;20:677. Erratum in: J Cataract Refract Surg 2007;33:2-3.
- Holladay J.T., Prager T.C., Chandler T.Y., Musgrove K.H., Lewis J.W. and Ruiz R.S. "A three-part system for refining intraocular lens power calculations". J Cataract Refract Surg. 14;17-24 (1988).
- Retzlaff J.A., Sanders D.R. and Kraft M.D. "Development of the SRK/T intraocular lens implant power calculation formula". J Cataract Refract Surg. 16; 339-340 (1990). Erratum in: J Cataract Refract Surg. 1990;16:528.
- Olsen T. "The Olsen formula". In: Shammas HJ, ed. Intraocular Lens Power Calculations. Thorofare, NJ: Slack; 2004; 27-40.
- Canovas C., Arias P. "Customized eye models for determining optimized intraocular lenses power". Biomed. Opt. Express 2011;2:1649-1662.

### Patient Card

An implant identification card, to be supplied to the patient, is included in the package. The patient should be instructed to keep the card as a permanent record of his/her implant and to show the card to any eye care practitioner he/she may see in the future.

### Reporting

Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as lens-related and that were not previously expected in nature, severity or rate of occurrence must be reported to Johnson & Johnson Surgical Vision, Inc. This information is being requested from all surgeons in order to document product usage.

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ICB00 +8.0D 3223952001

### Patient Lens Implant Identification Label

(Affix to Patient Lens Implant I.D. Card)

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# Intraocular Lens

*Patient Lens Implant Identification Card*

Johnson & Johnson Surgical Vision, Inc.  
1700 E. St. Andrew Place  
Santa Ana, CA 92705 USA  
1-877-266-4543

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