Instructions for Use

Astra Tech Implant System Precision Drill EV

Astra Tech Implant System EV may only be used by dentists/implantologists who have received the appropriate education and training.

Contraindications

The system comprises dental implants, cover screws, semi-retained abutments, other prosthetic components, as well as surgical and prosthetic instruments. Astra Tech Implant System EV components may only be used together with original products and instruments. The Astra Tech Implant System EV assessment is designed to simplify the use of control parameters by means of markings, color coding and/or geometrical design.

Contraindications of the product

The drill is used during implant surgery according to the Timing Protocol. The drill has R4-connection to R into a contra angle. The drill is single use.

Order information and Dimensions

See current Astra Tech Implant System EV product catalog for updated and detailed order information of the various drills.

Material

Stainless Steel

Indications for Use

This drill is intended to be used for implant treatment with Astra Tech Implant System EV.

Contraindications

Common 3.0 implants are not indicated to replace a central incisor, a canine, a premolar or a molar in the maxilla, nor to replace a canine, a premolar or a molar in the mandible. There are no other absolute contraindications for treatment with Astra Tech Implant System EV implants besides the one above and for implant treatment in general.

Clinical Application

The drill is used for starting the preparation of an appropriate osteotomy for the implant.

Warnings

Surgery is an oral cavity and oral rehabilitation includes general risks for complications which should be known by the clinician. To be mentioned, but not limited to, the following:

- Components that accidentally are dropped in the patient's mouth may be swallowed or aspirated which can result in suffocation or physical injury. Care must be taken to have control over small devices.
- Allergy or hypersensitivity to chemical ingredients of material used. Patient health history must be carefully assessed.
- Failure to start and prepare implant site position and depth correctly can result in harm to nerves, blood vessels, teeth or other sensitive structures. It is strongly recommended that clinicians always undergo special training before starting up new treatment methods.

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Fracture of mandible may occur in a very resorbable mandible, during surgery or during routine oral function after surgery. Highly angled abutments on small diameter implants are not recommended for use in the maxillary region.

When treatment planning for implant with 6 mm length consider the widest possible implant, a two-stage surgical approach and spinning of implants. Closely monitor the patient for peri-implant bone loss or changes in the implant's response to percussion. If the implant shows greater than 50 % bone loss or mobility, consider possible removal of implant.

Routine implant treatment is not recommended until the end of the jaw bone growth. Implants, prosthetic component or prosthesis fracture may occur due to unfavorable loading conditions. Mastication forces and retention design have to be evaluated carefully.

Inadequate drilling, incorrect implant site preparation, threaded instrument

November, 14th October 2019

Dr. ASTRID FRIEDRICH
Manager Regulatory Affairs

**Dentsply
Sirona**
Dentsply Sirona
Manufacturing GmbH
Bismarckstrasse 4
53157 Mainz, Germany

**Dentsply
Sirona**
Implants

For instructions for use
and technical drawings refer to:
www.dentsplysirona.com

Settings or excessive insertion on torque can result in failed implant. Bone quality and quantity, local infections or disturbed healing are all factors, which may influence implant survival. Excessive heating and overcompression may be negatively affected by the following factors:

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Side Effects/Adverse Reactions

Not applicable for Precision Drill EV

Adverse Reactions

During implant treatment, possible adverse reactions may include but are not limited to:

- Allergy (e.g., anaphylactic shock)
- Uncontrolled post-operative bleeding (e.g., damage to vessels/ligaments)
- Aspirated or swallowed components
- Permanent nerve injury
- Fractured mandibles
- Implant loss
- Multi-implantitis
- Component or prosthesis failure/fracture

Step-by-step instructions

- Tear out the blister onto a sterile area.
- Secure the drill by acquiring the blister.
- Separate the drill from the blister by pulling back the top of the blister.
- Engage the drill with the contra angle. This is recommended due to its sharpness.
- Drill at a maximum speed of 1500 rpm with profuse irrigation.

NOTE: The Precision Drill EV is an extremely sharp, sharp-use only drill and shall never be handled manually once out of its package. For detailed information, see latest version of appropriate manuals.

For detailed information see latest version of appropriate procedure manuals.

Cleaning and Sterilization

Not applicable.

The product is sterile.

The product is for single use only.

Packaging and Storage

The product is sterilized by irradiation. The component must not be used beyond the expiry date. The package containing a plastic blister which is sealed by a lid. The components must be stored in their sealed package, in a dry place, at normal room temperature, i.e. 18-25 °C.

Performance Data

Not applicable.

Transportation Conditions

Product should be transported and carried carefully in a manner that corresponds to the special transport precautions for each medical device category at a temperature range from -20 to +40 °C, with a humidity of up to 85 %.

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The shelf life is set up as 5 years.

Contraindications

For full implant systems warranty information, refer to the detailed terms and conditions in the Implant System Warranty booklet.

ATTENTION: Not all products described in the IFU may be regulatory cleared/registered in all markets. See current product labeling and brochures for product assessment and representative product photos or contact the local Dentsply Sirona sales office for product availability.

ALL trademarks and company names are the property of their respective owners. To improve readability for our customers, Dentsply Sirona Implants does not use "®" or "™" in body copy. However, Dentsply Sirona Implants does not waive any right to the trademark and nothing herein shall be interpreted to the contrary.

Authorized representative in the territory of the Russian Federation:

Dentsply Sirona LLC
115432, Russia, Moscow, Andropov str., 4, bld. 8
Tel. +7 (495) 723 10 17
sales@dentimplants.ru

Valid for registration (Russia)

Surgical instruments for Astra Tech Implant System EV single use

Explanation of Symbols

Symbol	Title
	Manufacturer
	Date of manufacture
	Expiry date
	Batch Code
	Dating Number
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Consult Instructions for Use
	European conformity (CE) mark
	Prescription only U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician
	Russian certification symbol according to GOST

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Instructions for Use

Astra Tech Implant System Guiding Cylinder/Fragment Fork/Screw Extractor/ Thread Cleaner

Astra Tech Implant System may only be used by dentists/physicians who have had the appropriate education and training.

Warnings and Precautions

The system comprises dental implants, cover screws, transmucosal abutments, other prosthetic components, as well as surgical and prosthetic instruments. Astra Tech Implant System components may only be used together with original products and instruments.

Intended Use

Astra Tech Implant System comprises the required instruments/tools for use during repair and repair process in treatment of wear failure. Guiding Cylinder is used to guide the drill bit to the center of the implant. Fragment Fork is used for removing fragments. Screw Extractor is used as a next step if the fragment cannot be removed with the fragment fork. Thread Cleaner is used for cleaning the inner threads of the implant or equipment after removing fragments. The products/instruments have an R&D certification. The products/instruments are single use.

Other Information and Dimensions

See current Astra Tech Implant System Repair procedure manual and product catalog for updated and detailed order information on the various products/instruments.

Materials

Stainless steel
Titanium Carbide - Fragment Fork

Intended Use

The products are intended to be used during repair and repair process of treatment of wear failure and to remove fragments in the implant and to remove the implant.

Contraindications

Not applicable.

Indications for Use

See above - Indications for Use.

Warnings

Surgery in the oral cavity and oral rehabilitation include general risks for complications which should be known by the clinician. To be mentioned, but not limited to, the following:

- Components that accidentally are trapped in the patient's mouth may be swallowed or aspirated which can result in suffocation or physical injury. Care must be taken to have control over oral devices.
- Allergy or hypersensitivity to chemical ingredients of material used. Patient health history must be strictly assessed.

Precautions

NOTE: Products intended for manual use shall not be used together with a contra angle.

RFI Safety Information

Not applicable for the products (Guiding Cylinder/Fragment Fork/Screw Extractor/Thread Cleaner).

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- Allergy (e.g. anaphylactic shock)
- Aspirated or swallowed components

Step-by-Step Instructions

- Clean and sterilize the products/instruments before use (see section Cleaning and Sterilization).

Dr. Astrid Evers-Meier
Manager Regulatory Affairs

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Dentsply
Sirona
Implants

For instructions for use
and symbols, please refer to
the product manual.

Instructions for Use

Storage and Transport

The products should be stored in their package, in a dry place, at normal room temperature, i.e. 15-25 °C.

Indications for Use

Not applicable.

Intended Use

Products should be transported and stored carefully in a sealed box, according to the special transport provisions for each medical device's nature at a temperature range from -30 to +50 °C with a humidity of up to 90 %.

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The shelf life is set up as 10 years.

Approved System Warranty

For full implant system warranty information, refer to the detailed terms and conditions in the Astra Tech Implant System Warranty booklet.

Product Availability

ATTENTION: Not all products described in this IFU may be regularly cleared/registered/used in all markets. See current product catalog and brochures for product assortment and representative product photos or contact the local Dentsply Sirona sales office for product availability.

Trademarks

All trademarks and company names are the property of their respective owners. To improve readability for our customers, Dentsply Sirona Implants does not use "®" or "™" in every copy. However, Dentsply Sirona Implants does not waive any right to use trademark and naming herein shall be interpreted to the contrary.

Dentsply Sirona LLC
115452, Russia, Moscow, Andropovov, 18, 5th fl.
Tel: +7 (495) 725 10 87
info@dentsplysirona.com

Surgeon documents for Astra Tech Implant System EV single use.

Explanation of Symbols

Symbol	Title
	Manufacture
	Date of manufacture
	Batch Code
	Catalog Number
	Item Status
	Do not reuse
	Consult instructions for use
	European conformity (CE) mark
	Prescription only
	U.S. Federal law restricts this device to sale by or on the order of a licensed doctor or physician
	Russian certification symbol according to GOST

CE

Manufacturer:
Dentsply Implants Manufacturing GmbH
Röntgenstr. 10, 4 65927 Hanau, Germany
www.dentsplysirona.com
www.dentsplysirona.com

Subject to change without notice.

For all services are available in all countries.

Subject to change:
IFU 585100, RU Rev. 001 8 2018/04

Dentsply
Sirona
Implants

CE

Manufacturer:
Dentsply Implants Manufacturing GmbH
Röntgenstr. 10, 4 65927 Hanau, Germany
www.dentsplysirona.com
www.dentsplysirona.com

Subject to change without notice.

For all products are available in all countries.

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(Handwritten signature of Alexander V. Dudarev)

Instructions for Use

Astra Tech Implant System Fragment Drills/Trephine Drills

Berlin, October 2019

Dr. Astrid Tschöke
Maximal Regeneration Aachen

Dentsply
Sirona
Dental
Manufacturing GmbH
Röntgenstr. 10
42699 Solingen, Germany

Dentsply
Sirona
Implants

For instructions for use
and general information visit
www.dentsplysirona.com

Astra Tech Implant System may only be used by dental professionals who have had the appropriate education and training.

Intended Use of the System

The system comprises dental implants, cover screws, transmucosal abutments, other prosthetic components, as well as surgical and prosthetic instruments. Astra Tech Implant System EV components may only be used together with original products and instruments.

The Astra Tech Implant System EV assessment is designed to simplify the use of dental components by means of markings, color coding and/or geometrical design.

Description of the Product

Astra Tech Implant System comprises the required instruments/products for use during repair and rescue process in treatment of maxillofacial.

The fragment drill is:

- Used during rescue process in case of screw fracture
- Marked with diameter

The trephine drill is:

- Used during rescue process in case of removal of implant
- Marked according to implant size (EV assessment)
- Marked according to diameter (TX assessment)
- Marked with length
- Marked with depth markings

The drills are single use and have RV-connection to fit into a contra angle.

Storage, Transportation and Disinfection
See current Astra Tech Implant System Repair procedure manual and product catalog for updates and specific order information on the various drills.

Warnings

No specific warnings.

Intended Use

The products are intended to be used during repair and rescue process of treatment of maxillofacial and to remove old parts in the implant or to remove the implant.

Contraindications

No contraindications.

Contraindications

No contraindications.

Contraindications

See above - Indications of use.

For details see Repair procedure manual.

Warnings

Surgery is the oral cavity and oral rehabilitation include general risks for complications which should be known by the dentist. To be mentioned, but not limited to, the following:

- Components that accidentally are dropped in the patient's mouth may be swallowed or aspirated which can result in suffocation or physical injury. Care must be taken to have control over small devices.
- Allergy or hypersensitivity to chemical ingredients in material used. Patient health history must be carefully assessed.
- Failure to plan and prepare implant site position and depth correctly can result in harm to nerves, blood vessels, teeth or other sensitive structures. It is strongly recommended that dentists always undergo special training before starting up new treatment methods.

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Fracture of mandible may occur in a very atrophic mandible during surgery or during routine oral function after surgery.

Highly angled abutments on small diameter implants are not recommended for use in the maxilla region.

Excessive torqueing, incorrect implant site preparation, incorrect instrument settings or excessive installation torque can result in failed implant.

Not for use information

Not applicable for Fragment Drill/Trephine Drills

Manufacturer:
Dentsply Sirona Manufacturing GmbH
Röntgenstr. 10, 42699 Solingen, Germany
www.dentsplysirona.com
info@dentsplysirona.com

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Instructions for Use

Storage and Transport

The products should be stored in their original packaging in a dry place at normal room temperature, i.e. 18-25 °C.

Disinfection

Not applicable.

Intended Use of the System

Products should be transported and carried carefully in a manner that corresponds to the special transport precautions for each medical device's section at a temperature range from -20 to +40 °C, with a humidity of up to 80 %.

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The shelf life is set to 10 years.

Intended Use of the System

For all implant systems warranty information, refer to the detailed terms and conditions in the Astra Tech System Warranty booklet.

Product Availability

ATTENTION: Not all products described in this IFU may be regularly distributed/available in all markets. See current product catalog and brochures for product availability and representative product photos or contact the local Dentsply Sirona sales office for product availability.

Trademarks

All trademarks and company names are the property of their respective owners. To improve readability for our customers, Dentsply Sirona Implants does not use "®" or "TM" to denote copyright. However, Dentsply Sirona Implants does not waive any right to the trademark and nothing herein shall be interpreted to the contrary.

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Explanation of Symbols

Symbol	Title
	Manufacturer
	Date of manufacture
	Batch Code
	Catalog Number
	Non-Sterile
	Do not reuse
	Consult Instructions for Use
	European conformity (CE) mark
	Prescription only
	Revision certification symbol according to GOST

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Instructions for Use

Astra Tech Implant System Guided Surgery with Astra Tech Implant System EV Single-Use Products EV-GS EV-GS Punch/Sleeve for Drill EV-GS/Guide fixation screw

Astra Tech Implant System EV may only be used by dentists/physicians who have had the appropriate education and training.

Description of the System

The guided surgery assessment set Astra Tech Implant System EV consists of all components which are necessary to treat Astra Tech Implants EV implants using a computerized case-specific drilling guide (e.g. Stereolitho Guide) which is designed by a compatible planning software (e.g. Simplant Software). The components of a drilling guide and a planning software are given when the design features coincide with those of the guided surgery assessment set Astra Tech Implant System EV.

The Astra Tech Implant System EV guided surgery assessment set is designed to simplify the use of correct components by means of markings, color coding and/or geometrical design.

Never change, wear, scratch, or alter, at all times, take the necessary precautions to avoid contamination of components during the implant treatment.

Description of the Products

The punch is used for soft tissue penetration during guided implant surgery according to the Drilling Protocol and has a RA-connection to fit into a cone angle and/or for single use.

The sleeve fits on the corresponding Drill EV-GS and guides the drill in the guiding section (e.g. a tube) of the guide, it is for single use.

The Guide fixation screw is intended to restrain movements of the guide during surgery and is for single use.

Order Information and Dimensions

See current Guided Surgery manual/product catalog for the computer guided implant treatment with the Astra Tech Implant System EV for updated and detailed order information of the various drills.

Material

Stainless steel
Titanium

Indications for Use

The components are intended to be used for guided implant treatment with Astra Tech EV Implants using a drilling guide designed by a 3D planning software.

Contraindications

Class II, III implants are not indicated to restore a central incisor, a canine, a premolar or a molar in the maxilla, nor to replace a canine, a premolar or a molar in the mandible.

There are no other absolute contraindications to treatment with Astra Tech Implant System Implants besides the one above and for implant treatment in general.

Clinical Application

- The Punch is an optional component used to perforate the soft tissue through the guide.
- The sleeve is necessary to guide the Drill EV-GS through the guide to the appropriate position.
- The Guide fixation screws to fix and secure the Implant Site Guide against relocating.

Warnings

Surgery in the oral cavity and oral rehabilitation include general rules for complications which should be taken by the dentist. It is mentioned, but not limited to, the following:

- Components that accidentally are dropped in the patient's mouth may be swallowed or aspirated which can result in suffocation or physical injury.
- Allergy or hypersensitivity to chemical ingredients of material used. Patient health history must be carefully assessed.
- Failure to plan and prepare implant site precisely and check correctly can result in harm to nerves, blood vessels, tooth or other sensitive structures. It is strongly recommended that clinicians always undergo special training before starting up new treatment methods.

Dentistry
to Strong
Dentistry Implants
Manufacturing GmbH
Rothschießer Chaussee 4
53157 Hürth, Germany

Dentistry
to Strong
Dentistry Implants
Manufacturing GmbH
Rothschießer Chaussee 4
53157 Hürth, Germany

For instructions for use
and available accessory visit
the website www.dentistrytostrong.com

Precautions

Components and instruments for Astra Tech Implant System are not compatible with Astra Tech Implant System EV.

Position the patient so that the danger of aspiration of components is minimized. All components that are used intra-orally must be secured to prevent aspiration or swallowing.

Fracture of mandible may occur in a very atrophic mandible, during surgery or during routine and function after surgery.

Highly angled abutments on small diameter implants are not recommended for use in the upper region.

When treatment planning for implants with a long length, consider the widest possible implant, a two-stage surgical approach and spinning of implants.

Closely monitor the patient for post-implant bone loss or change in the implant's response to percussion. If the implant shows greater than 50% bone loss or mobility, consider possible removal of implant.

Routine implant treatment is not recommended until the end of age jaw bone growth.

Implant, postoperative component or prosthesis fracture may occur due to unfavorable loading conditions. Fracture forces and restoration design have to be evaluated carefully.

Improper coating, incorrect implant site preparation, unwanted instrument settings or excessive insertion torque can result in failed implant.

Bone quality and quantity, local infections or disturbed initial healing are all factors, which may influence implant survival. Successful healing and osseointegration may be negatively affected by the following factors:

- compromised general health status
- compromised local health status in the jaw
- history of therapeutic radiation in the area
- alcohol or drug abuse
- smoking

Step-by-Step Instructions

Punch

- Open the package.
- Place the punch into a sterile area.
- Secure the punch by squeezing the holder.
- Due to the sharpness of the punch, expose the shaft of the punch by bending back the top of the holder and engage it with the screw and.
- Start rotating.
- Drill at a maximum speed of 800 rpm.

Sleeve for Drill

- Place the sleeve into a sterile area.
- Assemble sleeve and drill initially to the top of the guide and check the sleeve's smooth.
- Stop the sleeve with the widest part first on the corresponding drill EV-GS.
- The sleeve guides the drill into the guide.
- Remove the sleeve after the single patient surgical treatment.

Instructions for Use

Guide Fixation Screw

- The Guide fixation screw has to be stored in the planning software in a sterile package for implant.
- Place the guide in the mouth in its correct position.
- Place the punch into a sterile area.
- Prepare the bone cavity with the drill for Guide fixation screw. See IFU: 5325 - Resorbable Products EV-GS.
- Use the Screwdriver EV to insert the screw through the guide into the bone mandible.

NOTE: Be careful not to over-tighten.

For detailed information regarding the guided surgery procedure and the design features of the guided surgery assessment set Astra Tech Implant System EV, see latest version of the Guided surgery manual.

Cleaning and Sterilization

Not applicable.
The products are sterile.
This product is for single use only.

Packaging and Storage

The products are stored by radiation. The components must not be used beyond the expiry date. The package contains a probe which is secured by a lid.
The components must be stored in their sealed package, in a dry place, at ambient room temperature, i.e. 15-25 °C.

Documentation

Not applicable.

Transportation condition

Products should be transported and stored carefully in a manner that corresponds to the special transportation conditions for each medical device's nature at a temperature range from -30 to +60 °C, with a humidity of up to 80 %.

Terms of Use

Dental clinics and laboratories.

Application Area and Operating Conditions

Stomatology, Dental implantology.

Disposal of Products

For disposal of the packaging and components comply with the currently applicable national waste disposal regulations in your country.

Shelf Life

The shelf life is set up as 5 years.

Implant System Warranty

For full implant system warranty information, refer to the detailed terms and conditions in the Implant System Warranty booklet.

Product Availability

ATTENTION: Not all products discussed in this IFU may be regulatory cleared/CE marked in all markets. See current product catalog and brochures for product descriptions and representative product photos or contact the local Dentistry to Strong sales office for product availability.

Trademarks

All trademarks and company names are the property of their respective owners. To improve readability for our customers, Dentistry to Strong Implants does not use "®" or "™" as body copy.

However, Dentistry to Strong Implants does not waive any right to the trademark and nothing herein shall be interpreted to the contrary.

Dentistry to Strong Implants is a registered trademark of Dentistry to Strong Implants.
Dentistry to Strong LLC
11602, Vienna, Austria, Andrejovskaya, 10, 11602
Tel.: +7 (495) 725 10 17
www.dentistrytostrong.com

Symbol for Registration (Russia)
Surgical Instruments for Astra Tech Implant System EV single use

Explanation of Symbols

Symbol	Title
	Manufacturer
	Date of manufacture
	Use-by date
	Batch code
	Catalogue number
	Sterilized using radiation
	Do not reuse
	Do not use if package is damaged
	Do not flame
	Consult instructions for use
	European conformity (CE) mark
	Freezer storage only
	U.S. Federal law restriction: for use by or on the order of a licensed dentist or physician
	Revision certification symbol according to GOST

Manufacturer:
Dentistry to Strong Implants Manufacturing GmbH
Rothschießer Chaussee 4, 53157 Hürth, Germany
www.dentistrytostrong.com

Subject to change without prior notice.

Use of procedure and device is at user's risk.

Printed information:
IFU 3035-1_RU Rev. 002 B, 2018-04

Dentistry
to Strong
Dentistry Implants
Manufacturing GmbH
Rothschießer Chaussee 4, 53157 Hürth, Germany
www.dentistrytostrong.com

Manufacturer:
Dentistry to Strong Implants Manufacturing GmbH
Rothschießer Chaussee 4, 53157 Hürth, Germany
www.dentistrytostrong.com

Subject to change without prior notice.

Not all products are available in all countries.

Printed information:
IFU 3035-1_RU Rev. 002 B, 2018-04

Dentistry
to Strong
Dentistry Implants
Manufacturing GmbH
Rothschießer Chaussee 4, 53157 Hürth, Germany
www.dentistrytostrong.com

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- ~~YOUTH ORGANIZATIONS~~
- ~~REGISTRATION~~

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NOTARIUS PUBLICUS I GÖTEBORG

I, the undersigned, **PER BÄCKLUND**, Notary Public of Göteborg, Sweden, do hereby certify,

that **Göran Persmalm** –
with Swedish personal No. 590729-5058,
has personally signed to the hereto attached document and

that he is duly authorized to sign alone for and on behalf of
DENTSPLY Implants Manufacturing GmbH,
with address Steinzeugstrasse 50, 68229 Mannheim, Germany,
within the limits of a Power of Attorney, dated the 21st day of March 2017,
signed by **Norbert Rabenstein**, who, according to a Certificate issued by
Walter Büttner, Notary Public of the City of Schwetzingen, Germany, dated
the 21st day of March 2017, is duly authorized to sign alone for and on behalf
of the company.

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
e-mail:
notarius.publicus@lin.kil.se

Simplicity without compromise

The design philosophy of the Astra Tech Implant System EV is based on the natural dentition utilizing a site-specific, crown-down approach supported by an intuitive surgical protocol and a simple prosthetic workflow, for increased confidence and satisfaction for all members of the treatment team.

- Versatile implant designs, including straight, conical, sloped, short, narrow and wide, using only one surgical tray
- Flexible drilling protocol allows for preferred primary stability
- Restorative components including round and triangular options supporting soft tissue sculpturing
- Unique interface with one-position-only placement for
 - Atlantis patient-specific abutments
 - Self-guiding impression components that require only one hand for precise seating

The foundation of this evolutionary step remains the proven Astra Tech Implant System BioManagement Complex, well-documented for its long-term marginal bone maintenance and esthetic results.

Astra Tech Implant System[®]

CONTENTS

Introduction – Astra Tech Implant System[®] EV

Drilling protocol and bone classification	4
Implant assortment	4
Implant size/tooth position	5
Interface connection	5
Tray concept	6
Color coding	8

Treatment planning

Pre-operative procedures	9
Clinical application	10

Implant site preparation

Surgical components and instruments overview	12
Implant packaging and handling	16
Drilling protocol	18
Expanded drilling protocol	20
Osteotomy preparation	21
Step-by-step implant placement	22
One-stage surgical protocol	26
Two-stage surgical protocol	27
Step-by-step implant placement, 6 mm	28
Torque Wrench – surgical handling	30
Torque guide	32
Cleaning and sterilization	33

This manual is designed for qualified clinicians who have undergone at least basic surgical and in-clinic implant training. Staying current on the latest trends and treatment techniques in implant dentistry through continued education is the responsibility of the clinician.

All products may not be regulatory cleared/released/licensed in all markets. Please, contact the local Dentsply Sirona sales office for current product assortment and availability.

To improve readability for our customers, Dentsply Sirona does not use "or" in body copy. However, Dentsply Sirona does not waive any right to the trademark and nothing herein shall be interpreted to the contrary.

Product illustrations are not to scale.

Drilling protocol

The density and orientation of trabeculae in spongy bone differs from patient to patient and from location to location. These variations can create differences in perceived resistance when preparing the implant site and installing the implant. The marginal cortical bone, however, is most often homogeneous in quality but may vary in thickness from case to case, which requires a mandatory adaptive preparation method.

A unique drilling protocol has been developed to allow for preferred primary implant stability. The strategy with the drilling protocol is to provide apical bone support to the implant when indicated, but relieve the apex from bone contact when this support is not indicated.

The stepped osteotomy, providing apical bone support for the implant, is indicated in soft bone situations. In other situations the apical portion of the osteotomy can be widened using the V- or X-drill. The X-drill also widening the body portion of the osteotomy and is indicated in more dense bone situations e.g. in the lower jaw.




Implant assortment

OsseoSpeed EV implants are available in a versatile range of shapes, diameters and lengths for all indications, including situations with limited space and/or bone quantity.

Specific colors have been assigned to the different implant-abutment connection sizes, which are consistently used throughout the system and identified by symbols and colors.

OsseoSpeed® Profile EV

For additional information on OsseoSpeed Profile EV, please refer to the OsseoSpeed Profile EV manual/product catalog.


Int

The
Sys
a si
an
wo

Implant size/tooth position

The design philosophy of the Astra Tech Implant System EV is based on the natural dentition utilizing a site-specific, crown-down approach supported by an intuitive surgical protocol and a simple prosthetic workflow.

Multiple considerations are required for each tooth replacement, the support needed for the final restoration, the occlusal position, soft tissue healing, and implant design and size. The illustration indicates the recommended implant sizes in relation to the natural dentition, provided there is sufficient bone volume and space in relation to adjacent dentition.

Implant-abutment interface connection

The OsseoSpeed EV implant has a unique interface with a one-position-only placement for restorative procedures and components, e.g. the Atlantis patient-specific abutments. The interface also allows for the flexibility of six-position indexing of prefabricated abutments, while index-free abutments can be seated in any rotational position.

Abutment placement option

One-position-only
Atlantis patient-specific abutments will seat in one position only.

Six positions
Indexed abutments will seat in six available positions.

Index free
Index-free abutments will be seated in any rotational position.

Tray concept



The tray layout and components are organized to support the user throughout the entire surgical procedure. The tray design eliminates the need for rubber stoppers for holding drills and instruments, which simplifies the cleaning process.

The layout is printed on the overlay, which is snapped onto the tray base. This solution offers the possibility of adapting the tray's contents according to individual preferences.

Tray logics

The color-coded, large tray has a drill marking system for ease of use and effective handling throughout the procedure, based on the following principles:


- Drills for the spongy bone preparation are color-coded white and marked with drill numbers 1-6 on the drill shaft.
- Drills for the mandatory cortical bone preparation are color-coded according to the implant and marked with either an A or B for straight implants or A/B for conical implants.
- Drills for relieving the apical bone support and widening of the osteotomy are color-coded according to the implant and marked with V or X.
- In addition to the diameter, all drills shafts are marked with a number or letter for easy identification and reference.

Overlay 1

Accommodates all implant lengths and designs for OsseoSpeed EV straight and conical implants including OsseoSpeed Profile EV for the most commonly used implant diameters.

Overlay 2


Accommodates the complete range of diameters for OsseoSpeed EV straight and conical implants, including OsseoSpeed Profile EV, lengths 2-17 mm.

Overlay 3

Supports the full implant assortment.

Small Tray EV

A flexible, compact and convenient small tray for restorative and surgical procedures and a storage option.

Clinical application

OsseoSpeed EV implants are indicated for immediate placement in extraction sites and/or in partially or completely healed alveolar ridges using a one- or two-stage surgical procedure.

Implants with the OsseoSpeed surface are especially indicated for use in soft bone situations where implants with other surface treatments may be less effective.

OsseoSpeed EV implants can be used in an immediate loading protocol. However, for single-tooth replacement in soft bone or when using a 6 mm implant, where primary implant stability may be difficult to maintain, immediate loading may not be appropriate and thus not recommended.

Best
recent
pos
imp
load
mo

Based on mechanical strength considerations, it is recommended to always place the widest implant possible for the edentulous space. This is particularly important in the posterior regions of the jaws where loading forces are high and considerable bending moments could be generated.

In all cases, it is important to consider loading conditions when determining the number and spacing of implants.

Implant shape	General clinical application	3.0	3.6	4.2	4.8	5.4
	Suitable in the majority of situations.	For replacement of maxillary lateral and mandibular central and lateral incisors when there is not enough space for a wider implant.	Used in situations with limited bone volume or space between adjacent teeth, where a 4.2 mm implant is judged to be too wide.	Suitable in the majority of situations.	Used in situations with enough bone volume.	Suitable in situations with wide ridges and large edentulous spaces.
	In situations with limited bone volume where a larger prosthetic platform is preferred.	N/A	N/A	Suitable in situations with limited bone volume where a 3.6 mm implant diameter is the choice but where a larger prosthetic platform is preferred.	Suitable in situations with bone volume where a 4.2 mm implant diameter is the choice but where a larger prosthetic platform is preferred.	N/A

Precaution: When treatment planning for implant with 6 mm length consider the widest possible implant, a two-stage surgical approach and spacing of implants. Closely monitor the patient for peri-implant bone loss or change in the implant's response to percussion. If the patient shows greater than 90% bone loss or mobility, consider possible removal of implant.

Color coding


Throughout the Astra Tech Implant System EV, markings, color coding and geometrical designs simplify the correct identification of corresponding components.

Each implant-abutment connection size is identified by a specific color, which is used consistently throughout the system. The color is applied directly to components and instruments, as well as on packaging and informational material, where appropriate.

The following components and/or packaging are color-coded:

Healing components


Packaging for components on implant-level

Impression components on implant-level

CrestSpeed Profile EV

Abutment screws for all two-piece abutments

Laboratory components on implant-level

Pre

Pre

An

Wk

Se

to

Pre

Re

to

ra

or

A

It

of

h

I

I

Pre-operative procedures

Pre-operative examination

An evaluation of the patient's general and oral health, with clinical and radiographic examinations, must be performed. Particular attention should be given to mucous membranes, jaw morphology, dental and prosthetic history, and signs of oral dysfunction.

Radiographic analysis should be used to evaluate bone topography of the residual alveolar process. The initial radiographic evaluation, together with the clinical examination, is the basis for determining whether or not a patient is a candidate for implant treatment.

If the patient is found to be suitable, a more thorough clinical examination of the area for treatment and the opposing jaw should be performed. Any local pathology in the jaws should be treated before implant placement.

Pre-operative planning

Pre-operative planning should be based on the expected restorative treatment outcome. The Astra Tech Implant System EV assortment is designed to meet the prosthetic needs for the tooth replacement planned. The prosthetic versatility in materials, designs and sizes is aligned with the implant for support of the tooth replacements in the different positions in the jaw.

To achieve the expected outcome, treatment planning should include all stages of the procedure, from healing time and components to provisional and final restorations.

Today, digital processes with CBCT scans, together with optical surface scans, can replace or complement models mounted on an articulator and provide (analog or virtual) information of the relationship between jaws and teeth. A diagnostic wax-up with the missing teeth replaced provides important information in the planning phase. Based on analysis and evaluation of the occlusal table, force distribution and preferred sites for the implants, an optimal plan can be achieved. The transparent Radiographic Implant Guides displaying implants in different magnifications are helpful in planning optimal position, direction and implant size. When working in a digital environment, the planning software provides a library of the different implants.

Simply, the computer guided implant treatment software, can be used for the Astra Tech Implant System EV to ensure accurate planning for optimized implant position and placement.


Even though the final treatment approach may be determined at the time of surgery, consider the following based on the quality of supporting bone and expected initial stability of the implant(s):

- One- or two-stage surgical procedure
- Immediate or early loading protocol
- Expected healing time before loading

When the prerequisites for immediate loading cannot be met, an early loading protocol (at least six weeks healing period) may be considered.

In all situations, bone quality and quantity, primary stability achieved, design of restoration, and loading conditions should be carefully examined and assessed by the clinician when determining time to loading of implants for each individual case.

Before treatment begins, the patient should be informed about the results of the pre-operative examination and given a clear explanation of what is entailed by the planned treatment, including the expected outcome, maintenance requirements and risks involved.


Surgical components and instruments overview

Implant sites are prepared in a step-by-step procedure using different diameter drills, instruments and verification tools, ensuring an efficient and atraumatic preparation. All drilling in the bone should be performed at a maximum of 1500 rpm using profuse external irrigation with a saline solution. An intermittent drilling technique will help prevent heating of the bone and create a pumping effect for efficient removal of bone tissue.

Astra Tech Implant System EV drills:

- Excellent cutting properties
- Laser-etched depth indication lines
- Sterile packaging
- Multiple-use* with option for single-use
- In addition to the diameter, all drills shafts are marked with a number/letter for easy identification and reference
- Color-coded

Packaging

- Open the package.
- Pour the blister onto a sterile area.
- Secure the drill by squeezing the blister.
- Expose the drill shaft by bending back the top of the blister.

Pick-up

- Engage the drill with the control ring.

Guide Drill EV / Precision Drill EV


Used for marking and creating a starting point.

Note: The Precision Drill EV is an extremely sharp, single-use-only drill and should never be handled manually once out of its package.

* All drills except Precision Drill EV can be used for approximately ten cases. They should be carefully cleaned and sterilized after each surgery and replaced as soon as a decrease in their cutting efficiency is observed.

*Patent pending

Spongy bone preparation

Twist Drill Used for initial preparation and evaluation of the bone (1).

Step Drill Used for site preparation, resulting in a stepped osteotomy with apical bone support for the implant. The drills (2)-(6) provide guidance during the drilling process.

- Color: white
- Markings: diameter and drill number
- Lengths available in short (8-15 mm) and long (16-17 mm)

Cortical bone preparation - straight implants

Mandatory preparation of the cortical layer to reduce pressure in the bone around the implant neck.

- Color: corresponds to implant
- Markings: diameter and drill letter
- (A) - thin cortical bone < 2 mm
- (B) - thick cortical bone > 2 mm
- Length: one option

Note: There are separate cortical drills specific for the 5 mm implant.

Cortical bone preparation - conical implants

Mandatory preparation of the cortical layer to reduce pressure in the bone around the implant neck.


Apical border of the indication line indicates the minimum depth needed to fit the implant.

(A) In thin cortex, drill to the apical border of the depth indication line.




(B) In thick cortex, drill to the full depth of the depth indication line.

Make sure enough depth is provided for the entire implant.

- Color: corresponds to implant
- Markings: diameter and drill letter
- Length: one option

Surgical components and instruments overview

Alternative spongy bone preparation

V-Twist Drill

- extra apical preparation

Following the opening of the marginal cortical layer with apical drill ① or conical drill ②, the ③ drill is used to relieve the apical bone support when this support is not indicated.

- Color: corresponds to implant. Note: for conical implants, this color refers to the implant body diameter.
- Markings: diameter and ③
- Length: available in short (6-13 mm) and long (6-17 mm).

X-Stop Drill

- extra body preparation

Following the opening of the marginal cortical layer with apical drill ① or conical drill ②, the ③ drill, in addition to relieving the apical bone support, is used to widen the body portion of the osteotomy in situations with more dense bone e.g. the lower jaw.

- Color: corresponds to implant. Note: for conical implants, this color refers to the implant body diameter.
- Markings: diameter and ③
- Length: available in short (6-13 mm) and long (6-17 mm).

Direction Indicator EV

Used for visualizing the position and direction of the prepared osteotomy.

The narrow end is used after drill ① and the wider end is used after drill ③.

A laser marking indicates the 8 mm depth. The range (collar) of the instrument indicates the smallest (3.0) and greatest (5.4) implant platform size.

Direction Indicator EV is equipped with a hole for attaching a safety thread.

Implant size preparation

Measuring osteotomy depth

Carefully measuring the depth of the osteotomy. Use the same clinical reference point for the depth as for the planned implant position.

The depth should allow the implant to be level with or slightly submerged in relation to the adjacent marginal bone.

Drilling depth

The drilling depth is measured from the widest part of the drill up to the indication line. Lengths available in short (6-13 mm) and long (6-17 mm).

Additional tip depth is maximum 10 mm regardless of the drill diameter.

Implant Driver Extender EV-GS

Used for extending the length of a drill. Implant Driver EV or Implant Driver EV-GS.

Note: Implant Driver Extender EV-GS should only be used together with drills. Implant Driver EV or Implant Driver EV-GS. Ensure sufficient irrigation when using the extender.

Implant Depth Gauge EV

Used for measuring the depth of the implant site.

• Markings: corresponds to the implant lengths 6-17 mm.

The other end of the gauge can be used as a measuring probe.

• Markings: 0-25 mm

Note: The depth gauge is designed with a waist to facilitate the identification of the 13-15 mm markings.


Implant Driver EV

Used for picking up and inserting the implant in the prepared osteotomy. The Implant Driver EV can engage into any of the six positions of the implant. The body of the driver has a hexagonal geometry with splines.

Note: To facilitate optimal placement of pre-designed abutments, align one of the splines buccally.

The reference point (0°) of the depth markings is the intended bone level, i.e. the lowest point of the bevel.

- Color: corresponds to implant
- Markings: depth and diameter
- Length: two options

Implant packaging and handling

The packaging includes a 2D barcode for simplified inventory control

Packaging

- Packages contain protective blisters that hold the components.

Labels

- Three peel-off labels are provided for the patient's treatment record and for communicating with the restorative team.

Implant container

- Open the blister package.
- Pour the sterile inner container onto a sterile area.
- Remove the cap from the container using a twisting motion to expose the top of the implant.

Note: Marked with implant size and length

Implant pick up

- Attach the appropriate Implant Driver EV to the control angle or Surgical Driver Handle.

Implant pick up


- Make sure that the implant driver is fully seated into the implant (a).
- Press downwards to activate the carrying function before picking up the implant (b).
- When picking up the implant from the inner container, do not use excessive pressure.

Attach the Implant Driver EV


- Attach the Implant Driver EV by pressing it firmly into the Surgical Driver Handle EV. The driver is correctly seated when the outer conical marking is just in contact with the handle.

Manual implant pick up

- Use the Implant Driver EV together with Surgical Driver Handle EV to pick up the implant.

Drilling protocol for OsseoSpeed® EV - straight


Note: There are separate cortical drills specific for the 6 mm implant.

Drilling protocol for OsseoSpeed® EV - conical

Conical implants are available for situations with limited bone volume where a diameter 3.6 or 4.2 mm implant body is the choice, but where a larger prosthetic platform is preferred.


Note: If an additional spongy bone preparation with V or X drill is performed, use the color that refers to the implant body diameter.


Expanded drilling protocol for OsseoSpeed® EV – straight OsseoSpeed® EV – conical

In sites with compromised bone, the drilling protocol can be expanded to provide further guidance with more drilling steps. This approach is particularly useful in situations with an extremely narrow alveolar ridge (knife-edge).

OsseoSpeed® EV – straight


OsseoSpeed® EV – conical

Note: for conical implants, color refers to the implant body diameter

Osteotomy preparation


Mandatory

The spongius bone preparation results in an under-preparation compared to the implant diameter of 0.5 mm at the body portion. Consequent under-preparation at the apex is, from the widest the apex begins (towards the tip of the implant), 0.3, 0.4 and 0 mm respectively.

After completion of the spongius bone preparation, mandatory use of drill bit (A) or (B) will result in an under-preparation at the margin, of 0.5 mm or 0 mm respectively compared to the implant's nominal diameter.

Alternative*

The suggested osteotomy, which provides optimal bone support, is indicated in both bone situations. In other situations the osteotomy can be achieved using (A) or the (B) drill.

*Note: For cases where using the (A) drill results in an apical under-preparation, i.e. where the apex begins towards the tip of the implant, of 0.5 mm or 0 mm, respectively at the apex portion.


*Note: For edge and small under-ty use the (B) drill. This results in an under-preparation at the body portion of 0.25 mm. At the apical portion the under-preparation is 0.5 mm or 0 mm immediately from where the apex begins towards the tip of the implant.

*Note: BSA never allowed in soft bone situations.

Step-by-step implant placement: OsseoSpeed® EV

Below is a step-by-step protocol for the preparation and installation of an OsseoSpeed EV 4.2 S, 13 mm.


Note: All drilling should be performed at a maximum speed of 1500 rpm with profuse irrigation.

Twist Drill EV 1.9 ①

- An apical stop can be useful during the drilling procedure. If preferred, Guide Drill EV or Precision Drill EV can be used to create a pilot hole prior to using the Twist Drill ①.
- Drill in the planned direction to the appropriate depth.
- The drilling will provide valuable information about the cortical and spongy bone.
- Insert the smaller end of the Direction Indicator EV into the site to visualize/verify the direction.

Note: The depth should allow the implant to be level with or slightly submerged in relation to the adjacent marginal bone.

Step Drill EV 2.5/3.1 ③

- Drill in the planned direction to the appropriate depth.
- Insert the larger end of the Direction Indicator EV into the site to visualize/verify the direction.



Step Drill EV 3.1/3.7 ④

- Drill the implant site to the appropriate depth.

Cortical Drill EV

- Finalize the osteotomy by drilling to the full depth as indicated by the marked line. Use the cortical drill based on the cortical bone thickness:

Additional spongy bone preparation

- When deemed necessary, the preparation can be performed with either of the following drills:

V Twist Drill EV

Following the opening of the marginal cortical layer with cortical drill (A) or (B), the (V) drill can be used to relieve the apical bone support when this support is not indicated.

X Step Drill EV

Following the opening of the marginal cortical layer with cortical drill (B), the (X) drill, in addition to relieving the apical bone support, is used to widen the body portion of the osteotomy. In situations with more dense bone e.g. the lower jaw.

Measuring the osteotomy depth

- Carefully measure the depth of the osteotomy. Use the same clinical reference point for the depth as for the planned implant position.
- Verify the implant site depth after drilling using the Implant Depth Gauge EV.

Note: The depth should allow the implant to be level with or slightly submerged in relation to the adjacent marginal bone.

Step-by-step implant placement: OsseoSpeed® EV

Implant pick up

- Attach the appropriate Implant Driver EV to the contra angle.
- Make sure that the Implant driver is fully seated into the Implant (a).
- Press downwards to activate carrying function before picking up the Implant (b). Do not use excessive pressure.


Machine implant installation

- Install the Implant with the contra angle at low speed (25 rpm) under profuse irrigation and the maximum torque set to 45 Ncm.
- Allow the implant to work its way into the osteotomy. Avoid applying unnecessary pressure.

Note: Do not exceed 45 Ncm when installing the implant. If not completely seated before reaching 45 Ncm, reverse/ remove the implant and widen the osteotomy appropriately (see Additional osteotomy preparation).

It is recommended to have a titanium forceps available in case the Implant driver does not provide sufficient carrying function during the removal procedure.

Manual implant installation

- Attach the Implant Driver EV by pressing it firmly into the Surgical Driver Handle EV.
- The driver is correctly seated when the color coded marking is just in contact with the handle.
- Pick up and install the implant.

Attach

- Attach the Implant Driver EV and Surgical Driver Handle into the wrench until there is an audible click.


Positioning

- Position one of the diodes on the implant driver exactly to facilitate optimal placement of pre-designed abutments using the Torque Wrench EV together with Surgical Driver Handle.
- Release the implant driver by lifting it gently from the implant.

Finalizing implant installation

- Continue the procedure for finalizing the installation of the implant according to a one- or two-stage approach.

The reference point ("0") for the depth markings is the lowest point of the handle.

One-stage surgical protocol


HealDesign EV

- Round shapes are indicated for all positions in the mouth
- Triangular design for anterior implant sites and to mimic the shape of incisors and canines
- Heights and diameters are designed to correspond with final abutments to provide soft tissue sculpturing


Installation

- Select a HealDesign EV to develop the desired soft tissue anatomy.
- Pick up and install the sterile HealDesign EV directly from the blister package, using the Hex Screwdriver EV.

HealDesign EV, round

- Manually secure the healing abutment using light finger force (5-10 Ncm).

HealDesign EV, tri

- When using the two-piece triangle or HealDesign EV, first seat it with a Hex Screwdriver EV.
- Rotate the abutment sleeve until the desired index position is reached.
- Manually tighten the healing abutment using light finger force (5-10 Ncm).

T

O

■

Two-stage surgical protocol


Cover Screw EV

- One option for each interface connection

Cover Screw EV

- Insert the Cover Screw EV using the Hex Screwdriver EV.
- Tighten with light finger force (5-10 Ncm).
- Reposition the mucoperiosteal flaps carefully and suture tightly together.

Abutment installation

- After an appropriate healing phase, expose and reseat the Cover Screw EV with a Hex Driver EV.
- Install the selected healing abutment or restorative component/final abutment.




Note: When removing a two-piece component, keep the sleeve and screw assembled.

Step-by-step implant placement: 4.8



OsseoSpeed® EV, 6 mm

The drilling protocol for 6 mm implants is the same as for other implant lengths.


Direction Indicator EV

- Drill in the planned direction to the appropriate depth.
- When the drill (1) has been used, insert the smaller end of the Direction Indicator EV into the osteotomy to visualize/verify the direction.
- See the 6 mm markings on Direction Indicator EV.
- After drill (3), use the larger end of the indicator to verify the direction.
- Drill the implant site to the appropriate depth with drill (5).

6 mm – Cortical Drill EV

- Finalize the osteotomy by drilling to the marginal end of the depth indication line. Based on cortical bone thickness, use either the cortical drill (A) or (B) specific for the 6 mm implant.

Additional spongy bone preparation


- When deemed necessary, the preparation can be performed with either of the following drills:

Twist Drill EV

Following the opening of the marginal cortical layer with cortical drill (A) or (B), the Twist Drill EV can be used to relieve the apical bone support when this support is not indicated.


Step Drill EV

Following the opening of the marginal cortical layer with cortical drill (A) or (B), the Step Drill EV is used to widen the body portion of the osteotomy in situations with more dense bone e.g. the lower jaw.

Implant Depth Gauge EV

- After drilling, use the Implant Depth Gauge EV to verify the depth of the osteotomy.

Installation


- Insert the implant with the correct angle at low speed (25 rpm) and the maximum torque set to 45 Ncm.

Torque Wrench EV – surgical handling

A torque wrench together with the Surgical Driver Handle are used for implant installation and adjustment.

When used together with the Restorative Driver Handle the torque wrench can also be used for tightening abutment screws and bridge screws.

Assembly

- Assemble the head of the wrench and the body by pushing the components together and turning until there is an audible click.

Attach

- Attach the Implant Driver EV by pressing it firmly into the Surgical Driver Handle EV.
- The driver is correctly seated when the center contact marking is just in contact with the handle.
- Pick up and install the implant.
- Insert the driver handle and implant driver into the wrench until there is an audible click.


Handling

- Use a finger on the top of the driver handle to keep it steady and in place. Then gently pull the arm of the torque wrench in the direction of the arrow until the desired torque is achieved.
- Note:** The arm of the torque wrench must not go beyond the end of the scale, as this could result in inaccurate torque readings.
- Note:** The arrow on the head of the wrench shows the direction in which the wrench is functioning.

Disassembly

- Remove the driver handle from the wrench.
- Remove the head by pressing a finger into the recess (1) and gently pulling the head (2).

Cleaning and drying


- The three separated parts are now ready for cleaning using water and a brush. Let the parts dry.

Sterilization

- Follow the manufacturer's instructions for use.

Torque guide for Astra Tech Implant System® EV

Installation procedures	Recommended torque
Implant installation	44.5 Ncm 
Healing components	Manual/light finger force (5-10 Ncm)
Temporary restorations on all levels	15 Ncm 
Final restorations on implant level	25 Ncm 
Final restorations on abutment level	15 Ncm 

C

Cleaning and sterilization

All drills, except the single use Precision Drill EV, can be used approximately ten times. If drills are not reused, dispose them in a sharps container immediately after the implant procedure is completed.

Note: Single use products should not be reused.

Remove residual tissue or bone debris by immersing the used products in lukewarm water (40°C/104°F). Do not use fixation agents or hot water as this could influence subsequent cleaning results. Products should be kept in a wet environment until the next step is initiated. For Direct Driver EV Ø 3.3, Ø 4 and Ball Abutment Driver EV storage in a wet environment is mandatory.


If cleaning is delayed more than 120 minutes, place the devices in a bath of a cleaning and disinfection solution to avoid drying of soil and/or debris, blood and other contaminations.

Preparation for cleaning


Disassembly is required for the following products:

- Large Tray EV and Small Tray EV
- Impression components (plate-ups/transfers)

Large Tray EV

Small Tray EV

Manual procedure

Apply detergent, Neodisher MediClean-Forte (Dr. Weigert, Hamburg) or similar solution to all surfaces. Scrub the outer and, if applicable, the inner side of the product with a soft bristled nylon brush until all visible soil and/or debris are removed. Flush the inner channels/lumen with cleaning solution using an irrigation needle connected to a syringe. Check channels/lumen for residual soil and/or debris. Run the products in an ultrasonic bath with cleaning solution for minimum ten minutes. Drills and trays excluded. Rinse under clear running water until all trace of cleaning solution is removed. Flush the inner channels/lumen with water using an irrigation needle.

Prepare a bath with a disinfection solution, D122 Instrument disinfection (DORR SYSTEM-HYGIENE) or similar, according to the detergent manufacturer's instructions. Immerse the products completely for the time specified by the manufacturer. Flush the internal channels/lumen using an irrigation needle for a minimum of 3 times. Rinse under clear running water until all trace of disinfection solution is removed. Flush the inner channels/lumen with water using an irrigation needle.

Dry the products using medical compressed air and clean lint-free single-use wipes.

Automated procedure

Place instruments in a washer-disinfector, Vario TD or similar, according to recommendations from the supplier. Example of Vario TD washing programs:

- Pre-wash, 20°C
- Cleaning with detergent, Neodisher MediClean-Forte (Dr. Weigert, Hamburg) or similar solution 45-55°C
- Neutralization
- Intermediate rinse
- Disinfection, >90°C (preferable 93°C), 5 min
- Drying

Inspection and function testing

Drills shall be replaced as soon as their cutting capability diminishes. Discard blunt or damaged products.

Packaging pre-sterilization

Thoroughly dry everything prior to the sterilization process to prevent the risk of corrosion. Assemble the tray and re-position the drills and instruments using drill/letter numbers, where applicable. It is recommended to wrap the instruments and tray according to the sterilization wrap manufacturer's instructions. It is recommended to place the abutments, screws, and applicable products in a sterilization bag.

Note: For US: Use FDA cleared sterilization bag and 15 minutes dry time at the end of the steam sterilization cycle.

Sterilization

Steam sterilization with a pre-vacuum cycle (34°C/275°F for 3 minutes).

Sterilization procedure for zirconia products

The products should not be sterilized in a steam autoclave. The process can affect the mechanical properties of the material.

Statement Cleaning and sterilization of Astra Tech Implant System® EV products

The cleaning and sterilization instructions for Astra Tech Implant System EV assortment has been developed and validated by Dentsply Sirona. The instructions have been developed in accordance with the standards stated, please see below.

Both the VarioTD program (recommended for automated reprocessing) and the Neodisher Mediclean Forte (Dr. Weigert) detergent are recommendations and can be substituted with similar programs and detergents. For more information, please see <http://www.miele-professional.com> and/or www.drweigert.com.

Large Tray EV and Small Tray EV are made of PPSU (Polyphenylsulfone) material which may be sensitive for some chemicals containing acetate e.g. ethyl acetate. Consult your detergent manufacturer for compatibility of used detergent with PPSU if Neodisher Mediclean Forte is not used.

For ZrDesign abutment Liquid Chemical Sterilization/ High Level Disinfection is recommended.

Note: For US: Dry heat (160°C/320°F for 4 hours) is the recommended procedure.

For Atlantis abutment in zirconia Dry heat (160°C/320°F for 4 hours) is the recommended procedure.

Storage

The products should be stored, in their package, in a dry place, at normal temperature (18-25°C/64-77°F). Use the sterilized components within the stated time period from the sterile bag manufacturer.

Note: For maintenance and cleaning of Contra Angles and Torque Wrench EV, follow the manufacturer's instructions.

If alternative procedures are used it is the responsibility of the user to ensure that the cleaning and sterilization procedure chosen achieves the desired results.

- ANSI/ASME A17.9-2010 & A17.310 Comprehensive guide to steam sterilization and health care facilities in health care facilities.
- ANSI/AAMI ST89-2004/ISO 13088 Sterilization of medical devices - Information to be provided by the manufacturer for the processing of re-usable medical devices.
- AAMI TIR122:2010 Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers.
- EN ISO 17664:2004 Sterilization of medical devices - Information to be provided by the manufacturer for the processing of reusable medical devices.
- EN ISO 15883-2:2006 Washer-disinfectors - Part 2: General requirements, terms and definitions and tests.
- EN ISO 15883-3:2006 Washer-disinfectors - Part 3: Requirements and tests for washer-disinfectors employing thermal disinfection for surgical instruments, ancillary equipment, dental chairs, rooters, chisels, glassware, etc.
- SORTS 8283-5:2005 Washer-disinfectors - Part 5: Testable and methods for demonstrating cleaning efficacy.
- EN ISO 17663:2006 Sterilization of health care products - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.

About Dentsply Sirona Implants

Dentsply Sirona implants offers comprehensive solutions for all phases of implant therapy, including Ankylos, Astra Tech Implant System® and Xive® Implant Inset, digital technologies, such as Atlantis® patient-specific solutions and Simplant® guided surgery, Symbrics® regenerative solutions, and professional and business development programs, such as STEP-PS®. Dentsply Sirona implants creates value for dental professionals and allows for predictable and lasting implant treatment outcomes, resulting in enhanced quality of life for patients.

About Dentsply Sirona

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies, with a 150-year history of innovation and service to the dental industry and patients worldwide. Dentsply Sirona develops, manufactures, and markets a comprehensive solutions offering including dental and oral health products as well as other consumable medical devices under a strong portfolio of world class brands. As The Dental Solutions Company™, Dentsply Sirona's products provide innovative, high-quality and effective solutions to advance patient care and deliver better, safer and faster dentistry. Dentsply Sirona's global headquarters is located in York, Pennsylvania, and the international headquarters is based in Salzburg, Austria. The company's shares are listed in the United States on NASDAQ under the symbol XRAY.

Visit www.dentsplysirona.com for more information about Dentsply Sirona and its products.

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- ANSI/ASME STS-2010 & ASQ2010 Comprehensive guide to implant sterilization and quality assurance in health care facilities
- ANSI/ASME STS-2004/2010 Sterilization of medical devices. Information to be provided by the manufacturer for the sterilization of reusable medical devices
- AAMI TR-122010 Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities. A guide for medical device manufacturers
- EN ISO 17664:2004 Sterilization of medical devices - Information to be provided by the manufacturer for the processing of reusable medical devices
- EN ISO 13684-1:2005 Washer-disinfectors - Part 1: General requirements, terms and definitions and tests
- EN ISO 15695-2:2005 Washer-disinfectors - Part 2: Requirements and tests for washer-disinfectors equipped with ultrasonic and/or chemical disinfection systems, including equipment, tests, data, records, validation, etc.
- ISO 15855-2:2005 Washer-disinfectors - Part 2: Test rules and methods for determining cleaning efficacy
- EN ISO 17365:2006 Sterilization of health care products - Metal heat - Part 1: Requirements for the development, validation and control of a sterilization process for medical devices

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NOTARIUS PUBLICUS I GÖTEBORG

I, the undersigned, PER BÄCKLUND, Notary Public of Göteborg, Sweden,
do hereby certify:

that Göran Forsmalm –
with Swedish personal No. 590729-5058,
has personally signed to the hereto attached document, and



that he is duly authorized to sign alone for and on behalf of
DENTSPLY Implants Manufacturing GmbH,
with address Steinzeugstrasse 50, 68229 Mannheim, Germany,
within the limits of a Power of Attorney, dated the 21st day of March 2017,
signed by Norbert Rabenstein, who, according to a Certificate issued by
Walter Büttner, Notary Public of the City of Schwetzingen, Germany, dated
the 21st day of March 2017, is duly authorized to sign alone for and on behalf
of the company.

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DENTSPLY Implants
Manufacturing GmbH

Göran Forsmalm
Malmö 2018-06-18

Göran Forsmalm
Manager Regulatory Affairs

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- 1) EN ISO 17664 Sterilization of medical devices - Information to be provided by the manufacturer for the preparation of re-sterilizable medical devices (ISO 17664)
- 2) RQ Guideline 2012 Hygiene requirements for the reprocessing of medical devices (German Federal Health Gazette 2012; 56: 504-1310)
- 3) Guidelines of the DGKH (German Association for Hospital Hygiene), DGHV (German Association for Sterile Services) and AIG (Working Party Instrument Preparation) for the validation and routine monitoring of mechanical cleaning and thermal disinfection processes for medical devices, 4th edition 2014
- 4) Guidance for industry and FDA Staff - Reprocessing of Medical Devices in Health Care Settings: Validation Methods and Labeling, March 12, 2015
- 5) AAMI TIR12:2010, Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers
- 6) EN ISO 15883-1 Cleaning and disinfection appliances - Part 1: General requirements, definitions and test methods (ISO 15883-1)

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