To Whom it may concern

Ref.: Russian label on products to be distributed in Russia Federation


We, DENTSPLY Implants Manufacturing GmbH, hereby declare that products intended for distribution in the Russian Federation will be accompanied by an additional label available in Russian, as shown below.

Please note that following global product labels contain all product relevant information and standardized symbols according to ISO 15223, which are described in the Instructions for Use. This additional label in Russian was designed to do not duplicate information due to space restriction considering the size of label, and to indicate main and specific information for the country as the Authorized Representative and others shown in the Appendix.

Additionally, we would like to clarify that process of Country Specific Labeling is currently in implementation phase, thus the final labeling solution might slightly differ from the pictures submitted. We emphasize that in case of any change, it will not be content related.

The overview of final product to be delivered to the customer is the product in its sealed carton box provided in a sealed plastic bag, bearing the Instructions for Use and the Russian label that will be either printed on a paper placed in the plastic bag or affixed to the plastic bag as sticker. The global product label will not be overlapped.

Sincerely,

Dr. Aline Frare Mocruha
Regulatory Affairs Manager
May 27th 2020

Dentsply Implants
Manufacturing GmbH
Rodenbacher Chaussee 4
63457 Hanau, Germany

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


**Dentsply
Sirona**
Implants

Astra Tech Implant System® EV, Astra Tech Implant EV

Dr. Aline Frare Mocruha
Regulatory Affairs Manager

For Instructions for Use
and symbols glossary refer to
ifu.dentsplysirona.com

CAUTION: US federal law restricts this product to sale by or on order of a dentist or a physician.

Astra Tech Implant System EV may only be used by dentists/physicians who have had the appropriate education and training.

Description of the System

The Astra Tech Implant System EV is designed for replacing one or more missing teeth in patients to restore chewing function. The system comprises dental implants, cover screws, transmucosal abutments, other prosthetic components, as well as surgical and prosthetic instruments.

Astra Tech Implant System EV components may only be used together with original products and instruments.

The Astra Tech Implant System EV assortment is designed to simplify the use of correct components by means of markings, color coding and/or geometrical design.

Description of the Product

Astra Tech Implant EV is a screw-shaped dental implant with a defined surface roughness, achieved through a patented grit blasting technique, followed by a process involving a treatment in diluted hydrofluoric acid. The external coronal portion has minute threading – MicroThread.

There is a conical connection between the implants and the abutments giving a tight and stable connection – Conical Seal Design.

Astra Tech Implant EV has an inner geometry with a conical and indexing configuration.

Different profiles, diameters and lengths of the implant are available providing versatility according to the thickness and shape of the bone.

Order Information and Dimensions

See current Astra Tech Implant System EV product catalog for updated and detailed order information of the various implant products.

Material

Titanium – Grade 4

Indications for Use

The implants are intended for both one- and two-stage surgical procedures in the following situations and with the following clinical protocols:

- Replacing missing teeth in single or multiple unit applications in the mandible or maxilla.
- Immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- Especially indicated for use in soft bone applications where implants with other implant surface treatments may be less effective.
- Immediate and early loading for all indications, except in single tooth situations on implant shorter than 8 mm or in soft bone (type IV) where implant stability may be difficult to obtain and immediate loading may not be appropriate.

The intended use for Astra Tech Implant EV 3.0 S is limited to placement of maxillary lateral incisors and mandibular incisors.

Contraindications

Diameter 3.0 implants are not indicated to replace a central incisor, a canine, a premolar or a molar in the maxilla, nor to replace a canine, a premolar or a molar in the mandible.

Allergy or hypersensitivity to commercially pure titanium.

Clinical Application

See above – Indications for use.

Warnings

Surgery in the oral cavity and oral rehabilitation include general risks for complications which should be known by the clinician. To be mentioned, but not limited to, the following:

- Components that accidentally are dropped in the patient's mouth may be swallowed or aspirated which can result in suffocation or physical injury. Care must be taken to have control over small devices.
- Allergy or hypersensitivity to chemical ingredients of material used. Patient health history must be carefully assessed.

- Failure to plan and prepare implant site position and depth correctly can result in harm to nerves, blood vessels, teeth or other sensitive structures. It is strongly recommended that clinicians always undergo special training before starting up new treatment methods.

Inform us and your competent authority in the case of notice of life-threatening incidents or a severe deterioration in health status related to one of our products.

Precautions

- Components and instruments for Astra Tech Implant System are not compatible with Astra Tech Implant System EV.
- Fracture of mandible may occur in a very atrophic mandible, during surgery or during routine oral function after surgery.
- Highly angled abutments on small diameter implants are not recommended for use in the molar region.
- When treatment planning for implant with 6 mm length consider the widest possible implant, a two-stage surgical approach and splinting of implants. Closely monitor the patient for peri-implant bone loss or change in the implants' response to percussion. If the implant shows greater than 50 % bone loss or mobility, consider possible removal of implant.
- Routine implant treatment is not recommended until the end of the jaw bone growth.
- Implant, prosthetic component or prosthesis fracture may occur due to unfavorable loading conditions. Mastication forces and restoration design have to be evaluated carefully.
- Inadequate cooling, incorrect implant site preparation, incorrect instrument settings or excessive installation torque can result in failed implant.
- All products intended for single use must not be re-used. Failure to observe this precaution may prevent osseointegration of an implant or cause infection.
- Bone quality and quantity, local infections or disturbed initial healing are all factors, which may influence implant survival. Successful healing and osseointegration may be negatively affected by the following factors:
 - Compromised general health status
 - Compromised local health status in the jaws
 - History of therapeutic radiation in the area
 - Alcohol or drug abuse
 - Smoking

MRI Safety Information

Astra Tech Implant EV has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Astra Tech Implant EV in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Adverse Reactions

During implant treatment, possible adverse reactions may include but are not limited to:

- Allergy (e.g. anaphylactic shock)
- Uncontrolled post-operative bleeding (e.g. damage to arteria lingualis)
- Aspirated or swallowed components
- Permanent nerve injury
- Fractured mandible
- Implant loss
- Peri-implantitis
- Component or prosthesis failure/fracture

Step-by-Step Instructions

Clinical, radiographic as well as model analysis, are all essential prerequisites for the successful implant treatment. Pre-operative planning should be based on the expected restorative treatment outcome.

Each patient must be carefully evaluated from a general medical and from a dental/oral point of view.

Before any treatment procedure begins, all pathology in the jaws and remaining teeth must be resolved or be in the process of healing.

- Prepare the implant site according to chosen surgical procedure and type of implant.
Follow the recommended drilling and installation protocol.
- Open the blister package.
- Pour the sterile inner container onto a sterile area.

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.
Other languages are available upon request.

Not all products are available in all countries

Status of information
IFU 5698-en_RU Rev 002 A 2019-09

**Dentsply
Sirona**
Implants

Instructions for Use

- Remove the cap from the container using a twisting motion to expose the top of the implant.
- Attach the appropriate Implant Driver EV to the contra angle.
- Make sure that the implant driver is fully seated into the implant.
- Press downwards to activate carrying function before picking up the implant. Install the implant with a contra angle at low speed, 25 rpm under profuse irrigation. Set the maximum torque to 45 Ncm.
- Allow the implant to work its way into the osteotomy. Avoid applying unnecessary pressure.
- For the final seating use Torque Wrench EV with the Torque Wrench EV Surgical Driver Handle.

NOTE: Do not exceed 45 Ncm when installing the implant.
If not completely seated before reaching 45 Ncm, reverse/ remove the implant and widen the osteotomy appropriately.

- It is recommended to have a titanium forceps available in case the implant driver does not provide sufficient carrying function during the removal procedure.

For detailed information see latest version of appropriate procedure manuals.

Cleaning and Sterilization

Not applicable.

The product is sterile.

The product is for single use only.

Delivery – Storage – Disposal

The package comprises of a cardboard outer box containing a plastic blister which is sealed by a lid. Inside the blister there is an inner package containing the implant. The lid of the blister serves as a sterile barrier and the content of the blister is sterile. The product is sterilized by irradiation. The components must not be used beyond the expiry date.

The components must be stored, in their sealed package, in a dry place, at normal room temperature, i.e. 18–25 °C.

Do not use the product if the sterile package is broken or damaged.

For disposal of the packaging and components comply with the currently applicable national waste disposal regulations in your country.

Documentation

To ensure traceability, we recommend that the removable label on the package is attached to the patient record.

Transportation Condition

Products should be transported and carried carefully in a manner that corresponds to the special transport precautions for each medical device's nature. Transport and store the components in their original packaging in a dry place at room temperature. Do not expose to direct sunlight.

Terms of Use

Dental clinics and laboratories.

Application area and operating conditions

Stomatology. Dental implantology.

Shelf life

The shelf life is set up as 5 years.

Implant System warranty

For full implant systems warranty information, refer to the detailed terms and conditions in the Implant Systems Warranty leaflet.

Trademarks

All trademarks and company names are the property of their respective owners. To improve readability for our customers, Dentsply Sirona Implants does not use ® or ™ in body copy.

However, Dentsply Sirona Implants does not waive any right to the trademark and nothing herein shall be interpreted to the contrary.

Authorized Representative on the Territory of the Russian Federation

Dentsply Sirona LLC

115432, Russia, Moscow, Andropov av. 18, bld. 6

Tel. +7 (495) 725 10 87

cis@dentsplysirona.com

Valid for Registration (Russia)

Astra Tech Implant System® OsseoSpeed EV, OsseoSpeed Profile EV, Astra Tech Implant EV OsseoSpeed®, Astra Tech Implant EV Profile OsseoSpeed®

Explanation of Symbols

Symbol	Title
	Manufacturer
	Date of manufacture
	Use-by date
	Batch Code
	Catalog Number
	Sterilized using irradiation
	Do not resterilize
	Do not use if package is damaged
	Do not reuse
	Consult Instructions for Use
	European conformity (CE) mark
	Prescription only. U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.
	Russian certification symbol according to GOST

Recommended maximum installation torque for Astra Tech Implant EV is 45 Ncm.

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.
Other languages are available upon request.

Not all products are available in all countries.

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Description of the Product

Astra Tech Implant EV Profile is a screw-shaped dental implant with a defined surface roughness, achieved through a patented grit blasting technique, followed by a process involving a treatment in diluted hydrofluoric acid.

The external coronal portion has minute threading – MicroThread.

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The implants have an inner geometry with a conical and indexing configuration and One-position-only placement of all indexed restorative components.

Astra Tech Implant EV Profile has a sloped design to match a sloped bone ridge.

Different profiles, diameters and lengths of the implants are available providing versatility according to the thickness and shape of the bone.

The package for Astra Tech Implant EV Profile is marked with "P".

Order Information and Dimensions

See current Astra Tech Implant System EV Profile product catalog for updated and detailed order information of the various implant products.

Material

Titanium – Grade 4

Indications for Use

Astra Tech Implant EV Profile is intended for both one- and two-stage surgical procedures in the following situations and with the following clinical protocols.

- Replacing missing teeth in single or multiple unit applications in the mandible or maxilla.
- Immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- Especially indicated for use in soft bone applications where implants with other implant surface treatments may be less effective.
- Immediate and early loading for all indications except in single tooth situations in soft bone (type IV) where implant stability may be difficult to obtain and immediate loading may not be appropriate.
- Only together with Profile EV restorative components, Implant Driver Profile EV, Radiographic Implant Guides Profile EV and non-indexed prosthetic components within Astra Tech Implant System EV.

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See above – Indications for use.

Warnings

Surgery in the oral cavity and oral rehabilitation include general risks for complications which should be known by the clinician. To be mentioned, but not limited to, the following

- Components that accidentally are dropped in the patient's mouth may be swallowed or aspirated which can result in suffocation or physical injury. Care must be taken to have control over small devices.
- Allergy or hypersensitivity to chemical ingredients of material used. Patient health history must be carefully assessed.
- Failure to plan and prepare implant site position and depth correctly can result in harm to nerves, blood vessels, teeth or other sensitive structures.

It is strongly recommended that clinicians always undergo special training before starting up new treatment methods.

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Precautions

- Components and instruments for Astra Tech Implant System are not compatible with Astra Tech Implant System EV.
- Fracture of mandible may occur in a very atrophic mandible, during surgery or during routine oral function after surgery.
- Highly angled abutments on small diameter implants are not recommended for use in the molar region.
- Routine implant treatment is not recommended until the end of the jaw bone growth.
- Implant, prosthetic component or prosthesis fracture may occur due to unfavorable loading conditions. Mastication forces and restoration design have to be evaluated carefully.
- Inadequate cooling, incorrect implant site preparation, incorrect instrument settings or excessive installation torque can result in failed implant.
- All products intended for single use must not be re-used. Failure to observe this precaution may prevent osseointegration of an implant or cause infection.
- Bone quality and quantity, local infections or disturbed initial healing are all factors, which may influence implant survival. Successful healing and osseointegration may be negatively affected by the following factors
 - Compromised general health status
 - Compromised local health status in the jaws
 - History of therapeutic radiation in the area
 - Alcohol or drug abuse
 - Smoking

MRI Safety Information

Astra Tech Implant EV Profile has not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of Astra Tech Implant EV Profile in the MR environment is unknown.

Scanning a patient who has this device may result in patient injury.

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- Prepare the implant site according to chosen surgical procedure and type of implant.
Follow the recommended drilling and installation protocol.
- Open the blister package.
- Pour the sterile inner container onto a sterile area.
- Remove the cap from the container using a twisting motion to expose the top of the implant.
- Attach the appropriate Implant Driver EV to the contra angle.
- Make sure that the implant driver is fully seated into the implant.

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

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Status of information:
IFU 5699-en_RU Rev. 002 A 2019-09

Instructions for Use

- Press downwards to activate carrying function before picking up the implant. Install the implant with a contra angle at low speed, 25 rpm under profuse irrigation. Set the maximum torque to 45 Ncm.
- Allow the implant to work its way into the osteotomy. Avoid applying unnecessary pressure.
- For the final seating use Torque Wrench EV with the Torque Wrench EV Surgical Driver Handle.

NOTE: Do not exceed 45 Ncm when installing the implant. If not completely seated before reaching 45 Ncm, reverse/remove the implant and widen the osteotomy appropriately.

- Positioning: align the dimple/flat surface on the implant driver with the most apical point of the slope to facilitate optimal placement of the implant.
- It is recommended to have a titanium forceps available in case the implant driver does not provide sufficient carrying function during the removal procedure.

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The product is for single use only.

Delivery – Storage – Disposal

The package comprises of a cardboard outer box containing a plastic blister which is sealed by a lid. Inside the blister there is an inner package containing the implant. The lid of the blister serves as a sterile barrier and the content of the blister is sterile. The product is sterilized by irradiation.

The components must not be used beyond the expiry date.

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Documentation

To ensure traceability, we recommend that the removable label on the package should be attached to the patient record.

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Explanation of Symbols

Symbol	Title
	Manufacturer
	Date of manufacture
	Use-by date
	Batch Code
	Catalog Number
	Sterilized using irradiation
	Do not restenize
	Do not use if package is damaged
	Do not reuse
	Consult Instructions for Use
	European conformity (CE) mark
	Prescription only U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.
	Russian certification symbol according to GOST

Recommended maximum installation torque for Astra Tech Implant EV Profile is 45 Ncm.

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

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