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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of ZHEJIANG GONGDONG MEDICAL TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

Su Jing

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(Date: Nov. 27, 2023)

书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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
Option No. 1:

1) Vacuum sterile plastic tube LIQUIDPRF – 100 pcs. in individual packaging

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Option No. 2:

1) Vacuum sterile glass tube CLOTPRF – 100 pcs. in individual packaging

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2. Information about the manufacturer and site of production

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Zhejiang Gongdong Medical Technology Co., Ltd.

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Address of the location of the manufacturer and production:

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No.10 Beiyuan Ave., Huangyan 318020 Taizhou, Zhejiang, People's Republic of China

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Venous blood collection and autologous platelet-rich fibrin (PRF) creation for mixing with a bone graft to correct bone defects in donor patients, in the dental field.

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Vacuum Sterile Plastic Tube LIQUIDPRF is used to prepare liquid autologous platelet-rich fibrin (PRF) in the dental field.

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Vacuum Sterile Glass Tube CLOTPRF is used to prepare autologous platelet-rich fibrin (PRF) in solid form, in the dental field.

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5. Application environment and would-be users

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The medical device is used exclusively in the conditions of dental institutions by qualified medical personnel (doctor, nurse).

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There is no direct contact with patients' (human) body.

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Vacuum Sterile Plastic Tube LIQUIDPRF is used to prepare autologous platelet-rich fibrin (PRF) in liquid form.

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Vacuum Sterile Glass Tube CLOTPRF is used to prepare autologous platelet-rich fibrin (PRF) in solid form.

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There are no side effects subject to the storage, transportation and use requirements for the medical device set forth in the instructions for use.

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The principle of operation of the product is based on the presence of precisely dosed vacuum in hermetically sealed test tubes. When a vein is needled with one end of a special double-sided needle and a rubber stopper of a vacuum tube with the other, a required volume of blood is drawn directly from the vein into the test tube under vacuum.

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11. Description

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Golden Vac blood collection and processing medical vacuum sterile tube is sterile single-use vacuum tubes intended for collection of venous blood and creation of autologous platelet-rich fibrin for mixing with a bone graft to correct bone defects in donor patients, in the dental field. The tubes can be used in a centrifuge.

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There are following design versions of the Golden Vac blood collection and processing medical vacuum sterile tube: Vacuum Sterile Plastic Tube LIQUIDPRF and Vacuum Sterile Glass Tube CLOTPRF with their catalogue (REF) numbers: GD0100A1 and GD0100AG1, respectively.

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Vacuum Sterile Plastic Tube LIQUIDPRF is used to prepare liquid autologous platelet-rich fibrin (PRF), in the dental field. The tube is made of polyethylene terephthalate (PET), a clinically tested safe alternative to glass, with a butyl rubber stopper and an orange polyethylene cap.

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
Vacuum Sterile Glass Tube CLOTPRF is used to prepare solid autologous platelet-rich fibrin (PRF) in the dental field. The test tube is made of medical glass and has a butyl rubber stopper and a red polyethylene cap.

The cap can be easily removed and put on the test tube and acts as a protective mean, which prevents blood from contacting with medical personnel.

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The tube is labelled. Each tube is packed in its primary package - an individual bag that provides protection and tightness. A transport cardboard box contains 100 primary packages, placed in a plastic bag. The instructions for use are enclosed in the transportation box.

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The tolerance of overall dimensions, nominal capacity and weight of the product should be no more than $\pm 10\%$.

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13. Materials used to manufacture the medical device

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- 1) Choose the type of tube you need (Vacuum sterile glass tube CLOTPRF and Vacuum sterile plastic tube LIQUIDPRF)

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Secure the blood sampling needle in the holder. Make sure the needle fits snugly in the holder.

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- 3) Ask the patient to move the arm down. Select an area of venepuncture. Apply a tourniquet and treat the venepuncture site with an antiseptic.

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- 4) Remove the needle cap. When the vein is punctured, the patient's hand should be down and the tube stopper up.

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- 5) When piercing the stopper, place the tube exactly in the centre of the holder, otherwise there is a risk of puncturing the wall of the stopper, which will lead to a rapid loss of vacuum. Press down on the tube so that the needle in the holder enters the stopper.

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- 6) Remove the tourniquet as soon as blood begins to flow into the tube. If the volume of blood is less than the specified volume in the tube or blood does not flow into the tube, it is recommended to do the following to complete the collection:

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When the tube is filled to the specified volume and blood has stopped flowing into the tube, remove the tube from the holder.

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If necessary, repeat steps 5 ~ 7 to collect further samples from the same patient.

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When the last tube is full, remove the needle from the vein. Apply a sterile tissue to the puncture site and ask the patient to hold it until the bleeding stops.

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A small amount of blood may remain on top of the stopper after being drawn. When handling tubes, take the necessary precautions to avoid contact with blood samples.

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Dispose of needles according to your institution's requirements and regulations.

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Any blood return can harm the patient. Take the following precautions to avoid the blood return:

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Centrifugation is an important step that is required to separate the liquid part of blood from the cells.

Centrifugation errors may lead to the adverse effects:

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Make sure that tube caps do not rest on the walls of the centrifuge beaker: otherwise they may open during centrifugation.

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Centrifugation conditions for CLOT-PRF are:

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1300 g 10 min, 25 °C.

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Centrifugation conditions for LIQUID-PRF are:

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Restrictions: The amount of blood drawn depends on the altitude, ambient temperature, atmospheric pressure, needle diameter, patient's blood pressure, and blood viscosity.

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Do not use tubes with no or damaged caps, as well as with cracks or contamination.

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Do not use tubes after the expiration date.

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The tubes are for single use only. They must be disposed of in biosecure containers approved for this purpose.

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Use gloves and personal protective equipment to protect against splashes and spills of blood and exposure to blood borne pathogens.

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Place the tubes in the centre of the holder while piercing the cap to prevent puncture of the side wall and premature loss of vacuum.

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One device should only be punctured once during blood collection. If there is more than one puncture, then partial or complete loss of vacuum is possible.

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To avoid unexpected detachment of the stopper from the tube, hold the tube in the holder during the entire blood collection process by pressing on the bottom of the tube with your thumb.

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16. Information on compatible medical devices

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With vacuum tubes, use disposable double-sided needles, a hemostatic tourniquet, a holder, and a disposable disinfectant wipe (or a wipe with an antiseptic), approved for use in accordance with the legislation of the country of use. There are no technical specific requirements for compatible products.

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Product in primary packaging (package): 100 pcs.

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Shipping package (cardboard box): 1 pcs.

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It is desirable to combine force with a slight rotation of the stopper relative to the tube to reduce the force when detaching the stopper from the tube.

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- The product must be used before the expiration date.

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- The product does not emit toxic substances into the environment during storage and operation, and does not adversely affect humans when it is in direct contact.

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No maintenance, this product is not repairable.

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20. Environmental protection requirements when using the medical device

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Vacuum tubes do not adversely affect humans and the environment during operation, transportation and storage, provided that consumers meet the requirements set forth in the instructions for use.

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Golden Vac blood collection and processing medical vacuum sterile tube is supplied sterile.

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The radiation (gamma radiation) sterilization method is used, the absorbed dose is from 14.8 to 20.0 kGy. No additional sterilization is required.

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Primary packaging is a hermetically sealed polyethylene and paper bag.

Первичная упаковка – герметично запечатанный пакет из полиэтилена и бумаги.
Primary packaging functions as individual packaging. One side of the package is opaque and marked. The other one is transparent for ease of product identification. Each product (test-tube, sealed with a stopper and closed with a lid) is placed in a primary package.

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Products in the amount of 100 pieces in primary packaging are placed in a cellophane bag with a lock, which is placed in a transport package - a corrugated cardboard box. The shipping container is marked. Instructions for use of 1 pc. included in the shipping package.

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Labels of the primary packages of the product contain the following information:

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Labels contain the following information:

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
The following information is placed on the cardboard box:

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23. Conditions of transportation, storage and operation

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Storage

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During the entire shelf life, the products should be stored in the manufacturer's packaging in dry ventilated rooms at temperatures from + 2 to + 25 °C and relative humidity not more than 80%, at least 1 m from heat-emitting devices.

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During storage, the products must be protected from sunlight and atmospheric influences, must not be exposed to oils, gasoline, kerosene and other substances.

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The packed product is transported by all types of covered vehicles, in accordance with the haulage rules in force for this type of transport, at a temperature from -15 to +40 °C and a relative humidity of no more than 80%.

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The product does not require any maintenance and is beyond repair.

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The manufacturer guarantees the compliance of the products with the requirements, subject to the conditions of operation, transportation and storage.

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The warranty covers the entire shelf life of the product.

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All claims and appeals of consumers are accepted at the following address in the territory of the Russian Federation: Implant Guru, LLC, 43, Pokrovka Street, building 2, 105062, Moscow, Russia.

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After use, vacuum tubes are class B "epidemiological hazardous waste", unused tubes are class A "epidemiologically safe waste, close to solid utility waste in composition" according to SanPiN 2.1.3684-21

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ISO 11607-1:2019 «Packaging for terminally sterilized medical devices. Part 1: Requirements for materials, sterile barrier systems and packaging systems»

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ISO 15223-1:2021 «Medical products. Symbols used for medical devices labeling, tags and accompanying documents. Part 1. General requirements»

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29. Authorized representative in the RF

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For questions regarding the quality of the medical device, please, contact:

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