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Lessa Defi-Pad

Disposable multifunction electrodes for external
defibrillation synchronized cardioversion, cardiac
stimulation and ECG monitoring

Instruction for Use

Revision of 03.11.2023

Contacts

Manufacturer:

AB MEDICA GROUP, S.A.

Llobateras, 14-16-18, Talleres 7, Naves 4 y 6. Pol. Ind. Santiga

08210 Barberà. Barcelona. SPAIN

Tlf.: +34 93 7299400 Fax: +34 93 7191656

Authorized deputy in Russia:

Limited Liability Company

Emergency Medical Systems (EMS LLC)

Office 1, 1 Sormovskoe shosse str., 603074 Nizhny Novgorod, Nizhny Novgorod Oblast, Russia

Phone: +7 (831) 217-00-03

E-mail: info@emrmedsys.com

CAUTION!

Please, inform the manufacturer and/or its authorized deputy about all adverse incidents involving the product.

This Instruction for Use is applicable to the Lessa Defi-Pad disposable multifunction electrodes for internal defibrillation, synchronized cardioversion, cardiac stimulation and ECG Monitoring (hereinafter called Lessa Defi-Pad electrodes).

Upon agreement with the customer, the Instruction for Use can be provided in a single copy for a group of products.

CAUTION!

Please, pay special attention to such warnings in the text hereof, in order to avoid potential dangerous situations or operating errors.

Content

1 Symbols	5
2 Device description.....	6
2.1 Intended use	6
2.2 Intended users	6
2.3 Intended patients.....	6
3 Intended use.....	7
3.1 Indications for use.....	7
3.2 Contraindications for use.....	7
3.3 Potential side effects	7
3.4 Precautions for handling	7
3.5 Positioning	9
3.6 Operating procedure	9
4 Characteristics.....	11
4.1 Specification	11
4.2 Materials	12
4.3 Compatibility with defibrillator models.....	12
4.4 Product parts.....	15
4.5 Packaging details	15
4.6 Product disposal.....	15

Symbols

The following symbols are used in the product labelling:

Table 1.1 Labeling symbols

Symbol	Meaning	Symbol	Meaning
	see Instruction for Use		non-sterile
	lot No.		avoid exposure to sunlight
	reference No.		transporting and storage temperature range
	unique device identifier		operating temperature limits
	manufacturer		relative humidity range
	expiry date		latex-free
	EC Certificate with the Notified Body number		recycling sign (composite material)
	do not reuse		packaging: two electrodes in one pouch
	do not use when the package is damaged		Medical Device
	warning		

2 Device description

2.1 Intended use

Lessa Defi-Pad electrodes are constituted by a pair of adhesive pads provided with gels. It may be used in place of the manual reusable paddles.

Lessa Defi-Pad electrodes are indicated for:

- transthoracic external defibrillation.
- transthoracic synchronized cardioversion.
- transthoracic electrocardiographic (ECG) monitoring.
- temporary transthoracic cardiac pacing (non-invasive).

Lessa Defi-Pad electrodes allow the user to effectively operate in the treatment of rhythm disorders related to the above-mentioned applications, without the risk of accidental electrocution related to the use of normally available reusable paddles.

2.2 Intended users

Lessa Defi-Pad electrodes are designed to be used by health care professionals and specially trained people.

2.3 Intended patients

Depending on the model, Lessa Defi-Pad electrodes are designed for various types of patients in accordance with the labelling on the primary package (Figure 2.1):

- for adults: patients of over 8 years old or weight of over 25 kg;
- for pediatrics: patients of less than 8 years old or weight of less than 25 kg.

Notes:

1 Lessa Defi-Pad electrodes for adults are generally contraindicated in patients younger than 8 years old (weighing less than 25kg), but can be used if the size of the chest allow it, taking care that one pad do not come in contact with the other one. Follow the operating instructions of the defibrillator for energy to be delivered.

2 The use of Lessa Defi-Pad electrodes for pediatrics is generally contraindicated in patients older than 8 years old (weighing more than 25kg); In fact, the power release above 100J may increase the risk of burns.

Figure 2.1 — Patient-type labelling

Note — Patient-type labelling of electrodes, designed both for adults and pediatrics (universal), includes both labels (a & b Figure 2.1).

3.1 Indications for use

Lessa Defi-Pad electrodes shall be used as a part of a defibrillator, for urgent and emergency medical aid in case of acute cardiovascular diseases (events), as medically indicated.

3.2 Contraindications for use

Lessa Defi-Pad electrodes are not allowed to apply on skin that shows signs of irritation or injury.

3.3 Potential side effects

Plate adhesive may cause light cutaneous irritation.

The prolonged transthoracic stimulation or the repeated administration of defibrillation shock may cause more or less noticeable skin reddening according to the supplied energy.

A lack in adhesion and/or air presence under the electrode may cause burnings.

3.4 Precautions for handling

When handling the product, observe the following precautions:

- Use Lessa Defi-Pad electrodes only with the specified defibrillator models provided by labelling on the primary package (see 4.3 « Compatibility with Defibrillator Models»).
- Before use check the expiration date on the package. Do not use Lessa Defi-Pad electrodes after this date.
- Check that the packaging is intact: Do not use Lessa Defi-Pad electrodes otherwise.
- Use Lessa Defi-Pad electrodes within 24 h after the package is opened.
- Lessa Defi-Pad electrodes are to be applied within 30 minutes after removal of the protective coating.
- Do not use Lessa Defi-Pad electrodes if the gel is removed from the support or if it is ripped, torn or dry. Any discoloration localized on gel or on conductive foil does not affect the functionality of the product.
- Do not use Lessa Defi-Pad electrodes with signs of visible damage (e.g. the insulating coating of the contact, connector or cable has detached or there are tears in the foam support and / or in the electrode).
- Do not bend, do not cut and do not squash Lessa Defi-Pad electrodes.
- Lessa Defi-Pad electrodes neither are sterile nor can be sterilized.
- Lessa Defi-Pad electrodes are disposable. For use on a single patient. Discard after use.
- After an extended period of transthoracic pacing the ability to detect the evoked ECG signal can be reduced. In this case it is necessary to provide for the collection of the evoked ECG by a separate set of electrode pads for ECG.

CAUTION!

When handling the product, remember that:

- *Keep Lessa Defi-Pad electrodes well separated one from the other and be careful not to overlap them with other objects (ECG electrode pads, cables, transdermal patches, clothing etc.)*
 - *Lessa Defi-Pad electrodes needed to be replaced after 50 defibrillation shocks, 8 hours of pacing (30 minutes of pacing with pulses longer than 20ms) or after 24 hours on the patient's skin.*
 - *Do not supply a shock with manual metal paddle above the disposable electrodes Lessa Defi-Pad.*
 - *In case of prolonged pacing (greater than 30 minutes), check the skin of the patient for signs of irritation.*
 - *If the tracing is not sufficiently clear or in case of demand pacing, use an ECG patient cable and a separate set of ECG electrodes (if it is present).*
 - *Avoid applying the adhesive pad on the nipple or breast tissue.*
 - *Do not replace the adhesive pads once applied. If the position must be changed, remove and replace with new multifunction electrode pads.*
 - *To prevent accidental damage from electric shock, ensure that during discharge operators are not in contact with the electrode pads, with the patient, or with conductive parts close to the patient.*
 - *Avoid contact between parts of the patient's body such as exposed skin of head or limbs, conductive fluids such as gel, blood or saline and metal objects such as a bed frame or a stretcher which may provide unwanted pathways for the defibrillating current.*
 - *When defibrillators are used near oxygen sources or other flammable gases, use extreme care to avoid risk of fire or explosion.*
 - *All other medical devices that haven't applied parts with defibrillation protection must be separated from the patient during defibrillation.*
-

For ease of placement and for training purposes, the anterior-lateral side (Figure 3.1 a) is preferred for defibrillation and cardioversion; the anterior-posterior side (Figure 3.1 b) is more common in dynamics and in transthoracic pacing.

a) Anterior-lateral

b) Anterior-posterior

Figure 3.1 — Positioning

3.6 Operating procedure

- 1** Select the electrode application area
- Select the method for electrode application on the patient's body according to the Defibrillator's Instruction for use.
- 2** Prepare the patient's skin (if required)
- Check the skin in the electrode application area: it shall be clean, dry, without hair. Shave any hair and wipe the skin, if required (Figure 3.2 a, b).
- 3** Unpack the electrodes
- Unpack the electrodes, remove the electrodes from the primary package (Figure 3.2 c) and remove protective film from the electrodes (Figure 3.2 d).
- 4** Apply the electrodes
- Apply the adhesive pads one-to-one starting with one side and pressing progressively over the entire surface to avoid the formation of air bubbles and ensure complete adhesion to the skin (Figure 3.2 e).
- 5** Connect the electrodes
- Connect the electrodes' connector to the Defibrillator used (or the adapter cable).

6

Use according to the instruction

Proceed in accordance with instructions of the Defibrillator's Instruction for use.

Note — Do not leave the product on the patient's skin for more than 24 h.


7

Remove the electrodes

Remove the electrodes by slowly peeling them from the patient's skin. Dispose of the electrodes in accordance with the applicable national regulations for disposal of medical waste.

a) Shave hair


b) Wipe the skin

c) Remove the electrodes

d) Remove protective film

e) Press the electrodes down

Figure 3.2 — Operating procedure

4 Characteristics

4.1 Specification

Lessa Defi-Pad electrodes specification is provided in Table 4.1.

Table 4.1 — Lessa Defi-Pad electrodes specification

Parameter		Value
Pad geometry		
1 Active surface area	for adults	$\geq 50 \text{ cm}^2$
	for pediatrics	$\geq 15 \text{ cm}^2$
2 Foam surface area		$> 78 \text{ cm}^2$
3 Foam thickness		$\geq 1 \text{ cm}$
4 Gel surface area		$> 60 \text{ cm}^2$
5 Liner surface area (plastic support of compatible shape)		$> 120 \text{ cm}^2$
6 Gel rim – electrode rim distance		$> 3 \text{ mm}$
7 Gel rim – foam rim distance		$> 5 \text{ mm}$
Pad appearance		
8 Shape		Circular or similar
9 Adhesive appearance		Uniform
10 Gel appearance		Uniform
11 Appearance 30 days after production (gel, electrode, adhesive) at 25°C, 75% Rh		Uniform and homogeneous
12 Concentricity (electrode – foam – gel)		$\pm 1 \text{ mm}$
Cable geometry and appearance		
13 Length		$< 2 \text{ m}$
14 Cable appearance		Depending on the connector used
Connector geometry and appearance		
15 Shape and size of insulating housing, diameter and length of contacts, type of crimp contacts		Depending on compatible defibrillator
Electric functionality		
16 Test waveform attenuation at first 360 J discharge		$< 5 \%$
17 Test waveform attenuation after 50 360 J discharge		$< 5 \%$
18 Gel temperature after 15 360 J discharge 3 discharge/min		$< 41 \text{ }^\circ\text{C}$
Dielectric strength		
19 Cable insulation, altering voltage AV		$> 4000 \text{ Veff}$
20 Cable insulation, continuous voltage CV		$> 7500 \text{ Vdc}$
21 Cable insulation, DC resistance		$> 500 \text{ M}\Omega$
22 Insulation of the upper pad protection, altering voltage AV		$> 4000 \text{ Veff}$
23 Insulation of the upper pad protection, continuous voltage CV		$> 7500 \text{ Vdc}$
24 Insulation of the upper pad protection, DC resistance		$> 500 \text{ M}\Omega$
Separation		
25 Residual voltage peak on the cable (and pad) after 360 J discharge		$< 1 \text{ V}$

Parameter	Value
Recovery	
26 Signal recovery time after 360 J discharge	< 10 s
Leakage	
27 Leakage current in patient with mains voltage (250 Veff, 50 Hz) on AP (applied part)	<100 µA
Further electric parameters	
28 Small signal impedance Z _{acp} at 10 Hz	<3000 Ω
29 Small signal impedance Z _{acp} at 30 Hz	<5 Ω
30 Continuous offset voltage at rest V _{os}	<100 mV
31 Noise voltage and offset instability V _n	<100 µV _{p-p}
32 Residual voltage after defibrillation V _r at 4 s	<300 mV
33 Residual voltage after defibrillation V _r at 60 s	<400 mV
34 Impedance to large signals (only for gel-to-gel coupled pads) Z _{acg}	<3 Ω
Mechanical functionality	
35 Pad flexibility, curvature radius	< 1,5 cm
36 Pad adhesiveness, force of detachment from test surface	> 1 kg
37 Tensile strength of the electrode-cable connection	> 5 kg
38 Mechanical compatibility with the defibrillator connector	Corresponds designated defibrillator
Environmental conditions	
39 Operating temperature limits	from 0 °C to plus 50 °C
40 Transporting temperature limits	from plus 5 °C to plus 35 °C
41 Storage temperature limits (long term)	from plus 5 °C to plus 35 °C
42 Storage temperature limits (short term until 24 h)	from minus 30 °C to plus 65 °C
43 Relative humidity range	from 20 % to 80 %
Shelf life	
44 Shelf life (storage)	2 years since manufacturing (see package)
Weight	
45 Total weight	< 170 g

4.2 Materials

The materials coming into short-term contact with the human's intact skin are not toxic.

4.3 Compatibility with defibrillator models


Lessa Defi-Pad electrodes belong to the product family: 99.835.XX. Reference No determines the type of connector and features of the electrodes. Compatibility with defibrillator models is provided in Table 4.3.

Note — The combination of Lessa Defi-Pad electrodes and a defibrillator forms a medical system.

The following types of electrodes are provided:



- radiotransparent;
- pre-connected (external wire);
- adult, pediatric and universal.









Some paediatric models to be used with automatic defibrillators (AEDs) provide a system for reducing

Table 4. Compatibility with defibrillator models

Manufacturer	Defibrillator models and brands	Patient type	Radio-transparent	External wire	Reference	Connector design
PHILIPS	HeartStart MRx (with cable M3508A), HeartStart XL+, HeartStart FRx, Efficia DFM100	Adult	-	-	99.835.00	
		Universal	-	+	99.835.70U	
	HeartStart XL+, Efficia DFM100	Pediatric	-	-	99.835.01	
ZOLL	M Series, M Series CCT, E Series, R Series, X Series	Adult	-	-	99.835.02	
		Pediatric	-	-	99.835.03	
		Adult	+	-	99.835.51	
		Pediatric	+	-	99.835.51P	
		Adult	-	+	99.835.71	
	AED Plus, AED Pro	Pediatric	-	+	99.835.97	
PHYSIO-CONTROL	LifePak 12, LifePak 20 Monitor, LifePak CR Plus, LifePak 15, LIFEPAK 20e	Adult	-	-	99.835.04	
		Pediatric	-	-	99.835.05	
		Adult	+	-	99.835.52	
		Pediatric	+	-	99.835.52P	
		Adult	-	+	99.835.72	
		Pediatric	-	+	99.835.72P	
	LifePak 1000, LifePak 500	Adult	-	-	99.835.04	
		Pediatric	-	-	99.835.87	
		Adult	+	-	99.835.52	
		Adult	-	+	99.835.72	
		Pediatric	+	-	99.835.52P	
		Pediatric	-	+	99.835.72P	
COMEN	S3, S5, S6, S8	Adult	-	-	99.835.04	
		Pediatric	-	-	99.835.05	
		Adult	+	-	99.835.52	
		Pediatric	+	-	99.835.52P	
		Adult	-	+	99.835.72	
		Pediatric	-	+	99.835.72P	

Manufacturer	Defibrillator models and brands	Patient type	Radio-transparent	External wire	Reference	Connector design
MINDRAY	BeneHeart C1 Fully Automatic, BeneHeart C1A Fully Automatic, BeneHeart C2 Fully Automatic, BeneHeart C2A Fully Automatic, BeneHeart C1, BeneHeart C1A, BeneHeart C2, BeneHeart C2A, BeneHeart D3, BeneHeart D6, uMED 20, uMED 20A	Adult	—	—	99.835.04	
		Pediatric	—	—	99.835.05	
		Adult	+	—	99.835.52	
		Pediatric	+	—	99.835.52P	
		Adult	—	+	99.835.72	
		Pediatric	—	+	99.835.72P	
SCHILLER	DEFIGARD 1002, DEFIGARD 3002	Adult	—	—	99.835.08	
		Pediatric	—	—	99.835.09	
	DEFIGARD 4000, DEFIGARD 5000, DEFIGARD 6002, FRED easy	Adult	—	—	99.835.12	
		Pediatric	—	—	99.835.13	
	Defigard HD-7, FRED PA-1, FRED easyport	Adult	—	—	99.835.62	
		Pediatric	—	—	99.835.63	
PROGETTI	RESCUE SAM	Adult	—	—	99.835.18	
		Pediatric	—	—	99.835.61	
AXION	DA-N-01, DA-N-02, DA-N-05, DKI-N-11	Adult	—	—	99.835.18	
		Pediatric	—	—	99.835.19	
CORPULS	Corpuls 3	Adult	—	—	99.835.36	
		Pediatric	—	—	99.835.37	
METRAX	PRIMEDIC HeartSave PAD (SN up to 73944xxxxxx), PRIMEDIC HeartSave AED (SN up to 73945xxxxxx), PRIMEDIC HeartSave AED-M (SN up to 73946xxxxxx), PRIMEDIC HeartSave 6 (SN up to 73947xxxxxx), PRIMEDIC HeartSave AS (SN up to 73943xxxxxx), PRIMEDIC DefiMonitor	Adult	—	+	99.835.79	
WEINMANN	MEDUCORE Standard	Adult	—	+	99.835.84	
MEDIANA	A15, A15NP	Universal	—	+	99.835.89	

4.3 Product parts
The product parts (delivery set) are given in the table 4.3.

Part name	Manufacturer	Country
Lessa Defi-Pad disposable multifunction electrodes for external defibrillation synchronized cardioversion, cardiac stimulation and ECG Monitoring	AB MEDICA GROUP, S.A.	Spain
Instruction for use *	AB MEDICA GROUP, S.A.	Spain

** By agreement with the customer, the Instruction for use can be provided in a single copy for a group of products.*

4.5 Packaging details

Each pair of Lessa Defi-Pad electrodes is packaged in sealed pouches of opaque material suitable to protect the gel from light and moisture (primary packaging).
The carton box is used as the secondary packaging.

4.6 Product disposal

Dispose of Lessa Defi-Pad electrodes in accordance with the national regulations for the disposal of medical waste in force at the time of disposal.

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This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

[This Apostille is not valid for use anywhere within Spain]

[To verify the issuance of this Apostille, see <https://eregister.justicia.es/>]

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