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
Ulrich

1st of March, 2023
«Approved»

(signature)

Ulrich A. Becht

Managing director, Alfred Becht GmbH.

Alfred Becht GmbH
Postfach 11 45
D-77601 Offenburg
Tel.: + 49 781 60586-0
FAX: +49 781 60586-40

**INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE**

Articulating paper CLEAN+SAFE® for interocclusal registration

Edition 2

1. Name of the medical device

Articulating paper CLEAN+SAFE® for interocclusal registration

Versions:

1. Articulating strips, double-sided, straight (blue/red), 80 µm
2. Articulating strips, double-sided, straight (blue/blue), 80 µm
3. Articulating strips, double-sided, horseshoe type (blue/red), 80 µm
4. Articulating strips, double-sided, horseshoe type (blue/blue), 80 µm
5. Articulating strips, double-sided, straight (blue/red), 40 µm
6. Articulating strips, double-sided, straight (blue/blue), 40 µm

Hereinafter referred to as the CLEAN+SAFE® articulating paper for interocclusal registration, articulating paper, double-sided straight strips, double-sided strips horseshoe type, strips, medical device, product.

2 Information about Manufacturer, Developer and Authorized Representative of Manufacturer

Manufacturer:

Alfred Becht GmbH
Carl-Zeiss-Strasse 16, 77656 Offenburg, Germany
Tel.: + 49 781 60586-0
Fax: + 49 781 60586-40
www.becht-online.de

Developer:

Alfred Becht GmbH
Carl-Zeiss-Strasse 16, 77656 Offenburg, Germany
Tel.: + 49 781 60586-0
Fax: + 49 781 60586-40
www.becht-online.de

Authorized representative of manufacturer in the Russian Federation

PROTECO, JSC
Legal address: 192102, St. Petersburg, Salova st. 57, build. 1, lit. I.
Tel.: +7 (812) 635-88-90
E-mail: info@protecodent.ru

3 Intended Purpose and Intended Users

The medical device "Articulating paper CLEAN+SAFE® for interocclusal registration" is designed to be used in dentistry for detection of occlusal contacts and distribution of occlusal forces of newly installed dental restorations and dental prostheses and bridges, as well as for detection of occlusal

INSTRUCTION FOR USE

interferences.

Intended Users

Articulating paper is recommended for use by dentists and chairside assistants.

4 Device Description and Design

Articulating paper CLEAN+SAFE® for interocclusal registration is thin and tear-resistant sheets of paper in the form of straight stripes or horseshoe type (U-shaped) stripes, impregnated with pigment that colors the surface of teeth, dental prostheses, etc. when pressed. The product is available in two color options: blue/blue sides and blue/red sides. Differently colored sides (blue and red) allow marking the places of occlusion simultaneously and accurately. The high color intensity ensures clear results. The small thickness makes it possible to avoid false or blurred contacts between teeth and paper, as well as provides clear results when used on wet occlusal surfaces (gold, ceramics, polished metal, acryl).

The applied part of the articulating paper is a part of the sheet impregnated with pigment for coloring the surface of teeth, dental prostheses, etc. when pressed.

Articulating strips, double-sided, straight (blue/red), 40 µm and Articulating strips, double-sided, straight (blue/blue), 40 µm have a design feature – the strips are not fastened into a block, they are placed in the primary packaging separately from each other. Each strip is delimited by a protective insert. The limiter holds the strips in one position inside the primary packaging.

5 Indications, Contraindications and Possible Side Effects Associated with Intended Use of Medical Device

Indications: The articulating paper is a diagnostic dental product used in the following cases:

- Dental treatment (placement of fillings/crowns and prosthetics).
- Dental implantation (determination of the stability of various abutments).
- Determination of the spatial relationship of dentitions in occlusion.
- Identification of disorders in persons with bruxism.

Contraindications: Not detected.

Possible side effects associated with intended use of medical device: Irritation or allergic reaction due to the individual intolerance of raw materials of the device by a user.

6 Precautions

Precautions:

- Do not use the product for any purpose other than those listed in these instructions.
- Read the instructions for use before use.
- Do not use the product after the expiration date.
- This is a single-use product. Reuse of the product is prohibited.
- Use of the product only with gloves.

7 Method of Use

1. Before use, remove the strip from the package with forceps.

2. Plastic box of articulating strips (blue/red), (blue/blue) 40 μm , allows to withdraw strips easily with one hand.
3. Place the articulating paper on the teeth, close the upper and lower dental arches with light pressure, making eccentric movements of the lower jaw.
4. Open the upper and lower dental arches and remove the articulating paper. Marks will remain on the teeth showing where the dental arches get in contact.
5. Remove the existing eminences at the contact points of the filling using finishing diamond burs.
6. After initial treatment of high contacts, repeat the procedure as many times as necessary to obtain anatomically correct results.

8 Technical features.

Table 1 Articulating Paper Specifications

Product Name	Thickness (with pigment), μm	Overall size of the product / Dimensions of the applied part, mm
Articulating strips, double-sided, straight (blue/red)	80 $\pm$ 5	20 x 105 ($\pm$ 2) / 20 x 85 ($\pm$ 2)
Articulating strips, double-sided, straight (blue/blue)	80 $\pm$ 5	20 x 105 ($\pm$ 2) / 20 x 85 ($\pm$ 2)
Articulating strips, double-sided, horseshoe type (blue/red)	80 $\pm$ 5	23 x 68 ($\pm$ 2) / 23 x 60 ($\pm$ 2)
Articulating strips, double-sided, horseshoe type (blue/blue)	80 $\pm$ 5	23 x 68 ($\pm$ 2) / 23 x 60 ($\pm$ 2)
Articulating strips, double-sided, straight (blue/red), 40 μm	40 $\pm$ 5	20 x 100 ($\pm$ 2)
Articulating strips, double-sided, straight (blue/blue), 40 μm	40 $\pm$ 5	20 x 100 ($\pm$ 2)

9. Description of symbols used on consumer and transport package labeling

Table 2. Symbols and description

Symbol	Description	Symbol	Description
	Manufacturer		Manufacturing date
	Lot number		Best before
	Non-sterile		Keep dry
	Temperature limits		Product must not be disposed of as household waste
	Keep away from sunlight		Do not reuse

INSTRUCTION FOR USE
FOR THE MEDICAL DEVICE

Symbol	Description	Symbol	Description
	Consult instructions for use	[REF]	Reference number (catalogue number)
	Do not use if the package is damaged	CE	CE conformity mark

10 Requirements for Maintenance and Repair of Medical Device

The devices do not require any maintenance and repair.

11 Sterilization, Cleaning and Disinfection

All articulating paper versions are for single use only and supplied non-sterile. Disinfection and cleaning are not required.

12 Conditions of Transportation, Storage and Use

Products packed in shipping box should be transported by all types of covered vehicles according to applicable regulations for the carriage of goods by this type of transport, observing the precautions indicated on the shipping box.

Transportation conditions:

temperature: 0 to +25 °C;
relative humidity: 45-55 %;
atmospheric pressure: 800-1060 hPa.

Storage. The packaging provides protection against mechanical and climatic factors during transportation and storage. Products should be stored in well-ventilated areas away from sources of heat, sunlight, direct exposure to light sources and electrical equipment. Packages must be protected against dirt, dust, moisture and rupture. Store the product in a dry, clean and safe place in tightly closed packaging.

Storage conditions:

temperature: 0 to +25 °C;
relative humidity: 45-55 %;
atmospheric pressure: 800-1060 hPa.

Use. Devices are used in healthcare facilities, in conditions requiring to protect the hands of medical stuff when using the product for its intended purpose.

Transportation conditions:

- temperature: 0 to +25 °C;
- relative humidity: 45-55 %;
- atmospheric pressure: 800-1060 hPa

13 Disposal

The medical device "Articulating paper CLEAN+SAFE® for interocclusal registration" is considered hazardous waste of epidemiologic potential (devices contaminated with blood and/or other biological fluids of patients) after their use as intended.

In this case, the medical organization carries out only collection, decontamination, and temporary storage of the abovementioned waste. Collection is carried out using bags attached to a trolley or disposable, color-coded containers for hazardous medical waste indoors or in places for temporary storage of medical waste before transporting for destruction or disposal.

Specialized organizations licensed to handle hazardous medical waste of epidemiologic potential transport medical waste from medical institutions and recycle (dispose of) and destroy it.

Consumer packaging, defective and unused products with expired shelf life are classified as epidemiologically safe wastes (non-hazardous waste that does not come into contact with biological fluids of patients), similar to municipal solid waste (MSW).

Non-hazardous waste is collected using reusable containers or disposable bags (excluding color-coded bags for hazardous waste). Disposable bags are placed inside reusable containers. Waste containers should be marked and placed at a special site (indoors). After emptying, reusable containers are subject to washing and disinfection. The procedure for washing and disinfection of reusable containers is determined in accordance with the organization-specific waste management flowchart.

Transportation and disposal are performed considering the requirements established for a specific territory in accordance with the country-specific sanitary laws concerning the maintenance of populated areas and handling of industrial and household waste.

14 Shelf Life of Medical Device

The shelf life of the medical device is 10 years.

15 List of International Standards Applicable to Medical Device

EN ISO 13485, EN ISO 14971, EN ISO 15223-1, EN 1041, Директива 93/42/EWG (2007/47/EG), EN ISO 62366, EN ISO 10993-1, EN ISO 10993-5, EN IEC 17050-1, ZLG 3.9 A4, MEDDEV 2.7/1, Regulation (EU) 2017/745 MDR.

Note: The standards are used in the latest edition: the organization maintains and updates the list of reference standards as a record of its own management system in accordance with ISO 13485.

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D-Offenburg, 26.07.2023

CERTIFICATION

I hereby certify, that the above is the true signature, subscribed in my presence of

Mister Ulrich Alfred BECHT,

born 29.04.1957,

residing in Offenburg,

business address

D-77656 Offenburg, Carl-Zeiss-Str. 16,

– personally known to me –

After having inspected today the

**Commercial Register of the local court of
Freiburg im Breisgau / Germany**

I certify the following:

Incorporated there is in section

– HRB 471050 –


the company

**Alfred Becht GmbH
with its seat in Offenburg**

and Mr Ulrich Becht as its managing director, who is authorised, to represent the company alone.

Offenburg, Germany, 26.07.2023

**Notar Reinhard Körner
Civil Law Notary Public**

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