
Meril


Endo Surgery

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APPROVED

Meril Endo Surgery Pvt. Ltd.  
Additional General Manager-  
Regulatory Affairs  
Mr. Pratik Vasani

 /signature/  
stamp

"22" January 2024

## INSTRUCTIONS FOR USE

FOR MEDICAL DEVICE

**LINEAR ENDOSCOPIC STAPLE-CUTTING STAPLER WITH  
REPLACEABLE CASSETTES WITH STAPLES, DISPOSABLE, STERILE  
MIRUS™, VARIANT EXECUTIONS**


**Manufacturer** Meril Endo Surgery Pvt. Ltd.

Third Floor, E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala, Vapi, Valsad,  
Gujarat, 396191, India.

Version 1.1

**TRUE COPY**

  
**Ramesh Kumar Singh**  
Notary Public & Member of India  
Greater Mumbai & Thane District  
Regn. No. 11212

**Name of medical product:**

" Linear endoscopic stapling and cutting stapler with replaceable cartridges with staples, disposable, sterile MIRUS™, variant exwcutions"

**Manufacturer's name:** Meril Endo Surgery Pvt . Ltd. \_

**Address Manufacturer :** Third Floor, E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala , Vapi , Valsad , Gujarat, 396191, India.

**Phone:** 022-0260-2408000

**Email:** pratik.vasani@merillife.com

**Website:** www.merillife.com

**Developer Name:** Meril Endo Surgery Pvt . Ltd. \_

**Address Manufacturer :** Third Floor, E1-E3, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala , Vapi , Valsad , Gujarat, 396191, India.

**Name of the official representative in the Russian Federation:**

Limited Liability Company "Meril Medical "

**Address:** Russia, 108811, Moscow, vn.ter.g. Settlement Moskovsky, Kievskoe Highway 22nd (Moskovsky village) km, 4/4, room . 2/1.

**Phone:** +7 495 772 7643.

**Email:** tapas.sen@merillife.com

## 1. PURPOSE AND SCOPE OF APPLICATION

Linear endoscopic stapler- cutting stapler with replaceable cartridges with staples, disposable, sterile MIRUS<sup>TM</sup>, variant excisions designed for quick cutting /excision and suturing of soft and vascular tissues during open or minimally invasive surgical interventions, in abdominal surgery, during gynecological operations in the thoracic surgery for adults and children (hereinafter referred to as stapler (s), product).

Application area – Used only in medical institutions.

## 2. EQUIPMENT

Linear endoscopic stapling and cutting stapler with replaceable cartridges with staples, disposable, sterile MIRUS<sup>TM</sup>, variant excisions:

1. Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS<sup>TM</sup>, consisting of:

1.1 Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS<sup>TM</sup>

1.2 Three-row replaceable cassettes with MIRUS<sup>TM</sup> brackets:

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs.

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

1.3 Instructions for use 1 pc.

2. Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS<sup>TM</sup>, consisting of:

2.1 Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS<sup>TM</sup>

2.2 Three-row replaceable cassettes with MIRUS<sup>TM</sup> brackets:

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs.

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

2.3 Instructions for use 1 pc.

3. Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS<sup>TM</sup>, consisting of:

3.1 Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS<sup>TM</sup>

3.2 Three-row replaceable cassettes with MIRUS<sup>TM</sup> brackets:

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs.

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
  - MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
  - MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
  - MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).
- 3.3 Instructions for use 1 pc.

**4. Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS™, consisting of:**

**4.1 Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS™**

**4.2 Three-row replaceable cassettes with MIRUS™ brackets :**

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs.
- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

**4.3 Instructions for use 1 pc.**

**5. Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS™, consisting of:**

**5.1 Linear endoscopic stapling and cutting stapler MEC065 (short), 65 mm MIRUS™**

**5.2 Three-row replaceable cassettes with MIRUS™ brackets :**

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs.
- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

**5.3 Instructions for use 1 pc.**

**6. Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS™, consisting of:**

**6.1 Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS™**

**6.2 Three-row replaceable cassettes with MIRUS™ brackets :**

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs.
- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

6.3 Instructions for use 1 pc.

7. Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS<sup>tm</sup>, consisting of:

7.1 Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS<sup>tm</sup>

7.2 Three-row replaceable cassettes with MIRUS<sup>tm</sup> brackets :

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs.

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

7.3 Instructions for use 1 pc.

8. Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS<sup>tm</sup>, consisting of:

8.1 Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS<sup>tm</sup>

8.2 Three-row replaceable cassettes with MIRUS<sup>tm</sup> brackets :

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs.

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

8.3 Instructions for use 1 pc.

9. Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS<sup>tm</sup>, consisting of:

9.1 Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS<sup>tm</sup>

9.2 Three-row replaceable cassettes with MIRUS<sup>tm</sup> brackets :

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs.

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

9.3 Instructions for use 1 pc.

10. Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS<sup>tm</sup>, consisting of:

10.1 Linear endoscopic stapling and cutting stapler MEC160 (standard), 160 mm MIRUS™

10.2 Three-row replaceable cassettes with MIRUS™ brackets :

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs.

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

10.3 Instructions for use 1 pc.

11. Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS™, consisting of:

11.1 Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS™

11.2 Three-row replaceable cassettes with MIRUS™ brackets :

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs.

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

11.3 Instructions for use 1 pc.

12. Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS™, consisting of:

12.1 Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS™

12.2 Three-row replaceable cassettes with MIRUS™ brackets :

- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs.

- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

12.3 Instructions for use 1 pc.

13. Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS™, consisting of:

13.1 Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS™

13.2 Three-row replaceable cassettes with MIRUS™ brackets :

- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs.
  - MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
  - MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
  - MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
  - MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).
- 13.3 Instructions for use 1 pc.

14. Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS<sup>™</sup>, consisting of:

14.1 Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS<sup>™</sup>

14.2 Three-row replaceable cassettes with MIRUS<sup>™</sup> brackets :

- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs.
- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs. (if necessary).

14.3 Instructions for use 1 pc.

15. Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS<sup>™</sup>, consisting of:

15.1 Linear endoscopic stapling and cutting stapler MEC255 (extended), 255 mm MIRUS<sup>™</sup>

15.2 Three-row replaceable cassettes with MIRUS<sup>™</sup> brackets :

- MECRB60 (black), for fabrics with a thickness of 2.0-2.5 mm, row length of staples 60 mm, height of open staples 4.5, 4.5, 4.8 mm, up to 12 pcs.
- MECRW45 (white), articulated-vascular, row length of staples 45 mm, height of unclosed staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRP45 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 45 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).
- MECRW60 (white), for fabrics with a thickness of 1.0-1.5 mm, row length of staples 60 mm, height of open staples 2.6, 2.6, 3.0 mm, up to 12 pcs. (if necessary).
- MECRP60 (purple color), for fabrics with a thickness of 1.5-2.0 mm, the length of a row of staples is 60 mm, the height of an open staple is 3.5, 3.5, 4.0 mm, up to 12 pcs. (if necessary).

15.3 Instructions for use 1 pc.

### 3. INDICATIONS FOR USE

Open or minimally invasive general, gynecological, urological, thoracic, abdominal, pediatric, surgical interventions for rapid tissue cutting /excision and anastomosis.

### 4. CONTRAINDICATIONS

1. The endoscopic linear cutting instrument should not be used on tissue such as the liver or spleen where the compressibility is such that closing the instrument may cause fracture.
  2. Do not use the stapler on ischemic or necrotic tissue.
  3. Do not use the endoscopic linear stapler if surgical staples are contraindicated.
  4. Do not use the endoscopic linear stapler on inflamed, swollen mucous or tumor tissue.
  5. Do not use an endoscopic linear stapler if there is a suspicion that altered malignant tissue remains in the incision area.
- This stapler is not used if there are contraindications to the application of surgical staples.

#### **5. SIDE EFFECTS :**

- May cause bleeding or tissue damage
- Incomplete stitching/stitching failure/stitching error.

#### **Possible complications**

- Bleeding
- Vascular injury/hematoma formation
- Poor wound regeneration
- Allergic reactions

#### **Possible post-operative complications (see : Possible Complications Associated with Surgery):**

- Tissue trauma
- Pain/swelling
- Intestinal obstruction
- Problems with urination (bladder injury, retention)
- Infection/dehiscence
- Recurrent prolapse
- Foreign material
- Strangulated hernia

**6. USER** - specially trained medical personnel (doctor). The product should only be used by qualified personnel.

#### **7. RISK ASSESSMENT.**

##### **Risk of reuse:**

The product is intended for one time use only. Repeated use may lead to cross-contamination/infection and may contribute to non-fit for purpose. Consequently, the label indicates a prescription for single use only, and this information is also reflected in the instructions for use, to eliminate the possibility of this risk.

**Assessing the acceptability of the overall residual risk:**

- As a result of the risk management process, the frequency of possible risks for each risk event is reduced as much as possible and it is verified that risk control measures did not create additional risk.
- Major side effects are uncommon and are associated with the surgical process during and after use of the product.
- The overall residual risk after our risk control measures can be considered acceptable.

As a result of the above assessment, since the amount of acceptable risks is minimal and there is no unacceptable risk, it is concluded that stable risk management has been achieved.

## **8. SAFETY REQUIREMENTS**

4.1 Staplers are safe if the instructions for use are followed and do not have a harmful effect on the human body.

Work with the product is carried out strictly in accordance with the instructions for use.

4.2 Staplers are biologically safe and comply with ISO 10993-1 in terms of performance

- sensitization;
- irritation or intradermal reactivity;
- general toxicity;
- pyrogenicity .

4.3 Staplers must be equipped with a locking device that allows them to be removed from an organ or tissue only after the stitching has been completed.

The locking device does not allow operation if there are no staples in the cassette.

4.4. It is not allowed to use products if their individual packaging is damaged.

### **UNSUITABLE CRITERIA**

- ✓ Individual packaging damaged
- ✓ The stapler is damaged
- ✓ Cassette is damaged

## **9. OPERATING INSTRUCTIONS**

Staplers and cassettes must be used for the purposes established by this technical document, in strict accordance with the manufacturer's instructions for use.

The stapler must be resistant to climatic conditions, with air temperatures from -45 to +45 °C, air humidity up to 100% and relative air pressure from 84.0 kPa to 106.7 kPa.

When delivered in winter, to equalize the temperature of the products with the room temperature, they must be kept in the packaging at a temperature not lower than 18 ° C and relative air humidity (65 ± 5)% for at least 12 hours.

**Staplers and cassettes cannot be reused . Single use products. Dispose of after use.**

**Repeated sterilization is not permitted.**

**The implanted parts are titanium brackets, which do not have magnetic properties. MRI, CT and other types of scanning devices do not in any way affect the installed braces inside the human body. In rare cases, an increase in the temperature of titanium staples by less than 2°C° has been observed. With these examinations, the brackets are easily detected in the resulting images.**

**In case of suppuration or failure of the sutures, repeated surgical treatment is indicated. The need to remove, replace the implant or choose other surgical tactics is decided individually by the operating surgeon.**

**Once titanium staples have been correctly applied and the implant has been formed, no further maintenance of the implant is required.**

**Once their service life is complete, the braces remain with the patient for life or can be removed after healing during normal cell turnover or if the doctor deems it necessary. Implanted staples can be removed either during surgery or surgically after surgery. Permanent titanium alloy surgical braces are a well-established technology with extensive lifetime safety data.**

**Do not use the MECRW45 and MECRW60 3-row cassettes in the endoscopic linear stapler for any tissue with a layer thickness of less than 1.0 mm and more than 1.5 mm.**

**Do not use the MECRP45 and MECRP60 3-row cassettes in the endoscopic linear stapler for any tissue thickness less than 1.5 mm and more than 2.0 mm.**

**Do not use MECRB60 3-row cassettes in the endoscopic linear stapler for any tissue whose layer thickness is less than 2.0 mm and more than 2.5 mm.**

**Do not use the endoscopic linear stapler where the adequacy of hemostasis after use cannot be verified visually .**

**Do not use knife staplers on large vessels without ensuring proximal and distal control.**


**Before stitching with a stapler, it is necessary to carefully evaluate the thickness of the fabric. See Replacement Cassette Code Chart for requirements . If the fabric cannot be easily compressed to the height of the closed staple given in the table, or is easily compressed to a thickness less than the height of the closed staple given in the table, such fabric is not indicated for stitching with staples of the selected size, since it may be too thick or thin.**

## **HOW TO USE**

(Note: The data provided should not be used as technical guidance for anastomosis surgery.) When selecting the appropriate staple cartridge, always consider the overall thickness of the tissue and any staple line reinforcement material used.

### CHARGER

1. To load the endoscopic linear stapler with the appropriate cassette, insert the pin located on the distal end of the instrument shaft into the cassette. Make sure the alignment indicator on the cassette matches the alignment indicator on the stem. (Picture 1)
2. Press down on the cassette and turn it clockwise relative to the stapler to secure it. The alignment indicator on the stapler shaft will align with the alignment indicator on the cassette (Figure 2).
3. To ensure proper charging, cycle the stapler after installing the cassette. Squeeze the handle once to close the cassette jaws. Pull back the black return handles and make sure the cassette jaws are fully open (Figure 3)

Picture 1.


Figure 2.


Figure 3.

## DISCHARGE

1. To unload a cassette from the stapler , the lock lever must be in the neutral position. Make sure the cassette jaws are open by pulling the black return knobs all the way back. (Figure 4).
2. Pull the UNLOAD/UNLOCK button (located on the underside of the stapler shaft) back toward the stapler , rotate the cassette counterclockwise 45°, and remove the cassette from the stapler shaft (Figure 5)


Figure 4.


Figure 5.

#### **DIRECTIONS FOR USE .**


Insert the endoscopic linear stapler into a trocar sleeve of the appropriate size or larger.

Note: The Endoscopic Linear Stapler with White Violet Cassette is designed to be inserted and used through a trocar sleeve with a diameter of 12mm or larger (the sleeve is not supplied with the stapler ).


When using an endoscopic linear stapler with a black cassette, it must be inserted into a 15mm trocar sleeve (sleeve not supplied with the stapler ). The stapler can be reloaded and fired up to 12 times in one operation.

Technical characteristics of trocar sleeves:


Trocar sleeve, with an internal diameter of  $12 \pm 1$  mm.

Trocar sleeve, with an internal diameter of  $15 \pm 1$  mm.

1. The cassette jaws must be closed before inserting the instrument into the trocar sleeve. To do this, squeeze the handle (Fig. 6). Pull the UNLOAD/UNLOCK button (located on the underside of the shaft) back toward the instrument, rotate the cassette counterclockwise  $45^\circ$ , and remove the cassette from the instrument barrel. (Figure 6)


**Figure 6.**

2. With the stapler inserted into the body cavity, open the cassette jaws by pulling the black return handles all the way back. (Figure 7)

**Attention:**

1. Do not press down on the stapler handle while pulling back the black return handle.
2. Place the endoscopic linear stapler across the tissue being treated.

**Caution:** Make sure there are no obstructions (such as staples) in the stapler clamps. Processing through an obstruction may result in incomplete cutting and/or incorrectly formed staples.

**Figure 7.**

3. Adjust the lower swivel ring so that the stapler shaft rotates 360° (Figure 8). Rotate the top rotary ring so that the rotators rotate 35° in both directions. (Figure 9)


Figure 8.


Figure 9.


4. Close the jaws of the instrument across the tissue to be cut, squeezing the handle completely. The stapler is equipped with a safety lock; the device will not start the stapler or cut the fabric until the green button is pressed. (Figure 10).

**Attention:**

1. A safety interlock is provided to prevent the empty disposable cassette from being restarted.

Do not attempt to disable the safety interlock.

2. The jaws of the tool can be moved onto the fabric before starting by pulling the black return handle all the way back, allowing the jaws to open.

**Figure 10.**

3. To turn on the stapler, press the green start button. Squeeze the handle progressively until the oval clamp cap reaches the distal end of the cassette slot and the handle locks into place. To fully activate the cassette, successive presses on the handle are required (Figure 11).



Figure 11.

4. Once the tool is completely unloaded, open the clamps by pulling the black reset handle all the way back. Carefully remove the instrument from the fabric. After removing the instrument, the site of hemostasis should be checked. Minor bleeding can be controlled with electrocautery or hand sutures (Figure 12).

**Note:**

1. Do not attempt to insert or remove the instrument from the trocar hub while the instrument is in the articulated position.
2. This tool can be recharged and fired up to 12 times in one operation.

**Figure 12.**

**! WARNINGS AND PRECAUTIONS!**

- procedures should only be performed by surgeons who have undergone appropriate training and have residual knowledge in this field. Before performing any surgical intervention, it is recommended to familiarize yourself with the medical literature regarding the technique of performing the upcoming intervention, possible risks and complications.
- installation of replacement cartridges of the wrong type or size into the stapler can lead to tissue cutting without stitching
- during the procedure, do not try to lean on the front part of the cassette jaws, as this can lead to tissue damage and injury .

- when dissecting large vessels, follow the basic surgical principles of hemostasis control proximal and distal to the end of the intervention.
- make sure that the fabrics lie flat and are placed correctly between the cassette jaws. Any folding of the fabric along the replacement cassette, especially at the junction of the jaws, can cause the mechanical seam to fail.
- when applying sponges to the selected area, make sure that there are no obstacles between the sponges in the form of clips, stents, conductors, etc. Stitching in the presence of obstacles can lead to incomplete tissue dissection, incorrect formation of the mechanical staple suture and/or impossibility of opening the sponges.
- if the stitching mechanism malfunctions, do not attempt to use the stapler again.
- If the closing mechanism of the jaws is faulty and they do not grip the fabric, do not attempt to sew the fabric. Set aside and do not use this stapler.

#### **Clinical subsection:**

Examples of procedures in which staplers with replaceable staple cartridges can be used:

Table 1.

| General surgical                                                                                                       | Gynecological                                                            | Urological                                                                                                     | Thoracic                                                                                                             |
|------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>- Anastomoses</li> <li>- Partial gastropasty</li> <li>- Appendectomy</li> </ul> | <ul style="list-style-type: none"> <li>- Vaginal hysterectomy</li> </ul> | <ul style="list-style-type: none"> <li>- Prostatectomy</li> <li>- Nephrectomy</li> <li>- Cystectomy</li> </ul> | <ul style="list-style-type: none"> <li>- Lung resection and biopsy</li> <li>- Wedge resection of the lung</li> </ul> |

## **10.STORAGE AND SHELF LIFE**

Storage conditions:

Temperature: from -45 to +45 °C

Humidity up to 100%

Relative air pressure from 84.0 kPa to 106.7 kPa.

At a distance of at least 1 m from heating devices, and must be protected from contamination and exposure to aggressive environments.

Shelf life (shelf life) – 3 years from the date of sterilization.

The use of staplers after the expiration date is not allowed.

## **11. ENVIRONMENTAL PROTECTION REQUIREMENTS**

The materials of staplers and replacement cassettes do not have the ability to form toxic compounds in the air and wastewater. Solid waste generated during production must be processed or disposed of.

The main type of possible hazardous impact on the environment is pollution of the atmospheric air of populated areas, soils and waters as a result of

- emergency leaks (spills) of production materials;
- unorganized disposal of waste on the territory of the manufacturer or outside it;
- their arbitrary dumping in places not intended for this purpose.

## 12.DESIGN AND TECHNICAL REQUIREMENTS.


Figure 13.


Figure 14.

- |                           |                                                |
|---------------------------|------------------------------------------------|
| 1. Staple cassette        | 6. Feed sleeve                                 |
| 2. Cassette support       | 7. Rod                                         |
| 3. Shabot                 | 8. Adjustment knob                             |
| 4. Clamp cover            | 9. Upper turning part of the turning ring      |
| 5. Cassette rod           | 10. Lower rotating part of the rotating ring   |
| 11. Start button ( green) | 12. Return handle (groove indicated in black ) |
| 13. Fixed handle          | 14. Start lever                                |
| 15. Guide wedge           | 16. Brackets.                                  |
| 17. Knife                 |                                                |


Figure 15

Table 2.

| Name                              | Technical parameters |                           |
|-----------------------------------|----------------------|---------------------------|
|                                   | MECRW45 and MECRP45  | MECRW60, MECRP60, MECRB60 |
| Cassette with staples (l*w*h), mm | 74*10*6 ±1           | 89*11*7 ±1                |
| Shabot (l*w*h), mm                | 51*9*3.2 ±1          | 67*10*3.2 ±1              |
| Shabot elevation angle, °         | 13° ±1               | 13° ±1                    |
| Cassette support (l*w*h), mm      | 73*10*5 ±1           | 89*11*8 ±1                |
| Cutting angle, °                  | 29° ±1               | 29° ±1                    |
| Clamp cover (l*w*h), mm           | 24*11*9 ±2           | 24*13*10 ±2               |

|                         |            |            |
|-------------------------|------------|------------|
| Guide wedge (l*w*h), mm | 6*4*4 ±0.1 | 6*4*4 ±0.1 |
|-------------------------|------------|------------|



Figure 16

**Operating principle :** makes seams by installing a stacked row of titanium staples in a checkerboard pattern and at the same time separating the fabric between the rows.

Staplers are classified as a hand-held disposable surgical instrument , stitching, cutting, multi-staple , for obtaining multi-row linear sutures formed in a checkerboard pattern, equipped with a blade (knife) and a manual wedge-shaped drive.

Depending on the length of the working part, staplers are available in the following models:

**Figure 16.1**

Table 3.1 to Figure 16.1

| MODE<br>L  | DESCRIPTIO<br>N | TOTAL<br>LENGTH<br>OF<br>STAPLER<br>(WITHOUT<br>CASSETTE<br>) | STAPLE<br>R<br>HEIGHT | MAXIMU<br>M<br>NUMBER<br>OF<br>STARTS | ROD<br>LENGT<br>H | ROD<br>DIAMETER   | WEIGHT<br>WITHOUT<br>PACKAGIN<br>G (g) | PACKAG<br>E<br>WEIGHT<br>(g) | PACKAGING SIZE<br>(LENGTH*WIDTH*HEIGH<br>T) (mm)                                                                                      |
|------------|-----------------|---------------------------------------------------------------|-----------------------|---------------------------------------|-------------------|-------------------|----------------------------------------|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| MEC<br>065 | shortened       | 290±5mm                                                       | 150 ± 5<br>mm         | 12                                    | 65 ± 5<br>mm      | 12.40mm±0.3m<br>m | 235± 30                                | 405± 30                      | Individual packaging (blister<br>tray):<br><br>520*216*70 (±15 mm)<br><br>Secondary consumer<br>packaging:<br><br>544*217*75 (±15 mm) |
| MEC<br>160 | standard        | 380±5mm                                                       | 150 ± 5<br>mm         | 12                                    | 160 ± 5<br>mm     | 12.40mm±0.3m<br>m | 263± 30                                | 433± 30                      | Individual packaging (blister<br>tray):<br><br>520*216*70 (±15 mm)<br><br>Secondary consumer<br>packaging:                            |

|            |          |         |               |    |             |              |         |         |                                                                                                                                       |
|------------|----------|---------|---------------|----|-------------|--------------|---------|---------|---------------------------------------------------------------------------------------------------------------------------------------|
|            |          |         |               |    |             |              |         |         | 544*217*75 (±15 mm)                                                                                                                   |
| MEC<br>255 | extended | 480±5mm | 150 ± 5<br>mm | 12 | 255±5<br>mm | 12.40mm±0.3m | 290± 30 | 460± 30 | Individual packaging (blister<br>tray):<br><br>520*216*70 (±15 mm)<br><br>Secondary consumer<br>packaging:<br><br>544*217*75 (±15 mm) |

Table 3.2 to Figure 16.1

| MODEL | OVERALL<br>DIMENSIONS<br>OF THE<br>ADJUSTING<br>HANDLE<br>(L*W*H) | ROD<br>INLET<br>DIAMETER | OVERALL<br>DIMENSIONS<br>OF THE<br>ROTARY<br>PART OF THE<br>ROTARY<br>RING,<br>(L*W*H) | TRIGGER<br>BUTTON<br>DIAMETER | OVERALL<br>DIMENSIONS<br>OF RETURN<br>HANDLE,<br>(L*W*H) | OVERALL<br>DIMENSIONS<br>OF THE<br>RETURN<br>HANDLE<br>GROOVE,<br>(L*W) | OVERALL<br>DIMENSIONS<br>OF THE<br>FIXED<br>HANDLE,<br>(L*W*H) | OVERALL<br>DIMENSIONS<br>OF TRIGGER<br>LEVER,<br>(TOTAL<br>L*W) | INTERNAL<br>L*W |
|-------|-------------------------------------------------------------------|--------------------------|----------------------------------------------------------------------------------------|-------------------------------|----------------------------------------------------------|-------------------------------------------------------------------------|----------------------------------------------------------------|-----------------------------------------------------------------|-----------------|
|-------|-------------------------------------------------------------------|--------------------------|----------------------------------------------------------------------------------------|-------------------------------|----------------------------------------------------------|-------------------------------------------------------------------------|----------------------------------------------------------------|-----------------------------------------------------------------|-----------------|

|            |                  |          |                  |                |                  |           |                   |            |               |
|------------|------------------|----------|------------------|----------------|------------------|-----------|-------------------|------------|---------------|
| MEC<br>065 | 60*7*12 ±2<br>mm | Ø11±1 mm | 86*47*28±6<br>mm | Ø4.5±0.1<br>mm | 23*15*15±1<br>mm | 85*2±1 mm | 150*35*27±2<br>mm | 87*37±2 mm | 57*20±2<br>mm |
| MEC<br>160 | 60*7*12 ±2<br>mm | Ø11±1 mm | 86*47*28±6<br>mm | Ø4.5±0.1<br>mm | 23*15*15±1<br>mm | 85*2±1 mm | 150*35*27±2<br>mm | 87*37±2 mm | 57*20±2<br>mm |
| MEC<br>255 | 60*7*12 ±2<br>mm | Ø11±1 mm | 86*47*28±6<br>mm | Ø4.5±0.1<br>mm | 23*15*15±1<br>mm | 85*2±1 mm | 150*35*27±2<br>mm | 87*37±2 mm | 57*20±2<br>mm |

To equip staplers, disposable cartridges with staples should be used:

Table 4.

#### THREE-ROW CASSETTES FOR ENDOSCOPIC LINEAR STAPLER

| Cod<br>e        | Color | Staple<br>row<br>length<br>(mm) and<br>descriptio<br>n | Cut<br>length(mm<br>) | Stitchi<br>ng line<br>length(<br>mm) | Numb<br>er of<br>staple<br>s | Height of<br>open<br>bracket | Length of<br>open bracket<br>(mm) | Height of<br>bracket<br>in closed<br>position<br>(mm) | Numb<br>er of<br>rows | Weig<br>ht<br>witho<br>ut<br>packa<br>ging<br>(g) | Pack<br>age<br>weig<br>ht (g) | Package size<br>(length*width*height)<br>(mm) |
|-----------------|-------|--------------------------------------------------------|-----------------------|--------------------------------------|------------------------------|------------------------------|-----------------------------------|-------------------------------------------------------|-----------------------|---------------------------------------------------|-------------------------------|-----------------------------------------------|
| ME<br>CR<br>W45 | White | Articular-<br>vascular                                 | 40 ±5 mm              | 45 ±2<br>mm                          | 66                           | 1st row:<br>2.6±0.1mm,       | 1st row:<br>2.6±0.2mm,            | 1.0±0.3m<br>m                                         | 2 x 3<br>rows         | 60 ±5                                             | 92 ±<br>10                    | Individual packaging<br>(blister tray):       |

|                 |        |                                                                    |          |             |    |                                                                           |                                                                           |               |               |       |            |                                                                                                                                    |
|-----------------|--------|--------------------------------------------------------------------|----------|-------------|----|---------------------------------------------------------------------------|---------------------------------------------------------------------------|---------------|---------------|-------|------------|------------------------------------------------------------------------------------------------------------------------------------|
|                 |        | Cassette<br>45±2 mm                                                |          |             |    | 2nd row:<br>2.6±0.1mm,<br>3rd row:<br>3.0±0.1mm                           | 2nd row:<br>2.6±0.2mm,<br>3rd row:<br>3.0±0.2mm                           |               |               |       |            | 265*72*20 (±5 mm)<br><br>Secondary consumer<br>packaging:<br><br>280*105*80 (±5 mm)                                                |
| ME<br>CRP<br>45 | Violet | Cassette<br>for<br>medium<br>and thick<br>fabrics 45<br>±2 mm      | 40 ±5 mm | 45 ±2<br>mm | 66 | 1st row:<br>3.5±0.1mm,<br>2nd row:<br>3.5±0.1mm,<br>3rd row:<br>4.0±0.1mm | 1st row:<br>3.5±0.2mm,<br>2nd row:<br>3.5±0.2mm,<br>3rd row:<br>4.0±0.2mm | 1.6±0.3m<br>m | 2 x 3<br>rows | 60 ±5 | 92 ±<br>10 | Individual packaging<br>(blister tray):<br><br>265*72*20 (±5 mm)<br><br>Secondary consumer<br>packaging:<br><br>280*105*80 (±5 mm) |
| ME<br>CR<br>W60 | White  | Cassette<br>for<br>medium<br>tissues<br>and<br>vessels 60<br>±2 mm | 55 ±5 mm | 60 ±2<br>mm | 90 | 1st row:<br>2.6±0.1mm,<br>2nd row:<br>2.6±0.1mm,<br>3rd row:<br>3.0±0.1mm | 1st row:<br>2.6±0.2mm,<br>2nd row:<br>2.6±0.2mm,<br>3rd row:<br>3.0±0.2mm | 1.0±0.3m<br>m | 2 x 3<br>rows | 62 ±5 | 98 ±<br>10 | Individual packaging<br>(blister tray):<br><br>265*72*20 (±5 mm)<br><br>Secondary consumer<br>packaging:<br><br>280*105*80 (±5 mm) |
| ME<br>CRP<br>60 | Violet | Cassette<br>for<br>medium<br>and thick                             | 55 ±5 mm | 60 ±2<br>mm | 90 | 1st row:<br>3.5±0.1mm,<br>2nd row:<br>3.5±0.1mm,                          | 1st row:<br>3.5±0.2mm,<br>2nd row:<br>3.5±0.2mm,                          | 1.6±0.3m<br>m | 2 x 3<br>rows | 62 ±5 | 98 ±<br>10 | Individual packaging<br>(blister tray):<br><br>265*72*20 (±5 mm)                                                                   |

|                 |       |                                              |          |             |    |                                                                           |                                                                           |               |               |       |            |                                                                                                                                      |
|-----------------|-------|----------------------------------------------|----------|-------------|----|---------------------------------------------------------------------------|---------------------------------------------------------------------------|---------------|---------------|-------|------------|--------------------------------------------------------------------------------------------------------------------------------------|
|                 |       | fabrics 60<br>±2 mm                          |          |             |    | 3rd row:<br>4.0±0.1mm                                                     | 3rd row:<br>4.0±0.2mm                                                     |               |               |       |            | Secondary consumer<br>packaging:<br><br>280*105*80 ( ±5 mm)                                                                          |
| ME<br>CRB<br>60 | Black | Cassette<br>for thick<br>fabrics 60<br>±2 mm | 55 ±5 mm | 60 ±2<br>mm | 90 | 1st row:<br>4.5±0.1mm,<br>2nd row:<br>4.5±0.1mm,<br>3rd row:<br>4.8±0.1mm | 1st row:<br>4.5±0.2mm,<br>2nd row:<br>4.5±0.2mm,<br>3rd row:<br>4.8±0.2mm | 2.0±0.3m<br>m | 2 x 3<br>rows | 82 ±5 | 98 ±<br>10 | Individual packaging<br>(blister tray):<br><br>265*72*20 ( ±5 mm)<br><br>Secondary consumer<br>packaging:<br><br>280*105*80 ( ±5 mm) |

The stapler places titanium staples in a staggered pattern and simultaneously separates the fabric between rows.

The width of the cutting edge of the cassette knife should be no more than 8 microns, the height of irregularities on the surface of the knife should be no more than 10 microns.

The radii of blunting of the working part of the brackets are at least 0.3 mm.

The surface of the staplers should be free of cracks, nicks, scratches, burrs, shells, pores, chipped areas, delaminations, burns and other contaminants (scale, particles of grinding, polishing materials and traces of lubricant) .

The knife blade should be sharp without nicks, cracks, crumbles or dents.

External surfaces should be shiny or matte .

Soldered and welded seams (if any) must be tight and free of cracks, holes, and traces of uncooked material . In welds, individual pores with a depth of up to 0.2 mm and an area of up to 0.08 mm<sup>2</sup> are allowed in an amount of no more than 3 in a section no more than 5 mm long .

There should be no non-functional sharp edges in the design.

The surfaces of polymer components should not have kinks, nicks, dents, scratches, burrs, crumbled areas, foreign inclusions, dirt and stains visible to the naked eye.

The finger rests must not have burrs or sharp edges.

The cutting edges must be sharp along their entire length.

The stapler must be resistant to climatic conditions, with air temperatures from -45 to +45 °C, air humidity up to 100% and relative air pressure from 84.0 kPa to 106.7 kPa.

The cartridge with staples and the stapler rod must be resistant to aggressive biological fluids (blood, urine, gastric juice).

The connection of the components must be strong.

The possibility of spontaneous separation of components is not allowed.

Surface roughness of cutting edges  $R_a$  should be in the range from 0.05 to 0.63 microns. The roughness parameter  $R_a$  of the remaining surfaces of the staplers should be no more than:

- 0.32 microns – external;

The force required to carry out stitching with simultaneous cutting should be no more than:

- 200 N (20 kgf) – transmitted by the hand;

- 100 N (10 kgf) – transmitted by fingers.

The force required to install the cassette should be no more than 35 N (3.57 kgf).

When stitching, the staples should be completely pushed out of the grooves.

The staples must not fall out of the staple slots under the influence of their own weight. The brackets should not fall out when swinging.

The brace should form a closed "B" loop shape when fired.

The cassette should not fall out of the socket under the influence of its own weight unless the need for falling out is due to the design and operational features of the stapler .

The cassette lock must ensure that the pusher enters the magazine grooves without jamming.

movement of the jaws should be easy and smooth.

The movement of the stapler should be gentle, without blocking.

The start button of staplers must be painted green.

Staplers must be supplied sterile (sterilization is carried out using the gas method using ethylene oxide in accordance with ISO 11135, with confirmation of a sterility level (SAL) of  $10^{-6}$ ).

The design of the cassettes must ensure the interchangeability of their models within the appropriate standard size and with a standard feed sleeve diameter of 12 mm.

The stapler must be made of the following materials:

Table 5.

| Part name          | Material         | Brand or manufacturer of material | Type of contact with the human body                 |
|--------------------|------------------|-----------------------------------|-----------------------------------------------------|
| Lock lever         | Titanium,        | 1 UNSR50250                       | No contact                                          |
| Stock              | Titanium         | 1 UNSR50250                       | Brief contact (less than 24 hours) with soft tissue |
| Adjustment knob    | ABS+SS plastic   | BM Precision Co., Ltd.            | No contact                                          |
| Top swivel ring    | ABS+SS plastic   | BM Precision Co., Ltd.            | No contact                                          |
| Bottom swivel ring | ABS+SS plastic   | BM Precision Co., Ltd.            | No contact                                          |
| Start button       | ABS+SS plastic   | BM Precision Co., Ltd.            | No contact                                          |
| Return handle      | ABS+SS plastic   | BM Precision Co., Ltd.            | No contact                                          |
| Fixed handle       | ABS+SS plastic   | BM Precision Co., Ltd.            | No contact                                          |
| Start button       | ABS+SS plastic   | BM Precision Co., Ltd.            | No contact                                          |
| Start lever        | ABS+SS plastic . | BM Precision Co., Ltd.            | No contact                                          |

Dyes used in the manufacture of staplers : black, white, green, manufactured by Waston International Co.Ltd .

The cassette and brackets must be made of the following materials:

Table 6.

| Part name | Material and grade of material              | Type of contact with the human body                                                  |
|-----------|---------------------------------------------|--------------------------------------------------------------------------------------|
| Brackets  | Titanium grade 1 UNSR50250                  | Implanted products that have long-term contact (more than 30 days) with soft tissues |
| Cassette  | Titanium, grade 1 UNSR50250, ABS+SS plastic | Brief contact (less than 24 hours) with soft tissue                                  |

Dyes used in the manufacture of plastic parts of cassettes: black, white, purple produced by Waston International Co.Ltd . ).

#### MATERIAL CHARACTERISTICS ( Titanium , grade 1 UNSR50250):

Chemical composition (with percentage of detailed chemical composition)

Main material: Titanium grade 1 UNSR50250

Table 7.

| TYPE AND HARDNESS OF METAL |          |                      |               |
|----------------------------|----------|----------------------|---------------|
| Brand 1 UNSR50250          |          | Description of metal | Hardness , HB |
| 1 UNSR50250                | Titanium | Pure titanium        | 140 (± 10 )   |

#### MECHANICAL PROPERTIES OF MATERIAL Titanium , grade 1 UNSR50250:

Yield strength 0.2%, MPa – 170...310

Tensile strength, MPa – 240

Relative elongation, % – 24

**PHYSICAL CHARACTERISTICS:**

- All surfaces of the product must be free of sharp edges and burrs.
- All molded parts must not have an uneven surface.
- The entire surface should be free of pitting, cracks, etc.

|                                                      |       |
|------------------------------------------------------|-------|
| Metal density, g/cm <sup>3</sup>                     | 4.51  |
| Young's modulus, GPa                                 | 146   |
| Specific heat capacity, J/kg*K                       | 520   |
| Linear expansion coefficient, 10 <sup>-6</sup> m/ mK | 8.15  |
| Thermal conductivity, W/m*K                          | 18.85 |
| Electrical resistance, Ω* s m*10 <sup>-6</sup>       | 45    |

Properties at elevated temperatures: This metal has good resistance to scale formation up to 700°C.

Individual packaging of the stapler (subject to sterilization) must be made of the following materials:

Table 8.

| Package name | Material and grade of material                            | Technical specifications                              |
|--------------|-----------------------------------------------------------|-------------------------------------------------------|
| Blister tray | paper brand TYVEK 1059 B                                  | Thickness 750 microns                                 |
|              | PET produced by Waston International Co., Ltd.            | Thickness deviation 33%                               |
|              | Internal label: MIRRORCOAT self-adhesive paper            | Punching strength - 360±100 kPa                       |
|              | Label siliconized and mirror self-adhesive paper SILVASSA | Bonding strength 4.4<br>N/2.54 cm - 19.0<br>N/2.54 cm |

|  |                                      |  |
|--|--------------------------------------|--|
|  | <b>FLEXIBLE<br/>PRINTERS PVT LTD</b> |  |
|--|--------------------------------------|--|

Secondary consumer packaging of the stapler is made of 5- layer corrugated box , produced by Pekij Industries .

Individual packaging of staple cartridges (subject to sterilization) must be made of the following materials:

Table 9.

| <b>Package name</b> | <b>Material and grade of material</b>                                                      | <b>Technical specifications</b>                        |
|---------------------|--------------------------------------------------------------------------------------------|--------------------------------------------------------|
| <b>Blister tray</b> | <b>paper brand TYVEK 1059 B</b>                                                            | <b>Thickness 750 microns</b>                           |
|                     | <b>PET produced by Waston International Co., Ltd.</b>                                      | <b>Thickness deviation 33%</b>                         |
|                     | <b>Internal label: MIRRORCOAT self-adhesive paper</b>                                      | <b>Punching strength - 360±100 kPa</b>                 |
|                     | <b>Label siliconized and mirror self-adhesive paper SILVASSA FLEXIBLE PRINTERS PVT LTD</b> | <b>Bonding strength 4.4 N/2.54 cm - 19.0 N/2.54 cm</b> |

Secondary consumer packaging with a staple cartridge is a 5-ply cardboard box PRODUCED by Pekij Industries .

Materials must retain their properties after sterilization and must be resistant to the biological environment of the body.

The product does not contain materials of human or animal origin.

Titanium staples have high wear resistance and biological inertness. Materials that come into contact with the human body do not have a negative effect on the human body.

The design of consumer packaging must ensure tightness, prevent the penetration of microorganisms and changes in the toxicological properties of products, and also ensure protection of products from external mechanical damage.

Individual packaging should not change the toxicological properties of products.

### **13. MARKING**

on the individual packaging of the stapler includes:

- name and address of the manufacturer indicating the symbol according to GOST R ISO 15223-1;
- manufacturer's trademark;
- name of the stapler ;
- graphic image of a stapler ;
- symbol "batch code" ( LOT );
- catalog number ( REF );
- execution option;
- production date and corresponding symbol;
- the inscription "STERILE" and the sterilization method;
- the inscription "Non-toxic";
- the inscription "Pyrogen-free" and the corresponding symbol;
- symbol "Prohibition of re-use";
- expiration date and "use by" symbol ;
- symbol "refer to instructions for use";
- "do not re-sterilize" symbol;
- symbol "keep away from moisture";
- symbol "do not expose to sunlight";

- “do not use if packaging is damaged” symbol;
- sign “EU Marking and identification number of the notified body”;
- “Caution!” symbol;
- “temperature limit” symbol;
- operating, storage and transportation conditions;
- manipulation sign “Top”;
- QR code;
- license number ( Lic . No );
- sign and data “Official Representative in the European Union”;
- the inscription “for export only”;
- number of products in the package;
- information about disposal;
- number and date of issue of the registration certificate of the Russian Federation.

The additional label on the stapler contains information including:

- - name of the manufacturer and address of the manufacturer indicating the symbol in accordance with GOST R ISO 15223-1;
- manufacturer's trademark;
- name of the stapler in accordance with the registration certificate;
- execution option;
- symbol “batch code” ( LOT );
- catalog number ( REF );
- production date and corresponding symbol;

- expiration date and “best before” symbol ;
- the inscription “STERILE” and the sterilization method.

on the secondary consumer packaging of the stapler includes:

- name and address of the manufacturer indicating the symbol in accordance with GOST R ISO 15223-1;
- manufacturer's trademark;
- name of the stapler ;
- graphic image of a stapler ;
- symbol “batch code” ( LOT );
- catalog number ( REF );
- execution option;
- production date and corresponding symbol;
- the inscription “STERILE” and the sterilization method;
- the inscription “Non-toxic”;
- the inscription “Pyrogen-free” and the corresponding symbol;
- symbol “Prohibition of re-use”;
- expiration date and “use by” symbol ;
- symbol “refer to instructions for use”;
- “do not re-sterilize” symbol;
- symbol “keep away from moisture”;
- symbol “do not expose to sunlight”;
- symbol “number of products per package”;

- “do not use if packaging is damaged” symbol;
- sign “EU Marking and identification number of the notified body”;
- “Caution!” symbol;
- manipulation sign “Top”;
- “temperature limit” symbol;
- operating, storage and transportation conditions;
- QR code;
- license number ( Lic . No );
- sign and data “Official Representative in the European Union”.
- the inscription “for export only”;
- number of products in the package;
- information about disposal;
- number and date of issue of the registration certificate of the Russian Federation.

Individual packaging of staple cartridges includes information including:

- name and address of the manufacturer indicating the symbol in accordance with GOST R ISO 15223-1;
- manufacturer’s trademark;
- name of the medical product;
- name of the cassette with brackets and version;
- symbol “batch code” ( LOT );
- catalog number ( REF );
- production date and corresponding symbol;

- the inscription "STERILE" and the sterilization method;
- the inscription "Non-toxic";
- the inscription "Pyrogen-free" and the corresponding symbol;
- symbol "Prohibition of re-use";
- expiration date and "use by" symbol ;
- symbol "refer to instructions for use";
- "do not re-sterilize" symbol;
- symbol "keep away from moisture";
- symbol "do not expose to sunlight";
- "do not use if packaging is damaged" symbol;
- sign "EU Marking and identification number of the notified body";
- "Caution!" symbol;
- "temperature limit" symbol;
- operating, storage and transportation conditions;
- QR code;
- license number ( Lic . No );
- sign and data "Official Representative in the European Union".
- the inscription "for export only";
- number of products in the package;
- information about disposal;
- number and date of issue of the registration certificate of the Russian Federation.

**The additional label on the staple cartridges includes information including:**

- name of the manufacturer and address of the manufacturer indicating the symbol in accordance with GOST R ISO 15223-1;
- manufacturer's trademark;
- name of the medical product;
- name of the cassette with brackets and version;
- symbol "batch code" ( LOT );
- catalog number ( REF );
- production date and corresponding symbol;
- expiration date and "best before" symbol ;
- the inscription "STERILE" and the sterilization method;

**Information on the secondary consumer packaging of staple cartridges includes:**

- name and address of the manufacturer indicating the symbol in accordance with GOST R ISO 15223-1;
- manufacturer's trademark;
- name of the medical product;
- name of the cassette with brackets and version;
- symbol "batch code" ( LOT );
- catalog number ( REF );
- production date and corresponding symbol;
- the inscription "STERILE" and the sterilization method;
- the inscription "Non-toxic";
- the inscription "Pyrogen-free" and the corresponding symbol;

- symbol “Prohibition of re-use”;
- expiration date and “use by” symbol ;
- symbol “refer to instructions for use”;
- “do not re-sterilize” symbol;
- symbol “keep away from moisture”;
- symbol “do not expose to sunlight”;
- symbol “number of products per package”;
- “do not use if packaging is damaged” symbol;
- sign “EU Marking and identification number of the notified body”;
- “Caution!” symbol;
- “temperature limit” symbol;
- operating, storage and transportation conditions;
- QR code;
- license number ( Lic . No );
- sign and data “Official Representative in the European Union”.
- the inscription “for export only”;
- number of products in the package;
- information about disposal;
- number and date of issue of the registration certificate of the Russian Federation.

Explanation of graphic symbols:

Table 9.

| Symbol | Symbol meaning |
|--------|----------------|
|--------|----------------|

|                                                                                     |                                                    |
|-------------------------------------------------------------------------------------|----------------------------------------------------|
|    | Carefully!                                         |
|    | Do not reuse                                       |
|    | Manufacturer                                       |
|    | Keep away from moisture                            |
|    | Do not re-sterilize                                |
|    | date of manufacture                                |
|    | Do not expose to sunlight                          |
|    | Batch code                                         |
|    | Catalog number                                     |
|    | Refer to instructions for use                      |
|    | Do not use if packaging is damaged                 |
|    | Number of products in a stapler package            |
|   | Number of products in a cassette package           |
|  | EC marking and notified body identification number |
|  | Sell by                                            |
|  | Temperature limit                                  |
|  | Sterilized using ethylene oxide                    |

|                                                                                   |                                               |
|-----------------------------------------------------------------------------------|-----------------------------------------------|
|  | Official representative in the European Union |
|  | Non-pyrogenic                                 |
|  | Top                                           |

#### 14. PACKAGING , STERILIZATION AND PRODUCTION PROCESS

Each stapler is packed in a sterile blister (individual packaging), which protects the product from any mechanical influences.

The integrity of the packaging is ensured by durable thermal sealing - weld seam width  $15 \pm 2$  mm.

Sealed with heat seal at  $105 \pm 2^\circ\text{C}$ .

The product must be placed in the appropriate packaging slot and secured in it during transportation.

Each cassette with staples is packed in a sterile blister (individual packaging), which protects the product from any mechanical influences.

The integrity of the packaging is ensured by durable thermal sealing - weld seam width  $15 \pm 2$  mm.

The product must be placed in the appropriate packaging slot and secured in it during transportation.

The package with a stapler and the package with a cartridge with staples, together with instructions for use, must be placed in a cardboard box ( Document No. QC/FMT/STA/004) .

For transportation, preserved products in boxes (Document No. QC/FMT/STA/004) with operational documentation must be placed in boxes made of sheet wood materials or boards. Wooden boxes must be lined inside with waterproof material.


Packaging must provide:

- tightness;
- maintaining the sterility of the contents when stored in dry, clean and well-ventilated conditions;
- exclusion of penetration of microorganisms;
- protection of staplers from external mechanical damage;
- no risk of contamination of the contents during opening and removal from the package;
- reliable protection of the contents under normal conditions of transportation and storage;
- impossibility of resealing an opened package: the fact that the package has been opened must be obvious.


Individual packaging should not change the toxicological properties of packaged staplers . Migration of paints in individual packaging is not allowed.

**Stages of the production process.**

# Technological diagram ( Stapler ):

## Technological diagram (replaceable cassettes):

## **15. REPAIR INFORMATION**

The product cannot be repaired and must be disposed of after use.

## **16. STORAGE AND TRANSPORTATION IN THE TERRITORY OF THE RF**

Transportation is carried out by any type of covered transport, provided they are protected from pollution and mechanical damage, in accordance with the rules for the carriage of goods in force for this type of transport.

Staplers in the manufacturer's packaging must have

- vibration resistance when exposed to a vibration load with an amplitude of 0.26 mm and a frequency of 40 Hz;
- impact resistance when exposed to a shock load with an acceleration of 50 m/s<sup>2</sup>, a shock pulse duration of 10-12.5 ms and a blow frequency of 10-50 per minute;
- resistance to mechanical stress during transportation and withstand shaking tests for 1 hour with a rooting of 30 m/s<sup>2</sup> with a number of vibrations of 2-3 per second.

Transportation conditions:

Temperature: from -45 to +45 °C

Humidity up to 100%

Relative air pressure from 84.0 kPa to 106.7 kPa.

Storage conditions:

Temperature: from -45 to +45 °C

Humidity up to 100%

Relative air pressure from 84.0 kPa to 106.7 kPa.

Staplers must be stored at least at a distance 1 m from heating devices , and must be protected from contamination and exposure to aggressive environments.

## **17.INFORMATION ABOUT DISPOSAL IN THE TERRITORY OF THE RF**

After use, the products are disposed of in medical institutions as waste class B , and the products With expired for a period storage And Not former V use – By class And SanPiN 2.1.3684-21.

## **18.MANUFACTURER WARRANTY**

The manufacturer guarantees that the products comply with the requirements of this technical document, subject to the rules of use, transportation and storage.

Guaranteed shelf life ( shelf life ) is 3 years from the date of sterilization.

The use of staplers after the expiration date is not allowed.

Staplers in undamaged consumer containers must be pyrogen-free , non-toxic and sterile during the guaranteed shelf life.

If a malfunction is detected in the product that arose during the warranty period and led to a violation of their characteristics, the product will be replaced by the manufacturer free of charge.

### **Information about complaints**

All requests and complaints should be sent to:

**Name of the official representative in the Russian Federation:**

Limited Liability Company "Meril Medical "

**Address:** Russia, 108811 , Moscow, vn.ter.g. Settlement Moskovsky, Kievskoe Highway 22nd (Moskovsky village) km, 4/4, room . 2/1.

**Phone:** +7 495 772 7643.

**Email:** [tapas.sen@merillife.com](mailto:tapas.sen@merillife.com)

## **19. Brief review of clinical data.**

An assessment of all clinical literature has demonstrated the effectiveness of endoscopic staplers in various operations in both adults and children.

As stated in the clinical evaluation report ( CLE / REP /015 Rev.06 ), the device meets the essential safety requirements according to the essential requirements of ER1. Clinical literature has shown that staplers have proven to be safe in various surgeries such as abdominal, gynecological, gastrointestinal, pulmonary and thoracic surgery for resection, dissection and anastomosis, etc. in both children and adults patients. In addition, various preclinical physical, performance, and biocompatibility tests were conducted , confirming the safety and effectiveness of the device.

The risk-benefit profile was acceptable according to the core requirement of ER1. The risk analysis was carried out in accordance with the EN ISO 14971:2012 standard. A risk-benefit analysis shows that the benefits outweigh the risks.

## **Example of a warranty card**

### **WARRANTY CARD**

Product category \_\_\_\_\_

Manufacturer \_\_\_\_\_

Execution \_\_\_\_\_ option \_\_\_\_\_ Stamp \_\_\_\_\_

Serial number/art. \_\_\_\_\_ about sale \_\_\_\_\_

Date of sale \_\_\_\_\_

Seller organization \_\_\_\_\_

Seller's signature \_\_\_\_\_

Individual packaging has been checked and is not damaged. I have read and agree with the instructions for use.

Buyer's signature \_\_\_\_\_

Guaranteed shelf life ( shelf life ) is 3 years from the date of sterilization.

By issues guarantees Contact : Limited Liability Company "Meril Medical "

**Address:** Russia, 108811, Moscow, vn.ter.g. Settlement Moskovsky, Kievskoe Highway 22nd (Moskovsky village) km, 4/4, room . 2/1.

**Phone:** +7 495 772 7643.

**Email:** [tapas.sen@merillife.com](mailto:tapas.sen@merillife.com)

***ATTENTION! Make sure that the warranty card is filled out correctly. The serial number of the product must exactly match the number on the coupon. If the coupon is filled out inaccurately or incompletely, the product warranty is considered lost.***

**Staplers and replacement cartridges with staples must comply with the requirements established by the manufacturer in conjunction with the applicable standards:**

| <b>No.</b> | <b>Standard name</b>                                                                                                                                     |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1          | EN ISO 13485: Medical Devices – Quality Management System Requirements for Regulatory Purpose                                                            |
| 2          | EN 1041 - Information supplied by the manufacturer of medical devices                                                                                    |
| 3          | EN ISO 15223-1 : Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements |
| 4          | EN ISO 10993-1 : Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process                              |
| 5          | EN ISO 10993-3 : Biological evaluation of medical devices – Part 3: Test for Genotoxicity , Carcinotoxicity and Reproductive toxicity                    |
| 6          | EN ISO 10993-4 : Biological evaluation of medical devices – Part 4: Selection of Test for Interactions with Blood                                        |
| 7          | EN ISO 10993-5: Biological evaluation of medical devices – Part 5: Tests for In-Vitro Cytotoxicity                                                       |
| 8          | EN ISO 10993-6: Biological evaluation of medical devices – Part 6: Tests for local effects after implantation                                            |
| 9          | EN ISO 10993-7: Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residues                                                 |
| 10         | ISO 10993-10: Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization                                            |
| eleven     | EN ISO 10993-11: Biological evaluation of medical devices – Part 11: Tests for systemic toxicity                                                         |
| 12         | EN ISO 14971: Medical devices - Application of risk management to medical devices                                                                        |
| 13         | EN ISO 11607-1: Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems  |

|        |                                                                                                                                                                                             |
|--------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 14     | EN ISO 11607-2: Packaging for terminally sterilized medical devices -Part 2: Validation requirements for forming, sealing and assembly processes                                            |
| 15     | EN ISO 11135: "Medical Devices – Validation & Routine Control of Ethylene Oxide Sterilization"                                                                                              |
| 16     | EN ISO 5832-2: Implants for surgery -- Metallic materials -- Part 2: Unalloyed titanium                                                                                                     |
| 17     | EN ISO 5832-1: Implants for surgery — Metallic materials — Part 1: Wrought stainless steel                                                                                                  |
| 18     | ASTM F67: Standard Specification for Unalloyed Titanium, for Surgical Implant Applications                                                                                                  |
| 19     | ASTM F138: Standard Specification for Wrought 18Chromium-14 Nickel - 2.5 Molybdenum Stainless Steel Bar and Wire for Surgical Implants (UNS S31673)                                         |
| 20     | EN ISO 14644-1: Clean rooms and associated controlled environments -- Part 1: Classification of air cleanliness by particle concentration                                                   |
| 21     | EN ISO 14644-2: Clean rooms and associated controlled environments - Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration   |
| 22     | EN ISO 11737-1 : Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products (ISO 11737-1:2006)                        |
| 23     | EN ISO 11737-2 : Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process |
| 24     | MEDDEV 2.7/1 , revision 4 Evaluation of Clinical Data, A guide for Manufacturers and Notified Bodies.                                                                                       |
| 25     | Medical Device Directive 93/42/EEC: Medical devices directive as amended by DIRECTIVE 2007/47/EC                                                                                            |
| 26     | ASTM F 1980 : Standard Guide for Accelerated aging of sterile barrier systems for Medical Devices                                                                                           |
| 27     | ASTM F 88 : Standard test methods for seal strength of flexible barrier materials                                                                                                           |
| 28     | ICH Q1A(R2) Stability Testing of New Drug Substances and Products                                                                                                                           |
| 29     | ASTM D4169: Standard Practice for Performance Testing of Shipping Containers and Systems.                                                                                                   |
| thirty | ASTM D5276: Standard Test Method for Drop Test of Loaded Containers by Free Fall                                                                                                            |

|    |                                                                                |
|----|--------------------------------------------------------------------------------|
| 31 | ASTM D999: Standard Test Methods for Vibration Testing of Shipping Containers. |
| 32 | EN ISO 14630 - Inactive surgical implants - General requirements               |

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| <b>№</b> | <b>Наименование стандарта</b>                                                                                                                           |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1        | EN ISO 13485: Medical Devices – Quality Management System Requirements for Regulatory Purpose                                                           |
| 2        | EN 1041 - Information supplied by the manufacturer of medical devices                                                                                   |
| 3        | EN ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements |
| 4        | EN ISO 10993-1: Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process                              |
| 5        | EN ISO 10993-3: Biological evaluation of medical devices – Part 3: Test for Genotoxicity, Carcinotoxicity and Reproductive toxicity                     |
| 6        | EN ISO 10993-4: Biological evaluation of medical devices – Part 4: Selection of Test for Interactions with Blood                                        |
| 7        | EN ISO 10993-5: Biological evaluation of medical devices – Part 5: Tests for In-Vitro Cytotoxicity                                                      |
| 8        | EN ISO 10993-6: Biological evaluation of medical devices – Part 6: Tests for local effects after implantation                                           |
| 9        | EN ISO 10993-7: Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals                                               |
| 10       | ISO 10993-10: Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization                                           |
| 11       | EN ISO 10993-11: Biological evaluation of medical devices – Part 11: Tests for systemic toxicity                                                        |
| 12       | EN ISO 14971: Medical devices - Application of risk management to medical devices                                                                       |
| 13       |                                                                                                                                                         |

|    |                                                                                                                                                                                                   |
|----|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <b>EN ISO 11607-1: Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems</b>                                    |
| 14 | <b>EN ISO 11607-2: Packaging for terminally sterilized medical devices -Part 2: Validation requirements for forming, sealing and assembly processes</b>                                           |
| 15 | <b>EN ISO 11135: “Medical Devices – Validation &amp; Routine Control of Ethylene Oxide Sterilization”</b>                                                                                         |
| 16 | <b>EN ISO 5832-2: Implants for surgery -- Metallic materials -- Part 2: Unalloyed titanium</b>                                                                                                    |
| 17 | <b>EN ISO 5832-1: Implants for surgery — Metallic materials — Part 1: Wrought stainless steel</b>                                                                                                 |
| 18 | <b>ASTM F67: Standard Specification for Unalloyed Titanium, for Surgical Implant Applications</b>                                                                                                 |
| 19 | <b>ASTM F138: Standard Specification for Wrought 18Chromium-14 Nickel - 2.5 Molybdenum Stainless Steel Bar and Wire for Surgical Implants (UNS S31673)</b>                                        |
| 20 | <b>EN ISO 14644-1: Clean rooms and associated controlled environments -- Part 1: Classification of air cleanliness by particle concentration</b>                                                  |
| 21 | <b>EN ISO 14644-2: Clean rooms and associated controlled environments - Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration</b>  |
| 22 | <b>EN ISO 11737-1: Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products (ISO 11737-1:2006)</b>                        |
| 23 | <b>EN ISO 11737-2: Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process</b> |
| 24 | <b>MEDDEV 2.7/1, revision 4 Evaluation of Clinical Data, A guide for Manufacturers and Notified Bodies.</b>                                                                                       |
| 25 | <b>Medical Device Directive 93/42/EEC: Medical devices directive as amended by DIRECTIVE 2007/47/EC</b>                                                                                           |
| 26 | <b>ASTM F 1980: Standard Guide for Accelerated aging of sterile barrier system for Medical Devices</b>                                                                                            |
| 27 | <b>ASTM F 88: Standard test methods for seal strength of flexible barrier materials</b>                                                                                                           |
| 28 | <b>ICH Q1A(R2) Stability Testing of New Drug Substances and Products</b>                                                                                                                          |
| 29 | <b>ASTM D4169: Standard Practice for Performance Testing of Shipping Containers and Systems.</b>                                                                                                  |

|    |                                                                                  |
|----|----------------------------------------------------------------------------------|
| 30 | ASTM D5276: Standard Test Method for Drop Test of Loaded Containers by Free Fall |
| 31 | ASTM D999: Standard Test Methods for Vibration Testing of Shipping Containers.   |
| 32 | EN ISO 14630 - Inactive surgical implants - General requirements                 |

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