71e9-5421

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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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To see the Russian witness on the next page.

China Council for the Promotion
of International Trade

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证明 网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Manufacturer:

Baisheng Medical Co., Ltd

No. 11 Fusheng Road, Xinhui District, Jiangmen City, Guandong, China

Phone: +86-750-6628113

Fax: +86-7506616122

www.obs-medical.com

Authorized deputy in Russia:

Limited Liability Company

Emergency Medical Systems (EMS LLC)

Office 1, 1 Sormovskoe shosse str., 603074 Nizhny Novgorod, Nizhny Novgorod Oblast, Russia

Phone: +7 (831) 217-00-03

E-mail: info@emrmedsys.com

CAUTION!

Please, inform the manufacturer and/or its authorized deputy about all adverse incidents involving the product.

This Instruction for Use is applicable to Defibrillation Electrodes OBS-DE (hereinafter referred to as the product) of the following models: Defibrillation Electrode OBS-DE/P (electrodes for adults) and Defibrillation Electrode OBS-DE/W (electrodes for children).

CAUTION!

Please, pay special attention to such warnings in the text hereof, in order to avoid potential dangerous situations or operating errors.

Content

1 Symbols	5
2 Description	6
2.1 Intended Use	6
2.2 Intended Users	6
2.3 Intended Patients	6
3 Intended Use	7
3.1 Indications for Use	7
3.2 Contraindications for Use	7
3.3 Potential Side Effects	7
3.4 Precautions for Handling	7
3.5 Operating Procedure	8
4 Characteristics	10
4.1 Specification	10
4.2 Materials	11
4.3 Compatibility with Defibrillator Models	11
4.4 Product parts	13
4.5 Packaging Details	13
4.6 Product disposal	13

1 Symbols

The following symbols are used in the product labelling:

— see Instruction for Use

— Lot No.

— Reference No.

— Manufacturer

— Date of manufacture

— Best before

— Do not reuse

— Do not use when the package is damaged

— Non-sterile

— Avoid exposure to sunlight

— Temperature range

— Latex-free

— Warning

2 Description

2.1 Intended Use

Defibrillation Electrodes are also called Multifunction Electrode Pads, which are accessories for defibrillator devices. It is applied on the patient's chest and will be used for administering cardiac defibrillation shock (defibrillation and cardioversion), monitoring and for pacing when medically indicated.

2.2 Intended Users

The product is designed to be used by health care professionals and specially trained people.

2.3 Intended Patients

Depending on the model, the product is designed for various types of patients in accordance with the labelling on the primary package (Figure 2.1):

- Defibrillation Electrode OBS-DE/P is designed for patients of over 8 years old and weight of over 25 kg.
- Defibrillation Electrode OBS-DE/W is designed for patients of less than 8 years old and weight of less than 25 kg.


Figure 2.1 — Patient-type labelling

3 Intended Use

3.1 Indications for Use

The product shall be used as a part of a defibrillator, for urgent and emergency medical aid in case of acute cardiovascular diseases (events), as medically indicated.

3.2 Contraindications for Use

There are no contraindications to the product application.

3.3 Potential Side Effects

Potential side effects include skin reddening (which is normal) or burns (as a result of improper use).

3.4 Precautions for Handling

When handling the product, observe the following precautions:

- Do not use after expiration date.
- Do not crush, bend or fold the electrodes. Do not store them under heavy objects.
- Ensure that the patient's chest is clean and dry.
- Do not open the primary package until immediately prior to use.
- Before use, inspect the package, the electrodes, the cable, connectors and check that the product storage conditions are met. Do not use the product in case of any non-compliances.
- Use the product within 24 h after the package is opened.
- Do not use the product if gel is dry.
- Do not use additional gel on the electrodes.
- Do not re-use the product. If re-used, the electrical properties might be insufficient, which could lead to the patient's injury/hazard. Re-use also poses cross-infection risk for the patients.
- If long-term monitoring is required, use separate ECG electrodes and cable.
- Always keep a second pack of defibrillation electrodes with the defibrillator.
- Do not use the product for magnetic-resonance imaging (MRI).
- The operators shall avoid contact between parts of the patient's body, such as exposed skin of head or limbs, with conductive fluids, such as gel, blood or saline, and metal objects, such as a bed frame or stretchers, which may provide undesirable defibrillating-current conduction pathways.
- The electrodes used shall be replaced if any of the following conditions are met:
 - a) 50 defibrillation shocks are given;
 - b) pacing for 8 hours: 100 mA/80 bpm/40 ms;
 - c) the electrode is attached to the patient's skin for 24 hours.

CAUTION!

When handling the product, remember that:

Failure to comply with the instructions might result in the patient's burns or ineffective treatment.

- *Strangulation hazard: route the electrode cable away from the patient's throat.*
 - *The product shall be kept well clear of other electrodes or metal parts in contact with the patient.*
 - *Extensive pacing may lead to air bubbles. To ensure optimum contact with the patient's skin, smooth out any bubbles by hand.*
 - *Other medical devices that have no defibrillation-proof applied parts (without a relevant labelling) shall be disconnected from the patient during defibrillation.*
 - *Do not reposition the electrodes and do not touch the patient during defibrillation.*
 - *Keep the product clear of any other electrodes or metal parts in contact with the patient. Do not overlap the electrodes.*
-


3.5 Operating Procedure

- 1** **Select the electrode application area**
Select the method for electrode application on the patient's body according to the Instruction for Use for the defibrillator.
- 2** **Prepare the patient's skin (if required)**
Check the skin in the electrode application area: it shall be clean, dry, without hair. Shave any hair and wipe the skin, if required (Figure 3.1 a, b).
- 3** **Unpack the electrodes**
Remove the electrodes from the primary package (Figure 3.1 c) and remove protective film from the electrodes (Figure 3.1 d).
- 4** **Apply the electrodes**
Apply the electrodes on the patient's body and press them down firmly to ensure proper contact of the entire adhesive surface with the skin (Figure 3.1 e).
- 5** **Connect the electrodes**
Connect the electrodes' connector to the defibrillator used (Figure 3.1 f).
- 6** **Use according to the instruction**
Proceed in accordance with instructions of the Instruction for Use for the defibrillator used.
Note — Do not leave the product on the patient's skin for more than 24 h (Figure 3.1 g).

7

Remove the electrodes

Remove the electrodes by slowly peeling them from the patient's skin. Dispose of the electrodes in accordance with the applicable national regulations for disposal of medical waste.


a) Shave hair


b) Wipe the skin


c) Remove the electrodes


d) Remove protective film

e) Press the electrodes down

f) Connect the electrodes

g) Apply for no more than 24 h

Figure 3.1 — Operating procedure

4 Characteristics

4.1 Specification

The product specification is provided in the Table 4.1.

Table 4.1 — Product specification

Parameter	Value
1 Models	Defibrillation Electrode OBS-DE/P; Defibrillation Electrode OBS-DE/W.
2 Type	Non-reusable.
3 Reference Nos.	Determine the connector type and compatibility (refer to Table 4.3).
4 Composition	The product consists of five layers: - a protective film; - conductive gel; - a conductive metal foil; - a polyethylene foam backing (with acrylic adhesive); - a PET foil (with acrylic adhesive). The product materials coming into contact with the human body are not toxic (refer to Table 4.2).
5 Biocompatibility	The product is biocompatible and successfully passed cytotoxicity, skin irritation and sensitization evaluation according to GOST R IEC 60601-2-4-2013.
6 Weight	no more than 0.145 kg.
7 Dimensions	(165 × 112) mm — Defibrillation Electrode OBS-DE/P; (130 × 100) mm — Defibrillation Electrode OBS-DE/W. Cable length: (1.8±0.2) m; Cable outer diameter: (3±0.1) mm. Active (gel) area of the electrodes complies with requirements of GOST R IEC 60601-2-4-2013.
8 Potential risk of use class	Class I According to the order of the Ministry of Health of the Russian Federation of June 6, 2012 N 4n.
9 Operating conditions	Range of ambient temperatures: 0 °C to 35 °C; Relative air humidity: no more than 80 %; Atmospheric pressure: 50 kPa to 106 kPa.
10 Transportation and storage conditions	Range of ambient temperatures: 0 °C to 35 °C.
11 Electrical parameters	The products comply with requirements of GOST R IEC 60601-2-4-2013 and were tested in terms of - AC small signal impedance; - AC large signal impedance; - Offset instability and internal noise; - Defibrillation overload recovery; - DC offset voltage; - Insulation strength; - ECG signal recovery.
12 Mechanical parameters	The products comply with requirements of GOST R IEC 60601-2-4-2013 and were tested in terms of the electrode adhesion and contact with the patient.
13 Shelf life	2 years from the date of manufacture (refer to the package).

4.2 Materials





The materials coming into short-term contact with the human's intact skin are not toxic.




4.3 Compatibility with Defibrillator Models

The product is designed for use with the following defibrillator models (Table 4.3):

Note — The combination of OBS-DE defibrillation electrodes and a defibrillator forms a medical system.

Table 4.2 — Compatibility with defibrillator models

Reference	Connector design	Defibrillator model or brand
303A5030		PHILIPS HeartStart MRx (with cable M3508A), HeartStart XL+, HeartStart FRx, Efflicia DFM100
303A5035		MINDRAY BeneHeart C1 Fully Automatic, BeneHeart C1A Fully Automatic, BeneHeart C2 Fully Automatic, BeneHeart C2A Fully Automatic, BeneHeart C1, BeneHeart C1A, BeneHeart C2, BeneHeart C2A, BeneHeart D3, BeneHeart D6, uMED 20, uMED 20A PHYSIO-CONTROL LifePak 12, LifePak 20 Monitor, LifePak CR Plus, LifePak 15, LIFEPAK 20e, LifePak 1000, LifePak 500 COMEN S3, S5, S6, S8
303A4530		NIHON KOHDEN TEC-5521K (with adapter JC-755V), TEC-5531K (with adapter JC-755V), TEC-5621 (with adapter JC-755V), TEC-5631 (with adapter JC-755V)
303A2545		ZOLL M Series, M Series CCT, E Series, R Series, X Series
303A6030		ZOLL AED Plus, AED Pro
303A1525		CARDIAC SCIENCE PowerHeart AED G3, PowerHeart AED G3 Pro, PowerHeart AED G3 Automatic, FirstSave AED G3

Reference	Connector design	Defibrillator model or brand
303A2060		METRAX PRIMEDIC HeartSave PAD (SN up to 73944xxxxxx), PRIMEDIC HeartSave AED (SN up to 73945xxxxxx), PRIMEDIC HeartSave AED-M (SN up to 73946xxxxxx), PRIMEDIC HeartSave 6 (SN up to 73947xxxxxx), PRIMEDIC HeartSave AS (SN up to 73943xxxxxx), PRIMEDIC DefiMonitor
303A4531		NIHON KOHDEN TEC-5521K (with adapter JC-765V), TEC-5531K (with adapter JC-765V), TEC-5621 (with adapter JC-765V), TEC-5631 (with adapter JC-765V), TEC-7621K (with adapter JC-765V), TEC-7631 (with adapter JC-765V)
303A7035		ALTOMEDICA ALTDEF BIOSS AVD-1
303A4525		SCHILLER DEFIGARD 4000, DEFIGARD 5000, DEFIGARD 6002, FRED easy
303B1525		PROGETTI RESCUE SAM AXION DA-N-01, DA-N-02, DA-N-05, DKI-N-11

CAUTION!

Use the product only with the specified defibrillator models provided herein. Use of the product with any devices other than listed may result in potentially dangerous connection.

4.4 Product parts

The product parts (delivery set) are given in the tables 4.3 and 4.4.

Table 4.3 — Product parts for the Defibrillation Electrode OBS-DE/P

Part name	Manufacturer	Country
Defibrillation Electrode OBS-DE/P	Baisheng Medical Co., Ltd	China
Instruction for use *	Baisheng Medical Co., Ltd	China
* By agreement with the customer, the Instruction for use can be provided in a single copy for a group of products.		

Table 4.4 — Product parts for the Defibrillation Electrode OBS-DE/W

Part name	Manufacturer	Country
Defibrillation Electrode OBS-DE/W	Baisheng Medical Co., Ltd	China
Instruction for use *	Baisheng Medical Co., Ltd	China
* By agreement with the customer, the Instruction for use can be provided in a single copy for a group of products.		

4.5 Packaging Details

A pouch is used as a primary package of the product.

A box is used as the secondary package of the product.

A carton is used as a tertiary package of the product.

The product is packed as follows: 1 pc/pouch, 25 pouches/box, 4 boxes/carton (100 pouches/carton).

4.6 Product disposal

Dispose of the product in accordance with the national regulations for the disposal of medical waste in force at the time of disposal.

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FLINGU-YOU TRADING CO., LTD.
翻 译 公 司
FOR TRANSLATION ONLY
This is a true and correct translation
of the original document

本附加证明书证明文书上的署名、签署人职务的真实性。请核对文书上的印章属实。附加证明书不对文书内容予以证明。
This Apostille only certifies the authenticity of the signature, the capacity of the person who has signed the public document and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued.
Visit <https://consular.mfa.gov.cn/VERIFI> and scan the QR code to verify the issuance of this Apostille.

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本公文书 This public document

2. 签署人
has been
signed by 孙嘉
Sun Jia

3. 签署人身份
acting in the
capacity of 签证员
Authorized Official

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7. 签发人
by 李勇, 参赞/Li Yong, Counsellor
外交部领事司
Department of Consular Affairs of the Ministry of Foreign Affairs

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