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此公文書
This public document

2. 簽署人
has been signed by 沈怡君
Shen, Yi-Chun

3. 簽署人職務
acting in the capacity of 民間公證人
Notary public

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9. 章戳
Seal/stamp: 10. 簽署
Signature:

Shih, I-Chen

Secretary, Bureau of Consular Affairs

For The Minister of Foreign Affairs


11. 附註：
remarks:

本文件證明書僅證明簽名、章戳之真實性，並不證明文件內容之真實性。
This document authentication only certifies the authenticity of the signature, seal or stamp
and the capacity of the person who has signed the attached document. It does not validate the
contents of the document for which it was issued.

本文件證明書之驗證網址：
To verify the issuance of this authentication,
please peruse the following website:
<https://docauth.boca.gov.tw/BCCAWeb/index4/grv>

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


INSTRUCTION FOR USE

Bioport Crosslink Intra-Articular Injection Prothesis

Legal manufacturer
Maxigen Biotech Inc., Taiwan

Approved by :

Po-Hong Lai
Vice President/ Management Representative
Maxigen Biotech Inc.

1. Name of the Medical Device

Bioport Crosslink Intra-Articular Injection Prothesis, including:

- one disposable syringe without a needle containing a 2% solution of sodium hyaluronate, volume 3.0 ml – 1 pc.; - instructions for use – 1 pc.; - patient stickers – 2 pcs.

2. Legal Manufacturer and Authorized Representative Details

Legal manufacturer:

Full name: Maxigen Biotech Inc. Abbreviated name: MBI

Address (location): No. 88, Keji 1st Rd., Guishan Dist., Taoyuan City 33383, Taiwan

Phone: +88633287222, Fax: +88633287333

E-mail: contact.us@mbi.com.tw

Place of manufacture of the medical device:

Maxigen Biotech Inc., No. 88, Keji 1st Rd., Guishan Dist., Taoyuan City 33383, Taiwan

Authorized representative of the manufacturer in the Russian Federation:

Full name: Bioport Limited Liability Company Abbreviated name: Bioport LLC

Address (location): Office 408-1, Floor 4, 26V bldg 1, Bolshaya Pochtovaya Str., Basmanny Municipal Dist. (intracity territory of the federal city), 105082, Moscow

Phone: 8 (800) 222-38-02

E-mail: bioport@inbox.ru

Holder (owner) of the registration certificate:

Name: Nikovari Bio Tech Oy

Address: Retkeilijänkatu 5, 00980 Helsinki, Finland

Phone: +358449213972

E-mail: info@bioport.fi

3. Medical Device Description and Mechanism of Action

3.1 General Description

Bioport Crosslink Intra-Articular Injection Prothesis (hereinafter, the following names of the medical device can be used: Bioport Crosslink, Bioport Crosslink device, the device) is a sterile non-pyrogenic device for restoration of viscosity of synovial fluid. The device is based on viscose water solution of a highly purified crosslinked sodium hyaluronate fraction.

Hyaluronate is a natural polymer of the glucosaminoglycan family, a mucopolysaccharide consisting of repeating disaccharide units of N-acetylmannosamine and glucuronate interlinked by β 1- 3 and β 1-4 glycosidic bonds. Each pre-filled syringe contains 20 mg/ml (2.0%) of sodium hyaluronate dissolved in buffered physiological saline solution. The viscoelastic solution has a pH of 6.8-7.5, osmolarity of 270-330 mOsm/kg and viscosity is $> 126 \pm 16.19$ Pa, the degree of reticulation is $3.67 \pm 1.15\%$, the molecular weight of hualuronic acid is 30MDa, the amount of bacterial endotoxins is < 0.05 EE/mg.

3.2 Intended Use. Indication

Bioport Intra-Articular Injection Prothesis is intended for the treatment of pain caused by knee osteoarthritis in patients having not shown adequate response to a conservative non-pharmacological treatment, as well as to mild analgesics, such as acetaminophen. The effect of Bioport involves increasing the viscoelasticity of synovial fluid, reducing inflammatory response, protecting the joint tissues, as well as reducing osteoarthritis-caused pain when administered by intra-articular injections.

3.3 Method of Administration

Bioport Crosslink is administered by a single intra-articular injection using a standard procedure. The maximum possible volume the procedure is 3.0 ml.

All aseptic regulations must be followed. Joint effusion (if any) must be removed before the administration of Bioport Crosslink. The depth of injection into the joint cavity is determined

by the doctor individually. The biodegradation period of the product is 6 months. The duration of the clinical effect is 6 months.

3.4 Potential Users and Application

The device is intended for use only by qualified traumatologists or surgeons.

Application area: Traumatology and surgery.

4. Medical Device Application Risks, Contraindications, Expected and Predictable Side Effects Related to the Intended Use of the Medical Device

- Do not prescribe for patients with known hypersensitivity to hyaluronate preparations.

Contraindications

- Intra-articular injections are contraindicated in case of skin pathologies in the injection site.

No other contraindications for the intended use of Bioport Crosslink Intra-Articular Injection Prothesis found thus far.

Warnings and precautions

- The Bioport Crosslink device is intended only for intra-articular administration.
- Do not use in combination with disinfecting agents containing quaternary ammonium salts, as reaction with them may cause settling of hyaluronate.
- Do not use Bioport Crosslink if the packaging is opened or damaged.
- The device content is sterile if the packaging is not opened or damaged. Bioport Crosslink must be used immediately after the package opening. The content of the syringe is manufactured with the use of the aseptic filling procedure. Do not re-sterilize. The syringe is intended for single use only. Any Bioport Crosslink devices not used properly must be disposed. Do not reuse.

- Strictly comply with the aseptic injection procedure.
- Use Bioport Crosslink with caution in patients with severe forms of allergy.
- Use Bioport Crosslink with caution in case of venous or lymphatic stasis in the leg in which the injection is to be made.
- The intra-articular injection of Bioport Crosslink may be followed by transient and swelling of the relevant joint. The physician may prescribe the relevant pharmacological treatment.
- Joint effusion (if any) must be removed before the administration of Bioport Crosslink. Inject Bioport Crosslink into the joint through a sterile needle of gauge 20 (20G). Do not use the same syringe for removal of joint effusion and injection of Bioport Crosslink. Carefully remove the syringe plug and needle cap in compliance with the aseptic procedure.
- As it is the case for any other invasive procedure on the joint, it is recommended to avoid overloading the joint immediately after the intra-articular injection.
- It is also advisable to avoid the procedure when the patient is unable to undergo a medical follow-up during the next 48 hours.
- The safety and performance of repeated Bioport Crosslink treatment courses are not established.
- The safety and performance of the use of Bioport Crosslink on joints other than the knee, as well as for the treatment of pathologies other than the osteoarthritis, are not established.
- The safety and performance of the use of Bioport Crosslink in combination with other intra-articular injections are not established.
- It is prohibited to use Bioport Crosslink in pregnant or lactating women, as well as in children.
- Strictly comply with the procedures for maintaining the injection site sterility and use the proper injection method.
- The in vitro laboratory tests showed opalescence of sodium hyaluronate when mixing with other solutions containing quaternary ammonium compounds. Therefore, sodium hyaluronate should not be injected through a device previously used for injection of the said solutions.

Adverse reaction

Sodium hyaluronate is a natural substance present in human tissues. It is usually well-tolerated by the human. Some reports indicate transient post-intervention inflammatory response. Additionally, some reports indicate development of adverse reactions in some patients, such as pain in the injection site, knee swelling and/or joint effusion, local cutaneous reactions (rash, ecchymosis) and itching. A number of reports indicate temporary aggravation of inflammatory response in the treated knee in some patients after the injection in case of inflammatory arthritis, such as rheumatoid or gouty arthritis.

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TPD

5. Medical Device Specifications


The device is a syringe consisting of the barrel, hub (Luer Lock type), plug, plunger and plunger handle and prefilled with sodium hyaluronate solution. The syringe is used for injecting the sodium hyaluronate solution into the joint tissues.

Composition:

Each ml contains:

Sodium Hyaluronate - 20 mg;

Phys. Sodium phosphate buffer - q.s.


AB-N305	Bioport Crosslink Intra-Articular Injection Prothesis, including: - one disposable syringe without a needle containing a 2% solution of sodium hyaluronate, volume 3.0 ml – 1 pc.; - instructions for use – 1 pc.; - patient stickers – 2 pcs.
---------	--

Technical specifications and drawings:


Syringe:

- Barrel
3.0 ml glass barrel, ungraduated

Drawing

Description	Translation
Luer Lok taper 6%	Luer Lok taper 6%
1	Barrel outer wall
2	Barrel inner wall
3	Luer Lok taper
4	Plug

*The provided drawings indicate the dimensions in mm.

Technical specifications:

Total length from the plug to the flange: 91.10 mm (+1.20 mm / -2.10 mm)
Length from the taper (Luer Lock type) to the flange: 84.30 mm (± 1.10 mm)
Length from the barrel base to the flange: 72.20 mm (± 0.50 mm)
Taper diameter: 11.10 mm maximum
Male Luer Lock taper type

Attachment dimensions:

Outer diameter: 8.0 mm ± 0.1 mm
Inner diameter: 3.925 mm (+0 / -0.005) mm
Barrel outer diameter: 10.85 mm (± 0.10 mm)
Barrel inner diameter: 8.65 mm (± 0.20 mm)
Flange "wing" diameter: 8.25 mm minimum
Thumb rest thickness: 1.50 mm minimum
Thumb rest diameter 1: 14.70 mm (± 0.65 mm)
Thumb rest diameter 2: 17.75 mm (± 0.75 mm)

Plug:

Plug diameter d1: 4.0 mm maximum
Plug diameter d2: 6.5 mm maximum
Plug diameter d3: 8.0 mm maximum
Plug flange diameter: 2.0 mm maximum
Plug length: 14.0 mm maximum
Maximum "dead" space: 0.07 mL

• **Plunger handle**

Technical specifications:

Total plunger handle length: 74.50 mm (± 0.80 mm)
Plunger handle diameter: 8.0 mm (± 0.2 mm)
Plunger length: 7.70 mm (± 0.40 mm)
Plunger diameter: 9.15 mm (± 0.15 mm)
Plunger thread turns: 3 turns
Plunger flange diameter: 12.70 mm (± 0.40 mm)
Plunger flange thickness: 1.5 mm (± 0.2 mm)
Seal thickness: 1.2 mm (± 0.2 mm)
Labeling: BD (Becton Dickinson)

Syringe activation force:

Bioport Crosslink Intra-Articular Injection Prothesis: ≤ 25 N.

Syringe filling:

Sodium hyaluronate content: 20 mg/mL.

Recoverable volume: ≥ 3.0 mL/syringe.

6. Transportation, Storage and Application Conditions

Transportation conditions

Transportation temperature: from 15°C to 25°C

Relative humidity: 65% maximum

Atmospheric pressure: 70 – 106 kPa

Storage conditions

Storage temperature: from 15°C to 25°C

Relative humidity: 65% maximum

Atmospheric pressure: 70 – 106 kPa

Store in a cool, dry place, protected from direct sunlight.

Use environment

The device is intended to be used indoors at the ambient temperature: from 15°C to 25°C

Relative humidity: 65% maximum.

Atmospheric pressure: 86 – 106 kPa

7. Listing and Description of the Medical Device Materials to Come into Direct or Indirect Contact with the Human (Patient's) Body

Name	Type of human body contact	Manufacturer	Number of registration certificate/ regulatory document/ CAS number
BD Hypak glass syringes, sterile and non-sterile, with and without a needle, volume from 0.5 mL to 20 mL	Short-term contact (<24 h) with intact skin	Becton Dickinson France S.A.S, France	Federal Health Care Reference Book entry 2011/11237 dated 26.12.2011
Sodium Hyaluronate	Permanent (> 30 days) contact with bone tissue and internal	Bloomage Freda Biopharm Co., Ltd., China	Federal Reference Book entry 001305 dated 11.01.2016

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TPO

	environment		
Sodium chloride	Permanent (> 30 days) contact with bone tissue and internal environment	TAIWAN BIOTECH CO.,LTD	7647-14-5
Dibasic sodium phosphate (heptahydrate)	Permanent (> 30 days) contact with bone tissue and internal environment	Sigma-Aldrich, USA	7782-85-6
Sodium dihydrogen phosphate (dihydrate)	Permanent (> 30 days) contact with bone tissue and internal environment	Sigma-Aldrich, USA	10049-21-5
Water for injection	Permanent (> 30 days) contact with bone tissue and internal environment	TAIWAN BIOTECH CO.,LTD	7732-18-5


Bioport Crosslink Intra-Articular Injection Prothesis is a sterile medical device.

Handling of the device is allowed only with the use of personal protective equipment (gloves).

8. Medical Device Labeling and Packaging Information

8.1 Labeling


The labels are attached to the syringe, individual packaging (blister), patient's stickers, device cardboard box and transport container.


Example of the labeling of the Bioport Crosslink Intra-Articular Injection Prothesis cardboard box (side 1)

Example of the labeling of the Bioport Crosslink Intra-Articular Injection Prothesis cardboard box (side 2)


Example of the labeling of the Bioport Crosslink Intra-Articular Injection Prothesis cardboard box (lateral side 1)

Example of the labeling of the Bioport Crosslink Intra-Articular Injection Prothesis

MBI 和康生物科技股份有限公司
Maxigen Biotech Inc.
cardboard box (lateral side 2)

Example of the labeling of the Bioport Crosslink Intra-Articular Injection Prosthesis cardboard box (end face)


Translation of labels and description of symbols:

Text / Symbol	Translation
Nikovari Bio Tech Oy, Finland	Nikovari Bio Tech Oy, Finland
Bioport Crosslink	Bioport Crosslink
Intra-Articular Injection Prosthesis	Intra-Articular Injection Prosthesis
Ingredient: one single syringe without needle with 2% sodium hyaluronate solution	Composition: one single syringe without needle with 2% sodium hyaluronate solution.
3.0 ml	3.0 mL
	Catalogue number
	Batch code
	Manufacture date
	Expiry date
	Composition
	Temperature range from 15°C to 25°C
	Avoid exposure to sunlight
	Prohibition on reuse
	Consult instructions for use
	Sterilization using aseptic processing techniques

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
Example of the labeling of the Bioport Crosslink Intra-Articular Injection Prosthesis individual packaging

Translation of labels and description of symbols:

Text / Symbol	Translation
Bipoort Crosslink	Bipoort Crosslink
Intra-Articular Injection Prosthesis 2.0%	Intra-Articular Injection Prosthesis 2.0%
	Prohibition on reuse
	Consult instructions for use
	Sterilization using aseptic processing techniques

	Do not use in case of packaging damage
Aseptic filling Non-pyrogenic	Aseptic filling Non-pyrogenic
REF	Catalogue number
LOT	Batch code
	Expiry date
contents	Composition
Do not freeze	Do not freeze
Store at room temperature	Store at room temperature
Protected from light	Protected from light
Open	Open

Example of the Bioport Crosslink Intra-Articular Injection Prosthesis patient's stickers


Translation of labels and description of symbols:

Text / Symbol	Translation
Bioport Crosslink	Bioport Crosslink
Intra-Articular Injection Prosthesis 2.0%	Intra-Articular Injection Prosthesis 2.0%
REF	Catalogue number
LOT	Batch code
	Expiry date
contents	Composition

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3.0 ml


3.0 mL

Example of the labeling of the Bioport Crosslink Intra-Articular Injection Prothesis

Translation of labels and description of symbols:

Text / Symbol	Translation
Bioport Crosslink	Bioport Crosslink
Intra-Articular Injection Prothesis 2.0%	Intra-Articular Injection Prothesis 2.0%
	Catalogue number
	Batch code
	Expiry date
	Composition
3.0 ml	3.0 mL

Example of the labeling of the transport container

Symbols marked on the transport container:


Symbol	Translation
	Recyclable packaging

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TPO

	Fragile, handle with care
	Protect from moisture

Labeling in Russian provided for the circulation of the medical device in the Russian Federation:


8.2 Packaging Information

The individual package contains one single-use syringe without a needle, containing sodium hyaluronate solution. A locking plate (TAIRIPRO Homo grade) is placed on the thumb rest to prevent accidental pushing of the plunger.


The syringe is placed into the Tyvek-sealed blister tray (PET-GAG, CAS number 25038-59-9). The blister tray, along with the instructions for use (Simili paper, 220 g/cm) and patient's stickers (2 pcs, paper 0.1 IT) to be affixed in the patient's medical record, is placed into the cardboard box (Carta solida 45T cardboard).



200 (Two hundred) cardboard boxes are placed into the transport container (Corrugated cardboard BF K5xK5).

General view of the Bioport Crosslink Intra-Articular Injection Prothesis packaging contents

Packaging contents specifications:

Packaging contents	Material	Photo	Dimensions
Thumb rest lock	TAIRIPRO Homo grade		3.8 cm (± 0.2 cm) x 2.0 cm (± 0.2 cm) Weight: 1.23 g (± 0.10 g)
Individual packaging	Tyvek-sealed blister tray (PET-GAG, CAS number 25038-59-9)	 	18.0 cm (± 0.5 cm) x 5.4 cm (± 0.5 cm) x 2.3 cm (± 0.5 cm) Weight: 7.60 g (± 0.50 g)
Cardboard box	Carta solida cardboard 45T		18.0 cm (± 0.5 cm) x 2.5 cm (± 0.5 cm) x 6.0 cm (± 0.5 cm) Weight: 14.51 g (± 0.4 g)

Patient's stickers	Paper 0.11T		Size: 6.4 cm ($\pm$ 0.5 cm) x 2.5 cm ($\pm$ 0.5 cm) Weight: 0.38 g ($\pm$ 0.1 g)
Instruction for use	Simili paper, 220 g/cm		Size per item: 30 cm x 21 cm x 2 cm maximum Weight per item: 250 g maximum

Transport container specifications

Transport container	Corrugated cardboard BF K5xK5		Dimensions (L x W x H): (61 cm x 40 cm x 33 cm) maximum Weight: 1 kg maximum
---------------------	-------------------------------	--	---

Tyvek 1059B coating specification:

PROPERTIES	Standard	UOM	VALUE
Density	ISO 536	g/m ²	71.6-80.4
Thickness	ASTM F2251	μ m	163
Average porosity	ISO 5636-3	mL/min	>152
Tensile load	ASTM D5035	lb/inch (N/cm ²)	36.6 (25.23)
			39.2 (27.03)
Breaking strength	ASTM D1424	Lb (N)	0.67 (2.98)
			0.72 (3.20)
Elongation	ASTM D5035	%	21
			28
Bursting strength	ISO 2758	kPa	1,055

9. Supply Package and Compatible Devices Information

9.1 Supply Package

Each cardboard box of the Bioport Crosslink Intra-Articular Injection Prosthesis medical device contains:

- 1) Individual medical device packaging - 1 pc
- 2) Instructions for use - 1 pc
- 3) Patient's stickers - 2 pcs

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The transport container contains 200 (Two hundred) cardboard boxes.

9.2 Compatible devices

Prior to the injection of Bioport Crosslink into the joint, it is required to attach any sterile needle of gauge 20 (20G) to the syringe.

9.3 Accessories Information

Bioport Crosslink Intra-Articular Injection Prothesis is supplied without any accessories. The injection needle must be acquired separately by the user.

10. Description of Accessories, Medical Devices or Non-medical Devices Intended for Use in Combination with the Medical Device (if any)

See cl. 9.2, 9.3.

11. Information on the Presence of Medicinal Products for Human Use or Materials of Animal and (or) Human Origin in the Medical Device

Medical products information

The Bioport Crosslink Intra-Articular Injection Prothesis contains sodium hyaluronate

Manufacturer: Bloomage Freda Biopharm Co., Ltd. (China)

Registration in the Russian Federation: Federal Reference Book entry number FS-001305 dated 11.01.2016

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Information on the presence of materials of animal and (or) human origin in the medical device

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Not applicable. Bioport Crosslink does not contain any materials of animal and (or) human origin.

12. Information on the Procedure for Assembly, Installation, Adjustment, Calibration and Other Activities Required for Putting the Medical Device into Operation

Not applicable.

13. Requirements on Premises in which the Medical Device is to be Installed, as well as Requirements on Training or Qualification of People Involved in Installation of the Medical Device (if any)

Not applicable.

14. Information to Verify the Installation of the Medical Device and its Readiness for Safe Work/Operation

a) Scope and frequency of maintenance, including cleaning and disinfection of the medical device:

Not applicable. Bioport Crosslink is a single-use device.

b) list of information, keys, passwords and software provided by the manufacturer of the medical device and needed for installation, setup, operation and maintenance of the medical device:

Not applicable.

c) list of consumables (components, chemical agents), and procedure for their application and replacement:

Not applicable.

d) need for calibration in order to ensure the proper, safe operation of the medical device during its service life:

Not applicable.

e) methods for mitigation of risks related to installation, calibration or maintenance of the medical device:

Not applicable.

f) information on installation, setup, adjustment, calibration and other activities required for putting the medical device into operation and its proper use (application):

Not applicable.

g) information on the list of the main specifications related to the use (application) of the medical device, transportation and storage conditions (e.g. ambient air temperature and humidity, illumination and other characteristics):

see Section 6.

h) list of national standards applied by the manufacturer of the medical device:

Standards	Title
EN ISO 14630:2012 (GOST R ISO 14630-2017)	Non-active surgical implants - General requirements (ISO 14630:2012)
ASTM F2347-15	Standard Guide for Characterization and Testing of Hyaluronan as Starting Materials Intended for Use in Biomedical and Tissue Engineered Medical Product Applications
European Pharmacopoeia, Article 1472	Sodium hyaluronate
EN ISO 10993-1:2020 (GOST ISO 10993-1-2021)	Biological evaluation of medical devices - Part 1: Evaluation and Testing Within a Risk Management Process (ISO 10993-1:2018)
EN ISO 10993-3:2014	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity (ISO 10993-3:2014)
EN ISO 10993-5:2009 (GOST ISO 10993-5-2011)	Biological evaluation of medical devices - Part 5: Tests for <i>in vitro</i> cytotoxicity (ISO 10993-5:2009)
EN ISO 10993-6:2016 (GOST ISO 10993-6-2021)	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation (ISO 10993-6:2016)
ISO 10993-10:2021 (GOST ISO 10993-10-2011)	Biological evaluation of medical devices - Part 10: Tests for skin sensitization (ISO 10993-10:2021)
EN ISO 10993-11:2018 (GOST ISO 10993-11-2021)	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity (ISO 10993-11:2017)
EN ISO 10993-12:2012 (GOST ISO 10993-12-2015)	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials (ISO 10993-12:2021)
EN ISO 10993-16:2017	Biological evaluation of medical devices - Part 16: Toxicokinetic study design for degradation products and leachables (ISO 10993-16:2017)
EN ISO 10993-18:2020	Biological evaluation of medical devices - Part 18: Chemical characterization of medical device materials within a risk management process (ISO 10993-18:2020)
EN 556-2: 2015	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 2: Requirements for aseptically processed medical devices
EN ISO 11737-1:2018	Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on products. (ISO 11737-1:2018)
EN ISO 11607-1:2020 (GOST ISO 11607-1-2018)	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems. (ISO 11607-1:2019)
EN ISO 13408-1:2015	Aseptic processing of health care products - Part 1: General requirements (ISO 13408-1:2008/Amd 1:2013)

ISO 13408-7:2012	Aseptic processing of health care products - Part 7: Alternative processes for medical devices and combination products.
EN ISO 15223-1:2021 (GOST R ISO 15223-1-2020)	Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements (ISO 15223-1:2021)
EN ISO 20417:2021	Medical devices - Information to be supplied by the manufacturer (ISO 20417:2021)
EN ISO 14971: 2019	Medical devices - Application of risk management to medical devices (ISO 14971: 2019)
SanPiN 2.1.3684-21	"Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures"

15. Information on Medical Device Sterility, Method of Sterilization and Operational Procedure in Case of Affected Packaging Sterility

Bioport Crosslink Intra-Articular Injection Prothesis is a sterile single-use device. Sterilization method: sterilization involving the aseptic filling techniques. And the terminal sterilization method of Bioport Crosslink is Moist Heat sterilization.
Do not use if the package is damaged.

16. Information on the Procedure for Medical Device Treatment for its Reuse, Including Cleaning, Disinfection, Packaging and, if Required, the Method of Resterilization (if Reusable), as Well as Criteria for Unsuitability of the Medical Device for its Intended Use

Not applicable to this medical device, as Bioport Crosslink Intra-Articular Injection Prothesis is supplied sterile and is intended for single use only. Do not reuse. Do not resterilize. Do not recycle. Recycling may affect sterility, biocompatibility and functional integrity of the device.

17. Information Required for Identifying Medical Devices to Obtain a Safe Combination and Information on the Known Limitations of Combined Use of Medical Devices (for Medical Devices to be Used in Combination with Other Medical Devices and (or) Accessories)

Prior to the injection of Bioport Crosslink into the joint, it is required to attach any sterile needle of gauge 20 (20G) to the syringe.

The safety and performance of the use of Bioport Crosslink in combination with other intra-articular injections are not established.

18. Information on the Nature, Type, as well as (if Required) Intensity and Distribution of Radiation (Electromagnetic, Ionizing, Other) from the Medical Device and Methods of Protection of Consumers and Third Parties from Unintended Radiation during the Medical Device Operation (if the Medical Device Creates a Hazardous or Potentially Hazardous Radiation Level when Used as Intended)

Not applicable.

19. Information on Precautions to be Taken in Case of

a) non-serviceability of the medical device, its operation faults or functional deviations which may affect the safety of the medical device, including ones identified by external signs:

Do not use the device if package is damaged

b) impact on the function of the medical device due to external factors related to the use of the medical device in combination with other medical devices and (or) equipment or due to predictable factors such as external electromagnetic fields, electrostatic discharges, radiation (electromagnetic, ionizing, other), atmospheric pressure and its fluctuations, ambient air temperature and humidity:

For the operating conditions, see Section 6.

c) risk of electromagnetic interference caused by the medical device and affecting other medical devices, equipment and communication means during the diagnostics and assessment of its results or treatment in accordance with the intended use of the device (e.g. electromagnetic radiation produced by the medical device and affecting other equipment):

Not applicable.

20. Warnings and (or) Precautions to be Taken by the Consumer when Using a Medical Device Containing Medicinal Products for Human Use or Materials of Animal and (or) Human Origin or Carcinogenic, Mutagenic or Toxic Materials which, if Released or Leached, Cause Sensitization, Allergic Reaction or Reproductive Function Impairment

See Section 4.

21. Shelf-life

**The shelf life of Bioport Crosslink Intra-Articular Injection Prosthesis is 3 years.
The guaranteed storage life is equal to the shelf life.**

Service life: not applicable, as Bioport Crosslink is a single-use device.

22. Disposal

**Bioprot Crosslink Intra-Articular Injection Prosthesis must be used immediately after the package opening. The disposal must be performed by a safe and legally permitted method in compliance with the requirements of local regulations regarding medical waste management.
The disposal must be performed in compliance with the law in force in the Russian Federation. Waste class according to SanPIN 2.1.3684-21 «B».
Non-contaminated cardboard packaging must be disposed as domestic waste.**

23. Information on Circumstances Requiring the Consumer to Consult a Medical Officer

Not applicable. The device is intended to be used only by qualified medical specialists in the clinical environment.

24. Information on the First Issue and Latest Revision of the Operating Documentation

Version 3.

25. Procedure and Conditions for Disposal or Elimination of the Medical Device

Bioport Crosslink Intra-Articular Injection Prothesis must be used immediately after the package opening. The disposal must be performed by a safe and legally permitted method in compliance with the requirements of local regulations regarding medical waste management. The disposal must be performed in compliance with the law in force in the Russian Federation. Waste class according to SanPiN 2.1.3684-21 «5».

Non-contaminated cardboard packaging must be disposed as domestic waste.

26. Warranty

Maxigen Biotech Inc. guarantees that due attention was paid during the development and manufacturing of this device.

This warranty is provided in lieu of and excludes any other warranties not expressly set forth herein, whether express or implied warranties statutory or otherwise, including any implied warranties of merchantability or fitness for a particular purpose.

The processing, storage, cleaning and sterilization of this device, as well as factors related to the patient, diagnosis, treatment, surgical procedures, including those performed by specialists and doctors, as well as other matters beyond the control of Maxigen Biotech Inc., have a direct impact on the device itself and the results of its use.

The obligations of Maxigen Biotech Inc. under this warranty are confined to replacement of this device. Maxigen Biotech Inc. shall not be held liable for any accidental or indirect losses, damages or expenses directly or indirectly associated with the use of this device. Maxigen Biotech Inc. shall not take, as well as shall not empower any other parties to take on its behalf, any other or additional obligations or liability with respect to this device.

Maxigen Biotech Inc. shall not be held liable with respect to reused, recycle or sterilized devices and shall not provide any express or implied warranties with respect to such devices, including any warranties of merchantability or fitness for a particular purpose.

All inquiries related to the warranty in the Russian Federation shall be submitted to the authorized representative appointed by Maxigen Biotech Inc.

27. Claim Filing Procedure

Should any issues arise with respect to the quality of the Bioport Crosslink Intra-Articular Injection Prothesis medical device manufactured by Maxigen Biotech Inc., they must be communicated to the Authorized Representative of the manufacturer in the Russian Federation: Bioport Limited Liability Company (Bioport LLC), Office 408-1, Floor 4, 26V bldg 1, Bolshaya Pochtovaya Str., Basmanny Municipal Dist. (intracity territory of the federal city), 105082, Moscow
Phone: 8 (800) 222-38-02
E-mail: bioport@inbox.ru

047

No Responsibility For The Contents of This Document

Case No.	000608	Date	NOV. 15 2023
本文件由康生生物科技股份有限公司之董事、於臺灣臺北地方法院 附屬民間公證人以表轉合事務所認證。公證人 沈怡君			
Attested at the Xinyi Associated Notary Public Office of Taiwan Taipei District Court, R.O.C., that the seal(s) of MAXIGEN BIOTECH INC. in this document is/are authentic.			
Notary Public			
 Shen Yi-Chun			

ADD: 4F.-3, No.415, Sec. 4, Xinyi Rd., Taipei, Taiwan
TEL: 886-2-2729-3139 FAX: 886-2-2729-0509

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