SP x secret-41

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Addition to the Annex to IFU Microcatheters
Description of documents

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h h h h h h
Signature of document custodian (affiant)

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
State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 29th day of MAY,
2018 by Ashni De Alwis, proved

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Name of the medical device	Micro Catheter, in various versions Echelon, Rebar, Apollo, Marathon, Orion 21:
	I. Micro Catheter, in various versions Echelon, Rebar, Apollo, Marathon, Orion 21: 1. 1. Micro Catheter Echelon, versions: 1.1.1. Micro Catheter Echelon-10, tip shape 0°, consisting of: - Micro Catheter Echelon-10, tip shape 0° - Mandrel with pointed end 1.1.2. Micro Catheter Echelon-10, tip shape 45°, consisting of: - Microcatheter Echelon-10, tip shape 45° - Mandrel with pointed end - Split introducer 1.1.3. Micro Catheter Echelon-10, tip shape 90°, consisting of: - Micro Catheter Echelon-10, tip shape 90° - Mandrel with pointed end - Split introducer 1.1.4. Micro Catheter Echelon-14, tip shape 0°, consisting of: - Micro Catheter Echelon-14, tip shape 0° - Mandrel with pointed end 1.1.5. Micro Catheter Echelon-14, tip shape 45°, consisting of: - Micro Catheter Echelon-14, tip shape 45° - Mandrel with pointed end - Split introducer 1.1.6. Micro Catheter Echelon-14, tip shape 90°, consisting of: - Micro Catheter Echelon-14, tip shape 90° - Mandrel with pointed end - Split introducer 1.2. Accessories: - Adapter to the Micro Catheter Echelon 2. Micro Catheter Rebar, versions: 2.1. Micro Catheter Rebar-10, consisting of: - Micro Catheter Rebar-10 - Mandrel with pointed end 2.2. Micro Catheter Rebar-14, consisting of: - Micro Catheter Rebar -14 - Mandrel with pointed end 2.3. Micro Catheter Rebar -18, length 130 cm, consisting of: - Micro Catheter Rebar -18, length 130 cm

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
State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 31st day of Oct.,
20 17 by Ashni De Alwis —, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

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This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.

CONTENTS

One Rebar™ Micro Catheter with tip shaping mandrel

DESCRIPTION

The Rebar™ Micro Catheter is an endhole, single-lumen catheter designed to be introduced over a steerable guidewire into the vasculature. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a semi-rigid proximal shaft which transitions into the flexible distal shaft to facilitate the advancement of the catheter in the anatomy. Single or dual radiopaque markers at the distal end facilitate fluoroscopic visualization. The outer surface of the catheter is coated to increase lubricity.

INDICATIONS FOR USE

The Rebar™ Micro Catheter is intended for the controlled selective infusion of physician-specified therapeutic agents or contrast media into the vasculature of the peripheral and neuro anatomy.

CONTRAINDICATIONS

- It is not intended for use in the coronary vasculature.
- The Rebar™ Micro Catheter is contraindicated when, in the medical judgment of the physician, such a procedure may compromise the patient's condition.
- The Rebar™ Micro Catheter is contraindicated for neonatal and pediatric use.

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- Hematoma at the site of entry
- Vascular thrombosis
- Vessel perforation
- Thrombotic episodes
- Vessel spasm
- Neurological deficits including stroke and death
- Hemorrhage
- Ischemia
- Pain and tenderness

WARNINGS

- Infusion pressure with this device should not exceed pressure as specified in the package insert. Excess pressure may result in catheter rupture, which may result in patient injury.
- If flow through the micro catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the micro catheter to determine the cause of the obstruction or replace it with a new micro catheter. Excessive pressure may cause catheter rupture, which may result in patient injury.
- Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

PRECAUTIONS

- Prior to use, carefully examine the Rebar™ Micro Catheter and its packaging to verify that it has not been damaged during shipment.
- Emboloc agents must be used in accordance with manufacturer's instructions.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- Always monitor infusion rates when using this micro catheter.
- The micro catheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated in order to be lubricious. This can be accomplished by attaching a Y-connector to a continuous saline drip.
- When injecting contrast for angiography, ensure that the catheter is not kinked or occluded.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.
- Because the Rebar™ may be easily advanced into narrow, selective vasculature, repeatedly verify (for example by slightly withdrawing the catheter) that the Rebar™ has not been advanced so far as to interfere with its removal.

STORAGE

This device should be stored in a dry place at between 50°F (10°C) and 90°F (32°C).

DIRECTIONS FOR USE

NOTE: See package insert for Rebar™ functional characteristics.

- Place the appropriate guide catheter following recommended procedures. Connect a one-way stop-cock to the guide catheter to prevent backflow of blood during insertion of catheter. Connect a hemostatic side arm adapter one way stop-cock to allow for insertion of catheter and to continuously flush the guide catheter with saline.

NOTE: It is recommended that a guide catheter with a minimum internal diameter of 0.053" be used with the Rebar™ Micro Catheter.

- Gently remove the packaging coil for the Rebar™ Micro Catheter from the pouch.

- Flush the coil with heparinized saline through the female luer attached to the coil.
 - Flush the ID of the catheter with heparinized saline by attaching a saline filled syringe to the catheter hub.
 - Remove the appropriate steerable guidewire from its packaging and inspect for damage. Follow manufacturer's instructions for preparing and using the guidewire.
 - Carefully insert guidewire into the hub of the Rebar™ Micro Catheter and advance into the catheter lumen.
 - Remove the micro catheter by removing the hub from the clip and gently pulling on the hub.
 - Inspect the catheter prior to use to verify that it is undamaged. Retain the packaging coil for storage of the catheter when it is not in use during the procedure.
 - Close the one-way stop-cock.
 - Loosen the hemostatic valve.
 - Introduce the guidewire and micro catheter as a unit through the hemostatic sidearm adapter into the lumen of the guiding catheter. Advance guidewire/catheter assembly to the distal tip of the guiding catheter.
 - Tighten the valve around the micro catheter to prevent backflow, but allowing some movement through the valve by the micro catheter.
 - Open the one-way stop-cock.
 - Alternatively advance the guidewire and micro catheter until the desired site has been reached.
- NOTE:** To facilitate catheter handling the proximal portion of the micro catheter is uncoated to ensure a non-slip grip.
- When ready to infuse, withdraw the guidewire completely from the Rebar™ Micro Catheter. Connect a syringe containing infusate to the micro catheter hub, and infuse as required.
 - After completing the procedure, withdraw the micro catheter and discard.

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- Remove shaping mandrel from card and insert into distal tip of the catheter.
- Carefully bend catheter tip and shaping mandrel to desired shape. A slight over exaggeration of the shape may be required to accommodate for catheter relaxation.
- Shape the catheter by holding the shaped portion approximately 1 inch (2.5 cm not less than 1 cm) from the steam source for about 20 seconds (DO NOT EXCEED 30 SECONDS.)
- Allow catheter tip to cool in air or saline prior to removing mandrel.
- Remove mandrel from catheter and discard. Multiple shaping is not recommended.
- Inspect the tip for any damage that may have resulted from steam shaping the catheter tip. If any damage is found, do not use the catheter.

FUNCTIONAL CHARACTERISTICS

Flow Rates @ 100 PSI (690 kPa) and 300 PSI (2,070 kPa) – ml/sec

	Saline		76% Contrast/ Saline (50/50)		60% Ionic Contrast		76% Ionic Contrast	
Pressure	100	300	100	300	100	300	100	300
Rebar™-10								
105-5078-153	0.26	0.68	0.10	0.38	0.04	0.16	0.03	0.09
105-5078-170	0.26	0.59	0.10	0.38	0.04	0.15	0.02	0.08
Rebar™-14								
105-5080-143	0.4	0.8	0.19	0.6	0.07	0.30	0.03	0.11
105-5080-153	0.3	0.6	0.12	0.4	0.06	0.18	0.03	0.09
105-5080-170	0.2	0.6	0.08	0.3	0.05	0.15	0.01	0.07
Rebar™-18								
105-5083-153	0.9	1.5	0.4	1.0	0.2	0.6	0.1	0.3
105-5081-153	0.9	1.5	0.4	1.0	0.2	0.6	0.1	0.3
105-5081-130	0.9	1.7	0.5	1.1	0.2	0.7	0.1	0.4
105-5081-110	1.0	1.8	0.6	1.3	0.3	0.7	0.1	0.4
Rebar™-27								
105-5082-145	1.4	2.5	0.6	1.8	0.3	1.3	0.2	0.7
105-5082-130	1.8	3.0	1.1	2.6	0.7	1.7	0.3	0.9
105-5082-110	2.0	3.4	1.3	2.7	0.9	1.9	0.3	0.9

All Rebar™ micro catheters have a maximum allowable dynamic pressure of 700 PSI (4,830 kPa) and static pressure of 300 PSI (2,070 kPa).

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Ashni De Alwis
Signature of document custodian (affiant)

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
State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 31st day of Oct.,
2017 by Ashni De Alwis, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Steve Oberlander
Signature (Seal)

INSTRUCTIONS FOR USE

Apollo™ Onyx Delivery Micro Catheter

[Handwritten signature]

70728-001, Rev. 05/13

TABLE OF CONTENTS

Apollo™ Onyx Delivery Micro Catheter

English	3
Français	4
Deutsch	5
Italiano	6
Español	8
Svenska	9
Nederlands	10
Português	12
Suomi	13
Dansk	14
Ελληνικά	15
Čeština	17
Magyar	18
Русский	19
Polski	21
Türkçe	22
Norsk	23
Slovenčina	24
Română	25
Български	27
한국어	28
عربي	30
Symbol Glossary	31

English Instructions for Use

EN

Apollo™ Onyx Delivery Micro Catheter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.

This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures.

CONTENTS

One Apollo Onyx Delivery Micro Catheter with Syringe Adapter (50414).


DESCRIPTION

The Apollo Onyx Delivery Micro Catheter is a single-lumen, endhole catheter designed for the superselective infusion of physician-specified therapeutic agents such as embolization materials and diagnostic materials such as contrast media in tortuous, distal vessels. The catheter has a semi-rigid proximal shaft and a highly flexible distal shaft to facilitate the advancement of the catheter in the anatomy. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The outer surfaces of the catheter are coated to increase lubricity.

The Apollo catheter is designed to facilitate catheter retrieval in the event the catheter becomes entrapped within the vasculature. The distal section of the catheter incorporates a detachment zone that allows detachment of the distal tip when the force required to extract the catheter exceeds the force to detach the tip. The catheter has 2 radiopaque marker bands to visualize the position of the catheter and the detach zone area:

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The Apollo is packaged with a Syringe Adapter (50414). This device, attached to a 1 mL syringe filled with Onyx®, will reduce the dead space within the micro catheter luer hub. Reducing the dead space within the hub is intended to minimize the potential mixing of Onyx® and DMSO in the hub of the catheter during connection and injection.

INDICATIONS FOR USE

The Apollo Onyx Delivery Micro Catheter is intended to access the neuro vasculature for the controlled selective infusion of physician-specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media.

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- Puncture site hematoma
- Pain and tenderness
- Vessel perforation
- Thrombotic episodes
- Vessel spasm
- Neurological deficits including stroke and death
- Hemorrhage
- Vascular thrombosis

WARNINGS

- Not intended for use with embolization particles, detachable coils or Onyx HD500.
- Do not steam shape the tip of the microcatheter. Steam shaping the catheter tip may cause damage to the detachment zone and unintended detachment.

WARNINGS

- Do not use a cannula or needle to introduce the end of the catheter into a connector. Insertion of a cannula or needle may cause damage to the detachment zone and result in unintended detachment.
- Always handle the distal end of the catheter with care to avoid damage to the detachment zone and unintended detachment.
- Infusion pressure with this device should not exceed 690 kPa/100 psi. Pressure in excess of 690 kPa/100 psi may result in catheter rupture, possibly resulting in patient injury.
- If flow through the catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Remove the catheter and replace it with a new catheter. Excessive pressure may cause catheter rupture, possibly resulting in patient injury.
- Never advance or withdraw an intraluminal device against resistance. Excessive force against resistance may result in damage to the device or vessel perforation.
- When injecting contrast for angiography, ensure that the catheter is not kinked, prolapsed or occluded. Remove excess slack in the catheter to reduce the potential for catheter kink or prolapse.
- Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damage or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.
- Regardless of the liquid embolic being used, leave a gap between the reflux and the proximal marker band. Excessive reflux may result in difficult catheter removal.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

PRECAUTIONS

- Select tip size based on angiography. The detachment zone should never be distal to the last tortuous curve of the vessel. Refluxing over the detachment zone distal to the last tortuous curve may result in catheter entrapment. Do not place catheter such that the detached tip could interfere with patent vessels.
- Prior to use, carefully examine the Apollo and its packaging to verify that it has not been damaged during shipment. Do not touch or manipulate the catheter tip prior to use.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- During navigation, check that the distal tip of the catheter is not kinked before passing the guidewire through it. Kinking or prolapsing of the catheter may result in unintended rupture of the catheter.
- Always monitor infusion rates when using the catheter.
- The Apollo has a hydrophilic coating on the outside of the catheter which must be kept hydrated.
- This catheter is not intended for use with chemotherapy agents.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.
- Navigating or repositioning the catheter while it is in a wedged position or with vessels that are in vasospasm may cause premature tip detachment.
- Do not reposition catheter after start of Onyx injections.
- When performing angiography, it is recommended to use a 3cc syringe rather than a 1cc syringe to reduce the risk of catheter over-pressurization.
- The Apollo is a flow directed micro catheter that can optionally be used with hydrophilic, 0.010" or less sized guidewires. The Apollo is not compatible with non-hydrophilically coated guidewires or guidewires greater than 0.010" in diameter. It is recommended that the Apollo be used with an appropriately sized guiding catheter which allows adequate clearance (minimum internal diameter of 0.053" or 1.35mm).
- When withdrawing the catheter, monitor the distal tip under angiography. Pulling the catheter against significant resistance can cause patient injury. If significant catheter resistance is felt, refer to the precaution in the procedure section below for guidance.
- If catheter entrapment is suspected (with any embolic agent), fast catheter retrieval technique may result in catheter shaft separation and potential vascular damage. Follow catheter retrieval instructions at the end of instructions for use.

STORAGE

Store the Apollo Onyx Delivery Micro Catheter in a dry place between 50°F (10°C) and 90°F (32°C).

DIRECTIONS FOR USE

Procedure

1. Place the appropriate guiding catheter following recommended procedures. Connect a one-way stop-cock to the guiding catheter to prevent backflow of blood during insertion of catheter. Connect a hemostatic side arm adapter and a one-way stop-cock to allow for insertion of catheter and to facilitate continuous flush of the guiding catheter with saline.
NOTE: It is recommended that a guiding catheter with a minimum internal diameter of 0.053" (1.35 mm) be used with the Apollo Onyx Delivery Micro Catheter.
2. Gently remove the packaging coil for the Apollo Onyx Delivery Micro Catheter from the package. Flush the coil with heparinized saline through the female luer attached to the coil.

3. Flush the coil with saline through the female luer attached to the coil.
4. Remove the catheter by removing the hub from the clip and gently pulling on the hub.
5. Prior to use, flush the catheter lumen with heparinized saline by attaching a saline filled syringe to the catheter hub.
6. Remove the appropriate steerable guidewire from its packaging (the Silverspeed .010 and the Mirage .008 have been qualified for use with the Apollo) and check for damage. Follow the manufacturer's instructions for preparing and using the guidewire.
7. Insert the catheter, taking care to not touch or manipulate the tip prior to use to verify that it is undamaged. Retain the packaging coil for storage of the catheter when it is not in use during the procedure.
8. Carefully insert the guidewire into the hub of the micro catheter and advance the guidewire into the catheter lumen.
9. A split introducer is loaded on the proximal end of the catheter to aid in insertion into the hemostatic side-arm adapter. To use, slide the split introducer from the proximal end until it covers the distal end of the catheter.
10. Close the one-way stop-cock.
11. Loosen the hemostatic valve.
12. Introduce the guidewire, micro catheter as a unit through the hemostatic valve. Advance the guidewire/catheter assembly to the distal tip of the guiding catheter. Slide the introducer back to the proximal end of the catheter next to the hub.
13. Tighten the valve around the catheter to prevent backflow, but allowing some movement through the valve by the catheter.
14. Open the one-way stop-cock.
15. The catheter may be advanced through the guidecatheter and the vasculature by gently pushing the proximal shaft. It is recommended to use a guidewire during navigation to reduce the risk of catheter kink or prolapse.

WARNING

- Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damage or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.
- 16. To infuse, connect a syringe with infusate to the catheter luer, and infuse as required.

SYRINGE ADAPTER (50414) DIRECTIONS FOR USE:

Usable Length	Apollo Minimum Dead Space Volume	Catheter Dead Space with Syringe Adapter (50414)	Approximate Infusion Rate at 100 psi (690 kPa)	
			Water	Contrast (76% Renografin)
165 cm	>.23 ml	>.20ml	12 ml/min	3 ml/min

1. Prepare the Onyx according to the Onyx Instructions for Use.
2. Following aspiration of Onyx into syringe, detach the needle and attach the Syringe Adapter to the syringe. Purge air from the Syringe Adapter.
3. Immediately connect the Onyx syringe with the Syringe Adapter to the hub of the previously DMSO primed micro catheter according to the Onyx Instructions for Use. Make sure the connection is tight and there is no air in the hub following connection.
4. Proceed with the Onyx injection according to the Onyx Instructions for Use.
5. If a second Onyx syringe is required, **do not remove the Syringe Adapter from the catheter**. When the new Onyx syringe is ready, simply remove the empty syringe from the interface device's proximal end and connect the next syringe to the Syringe Adapter being careful to avoid introducing air.

CATHETER RETRIEVAL INSTRUCTIONS:

PRECAUTION: If catheter entrapment is suspected (with any embolic agent), fast catheter retrieval technique may result in catheter shaft separation and potential vascular damage.

1. Slowly remove any 'slack' in the distal catheter shaft.
2. Gently and slowly apply 3-5 cm traction to the catheter to begin catheter retrieval.
3. Detachment of the catheter may be observed by visualizing the separation of the catheter between the distal and proximal marker bands.
4. In the event that the catheter does not detach, assess the following parameters by observing the distal shaft of the catheter:

- Vessel straightening
- Catheter tip releasing from embolic cast
- Traction on embolic cast

NOTE: Do not apply more than 20 cm of traction to catheter to minimize risk of catheter separation, proximal to the detachment zone.

5. Under some difficult clinical situations, it may be safer to leave a flow-directed catheter in the vascular system rather than risk rupturing the malformation and, consequently a hemorrhage, by exercising too much traction on an entrapped catheter.

This is accomplished by stretching the catheter and cutting the shaft near the entry point of vascular access allowing the catheter to remain in the artery.

If catheter breaks during removal, distal migration or coiling of the catheter may occur. Same day surgical resection should be considered to minimize risk of thrombosis.

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Echelon™ Micro Catheter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician. This device should be used only by physicians with a thorough understanding of angiography and/or percutaneous neurointerventional procedures.

CONTENTS

One Echelon™ Micro Catheter with tip shaping mandrel.

DESCRIPTION

The Echelon™ Micro Catheter is an endhole, single-lumen catheter designed to be introduced over a steerable guidewire into the vasculature. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a semi-rigid proximal shaft which transitions into the flexible distal shaft to facilitate the advancement of the catheter in the anatomy. Dual radiopaque markers at the distal end facilitate fluoroscopic visualization. The outer surface of the catheter is coated to increase lubricity.

INDICATIONS FOR USE

The Echelon™ Micro Catheter is intended to access peripheral and neuro vasculature for the controlled selective infusion of physician-specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media.

CONTRAINDICATIONS

- Not intended for use in the coronary vasculature.
- Contraindicated when, in the medical judgment of the physician, such a procedure may compromise the patient's condition.
- Contraindicated for neonatal and pediatric use.

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- Hematoma at the site of entry
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Vascular thrombosis
- Thrombotic episodes
- Neurological deficits including stroke and death
- Ischemia

WARNINGS

- Infusion pressure with this device should not exceed 600 PSI (4137 kPa). Pressures in excess of 600 PSI (4137 kPa) may result in catheter rupture, which may result in patient injury.
- If flow through the micro catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the micro catheter to determine the cause of the obstruction or replace it with a new micro catheter. Excessive pressure may cause catheter rupture, which may result in patient injury.
- Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation or other patient injuries.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

PRECAUTIONS

- Ensure embolic material compatibility with catheter prior to use.
- Prior to use, carefully examine the micro catheter and packaging to verify that they have not been damaged during shipment.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- Always monitor infusion rates when using the micro catheter.
- When injecting contrast for angiography, ensure that the micro catheter is not kinked or occluded.
- The micro catheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated in order to be lubricious. This can be accomplished by attaching the Y-connector to a continuous saline drip.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the micro catheter without observing the resultant tip response.
- Because the micro catheter may be easily advanced into narrow, selective vasculature, repeatedly verify (for example by slightly withdrawing the catheter) that the micro catheter has not been advanced so far as to interfere with its removal.

STORAGE

The micro catheter should be stored in a dry place at between 50°F (10°C) and 90°F (32°C).

DIRECTIONS FOR USE

Flow Rates @ 100 PSI (690 kPa) and 300 PSI (2,070 kPa) – mL/sec

	Saline		76% Contrast/ Saline (50/50)		60% Ionic Contrast		76% Ionic Contrast	
Pressure	100	300	100	300	100	300	100	300
Echelon 14 (105-5092-150, 145-5092-150, 190-5092-150)								
Flowrates	0.25	0.85	0.33	0.46	0.16	0.27	0.06	0.09
PVA Particles	Maximum 500µm							
Guide Catheter Minimum ID	0.050" (1.28 mm)							
Echelon 10 (105-5091-150, 145-5091-150, 190-5091-150)								
Flowrates	0.25	0.85	0.33	0.46	0.16	0.27	0.06	0.09
PVA Particles	Maximum 500µm							
Guide Catheter Minimum ID	0.035" (0.89mm)							

1. Place the appropriate guide catheter following recommended procedures. Connect a hemostatic side arm adapter to the guide catheter to prevent backflow of blood during insertion of micro catheter. Connect a one-way stopcock to the hemostatic side arm adapter to continuously flush the guide catheter with saline.
2. Model numbers 145-5091-150, 145-5092-150, 190-5091-150, and 190-5092-150 are provided in a tray with lid, and the catheter tip is held within the groove of a curve retainer. To open the tray, lift the lid at the corner of the tray with the catheter hub. The catheter and curve retainer may be removed from the tray together.
3. To free the catheter tip from the curve retainer, pinch two opposite corners back. This opens the groove. Lift the catheter tip out of the groove.
4. Hydrate the catheter outer surface.
 - a. If the catheter is provided in a coil dispenser, flush the coil with heparinized saline through the female luer attached to the coil.
 - b. Otherwise, the catheter may be hydrated in the tray.
5. Prior to use, flush the catheter lumen with heparinized saline by attaching a saline filled syringe to the catheter hub.
6. Remove the appropriate steerable guidewire from its packaging and inspect for damage. Follow manufacturer's instructions for preparing and using the guidewire.
7. Inspect the catheter prior to use to verify that it is undamaged. Retain the packaging for storage of the catheter if it is not used immediately in the procedure.
8. Carefully insert guidewire into the hub of the micro catheter and advance guidewire into the catheter lumen.
9. Model numbers 145-5091-150, 145-5092-150, 190-5091-150, and 190-5092-150 include a split introducer to aid in insertion into the hemostatic side arm adapter. To use, slide the Split Introducer from the proximal end until it covers the distal end of catheter.
10. Close the one-way stopcock.
11. Loosen the hemostatic valve.
12. Introduce the guidewire and micro catheter as a unit through the hemostatic port of the hemostatic sidearm adapter into the lumen of the guiding catheter. Advance guidewire/catheter assembly to the distal tip of the guiding catheter.
 - a. If using the Split Introducer slide the introducer back to the proximal end of the catheter next to the hub.
 - b. Or peel introducer away from the microcatheter.
13. Tighten the valve around the micro catheter to prevent backflow, but allowing some movement through the valve by the micro catheter.
14. Open the one-way stopcock.
15. Alternatively advance the guidewire and micro catheter until the desired site has been reached.
16. When ready to infuse, withdraw the guidewire completely from the micro catheter. Connect a syringe containing infusate to the hub of the micro catheter, and infuse as required.
17. After completing the procedure, withdraw the micro catheter and discard.

NOTE: To facilitate catheter handling, the proximal portion of the micro catheter is uncoated to ensure a non-slip grip.

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Remove shaping mandrel from card and insert into distal tip of the catheter.

18. Carefully bend catheter tip and shaping mandrel to desired shape. A slight over exaggeration of the shape may be required to accommodate for catheter relaxation.
19. Shape the catheter by holding the shaped portion approximately 1 inch (2.5 cm not less than 1 cm) from the steam source for about 20 seconds (DO NOT EXCEED 30 SECONDS.)
20. Allow catheter tip to cool in air or saline prior to removing mandrel.
21. Remove mandrel from catheter and discard. Multiple shaping is not recommended.
22. Inspect the tip for any damage that may have resulted from steam shaping the catheter tip. If any damage is found, do not use the catheter.

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[Signature]
Signature of document custodian (affiant)

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
State of California

County of Orange

Subscribed and sworn to (or affirmed) before me on this 31st day of Oct.
20 17 by Ashni De Alwis, proved
 to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

[Signature]
Signature

INSTRUCTIONS FOR USE

Marathon™ Flow Directed Micro Catheter

Handwritten signature

70879-001 Rev. 04/14

This page is intentionally left blank.

TABLE OF CONTENTS

Marathon™ Flow Directed Micro Catheter

English	5
Français	6
Deutsch	8
Italiano	9
Español	11
Svenska	12
Nederlands	14
Português	15
Suomi	17
Dansk	18
Ελληνικά	20
Česky	21
Magyar	23
Русский	24
Polski	26
Türkçe	27
Norsk	29
Slovenčina	30
Română	32
Български	33
한국어	35
عربي	38
Symbol Glossary	39

This page is intentionally left blank.

English Instructions for Use

EN

Marathon™ Flow Directed Micro Catheter

CAUTION

Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician. This device should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures.

CONTENTS

- One Marathon™ Micro Catheter
- Shaping Mandrel

ACCESSORIES: (May be included with device model.)

- 0.007" Stylet or smaller
- 0.010" Micro Guidewire or smaller
- 1.5F Introducer Sheath or larger

DESCRIPTION

The Marathon™ Flow Directed Micro Catheter is a single-lumen, endhole catheter designed for the subselective infusion of physician-specified therapeutic agents such as embolization materials and diagnostic materials such as contrast media in tortuous, distal vessels. The catheter has a semi-rigid proximal shaft and a highly flexible distal shaft to facilitate the advancement of the catheter in the anatomy. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a radiopaque marker at the distal end to facilitate fluoroscopic visualization. The outer surfaces of the catheter are coated to increase lubricity. Microcatheter may be used with stylet, guidewire or introducer sheath to increase the rigidity of the distal section during introduction into the guiding catheter.

INDICATIONS FOR USE

The Marathon™ Flow Directed Micro Catheter is intended to access the peripheral and neuro vasculature for the controlled selective infusion of physician-specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media.

CONTRAINDICATIONS

- The Marathon™ Flow Directed Micro Catheter is contraindicated when, in the medical judgment of the physician, use of such product may compromise the patient's condition.
- Not intended for use in the coronary vasculature.
- The Marathon™ Flow Directed Micro Catheter is contraindicated for neonatal and pediatric use.

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- | | |
|--------------------------|--|
| • Puncture site hematoma | • Pain and tenderness |
| • Vessel perforation | • Thrombotic episodes |
| • Vessel spasm | • Neurological deficits including stroke and death |
| • Hemorrhage | • Vascular thrombosis |

WARNINGS

- Infusion pressure with this device should not exceed 690 kPa/100 psi. Pressure in excess of 690 kPa/100 psi may result in catheter rupture, possibly resulting in patient injury.
- If flow through the catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the catheter to determine the cause of the obstruction or replace it with a new catheter. Excessive pressure may cause catheter rupture, possibly resulting in patient injury.
- Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation.
- When injecting contrast for angiography, ensure that the catheter is not kinked, prolapsed or occluded. Remove excess slack in the catheter to reduce the potential for catheter kink or prolapse.
- Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damage or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.
- Do not "over-shape" the catheter to achieve the desired shape angle. Over-shaping can lead to catheter kink or prolapse.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

- The stylet (supplied with some models) is not to be used as a guidewire. The stylet should never be manipulated within the catheter. **Do not advance the stylet beyond the tip of the catheter.** Use of the stylet as a guidewire could cause damage to the catheter and/or injure the patient.

PRECAUTIONS

- Prior to use, carefully examine the Marathon™ and its packaging to verify that they have not been damaged during shipment.
- Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
- Always monitor infusion rates when using this catheter.
- The Marathon™ has a hydrophilic coating on the outside of the catheter which must be kept hydrated.
- This catheter is not intended for use with chemotherapy agents.
- When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.
- Because the Marathon™ may be easily advanced into narrow, selective vasculature, repeatedly verify that the Marathon™ has not been advanced so far as to interfere with its removal by slightly withdrawing the catheter before each infusion.
- When performing angiography, it is recommended to use a 3cc syringe rather than a 1cc syringe to reduce the risk of catheter over-pressurization.
- The Marathon™ is a flow directed micro catheter that can optionally be used with hydrophilic, 0.010" or less sized guidewires. The Marathon™ is not compatible with non-hydrophilically coated guidewires greater than 0.010" in diameter.
- It is recommended that the Marathon™ be used with an appropriately sized guiding catheter which allows adequate clearance (minimum internal diameter of 0.053" or 1.35mm).

STORAGE

Store the Marathon™ Flow Directed Micro Catheter in a dry place between 50°F (10°C) and 90°F (32°C).

DIRECTIONS FOR USE

Marathon™ Flow Directed Micro Catheter

1. Place the appropriate guiding catheter following recommended procedures. Connect a one-way stop-cock to the guiding catheter to prevent backflow of blood during insertion of catheter. Connect a hemostatic side arm adapter one-way stop-cock to allow for insertion of catheter and to continuously flush the guiding catheter with saline.
NOTE: It is recommended that a guiding catheter with minimum internal diameter of 0.053" (1.35mm) be used with the Marathon™ Flow Directed Micro Catheter.
2. Gently remove the packaging coil for the Marathon™ Flow Directed Micro Catheter from the package.
3. Flush the coil with saline through the female luer attached to the coil.
4. Remove the catheter by removing the hub from the clip and gently pulling on the hub.
5. Inspect the catheter prior to use to verify that it is undamaged. Retain the packaging coil for storage of the catheter when it is not in use during the procedure.
6. Flush the ID of the micro catheter coil with saline.
7. To assist with insertion of the microcatheter in to the hub, the following accessories may be used: a guidewire, stylet, and/or introducer sheath.
8. If using a stylet and/or guidewire:
 - Prepare stylet or guidewire per manufacture's instructions; product should be hydrated.
 - Remove the stylet or guidewire from its packaging coil
 - Insert the stylet or guidewire through the hub of the Micro Catheter.
 - Advance the stylet or guidewire through the catheter, keeping the catheter body straight. **Do not advance the stylet past the tip of the catheter.**
 - Close the one-way stop-cock.
 - Loosen the hemostatic valve.
 - Introduce the catheter through the hemostatic adapter.
 - Tighten the valve around the catheter to prevent backflow, but allowing some movement through the valve by the catheter.
 - Open the one-way stop-cock.
 - Advance the catheter such that the stylet or guidewire remains in the guiding catheter. **Do not advance the stylet past the tip of the guiding catheter.**
 - Once the catheter has been advanced to the end of the guiding catheter, if using a stylet withdraw the stylet from the catheter.
9. The catheter may be advanced through the vasculature by gently pushing the proximal shaft. It is recommended to use a guidewire during navigation to reduce the risk of catheter kink or prolapse.
WARNING: Verify catheter integrity prior to re-inserting guidewire or injecting embolic material to prevent vascular damaged or unintended embolization. Catheter integrity is verified by angiographically confirming that contrast agent is exiting only from the catheter tip while viewing the entire distal section of the catheter.

10. To infuse, connect a syringe with infusate to the catheter luer, and infuse as required.

Functional Characteristics

Usable Length	Minimum Dead Space Volume	Maximum Embolization Particle Size	Approximate Infusion Rate at 100 psi (690 kPa)	
			Water	Contrast (76% Renografin)
165 cm	0.23 ml	250 µ	12 ml/min	3 ml/min

Note: Embolization particle concentration and injection rate should be verified by the physician prior to use. High concentrations of particles have a tendency to aggregate in the catheter hub.

STEAM SHAPING MANDREL

In order to maintain catheter integrity and dimensional stability of the inner diameter, it is recommended that the user follow these instructions.

WARNINGS

Shaping Mandrel is not intended for use in the human body.

Use only a steam source to shape the catheter tip. Do not use other heat sources.

Prior to use, inspect the catheter tip for any damage that may have resulted from shaping. Do not use a catheter that has been damaged in any way. Damaged catheters may rupture causing vessel trauma or tip detachment during steering maneuvers.

1. Remove shaping mandrel from card and insert into distal tip of the catheter.
2. Carefully bend catheter tip and shaping mandrel to desired shape. A slight over exaggeration of the shape may be required to accommodate for catheter relaxation.
3. Shape the catheter by holding the shaped portion approximately 1 inch (2.5 cm not less than 1 cm) from the steam source for about 20 seconds (DO NOT EXCEED 30 SECONDS).
4. Allow catheter tip to cool in air or saline prior to removing mandrel.
5. Remove mandrel from catheter and discard. Multiple shaping is not recommended.
6. Inspect the tip for any damage that may have resulted from steam shaping the catheter. If any damage is found, do not use the catheter.

PRECAUTIONS

If catheter entrapment is suspected (with any embolic agent), fast catheter retrieval technique may result in catheter shaft separation and potential vascular damage.

1. Slowly remove any 'slack' in the distal catheter shaft.
 2. Gently and slowly apply 3-5 cm traction to the catheter to begin catheter retrieval.
 3. Should catheter removal become difficult, the following will assist in catheter retrieval. Assess the following parameters by observing the distal shaft of the catheter:
 - Vessel straightening
 - Traction on embolic cast
 - Catheter tip releasing from embolic cast
 4. Further traction (of 3-5 cm) may be applied gently if necessary to transfer distal catheter shaft.
 - Hold this traction for a few seconds and release
 - Assess traction of vasculature to minimize risk of hemorrhage
 5. This process can be repeated intermittently until catheter is retrieved.
- NOTE:** Do not apply more than 20 cm of traction to catheter to minimize risk of catheter separation.
6. Under some difficult clinical situations, it may be safer to leave a flow-directed catheter in the vascular system rather than risk rupturing the malformation and, consequently a hemorrhage, by exercising too much traction on an entrapped catheter.

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State of California


County of Orange

Subscribed and sworn to (or affirmed) before me on this 31st day of Oct.
20 17 by Ashni De Alwis, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

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DESCRIPTION

The Orion Micro Catheter is a variable stiffness, single lumen micro catheter designed to be introduced over a steerable guidewire for accessing small, tortuous vascular areas. The outer surface of the micro catheter's distal segment is coated with a hydrophilic material to provide lubricity during use. The proximal end of the micro catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The micro catheter also incorporates a PTFE liner to facilitate movement or introduction of devices passed through its lumen. The Orion Micro Catheter has dual radiopaque markers at the distal tip to facilitate fluoroscopic visualization. The distal tip of the micro catheter is shapeable. The Orion micro catheter has a Fluoro-Safe™ marker at 88 cm from the distal tip of the microcatheter.

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CONTRAINDICATIONS

Not intended for delivering liquid embolics (e.g. Onyx™) into the neurovasculature. If the micro catheter tip becomes entrapped within an embolic mass it may be difficult to remove the micro catheter from the patient. Entrapped devices should not be implanted and may need to be surgically resected.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to, the following:

- Distal emboli (air, foreign material, tissue, or thrombus)
- Hemorrhage
- Ischemia
- Neurological deficit that could lead to stroke and death
- Hematoma
- Infection
- Pain at insertion site
- Vessel spasm, dissection, perforation or injury

PRECAUTIONS

1. This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures.
2. Inspect the micro catheter prior to use for any kinks or bends. Any micro catheter damage may decrease the desired performance characteristics.
3. Always monitor infusion rates when using the micro catheter.
4. The micro catheter has a lubricious hydrophilic coating on the outside of the micro catheter. It must be kept hydrated in order to be lubricious.
5. Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions.
6. When the micro catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the micro catheter without observing the resultant tip response.

WARNINGS

- THIS DEVICE IS INTENDED FOR SINGLE USE ONLY. Discard after one procedure. Structural integrity and/or function may be impaired after reuse or cleaning.
- Never advance or withdraw an intraluminal device against resistance until the cause of resistance is determined by fluoroscopy. If the cause cannot be determined, withdraw the catheter. Movement of the micro catheter against resistance may result in damage to the micro catheter, or the vessel.
- Infusion pressure with this device should not exceed 300 psi (2,068 kPa). Pressures in excess of 300 psi (2,068 kPa) may result in micro catheter rupture, which may result in patient injury.

WARNINGS

- If flow through the micro catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the micro catheter to determine the cause of the obstruction or replace it with a new micro catheter. Excessive pressure may cause micro catheter rupture, which may result in patient injury.
- Excessive force or tortuosity may cause the proximal segment of micro catheter to kink permanently.
- This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization increase the risks of patient infection and compromised device performance.

STORAGE

- Do not use if package is opened or damaged.
- Use prior to the "Use By" date.
- This device should be stored in a dry place at between 50°F (10°C) and 90°F (32°C).

DIRECTIONS FOR USE

1. Place the appropriate guide catheter following recommended procedures. Connect a one-way stop-cock to the guide catheter to prevent backflow of blood during insertion of micro catheter. Connect a hemostatic side arm adapter to the one way stop-cock to allow for insertion of micro catheter and to continuously flush the guide catheter with saline. Flush the lumen of the catheter with heparinized saline.
NOTE: It is recommended that a guide catheter with a minimum internal diameter of 0.044" be used with the Orion Micro Catheter.
2. Gently remove the packaging coil for the Orion Micro Catheter from the pouch.
3. Flush the coil with heparinized saline through the female luer attached to the coil.
4. Flush the ID of the micro catheter with heparinized saline by attaching a saline filled syringe to the catheter hub.
5. Remove the appropriate steerable guidewire from its packaging and inspect for damage. Follow manufacturer's instructions for preparing and using the guidewire.
6. Carefully insert guidewire into the hub of the Orion Micro Catheter and advance into the micro catheter lumen.
7. Remove the micro catheter: by removing the hub from the clip and gently pulling on the hub.
8. Inspect the micro catheter prior to use to verify that it is undamaged. Retain the packaging coil for storage of the micro catheter when it is not in use during the procedure.
9. Close the one-way stop-cock.
10. Loosen the hemostatic valve.
11. Introduce the guidewire and micro catheter as a unit through the hemostatic sidearm adapter into the lumen of the guiding catheter. Advance guidewire/micro catheter assembly to the distal tip of the guiding catheter.
12. Tighten the valve around the micro catheter to prevent backflow, but allowing some movement through the valve by the micro catheter.
13. Open the one-way stop-cock.
14. Alternatively advance the guidewire and micro catheter until the desired site has been reached.
15. When ready to infuse, withdraw the guidewire completely from the Orion Micro Catheter. Connect a syringe containing infusate to the catheter hub, and infuse as required.
16. After completing the procedure, withdraw the catheter and discard.

STEAM SHAPING MANDREL

In order to maintain micro catheter integrity and dimensional stability of the inner diameter, it is recommended that the user follow these instructions:

WARNINGS

- Shaping Mandrel is not intended for use in the human body.
 - Use only a steam source to shape the micro catheter tip. Do not use other heat sources.
 - Prior to use, inspect the micro catheter tip for any damage that may have resulted from shaping. Do not use a micro catheter that has been damaged in any way. Damaged micro catheters may rupture causing vessel trauma or tip detachment during steering maneuvers.
1. Remove shaping mandrel from card and insert into distal tip of the catheter.
 2. Carefully bend micro catheter tip and shaping mandrel to desired shape. A slight over exaggeration of the shape may be required to accommodate for catheter relaxation.
 3. Shape the micro catheter by holding the shaped portion approximately 1 inch (2.5 cm) not less than 1 cm) from the steam source for about 20 seconds (DO NOT EXCEED 30 SECONDS.)

4. Allow micro catheter tip to cool in air or saline prior to removing mandrel.
5. Remove mandrel from micro catheter and discard. Multiple shaping is not recommended.
6. Inspect the tip for any damage that may have resulted from steam shaping the micro catheter tip. If any damage is found, do not use the micro catheter.

Useable Length (cm)	Dead Space Volume (ml)	Maximum Injection Pressure [psi(kPa)]	100% Saline (ml/sec)	50% / 50% Contrast-76 / Saline (ml/sec)	100% Contrast-76 (ml/sec)
150	0.48	100 (689)	0.8	0.4	0.08
		300 (2,068)	1.6	1.1	0.3

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COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Ashni De Alwis, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

Local Appendix to IFU Russia Reg.
Description of documents

Is a true and exact copy of the correct and complete original document.

M. M. H.

Signature of document custodian (affiant)

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California

County of Orange


Subscribed and sworn to (or affirmed) before me on this 2th day of November

20 17 by Ashni De Alwis, proved

to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Steve Oberlander
Signature

(Seal)

Appendix to the Instruction for use

1. Name of the medical device	Micro Catheter, in various versions Echelon, Rebar, Apollo, Marathon, Orion 21:
	I. Micro Catheter, in various versions Echelon, Rebar, Apollo, Marathon, Orion 21: 1. Micro Catheter Echelon, versions: 1.1. Micro Catheter Echelon-10, tip shape 0°, consisting of: - Micro Catheter Echelon-10, tip shape 0° - Mandrel with pointed end 1.2. Micro Catheter Echelon-10, tip shape 45°, consisting of: - Microcatheter Echelon-10, tip shape 45° - Mandrel with pointed end - Split introducer 1.3. Micro Catheter Echelon-10, tip shape 90°, consisting of: - Micro Catheter Echelon-10, tip shape 90° - Mandrel with pointed end - Split introducer 1.4. Micro Catheter Echelon-14, tip shape 0°, consisting of: - Micro Catheter Echelon-14, tip shape 0° - Mandrel with pointed end 1.5. Micro Catheter Echelon-14, tip shape 45°, consisting of: - Micro Catheter Echelon-14, tip shape 45° - Mandrel with pointed end - Split introducer 1.6. Micro Catheter Echelon-14, tip shape 90°, consisting of: - Micro Catheter Echelon-14, tip shape 90° - Mandrel with pointed end - Split introducer 2. Micro Catheter Rebar, versions: 2.1. Micro Catheter Rebar-10, consisting of: - Micro Catheter Rebar-10 - Mandrel with pointed end 2.2. Micro Catheter Rebar-14, consisting of: - Micro Catheter Rebar -14

Available upon request

Appendix to the Instruction for use
Name of Medical Device

- Mandrel with pointed end
- 2.3. Micro Catheter Rebar -18, length 130 cm, consisting of:
 - Micro Catheter Rebar -18, length 130 cm
 - Mandrel with pointed end
- 2.4. Micro Catheter Rebar-18, length 153 cm, consisting of:
 - Micro Catheter Rebar -18, length: 153 cm
 - Mandrel with pointed end
- 2.5. Micro Catheter Rebar -027, length 130 cm, consisting of:
 - Micro Catheter Rebar -027, length 130 cm
 - Mandrel with pointed end
- 2.6. Micro Catheter Rebar -027, length 145 cm, consisting of:
 - Micro Catheter Rebar -027, length 145 cm
 - Mandrel with pointed end
- 3. Micro Catheter Apollo Onyx Delivery, versions:
 - 3.1. Micro Catheter Apollo Onyx Delivery, the length of the detachable tip is 1.5 cm, consisting of:
 - Microcatheter Apollo, length of detachable tip 1.5 cm
 - adapter for syringe
 - divided introducer
 - 3.2. Micro Catheter Apollo Onyx Delivery, the length of the detachable tip is 3.0 cm, consisting of:
 - Micro Catheter Apollo, length of detachable tip 3,0 cm
 - adapter for syringe
 - divided introducer
 - 3.3. Micro Catheter Apollo Onyx Delivery, the detachable tip length of 5.0 cm, consisting of:
 - Micro Catheter Apollo, length of detachable tip 5,0 cm
 - adapter for syringe
 - divided introducer
- 4. Micro Catheter Marathon Flow Directed without probe, consisting of:
 - Micro Catheter Marathon
 - Mandrel
- 5. Micro Catheter Orion 21, consisting of:
 - Micro Catheter Orion
 - Mandrel for form setting
- II. ACCESSORIES:

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	1. Adapter to the Micro Catheter Echelon 2. Adapter for syringe to Micro Catheter Apollo
2.1 With regards to the manufacturer of the medical device - full name and, if available, abbreviated name, including company name, legal form of the legal entity, location address or last name, first name and (if available) patronymic, details of the identity document, address of residence of the individual entrepreneur, as well as phone numbers and (if available) e-mail address of the legal entity or individual entrepreneur;	Developer: Micro Therapeutics, Inc. Dba Ev3 Neurovascular 9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700 Manufacturer: Micro Therapeutics, Inc. Dba Ev3 Neurovascular 9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700 Holder of Registration Certificate in Russia (name, address, tel, fax, e-mail, website) Micro Therapeutics, Inc. Dba Ev3 Neurovascular 9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700 Addresses of Manufacturing Sites for medical device: Micro Therapeutics, Inc. Dba Ev3 Neurovascular 9775 Toledo Way, Irvine, CA 92618, USA +1-949-837-3700 Medtronic Mexico S.de R.L.de CV, Av. Paseo Cucapah, 10510 El Lago, C.P. 22210, Tijuana, Baja California, MEXICO
2.2 With respect to authorized representative of the manufacturer	Authorized Representative of Manufacturer: “Medtronic” LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya., 10 Tel: (495) 580-73-77, Fax: (495) 580-73-78 www.medtronic.ru

Available upon request

Appendix to the Instruction for use
Name of Medical Device

of the medical device - full and, if available, abbreviated name, including company name, legal form of the legal entity, address (location), telephone numbers and (if available) e-mail address of the legal entity	
3. Indication of the medical device with potential customer (for example, a health professional)	For selective infusion physician-specified therapeutic agents or contrast media into peripheral or neuro vasculature
4. Functional characteristics and indication of the medical device	Indications: Echelon The Echelon Micro Catheter is intended to access peripheral and neuro vasculature for the controlled selective infusion of physician-specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media. Rebar The Rebar Micro Catheter is intended for the controlled selective infusion of physician-specified therapeutic agents or contrast media into the vasculature of the peripheral and neuro anatomy. Apollo The Apollo Onyx Delivery Micro Catheter is intended to access the neuro vasculature for controlled selective infusion of physician-specified therapeutic agents such as embolization materials and diagnostic materials such as contrast media. Marathon

Available upon request

The Marathon Micro Catheter is intended to access the peripheral and neuro vasculature for the controlled selective infusion of physician- specified therapeutic agents such as embolization materials and of diagnostic materials such as contrast media.

Orion 21

This device is intended for the controlled selective infusion of physician-specified therapeutic agents or contrast media into the vasculature of the peripheral and neuro anatomy.

Echelon

The Echelon Micro Catheter is an endhole, single-lumen catheter designed to be introduced over a steerable guidewire into the vasculature. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a semi-rigid proximal shaft which transitions into the flexible distal shaft to facilitate the advancement of the catheter in the anatomy. Dual radiopaque markers at the distal end facilitate fluoroscopic visualization. The outer surface of the catheter is coated to increase lubricity.

Rebar

The Rebar Micro Catheter is an endhole, single-lumen catheter designed to be introduced over a steerable guidewire into the vasculature. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a semi-rigid proximal shaft which transitions into the flexible distal shaft to facilitate the advancement of the catheter in the anatomy. Single or dual radiopaque markers at the distal end facilitate fluoroscopic visualization. The outer surface of the catheter is coated to increase lubricity.

Apollo


The Apollo Onyx Delivery Catheter is a single-lumen, endhole catheter designed for the superselective infusion of physician-specified therapeutic agents such as embolization materials and diagnostic materials such as contrast media in tortuous, distal vessels. The catheter has a semi-rigid proximal shaft and a highly flexible distal shaft to facilitate the advancement of the catheter in the anatomy. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The outer surfaces of the catheter are coated to increase lubricity. The Apollo catheter is designed to facilitate catheter retrieval in the event the catheter becomes entrapped within the vasculature. The distal section of the catheter incorporates a detachment zone that allows detachment of the distal tip when the force required to extract the catheter exceeds the force to detach the tip. The catheter has 2 radiopaque marker bands to visualize the position of the catheter and the detach zone area:

Available upon request

Appendix to the Instruction for use
Name of Medical Device

- Proximal to the detach zone
- At the distal end of the catheter

A	B Detachable Tip Length	C	D
0.5 mm	1.5 cm / 15 mm	1.25 mm	1.9 F/0.63 mm
0.5 mm	3 cm / 30 mm	1.25 mm	1.9 F/0.63 mm
0.5 mm	5 cm / 50 mm	1.25 mm	1.9 F/0.63 mm

1. Proximal
2. Detach Point
3. Distal

The Apollo is packaged with a Syringe Adapter. This device, attached to a 1 mL syringe filled with Onyx, will reduce the dead space within the micro catheter luer hub. Reducing the dead space within the hub is intended to minimize the potential mixing of Onyx and DMSO in the hub of the catheter during connection and injection.

Marathon

The Marathon Micro Catheters are single-lumen, endhole catheters designed for the subselective infusion of physician specified therapeutic agents or contrast media in tortuous, distal vessels. The catheter has a semi-rigid proximal shaft and a highly flexible distal shaft to facilitate the advancement of the catheter in the anatomy. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The catheter has a radiopaque marker at the distal end to facilitate fluoroscopic visualization. The outer surfaces of the catheter are coated to increased lubricity. The stylet accompanying the catheter is used to increase the rigidity of the distal section during introduction into the guiding catheter.

Available upon request

Appendix to the Instruction for use
Name of Medical Device

Orion 21

The Orion Micro Catheter is a variable stiffness, single lumen micro catheter designed to be introduced over a steerable guidewire for accessing small, tortuous vascular areas. The outer surface of the micro catheter's distal segment is coated with a hydrophilic material to provide lubricity during use. The proximal end of the micro catheter incorporates a standard luer adapter to facilitate the attachment of accessories. The micro catheter also incorporates a PTFE liner to facilitate movement or introduction of devices passed through its lumen. The Orion Micro Catheter has dual radiopaque markers at the distal tip to facilitate fluoroscopic visualization. The distal tip of the micro catheter is shapeable. The Orion micro catheter has a Fluoro-Safe marker at 88 cm from the distal tip of the microcatheter.

The Orion Micro Catheter is compatible with DMSO.


Key functional parts and components are summarized in Figure below and the materials section listed below. The medical device does not contain any software or electronic components.

Echelon

Echelon

Available upon request


Appendix to the Instruction for use
Name of Medical Device

Micro Catheter Echelon-10, tip shape 0°

Available upon request


Appendix to the Instruction for use
Name of Medical Device

Micro Catheter Echelon-10, tip shape 45°

Available upon request


Appendix to the Instruction for use
Name of Medical Device

Micro Catheter Echelon-10, tip shape 90°

Available upon request


Appendix to the Instruction for use
Name of Medical Device

Micro Catheter Echelon-14, tip shape 0°

Available upon request


Appendix to the Instruction for use
Name of Medical Device

Micro Catheter Echelon-14, tip shape 45°

Available upon request

Appendix to the Instruction for use
Name of Medical Device


Micro Catheter Echelon-14, tip shape 90°


Rebar

Available upon request

Appendix to the Instruction for use
Name of Medical Device


Apollo

Marathon

Available upon request

Appendix to the Instruction for use
Name of Medical Device

Orion 21

Available upon request

Appendix to the Instruction for use
Name of Medical Device

5. Risks of using the medical device, contraindications, expected adverse effects associated with the use of the medical device for intended purpose	Contraindications: Echelon • Not intended for use in the coronary vasculature. • Contraindicated when, in the medical judgment of the physician, such a procedure may compromise the patient's condition. • Contraindicated for neonatal and pediatric use. Rebar • It is not intended for use in the coronary vasculature. • The Rebar Micro Catheter is contraindicated when, in the medical judgment of the physician, such a procedure may compromise the patient's condition. • The Rebar Micro Catheter is contraindicated for neonatal and pediatric use. Apollo

Available upon request

Appendix to the Instruction for use
Name of Medical Device

- The Apollo Onyx Delivery Micro Catheter is contraindicated when, in the medical judgment of the physician, use of such product may compromise the patient's condition.
- Not intended for use in the coronary vasculature.
- The Apollo Onyx Delivery Micro Catheter is contraindicated for neonatal and pediatric use.

Marathon

- The Marathon Micro Catheter is contraindicated when, in the medical judgment of the physician, use of such product may compromise the patient's condition.
- Not intended for use in the coronary vasculature.
- The Marathon Micro Catheter is contraindicated for neonatal and pediatric use.

Orion 21

Not intended for delivering liquid embolics (e.g. Onyx) into the neurovasculature. If the micro catheter tip becomes entrapped within an embolic mass it may be difficult to remove the micro catheter from the patient. Entrapped devices should not be implanted and may need to be surgically resected.

Adverse events (Potential complications)

Potential complications include, but are not limited to:

Echelon

- Hematoma at the site of entry
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Vascular thrombosis
- Thrombolytic episodes
- Neurological deficits including stroke and death
- Ischemia

Rebar

- Hematoma at the site of entry
- Vessel perforation

Available upon request

Appendix to the Instruction for use
Name of Medical Device

- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Vascular thrombosis
- Thrombolytic episodes
- Neurological deficits including stroke and death
- Ischemia

Apollo

- Puncture site hematoma
- Vessel perforation
- Vessel spasm
- Pain and tenderness
- Thrombolytic episodes
- Neurological deficits including stroke and death

Marathon

- Puncture site hematoma
- Vessel perforation
- Vessel spasm
- Hemorrhage
- Pain and tenderness
- Thrombolytic episodes
- Neurological deficits including stroke and death
- Vascular thrombosis

Orion 21

- Distal emboli (air, foreign material, tissue, or thrombus)
- Hemorrhage
- Ischemia
- Neurological deficits including stroke and death
- Hematoma
- Infection
- Pain at insertion site

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	Vessel spasm, dissection, perforation or injury In accordance with BS EN ISO 14971:2012 and internal Quality procedures, a risk assessment of the Microcatheters was performed taking into consideration the possible hazards and associated risks related to the product design and clinical usage of the device. The results of the overall risk analysis demonstrate that the associated risks of using the Device are adequately mitigated through the device design and Instructions for Use.									
6. Technical characteristics of the medical device	General characteristics									
		Overall length, cm	Distal OD, mm	Proximal OD, mm	ID, mm	Marker band width, mm	Working length (effective length), cm	Distance from distal end to first radiopaque marker. cm	Distance Between Marker Bands, mm	Mass of microcatheter, g
	Echelon-10									
	tip shape 0°	155 ± 5%	<0,66	< 0,76	0,419 ± 5%	0,5 ± 0,01	150 ± 5%	0,5 ± 5%	29,65 ± 1	3 ± 0,5
	tip shape 45°	155 ± 5%	<0,66	< 0,76	0,419 ± 5%	0,5 ± 0,01	150 ± 5%	0,5 ± 5%	29,65 ± 1	3 ± 0,5
	tip shape 90°	155 ± 5%	<0,66	< 0,76	0,419 ± 5%	0,5 ± 0,01	150 ± 5%	0,5 ± 5%	29,65 ± 1	3 ± 0,5
	Echelon-14									
	tip shape 0°	155 ± 5%	<0,69	< 0,86	0,419 ± 5%	0,5 ± 0,01	150 ± 5%	0,5 ± 5%	29,65 ± 1	3 ± 0,5
	tip shape 45°	155 ± 5%	<0,69	< 0,86	0,419 ± 5%	0,5 ± 0,01	150 ± 5%	0,5 ± 5%	29,65 ± 1	3 ± 0,5
	tip shape 90°	155 ± 5%	<0,69	< 0,86	0,419 ± 5%	0,5 ± 0,01	150 ± 5%	0,5 ± 5%	29,65 ± 1	3 ± 0,5
	Rebar									

Available upon request

Appendix to the Instruction for use
Name of Medical Device

Rebar -10	159 ± 5%	0,57 ± 0,05	0,77 ± 0,05	0,40 ± 0,05	0,5 ± 0,01	153 ± 2,5	0,5 ± 5%	29,5 ± 1	2 ± 0,5
Rebar -14	159 ± 5%	0,63 ± 0,05	0,8 ± 0,05	0,40 ± 0,05	0,5 ± 0,01	130 ± 2,5	0,5 ± 5%	29,5 ± 1	2 ± 0,5
Rebar -18									
length (cm) 130	136,5 ± 5%	0,93 ± 0,05	0,99 ± 0,05	0,50 ± 0,05	0,5 ± 0,01	130 ± 2,5	0,5 ± 5%	N/A	2 ± 0,5
length (cm) 153	159 ± 5%	0,93 ± 0,05	0,99 ± 0,05	0,50 ± 0,05	0,5 ± 0,01	153 ± 2,5	0,5 ± 5%	29,5 ± 1	2 ± 0,5
Rebar - 027									
length (cm) 130	135 ± 5%	≤ 1,09	≤ 1,04	≥ 0,66	0,5 ± 0,01	130 ± 2,5	0,5 ± 5%	N/A	2 ± 0,5
length (cm) 145	150 ± 5%	≤ 1,09	≤ 1,04	≥ 0,66	0,5 ± 0,01	145 ± 2,5	0,5 ± 5%	N/A	2 ± 0,5
Apollo									
length of the tip (cm) 1,5	172 ± 3	0,51 ± 0,07	0,91 ± 0,05	0,38 ± 0,02	0,5 ± 0,01	165 ± 0,4	0,5 ± 5%	N/A	2 ± 0,5
length of the tip (cm) 3,0	172 ± 3	0,51 ± 0,07	0,91 ± 0,05	0,38 ± 0,02	0,5 ± 0,01	165 ± 0,4	0,5 ± 5%	N/A	2 ± 0,5
length of the tip (cm) 5,0	172 ± 3	0,51 ± 0,07	0,91 ± 0,05	0,38 ± 0,02	0,5 ± 0,01	165 ± 0,4	0,5 ± 5%	N/A	2 ± 0,5
Marathon,									
With mandrel	170 ± 3	< 0,55	< 0,99	Prox > 0,4 Dist > 0,3	0,3 ± 0,01	165 ± 2,5	0,5 ± 5%	N/A	3 ± 0,5
Orion 21									
Orion 21,	157 ± 3	< 0,93	< 0,88	> 0,55	0,5 ± 0,01	150,0 ± 2,5	0,5 ± 5%	30 ± 1	3 ± 0,5

Available upon request

Appendix to the Instruction for use
Name of Medical Device

7. Description of accessories, medical devices or non-medical devices that are intended for use in combination with the medical device (if applicable)	Adapter to the microcatheter Echelon		Overall length, mm: 37 Thread length, mm: 5 Inner diameter of thread, mm: 4 Diameter of the expanding part, mm: 10 Inner diameter of the metal part, mm: 0,6 Outer diameter of the metal part, mm: 0,9 Hardness of the metal part, HB: not less 50 The roughness of the metal part Ra: no more 0,16 mkm Weight, g: 0,69
	Adapter for syringe to Apollo microcatheter		Overall length, mm: 25 Thread length, mm: 5 Inner diameter of thread, mm: 4 Diameter of the expanding part, mm: 10 Diameter of spout inner, mm: 0,4 Weight, g: 0,63
8. Information on the presence in the medical device of a medicinal product for human use, human and/or animal materials	Not applicable		
9. Information on the procedure of installation, adjustment, calibration and	Not applicable		

Available upon request

Appendix to the Instruction for use
Name of Medical Device

other actions necessary to put the medical device into operation	
10. Requirements for the premises in which the installation of the medical device shall be carried out, as well as requirements for training or qualification of persons who install the medical device (if applicable)	Not applicable
11. Information necessary to verify the correct installation of the medical device and its readiness for safe operation, including:	Not applicable
a. Maintenance and its frequency, including cleaning and disinfection of the medical device;	Not applicable
b. List of information	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

provided by the manufacturer of the medical device, keys, access passwords, programs required for installation, adjustment, operation and maintenance of the medical device;	
c. List of consumables (components, reagents), as well as the procedure for their use and replacement;	Not applicable
d. Need for calibration to ensure the proper and safe operation of the medical device during its lifetime;	Not applicable
e. Risk minimization methods associated with the installation, calibration or maintenance of the	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

medical device;				
f. Information on installation, adjustment, calibration and other actions necessary for putting the medical device into operation and its correct operation (application);	Not applicable			
g. Information on the list of the main characteristics for the operation (application) of the medical device, transportation and storage conditions (for example, temperature and humidity, lighting and other characteristics)	Individual characteristics of catheters:			
	Echelon-10			
		tip shape 0°	tip shape 45°	tip shape 90°
	Length of outer coating, cm	100 ± 5	100 ± 5	100 ± 5
	Tensile, N	≥ 3	≥ 3	≥ 3
	Diameter of marker band, mm	0,5 ± 0,01	0,5 ± 0,01	0,5 ± 0,01

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	The nominal size of the catheter, cm	$150 \pm 5\%$	$150 \pm 5\%$	$150 \pm 5\%$
	The number of revolutions per 360° without kinking	> 5	> 5	> 5
	The maximum guide wire diameter, mm	$< 0,35$	$< 0,35$	$< 0,35$
	Distance from the distal end to the marker on the shaft, cm	88 ± 2	88 ± 2	88 ± 2
Echelon-14				
	tip angle 0°	tip angle 45°	tip angle 90°	

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	Length of outer coating, cm	100 ± 5	100 ± 5	100 ± 5
	Tensile, N	≥ 3	≥ 3	≥ 3
	Diameter of marker band, mm	$0,5 \pm 0,01$	$0,5 \pm 0,01$	$0,5 \pm 0,01$
	The nominal size of the catheter, cm	$150 \pm 5\%$	$150 \pm 5\%$	$150 \pm 5\%$
	The number of revolutions per 360° without kinking	> 5	> 5	> 5
	Distance from the distal end to the marker on the shaft, cm	88 ± 2	88 ± 2	88 ± 2
Rebar				
	Rebar -10	Rebar -14		
Tensile, N	≥ 3	≥ 3		
Length of outer coating, cm	100 ± 5	100 ± 5		

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	The nominal size of the catheter, cm	$153 \pm 5\%$	$130 \pm 5\%$
	Diameter of marker band, mm	$0,5 \pm 0,01$	$0,5 \pm 0,01$
	The number of revolutions per 360° without kinking	> 5	> 5
	Rebar -18		
		Length 130 cm	Length 153 cm
	Tensile, N	$\geq 4,5$	$\geq 4,5$
	Length of outer coating, cm	$130 \pm 2,5$	$153 \pm 2,5$
	The nominal size of the catheter, cm	$130 \pm 5\%$	$153 \pm 5\%$
	Diameter of marker band, mm	$0,5 \pm 0,01$	$0,5 \pm 0,01$
	The number of revolutions per 360° without kinking	> 5	> 5
	Rebar -27		
		Length 130 cm	Length 145 cm

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	Tensile, N	$\geq 4,5$	$\geq 4,5$
	Length of outer coating, cm	$130 \pm 2,5$	145 ± 2.5
	The nominal size of the catheter, cm	$130 \pm 5\%$	$145 \pm 5\%$
	Diameter of marker band, mm	$0,5 \pm 0,01$	$0,5 \pm 0,01$
	The number of revolutions per 360° without kinking	> 5	> 5

Apollo			
	Tip length (cm) 1,5	Tip length (cm) 3,0	Tip length (cm) 5,0
Length of outer coating, cm	100 ± 5	100 ± 5	100 ± 5
Tensile, N	$\geq 4,5$	$\geq 4,5$	$\geq 4,5$
Diameter of fluorescence marker, mm	$\geq 0,99$	$\geq 0,99$	$\geq 0,99$
Diameter of marker band, mm	$0,5 \pm 0,01$	$0,5 \pm 0,01$	$0,5 \pm 0,01$

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	The nominal size of the catheter, cm	$170 \pm 5\%$	$170 \pm 5\%$	$170 \pm 5\%$
	The number of revolutions per 360° without kinking	> 5	> 5	> 5
	Location of fluorescence marker	Distal edge of marker is $88 \pm 2\text{cm}$ from of the proximal edge of the distal marker band of the catheter.		
Marathon				
	Length of outer coating, cm	100 ± 5		
	Tensile, N	$\geq 4,5$		
	Diameter of fluorescence marker, mm	$\leq 0,9$		
	Diameter of marker band, mm	$\leq 0,5$		
	The nominal size of the catheter, cm	170 ± 3		
	The number of revolutions per 360° without kinking	> 5		
	Location of fluorescence marker	Distal edge of marker is $88 \pm 2\text{cm}$ from of the proximal edge of the distal marker band of the catheter.		
Orion 21				

Available upon request

Appendix to the Instruction for use
Name of Medical Device

Length of outer coating, cm	100 ± 5
Tensile, N	≥ 4,5
Diameter of fluorescence marker, mm	≤ 0,9
Diameter of marker band, mm	≤ 0,5
The nominal size of the catheter, cm	160 ± 3
The number of revolutions per 360° without kinking	> 5
Location of fluorescence marker	Distal edge of marker is 88 ± 2cm from of the proximal edge of the distal marker band of the catheter.

The medical device is intended for use under standard conditions. A medical device is allowed to be used only by specialists who have passed possession and have relevant work experience.

Operating conditions in the health facility:

t°	10 to 40 C°
Pressure	700-1060 gPa
Humidity	30 -75 %

Operating conditions in the internal environment of the body:

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	<table><tr><td>t°</td><td>32 to 42 C°</td></tr><tr><td>Pressure</td><td>700-1060 gPa</td></tr><tr><td>Humidity</td><td>100 %</td></tr></table>	t°	32 to 42 C°	Pressure	700-1060 gPa	Humidity	100 %					
	t°	32 to 42 C°										
	Pressure	700-1060 gPa										
	Humidity	100 %										
	Storage Conditions Store with other sterile devices in normal operating room temperatures in a dry place and in a manner that protects the integrity of the package and the sterile barrier. Rotate inventory so that device is used before the expiry date on the package.											
	<table><tr><td>Air temperature, °C</td><td>+10 - +32</td></tr><tr><td>Air humidity (relative), non-condensing, %</td><td>< 75</td></tr><tr><td>Atmosphere pressure, hPa</td><td>860 - 1060</td></tr></table>	Air temperature, °C	+10 - +32	Air humidity (relative), non-condensing, %	< 75	Atmosphere pressure, hPa	860 - 1060					
	Air temperature, °C	+10 - +32										
	Air humidity (relative), non-condensing, %	< 75										
	Atmosphere pressure, hPa	860 - 1060										
	Transport Conditions											
<table><tr><td>Air temperature, °C</td><td>-35 +40</td></tr><tr><td>Air humidity (relative), non-condensing, %</td><td>30 - 95</td></tr><tr><td>Atmosphere pressure, hPa</td><td>700 - 1060</td></tr></table>	Air temperature, °C	-35 +40	Air humidity (relative), non-condensing, %	30 - 95	Atmosphere pressure, hPa	700 - 1060						
Air temperature, °C	-35 +40											
Air humidity (relative), non-condensing, %	30 - 95											
Atmosphere pressure, hPa	700 - 1060											
h. List of national standards applied by the manufacturer												
International Standards:												
<table><tr><th>Number</th><th>Standard Title</th></tr><tr><td>ISO 594-1:1986E</td><td>Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment -- Part 1: General requirements</td></tr><tr><td>ASTM F 2052-06</td><td>Standard Test Method for Measurement of Magnetically Induced Displacement Force on Passive Implants in the Magnetic Resonance Environment</td></tr><tr><td>ASTM F 2096-04</td><td>Standard Test Method for Detecting Gross Leaks in Porous Medical Packaging by Internal Pressurization (Bubble Test)</td></tr><tr><td>ASTM F 2182-02a</td><td>Standard Test Method for Measurement of Radio Frequency Induced Heating Near Passive Implants During MR Imaging</td></tr><tr><td>ASTM F 2119-07</td><td>Standard Test Method for Evaluation of MR Imaging Artifacts from Passive Implants</td></tr></table>	Number	Standard Title	ISO 594-1:1986E	Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment -- Part 1: General requirements	ASTM F 2052-06	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Passive Implants in the Magnetic Resonance Environment	ASTM F 2096-04	Standard Test Method for Detecting Gross Leaks in Porous Medical Packaging by Internal Pressurization (Bubble Test)	ASTM F 2182-02a	Standard Test Method for Measurement of Radio Frequency Induced Heating Near Passive Implants During MR Imaging	ASTM F 2119-07	Standard Test Method for Evaluation of MR Imaging Artifacts from Passive Implants
Number	Standard Title											
ISO 594-1:1986E	Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment -- Part 1: General requirements											
ASTM F 2052-06	Standard Test Method for Measurement of Magnetically Induced Displacement Force on Passive Implants in the Magnetic Resonance Environment											
ASTM F 2096-04	Standard Test Method for Detecting Gross Leaks in Porous Medical Packaging by Internal Pressurization (Bubble Test)											
ASTM F 2182-02a	Standard Test Method for Measurement of Radio Frequency Induced Heating Near Passive Implants During MR Imaging											
ASTM F 2119-07	Standard Test Method for Evaluation of MR Imaging Artifacts from Passive Implants											

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	ASTM F 2503-05	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment
	ISO 15223-1:2012	Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
	EN 1041:2008	Information supplied by the manufacturer of medical devices
	EN ISO 14971:2012	Medical devices – Application of risk management to medical devices
	EN ISO 10993-1 :2009	Biological evaluation of medical devices – Part 1: Evaluation and testing
	EN ISO 10993-3:2003	Biological Evaluation of Medical Devices Part 3: Tests for Genotoxicity, Carcinogenicity and Reproductive Toxicity
	EN ISO 10993-4:2002; A1:2006	Biological Evaluation of Medical Devices-Part 4: Hemocompatibility
	EN ISO 10993-5:2009	Biological Evaluation of Medical Devices -Part 5: Tests for In Vitro Cytotoxicity
	EN ISO 10993-6:2007	Biological Evaluation of Medical Devices -Part 6: Implantation
	EN ISO 10993-7:2008/Cor. 1:2009	Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals
	EN ISO 10993-10:2002; A1:2006	Biological Evaluation of Medical Devices -Part 10: Tests for Irritation and Delayed-type Hypersensitivity
	EN ISO 10993-11:2006	Biological Evaluation of Medical Devices -Part 11: Test for Systemic Toxicity
	EN ISO 10993-12:2007	Biological Evaluation of Medical Devices -Part 12: Sample Preparation and Reference Materials
	ISO 10993-13:2010	Biological evaluation of medical devices – Part 13: Identification and quantification of degradation products from polymeric medical devices
	ISO 10993-17:2002	Biological evaluation of medical devices – Part 17: Establishment of allowable limits for leachable substances
	ISO 10993-18:2005	Biological evaluation of medical devices – Part 18: Chemical characterization of materials

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	ISO 15223-1:2012	Medical devices - Symbols to be used with Medical Device Labels, Labeling and Information to be Supplied –Part 1: General Requirements												
12. Information on the sterile condition of the medical device, method of its sterilization and the procedure for dealing with a sterile package violation (if the medical device is delivered sterile), or if the medical device is not sterile, indication of the need for sterilization before use	Microcatheters are sterilized using medical equipment, using the standard ethylene oxide sterilization cycle (EtO), which results in a guaranteed level of sterility (SAL) of 10⁻⁶. ANSI / AAMI / ISO 11135-1: 2007 and ANSI / AAMI / ISO 11135-2: 2008. The microcatheters are sterilized according to the parameters listed in the table below: <table><tr><td>Sterility Validation Method.</td><td>ANSI/AAMI/ISO 11135-1 2007, Sterilization of health care products – Ethylene oxide – Requirements for development, validation and routine control of a sterilization process for medical devices. Overkill approach.</td></tr><tr><td>EO Residuals</td><td>ISO 10993-7 2008 Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals.</td></tr><tr><td>Sterilization Method</td><td>100% ethylene oxide gas</td></tr><tr><td>Sterility Assurance</td><td>10⁻⁶</td></tr><tr><td>Pyrogen Tests</td><td>The device is non pyrogenic. Pyrogen testing is conducted on a lot-to-lot basis using the “Guideline for Validation of LAL Test as an End-Product Endotoxin Test for Medical Devices, (FDA 1978).”</td></tr><tr><td>Contract Sterilization</td><td>Sterigenics, Inc., Los Angeles, CA</td></tr></table>		Sterility Validation Method.	ANSI/AAMI/ISO 11135-1 2007, Sterilization of health care products – Ethylene oxide – Requirements for development, validation and routine control of a sterilization process for medical devices. Overkill approach.	EO Residuals	ISO 10993-7 2008 Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals.	Sterilization Method	100% ethylene oxide gas	Sterility Assurance	10 ⁻⁶	Pyrogen Tests	The device is non pyrogenic. Pyrogen testing is conducted on a lot-to-lot basis using the “Guideline for Validation of LAL Test as an End-Product Endotoxin Test for Medical Devices, (FDA 1978).”	Contract Sterilization	Sterigenics, Inc., Los Angeles, CA
Sterility Validation Method.	ANSI/AAMI/ISO 11135-1 2007, Sterilization of health care products – Ethylene oxide – Requirements for development, validation and routine control of a sterilization process for medical devices. Overkill approach.													
EO Residuals	ISO 10993-7 2008 Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals.													
Sterilization Method	100% ethylene oxide gas													
Sterility Assurance	10 ⁻⁶													
Pyrogen Tests	The device is non pyrogenic. Pyrogen testing is conducted on a lot-to-lot basis using the “Guideline for Validation of LAL Test as an End-Product Endotoxin Test for Medical Devices, (FDA 1978).”													
Contract Sterilization	Sterigenics, Inc., Los Angeles, CA													
13. Information	Not applicable													

Available upon request

Appendix to the Instruction for use
Name of Medical Device

on the procedure for the preparation of the medical device for its reuse, including cleaning, disinfection, packaging and, if necessary, method of re-sterilization (if the medical device is reusable), as well as criteria for the unfitness of the medical device for use	
14. Information necessary for the identification of medical devices in order to obtain a safe combination, and information on known restrictions on concomitant use of medical devices (for medical devices intended for use with other medical devices and/or supplies)	Not applicable
15. Information on the nature, type,	Not applicable.

Available upon request

Appendix to the Instruction for use
Name of Medical Device

and (if necessary) the intensity and distribution of radiation (electromagnetic, ionizing, other) of the medical device and methods for the protection of consumers and third parties from unintentional radiation during the operation of the medical device (if the medical device creates a dangerous or potentially hazardous radiation level when used as intended)	
16. Information on precautions taken in case of:	WARNINGS • Infusion pressure with this device should not exceed 600 PSI (4137 kPa). Pressures in excess of 600 PSI (4137 kPa) may result in catheter rupture, which may result in patient injury. • If flow through the micro catheter becomes restricted, do not attempt to clear the device by high pressure infusion. Either remove the micro catheter to determine the cause of the obstruction or replace it with a new micro catheter. Excessive pressure may cause catheter rupture, which may result in patient injury. • Never advance or withdraw an intraluminal device against resistance until the cause of the resistance is determined by fluoroscopy. Excessive force against resistance may result in damage to the device or vessel perforation or other patient injuries. • This device is supplied STERILE for single use only. Do not reprocess or re-sterilize. Reprocessing and re-sterilization

Available upon request

Appendix to the Instruction for use
Name of Medical Device

	increase the risks of patient infection and compromised device performance. Warnings Shaping Mandrel is not intended for use in the human body. Use only a steam source to shape the catheter tip. Do not use other heat sources. Prior to use, inspect the catheter tip for any damage that may have resulted from shaping. Do not use a catheter that has been damaged in any way. Damaged catheters may rupture causing vessel trauma or tip detachment during steering maneuvers. PRECAUTIONS • Ensure embolic material compatibility with catheter prior to use. • Prior to use, carefully examine the micro catheter and packaging to verify that they have not been damaged during shipment. • Prior to use, all accessory devices and agents should be fully prepared according to the manufacturer's instructions. • Always monitor infusion rates when using the micro catheter. • When injecting contrast for angiography, ensure that the micro catheter is not kinked or occluded. • The micro catheter has a lubricious hydrophilic coating on the outside of the catheter. It must be kept hydrated in order to be lubricious. This can be accomplished by attaching the Y-connector to a continuous saline drip. • When the infusion catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the micro catheter without observing the resultant tip response. • Because the micro catheter may be easily advanced into narrow, selective vasculature, repeatedly verify (for example by slightly withdrawing the catheter) that the micro catheter has not been advanced so far as to interfere with its removal.
a. Failures of the medical device, malfunction or deviations in the functioning that may affect the safety of the medical device, including those determined by external signs;	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

b. Effects on the operation of the medical device of external factors associated with the use of the medical device in combination with other medical devices and/or equipment, or such predictable factors as external electromagnetic fields, electrostatic discharges, radiation (electromagnetic, ionizing, other), atmospheric pressure and its differences, humidity and air temperature;	Not applicable
c. The risk of electromagnetic interference created by the medical device for other medical devices, equipment and communication	Not applicable

Available upon request

Appendix to the Instruction for use
Name of Medical Device

equipment during the conduct and evaluation of the results of diagnosis, treatment or its intended use (for example, electromagnetic radiation from a medical device that affects other equipment);	
17. Warnings and/or precautions taken by the consumer in the application of the medical device containing a medicinal product for human use, human and/or animal material, materials that are carcinogenic, mutagenic or toxic, whose possible release or washing-out leads to sensitization, allergic reaction or negatively affects the reproductive	Not applicable. The medical device does not contain materials of animal or human.

Available upon request

Appendix to the Instruction for use
Name of Medical Device

function	
18. Warnings and/or precautions taken by the consumer when disposing of the medical device, accessories and consumables used together with it (if any), including information about the infectious, microbial, environmental or physical hazard of the medical device	For disposing of the microcatheters is required to provide the mandatory implementation of all the necessary procedures in accordance with local applicable laws and regulations.
19. Information on the circumstances which should prompt the consumer to consult a health professional	A medical device used only by qualified medical personnel who have undergone specialized training and have an appropriate level of surgical experience.
20. Information on the initial release or the latest revision of the operational documentation;	Not applicable
21. Procedure	<u>After use, this product may be a potential biohazard.</u> Dispose the medical device in accordance with the requirements

Available upon request

Appendix to the Instruction for use
Name of Medical Device

and conditions for the disposal or destruction of a medical device.	of the national legislation on environment and health and/or hospital Unused, expired devices are classified as medical waste class «A» Catheters extracted from the patient's body are medical waste Class B, that is epidemiologically hazardous waste. Waste Class B subjects to mandatory disinfection (disinfection)/ disposal. The choice of the method of disinfection/ neutralization is determined by the capabilities an organization engaged in medical and/ or pharmaceutical activity and performed in the development of the treatment scheme of medical waste.
--	--

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