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
Authorized


Signature

**INSTRUCTION FOR USE
OF THE MEDICAL DEVICE**

**STERILE ANTI-ADHESIVE ABSORBABLE GEL
STERILE ANTIADHESIN[®],
1,5 g, 3 g, 4 g, 5 g, 6 g, 10 g SYRINGES.**

Manufactured by Genewel Co., Ltd., Republic of Korea

1. Name of the medical device

Antiadhesin[®] anti-adhesive, absorbable, sterile gel, 1,5 g, 3 g, 4 g, 5 g, 6 g, 10 g syringes.

2. Purpose

Antiadhesin[®], anti-adhesive, absorbable, sterile gel, 1,5 g, 3 g, 4 g, 5 g, 6 g, 10 g syringes (hereinafter referred to as Antiadhesin[®] gel, medical device) is intended for the prevention of adhesions that can occur after any types of surgery and can be manifested as the fusion of organs or tissues between each other. The therapeutic effect of the gel is related to the formation of a temporary artificial barrier between the injured tissues, ensuring effective separation of surfaces during healing. Antiadhesin[®] gel is indicated for use as an adjuvant in open and laparoscopic surgeries to reduce adhesions, reduce the spread and severity of the postoperative adhesions in the surgical area, as well as the associated serious complications such as intestinal obstruction, pain, and infertility, occurrence of which increases the risk of repeated surgical interventions.

The medical device is intended for use in healthcare facilities by qualified healthcare professionals.


3. Functional characteristics

Antiadhesin[®] gel is intended for the prevention of adhesions after surgeries on organs and tissues with serous coating: abdominal organs (stomach, intestines, spleen, liver); pelvic organs (uterus, ovaries, fallopian tubes), pleural cavity, pericardial cavity, testicular membrane, as well as during operations in the uterine cavity, on the bladder, mammary glands, in spine surgery. Antiadhesin[®] gel can be used in operations where there is a risk of soft tissue adhesion. For example, during neck surgery, when there is a high risk of scar formation throughout the neck layers. Antiadhesin[®] gel can be used as an installation gel in plastic surgery. The gel facilitates the installation of breast implants. Antiadhesin[®] gel can be used in such ophthalmic operations as strabismus surgery, endonasal dacryocystorhinostomy, as well as surgical restoration of the orbit wall using transconjunctival access to prevent the development of adhesions and ectropion of the lower eyelid in the early postoperative period.

After application of Antiadhesin[®] gel to the surgical site, it adheres to the anatomical surfaces without flowing and forms a viscous, lubricant coating or anti-adhesive film, which ensures the sliding of adjacent injured surfaces, preventing them from sticking together. Antiadhesin[®] gel is a biodegradable coating that separates the mating surfaces only during the critical phase of wound healing and postoperative adhesion, lasting for 7 days, without affecting the normal regeneration processes.

4. Main characteristics

Appearance	The medical device is a colorless clear gel in a pre-filled disposable syringe with a Luer lock, an application catheter.	
Composition	Composition of Antiadhesin [®] gel per 1 ml:	
	Sodium hyaluronate	2,50 mg
	Sodium carboxymethyl cellulose	5,00 mg
	NaCl	8,00 mg
	KCl	0,20 mg
	KH ₂ PO ₄	0,20 mg
	Na ₂ HPO ₄ · 12H ₂ O	1,15 mg
	Water for injections	1 ml
The gel is a mixture of two anionic polysaccharides: sodium hyaluronate (SH) and sodium carboxymethyl cellulose (SCMC).		

	Sodium hyaluronate (SH) Sodium hyaluronate is a natural anionic polysaccharide composed of D-glucuronic acid and D-N-acetylglucosamine linked via alternating β-1,4- and β-1,3-glycosidic bonds. SH may be detected in various tissues of the body such as connective and cartilaginous tissue, as well as in the skin, vitreous humor, and synovial fluid. SH is one of the main ingredients of the extracellular matrix; it is a high molecular viscoelastic substance with hydrophilic and non-immunogenic properties. Due to these features, SH can act as a lubricant when applied onto the injured surface. This was the basis for its use to reduce or prevent injury of the surgical area. Sodium carboxymethyl cellulose (SCMC) Sodium carboxymethyl cellulose is also an anionic polysaccharide. By moving the glycosidic hydroxyl group into the carboxymethyl group, SCMC becomes a more hydrophilic compound. Just as SH, SCMC has hygroscopic properties and can be used for pharmaceutical, cosmetic, and food purposes as an adjuvant, lubricant, and stabilizing agent.																															
Technical and functional characteristics	Antiadhesin[®] gel contains highly purified sodium salt of hyaluronic acid obtained without the use of genetically modified microorganisms. It is a sterile, biocompatible, biodegradable, bioabsorbable, non-pyrogenic medical device. It does not contain antigens. The medical device is biologically safe: non-pyrogenic, has no cytotoxic, genotoxic, irritant, and sensitizing effects, acute and subacute toxicity.  <table><tr><th rowspan="2">Name of the device</th><th rowspan="2">Weight of content Deviation: $\pm 7\%$</th><th colspan="2">Catheter tube</th><th colspan="3">Syringe</th></tr><tr><th>Diameter (a, mm)</th><th>Length (b, mm)</th><th>Volume</th><th>Diameter (c, mm)</th><th>Length (d, mm)</th></tr><tr><td rowspan="6">Antiadhesin[®] gel</td><td>1,5 g</td><td rowspan="6">1,7$\pm$0,085</td><td>80$\pm$4,0</td><td rowspan="3">5 ml</td><td rowspan="3">15$\pm$0,10</td><td rowspan="3">76,90$\pm$0,2</td></tr><tr><td>3 g</td><td colspan="2" rowspan="2">70$\pm$3,5</td></tr><tr><td>4 g</td></tr><tr><td>5 g</td><td rowspan="3">10 ml</td><td rowspan="3">18$\pm$0,10</td><td rowspan="3">97,70$\pm$0,2</td></tr><tr><td>6 g</td><td colspan="2" rowspan="2">70$\pm$3,5</td></tr><tr><td>10 g</td></tr></table>	Name of the device	Weight of content Deviation: $\pm 7\%$	Catheter tube		Syringe			Diameter (a, mm)	Length (b, mm)	Volume	Diameter (c, mm)	Length (d, mm)	Antiadhesin [®] gel	1,5 g	1,7 $\pm$ 0,085	80 $\pm$ 4,0	5 ml	15 $\pm$ 0,10	76,90 $\pm$ 0,2	3 g	70 $\pm$ 3,5		4 g	5 g	10 ml	18 $\pm$ 0,10	97,70 $\pm$ 0,2	6 g	70 $\pm$ 3,5		10 g
Name of the device	Weight of content Deviation: $\pm 7\%$			Catheter tube		Syringe																										
		Diameter (a, mm)	Length (b, mm)	Volume	Diameter (c, mm)	Length (d, mm)																										
Antiadhesin [®] gel	1,5 g	1,7 $\pm$ 0,085	80 $\pm$ 4,0	5 ml	15 $\pm$ 0,10	76,90 $\pm$ 0,2																										
	3 g		70 $\pm$ 3,5																													
	4 g																															
	5 g		10 ml	18 $\pm$ 0,10	97,70 $\pm$ 0,2																											
	6 g					70 $\pm$ 3,5																										
	10 g																															

Areas of application	Antiadhesin® gel is intended for the prevention of adhesions after all types of surgical interventions (either open or laparoscopic) on organs and tissues where there is a risk of postoperative adhesions: in general surgery, neurosurgery (spinal cord and brain membranes), otolaryngology (middle ear cavity, nasal cavity, and paranasal sinuses), ophthalmology, thoracic surgery, cardiac surgery, abdominal surgery, orthopedic surgery, obstetrics, and gynecology (including the uterine cavity), urology, plastic surgery.
Toxicity	Non-toxic
Possible application-related risks	Since safety and efficacy of the medical device use in combination with other anti-adhesive barriers, hemostatic agents, medicinal products; in patients with lymphatic and coagulation disorders; in patients with severe hepatic and renal impairment; in patients with malignant tumors; in patients with autoimmune or immune disorders, in pediatric patients, in pregnant and breastfeeding women have not been established, it should be used with caution. It is recommended to avoid pregnancy until the end of the first menstrual cycle after using the device. When using Antiadhesin® gel, reactions to foreign bodies, as well as to any implanted material, may occur.
Sterilization	The gel is sterilized by membrane filtration. The packaged medical device is sterilized using ethylene oxide.
Compliance with standards	Applied standards EN 556-1:2001, EN 556-2:2015, EN ISO 20417:2021, EN ISO 10993-1:2018, EN ISO 10993-3:2014, EN ISO 10993-5:2009, EN ISO 10993-6:2016, EN ISO 10993-7:2008, ISO 10993-10:2021, EN ISO 10993-11:2018, EN ISO 11135-1:2014, ISO 11135:2014 (equivalent to GOST ISO 11135:2017), EN ISO 11607-1:2019, EN ISO 11607-2:2019, EN ISO 13485:2016, EN ISO 14155:2020, MEDDEV 2.7.1 version 4 (July 2016), EN ISO 13408-1:2015, EN ISO 13408-2:2018, ISO 14644-1:2015, ISO 14644-2:2015, ISO 14644-3:2019, EN ISO 14971:2019, ICH Q1A (R2), ASTM F88/F88M:21, ASTM F1980-21, ASTM F838-20, ISTA 2A: 2012.
Packaging	Antiadhesin® gel is packed into a pre-filled disposable 1.5 g, 3 g, 4 g, 5 g, 6 g, 10 g syringes; it is sterile and ready for use. The syringe and applicator catheter are sealed in a blister made of a transparent polymer film; the blister with instruction for use and an adhesive label are placed in a carton.

5. Therapeutic indications

Prevention of adhesions after any type of surgeries on organs and tissues where there is a risk of postoperative adhesions.

6. Contraindications

Do not prescribe to patients with established hypersensitivity reactions to hyaluronic acid, carboxymethyl cellulose, and salts thereof.

Use of Antiadhesin® gel is contraindicated in case of an overt infection or contamination at the surgical site.

7. Undesirable effects

Allergic reactions, infection, pain, inflammation.


8. Method of administration

The medical device should be kept at room temperature before use.


Device preparation for operation

- Open the package in sterile conditions *immediately before use*, remove the syringe and the application catheter.

- Remove the rubber cap from the Luer lock by turning the syringe body.
Warning: avoid contact of the syringe tip with foreign objects after removal of the rubber cap.

- Connect the application catheter in the sheath to the Luer lock tip of the syringe and turn it clockwise until a fixed connection is achieved. Remove the sheath immediately before use.

Procedure for device operation

- The target management area should be as dry as possible. It is necessary to remove excess liquid using an aspirator, as well as all hemostatic gauze sponges and swabs inserted after the operation.
- Ensure proper hemostasis: make sure that bleeding in the surgical site is completely stopped.
- Evenly cover all the anatomical surfaces, where adhesion should be prevented, with the gel. The amount of gel to be injected depends on the specific case and area of application. One should use as much gel as it is necessary to apply it in one thin layer to the areas that require prevention of adhesions.
- Dispose of the remaining gel, used syringe and catheter after their use.

9. Storage and safe use conditions

Keep in a place protected from direct sunlight at a temperature of +1 °C to +30 °C and a humidity not greater than 65%.

Shelf life: 36 months from the date of manufacture.

Keep out of reach of children.

Do not use the device after the expiration date!

Use the device in accordance with the Instruction for Use.

Do not use the device with the damaged package and with traces of moisture inside the package! The product is sterile!

The product is disposable. Do not expose the product to re-sterilization.

10. Transportation conditions

The devices should be transported at +1 °C to +30 °C and a humidity not greater than 65% by any means of transport in covered vehicles while observing temperature conditions in accordance with the local laws and rules for the transportation of goods applied for each type of transport. Where applicable, use thermal vehicles for transportation of medical devices under climatic conditions specified by the manufacturer.

11. Components of the medical device

Antiadhesin® anti-adhesive, absorbable, sterile gel in a syringe, 1 pce.

Application catheter, 1 pce.

Instruction for Use, 1 pce.

Adhesive label, 1 pce.

12. Disposal procedure

In accordance with Sanitary and Epidemiological Rules and Norms SanPiN 2.1.3684-21, a used device (i.e., a device that has been in contact with the biological fluids of patients) is disposed of as Class B medical waste, and a damaged or expired device that has not been in contact with the biological fluids of patients is disposed of as household waste of Class A.

Manufacturer

Genewel Co., Ltd., Republic of Korea

6F, 37, Sagimakgol-ro 62 beon-gil, Jungwon-gu, Seongnam-si, Gyeonggi-do, Republic of Korea, Postal code: 13211.

Authorized representative of the manufacturer in the Russian Federation

Nizhpharm JSC, Russia

7 Salganskaya St., 603105 Nizhny Novgorod

Tel.: +7 (831) 278-80-88, fax: +7 (831) 430-72-28

E-mail: med@stada.ru

STADA

Warranty obligations


The manufacturer guarantees the quality of the device throughout the entire shelf-life. The device does not require maintenance and repair.

The quality of this medicinal product and manufacturing facilities of Genewel Co., Ltd., Korea meet all international quality standards, which was confirmed by certification under ISO 13485:2016.

Explanation of symbols on the package

	Manufacturer
	Date of manufacture
	Temperature range +1 °C to +30 °C
	Batch code

		Use before ...
		See the Instruction for Use
		With caution! See the Instruction for Use
		Do not use if the package is damaged
STERILE	EO	Sterilization with ethylene oxide
		Repeated use is prohibited
		Do not re-sterilize

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Lim, Youngill -----

attorney-in-fact of
Sang-Deok Han -----
CEO -----
GENEWEL CO., LTD. -----

appeared before me and admitted
said principal`s subscription to
the attached -----
Certificate of approval -----

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This is hereby attested on this
26th day of Jan. 2024 at this office.

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Name of the office
HAHNBAEK LAW & NOTARY OFFICE
Belong to Seoul Central
District Prosecutor's Office
Address of the office
138,BANPO-DAERO,SEOCHO-GU,SEOUL,
KOREA 4F (SEOCHO-DONG)

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Signature of the Notary Public
SEO JAE JONG

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.2.

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This public document

2. has been signed by SEO JAE JONG

3. acting in the capacity of Notary Public

4. bears the seal/stamp of HAHNBAEK LAW and NOTARY OFFICE

Certified

To verify the Apostille, please refer to the website below.
<https://www.apostille.go.kr>

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Authorized


Signature

**INSTRUCTION FOR USE
OF THE MEDICAL PRODUCT**

**STERILE ANTI-ADHESIVE ABSORBABLE GEL
STERILE ANTIADHESIN[®],
1,5 g, 3 g, 4 g, 5 g, 6 g, 10 g SYRINGES.**

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2. Purpose

Antiadhesin[®] anti-adhesive, absorbable, sterile gel,
1,5 g, 3 g, 4 g, 5 g, 6 g, 10 g syringes (hereinafter referred to as Antiadhesin[®] gel, medical device) is intended for the prevention of adhesions that may occur after any types of surgeries and can be manifested as the fusion of organs or tissues between each other. The therapeutic effect of the gel is related to the formation of a temporary artificial barrier between the injured tissues, ensuring effective separation of surfaces during healing. After application of Antiadhesin[®] gel onto the surgical site, it adheres to the anatomical surfaces without flowing, and forms a viscous, lubricant coating or anti-adhesive film, which ensures the sliding of adjacent injured surfaces, preventing them from sticking together. Antiadhesin[®] gel is a biodegradable coating that separates the mating surfaces only during the critical phase of wound healing and postoperative adhesion, lasting 7 days, without affecting the normal regeneration processes.

Antiadhesin[®] gel is indicated for use as an adjuvant in open and laparoscopic surgeries to reduce adhesions, reduce the spread and severity of postoperative adhesions in the surgical area, as well as the associated serious complications such as intestinal obstruction, pain and infertility, occurrence of which increases the risk of repeated surgical interventions.

Antiadhesin[®] gel is intended for the prevention of adhesions after surgeries on organs and tissues with serous coating: abdominal organs (stomach, intestines, spleen, liver); pelvic organs (uterus, ovaries, fallopian tubes), pleural cavity, pericardial cavity, testicular membrane, as well as during surgeries in the uterine cavity, on the bladder, mammary glands, in spine surgery. Antiadhesin[®] gel can be used for operations where there is a risk of soft tissue adhesion. For example, during neck surgeries, where there is a high risk of scar formation throughout the neck layers. Antiadhesin[®] gel may be used as an installation gel in plastic surgery. The gel facilitates installation of breast implants. Antiadhesin[®] gel can be used in ophthalmic operations.

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3. Main characteristics

Appearance	The medical device is a colorless clear gel in a pre-filled disposable syringe with a Luer lock, an application catheter.														
Composition	Composition of Antiadhesin[®] gel per 1 mL: <table><tr><td>Sodium hyaluronate</td><td>2.50 mg</td></tr><tr><td>Sodium carboxymethyl cellulose</td><td>5.00 mg</td></tr><tr><td>NaCl</td><td>8.00 mg</td></tr><tr><td>KCl</td><td>0.20 mg</td></tr><tr><td>KH₂PO₄</td><td>0.20 mg</td></tr><tr><td>Na₂HPO₄ · 12H₂O</td><td>1.15 mg</td></tr><tr><td>Water for injections</td><td>1 mL</td></tr></table> The gel is a mixture of two anionic polysaccharides: sodium hyaluronate (SH) and sodium carboxymethyl cellulose (SCMC). Sodium hyaluronate (SH) Sodium hyaluronate is a natural anionic polysaccharide composed of D-glucuronic acid and D-N-acetylglucosamine linked via alternating β-1,4- and β-1,3-glycosidic bonds. SH may be detected in various tissues of the body such as connective and cartilaginous tissue, as well as in the skin, vitreous	Sodium hyaluronate	2.50 mg	Sodium carboxymethyl cellulose	5.00 mg	NaCl	8.00 mg	KCl	0.20 mg	KH ₂ PO ₄	0.20 mg	Na ₂ HPO ₄ · 12H ₂ O	1.15 mg	Water for injections	1 mL
Sodium hyaluronate	2.50 mg														
Sodium carboxymethyl cellulose	5.00 mg														
NaCl	8.00 mg														
KCl	0.20 mg														
KH ₂ PO ₄	0.20 mg														
Na ₂ HPO ₄ · 12H ₂ O	1.15 mg														
Water for injections	1 mL														

humor and synovial fluid. SH is one of the main ingredients of the extracellular matrix; it is a high molecular viscoelastic substance with hydrophilic and non-immunogenic properties. Due to these features, SH can act as a lubricant when applied onto the injured surface. This was the basis for its use to reduce or prevent injury of the surgical area.


Sodium carboxymethyl cellulose (SCMC)

Sodium carboxymethyl cellulose is also an anionic polysaccharide. By moving the glycosidic hydroxyl group into the carboxymethyl group, SCMC becomes a more hydrophilic compound. Just as SH, SCMC has hygroscopic properties and may be used for pharmaceutical, cosmetic and food purposes as an adjuvant, lubricant, and stabilizing agent.

Technical and functional characteristics

Antiadhesin[®] gel contains highly purified sodium salt of hyaluronic acid obtained without the use of genetically modified microorganisms. It is a sterile, biocompatible, biodegradable, bioabsorbable, non-pyrogenic medical device. It does not contain antigens.

The medical device is biologically safe: non-pyrogenic, has no cytotoxic, genotoxic, irritant and sensitizing effects, acute and subacute toxicity.

Name of the device	Weight of content	Catheter tube		Syringe		
	Deviation: ±7%	Diameter (a, mm)	Length (b, mm)	Volume	Diameter (c, mm)	Length (d, mm)
Antiaghesin® gel	1,5 g	1.7 ± 0.085	80 ± 4,0	5 ml	15 ± 0,10	76,90 ± 0,2
	3 g		70 ± 3,5			
	4 g					
	5 g		70 ± 3,5	10 ml	18 ± 0,10	97,70 ± 0,2
	6 g					
	10 g					

Areas of application

Antiadhesin[®] gel is intended for the prevention of adhesions after all types of surgical interventions (either open or laparoscopic) on organs and tissues where there is a risk of postoperative adhesions: in general surgery, neurosurgery (spinal cord and brain membranes), otolaryngology (middle ear cavity, nasal cavity and paranasal sinuses), ophthalmology, thoracic surgery,

	cardiac surgery, abdominal surgery, orthopedic surgery, obstetrics and gynecology (including the uterine cavity), urology, plastic surgery.
Toxicity	Non-toxic
Possible application-related risks	Since safety and efficacy of the medical device use in combination with other anti-adhesive barriers, hemostatic agents, medicinal products; in patients with lymphatic and coagulation disorders; in patients with severe hepatic and renal impairment; in patients with malignant tumors; in patients with autoimmune or immune disorders, in pediatric patients, in pregnant and breastfeeding women, have not been established, it should be used with caution. It is recommended to avoid pregnancy until the end of the first menstrual cycle after using the device. When using Antiadhesin® gel, reactions to foreign bodies, as well as to any implanted material, may occur.
Sterilization	The gel is sterilized by membrane filtration. The packaged medical device is sterilized using ethylene oxide.
Compliance with standards	Applied standards EN 556-1:2001, EN 556-2:2015, EN ISO 20417:2021, EN ISO 10993-1:2018, EN ISO 10993-3:2014, EN ISO 10993-5:2009, EN ISO 10993-6:2016, EN ISO 10993-7:2008, ISO 10993-10:2021, EN ISO 10993-11:2018, EN ISO 11135-1: 2014, ISO 11135:2014, EN ISO 11607-1:2019, EN ISO 11607-2:2019, EN ISO 13485:2016, EN ISO 14155:2020, MEDDEV 2.7.1 version 4 (July 2016), EN ISO 13408-1:2015, EN ISO 13408-2:2018, ISO 14644-1:2015, ISO 14644-2:2015, ISO 14644-3:2019, EN ISO 14971:2019, ICH Q1A (R2), ASTM F88/F88M:21, ASTM F1980-21, ASTM F838-20, ISTA 2A: 2012.
Packaging	Antiadhesin® gel is packed into a pre-filled disposable 1,5 g, 3 g, 4 g, 5 g, 6 g, 10 g syringes; it is sterile and ready for use. The syringe and applicator catheter are sealed in a blister made of transparent polymer film, blister with instructions for use and an adhesive label are placed in a carton.

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5. Contraindications

Do not prescribe to patients with established hypersensitivity reactions to hyaluronic acid, carboxymethyl cellulose and salts thereof.

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
Device preparation for operation


- Open the package in sterile conditions **immediately before use**, remove the syringe and the application catheter.

- Remove the rubber cap from the Luer lock by turning the syringe body.
Warning: avoid contact of the syringe tip with foreign objects after removal of the rubber cap

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Procedure for device operation

- The target management area should be as dry as possible. It is necessary to remove excess liquid using an aspirator, as well as all hemostatic gauze sponges and swabs inserted after the operation.
- Ensure proper hemostasis: make sure that bleeding in the surgical site is completely stopped.
- Evenly cover all the anatomical surfaces, where adhesion should be prevented, with the gel.
- Dispose of the remaining gel, used syringe and catheter after their use.

8. Storage and safe use conditions, shelf life

Keep in a place protected from direct sunlight at a temperature of +1°C to +30°C and a humidity not greater than 65%.

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Use the device in accordance with the Instructions for use.

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9. Transport conditions

Products are transported at a temperature of +1°C to +30°C at a humidity of not more than 65% by all types of covered vehicles in compliance with the temperature regime in accordance with local legislation and the rules for the transportation of goods applied for each mode of transport. If necessary, use thermal machines for transporting medical devices under the climatic conditions specified by the manufacturer.

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Adhesive label, 1 pce.

11. Disposal procedure

In accordance with SanPiN 2.1.3684-21, a used product (i.e., one that has had contact with patients' biological fluids) is disposed of as class B medical waste, and a damaged or expired product that has not been in contact with patients' biological fluids is disposed of as household class A waste.

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
STADA

Warranty obligations

The manufacturer guarantees the quality of the device throughout the entire shelf-life. The device does not require maintenance or repair.

The quality of this medicinal product and manufacturing facilities of Genewel Co., Ltd., Korea, meet all international quality standards that has been certified by ISO 13485:2016.

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	Manufacturer
	Date of manufacture
	Temperature range +1°C to +30°C
	Lot number
	Use until...

	Refer to instructions for use	
	Carefully! Refer to instructions for use	
	Do not use if the package is damaged	
STERILE	EO	Ethylene oxide sterilization
	Reuse prohibition	
	Do not re-sterilize	

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Young Ok Lim -----

attorney-in-fact
Sang-Deok Han -----
CEO -----
GENEWEL CO., LTD. -----

appeared before me and admitted
said principal`s subscription to
the attached -----

INSTRUCTION FOR USE OF -----
THE MEDICAL PRODUCT -----

This is hereby attested on this
31st day of Jul, 2023 at this office.

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