To Whom It May Concern

**Registration documents for RelyX U200 and RelyX U200 Automix
- Summary of Technical Documentation (STED) -**

We, 3M Deutschland GmbH, Carl-Schurz-Str. 1, 41453 Neuss, herewith confirm that the following pages are an excerpt of the original Summary of Technical Documentation.

Seefeld, April 7, 2017

Yours Sincerely,

3M Deutschland GmbH
3M Oral Care

Elfriede Zedler
Regulatory Affairs

Summary Technical Document

**RelyX U200
RelyX U200 Automix**

3M Deutschland GmbH, Germany

CONFIDENTIAL!

The information contained in this dossier is 3M confidential. 3M subsidiaries are advised not to disclose information contained in this Registration Dossier to third parties other than local Health Authorities unless confidentiality agreements have been established with these third parties. The consolidated technical document given below is settled for the purposes of medical device state registration in the territory of the Russian Federation.

The details are provided to the extent required and in line with Medical Device state registration regulations (Resolution of the Government of Russian Federation № 1416 of 27 December 2012).

1 Table of Contents

1	TABLE OF CONTENTS	1
2	DEVICE DESCRIPTION.....	2
2.1	EXECUTIVE SUMMARY	2
2.2	PRODUCT NAME.....	2
2.3	PRODUCT DESCRIPTION AND MAIN CHARACTERISTICS.....	2
2.4	INDICATIONS/ CONTRAINDICATIONS.....	3
2.5	RESTRICTIONS ON USE	3
2.6	PRECAUTIONS.....	3
2.7	POTENTIAL COMPLICATIONS	4
2.8	PRINCIPLES OF OPERATION	4
2.9	INTERNATIONAL CODES AND COUNTRY REGULATORY CLASSIFICATIONS.....	7
2.10	COMMERCIAL PRESENTATION.....	8
2.11	CONSTRUCTION	8
2.12	PRODUCT PHOTO	12
2.13	PACKAGING DESCRIPTION.....	13
2.14	OPERATING ENVIRONMENT.....	14
2.15	STORAGE, USE AND SHELFIFE	14
3	FINISHED PRODUCT SPECIFICATIONS.....	14
4	INSTRUCTIONS FOR USE. GENERAL INFORMATION.	15
4.1	RELYX U200.....	15
4.2	RELYX U200 AUTOMIX.....	19
5	LABELS.....	25
6	TECHNICAL PROFILE.....	36
6.1	PRODUCT DESIGN REQUIREMENTS AND TEST RESULTS	38
7	MANUFACTURING INFORMATION.....	43
7.1	NAME AND ADDRESS OF THE MANUFACTURER.....	43
7.2	MANUFACTURING PROCESS	44
8	PRODUCT STANDARDS COMPLIANCE	45
9	BIOCOMPATIBILITY STATEMENT	45
10	STERILIZATION.....	46
11	PROTECTION OF THE ENVIRONMENT	46
12	MEDICINAL SUBSTANCES.....	46
13	HUMAN BLOOD DERIVATIVES	46
14	ANIMAL-DERIVED SUBSTANCES.....	46
15	UTILIZATION/LIQUIDATION OF MD	46
16	TRANSPORTATION INFORMATION	47
17	WARRANTY	47
18	RISK MANAGEMENT SUMMARY.....	47

2.4 Indications/ Contraindications

Indications

- Final cementing of all-ceramic, composite, or metal inlays, onlays, crowns and bridges; 2-3-unit Maryland bridges and 3-unit inlay/onlay bridges (contraindicated for patients with bruxism or periodontitis)
- Final cementing of posts and screws
- Final cementation of all-ceramic, composite, or metal restorations on implant abutments.
- Final cementation of Lava zirconia build-ups for two-piece abutments - only in accordance with the Instructions for Use for Lava Frame or Lava Plus respectively.

Contraindications

None

2.5 Restrictions on use

For use only by dental professionals.

2.6 Precautions

For Patients and Dental Personnel

- Base paste: contact with eyes may cause severe eye damage. Wear eye protection to prevent injury. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Base paste: contains sodium persulfate that may trigger an allergic respiratory reaction in certain individuals. This product may not be appropriate for use in those individuals with known sensitivity to sulfites since a cross-reaction may occur with sodium persulfate.

For Patients

- This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate and/or peroxide (especially perester) allergies.
- If prolonged contact with oral tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use.

For Dental Personnel

- This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water.
- The use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If the product contacts the glove, remove and discard the glove, wash the hands immediately with soap and water and then re-glove.
- If allergic reaction occurs, seek medical attention as needed. 3M ESPE MSDSs can be obtained from www.mmm.com or contact your local subsidiary. If allergic reaction occurs, seek medical attention as needed.

2 Device Description

2.1 Executive summary

Dental cements are used in restorative dentistry for the final cementation of dental restorations. become a permanent part of the organism and have to meet a large number of physical and mechanical as well as safety requirements. All dental materials are exposed to the constantly changing environment of the oral cavity - a supplementary stress factor for the materials used.

Cements are classified according to their chemical composition, which define the performance properties into Zinc phosphate cements, Carboxylate cements, Glass ionomer cements, Hybrid cements (non-modified glass ionomer cements and compomers), resin cements and self-adhesive resin cements. The class of the self-adhesive resin cements was created in 2002 by 3M ESPE, introduced RelyX Unicem.

RelyX U200 Automix represents the "next generation" of self-adhesive resin cement, based on RelyX Unicem technology, displaying additional features and new automix syringe delivery. It can be used for all type of restorative materials and does not need any pre-treatment of the tooth like etching, priming and bonding for standard indications. It is dual-curing and requires no absolute moisture exclusion during manipulation. RelyX Unicem technology. It can be used for all type of restorative materials and does not need any pre-treatment of the tooth like etching, priming and bonding for standard indications. It is dual-curing and requires no absolute moisture exclusion during manipulation. RelyX U200 is available as RelyX U200 Automix (represents the new automix syringe delivery system) and RelyX U200 Clicker (represents the clicker delivery system as hand mix system).

2.2 Product name

Authority Brand*: 3M ESPE

Brand/Product Name*: RelyX U200 Automix

Generic Descriptor: Self-Adhesive Resin Cement

Authority Brand*: 3M ESPE

Brand/Product Name*: RelyX U200

Generic Descriptor: Self-Adhesive Resin Cement

*3M, ESPE and RelyX are trademarks of 3M Company or 3M Deutschland GmbH. These trademarks may appear in conjunction with the product name(s) listed below on product packaging, inserts and other labeling.

2.3 Product description and main characteristics

RelyX U200 Automix a dual-curing, self-adhesive resin cement supplied in an automix syringe. It is used for the adhesive cementation of indirect all-ceramic, composite, or metal restorations and posts and screws. When RelyX U200 Automix is used, bonding and conditioning of the tooth structure are not necessary. Unlike other cementation materials, RelyX U200 Automix is characterized by its stability in combination with good flowability under pressure (structural viscosity).

RelyX U200 is a dual-curing, self-adhesive resin cement supplied in a clicker for manual mixing. It is used for the adhesive cementation of indirect all-ceramic, composite, or metal restorations and for posts and screws. When RelyX U200 is used, bonding and conditioning of the tooth structure are not necessary.

The cement is available in various shades. RelyX U200 Automix and Clicker contain a bi-functional (meth)acrylate. The proportion of inorganic fillers is about 43% by volume; the particle size (D 90%) is about 12.5 µm. The mixing ratio, based on volume, is 1 part base paste : 1 part catalyst.

Description of first aid measures

Inhalation:

Remove person to fresh air. If you feel unwell, get medical attention.

Skin Contact:

Immediately wash with soap and water. Remove contaminated clothing and wash before reuse. If signs/symptoms develop, get medical attention.

Eye Contact:

Flush with large amounts of water. Remove contact lenses if easy to do. Continue rinsing. If signs/symptoms persist, get medical attention.

If Swallowed:

Rinse mouth. If you feel unwell, get medical attention.

2.7 Potential complications

No known side effects are reported on.

2.8 Principles of operation

RelyX U200/ RelyX U200 Automix

Initiated by exposure of the mixed RelyX U200/RelyX U200 Automix to a dental curing light or by redox reaction, the methacrylate functions of Glypho and other Methacrylates undergo a radicalic polymerization. Basically, the chemistry of RelyX U200/RelyX U200 Automix comprises four different curing mechanisms:

Figure 1 Light induced polymerization


Figure 2 Chemical curing I




Figure 3 Chemical curing II

Chemical curing reaction I (see figure 2) provides for increased mechanical stability, chemical curing reaction II (see figure 3) provides for adhesive strength. The radicals, regardless by which mechanism they were developed, initiate a radicalic polymerization of the methacrylate functions of the monomers (see figure 4).

Figure 4 Radicalic Polymerisation

In a first step the acidic functions deprotonize (see figure 5). Calcium and aluminum from the reactive fillers react with protons by ion interchange and are dissolved (see figure 6).

Figure 5 Acid-base reaction




Figure 6 Release of Ions

The dissolved ions react with the phosphate functions resulting in chelate complexes. Due to reactions, the acid groups are neutralized (see figure 7) and a long term stability of the material is obtained.

from the
)

Figure 7

axes. Due to
y of the mat

2.9 International codes and country regulatory classifications

RelyX U200/ RelyX U200 Automix

Global Nomenclature	
GMDN	35870 Resin-composite dental cement
UMDNS	16-707 Cement, Dental, Resin Composite
Country	Regulatory Classification
US	Class II, 21 CFR § 872.3275(b) Dental cement other than zinc oxide-eugenol
EU	Class IIa per Rule 8, Annex IX, of the Medical Device Directive 93/42/EEC
Canada	Class 3
Duration of contact in accordance to ISO 10993:	Intended to be in contact with the body for greater than 30 days
The rate of application	Single
Sterilization	Non-sterile

2.10 Commercial presentation**RelyX U200**

Item	Catalogue/Order Number
1 Clicker (11 g), colour shade: Translucent Accessories (1 mixing pad)	56877
1 Clicker (11 g), colour shade: A2 Universal Accessories (1 mixing pad)	56878
1 Clicker (11 g), colour shade: A3 Opaque Accessories (1 mixing pad)	56879

RelyX U200 Automix

Item	Catalogue/Order Number
1 x 8.5 g syringe, colour shade: Translucent 10 x RelyX U200 Automix Tips, Mixing Tips Regular 5 x RelyX U200 Automix Tips, Mixing Tips Wide 5 x RelyX U200 Automix Endo Tips	56897
1 x 8.5 g syringe, colour shade: A2 Universal 10 x RelyX U200 Automix Tips, Mixing Tips Regular 5 x RelyX U200 Automix Tips, Mixing Tips Wide 5 x RelyX U200 Automix Endo Tips	56895
1 x 8.5 g syringe, colour shade: A3 Opaque 10 x RelyX U200 Automix Tips, Mixing Tips Regular 5 x RelyX U200 Automix Tips, Mixing Tips Wide 5 x RelyX U200 Automix Endo Tips	56896
ACCESSORIES	
30 RelyX U200 Automix Tips, Mixing Tips Regular	56913
15 RelyX U200 Automix Tips, Mixing Tips Wide 15 RelyX U200 Automix Endo Tips	56914
15 RelyX U200 Automix Tips, Mixing Tips Wide 15 RelyX U200 Automix Intraoral Tips	56917

2.11 Construction**Composition of RelyX U200/ RelyX U200 Automix**

Base			
Function	Ingredient	CAS-Number	Quantity
Filler	Silane Treated Glass Filler	none	51 – 52
Monomer	Dimethacrylate	109-16-0	12.6 – 15
Monomer	Phosphoric acid methacrylates	1224866-76-5	22 – 23
Filler	Silane treated fumed silica	68909-20-6	6.10 – 7.5
Filler	Glass filler	65997-17-3	1.74 – 2.3
Initiator	Sodium Persulfate	7775-27-1	1.68 – 2.1
Initiator	Peroxide	13122-18-4	0.17 – 0.2
Stabilizer	Stabilizer	128-37-0	< 0.1
Initiator	Copper acetate	6046-93-1	< 0.1

logue/Order
Number

Catalyst			
Function	Ingredient	CAS-Number	Quantity %*
Filler	Silane treated glass filler	none	58 – 65
Monomer	Dimethacrylates	27689-12-9	23 – 26
		72829-09-5	3.24 – 4.40
Initiator	Salt of substituted barbituric acid	945012-02-2	2.43 – 3.30
Filler	Silane treated fumed silic	68909-20-6	1.94 – 2.93
Initiator	Calcium hydroxide	1305-62-0	1.15 – 1.75
Initiator	Sodium ptoluenesulfinate	824-79-3	1.21 – 1.73
Initiator	Methacrylated amine	93962-71-1	0.68 – 0.90
		93962-70-0	0.09 – 0.20
Initiator	Camphorquinone/ dl-Camphorquinone	10373-78-1	0.17 – 0.23
Stabilizer	Stablizer	128-37-0	< 0.1
		150-76-5	< 0.1
Pigment	Pigments	82-38-2	≤ 0.1
		13463-67-1	
		100280-62-6	

*The composition is presented as a target composition or 'mean composition'. Process optimization and variation in raw materials content of particular components may impact the actual released product formulation.

56896

Pigments in RelyX U200 Automix & RelyX U200	CAS-Number
Pigment White 6, CI 77891	13463-67-7
Quinoline Yellow Lake	100208-62-6
Solvent Red 111, CI 60505	82-38-2

56913

56914

56917

Mixing tips (regular)/ Mixing tips (wide)

Function	Ingredient	CAS-Number
Mixing tube/Cover plate	Polypropylene	9003-07-0
Rotor	Polyoxymethylene	9002-81-7, 24969-26-4

Quantity %

51 – 52

12.6 – 15.4

22 – 23

6.10 – 7.52

1.74 – 2.36

1.68 – 2.01

0.17 – 0.23

< 0.1

< 0.1

Intra-oral tips: Polypropylene

CAS 9003-07-0


Endo tips: Polypropylene

CAS 9003-07-0


Mix pad: Paper

CAS 9004-34-6


Drawings


Automix Endo Tips


Automix Intraoral Tips

Mixing Tips Regular


Mixing Tips Wide


2.12 Product photo

RelyX U200 Automix Syringe

RelyX U200 Clicker

RelyX U200

2.13 Packaging description

Materials

RelyX U200

Item	Part	Material	CAS	Dimensions, mm
Outer package*		Cardborad box	9004-34-56	225,0 x 142,0 x 45,0
Primary package*	Clicker (complete)	Polypropylene and Polyethylene	9003-07-0 9002-88-4	l=190 w=22,1 h=34/ Ø 9,7
	Catridge 	Polyethylene (black)	9002-88-4	l=140, w=22,1 h/Ø 9,7
	Plunger 	Glass fiber reinforced polypropylene (yellow)	9003-07-0	l=115 w=19,4 h=5
	Cap 	Polypropylene (yellow)	9003-07-0	l=27,2 w=24 h=11,5 h with nose: 34
Secondary package		Polyethylene laminated aluminum foil pouch with printed label on outer surface	9002-88-4	l: 270 ±3 x w: 120 ± 2

RelyX U200 Automix

Item	Part	Material	CAS	Dimensions, mm
Outer package*		Cardborad box	9004-34-56	225,0 x 142,0 x 45,0
Primary package*	Syringe (complete)	Polypropylene and Polyethylene	9003-07-0 9002-88-4	l=160 w=26 h=20.6/Ø=9 with cap=171.8
	Catridge 	Polyethylene (black)	9002-88-4	l=93,3 w=26 h=20.6/Ø=9
	Plunger 	Glass fiber reinforced polypropylene (yellow)	9003-07-0	l=86,4 w=7.2 h=20.6/Ø=7,1
	Cap 	Polypropylene (yellow)	9003-07-0	Ø 11,8 x h 14,5
Secondary package		Polyethylene laminated aluminum foil pouch with printed label on outer surface	9002-88-4	l=270 ±3 w=120 ± 2

* ranges not defined

Accessories

Item	Material	CAS	Dimensions, mm
Primary package	Polyethylene	9002-88-4	L: 150 $\pm$ 2 x W: 175 $\pm$ 2

2.14 Operating Environment

The products are designed to be used at room temperature. If stored in cooler allow product reach room temperature prior to use. Ambient temperatures routinely higher than 27°C/80°F reduce shelf life.

2.15 Storage, Use and Shelflife

RelyX U200/ RelyX U200 Automix

The shelf-life of the product is 18 months when stored in the pouch at 15°-25 °C/59°- 77 °F. Do not use the product after expiration date.

After removal from the pouch, use the product within 6 months and before the lapse of the expiration date. Constantly high humidity accelerates the setting and must be avoided.

Store away from heat.

3 Finished Product Specifications

RelyX U200/ RelyX U200 Automix

Methods	Limit	Reference
Setting time (start)	> 0.5 min	Inhouse -method
Appearance base, Color, Odor	Solid, tooth colored paste with slight acrylic odor	Organoleptic test
Appearance catalyst, Color, Odor	Solid, tooth colored paste with slight acrylic odor	

4 Instructions for Use. General information.

4.1 RelyX U200

Tooth Preparation

Preparation of the Cavity/Tooth Stump

- Prior to final cementation, clean the prepared stump or the cavity thoroughly with pumice slurry, rinse with a water spray, and lightly air dry in only 2-3 bursts of water-free and oil-free air, or use cotton pellets to dry it off. Do not overdry!
 - The cavity should be just dry enough that the surface has a slightly glossy appearance. As is the case with any permanent cement, over-drying can lead to post-operative sensitivity.
- Do not use substances such as desensitizers, disinfectants, astringents, dentin sealants, rinsing solutions containing EDTA, etc., after the final cleaning with pumice slurry and water. Their residues may have a detrimental effect on the bonding strength and setting reaction of the cement.

Preparation of Root Canals

- Treat the root canal endodontically as usual (root canal filling with guttapercha and removal of the root canal filling, leaving 4 mm of gutta-percha apically).
- Clean the root canal with a 2.5%-5.25% solution of sodium hypochlorite (NaOCl).
- Rinse immediately with water and dry with paper points.
- We recommend the use of a rubber dam during the cementation of posts.

Preparation for Maryland and Inlay/Onlay Bridges

Restored teeth must have an adequate enamel surface for bonding. They should be healthy or only slightly restored and the periodontal conditions should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- Prepare retentive elements such as cingular rests and/or approximal grooves.
- Use a rubber dam and keep the cavity free of any contamination during cementation.
- Etch the enamel surface of the cavity with 37% phosphoric acid for 15-20 sec. Then rinse thoroughly with water and dry with water-free and oil-free air. In the case of exposed dentin, make sure to selectively etch the enamel to avoid post-operative sensitivity.

Trial Fitting of the Restoration

- Remove the temporary restoration and thoroughly clean the tooth surface of any residue from the temporary cement.
- Make a trial insertion of the final restoration and check its fit and contact points.
 - If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely.
 - For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.
- Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment and until the final cementation.

Pre-treatment of Restorations**Pre-treatment of Metal Restorations**

Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- Blast the restoration surface to be luted with aluminum oxide <40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".

Pre-treatment of Etchable Glass Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.
- Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use, e.g.:
 - RelyX Ceramic Primer: allow to react for 5 sec, then air dry so that the solvent evaporates completely.

Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

Alternative 1:

- Blast the restoration surface to be luted with aluminum oxide <40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Alternative 2:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet Prep the blast-coating agent CoJet Sand from a distance of 2-10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".

Pre-treatment of Composite Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Blast the composite restoration surface to be luted with aluminum oxide <40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Pre-treatment of Glass Fiber-reinforced Posts

Please follow the instructions for use of the post to be used. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Clean the post with alcohol and dry it with water-free and oil-free air.
- When using RelyX Fiber Post application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with the relevant instructions for use.

Pre-treatment of Implant Abutments

Please follow the recommendations for the relevant restorative material (e.g., metal, zirconia).

Pre-treatment of Maryland and Inlay/Onlay Bridges Metal, Zirconia and Aluminium Oxide

Surface:

- **Coat (silicate)** the restoration surface to be luted with the micro-blasting device CoJet Prep and the blast-coating agent CoJet Sand from a distance of 2-10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand.
- **Blow away** any residues of the blasting agent with water-free and oil-free air.
- **Apply** an appropriate silane in accordance with the instructions for use.

For etchable glass ceramics, please refer to "Pre-treatment of Etchable Glass Ceramic Restorations".

Dosing

When dosing the first time:

- **Remove** the clicker dispenser from the foil pouch and note the removal date on the clicker label. Discard the pouch.
- **Remove** the cap from the clicker dispenser.
- **Dispense** a small quantity of paste on a mixing pad by pressing down the clicker lever. This levels the filling quantity of base paste and catalyst paste in the clicker. Discard the paste which has been dispensed.

When dosing, always:

- **Remove** the cap.
- **Dispense** the paste on a mixing pad:
 - Press the clicker lever and hold it down.
 - When the paste stops flowing, wipe off the clicker openings on the mixing pad. Doing so dispenses the correct volume of the base paste and catalyst paste.
 - Dispense only as much paste on the mixing pad (number of clicks) as required for cementing.
- **Cleaning** the clicker openings:
 - Wipe the clicker openings with a clean cloth to prevent blockage from cross-contamination of the pastes.
- **Replace** the cap by pressing until it locks into place noticeably.

Times

The processing and setting times depend on the ambient and oral temperature. The times shown are based on conditions relevant for practice. As is the case with every composite cement, the setting of RelyX U200 cement slows down significantly at room temperature.

RelyX U200 is a dual-curing cement and therefore also sensitive to natural or artificial light. The working time is significantly reduced during application under operating lights!

	min:sec
Working time from start of mixing	02:30
Light-curing in a lightwave range of 400-500 nm	
- Single surface, from occlusal	00:20
- Any other surface additional	00:20
- RelyX Fiber Post posts, from occlusal	00:40
Self-curing	
- Setting time after start of mixing	06:00

Mixing

- Mix base paste and catalyst paste into a homogenous paste within 20 seconds using a spatula. Avoid incorporating air bubbles.
- Then cover the cement with a light protection cover to protect it from intensive light.

Application

- Protect the working area from contamination with water, blood, saliva and sulcal fluid during the application and setting phase.
- Apply the mixed cement evenly to the entire cavity and as appropriate to the bottom side of the inlay/onlay, or fill the crown with cement.
- Seat the restoration firmly and stabilize long enough for the cement to fully set.

Application in the Root Canal

- Do not use Lentulo-Spirals to insert the cement in the root canal as this can excessively accelerate setting.
- Spread cement to the post.
 - As an alternative, the cement can be inserted into the root canal with the aid of an applicator tip, e.g., Centrix® Needle Tubes (Centrix is a registered trademark of the U.S.-based company Centrix); this will minimize the entrapment of air bubbles.
- Immediately place the post in the pretreated root canal; apply moderate pressure to hold it in position. We recommend rotating the post slightly during insertion to avoid the inclusion of air bubbles.

Removal of Excess

Notes on oxygen inhibition: As is the case with all composite materials, RelyX U200 is subject to oxygen inhibition, i.e., the upper layer (about 50 µm) which is in contact with atmospheric oxygen during the polymerization does not harden. Leaving sufficient excess permits the removal of this uncured layer during shaping/polishing without leaving behind a deficit.

Using light-curing:

- After brief curing (about 2 sec using a standard polymerization device), use a suitable instrument (e.g., scaler) to remove the excess. Larger volume excess can be removed more easily.
- Immediately light-cure after removal (see Polymerization and Shaping).

Using self-curing:

- Use a suitable instrument (e.g., scaler) to remove the excess during the self-curing phase (about 3 min after beginning the mixing in the "gel phase").
- Immediately cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.

Alternatively:

- Remove the excess before polymerization, e.g., using a sponge pellet, while using a suitable instrument to hold the restoration in position.
- Cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.
- Immediately light-cure or wait for the self-curing to finish (see Polymerization and Shaping).

Polymerization and Shaping

- We recommend light-curing the cement through the restoration when doing ceramic and composite work. Select the exposure times appropriate for the number of surfaces (please refer to "Times"). Polymerize the cement through the post when using translucent posts. The exposure time depends on the translucency of the post being used; for RelyX Fiber Post, it is 40 sec.
- Shape any remaining uneven areas on the edges of restorations and polish the marginal areas with diamond devices, aluminum oxide coated discs (e.g., Sof-Lex™), and diamond polishing paste.
- Then check the occlusion.

After Shaping

- Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing on implant abutments.

Hygiene and Disinfection

- Clean a contaminated RelyX U200 with cleaning agents commonly used in the dental practice.
- Use a cloth soaked in a disinfectant commonly used in the dental practice to disinfect the clicker.

Notes

- RelyX U200 can prematurely polymerize when exposed to natural or artificial lighting. Therefore dispense the paste on a mixing pad only just before mixing and apply immediately after mixing; avoid intensive light exposure during application.
- Deviations from the specified mixing ratio (1:1) adversely affect the material properties.

4.2 RelyX U200 Automix**Tooth Preparation****Preparation of the Cavity/Tooth Stump**

- Prior to final cementation, clean the prepared stump or the cavity thoroughly with pumice slurry, rinse with a water spray, and lightly air dry in only 2-3 bursts of water-free and oil-free air, or use cotton pellets to dry it off. Do not overdry!
 - The cavity should be just dry enough that the surface has a slightly glossy appearance. As is the case with any permanent cement, overdrying can lead to post-operative sensitivity.

Do not use substances such as desensitizers, disinfectants, astringents, dentin sealants, rinsing solutions including EDTA, etc., after the final cleaning with pumice slurry and water. Their residues may have a detrimental effect on the bonding strength and setting reaction of the cement.

Preparation of Root Canals

- Treat the root canal endodontically as usual (root canal filling with guttapercha and removal of the root canal filling, leaving 4 mm of guttapercha apically).
- Clean the root canal with a 2.5%-5.25% solution of sodium hypochlorite (NaOCl).
- Rinse immediately with water and dry with paper points.
- We recommend the use of a rubber dam during the cementation of posts.

Preparation for Maryland and Inlay/Onlay Bridges

Abutment teeth must have an adequate enamel surface for bonding. They should be healthy or slightly restored and the periodontal conditions should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- Prepare retentive elements such as cingular rests and/or approximal grooves.
- Use a rubber dam and keep the cavity free of any contamination during cementation.
- Etch the enamel surface of the cavity with 37% phosphoric acid for 15-20 sec. Then rinse thoroughly with water and dry with water-free and oil-free air. In the case of exposed dentin, make sure to selectively etch the enamel to avoid post-operative sensitivity.

Preparatory Measures

- Remove the temporary restoration and thoroughly clean the tooth surface of any residue of the temporary cement.
- Make a trial insertion of the final restoration and check its fit and contact points.
 - If a low viscosity silicone is used for the trial fit, the silicone residue must then be removed completely.
 - For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.
- Avoid any contamination whatsoever of the surfaces to be treated during pre-treatment until the final cementation.

Pre-treatment of Restorations**Pre-treatment of Metal Restorations**

Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- Blast the restoration surface to be luted with aluminum oxide <40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay bridges".

Pre-treatment of Etchable Glass Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Use hydrofluoric acid to etch the glass ceramic restoration surface to be luted.
- Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use, e.g.:
 - RelyX[®] Ceramic Primer: allow to react for 5 sec, then air dry so that the solvent evaporates completely.

Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

Alternative 1:

- Blast the restoration surface to be luted with aluminum oxide <40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Alternative 2:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet® Prep and the blast-coating agent CoJet Sand from a distance of 2-10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Composite Restorations

Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Blast the composite restoration surface to be luted with aluminum oxide <40 µm.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Pre-treatment of Glass Fiber-reinforced Posts

Please follow the instructions for use of the post to be used. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Clean the post with alcohol and dry it with water-free and oil-free air.
- When using RelyX® Fiber Post application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with the relevant instructions for use.

Pre-treatment of Implant Abutments

Please follow the recommendations for the relevant restorative material (e.g., metal, zirconia).

Pre-treatment of Maryland and Inlay/Onlay Bridges

Metal, Zirconia and Aluminium Oxide Surfaces:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet Prep and the blast-coating agent CoJet Sand from a distance of 2-10 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prep and CoJet Sand.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For etchable glass ceramics, please refer to "Pre-treatment of Etchable Glass Ceramic Restorations".

Times

The processing and setting times depend on the ambient and oral temperature. The times shown are based on conditions relevant for practice. As is the case with every composite cement, the setting of RelyX U200 Automix cement slows down significantly at room temperature.

RelyX U200 Automix is a dual-curing cement and therefore also sensitive to natural or artificial light. The working time is significantly reduced during application under operating lights!

	min:sec
Working time from start of mixing	02:30
Light-curing in a lightwave range of 400-500 nm	
- Single surface, from occlusal	00:20
- Any other surface additional	00:20
- RelyX Fiber Post posts, from occlusal	00:40
Self-curing	
- Setting time after start of mixing	06:00

Application

➤ Select a RelyX U200 Automix syringe containing the desired color and place it ready for use.
When using a new RelyX U200 Automix syringe:

- Remove the automix syringe from the foil bag and discard the bag.
- Note down the date the syringe was removed on the syringe label.
- Remove and discard the sealing cap from the automix syringe. Check the syringe opening for blockage; remove any paste plugs.
 - **The sealing cap must not be reused to close the syringe** because doing so might introduce air bubbles into the material.
- Squeeze out a small quantity of paste onto a block to equal the base paste and catalyst pastes from the automix syringe. Discard the paste which has been squeezed out.

During every application:

Protect the working area from contamination with water, blood, saliva and sulcal fluid during application and setting phase.

- Remove and discard the mixing tip remaining on the automix syringe from the previous application.
- Attach one of the two types of mixing tips ("Standard" or "Wide") and secure it by turning the cap to the right.
 - The mixing tip "Standard" is used without an Intraoral Tip for application in easily accessible areas.
 - The mixing tip "Wide" is used for the application with the Intraoral Tip or the EndoTip.
- When using the mixing tip "Wide", attach an Intraoral Tip or an EndoTip.
- Squeeze out and discard a peppercorn-size quantity of RelyX U200 Automix until an evenly mixed paste in a homogeneous color flows out of the tip.
 - The paste requires a certain amount of time to flow through the mixing tip. The flow rate speed cannot be accelerated by increasing the pressure on the plunger.
 - As soon as the pressure on the plunger decreases, the material flow stops and the material begins to set up. Do not use force to press out paste that has set as this could cause damage to the mixing tip and the automix syringe.

- Apply RelyX U 200 Automix evenly to the entire cavity and as appropriate to the bottom side of the inlay/onlay or fill the crown with cement.
 - Keep the opening of the mixing tip "Standard" or of the Intraoral or Endo Tip immersed in the material during the entire application to prevent the inclusion of any air bubbles.
- Seat the restoration firmly and stabilize long enough for the cement to set fully.
- Leave the used mixing tip on the automix syringe as a cap until the next application.

Application in the Root Canal

- See instructions under "Application" for handling of the automix syringe.
- Do not use Lentulo-Spirals to insert the cement in the root canal as this can excessively accelerate setting.
- Attach an Endo Tip to the mixing tip "Wide" for application in the root canal
- Insert the EndoTip as deeply as possible in the root canal and apply RelyX U200 Automix, beginning apically. Keep the tip of the EndoTip immersed in the cement and slowly move the Endo Tip upwards as the level of the paste rises.
- Do not remove the Endo Tip from the cement until the root canal has been completely filled.
 - This so-called immersion filling should not be carried out in less than 5 sec.; this will minimize the entrapment of air bubbles.
- Place the post in the root canal filled with cement; apply moderate pressure to hold it in position. We recommend rotating the post slightly during insertion to avoid the inclusion of air bubbles.

Removal of Excess

Notes on oxygen inhibition: As is the case with all composite materials, RelyX U200 Automix is subject to oxygen inhibition, i.e., the upper layer (about 50 µm) which is in contact with atmospheric oxygen during the polymerization does not harden. Leaving sufficient excess permits the removal of the uncured layer during shaping/polishing without leaving behind a deficit.

From Restorations

Using light-curing:

- After brief curing (about 2 sec using a standard polymerization device), use a suitable instrument (e.g., scaler) to remove the excess. Larger volume excess can be removed more easily.
- Immediately light-cure after removal (see Polymerization and Shaping).

Using self-curing:

- Use a suitable instrument (e.g., scaler) to remove the excess during the self-curing phase (about 3 min after beginning the mixing in the "gel phase").
- Immediately cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.

Alternatively:

- Remove the excess before polymerization, e.g., using a sponge pellet, while using a suitable instrument to hold the restoration in position.
- Cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.
- Immediately light-cure or wait for the self-curing to finish (see Polymerization and Shaping).

From Posts

- Remove the cement with a suitable instrument or a cotton pellet.

Polymerization and Shaping

- We recommend light-curing the cement through the restoration when doing ceramic composite work. Select the exposure times appropriate for the number of surfaces (please refer to "Times"). Polymerize the cement through the post when using translucent posts; exposure time depends on the translucency of the post being used; for RelyX Fiber Post, it is 10 sec.
- Shape any remaining uneven areas on the edges of restorations and polish the margins with diamond devices, aluminum oxide coated discs (e.g., Sot-Lex™), and diamond polishing paste.
- Then check the occlusion.

After Shaping

- Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing on implant abutments.

Hygiene and Disinfection


- We recommend the use of commonly available hygienic protective covering to avoid contamination of the RelyX U200 Automix syringes during the treatment.
- Clean contaminated RelyX U200 Automix syringes with cleaning agents commonly used in dental practice.
- Use a cloth soaked in a disinfectant commonly used in the dental practice to disinfect the automix syringe.

Notes


- RelyX U200 Automix can prematurely polymerize when exposed to natural or artificial lighting. Therefore avoid intensive light exposure during application.

5 Labels

RelyX U200

Sample of inner label

3M ESPE

RelyX™ U200

- Self-Adhesive Resin Cement
- Cemento definitivo de resina autoadhesivo
- Cimento resinoso autoadesivo
- 자가접착형 컴포지트 레진시멘트
- 自黏性樹脂黏合劑
- செய்நிறைந்த பிளாஸ்டிக்
- Ciment de scellement composite auto-adhésif
- Цемент самоклеящийся универсальный композитный для фиксации RelyX U200 в наборе (в тозирующей системе Clicker)
- Самодететивен композитен цемент
- Kompozitni cement za adhezivno cementiranje protetskih nadomjestaka
- Önkondicionáló rezincement
- Samoprzylegający cement do osadzania prac protetycznych
- Ciment rășinic autoadeziv
- Samoadhezivný žltkový cement
- Samu-adhezivni kompozitni cement
- Samoadhezivni kompozitni pryskyfičný cement
- Self-Adesiv Resin Siman
- Isekinnituv komposittsement
- Palimėjods svečia cementis
- Savaiminio sukibimo dervinis cementas fiksavimui
- Самодететивний композитний цемент

CONTENTS

11 g Base Paste/
Catalyst Paste

1

Accessories

Contains: Disodium peroxodisulphate.
May produce an allergic reaction.

Made in Germany by
3M Deutschland GmbH
Dental Products
Carl-Schurz-Str. 1
41453 Neuss - Germany
(049 (0) 8152/700-0)

3M, ESPE and RelyX are trademarks of 3M or 3M Deutschland GmbH.
© 2016, 3M. All rights reserved.

Larger images of box label

RelyX U200 Automix

3M ESPE


44000188274/01

**RelyX™ U200
Automix****Self-Adhesive
Resin Cement****8,5 g Base Paste/Catalyst Paste****3M Deutschland GmbH
41453 Neuss - Germany**TR
TR
TR
TR
TR
TR
TR
TR
TR
TR
TR

Open Date:

Sample of inner label

Sample of box

3M ESPE

RelyX™ U200

Automix

- | | |
|---|--|
| (CS) Self-Adhesive Resin Cement | (HU) Önkondicionáló rezincement |
| (ES) Cemento definitivo de resina autoadhesivo | (PL) Samoprzylegający cement do osadzania prac protetycznych |
| (PT) Cimento resinoso autoadesivo | (RO) Ciment rășinic autoadeziv |
| (KO) 자가 접착 레진시멘트 | (SK) Samoadhezivný živcový cement |
| (ZH) 自黏性樹脂黏合劑 | (SI) Samo-adhezivni kompozitni cement |
| (TH) เรซินซีเมนต์ที่เป็นสารยึดติดในตัว | (CS) Samoadhezivní kompozitní pryskyřičný cement |
| (FR) Ciment de scellement composite auto-adhésif | (IT) Self-Adeziv Resin Siman |
| (RU) Цемент самоадгезивный универсальный композитный для фиксации RelyX U200 в наборе (в дозаторе для автоматического смешивания) | (LT) Isekinnituv komposittsement |
| (PT) Самоадхезивен композитен цемент | (LV) Pašīmējošs sveķu cements |
| (HR) Kompozitni cement za adhezivno cementiranje protetskih nadomjestaka | (NL) Savaiminio sukibimo dervinis cementas fiksavimui |
| | (UK) Самоадгезивний композитний цемент |

CONTENTS

8.5g Base Paste/
Catalyst Paste

1

RelyX™ U200
Automix Tips

10

5

RelyX™ U200
Automix Endo Tips

5

Contains/Enthält: Disodium peroxodisulphate/Dinatriumperoxodisulfat.
May produce an allergic reaction/Kann allergische Reaktionen hervorrufen.

Made in Germany by
3M Deutschland GmbH
Dental Products
Carl-Schurz-Str. 1
41453 Neuss - Germany
0049 (0) 8152/700-0

3M, ESPE and RelyX are trademarks of 3M or 3M Deutschland GmbH.
Used under license in Canada. © 2016, 3M. All rights reserved.

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3M Deutschland GmbH
Dental Products
Carl-Schurz-Str. 1
41453 Neuss - Germany

40001821506

3M, ESPE and RelyX
are trademarks of 3M or
3M Deutschland GmbH.
© 2016, 3M
All rights reserved.


Sample of Mixing Tips-Regular label

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
Endo Tips
15 

REF 56914

3M ESPE and RelyX
are trademarks of 3M or
3M Deutschland GmbH
© 2016 3M
All rights reserved

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
Intraoral Tips

15  **REF 56917**

3M, ESPE and RelyX are trademarks of 3M or 3M Deutschland GmbH
 © 2016, 3M.
 All rights reserved.

Mixing
Tips WideIntraoral
Tips

4400189231/06

0123**3M ESPE I.D. No. 70201139386****LOT 000000****56917****REF**

4400176201139386/50000000

Sample of Mixing Tips – Intraoral label

3M ESPE**RelyX™ Unicem 2 Automix****RelyX™ U200 Automix**

Mixing Tips Wide

5 

Endo Tips

5 Mixing
Tips WideEndo
Tips**3M Deutschland GmbH****41453 Neuss - Germany****LOT 000000**

44000 189199 02

Sample of inner label for accessories

Explanation of symbols on labels and packaging:

Symbol	Description
	Manufacturer
	Do not re-use
	Use-by Date
	Catalogue number
	Lot number
	Product meets the requirements of the applicable EC directives
	Non-corrugated fiberboard (paperboard)
	Caution
	Low-density polyethylene
	Green Dot
	Contents
	Limit of temperature

6 Technical profile

RelyX U200

Clicker (11 g)

Dimensions, mm: length: 190 ± 5 mm, width: 26 ± 3 mm, height: 10 ± 2 mm.

Extrusion force, N: 4.5

Flexural strength, MPa: 99

Compressive strength, MPa: 291

Modulus of elasticity, GPa: 6.6

Surface hardness, MPa: 190

Film thickness, μm : 13

Water sorption, $\mu\text{g}/\text{mm}^3$: 28

Solubility, $\mu\text{g}/\text{mm}^3$: 0

Expansion after 1 month, %: 0.63

Base paste

General Physical Form: Solid

Specific Physical Form: Paste

Odor, Color, Grade: tooth-colored paste with slight acrylic odor

Density, g/cm^3 : 2 - 2.2

Specific Gravity, [Ref Std: WATER=1]: 2 - 2.2

Solubility in Water: Negligible

Catalyst paste

General Physical Form: Solid


Specific Physical Form: Paste

Odor, Color, Grade: tooth-colored paste with slight acrylic odor

Density, g/cm^3 : 2 - 2.2

Specific Gravity, [Ref Std: WATER=1]: 2 - 2.2

Solubility in Water: Nil


Mixing pad

Dimensions, mm: length: 94 ± 1 , width: 70 ± 1 , height: 9 ± 0.3

Number of pages: 45

RelyX U200 Automix

8.5 g syringe

Dimensions, mm: length: 170 ± 5 mm, width: 18 ± 2 mm, height: 9 ± 2 .

Extrusion force, N : 40.9

Flexural strength, MPa: 99

Compressive strength, MPa: 291

Modulus of elasticity, GPa: 6.6

Surface hardness, MPa: 190

Film thickness, μm : 13

Water sorption, $\mu\text{g}/\text{mm}^3$: 28

Solubility, $\mu\text{g}/\text{mm}^3$: 0

Expansion after 1 month, %: 0.63

Base paste

General Physical Form: Solid

Specific Physical Form: Paste

Odor, Color, Grade: tooth-colored paste with slight acrylic odor

Density, g/cm^3 : 2 - 2.2

Specific Gravity, [Ref Std: WATER=1]: 2 - 2.2

Solubility in Water: Negligible

Catalyst paste

General Physical Form: Solid

Specific Physical Form: Paste


Odor, Color, Grade: tooth-colored paste with slight acrylic odor

Density, g/cm³: 2 - 2.2

Specific Gravity, [Ref Std: WATER=1]: 2 - 2.2


Solubility in Water: Nil

Mixing Tips Regular


Dimension, mm: length: 39.7 ±1, diameter: 11.8 ±0.1

Mixing Tips Wide


Dimension, mm: length: 33.4 ±1, diameter: 11.8 ±0.1

Automix Endo Tips

Dimension, mm: length: 15.9±0.2mm, diameter: 5±0.2mm

Automix Intraoral Tips

Dimension mm: length: 14.2±0.2 diameter: 5±0.2

6.1 Product design requirements and test results**RelyX U200**

Product Name: RelyX Unicem 2 Clicker / RelyX U200, Shade: A2

Batch/Lot number: 626921

Test method: ISO 4049:2009, Classification: Type 2, Class 3 acc. to ISO 4049 (dual cure materials)

Test Conditions: Temperature: 23°C (± 1 °C); Humidity: 50 % (± 10 %)

Test Results

Test: Test methods: Limit/Unit:	Film Thickness ISO 4049 < 50µm	Working Time ISO 4049 > 60 sec.	Setting Time ISO 4049 < 10:00 min:sec	Shade ISO 4049	Color Stability ISO 4049
	7	pass	03:23	pass	pass

Test: Test methods: Limit/Unit:	Water Sorption ISO 4049 < 40 jig/mm ³	Solubility ISO 4049 <7.5 µg/mm ³	Flexural Strength ISO 4049 > 50 Mpa	Radiopacity ISO 4049 > 1 mm
	27 + 3.5	-4.4 + 2.7	90 + 6.5	1.94

Test specimen light cured with the LED Curing Light Elipar DeepCure-S

Please note that though being in scope for ISO 4049:2000, polymer-based luting materials that have an adhesive component with their structure like RelyX Unicem 2 Clicker /RelyX U200 are out of scope in ISO 4049:2009.

Stability certificate

Product Name: RelyX U200

Manufacturing Site: Seefeld, Germany

Test Method: Inhouse method: Rheometer

Storage Conditions: Room temperature, 36°C

Study Initiation Date: 05.03.2013 / 13.06.2013 / 04.09.2013

Study End Date: 04.09.2014/ 15.12.2014/04.03.2015

Shelf Life: The given shelf life of 18 months can be concluded from the following stability data:

Specification:

Parameter	Limit	Method
Setting time (start)	>0.5 min	Inhouse - method Rheometer

	Batch #508609	Batch #518424	Batch #527849
Date of Manufacturing	03/2013	06/2013	09/2013

Storage Conditions: Room Temperature

Storage time [months]	Batch #508609	Batch #518424	Batch #527849
0	3.1	3.1	2.8
3	4.2	3.1	2.7
6	4.0	3.0	2.6
9	3.6	2.6	2.9
12	3.0	3.0	3.1
18	3.4	3.0	2.9

Storage Conditions: 36 °C

Storage time [months]	Batch #508609	Batch #518424	Batch #527849
0	3.1	3.1	2.8
3	3.4	2.9	3.1
6	4.6	3.2	3.3
9	polymerized	polymerized	polymerized

RelyX U200 Automix

Product Name: RelyX Unicem 2 Automix / RelyX U200 Automix, Shade A2

Batch/Lot number: 630695

Test method: ISO 4049:2009, Classification: Type 2, Class 3 acc. to ISO 4049 (dual cure materials)

Test Conditions: Temperature: 23°C (± 1 °C); Humidity: 50 % (± 10 %)

Test Results

Test: Test methods: Limit/Unit:	Film Thickness ISO 4049 < 50/ μ m	Working Time ISO 4049 > 60 sec.	Setting Time ISO 4049 < 10:00 min:sec	Shade ISO 4049	Color Stability ISO 4049
	10	pass	03:18	pass	pass

Test: Test methods: Limit/Unit:	Water Sorption ISO 4049 < 40 pg/mm ³	Solubility ISO 4049 <1.5 pg/mm ³	Flexural Strength ISO 4049 > 50 Mpa	Radiopacity ISO 4049 > 1 mm
	24 + 2.4	-5.7+ 0.5	88 + 12	2.07

Test specimen light cured with the LED Curing Light Elipar DeepCure-S

Please note that though being in scope for ISO 4049:2000, polymer-based luting materials that have an adhesive component with their structure like RelyX Unicem 2 Automix / RelyX U200 Automix are out of scope in ISO 4049:2009.

Stability certificate

Product Name: RelyX U200 Automix

Manufacturing Site: Seefeld, Germany

Test Method: Inhouse method: Rheometer

Storage Conditions: Room temperature, 36°C

Study Initiation Date: 10.07.09 /10.9.2009/9.10.2009

Study End Date: 12.1.2011/14.1.2011/28.3.2011

Shelf Life: The given shelf life of 18 months can be concluded from the following stability data:
Specification:

Parameter	Limit	Method
Setting time (start)	> 0.5 min	Inhouse - method Rheometer

	Batch MP0029	Batch MP0030	Batch MP0037
Date of Manufacturing	07/2009	07/2009	10/2009

Storage Conditions: Room Temperature

Storage time [months]	Batch MP0029	Batch MP0030	Batch MP0037
0	3.1	3.4	3.4
1	3.2	3.3	3.2
3	3.5	3.6	4.0
6	3.7	4.0	4.2
9	4.4	4.7	4.3
12	4.7	5.0	4.3
15	4.5	5.6	5.7
18	4.2	6.1	6.1

Storage Conditions: 36 °C

Storage time (months]	Batch MP0029	Batch MP0030	Batch MP0037
0	3.1	3.4	3.4
1	3.8	4.2	3.6
3	5.4	4.9	4.9
6	4.9	4.7	base polymerized
9	base polymerized	base polymerized	

7 Manufacturing Information**7.1 Name and address of the manufacturer**

The name and address of the manufacturer are provided below:

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
Flow chart:

Description:

1. Ingredients are weighed and added together as directed in the process standard
2. Mixture is blended together as directed in the process standard
3. Mixture is inspected for conformity to specifications
4. Product primary packaging is filled according to process standard
5. Labeling applied to packaging and inspected for conformity to specifications
6. Secondary packaging is filled according to process standard
7. Labeling applied to finished package
8. Final Product is inspected and tested for conformity to specifications
9. Product is placed in inventory for distribution to customer

8 Product Standards Compliance

(International - Horizontal & Vertical Standards)

RelyX U200 / RelyX U200 Automix

This product has been evaluated using the processes described in the standards below and found to be in compliance relevant requirements:

Standard	Type of Standard	Compliance
ISO 10993-1:2009 ISO 7405:2008/AMD 1:2013	Biocompatibility Assessment	Y
ISO 13485:2003 + Cor 1: 2009	Quality Management	Y
ISO 14971:2007	Risk Management	Y
ISO 4049:2009	Product Efficacy/Performance	Y
ISO 15223-1:2012 EN 980:2008 EN 1041:2008	Labeling	Y

9 Biocompatibility Statement

RelyX U200 / RelyX U200 Automix

A Diplomate of the American Board of Toxicology has assessed the safety of these products. Standard risk assessment techniques and consideration of internationally recognized guidelines were used in this evaluation.

RelyX U200 Automix and Clicker are safe for their intended use based on the following considerations:

- 1) Biocompatibility test results for toxicologically equivalent products;
- 2) A review of the intrinsic hazards of individual components in relation to the amount used in the product;
- 3) Material characterization studies of sterilized and unsterilized toxicologically equivalent materials; and
- 4) The favorable history of clinical biocompatibility for a toxicologically equivalent product.

The biocompatibility assessment for this product was conducted in accordance with the following standards:

- 1) Testing guidelines outlined in the FDA General Program Memorandum G95;
- 2) ISO 10993-1:2009(E) Biological evaluation of medical devices – Part 1: Evaluation and testing; in addition, relevant detailed guidance in ISO Standards 10993-3:2014 (Tests for genotoxicity, carcinogenicity and reproductive toxicity), 10993-5:2009 (Tests for in vitro cytotoxicity), 10993-10:2010 (Tests for irritation and delayed-type hypersensitivity); and 10993-11:2006 (Tests for systemic toxicity) was considered;
- 3) ISO 7405:2008/2013 Amd Dentistry- Evaluation of Biocompatibility of Medical Devices in Dentistry;

- 4) Japan: PFSB/ELD/OMDE Notification No 0301-1 Mar. 1, 2012; (as translated by 3M Health Care Japan, August 6, 2012);
- 5) 3M ESPE Standard Operating Procedure 04-200.

RelyX U200 Automix and Clicker are externally communicating products intended for more than 30 days contact with tissue, bone or dentin (ISO 10933 and ISO 7405, G95) and coupling instruments between the inside and outside of the body (PFSB).

10 Sterilization

Sterilization is not applicable. 3M ESPE's RelyX U200 / RelyX U200 Automix is not labelled nor otherwise represented as sterile, nor is it intended to be sterilized by the user.

11 Protection of the environment

Avoid release to the environment.

12 Medicinal substances

Medical devices does not contain or deliver medicinal substances.

13 Human Blood Derivatives

Medical devices does not contain or deliver human blood derivatives

14 Animal-Derived Substances

No materials of animal origin are intentionally being added to the products.

15 Utilization/liquidation of MD

Disposal of disposable products, products with expired shelf life, products with expired service life is release in accordance with local/regional/national/international regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" to the waste class B.

This medical device isn't subject to maintenance and repair.

16 Transportation information

During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

17 Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

18 Risk Management Summary

March 31st, 2017

ISO 14971 Risk Management Report

RelyX U200 Clicker / RelyX U200 Automix

The environmental, health and safety (EHS) risks for RelyX U200 Clicker / RelyX U 200 Automix were evaluated using a process compliant with ISO 14971:2007 (Risk Management for Medical Devices) and specific procedures and practices outlined by division Standard Operating Procedures. In accordance with ISO 14971, division management has established acceptable levels of EHS risk for the products marketed by 3M Oral Care and it is these criteria against which the risks identified for RelyX U200 Clicker / RelyX U 200 Automix were evaluated.

In the risk assessment process, several sources of information were reviewed for potential risks or EHS information. These include: the product biocompatibility assessment, the clinical assessment, the complaint history of similar products (including both quality and clinical complaints) and published clinical toxicity literature on this or similar products and product components. In addition, changes to the product composition or manufacturing process were reviewed with respect to any EHS impacts they could introduce individually or collectively.

For the periodic safety update review of RelyX U200 Clicker / RelyX U 200 Automix conducted on September 8, 2011, there were a total of 4 EHS risks identified. After complete analysis and resolution of risk mitigation activities, all risks relevant for the present project phase were deemed to be acceptable. In the second periodic safety update review of RelyX U200 Clicker / RelyX U 200 Automix conducted on November 3, 2014, there were no EHS risks identified.

The risk management report was reviewed by the 3M Oral Care management.

A review of the product chemistry and division standard operating practices resulted in the following: Product Attribute	Division Practice
As sold, product contains uncured methacrylates, known dermal sensitizers	Material Safety Data Sheet and Instructions for Use contain information regarding sensitization and proper handling
As sold, product contains sodium peroxodisulfate, known sensitizers	Material Safety Data Sheet and Instructions for Use contain information regarding sensitization and proper handling
As sold the product shows acidity during application that could lead to post-operative sensitivities	Instructions for use contain information regarding use of the product in deep cavities (use of pulp protection)
As sold the base paste contains an acidic component that could cause eye irritation	Material Safety Data Sheet and Instructions for Use contain information regarding eye protection and proper handling

Product risks related to ergonomics and human factors were considered during the risk assessment process and none were identified.

Based on current knowledge RelyX U200 Clicker / RelyX U 200 Automix don't contain any ingredients of animal or human origin.

I hereby certify that the document
affixed hereto is a true copy
of the original.

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3M, ESPE and RelyX are trademarks of 3M or 3M Deutschland GmbH.   2016, 3M. All rights reserved.