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Carl-Schurz-Straße 1
41453 Neuss
Germany

Postal Address:
ESPE Platz
82229 Seefeld
Germany

☎+49 (0)8152 700-0
☎+49 (0)8152 700-1366
Internet: www.3MESPE.com
E-Mail: info3MESPE@3M.com
WEEE-Reg.-Nr. DE 36963167
VAT-ID: DE 120679179


To whom it may concern

March 28, 2017

Registration documents for Impregum Penta H DuoSoft - IFU

We 3M Deutschland GmbH, Carl-Schurz-Str. 1, 41453 Neuss, herewith confirm that the following pages are an excerpt of the Information for Use.

3M Deutschland GmbH
3M Oral Care

Christine Schalk
Regulatory Affairs

General Information

Impregum™ Penta™ H DuoSoft™ is heavy-bodied and light-bodied consistency polyether impression materials, respectively, to be used with the mixing device, Pentamix™. The mixing ratio is 5 volumes base paste: 1 volume catalyst.

Due to its lower Shore hardness, the impression material is suitable for use with impressions with slight undercuts even without additional blocking out.

Impregum Penta H DuoSoft and Impregum Penta L DuoSoft are not available in all countries.

☛ These Instructions for Use should not be discarded for the duration of product use. For details on Pentamix and accessories, Polyether Adhesive, Penta™ Elastomer syringe, Penta™ Mixing Tips - Red, and Impregum™ Garant™ L DuoSoft™, please refer to the corresponding Instructions for Use.

Indications

- Impression of inlay, onlay, crown, bridge, and veneer preparations
- Functional impressions
- Fixation and Implant impressions

Precautionary Measures

For Patients

The base paste of this product contains substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an adverse reaction is observed in a patient discontinue use of the product in this patient and seek medical attention as needed.

Known allergens are 1-Dodecylimidazole and mint flavors.

For Dental Personnel

The base paste of this product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended.

3M MSDSs can be obtained from www.mmm.com or contact your local subsidiary.

Preparation

Blocking out

- Block out undercuts or areas where gingival recession exists to prevent the impression material from "locking" onto the tooth structure. Failure to block out may make tray removal difficult or cause extraction of natural teeth or prostheses. Blocking out is especially important when using impression materials with high Shore hardness.

Prepare impression trays

All impression trays generally used for precision impressions are suitable.

- For sufficient adhesion, apply a thin layer of Polyether Adhesive to the tray and allow to dry completely according to the Polyether Adhesive Instructions for Use.

Prepare Pentamix/Penta DuoSoft cartridge/poly bag

- Place the Impregum Penta H and L DuoSoft poly bags in the corresponding designated 3M ESPE Penta cartridges only!
- Place the cartridge in the Pentamix mixing unit and attach a new Penta Mixing Tip - Red.
 - If a new mixing tip is installed when the cartridge is inserted, check whether the drive shaft is engaged before you start mixing.
- With newly filled cartridges, start mixing and extrude and discard the first unevenly mixed paste prior to the first use for impression taking. Do not use the paste to take an impression until the color of the paste is homogeneous.

Retraction

Retraction cords as well as solutions and pastes based on aluminum chloride, aluminum sulfate, or iron sulfate are suitable retraction agents. Recommendation: 3M ESPE Astringent Retraction Paste/Retraction Capsule. The paste containing aluminum chloride is directly applied from the retraction capsule into the sulcus. The paste can be used as an alternative to, or in combination with, retraction cords or other retraction methods, depending on the clinical situation and the practitioner's application technique.

- Areas from which impressions are to be taken should be kept dry.
- Prior to taking the impression, completely remove all residual retraction agents by rinsing and drying.

Dosing and Mixing

- Dosing and mixing are performed automatically in the Pentamix device.

Times

	Working Time from start of mixing* min:sec	Setting Time from start of mixing* min:sec	Intraoral Setting Time** min:sec
Symbols on product		-	
Impregum Penta H DuoSoft Penta Putty	2:30	6:00	3:30

The working times shown are in conformity with ISO 4823:2000. Higher temperatures will shorten the working time, while lower temperatures will prolong it. We recommend a working time of 1:45 min:sec and compliance with the indicated setting time from the start of mixing to assure an ideal impression result at all times under clinical conditions.

*Start of mixing = entry of paste into the mixing tip.

**For product combinations, longer times are valid.

Impression Taking

One-step technique

- For application around the preparation, attach the Penta Elastomer syringe to the mixing tip of the Pentamix device and fill with Impregum Penta L DuoSoft.
- Subsequently change cartridges and load Impregum Penta H DuoSoft onto the tray previously prepared with adhesive.

- While loading the tray, apply the mixture into the sulcus or cavity from bottom up, always **keeping** the tip of the application nozzle immersed in the material and in contact with the **surface** of the tooth.
 - Depending on the number of units to be syringed, initiate application around the preparation appropriately to finish syringing and tray filling at the same time. Failing this, the syringe material will set more rapidly than the tray material, leading to distortions of the impression.
- Immediately after application around the preparation, position the loaded tray in the mouth and hold it there without applying pressure until setting is complete.
- To remedy initial adhesion ("setting the valve"), especially with upper jaw impressions, remove the tray from the gingiva on one side in a posterior position. If this proves difficult, it may be necessary to carefully blow some air or water between the impression and the gingiva.

Functional Impression

- After preparing the custom tray with adhesive, load Impregum Penta H DuoSoft onto the tray, place the tray in the desired location in the mouth, and have the patient perform **functional** movements.
- If required, use Impregum Penta L DuoSoft or Impregum Garant L DuoSoft for precise impression taking (correction).

After Impression Taking

- Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual impression material from the mouth.

Hygiene

- Place the impression in a standard disinfectant solution commonly used for impressions for the period of time recommended by the manufacturer. **Excessive disinfection may damage the impression.**
- After disinfection, rinse the impression under running water for approx. 15 sec.

Model Preparation

- Prepare a cast from the impression with a commercial specialized stone plaster no earlier than 30 min and no later than 14 days after taking the impression.
- To avoid introducing bubbles into the model, briefly pre-rinse with water and dry with air before pouring.
- **Do not use surfactants** as these impair the quality of polyether impressions and are not required!
- Polyether impressions can be silver-coated, whereas copper-coating is not feasible.

Cleaning

- Paste that has not set may be removed with ethanol or by rinsing with water and soap.

Notes

- At temperatures below 18°C/64° F, the viscosity of the pastes may increase sufficiently to make mixing in the device difficult. After keeping the pastes at 18°C/64° F for at least one day, the processability is re-established without compromising quality.
- Direct exposure to sun and moist storage conditions damage the impression.
- Polyether impressions should not be exposed to solvent-containing liquids as this may result in swelling and imprecise modelling.
- 3M ESPE Polyether materials may only be combined with 3M ESPE polyether materials. When doing so, it is absolutely essential to observe the working and setting times of the individual products.

Storage and Stability

Store the product at 18-25°C/64-77° F. Do not refrigerate!

Do not use after the expiration date.

Store impressions dry and below 30°C/86° F in the dark.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

Information valid as of June 2016

RUS: Authorised representative of the manufacturer in RF:

ZAO "3M Russia", Russia 121614, Moscow, Krylatskaya str. bldg. 17/3 T.+74957847474

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
To whom it may concern

March 28, 2017

Registration documents for Impregum Garant L DuoSoft - IFU

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Product Description

Impregum™ Garant™ L DuoSoft™ is a light-bodied consistency polyether impression material supplied in the Garant cartridge. The mixing ratio is 2 volumes base paste: 1 volume catalyst. Due to its lower shore hardness, the impression material is suitable for use with impressions with slight undercuts even without additional blocking out.

For details on Garant™ Dispenser, Garant™ Mixing tips and Garant™ Intraoral Tips, Polyether Adhesive, Impregum™ Penta™ H DuoSoft™, please refer to the corresponding Instructions for Use

☛ These Instructions for Use should not be discarded for the duration of product use.

Indications

- Impression of inlay, onlay, crown, bridge and veneer preparations
- Functional impressions
- Fixation and Implant impressions

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Preparation

Garant dispenser/cartridge:

- Upon use of a new cartridge and with each new mixing process, start without attaching the mixing tip by extruding a small quantity of paste until both components emerge uniformly.
- Remove clogged paste, if necessary.
- Subsequently, install the Garant Mixing Tip White and Garant Intraoral Tip White.

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Retraction

Solutions based on aluminum hydroxide chloride or aluminum sulphate, are suitable retraction agents. Retraction with threads or rings containing epinephrine (adrenaline), 8-hydroxyquinoline sulphate or iron (III) sulphate may impair the setting of polyether impression materials.

- Thoroughly **examine** and explore the sulcus of the prepared teeth and surrounding dentition. Remove any **residual** cured impression material from the mouth.

Hygiene

- Place the impression in a standard disinfectant solution commonly used for impressions, for the period of time recommended by the manufacturer. Excessive disinfection may damage the impression.
- After disinfection, rinse the impression under running water for approx. 15 sec.

Model Preparation

- Prepare a cast from the impression with a commercial specialized stone plaster no earlier than 30 min and no later than 14 days after taking the impression.
- To avoid introducing bubbles into the model, briefly pre-rinse with water, and dry with air before pouring.
- Do not use surfactants as these **impair** the quality of polyether impressions and are not **necessary!**
- **Polyether** impressions can be silver-coated, whereas copper-coating is not feasible.

Cleaning

- **Garant Dispenser:** Remove paste that has not set using an alcohol-soaked cloth. The dispenser handle and plunger can be autoclaved up to a temperature of 135° C/275° F. Disassemble the dispenser before autoclaving. Glutaraldehyde-based solutions can be used for disinfection.
- **Impression Tray:** The adhesive can be removed from impression trays using acetone.

Notes

- Turning or reattaching the mixing tip without subsequent impression taking may lead to carry-over of paste and clogging.
- At storage temperatures below 18°C/64° F, the pastes may contract such way that precise dosing cannot be guaranteed even at room temperature.
- Direct exposure to sunlight and damp storage conditions may damage the impression.
- Polyether impressions should not be exposed to solvent-containing liquids, as this may result in swelling and imprecise modelling.
- 3M ESPE Polyether materials may only be combined with 3M ESPE polyether materials. Doing so the working and setting times of the individual products must necessarily be kept.

Storage and Stability

Store the product at 18-25° C/64-77° F. **Do not refrigerate!**

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Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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To whom it may concern

March 28, 2017

Registration documents for Intra-oral syringe - purple - IFU

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Product Description

Intra-oral Syringes - Purple are used to apply polyether impression materials around preparations. They have a maximum capacity of 1.5 ml of impression material and are filled individually with the aid of the hand dispenser (e.g., Garant™) right from 50 ml cartridges. The Intra-oral Syringes are available in the colors purple for 3M ESPE polyether impression materials (mixing ratio 2:1) and green for addition-type silicone impression materials (VPS) (mixing ratio 1:1).

- ☛ These Instructions for Use must be kept for reference for the duration of product use. Please refer to the relevant instructions for use for details about the impression materials used in the application syringe.

Indications

- Intra-oral Syringes - Purple: Syringing preparations with polyether precision impression materials in accordance with ISO 4823 Type 3 and Type 2, provided the polyether impression materials are intended for syringing.

The Intra-oral Syringes may not be used for impression materials with a higher consistency or for other products in paste form (e.g., impression materials for filling the tray in accordance with ISO 4823 Type 1 and Type 0 or pastes for bite registration).

Application

The mixing spirals and plungers are color-coded to ensure clear correlation of the Intra-oral Syringe to the correct material class (polyether or addition-type silicone impression materials): Purple = Mixing ratio 2:1 = 3M ESPE polyether impression materials (e.g., Impregum™, Permadyne™).

Green = Mixing ratio 1:1 = addition-type silicone impression materials (VPS).

The content of an Intra-oral Syringe is sufficient for syringing about 2-4 preparations, depending on the quantity applied around each unit. An Intra-oral Syringe must only be used on one patient!

- Have the hand dispenser with the inserted 50 ml cartridge ready for use (see Garant installation instructions, for example).
- Remove the original cap or the mixing tip left on the cartridge from the previous use from the 50 ml cartridge. Retain the original cap!
- Check the cartridge openings for possible blockages in accordance with the instructions for use for the impression material; press out a small quantity of base paste and catalyst to equal the pastes.
- Remove the plunger of the Intra-oral Syringe.
- Place the Intra-oral Syringe on the openings of the Garant cartridge and hold firmly in place while filling. Leave the mixing tip in an angled position; do not fold open.
- Use the hand dispenser to fill the Intra-oral Syringe with the required amount of impression material, but do not exceed the marking line (about 1 cm below the mixing tip joint).
- Remove the Intra-oral Syringe from the 50 ml cartridge and close the 50 ml cartridge, using the original cartridge cap.
 - When the cap is used more than once, make sure that the catalyst paste does not come into contact with the base paste to prevent paste carryover.

- Place the plunger in the Intra-oral Syringe, pushing the paste in the syringe barrel no further than to the joint of the mixing tip.
 - Any paste extruded on the sides of the mixing tip joint must be wiped off before the application in the mouth.
 - Once prepared as described above, the Intra-oral Syringe can be stored for as long as 12 hours. It may be necessary to prepare a number of Intra-oral Syringes in advance, depending on the number of preparations to be syringed, in the case of larger prosthetic work.
- Immediately before application, fold out the mixing tip of the Intra-oral Syringe at an angle of 180° (until resistance is felt).
 - The mixing tip must be folded out completely to prevent leakages at the mixing tip joint.
- Press out a small quantity of paste onto a mixing block until the base paste and catalyst flow out in an even, homogeneous mixture.
- Immediately apply the impression material around the preparation in accordance with the instructions for use.
 - **Refilling an Intra-oral Syringe which has been emptied or partially emptied is prohibited.**

Notes

- Intra-oral Syringes cannot be autoclaved.
- The Intra-oral Syringe - Purple with the mixing ratio 2:1 may be used solely for 3M ESPE polyether impression materials, and the Intra-oral Syringe - Green with the mixing ratio 1:1 may be used solely for addition-type silicone impression materials (VPS). **Filling the Intra-oral Syringe with the wrong paste will cause incorrect mixing results.**

Disposal

- Dispose of Intra-oral Syringes, just as other mixing tips, along with the other practice waste.

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Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

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