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APPROVE

Director

«ClaroNav Kolahi Inc.»

Ahmad Kolahi, CEO

«_01_»_December_2023

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SignNow e-signature ID: a008b866cf...
12/01/2023 19:59:27 UTC

User Manual

Surgical navigation equipment Navient, versions Otorhinolaryngological (95-NE-NE) and Universal neurosurgical and otorhinolaryngological (95-NC-NC)

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Emergo Europe
Prinsesgracht 20
2514 AP The Hague
The Netherlands

2797

Caution: U.S. Federal (and Canadian) law restricts this device to sale by or on the order of a licensed medical practitioner.

Contents

1	Indications, Safety and Warnings.....	8
1.1	Intended use.....	8
1.2	Contraindications	9
1.3	System Installation	9
1.4	Conventions used in this manual	9
1.5	General warnings	10
1.6	General cautions	11
1.7	Device and Label symbols	12
1.8	Safety warnings and compliance labels	14
1.9	Environmental standards	17
1.10	Data Security	17
1.11	General information	18
1.11.1	Transport	18
1.11.2	Laptop Battery Usage and Maintenance.....	19
1.11.3	Cleaning & Disinfection	19
1.11.4	Disposal	20
2	Basics	21
2.1	Principles of operations.....	21
2.2	Navient components and accessories.....	24
3	Safety Information.....	26
3.1	Preventing navigation inaccuracy	26

3.2	Cleaning and Sterilization prior to usage	27
3.2.1	Cleaning.....	27
3.2.2	Steam sterilization (Autoclave)	29
3.2.3	Disinfection.....	29
3.2.4	Cart cleaning.....	30
4	Hardware components.....	31
4.1	Navient cart	31
4.2	CameraBox (Tracking system)	32
4.3	Laptop computer	32
4.4	Accessory kit.....	33
4.4.1	Registrator	34
4.4.2	Maxillary Seeker	34
4.4.3	Frontal Seeker	35
4.4.4	Bayonet pointer.....	36
4.4.5	Instrument Tip Calibrator	36
4.4.6	Patient Tracker	37
4.4.7	Universal tag.....	37
4.4.8	Universal tag 360.....	38
4.4.9	Universal Adaptors	38
4.4.10	Silicone Tape.....	40
4.4.11	Cranial Articulated Arm	41
4.4.12	Dual Adaptor	42
4.4.13	Cranial Patient Reference.....	42
1.1.1	Biopsy Aiming Device	42

1.1.2	Reducing Tube	43
1.1.3	Biopsy Pointer	43
1.1.4	Cranial Pointer	44
5	Software Workflow	45
5.1	Workflow overview	45
5.2	Splash screen.....	45
5.3	Patient data loading and importing	46
5.3.1	Loading data	46
5.3.2	Importing Data	48
5.4	Stage 1: IMAGES.....	50
5.4.1	Image Fusion	52
5.4.2	Tool strips buttons.....	56
5.5	Stage 2: PLANNING.....	57
5.5.1	Landmarks	57
5.5.2	How to mark CT landmarks	57
5.5.3	How to modify marked landmarks.....	58
5.5.4	How to use fiducials as landmarks	58
5.5.5	3D Modeling	59
5.5.6	Biopsy	60
5.6	Stage 3: Set-Up	61
5.6.1	Tasks to be done in this stage	61
5.6.2	Calibrating an instrument tip	62
5.7	Stage 4: Registration	63
5.7.1	Landmark Based Registration.....	63

5.7.2	Trace Registration	64
5.7.3	Floating menu.....	65
5.7.4	Using Cranial Patient Reference Frame.....	66
5.8	Stage 5: Navigation.....	66
5.8.1	Accuracy Check.....	66
5.8.2	Navigation	67
5.9	Settings.....	69
6	PACS Configuration	72
6.1	Introduction.....	72
6.2	Navient DICOM Entity configuration.....	72
6.2.1	AET value	73
6.2.2	Host IP Address value	74
6.2.3	Port Number value	74
6.3	Remote DICOM device configuration.....	74
6.3.1	Add new DICOM remote device.....	74
6.3.2	Edit existing remote DICOM device	79
6.3.3	Delete existing remote DICOM device	80
6.3.4	Test remote device.....	81
6.4	Data Browser change after Remote DICOM device is configured	84
7	Service and Support	87
	Appendix A: The R-Fine Application	88
	R-Fine Quick Guide	88
	Appendix C: Troubleshooting Guide	91
	Appendix D: Adjusting the cart arm's joint friction.....	93

Appendix E: Electrical Safety and EMC Declarations94

1 Indications, Safety and Warnings

1.1 Intended use

Cranial edition	ENT edition
Navient is a computerized surgical navigation system intended as an aid for precisely locating anatomical structures in either open or percutaneous procedures. The Navient system is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure, such as the skull, long bone or vertebra, can be identified relative to a CT, MR based model of the anatomy. Example procedures include but are not limited to: Cranial procedures: • resection of brain tumors; intracerebral: supratentorial (frontal, parietal and temporal lobes, subcortical nuclei, lateral ventricles, 3rd ventricle of the brain), subtentorial (cerebellum, brain stem, including the midbrain and medulla oblongata, pons, 4th ventricle); - extracerebral: tumors of blood vessels, bones, membranes and neuromas of any location; • surgery on the base of the skull; • biopsybrain; • craniotomy; • puncture and drainage of the ventricles of the brain in adults and children;	Navient is a computerized surgical navigation system intended to guide sinus and trans-nasal skull base endoscopic surgery by dynamically presenting the location of the tip of a surgical instrument mapped to a corresponding location in a pre-acquired CT scan of the patient's head. The device is intended for use by a qualified ENT surgeon (otolaryngologist). Navient is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure, such as the skull, can be identified relative to images of the anatomy. Surgery procedures include but are not limited to the following: • Transphenoidal procedures • Maxillary antrostomies • Ethmoidectomies • Sphenoidectomies • Sphenoid explorations • Turbinate resections • Frontal sinusotomies • Intranasal procedures • Intranasal tumor resections • All ENT related skull base surgery

• implantation of a ventricular drainage system in adults and children; • thalamotomy. ENT procedures: • transsphenoidal procedures; • maxillary antrostomy; • ethmoidectomy; • sphenoidotomy; • resection of nasal turbinates; • frontal sinusotomy; • intralesional procedures.	
--	--

- **Note:** Navient Cranial edition is not available for the USA customers.

1.2 Contraindications

Medical conditions which contraindicate the use of Navient include any medical conditions which may contraindicate the medical procedure itself.

The use of Navient should be terminated when the accuracy check prior to the initiation of tool guidance has failed. For more details see section 5.8.1.

1.3 System Installation

Navient system is installed by qualified CKI personnel or distributors. Please contact CKI if you have any question about installation.

1.4 Conventions used in this manual

 CAUTION	Cautions are directions which, if they are not followed, can cause damage to the system described in this manual.
 WARNING	Warnings are directions which, if they are not followed, can cause fatal or serious injuries to a user, engineer, patient or any other person or can lead to a mistreatment.

1.5 General warnings

Caution: U.S. Federal (and Canadian) law restricts this device to sale by or on the order of a licensed medical practitioner.

The system and its associated applications should be used only by qualified medical professionals who are thoroughly trained and experienced in performing Navient-guided surgery.

Navient should be used only as an adjunct for surgical guidance. Navient is not a replacement for the surgeon's knowledge, expertise, or judgment.

The user is instructed to perform accuracy checks prior to navigation. If navigation guidance seems inaccurate and recommended steps to restore accuracy are not successful, abort use of the system.

Navient should not be used when there are noticeable patient motion artifacts in the CT image or if the CT image does not include any of the patient forehead, dorsum of nose, base of nose and eye laterals.

Inspect all system components before use. If visibly damaged, do not use the system.

Do not attempt to modify or alter the system or any of its components in any way, since doing so can impair its function.

All Navient cables listed below are in compliance with IEC 60601-1-2:2007 requirements. If spare cables are needed, order directly from the Navient manufacturer.

- 950-NE-LPPAE – Navient Laptop and IR LED Power Extension Cable
- 950-NE-FWCB – Navient Camera FireWire cable
- 950-NE-MPCNA – Navient Main power cable
- 950-NE-MPS – Navient Power Splitter

The use of accessories and cables other than those specified, with the exception of accessories and cables qualified and sold by the manufacturer of the Navient may result in increased emissions or decreased immunity of the Equipment and may cause the Navient to be non-compliant with the requirements of IEC 60601-1-2:2007

Navient should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Navient System should be observed to verify normal operation in the configuration in which it will be used.

1.6 General cautions

Do not install any other software on the laptop to prevent compromising the system performance.

Do not connect any other equipment to the laptop except the accessories provided (Thunderbolt connection for the tracking camera and DVD drive and an external memory USB drive)

While the laptop computer and the tracking camera may be operated on battery power alone for some period, the tracking system IR LED will remain turned off when operated on battery and therefore performing guidance on battery alone is not recommended. It is also possible that the battery will become fully discharged in the middle of the operation, causing the system to shut down.

When the system is connected to a main power source, make sure it is positioned in a way that would allow for a quick and easy disconnection from the mains, should the need arise.

Verify that all instruments (Registrator, Maxillary and Frontal Seekers, Bayonet Pointer, Calibrator, Silicon strap, Universal Tag, Clamp and thumb screw) have been steam sterilized (Autoclaving) according to the instructions provided in Section 0.

Patient Tracker must be disinfected before surgery following the instructions provided in Section 3.2. Do not heat the Patient Tracker above 80°C (180°F) or attempt to steam-sterilize it, since this may create geometrical distortions which will compromise the system's accuracy.

None of the electronic components of Navient is designed for sterilization, and may be damaged if sterilization is attempted. If contact of sterile glove surfaces with electronic components is needed during surgery, use a sterile barrier, such as a sterile disposable nylon sheet.

The system has been tested against the requirements of IEC 60601-1-2. However, RF interference could hamper its operation or the operation of other nearby electrical devices. If you suspect either of these conditions, move the interfering equipment more than 30cm (12in) away from Navient, separate the equipment with an RF barrier, or discontinue use of the system.


Avoid dripping any fluids into or on the laptop computer surface. Disconnect the power; turn the laptop computer off and allow the system to dry if you suspect fluids may have entered any part of the system.

During system operation, maintain the acceptable ranges for temperature and humidity (for details, see Section 1.9)

1.7 Device and Label symbols


Symbol	Definition
	Catalogue number
	System serial number
	Follow instructions for use
	Caution, consult documents
	Not suitable for regular trash disposal
	Certification mark

Symbol	Definition
	Manufacturer's name
	Batch code
	Upper temperature limit – Do not store in temperature above 45°C
	Non-sterile
	Re-use is not allowed
	Fragile item, handle carefully
	Do not expose to sunlight
	Humidity limitation
	Atmospheric pressure limitation
	Temperature limitation


Symbol	Definition
	Double insulation (class II device)
	Follow instructions for use
	Do not push the equipment while in use
	Power button


1.8 Safety warnings and compliance labels

The Navient system is fitted with a label. This is located at the rear panel of the system.


Navient ENT edition label (model 955-NE-NE) for Europe (top picture) and US & Canada (bottom picture)



Navient Cranial edition label (model 955-NC-NC) for Europe (top picture) and Canada (bottom picture)

1.9 Environmental standards

The following table details the recommended environmental standards required for correct operation of Navient.

	Operating temperature	15° to 30°C (59° to 86°F)
	Storage temperature	-25°C to 45°C (-13° to 113°F)
	Humidity	10% to 90% noncondensing
	Atmospheric pressure	700 - 1060 hPa (up to 2000 m)

1.10 Data Security

Navient loads patient datasets from external data sources using CD/DVD, USB or a network drive. Only use CD/DVD and USB drives created by secure sources. Do not connect them to the Navient laptop if they are not from well trusted sources as it may contain virus or malware and can damage Navient laptop. The DVD and USB drives should only include patient scan information generated by the scanning devices.

Navient comes with pre-installed Microsoft firewall and provides basic security when the system is connected to the network. Do not connect Navient laptop to an unsecure wireless network.

Note 1: Navient is provided with a secure password and user Identification number. It is recommended to change the laptop password at the first use (The default password for the first login is "Navient") and every 6 months afterward and to always lock the Navient laptop screen when left unattended.

Note 2: It is recommended that wireless access be turned off all the times except for loading data from the local network drives or online support. The Navient laptop is like any other computer on your network and should be protected according to your local IT security policy.

1.11 General information


Do not drop system components or mishandle them in such a manner as to cause them physical damage. This may impair device function. Even if the component works immediately after being dropped, operational damage may have occurred that may not be observed until some future time.

The system and its associated applications contain no user-repairable parts. For repair or replacement of any part of the system or application, contact a technical support representative.

1.11.1 Transport

When moving Navient ENT edition cart, fold the arm into the pole, and fold the tray with the laptop against the folded arm (see the picture). Hold the cart by its handle and wheel it to its new location.

To prevent the cart from rolling down a slope, place both wheels inside the two-wheel stoppers provided in the accessories kit, as shown in the picture below.

When moving Navient Cranial edition cart, push the extension pole down, dock the tracking camera and rotate the monitor in to 90deg angle. (see the picture). Hold the cart by its handle and wheel it to its new location.

To prevent the cart from rolling down a slope, lock the wheels by pushing the caster lock down.

1.11.2 Laptop Battery Usage and Maintenance

The laptop uses Lithium-Ion rechargeable batteries. For optimal performance, Navient should be plugged into a mains power outlet during surgical use, in which case the batteries will not be used to power the laptop.

1.11.3 Cleaning & Disinfection

When cleaning Navient, first shut down the computer and unplug the power cable. Use a damp, soft, lint-free cloth to clean the cart and the laptop exterior. Avoid getting moisture in any openings. Do not spray liquid directly on the laptop. Don't use aerosol sprays, solvents, abrasives, or cleaners containing hydrogen peroxide that might damage the surface finish of the parts.

Use a damp cloth soaked with alcohol solution (70% or higher) to disinfect surfaces that may come into occasional contact with liquids during surgery, such as the laptop's keyboard and screen, the camera box and the front pole holding them. In preparation for surgery, ensure that all surfaces that may need to be touched by sterile gloves, such as the tracking camera, the laptop trackpad and the front pole, are covered with sterile sheets.

Navient Patient Tracker is a reusable part. Between uses, it should be scrubbed with a wet and soapy brush and thoroughly rinsed with deionized (DI) water to remove any visible residual and dried with a clean, soft cloth. Examine the part for visible soil. If present, repeat cleaning. Once the part is cleaned it should be disinfected using a nylon compatible disinfectant, such as Metricide, Sporox or Cidex. Patient tracker tape strips are single use parts and should be disposed right after using. Place a new tape strip on the Patient Tracker after the Patient Tracker was disinfected.

Do not heat the Patient Tracker above 80°C (180°F). Steam-sterilizing it may damage the part and create geometrical distortions which will compromise the system accuracy in unpredictable ways.

Navient Universal Tag Adapter should be mounted on the surgical instrument and sterilized together with it using autoclaving.

1.11.4 Disposal


This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact CKI or your distributor for information concerning the decommissioning of your equipment.

If assistance is desired in the proper disposal of the Navient product (including unused disposable accessories and components), after its useful life, please return to CKI or your distributor.


2 Basics

2.1 Principles of operations

Navient's guidance function is based on the patient images acquired prior to the procedure, combined with optical measurements of the pose of navigated instruments relative to the patient's head. To enable navigation, the Patient reference is attached to the patient's forehead (on ENT edition) or skull holder (on Cranial edition) to enable tracking of the patient's head. The patient images are then spatially registered with the patient's head by matching landmark locations marked on both the image and the patient, followed by matching a path traced by the user on the patient's face with a model of patient's face surface automatically generated from the image data. Subsequently, Navient presents the location of a navigated instrument tip as mapped on the image data through the registration transform.

Patient Reference attached on patient forehead (ENT edition)
attached on skull holder


Patient

Reference


Navient provides the following instruments:

ENT Instruments	Cranial Instruments	Biopsy Instruments
1. Patient Tracker (2 sets) 2. Maxillary Seeker 3. Frontal Seeker 4. Bayonet Pointer 5. Registration Pointer 6. Calibrator 7. Universal Tag (2 sets) 8. Four C-Clamps	1. Articulated arm 2. Two cranial frames 3. Cranial Pointer 4. Registration Pointer 5. Calibrator 6. Universal Tag 7. Four C-Clamps	1. Articulated arms 2. Dual adaptor 3. Biopsy aiming device 4. Biopsy Aiming Pointer 5. Universal Tag 6. Reducing tube 2.2mm

Navient ENT Instruments

Optionally, a second set of Instrument kit, called Set B, is included in the Navient system. Set B would allow using navigation on back-to-back OR cases without having to wait for sterilization of the main set. Like the main set of instruments, Set B includes navigation instruments, adaptors and patient trackers (see the picture above). Set B part numbers are similar to the main set, except it is affixed with "-B". It is also marked with a different serial number.


Navient Cranial Instruments

All navigation instruments are pre-calibrated except the universal tags. Universal Tags can be attached to any surgical instrument using the C-Clamps. The location of the tip of the instrument relative to the tag is calibrated by placing the tip at a designated location on the Calibrator for a few seconds.


The Tracking CameraBox contains a stereoscopic Navient camera and an IR light source, both operating in the near infrared (NIR) light spectrum. The tracking software identifies and pin-points the center of checkered-pattern targets, called XPoints, marked on the instruments or the Patient Tracker. The constellation of XPoints marked on an object, such as an instrument or Patient Tracker, is called a "marker". Software running on the laptop computer processes the stereoscopic video images captured by Navient Camera to detect and triangulate the positions of XPoints. The relative positions and orientations of the detected XPoints are compared to marker templates stored in the computer memory to identify the instrument and its pose relative to the tracker. At least 3 or 4 XPoints need to be detected for a successful marker identification and pose measurement. The location and orientation of the tip of each instrument relative to its XPoints is also stored in the marker template, and is used in mapping the tip to the CT image.

2.2 Navient components and accessories

The picture below shows a fully assembled Navient ENT and Cranial systems. The ENT system includes a MacBook Pro laptop and a tracking camera mounted on the Navient cart. The cart arm can be folded to make it compact and easily portable.

The Cranial system includes a PC, 27" touch screen monitor and a tracking camera all mounted on a variable height cart. The height of the cart can be adjusted by stepping on a foot pedal between the two front casters.

The PC runs the Navient software which displays the position of pre-calibrated instruments superimposed on the pre-operative patient scans. The Tracking system contains an infra-red (IR) light emitting diodes (LEDs) lighting panel and the Navient Camera, a stereoscopic tracking camera, which feeds a video stream to the laptop for further processing.

3 Safety Information

3.1 Preventing navigation inaccuracy

While a properly used navigation system reduces the likelihood of surgical complications, the system is fallible and may occasionally misguide you. The system requires navigation accuracy verification at regular intervals.

Navient should not be used to make important surgical decision without first validating its correctness using recognizable surface and landmarks.

To minimize the risk of injury to the patient due to incorrect localization, navigation accuracy must be verified at regular intervals.

You can reduce the likelihood of navigation errors by being aware of the key factors degrading accuracy and acting to minimize or prevent them:

Registration error: Mis-registration is the most common reason for localization errors. Following the registration step, immediately perform an accuracy check, and, in case of doubt, redo the registration (see section 5.8).

Instrument tip calibration error: Navient tracks the position of markers (the checkerboard patterns on each tool), rather than the tips of instruments. To map the location of the instrument tip to the CT image volume, Navient needs to accurately know the position of the instrument tip relative to the relevant marker. For accurate operation, it is critical that that position would be validated, and, as needed, correctly calibrated.

Patient Tracker movement: Movement of the Patient Tracker relative to the skull will invalidate the registration. When mounting it, clean the patient forehead using alcohol wipe and ensure that the tag is well attached to the skin and avoid touching it for the duration of the surgery. In case you suspect it may have moved, perform an accuracy check and, as needed, repeat the registration.


Tracking camera mishandling: The tracking system is a high precision instrument and should be treated carefully. When handling its enclosure, avoid shaking, dropping or mishandling of the device. Do not open the enclosure. In case you suspect the tracking system may not be operating correctly (e.g., following an accidental shock or when marker detection becomes unreliable), you should check tracking system accuracy by performing the R-Fine procedure. The R-Fine application is installed on the Navient computer and is activated using a desktop icon (see Appendix A for details).

3.2 Cleaning and Sterilization prior to usage

3.2.1 Cleaning

The following instructions for cleaning must be followed prior to use for the following Navient parts:

- 950-NE-BP Bayonet Pointer
- 950-NE-RP Registration Pointer
- 950-NE-FS Frontal Seeker
- 950-NE-MS Maxillary Seeker
- 950-NE-UT1 Universal Tracker 1
- 950-NE-UT1 Universal Tracker 2
- 950-NE-TC Tip Calibrator
- 950-NE-PT Patient Tracker
- 950-NE-UTC Small C-Clamp
- 950-NE-UTCM C-Clamp Medium
- 950-NE-UTCL C-Clamp Large
- 950-NE-UTCXL C-Clamp X Large
- 950-NE-TS Thumb Screw
- 950-NE-TSL Thumb Screw Long
- 950-NC-CAA Cranial Articulated Arm
- 950-NC-CRF Cranial Reference Frame
- 950-NC-BAD Biopsy Aiming Device
- 950-NC-SP Straight pointer
- 950-NC-BP Biopsy probe
- 950-NC-UT1 Universal Tracker 1
- 950-NC-UT1 Universal Tracker 2
- 950-NC-TC Tip Calibrator

Warnings:

- Clean the instruments prior to first use and after each use.
- Before sterilization, disassemble the Universal Tag from its Adaptor or C-Clamp and the thumb screw(s).
- Remove and properly dispose of the single-use silicon tape used for Universal Tag and the Patient Tracker dual sided tape.

Manual Cleaning Instructions:

1. Promptly and thoroughly rinse the instrument under running tap water to remove any visible soil. This step should be done at the Point-of-Use to prevent the drying of soil on the device surface. A soft-bristled brush or clean cloth may be used to aid in soil removal.
2. Prepare an enzyme-based detergent according to the detergent manufacturer's recommendations, such as Medline Single Enzymatic Detergent and Presoak.
3. Immerse the instrument in the enzymatic detergent solution and let it soak for a minimum of two (2) minutes.
4. Thoroughly clean the instrument using a soft-bristled brush and remove all foreign matter.
5. Flush crevices, lumens, and blind holes in the instrument with fresh detergent solution until the solution runs clear.
6. Thoroughly rinse the instrument in deionized (DI) water for a minimum of one (1) minute to remove any residual detergent.
7. Flush crevices, lumens, and blind holes in the instrument with DI water.
8. Dry the instrument with a clean, soft cloth.
9. Examine the instrument for visible soil. If present, repeat cleaning.

Automated Cleaning Instructions:

1. Promptly and thoroughly rinse the instrument under running tap water to remove any visible soil. This step should be done at the Point-of-Use to prevent the drying of soil on the device surface. A soft-bristled brush or clean cloth may be used to aid in soil removal.
2. Prepare an enzyme-based detergent according to the detergent manufacturer's recommendations, such as Medline Single Enzymatic Detergent and Presoak.
3. Add cleaning solution to the tank of water (Allowing for the volume of the parts you will be cleaning and operation level line)
4. Place the instruments into a basket, slowly immerse it in the cleaner bath, activate ultrasonic machine and set the timer for 5 to 10 minutes and the temperature as per manufacturer's recommendations. (Do not allow instruments to contact the tank bottom).
5. When cleaning is finished, slowly remove them from the cleaner.
6. Thoroughly rinse the instrument in deionized (DI) water for a minimum of one (1) minute to remove any residual detergent.
7. Flush crevices, lumens, and blind holes in the instrument with fresh deionized (DI) water until the solution runs clear.
8. Dry the instrument with a clean and soft cloth and place it on clean bench.
9. Examine the instrument for visible soil. If present, repeat cleaning.

3.2.2 Steam sterilization (Autoclave)

The following instructions for sterilization must be followed prior to use of all instruments except the Patient Tracker.

Warnings:

- Instruments must be adequately cleaned prior to sterilization.
- Remove the plastic backing of the silicon tape prior to sterilization.
- Before sterilization, disassemble the Universal Tag from its Adaptor or C-Clamp and the thumb screw(s).

Instructions:

1. Instruments should be individually wrapped in two layers of 1-ply polypropylene wrap (Kinguard KC200 – 510(k) K082554) using sequential envelope folding techniques.
2. Cut silicon tape in pieces of 6inch (15cm) long, remove its plastic backing and wrap them individually.
3. The sterilization process should be done according to the following parameters:
 - Sterilizer Type: Prevacuum
 - Preconditioning Pulses: 4
 - Temperature: 132 °C
 - Exposure Time: 4 minutes
 - Dry Time: 30 minutes

Note: The Cranial Articulated Arm should be draped with thick plastic sterile drape.

3.2.3 Disinfection

The following instructions for disinfection must be followed prior to use for the Navient 950-NE-PT Patient Tracker.

Warnings:

- Remove and dispose the dual sided patient tracker tape
- Patient Tracker must be adequately cleaned prior to disinfection.
- Do not sterilize patient tracker. It will be damaged if heated above 80°C (176°F).

Instructions:

1. Follow the Cleaning instructions above and thoroughly clean the Patient Tracker.
2. Prepare FDA cleared/CE marked nylon compatible disinfectant (such as Metricide, Sporox or Cidex) according to the manufacturer recommendation.
3. Soak the patient tracker inside the disinfectant for the duration recommended by the manufacturer.

4. Remove the Patient Tracker from the solution and rinse thoroughly with minimum of 2 liters of sterile water for a minimum of one (1) minute to remove any residual disinfectant.
5. Thoroughly dry it using a soft, lint-free sterile cloth and remove any residual sterile water.
6. Wrap the patient tracker with a lint-free sterile cloth and enclosed in a sealable bag or container. If instrument will not be used immediately, it should be stored in a well-ventilated, dust-free closed cabinet in a manner that minimizes contamination or damages for up to 7 days, after which time the instrument should be disinfected again before use.

3.2.4 Cart cleaning

When cleaning the Navient cart, first shut down the computer and unplug the power cable. Use a damp, soft, lint-free cloth to clean the cart and the laptop exterior. Avoid getting moisture in any openings. Do not spray liquid directly on the laptop. Don't use aerosol sprays, solvents, abrasives, or cleaners containing hydrogen peroxide that might damage the surface finish of the parts.


Use a damp cloth soaked with alcohol solution (70% or higher) to disinfect surfaces that may come into occasional contact with liquids during surgery, such as the laptop's keyboard and screen, the camera box and the front pole holding them. In preparation for surgery, ensure that all surfaces that may need to be touched by sterile gloves, such as the tracking camera, the laptop trackpad and the front pole, are covered with sterile sheets.

4 Hardware components

4.1 Navient cart

Navient cart's main purpose is to hold the monitor and the tracking camera. The cart is light weight and the arm height can be adjusted. This allows the user to easily position the system in the operating room.

On the ENT cart, the arm is opened by placing a foot over a black tab near the front of the base while pulling the front pole away from the back post. The arm is folded back by pushing it downward. The arm's friction can be adjusted (see Appendix D).


On the Cranial cart, the extension pole raises by pushing down the black pedal between the front two casters. The pole is contracted by holding the foot pedal down and pushing the extension pole down.

4.2 CameraBox (Tracking system)


The CameraBox contains a stereoscopic Navient camera and an IR LED light source, both operating in the near infrared (NIR) light spectrum. The tracking system identifies and pin-points checkered marker patterns. It is important to have a clear line of sight between the tracking camera and the targets. There are two status indicator diodes in front of the CameraBox:

1. Green diode indicates the status of the Navient camera:
 - a. Dim (Normal): Camera not activated and not streaming images.
 - b. Bright (Normal): Camera is activated and streaming images.
 - c. Flashing (Abnormal): Camera has hardware or firmware issues. Contact support.
2. Blue diode indicates the status of the IR light:
 - a. ON (Normal): IR light is powered on and functioning.
 - b. OFF (Abnormal): means the IR light is not powered on. Make sure the main power cord is connected to the power outlet and the DC power socket is plugged in on the back of the CameraBox.


4.3 Laptop computer

The Navient software runs on a 15.4" Apple Mac Book Pro laptop. Windows 10 operating system is installed on the laptop. The tracking camera's live video stream is received by the laptop via a thunderbolt port and is processed by the Navient application.

4.4 Accessory kit

Navient system includes a set of accessories. All parts are steam sterilisable except the Patient Tracker. All instruments are marked with black and white patterns which make them trackable by the tracking system. It is important to have a clear line of sight between the tracking system and the targets.


Note: Second set of Navigation Instruments, called set B, can be ordered separately. Set B would allow using navigation on back-to-back OR cases without having to wait for sterilization of the main set.

4.4.1 Registrator


The picture below shows Navient registration pointer, called the Registrator. The Registrator is made of aluminium and is marked with black and white patterns that enable it to be detected by the tracking system. During registration, the surgeon touches patient landmarks with the Registrator's tip, and performs the trace registration by gliding the tip on the patient forehead.

Holding the Registrator above the patient's nose activates the floating menu.


4.4.2 Maxillary Seeker

Picture below shows Navient Maxillary Seeker instrument. Maxillary seeker is a common instrument to reach maxillary sinus area. The Navient maxillary seeker includes an ergonomic attachment made of aluminium with black and white pattern that enables it to be detected by the tracking system. It is important to have a clear line of sight between the tracking system and the black and white patterns. Once detected, the Navient software will display the real time location of the instrument tip super-imposed on the patient CT data.


4.4.3 Frontal Seeker

Picture below shows Navient Frontal Seeker instrument. Frontal seeker is a common instrument to reach frontal sinus area. The Navient frontal seeker includes an ergonomic attachment made of aluminium with black and white pattern that enables it to be detected by the tracking system. It is important to have a clear line of sight between the tracking system and the black and white patterns. Once detected, The Navient software will display the real time location of the instrument tip super-imposed on the patient CT data.

4.4.4 Bayonet pointer

Picture below shows the Bayonet pointer instrument. It is a bayonet shape pointer with a long tip made of stainless steel to allow surgeon to localize internal anatomies. The black and white pattern enables it to be detected by the tracking system. It is important to have a clear line of sight between the tracking system and the black and white patterns. Once detected, The Navient software will display the real time location of the instrument tip super-imposed on the patient CT data.

4.4.5 Instrument Tip Calibrator

Picture below shows the Instrument Tip Calibrator. Its main purpose is to calibrate instruments tip and make them trackable. A black and white pattern on the calibrator enables it to be trackable by the tracking system. To calibrate an instrument tip, the tip must be placed on the small dimple on the calibrator and then both the instrument and the calibrator must be presented to the camera. It is important to have a clear line of sight between the tracking system and the black and white patterns. Once instrument and calibrator are detected, the software will display a progress bar indicating the tip calibration process.

Dimple to place the instrument tip


Calibrator handle

progress bar indicating the tip

4.4.6 Patient Tracker

Picture below shows the Patient Tracker. It is made of plastic and is attached to the patient forehead with a double sided tape. It remains on the patient forehead for the entire duration of operation. The black and white pattern enables it to be trackable by the camera. It is important to have a clear line of sight between the tracking camera and the black and white patterns. Its main purpose is to locate the spatial location of the patient head.

Important Note: It is recommended to clean the patient forehead using alcohol wipe and remove any excess oil from the skin before placing the patient Tracker. It is important that the Patient Tracker shouldn't move or bent during and after registration.



4.4.7 Universal tag

Universal tag is designed to be attached to a shaver or endoscope instruments and enable them to be trackable by the tracking system. It does not contact patient. A silicon tape is used to attach the tag to the surgical instruments. Once attached, the user will need to go through the instrument tip calibration process. Once calibrated, the Navient software will display the real time location of the instrument tip super-imposed on the patient CT data.

4.4.8 Universal tag

360

Universal tag 360 is designed to be attached to a shaver or endoscope instruments and enable them to be trackable by the tracking system. It does not contact patient. A silicon tape is used to attach the tag to the surgical instruments. Once attached, the user will need to go through the instrument tip calibration process. Once calibrated, the Navient software will display the real time location of the instrument tip super-imposed on the patient CT data.

4.4.9 Universal

Adaptors

Navient Universal adaptors include four **C-clamps** and one **Doggie Adaptor**. They are intended to connect the universal tags to any surgical instruments and convert them to the navigated instruments.

Navient User Manual

C-clamps are used to clamp to any surgical instruments with 4mm to 25mm diameter. Each C-Clamp has two thumb screws:

- Small thumb screw to connect the universal tag to the C-clamp
- Long thumb screw to secure the C-clamp to the surgical instrument

C-Clamps come in 4 different sizes:

1. Small: for 4-10mm diameter instruments
2. Medium: for 10-15mm diameter instruments
3. Large: for 15-20mm diameter instruments
4. X Large: for 20-25mm diameter instruments

Below are the instructions for using the C-clamps:

1. Select the right size of C-clamp for your instrument.
2. Use the small thumb screw and attach the universal tag on the C-clamp.
3. Attach the large thumb screw.
4. Place the instrument on the C-clamp and tighten the thumb screw.

Important note: The instrument should be held firmly by the thumb screw pushing it against the "V" groove on the adaptor. Verify the final assembly and make sure the universal tag is tightly connected to the instrument.

If none of the C-clamps fit the surgical instrument, then a Doggie adaptor could be utilized. Doggie adaptor can be attached to any surgical instrument using self-adhesive silicon tape.

Below are the instructions for using the Doggie adaptor:

1. Cut silicon tape in 15-20cm (6-8 inch) long, remove the backing and have it sterilized.
2. Use the small thumb screw and attach the universal tag on the doggie.
3. Find the right location to attach the universal tag on the instrument. Make sure the tag will be always facing to the tracking camera.
4. Wrap the silicon tape around the adaptor. Pull the silicon tape and wrap it tightly. Cover all 4 legs of the "doggie" adaptor.
5. Verify the final assembly and make sure the universal tag is tightly connected to the instrument.

Doggie adaptor connected to a multi-debrider

4.4.10 Silicone Tape

It is a self-adhesive strong tape which is used to attach the universal tag to the surgical instruments like shaver and endoscope. Prior to use it should be steam sterilized. Cut silicon tape in pieces of 15-20cm (6-8inch) long, remove its plastic backing and wrap them individually and sterilize them according to section 3.2.

Silicon tape is a single use part and should be disposed after the use.

Silicon tape is steam sterilizable

Silicon tape warped for sterilization

4.4.11 Cranial Articulated Arm

It is a custom-made flexible arm to clamp the patient tracker to the Mayfield of the patient table. It is made of stainless steel and aluminium and is steam sterilizable.


4.4.12 Dual Adaptor

It connects to the skull holder and provides two ports to connect two articulated arms, one for the patient reference frame and the other for the aiming device. It is made of Al and is steam sterilizable.


4.4.13 Cranial Patient Reference

It is a custom-made patient reference made of Al and steam sterilizable.

1.1.1 Biopsy Aiming Device

It is used to aim the biopsy needle to the target point. It is made of aluminum and is sterilizable.

1.1.2 Reducing Tube

It is made of stainless steel and is steam sterilizable.

1.1.3 Biopsy Pointer

It is a custom-made short tip probe made of aluminum and is sterilizable. It is used to help with adjusting the biopsy aiming device.

1.1.4 Cranial Pointer

It is a custom-made straight probe made of aluminum and is sterilizable.

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Navient Repository is a folder where the Navient application stores all imported datasets locally. Section 5.2 explains importing patient scan data from CD/DVD to the repository and loading it to the application for guidance.

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
1. **Images:** Load the CT or MR images, review and explore images and perform image fusion.
2. **Planning:**
 - a. **Landmarks:** Mark easily identifiable anatomical landmarks on the patient scan.
 - b. **3D modeling:** Generate 3D model of the tumor or sinus.
 - c. **Biopsy:** plan biopsy trajectory or define target points for navigation.
3. **Set Up:** attach the Patient reference, set up the tracking system position and verify the instruments' tip calibration accuracy.
4. **Registration:** Using registrator, point to the marked landmarks on the actual patient's skin, then trace a path on the patient's face to refine the registration.
5. **Navigation:** Observe the position of the tip of various instrument in the image volume.

5.2 Splash screen

When started, Navient displays a "splash" screen that provides the software version number and license information. The splash screen requires the user to select the procedure type. Navient supports trans-cranial and trans-nasal procedures. If the system includes more than one set of

instruments, then a dropdown box will be displayed. The user must select the serial number of instruments in order to proceed.

5.3 Patient data loading and importing

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
5.3.1 Loading data

The **Data Browser** searches through the **Navient Repository** folder and shows the search results in two panes. On the top pane it lists the patient information for all patients for which some DICOM-formatted images were found. The browser auto selects the first discovered

patient, but the user can click on other rows to change that selection. On the pane below, the discovered image series (slice stack) for the selected patient are shown.

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1. Click on the Import button to open up the folder browser.
2. Select the device that you want to import the data and browse and select the patient data folder.
3. Click **Select Folder** to return to the Data Browser form. Data Browser will display all scans found under the selected folder.
4. Select the scan that you want to import and click on **Import & Load** (see picture below). It will import the patient data to the local **Navient Repository** and then will load the data from the repository.

The screenshot shows the 'Data Browser' window. At the top, there's a search bar and a status bar indicating 'Found 1 Patient'. Below this is a table with patient information:

Patient ID	Patient Name	Sex	Age	D.O.B.	Study Date	Description
Anonymized	Navient Fusion MR	M	?	?	0001-01-01	

Below the patient table is a section for 'Scans' with a table of scan details:

Modality	Images	Thickness	Spacing X/Y/Z	FOV HxW	Derived?	Date	Description
MR	100	2.00 mm	0.94x1.99 mm	20x24 cm	N	0001-01-01	

At the bottom of the window, there is a button labeled 'Import & Load' with a circular arrow icon. A text label 'Import selected patient scan to the local repository and load' points to this button. Other buttons at the bottom include 'Help', '?', and 'Cancel'.

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To avoid this from happening again, please do not remove the white adapter cable from the laptop.

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

- A. **Verify Patient info:** Make sure the patient info displayed on the side bar panel matches with the actual patient.
- B. **Adjust 2D and 3D images contrasts:** use the toolbar contrast button and adjust the windowing width and center of images. The software will preserve all changes from this stage and will display the exact same contrast during navigation stage.
 - C. Double clicking on the Contrast button resets the contrast to the default value.
 - D. Click on the Contrast button and drag up/down on the 2D or 3D images to change windowing center. Drag left/right to change windowing center.
 - E. Click on Presets button to open the presets drop-down menu. The menu lists available 3D presets options for rendering 3D model. Select an appropriate preset to visualize 3D model as desired.
- F. **Correct 3D image orientation:** The patient 3D image must be facing forward as displayed below. Use the Rotate toolbar to orient the 3D image. This initial orientation will be used as reference for proper presentations in the following stages.

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


It is important to review and make sure the two images are fused accurately. Select Axial, Coronal and Sagittal views to review the accuracy of fusion. If the images are not fused

accurately, then go back to "Merge" stage and repeat marking the landmarks and the image merging process.

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

Refer to Navient Scanning Instructions for a comprehensive Scanning instructions.

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Note: A double click (double tap on touch screen) on a toolstrip button will reset the parameters it controls to their default value.

5.5 Stage 2: PLANNING

Planning stage has 3 substages:

- Landmarks
- 3D modeling
- Biopsy


Landmarks are mandatory but the two are optional.

5.5.1 Landmarks

In this stage, the user must place 4 - 7 anatomical landmarks on the patient scan. Right click (Tap and hold on touch screen) to place a landmark. Once at least 4 landmarks were marked, the Set-Up button will be activated enabling to proceed to the next stage.

5.5.2 How to mark landmarks

1. Identify the landmark. Use the toolbar buttons (like scroll, rotate and contrast buttons) to identify the landmark on the patient scan.
2. Right click (Tap and hold on the touch screen) on the landmark, to mark the landmarks.
3. Review the location on 2D and 3D images to make sure the correct location is selected.
4. Click on "Save" button to save the landmark.

5.5.3 How to modify marked landmarks

1. From the Landmarks list select the landmark that you need to modify
2. Once selected, the corresponding cross bar on the main window will be highlighted
3. Right click on the new location and then click "Save" to overwrite the new location.
4. If necessary, use tool strip buttons to rotate, zoom, scroll, or pan the image.


5.5.4 How to use fiducials as landmarks

Using CT or MR fiducials for registration is highly recommended if the navigation is intended in the back of the cranium. The fiducials are recommended to be placed around the intended surgical field (i.e. back of the skull if the navigation is intended to be done on the back of the cranium). A minimum of 4 fiducials are need (7 landmarks recommended). The patient must be scanned with the fiducials and the fiducials location shouldn't be changed. The user needs to right click with a mouse to mark each landmark. Once all fiducials were marked, the **Set-up** button will be activated enabling to proceed to the next stage.

1. Identify which fiducials you are going to mark. Use Contrast toolstrip button to make fiducials detectable.
2. Place the cursor over the fiducial in the 3D image and right click.
3. Click "Save" button to store the landmark.

5.5.5 3D Modeling

This is an optional stage and allows the user to segment and generate 3D model of tumor or sinus. The software will display the 3D model on the navigation window.


Instructions for creating 3D model:

1. Scroll and find the tumor or sinus on any of the 2D viewports
2. Right click (Tap and hold on touch screen monitor) on the middle of tumor. The software will automatically segment the tumor and will display green contour around that.
3. Modify the contour by:
 - a. Right clicking on another place on the tumor on the same 2D viewport or any other 2D viewport
 - b. Use the slider and modify the size of the 2D segmented object. This will automatically update the 2D segmentation.
4. Click on "3D model" button to generate 3D model of the tumor. Tumor volume will be displayed.
5. To edit the model, click on the "Edit" button, and then drag the yellow contour on the 2D viewports. The 3D viewport will be updated automatically.
6. Click on "Smooth" button to dilute the 3D model and make it smoother.
7. Clicking on the "Clear" button will delete all the models.
8. Use alpha blending slider to control the contrast of the 3D model.

5.5.6 Biopsy

This is an optional stage and allows the user to define biopsy trajectory or define target points. The biopsy trajectory or target point will be displayed during the navigation stage.

Instruction for planning biopsy trajectory:

1. Use the scroll button and find the tumor location on 2D viewports.
2. Right click (Tap and hold on touch screen monitor) on the center of the tumor.
3. Click "Set target" button to store the target location.
4. Scroll on the 2D viewports or rotate 3D viewport and identify the entry point.
5. Right click on the entry point.
6. Click on "Set Entry" to store it as entry.
7. Once both target and entry are set, the software will display a trajectory line from entry to the target. The distance between the two points will be displayed on a slider.
8. Drag the slider to review the trajectory. Click "play" button to play it.
9. Enable "probe eye view" to see the trajectory from a perpendicular view.


5.6 Stage 3: SET-UP

This stage is designed to help the user with setting up the patient, positioning the tracking system and verifying/calibrating the instruments tips. The live video from tracking system is displayed in the side panel. It helps the user to orient and re-position it toward the field. A proper warning/instruction will be displayed on the live video if it is not positioned properly.

5.6.1 Tasks to be done in this stage

1. **Verify tracking system:** Make sure the tracking system is active and the live video is displayed. Follow the troubleshooting instruction if it is inactive.
2. **Attach the Patient Reference:** Stick the Patient Tracker to the patient's forehead (for trans-nasal procedure), or mount it on the articulated arm (for trans-nasal procedures). Follow the instruction displayed on the screen.
3. **Verify the Patient Tracker:** Position the tracking system so that it detects the Patient Reference. Once detected, patient reference detection indicator will turn to a green checkmark.
4. **Position and orient the tracking system:** Follow the instruction on the main UI and place the camera at correct distance from the Patient Reference. The tracking system must be facing to the Patient Reference. Follow any orientation adjustment instructions shown over the video image, if any.
5. **Verify Instruments tip (optional):** present the instruments to the tracking camera and verify the tip accuracy in the magnified window (displayed on top right corner of the camera video).


Click on **Registration** button or use the floating menu to proceed to the next stage.

5.6.2 Calibrating an instrument tip

In the Set Up stage a zoomed in image of the tip is displayed on the right sidebar panel. It provides visual feedback to verify the accuracy of the tip calibration. If the tip is not overlaid on the center of the hair cross then tip calibration is necessary.

Navient package includes a calibrator block that allows tip calibration. Simply place the tip of the instrument (or any surgical tool with a universal tracker clamped on it) on the dimple and present the calibrator and the instrument to the tracking camera. A progress bar will pop up showing the tip calibration progress. It typically takes a few seconds to complete the calibration. Once the tip calibration is complete verify the accuracy of the tip by looking at the tip zoomed in window.


5.7 Stage 4: REGISTRATION

To allow registration, you need to establish a frame of reference for the patient's head by attaching the Patient Tracker to the patient's forehead using a fresh dual-sided tape strip.

Navient registration process has two stages, landmark based registration followed by trace registration refinement. Both operations are designed to be conducted without touching the trackpad or keyboard. The Registrator pointer must be used for registration. Holding the Registrator pointer tip above the patient reference will activate the floating menu that enables switching to the next/previous stages. You could also reset the registration process using the floating menu (See section 5.7.3).

5.7.1 Landmark Based Registration

In this stage the user will be provided with visual and audio feedback to identify patient landmarks corresponding to the patient scan landmarks. Holding the Registrator pointer tip still on the patient landmarks will activate the sampling progress bar. It is very important to hold the Registrator tip as close as possible to the corresponding landmark shown on the 3D image.


Note 1: Landmarks registration using fiducials is similar to anatomical landmarks. Simply place the registrator pointer tip on the corresponding fiducial and wait for the progress bar to complete.

Note 2: Hide the registrator markers if you are not ready for registration. Place the registrator tip on the landmark location and then uncover registrator.

5.7.2 Trace Registration

Once the landmark registration is complete, the instruction for trace registration will be presented. Hold the Registrator tip next to the eye lateral right or ear tragus right (depending on which landmark was selected in the stage Mark Landmarks) and wait for the progress bar. Once progress bar is complete, the user needs to slowly glide the Registrator pointer around the patient forehead and nose area. It normally takes less than a minute to complete the process. Once it is done the application will proceed to the next stage.


Note: By default the software will go through trace based registration, the user can decide to skip it (by selecting **Navigation** circle in the floating menu or clicking on **Navigation** button in the side bar). A few reasons to skip trace based registration are:

- a) Patient scan doesn't cover the required tracing area
- b) Patient scan is low resolution or has large slice spacing
- c) Patient scan has significant motion artifacts
- d) Patient scan is very old and doesn't look like the real patient.

5.7.3 Floating menu

The user can interact with the software using gesture based commands. Simply hold the Registrator tip above the Patient Reference to enable the floating menu.

- **Activate Circle:** Using the Registrator move the cursor inside the **Activate** circle and wait for the progress bar. Once it is activated, you can proceed with selecting other circles.
- **Restart Landmarks:** This will restart Landmarks registration from beginning. Simply bring the cursor (Registrator tip) inside this circle and wait for the progress bar to complete
- **Restart Trace:** this will restart trace registration from beginning. Simply bring the cursor (Registrator tip) inside the circle and wait for the progress bar to complete.
- **Back Circle:** moves back to the **Set-Up** stage
- **Continue Circle:** this will be activated if the registration was completed before. If activated the software will proceed to the **Navigation** stage.

5.7.4 Using Cranial Patient Reference Frame

For cranial application, once registration is complete, remove the non-sterile Cranial Patient Reference and drape the patient.

Note the orientation of the Navient logo on the surface of the reference. When you re-attach the reference after sterilization, make sure the logo resides in the same orientation.

Match the mounting screw and the alignment post on the underside of the sterile reference frame with the corresponding screw hole and alignment notch on the attachment swivel of the Articulating Support Arm. Puncture the plastic drape with the mounting screw and secure the reference to the attachment swivel by tightening the mounting screw.


Check the navigation accuracy and make sure the draping process hasn't degraded the registration accuracy.

5.8 Stage 5: NAVIGATION

This stage starts with presenting accuracy check instructions. If the accuracy is not satisfactory then user can activate the floating menu to repeat the registration. If the accuracy is satisfactory then full screen navigation can be continued.

5.8.1 Accuracy Check

When entering the navigation stage, the distance between the tooltip to patient skin will be displayed. Use the Registrator and hold the instrument tip on various landmarks on the patients and check the distance from the skin. The navigation is accurate if the distance is less than 2mm.

Important Note: The registration must be repeated if the registration check shows an unacceptable registration error. Use the floating menu (or click on Registration button) to go back to the registration stage. Use only landmark based registration if the trace based registration fails. Landmark based registration should be used as the last option and should be followed by accuracy check. If the accuracy is still not as expected, consider terminating the use of Navient for guidance.

5.8.2 Navigation

In the Navigation stage, the user is presented with three viewports. By default, a 2D sagittal view is displayed in the larger viewport, and the axial and coronal viewports are displayed in two smaller viewports. The user can modify the viewports configuration in the settings. A cone representation of the pointer stem is shown when any trackable instrument is detected inside patient vicinity.


When a navigated tool and the reference marker are both being tracked, the sidebar control panel will be hidden to maximize the viewports size. The smart zoom feature will zoom up the region that the pointer tip is navigating. However if the tool or reference marker are not being tracked, the sidebar panel will be displayed along with the live video images. This would help user to identify why the navigation is inactive.

**Instrument tip location
and orientation**

If biopsy was planned in stage 2, then the software will display the biopsy planning tools as well as target navigation view on the side panel. It is allowed to modify the target and entry points for the biopsy.

If 3D model was generated in stage 2, then the software will display a toggle button (on-off) to enable or disable displaying of the 3D model on the patient scan. An Alpha blending slider allows controlling the opacity of the patient scan.

A tip extension toggle button (on-off button) allows virtual extension of the pointer tip. Use the slider to set the extension length. A warning message will be displayed on the screen if the extension is enabled. Extended tip will be displayed as red dashed line.


Important Note: Tip extension is primarily used during planning stage, like identifying the entry point. It is rarely used during navigation. Make sure the feature is disabled during the navigation.

5.9 Settings


Navient Settings form allows customizing and configuring Navient application settings. Following application properties are customizable:

- **Collage Layout:** Adjusts the navigation window layout. The default layout is set to **Coronal Axial Sagittal** with Coronal displayed in a large viewport and Axial and Sagittal displayed

in two smaller viewports. To change the layout, select the radio button next to the layout and click on "Save & Close".

- **Vertical splitter Fraction:** Sets the ratio of the width of the large window to the entire main window (x/y in the diagram below).

- **Auto Zoom:** Enables/disables automatic zooming feature. When the navigation instrument tip is still, it automatically zooms in to show magnified region surrounding the tip location.
- **Navigation volume margin:** a buffer zone around the patient scan. Navigation will be enabled if the navigation probe tip is inside the volume.
- **Restore Default:** Restores the default manufacturer settings.
- **Save & Close:** Stores all the changes permanently. It will restore the next time that application launches.
- **Cancel:** Ignores the changes and closes the settings window.

6 PACS Configuration

6.1 Introduction

This document explains how to configure a remote DICOM device, which is going to be used for querying and retrieving DICOM data as a result of a request from the Navient application, or for sending DICOM data to the Navient application with no request from the application.

DICOM files can be exchanged between two DICOM entities that are capable of receiving image and patient data in DICOM format. These DICOM devices must implement and support the DICOM classes/services needed for communication between devices over a network using TCP/IP protocol.

To be able to communicate between each other, each DICOM device needs to be aware about other device to establish a DICOM communication session. To achieve this condition an appropriate DICOM configuration must be done on both communicating devices. Each DICOM device has its own way for such configuration, therefore this document describes how to configure a remote DICOM device on Navient application side only.

6.2 Navient DICOM Entity configuration

Each DICOM device must have defined at least three following DICOM properties.

- Application Entity Title (AET)
- Host IP Address
- Port IP

When the Navient application is launched for the first time, the initial DICOM configuration file is created automatically. This file will contain the default values for two properties: AET and Port number, which can be replaced by a customer of Navient with their desired values. The Host IP Address of the PC where Navient is installed will be left empty and must be assigned by a customer.

To complete the Navient Local DICOM Entity configuration click on the button "Settings" on the Navient side panel, and then select "PACS" tab.

73

Navigation | Planning | Advanced | **PACS** ←

Local DICOM Application Entity

Application Entity Title:

Host IP Address: !

Port Number:

Remote DICOM Archives

No DICOM remote archive configured.
Use button "New" to add any available DICOM archive.

New Edit Delete Test

Restore Default Save & Close Cancel

And then define three following properties of Local DICOM Entity.

6.2.1 AET value

Default AET value, assigned by the application, is combined from the prefix "NAVI_" and the name of the computer where Navient is installed. If the computer name is longer than 11 characters, then it will be truncated to 11 first letters because the maximum allowed length of AET string is 16 characters.

Default AET is constructed in this way to make sure that each workstation with installed Navient will have unique AET value in case a customer will chose to leave the default AET for all workstations on their network.

Customers can replace the default AET value according to their policy for AET naming.

6.2.2 Host IP Address value

The IP address of the computer with Navient installed must be assigned by a customer before configuring a first remote DICOM device on that workstation.

6.2.3 Port Number value

Internet Assigned Numbers Authority (IANA) has reserved the following IP port numbers for DICOM devices.

- 104 well-known port for DICOM over Transmission Control Protocol (TCP). Since 104 is in the reserved subset, many operating systems require special privileges to use it.
- 2761 registered port for DICOM using Integrated Secure Communication Layer (ISCL) over TCP.
- 2762 registered port for DICOM using Transport Layer Security (TLS) over TCP.
- 11112 registered port for DICOM using standard, open communication over TCP.

Please note that DICOM standard recommends but does not require the use of these port numbers.

The port number 11112 has been selected as a default value for Navient, but the customers can replace this value according to their policy for assigning the port numbers to the DICOM devices.

IMPORTANT: If any changes are done on the form "Settings", then the button "Save and Close" must be used to close the form, otherwise the changes will be lost.

6.3 Remote DICOM device configuration

6.3.1 Add new DICOM remote device

Use the button “New” on the “PACS” tab to add a new remote DICOM device to the configuration.

Navigation Planning Advanced **PACS** Navigation stage co

Local DICOM Application Entity

Application Entity Title:

Host IP Address: !

Port Number:


Remote DICOM Archives

No DICOM remote archive configured.
Use button "New" to add any available DICOM archive.

Before adding the first remote DICOM device all three properties of Local DICOM Entity must be defined. Otherwise there will be the following error message popped up on clicking the button “New”.


After all properties of Local DICOM Entity are completed and the button "New" is clicked, the following dialog will be popped up.

"Displayed Name" can be any string and is used to show this DICOM device in the list of the remote DICOM devices configured in Navient application, which will be shown on the "DICOM Archives" tab after adding the first DICOM device.

Three other properties on that dialog are common DICOM device properties described in the section #2 “Navient DICOM Entity configuration” . Their values must be obtained from the configuration of the particular DICOM device to be added to the Navient configuration.

After entering the values for all properties the button “Add” should be used for adding this device or button “Cancel” to discard this dialog.

When the dialog is closed after clicking the button “Add” the “DICOM Archives” tab will have a modified GUI where the new data grid control with the DICOM devices list replaces two messages shown in previous tab’s GUI.

Navigation
Planning
Advanced
PACS

Local DICOM Application Entity

Application Entity Title: NAVI_AHMAD-PC

Host IP Address: 10.168.20.249

Port Number: 11112

Remote DICOM Archives

DICOM Archive Name	Ping Test Result	DICOM Echo Test Result
iSite		

☐ Enable DICOM data push from PACS - Only one time

Import Progress

New Edit Delete Test

Restore Default Save & Close Cancel

Now newly added DICOM device can be selected for editing, removing or testing DICOM connectivity with this device.

Navigation
Planning
Advanced
PACS

Local DICOM Application Entity

Application Entity Title: NAVI_AHMAD-PC

Host IP Address: 10.168.20.249

Port Number: 11112

Remote DICOM Archives

☒ DICOM Archive Name	Ping Test Result	DICOM Echo Test Result
☒ iSite		

☐ Enable DICOM data push from PACS - Only one time

Import Progress

New Edit Delete Test

Restore Default Save & Close Cancel

6.3.2 Edit existing remote DICOM device

The button "Edit" will be enabled if only one remote DICOM device is selected in case there is already more than one device configured.

On pressing the button "Edit" the same dialog as for adding a new device will be popped up. The only difference is that the button "Save" will replace the button "Add". The changes can be saved using this button or discarded using the button "Cancel".

30

Edit DICOM Archive

Displayed Name: iSite

Application Entity Title: STENTOR_SCU

Host IP Address: 10.168.10.26

Port Number: 107

Save Cancel

6.3.3 Delete existing remote DICOM device

The button "Delete" will be enabled if one or more devices is selected. All selected devices will be removed from configuration.

If all configured devices are removed, then the tab "PACS" GUI will be changed back to initial state with no device configured and showing two messages instead.

Navigation
Planning
Advanced
PACS
Navigation stage co

Local DICOM Application Entity

Application Entity Title: NAVI_AHMAD-PC

Host IP Address: Select one

Port Number: 11112

Remote DICOM Archives

No DICOM remote archive configured.
Use button "New" to add any available DICOM archive.

New
Edit
Delete
Test

Restore Default
Save & Close
Cancel

6.3.4 Test remote device

The button "Test" will be enabled if one or more devices is selected. All selected devices will be tested for DICOM connectivity.

Two following tests will be performed for each selected device.

- A regular network Ping command will be executed for the IP where the remote DICOM device is resided. So, if the IP of the DICOM device is correct and it is accessible on the network, then the test result should be successful.
- DICOM Echo request will be sent to the remote device. Each DICOM device must provide the service for that special DICOM request. So, if the DICOM device configuration is done correctly on both communicating sides, then the test result should be successful.

In case both tests are successful it will be shown in the following way.

Navigation | Planning | Advanced | **PACS**

Local DICOM Application Entity

Application Entity Title:

Host IP Address:

Port Number:

Remote DICOM Archives

☐ DICOM Archive Name	Ping Test Result	DICOM Echo Test Result
☒ iSite	Success	Success

☒ Enable DICOM data push from PACS - Only one time

If Ping test failed, then the results will be shown in the following way. Since the remote DICOM device is not reachable on the network, failure for the second test is of course expected.

Navigation
Planning
Advanced
PACS

Local DICOM Application Entity

Application Entity Title: NAVI_AHMAD-PC

Host IP Address: 10.168.20.249

Port Number: 11112

Remote DICOM Archives

☐	DICOM Archive Name	Ping Test Result	DICOM Echo Test Result
☐	iSite	Success	Success
☐	MDC	Success	Failed*

*See the log file for the details.

☐ Enable DICOM data push from PACS - Only one time

Import Progress

New Edit Delete Test

Restore Default Save & Close Cancel


Note: The reason for any test failure can be checked in the log file "Wif.log" under the Navient installation folder. Default installation folder is "C:\Navient

If the Ping test was successful but DICOM Echo test failed it means that while the remote DICOM device is reachable on the network, the DICOM configuration wasn't done correctly and must be verified on both communicating devices.


Note: If any changes are done on any tab of the form "Settings", then the button "Save and Close" must be used to close the form, otherwise the changes will be lost

6.4 Data Browser change after Remote DICOM device is configured


Before any remote DICOM device is configured in Navient application the "Data browser" window doesn't have "PACS" button, like in the screen below.

After adding the first remote DICOM device to the configuration, a new "PACS import" button will be automatically added to that context menu.

Clicking on “PACS Import” will open the list of the configured remote DICOM devices available for a selection.

Select the filter parameters and click on "OK". It will open a new window and display all datasets on the PACS.

Data Explorer

Patients

Free of all issues
23.3 GB

No of full studies in
C:\Users\user1\CLIN\RAW\Desktop\Navis21 Repository

Found 28 Patients

Patient ID	Patient Name	Sex	Age	D.O.B	Study Date	Description	Modality	Referring Physician
phantom	3000 Images	O	?	?	8/1/2005		CT	
1348370spine	3367 L-spine 1mm	M	?	?	6/20/2001		CT/KO	
23 CAROTIDES	AVA study2 Priors-Remote	O	98	9/28/1921	8/12/2002	RT	CT	DR TAPIAS
AvaFull 3	AvaFull Study Patient3	F	75	2/14/1944	3/21/2014	CT ANGIO ABDOMEN...	CT	BUTLER WILLIAM EDWA...
02100335	DEMO ABDOMEN-P	M	?	?	1/31/2000		CT	
02/01/1978	DEMO DENTACT	O	?	?	7/24/1996		CT	
02/01/1978	DEMO DENTACT Priors-Remote	M	?	?	7/24/1996		CT	
Dental Planning	Dental Routine	O	?	?	12/7/2000		CT	
LNA ICT	Follow-up 1 2Y Remote Prior	M	?	?	5/18/2009	La Pitie	CT	BR64
LNA ICT	Follow-up 1 2Y Remote Prior 2	M	?	?	5/18/2009	La Pitie	CT	BR64
LNA ICT	Follow-up 1 2Y Remote Prior 3	M	?	?	5/18/2009	La Pitie	CT	BR64
LNA ICT	Follow-up 1 2Y Remote Prior 4	M	?	?	5/18/2009	La Pitie	CT	BR64
1234567	IJ Dental Priors-Remote	M	78	9/8/1941	11/2/2007	CT Maxillofacial - Faci...	CT	WOWWHATANAbsound
FilterL_mAs200	LNA Phantom high dose Middle...	O	43	6/1/1976	11/11/2003		CT	Dr. Sus
Anonymized	Navient Fusion CT	M	?	?	1/1/2100		CT	
Pat7	NavisPine C10	M	?	?	2/22/2018		CT	

Scans

Series No	# Images	Modality	Date	Description	Matrix Size	Derived?
-----------	----------	----------	------	-------------	-------------	----------

7 Service and Support

If support is needed, contact us by:

Email: support@claronav.com

Website: <http://www.claronav.com/contact/>

Phone:

- (US & Canada): 1-647-951-1525
- Toll Free: 1 (844) 628-4336
- Contact your local representative

Fax: 1-647-951-1524

Regular Mail:

ClaroNav Kolahi Inc.
1140 Sheppard Avenue West, Unit 10
Toronto, Ontario, Canada, M3K 2A2

Authorized representative:

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Appendix A: The R-Fine Application

The tracking camera is specifically designed to preserve for a long time the high precision of its optical calibration performed at the factory. But its accuracy can suffer from undesired hardware changes for a variety of reasons, such as environmental stresses or rough handling. R-Fine provides a fast and intuitive method to evaluate the camera's calibration accuracy and the possibility to refine the calibration, if required. Navient application will remind the user to perform R-Fine every 12 months.


R-Fine operates by detecting a special AccBar marker (included in the Navient accessory kit). The AccBar must be placed by the user in a number of different regions and orientations, as guided by the software using semi-transparent image overlays. Once correctly positioned, the AccBar's coordinates are sampled accuracy analysis.

R-Fine Quick Guide


Upon opening the R-Fine application, the user is presented with a semi-transparent image of the AccBar tool. The user must position and orient the tool as shown on the screen.

- A *Yellow* guide is shown when the camera cannot detect the bar.
- A *Red* guide is shown when the bar is detected but not in the correct position.
- A *Green* guide is shown when the bar is positioned and oriented correctly.


Try to overlap AccBar with the guide


Sometimes it is required to move the bar towards or away from the camera

The bar must be oriented as depicted


In this orientation the bar is held pointing towards the camera

After all samples have been collected, the R-Fine application will determine whether the camera is out of alignment or not. The application will automatically update the calibration file if the

improvement is greater than 10%. In this scenario, the refinement process should be repeated until the improvement seen is less than 10%. If the improvement is greater than 10% after 3 consecutive attempts, please contact support. This is an indication that the camera is out of calibration

Improvement is greater than 10%. Repeat R-Fine procedure.

Improvement is less than 10%. Refinement is successful.

Improvement is greater than 10% for 3 consecutive tries. Contact support.

Appendix C: Troubleshooting Guide

The various error messages are described and explained here, as well as the possible causes and solutions.



Message	Explanation	Resolution
Another copy of Navient is running. Please wait for it to terminate, then restart	You are allowed to only execute a single instance of the Navient application at a time	Close the running application or wait for it to exit and then retry
Laptop power cable is disconnected. Make sure the Navient power cord is plugged in and both the laptop and the tracking system are powered on.	The laptop power cable is not connected or there is no power.	Make sure the laptop power cable is connected and Navient's main power cord is connected to a proper power outlet.
Camera may be out of calibration. Running R-Fine every 12 months is highly recommended	For optimal navigational accuracy the Navient camera needs to be calibrated every 12 months	Click "Yes" to run R-Fine or "No" to ignore (not recommended). See <i>Appendix A: The R-Fine Application</i>
Low resolution warning, spacing is more than 1.0mm. Navigation accuracy may be degraded	The spacing between voxels should be less than 1.0mm.	Use extra caution when navigating and perform accuracy checks by touching visible nearby structures prior to navigation.
Import failed, Invalid path	You may be trying to import data from the main Navient data repository (where it already resides), or from a directory that is not accessible, such as one on a network drive or portable medium that is not currently mapped.	Correct the problem with the path and try again.
Import failed. At least 100 slices are required	Navient is able to process datasets with 100 or more slices	Refer to Navient <i>CT Scanning Instructions</i>
Import failed. There is not enough disk space to import. Do you want to delete data from the repository	Disk is full	In order to continue existing datasets need to be deleted
Directory exists, overwrite patient data?	The data being imported already exists in the data repository	Select Yes to overwrite the existing data and plans file. Select No to cancel import.

Message	Explanation	Resolution
Error mapping drive	There is an issue with mapping the indicated drive	Please check the drive path and user credentials (if required)
Slice stack loading warning. More than one stack of slices found in the folder.	The dataset may be corrupted (e.g. missing or damaged files), or multiple scans are stored in the same folder.	The application will open the data and it is recommended that the user verify the patient data is complete by browsing through the slices
Tracking camera error. Please disconnect and reconnect the camera then wait a few seconds	The MT camera has been either disconnected or there is a communication error	Restart the application and if the application is still unable to connect to camera then restart the computer
Confirm delete, Permanently delete all data of the selected patient(s) from the repository?	Deleting a dataset from the Navient repository	Selecting "Yes" will remove the data from the repository
There is not enough disk space	The location that you are exporting the data to does not have enough disk space	Free some space and try again
Directory exists, Overwrite patient data?	The data being exported already exists in the destination folder	Selecting "Yes" will overwrite the destination folder
Adjust camera, move Right/Left/Up/Down/Closer/Further	In order to obtain optimal accuracy, Forehead Tag needs to be within the camera field of view	In case the markers are outside the bounds the application will warn the user to adjust the camera accordingly
Camera operating outside verified temperature range	The camera may be heating up	The cause for the rise in temperature should be investigated, once the issue is mitigated wait for the camera to cool down before continuing
Navigation inactive	During navigation if the forehead tag or the navigation pointer are not detected (for any reason e.g. line of sight interruption) the application will show the side bar with live camera video.	Check if something is obstructing the FOV. Make sure there is direct line of sight between the navigation pointers and the camera box.

Appendix D: Adjusting the cart arm's joint friction

When open, the cart's arm should allow for easy, yet stable, adjustment of the height of the front pole holding the laptop and camera box. The stability of the pole's height is controlled by friction in the joint connecting the arm to the back post. This friction may be adjusted by turning a screw located on the post's left side, which can be accessed through a hole in the plastic cover, as shown below.

Insert the supplied 1/8" Allen key through the hole and turn the screw clockwise to tighten it.

Appendix E: Electrical Safety and EMC Declarations

Manufacturer's Declaration for Navient (ME System) IEC 60601-1-2: 2014 (Ed 4.0)

Table 1 – Emission Class and Group Compliance

The Navient is intended for use in the electromagnetic environment specified below. The customer or the user of the Navient should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The Navient uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The Navient is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the Navient or shielding the location.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

Table 2 – Immunity Test Level Compliance

The Navient is intended for use in the electromagnetic environment specified below. The customer or the user of the Navient should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV Contact ±15 kV Air	±8 kV Contact ±15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-6	±0.5 kV, ±1 kV line to line & ±0.5 kV, ±1 kV, ±2 kV line to ground	±0.5 kV, ±1 kV line to line & ±0.5 kV, ±1 kV, ±2 kV line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T (100% dip in U_T) for 0,5 cycle 0% U_T (100% dip in U_T) for 1 cycles 70% U_T (30% dip in U_T) for 25 cycles 0% U_T (100% dip in U_T) for 5 sec	0% U_T (100% dip in U_T) for 0,5 cycle 0% U_T (100% dip in U_T) for 1 cycles 70% U_T (30% dip in U_T) for 25 cycles 0% U_T (100% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Navient requires continued operation during power mains interruptions, it is recommended that the Navient be powered from an uninterruptible power supply or a battery.
Power Frequency Magnetic Field (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Table 3 – Immunity Test Level Compliance

The Navient is intended for use in the electromagnetic environment specified below. The customer or the user of the Navient should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the Navient including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	6 Vrms ISM bands inside 150 kHz to 80 MHz	6 Vrms ISM bands inside 150 kHz to 80 MHz	
	3 V/m 80 MHz to 2,7 GHz	3 V/m 80 MHz to 2,7 GHz	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.7 GHz
	RF communication equipment inside 80 MHz to 6 GHz	RF communication equipment inside 80 MHz to 6 GHz	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Navient is used exceeds the applicable RF compliance level above, the Navient should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Navient

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4 – Recommended Separation Distances

Recommended separation distances between portable and mobile RF communications equipment and the NAVIENT			
The Navient is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Navient can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Navient as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.24
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Ahmad Kolahi, CEO

«_01_»_December_2023

Electronically signed & declared via online video before me in Toronto, ON while the Affiant / Declarant was located in Toronto on December 01, 2023, as per Ontario Reg. 431/20 & the Electronic Commerce Act.

N. Shum

SignNow e-signature ID: f488de7fc9...
12/01/2023 20:01:16 UTC

Narek Mirzoyan Licensed Paralegal & Notary Public Law Society of Ontario Licensee #P12490

A. Kolahi

SignNow e-signature ID: 62342ce89c...
12/01/2023 20:01:45 UTC

SignNow e-signature ID: 00de0b7e16...
12/01/2023 19:58:38 UTC

User Manual Appendix
Surgical navigation equipment Navient, versions
Otorhinolaryngological (95-NE-NE) and Universal neurosurgical
and otorhinolaryngological (95-NC-NC)

Manufacture: ClaroNav Kolahi Inc., Canada

Version 1.1

Appendix to the Operating Manual
Surgical navigation equipment NaviENT, Otorhinolaryngological (95-NE-NE) and Universal neurosurgical and otorhinolaryngological (95-NC-NC) versions
(hereinafter referred to as NaviENT equipment, equipment, product(s))

Before using this equipment, please read the instruction manual and this manual appendix carefully and retain it for future reference!

SAFETY INSTRUCTIONS

SAFETY INFORMATION:

- While properly used NaviENT equipment reduces the likelihood of surgical complications, the equipment is error-prone and can sometimes mislead you. NaviENT equipment requires regular verification of navigation accuracy;
- NaviENT equipment should not be used to make important surgical decisions without first verifying its correctness using a recognizable surface and landmarks;
- To minimize the risk of patient injury due to mislocalization, navigation accuracy should be checked at regular intervals;
- You can reduce the likelihood of navigation errors by being aware of the key factors that reduce accuracy and taking action to minimize or prevent them:
 - **Registration error:** Incorrect registration is the most common cause of localization errors. After the registration step, immediately perform an accuracy check and, if in doubt, repeat the registration;
 - **Tool tip calibration error:** NaviENT equipment tracks the position of markers (on each instrument) rather than the tips of the instruments. To match the location of the instrument tip to the CT volumetric image, the NaviENT equipment needs to know exactly the position of the instrument tip relative to the corresponding marker. For accurate operation it is essential that this position is confirmed and, if necessary, correctly calibrated;
 - **Moving the patient tracker:** Moving the patient tracker relative to the patient's head will invalidate the registration. When installing it, wipe the patient's forehead with an alcohol wipe and make sure the patient tracker is well attached to the skin and do not touch it during surgery. In case you suspect it may have moved, perform an accuracy check and re-register if necessary;
 - **Incorrect handling of the security camera unit:** The security camera unit is a precision component and should be handled with care. When handling the device body, avoid shaking, dropping or mishandling it. Do not open the case. In the event that you suspect that the security camera unit may not be operating properly (for example, after an accidental impact or when marker detection becomes unreliable), you should check the accuracy of the security camera unit by performing the R-Fine procedure. The R-Fine application is installed on the laptop and activated using the desktop icon.

PLEASE READ THE INSTRUCTION MANUAL AND THIS APPENDIX BEFORE USING THE EQUIPMENT

GENERAL WARNINGS:

- Do not install any other software on your laptop to prevent NaviENT hardware performance degradation.
- Do not connect any other equipment to the laptop other than the included accessories (Thunderbolt connection for security camera unit and DVD drive and USB drive with external memory).
- Although the laptop and security camera unit can operate on battery power alone for some time, the IR indicator on the security camera unit will remain off when running on battery power and therefore it is not recommended to operate on battery power alone. It is also possible that the battery will completely discharge in the middle of an operation, causing the NaviENT equipment to shut down.
- When NaviENT equipment is connected to the main power source, ensure that it is located in such a way that it can be quickly and easily disconnected from the mains if necessary.
- Ensure that all instruments (recorder, maxillary and anterior probes, bayonet pointer, calibrator, silicone tape, universal marker, clamps and large head screws) have been steam sterilized (autoclaving) according to the instructions in the instruction manual.
- The patient tracker must be disinfected before surgery according to the instructions in the instruction manual. Do not heat the patient tracker above 80°C (180°F) or attempt to steam sterilize it, as this may cause geometric distortion that will compromise the accuracy of the NaviENT equipment.
- None of the electronic components of NaviENT equipment are intended to be sterilized and may be damaged if sterilization is attempted. If contact between sterile glove surfaces and electronic components is required during surgery, use a sterile barrier such as a sterile, disposable nylon sheet.
- NaviENT equipment has been tested to comply with the requirements of the IEC 60601-1-2 standard. However, radio frequency interference may interfere with its operation or the operation of other nearby electrical devices. If you suspect one of these conditions, move the interfering equipment more than 30 cm (12 inches) from the NaviENT equipment, separate the equipment with a radio frequency barrier, or stop using the NaviENT equipment.
- Avoid contact of any liquids with the surface of the laptop or NaviENT equipment. Turn off the power; Turn off the laptop and allow the equipment to dry if you suspect that liquids may have entered any part of the equipment.
- Maintain within acceptable temperature and humidity ranges while operating NaviENT equipment.
- Do not drop or mishandle NaviENT equipment components in such a way as to cause physical damage. This may impair the performance of NaviENT equipment. Even if a component works immediately after being dropped, operational damage may have occurred that may not be seen until some time in the future.
- The laptop uses lithium-ion batteries. For optimal performance, NaviENT equipment should be connected to a power outlet during surgical use, in which case batteries will

not be used to power the laptop.

Do not use NaviENT equipment in any of the following locations:

- in places where the ambient temperature or humidity does not meet the requirements;
- in places where there may be strong shock or vibration;
- in places where dust accumulates.

When operating NaviENT equipment in the following locations, provide appropriate and sufficient protective measures:

- in areas exposed to static electricity or any other interference;
- in places exposed to intense electric or magnetic fields;
- in places located near power lines.

CONTENT

Introduction	7
1 Purpose and contraindications	8
2 General description of NaviENT surgical navigation equipment	10
3 Design features and appearance of NaviENT surgical navigation equipment	11
3.1 Surgical navigation equipment NaviENT, Otorhinolaryngological version (95-NE-NE)	11
3.2 Surgical navigation equipment NaviENT, version Universal neurosurgical and otorhinolaryngological (95-NC-NC)	23
4 Specifications	40
5 Materials used in manufacturing	46
6 Use of NaviENT surgical navigation equipment	49
6.1 Unpacking and installation	49
6.2 Preparing NaviENT surgical navigation equipment for operation	49
7 Rules for care, cleaning, disinfection and sterilization	51
8 Marking	54
9 Packaging	57
10 Safety requirements and environmental protection	58
11 Transportation and storage	59
12 Maintenance and repair	60
13 Software workflow	65
14 Disposal information	74
15 Manufacturer's warranty	75
16 Information about complaints	77
17 Electromagnetic compatibility	78
List of regulatory documentation used by the manufacturer to confirm compliance	82

INTRODUCTION

The main purpose of the manual and this appendix to the operating manual is to provide the user with instructions for working with NaviENT surgical navigation equipment, as well as information regarding its repair and maintenance. The manual provides technical characteristics and performance indicators of the NaviENT surgical navigation equipment, its installation, operation and maintenance, as well as safety information. Before use, the user must carefully read the manual and this appendix to the operating manual in order to correctly and accurately operate the NaviENT surgical navigation equipment, achieving the specified safety standards and performance indicators.

1 PURPOSE AND CONTRAINDICATIONS

Surgical navigation equipment NaviENT, Otorhinolaryngological (95-NE-NE) and Universal Neurosurgical and Otorhinolaryngological (95-NC-NC) versions, is intended for precise determination of the location of anatomical structures and dynamic representation of the location of the tip of the surgical instrument, matched with the corresponding location in the previously acquired CT scan of the patient's head, for open or percutaneous surgical procedures.

NaviENT equipment is suitable for professional use in operating rooms of hospitals and clinics by otolaryngologists and neurosurgeons.

Application area: in operating rooms of hospitals and clinics.

Indications for use:

Universal neurosurgical and otorhinolaryngological (95-NC-NC)	Otorhinolaryngological (95-NE-NE)
Cranial procedures: • resection of brain tumors; intracerebral: supratentorial (frontal, parietal and temporal lobes, subcortical nuclei, lateral ventricles, 3rd ventricle of the brain), subtentorial (cerebellum, brain stem, including the midbrain and medulla oblongata, pons, 4th ventricle); extracerebral: tumors of blood vessels, bones, membranes and neuromas of any location; • surgery on the base of the skull; • brain biopsy; • craniotomy; • puncture and drainage of the ventricles of the brain in adults and children; • implantation of a ventricular drainage system in adults and children; • thalamotomy. ENT procedures: • transsphenoidal procedures; • maxillary antrostomy; • ethmoidectomy; • sphenoidotomy; • resection of nasal turbinates; • frontal sinusotomy; • intralesional procedures.	• transsphenoidal procedures; • maxillary antrostomy; • ethmoidectomy; • sphenoidotomy; • resection of nasal turbinates; • frontal sinusotomy; • intralesional procedures; • intranasal tumor resections; • operations on the base of the skull.

Contraindications for use:

Medical conditions that contraindicate the use of Navient include any medical condition that may contraindicate the medical procedure itself.

Use of Navient should be discontinued if the accuracy check prior to tool guidance fails.

Possible complications:

The same complications as if the operation was performed without Navient. Navient does not add any new complications.

Possible side effects:

There are no known side effects from using the device. It is used as an aid in surgical procedures but does not diagnose disease or make treatment recommendations.

2 GENERAL DESCRIPTION OF NAVIGATION SURGICAL EQUIPMENT NAVIENT

NaviENT equipment, depending on the version, consists of a cart, a laptop with pre-installed software or a system unit with pre-installed software with a monitor, a security camera unit, a power cord and tools.

NaviENT equipment, depending on the design, functions and other characteristics, is manufactured in the following versions:

- Otorhinolaryngological (95-NE-NE);
- Universal neurosurgical and otorhinolaryngological (95-NC-NC).

3 DESIGN FEATURES AND APPEARANCE OF NAVIENT NAVIGATION SURGICAL EQUIPMENT

3.1 NAVIGATION SURGICAL EQUIPMENT NAVIENT, OTORINOLARYNGOLOGICAL OPTION (95-NE-NE)

DESIGN FEATURES AND APPEARANCE


Figure 1. Appearance

- Power cable – power cable;
- Thunderbolt adapter – adapter for connecting an Apple MacBook Pro laptop via a Thunderbolt connector to a security camera unit;
- DVD USB cable – DVD USB cable;
- Trackpad – trackpad (touchpad);
- DVD ROM Eject button – DVD ROM eject button;
- DVD ROM – DVD ROM;
- Camera LED indicator – LED indicator of the security camera unit;
- IR LED – infrared LEDs of the security camera unit;
- Tracking camera box – body of the tracking camera unit;
- Foldable arm – folding arm;
- Front wheels – front rollers;
- MacBook Pro laptop – Apple MacBook Pro laptop;
- Power button – power button;
- USB port – connector for USB;

Foldable arm locker – folding arm holder;
Navient cameras – navigation cameras;
Cart handle – handle of the cart of otorhinolaryngological equipment;
Foldable arm friction adjustment bolt – bolt for adjusting the folding arm;
Power cord wrap – power cord;
Side bins – side containers;
Rear wheels - rear rollers;
Stepper - stepper.

NaviENT equipment, Otorhinolaryngology version (95-NE-NE) includes an Apple MacBook Pro laptop with a charger for control optical navigation and planning “NaviENT”, a security camera unit with a power adapter, a power cord installed on the trolley of otorhinolaryngological equipment, as well as a set of tools for navigation during otorhinolaryngological interventions. The folding trolley arm can be folded to make the equipment compact and easy to move.

When moving an otolaryngology equipment cart, fold the folding arm into the rack and press the laptop tray against the folded arm (see Figure 2). Hold the otolaryngology equipment cart by the handle and roll it to a new location.


Figure 2. Appearance of otorhinolaryngological equipment trolleys

The main purpose of the otorhinolaryngology equipment cart is to hold a laptop and a security camera unit. The otorhinolaryngology equipment cart is lightweight, and the height of the folding arm can be adjusted. This allows the user to easily position NaviENT equipment in the operating room.

The equipment cart has a folding arm inside. To put the cart into working condition, pull the folding lever towards you; to move the cart back to the transport position, you need to press the folding lever from above.

The cart of otorhinolaryngological equipment is equipped with two stoppers for rollers. To prevent the cart from rolling downhill, place both casters in stoppers as shown in the illustration below. When moving equipment, remove the stoppers from the casters and place them in the side pockets of the equipment cart.

Figure 3. External view of the roller stopper

The security camera unit contains navigation cameras and infrared LEDs that operate in the near-infrared (NIR) light spectrum.


Figure 4. Block kUS tracking

The binding system identifies and captures checkered marker patterns. It is important to have a clear line of sight between the security camera unit and the markers on the instruments. There are two LED status indicators on the front panel of the security camera unit:

1. The green indicator shows the status of the navigation cameras:
 - dim (normal): cameras are not activated and do not broadcast images;
 - bright (normal): the camera is activated and broadcasts images;
 - Blinking (abnormal): Cameras have hardware or firmware problems. Contact Support.
2. The blue indicator shows the status of the infrared LEDs:
 - On (Normal): The infrared LEDs are on and functioning;
 - Off (abnormal): means the infrared LEDs are not turned on.

NaviENT equipment, Otorhinolaryngological version (95-NE-NE) has pre-installed NaviENT software (hereinafter referred to as NaviENT software).


Figure 5. Appearance of the Apple MacBook Pro laptop

NaviENT software runs on an Apple MacBook Pro 15.4-inch laptop. The laptop has the Windows 10 operating system installed. The live video stream of the security camera unit is received by the laptop via the Thunderbolt port and processed by NaviENT software.

NaviENT software displays the position of pre-calibrated instruments overlaid on the patient's preoperative scan.

NaviENT software is pre-installed on the laptop. A shortcut for the NaviENT Software application, the NaviENT Software Repository, and the NaviENT Software Record folders are placed on the desktop. Information about NaviENT software corresponds to Table 1.

Table 1: NaviENT Software Details

Software name	Software identification name	Software version number (identification number)	Digital software identifier (executable code checksum)	Algorithm for calculating the digital software identifier	Software release date
Navient	NES-R64-WIN	2.2.1	N/A	N/A	June 19, 2022

Surgical navigation equipment NaviENT, Otorhinolaryngological version (95-NE-NE) is equipped with a set of tools for navigation during otorhinolaryngological interventions.

All instruments must be steam sterilized before use, with the exception of the patient tracker.

All instruments are marked with black and white patterns (markers), which allow them to be tracked using security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the instruments.

A set of tools for navigation during otorhinolaryngological interventions includes:

- patient tracker (950-NE-PT);
- bayonet pointer (950-NE-BP);
- registrar (950-NE-RP);
- frontal seeker (950-NE-FS);
- maxillary seeker (950-NE-MS);
- universal tag (950-NE-UT1);
- universal tag 360 (950-NE-UT2);
- universal C-clamp small (950-NE-UTCS);
- universal C-clamp medium (950-NE-UTCM);
- universal C-shaped clamp large (c-clamp large) (950-NE-UTCL);
- universal clamp C-shaped X large (c-clamp X large) (950-NE-UTCXL);
- calibrator (950-NE-TC);
- tray for sterilizing instruments for navigation during otorhinolaryngological interventions (950-NE-ST).*


Figure 6. Appearance of tools for navigation during otorhinolaryngological interventions

***Note: all instruments, being reusable, are packaged inside a tray for sterilization upon delivery**

Also, NaviENT equipment, Otorhinolaryngological version (95-NE-NE), can additionally be equipped with a second set of instruments (identical to the first) for navigation during otorhinolaryngological interventions, each instrument of which is marked with the letter "B". Instruments for navigation during otorhinolaryngological interventions with the marking "B" allow you to use NaviENT equipment without waiting for sterilization of the main set. The additional set of tools for navigation during otorhinolaryngological interventions includes the same tools as the main set.


Figure 7. Appearance of tools for navigation during otorhinolaryngological interventions additional set

The recorder is made of aluminum, with black and white patterns (markers) applied, which allow it to be detected by security cameras. During recording, the otolaryngologist or surgeon touches the patient's landmarks with the tip of the recorder and performs a trace registration by sliding the tip of the recorder across the patient's forehead.


Figure 8. Appearance registrar

The maxillary probe is a common instrument for reaching the maxillary sinus area. The maxillary probe includes an ergonomic aluminum attachment with black and white patterns (markers) that allows it to be detected by security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the instrument. Once detected, NaviENT software will display the real-time location of the maxillary probe tip overlaid on the patient's CT scan data.

Figure 9. Appearance maxillary probe

The frontal probe is a common instrument for reaching the frontal sinus area. The front probe includes an ergonomic aluminum attachment with black and white patterns (markers) that allows it to be detected by security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the instrument. After software detection NaviENT will display the real-time location of the frontal probe tip overlaid on the patient's CT scan data.


Figure 10. Appearance of a frontal probe

Bayonet pointer is a bayonet-shaped pointer with a long stainless steel tip that allows the otolaryngologist or surgeon to localize internal anatomical changes. Black and white patterns (markers) allow it to be detected by security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the signpost. After discovery BY NaviENT will display the real-time location of the pointer tip overlaid on the patient's CT scan data.

Figure 11. Appearance of a bayonet pointer

The calibrator is designed to calibrate instrument tips to make them traceable. Black and white patterns (markers) on the calibrator allow it to be tracked using security cameras. To calibrate the instrument tip, the tip must be placed in a small recess on the calibrator, and then both the instrument and calibrator must be presented to the camera. It is important to have a clear line of sight between security cameras and markers. Once the instrument and calibrator are detected, NaviENT software will display a progress bar showing the tip calibration progress.

Figure 12. Calibrator appearance

Dimple to place the instrument tip - recess for placing the tip of the instrument
Calibrator handle – calibrator handle

The patient tracker is made of plastic and is attached to the patient's forehead with double-sided silicone tape. It remains on the patient's forehead throughout the operation. Black and white patterns (markers) allow you to track it using security cameras. It is important to have a clear line of sight between the security camera unit and the black and white patterns on the tracker. Its main purpose is to determine the spatial location of the patient's head.


Figure 13. Patient tracker appearance

Patient tracker marker – patient tracker marker
Double sided tape – double-sided silicone tape
Direction indicator icon – direction indicator icon

The universal marker is designed to be attached using a universal C-shaped clamp to endoscope instruments and allows them to be tracked using security cameras. He has no contact with the patient.

Once attached, the user will need to go through a tool tip calibration process. Once calibrated, NaviENT software will display the real-time location of the instrument tip overlaid on the patient's CT scan data.


Figure 14. Appearance of a universal marker

Markers – markers (black and white patterns)

Thumb Screw – screw with a large head

C-Clamp Adaptor – universal C-shaped clamp

The 360 Universal Marker is designed to be attached using a Universal C-clamp to endoscope instruments and allows them to be tracked using security cameras. He has no contact with the patient. Once attached, the user will need to go through a tool tip calibration process. Once calibrated, NaviENT software will display the real-time location of the instrument tip overlaid on the patient's CT scan data.


Figure 15. Appearance of the 360 universal marker

Universal C-clamps include four C-clamps of different sizes. They are designed for attaching universal markers to any surgical instruments.

Figure 16. Appearance of universal C-shaped clamps

C-Clamp Small – universal C-shaped small clamp

C-Clamp Medium – universal C-shaped medium clamp

C-Clamp Large – universal C-shaped large clamp

C-Clamp X Large – universal clamp C-shaped X large

Universal C-shaped clamps are used for attachment to any surgical instruments with a diameter from 4 mm to 25 mm. Each universal clamp has two large head screws:

- a small, large-head screw for attaching the multi-purpose marker to the clip;
- a long screw with a large head for attaching the clamp to a surgical instrument.

Universal C-shaped clamps are made in 4 different sizes:

- small: for instruments with a diameter of 4-10 mm;
- medium: for instruments with a diameter of 10-15 mm;
- large: for instruments with a diameter of 15-20 mm;
- X large: for instruments with a diameter of 20-25 mm.

To attach the patient tracker, NaviENT equipment is equipped with a silicone tape. This is a self-adhesive, durable tape. It can also be used to attach the Universal Marker and the 360 Universal Marker to a surgical instrument. Before use, the silicone tape is cut into 15-20 cm lengths and steam sterilized individually.

NaviENT equipment is also equipped with an ENT phantom.


Figure 16.1. ENT phantom

Principle of operation

The navigation function in the NaviENT Surgical Navigation Equipment, Otorhinolaryngological version (95-NE-NE) is based on images of the patient obtained before the procedure, combined with optical measurements of the position of the instrument tips relative to the patient's head. To enable navigation, the patient tracker is attached to the patient's forehead to provide tracking of the patient's head. Images of the patient are then registered to the position of the patient's head in space by matching the locations of markers marked both on the image and on the patient's head, followed by matching the path traced by the user on the patient's face to a model of the patient's facial surface automatically generated from image data. Subsequently, NaviENT Surgical Navigation Equipment, Otorhinolaryngological version (95-NE-NE) represents the location of the tip of the required surgical instrument, displaying it on these images using a tracking camera unit and NaviENT software.


Figure 17. Attaching the patient tracker to the patient's forehead

Figure 17.1. Storage case for ENT phantom

EQUIPMENT:

1) Surgical navigation equipment NaviENT, Otorhinolaryngological version (95-NE-NE) assembled:

1.1) Trolley for otorhinolaryngological equipment – 1 piece;

1.2) Security camera unit – 1 piece;

1.3) Laptop for controlling optical navigation and planning with pre-installed NaviENT software – 1 pc.;

1.4) Laptop charger with extended cable – 1 piece;

1.5) Power cord – 1 piece;

1.6) AccBar camera calibrator – 1 piece;

1.7) Stopper for the roller – 2 pcs;

2) Cover for the security camera unit – 1 piece;

3) A set of instruments for navigation during otorhinolaryngological interventions, up to 10 pieces (if necessary), consisting of:

patient tracker – up to 10 pcs (if necessary);

bayonet pointer – up to 10 pcs (if necessary);

recorder – up to 10 pcs (if necessary);

front probe – up to 10 pcs (if necessary);

maxillary probe – up to 10 pcs (if necessary);

universal marker – up to 10 pcs (if necessary);

universal marker 360 – up to 10 pcs (if necessary);

universal C-shaped small clamp – up to 10 pcs (if necessary);

universal C-shaped medium clamp – up to 10 pcs (if necessary);

universal C-shaped large clamp – up to 10 pcs (if necessary);

universal C-shaped clamp X large – up to 10 pcs (if necessary);

calibrator – up to 10 pcs (if necessary);

tray for sterilizing instruments for navigation during otorhinolaryngological interventions - up to 10 pcs (if necessary).

4) Silicone tape – up to 10 pcs (if necessary);

5) ENT phantom – up to 10 pieces (if necessary) for training;

6) Operation manual – 1 pc.

3.2 NAVIGATION SURGICAL EQUIPMENT NAVIENT, UNIVERSAL NEUROSURGICAL AND ENTORHOLARINGOLOGICAL OPTION (95-NC-NC)

DESIGN FEATURES AND APPEARANCE


Figure 18. Appearance

Tracking Camera Arm – lever of the tracking camera unit

Extension pole - extension lever

Cart handle – cart handle for universal neurosurgical and otorhinolaryngological equipment

Tracking camera – block of tracking cameras

Touchscreen Monitor – 27-inch touchscreen monitor

Extension pole lock / release pedal – pedal for locking/unlocking the extension lever

Caster lock - roller lock

NaviENT equipment, version Universal neurosurgical and otorhinolaryngological (95-NC-NC) includes a system unit, a 27-inch touch monitor, a power cord with a power adapter and a security camera unit installed on the trolley of universal neurosurgical and otorhinolaryngological equipment, as well as sets tools for navigation during otorhinolaryngological interventions, neurosurgical interventions and biopsies.

7-2

The height of the multipurpose neurosurgical and otolaryngology equipment cart can be adjusted by pressing the lock/unlock pedal on the extension arm between the two front casters (see Figure 19). This allows the user to easily position NaviENT equipment in the operating room.


Figure 19. Appearance of navigational surgical equipment NaviENT, version Universal neurosurgical and otorhinolaryngological (95-NC-NC) in working order and compact for moving

Figure 19.1. Appearance power adapter

When moving the multipurpose neurosurgical and otolaryngology equipment cart, lower the extension arm, secure the security camera unit, and rotate the touch monitor to a 90-degree angle. (See Figure 19 on the right). Hold the cart by the handle and roll it to a new location.

The main purpose of the universal neurosurgical and otorhinolaryngological equipment trolley is to hold the system unit, touch monitor and security camera unit.

The trolley of universal neurosurgical and otorhinolaryngological equipment is equipped with rollers with locks. To prevent the cart from rolling downhill, secure the casters by pushing the caster lock down.

The security camera unit contains navigation cameras and infrared LEDs that operate in the near-infrared (NIR) light spectrum.


Figure 20. Block kUS tracking

The binding system identifies and captures checkered marker patterns. It is important to have a clear line of sight between the security camera unit and the markers on the instruments. There are two LED status indicators on the front panel of the security camera unit:

1. The green indicator shows the status of the navigation cameras:
 - dim (normal): cameras are not activated and do not broadcast images;
 - bright (normal): the camera is activated and broadcasts images;
 - Blinking (abnormal): Cameras have hardware or firmware problems. Contact Support.
2. The blue indicator shows the status of the infrared LEDs:
 - On (Normal): The infrared LEDs are on and functioning;
 - Off (abnormal): means the infrared LEDs are not turned on.

Surgical navigation equipment NaviENT, version Universal neurosurgical and otorhinolaryngological (95-NC-NC) has pre-installed NaviENT software (hereinafter referred to as NaviENT software).

NaviENT software runs on an Intel NUC8I7BEH system unit with an Acer T272HUL 27-inch touchscreen monitor. The Windows 10 operating system is installed on the system unit. The live video stream of security cameras is received by the system unit and processed by NaviENT software.


Figure 21. External view of the system unit


Figure 22. Appearance of the touch monitor

NaviENT software displays the position of pre-calibrated instruments overlaid on the patient's preoperative scan.

NaviENT software is pre-installed on the system unit. A shortcut for the NaviENT Software application, the NaviENT Software Repository, and the NaviENT Software Record folders are placed on the desktop. Information about NaviENT software corresponds to Table 2.

Table 2: NaviENT Software Details

Software name	Software identification name	Software version number (identification number)	Digital software identifier (executable code checksum)	Algorithm for calculating the digital software identifier	Software release date
Navient	NES-R64-WIN	2.2.1	N/A	N/A	June 19, 2022

Surgical navigation equipment NaviENT, version Universal neurosurgical and otorhinolaryngological (95-NC-NC) is equipped with sets of tools for navigation presented in Table 3.

Table 3: Navigation Toolboxes**

A set of tools for navigation during otorhinolaryngological interventions	A set of tools for navigation during neurosurgical (cranial) interventions	Biopsy Navigation Toolkit
1. Patient tracker 2. Maxillary probe 3. Front probe 4. Bayonet pointer 5. Registrar 6. Calibrator 7. Universal marker 8. Universal marker 360 9. Universal clamp C-shaped small 10. Universal clamp C-shaped medium 11. Large universal C-shaped clamp 12. Universal clamp C-shaped X large 13. Tray for sterilization of instruments for navigation during otorhinolaryngological interventions	1. Articulated arm 2. Two cranial reference frames 3. Straight probe 4. Registrar 5. Calibrator 6. Universal marker 360 7. Universal clamp C-shaped small 8. Universal clamp C-shaped medium 9. Large universal C-shaped clamp 10. Universal clamp C-shaped X large 11. Tray for sterilization of instruments for navigation during neurosurgical (cranial) interventions	1. Articulated arm 2. Double adapter 3. Biopsy targeting device 4. Biopsy probe 5. Universal marker 360 6. Reducing tube 2.2 mm 7. Tray for sterilizing instruments for biopsy navigation


Figure 23. Appearance of tools for navigation during otorhinolaryngological interventions
 **Note: all instruments, being reusable, are packaged inside trays for sterilization upon delivery.


Figure 24. Appearance of tools for navigation during neurosurgical (cranial) interventions


Figure 25. Appearance of tools for biopsy navigation

All instruments must be steam sterilized before use, with the exception of the patient tracker.

All instruments are marked with black and white patterns (markers), which allow them to be tracked using security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the instruments.

Also, NaviENT equipment, version Universal Neurosurgical and Otorhinolaryngological (95-NC-NC), can additionally be equipped with a second set (identical to the first) of instruments for navigation during otorhinolaryngological interventions, each instrument of which is marked with the letter "B". Instruments for navigation during otorhinolaryngological interventions with the marking "B" allow you to use NaviENT equipment without waiting for sterilization of the main set. The additional set of tools for navigation during otorhinolaryngological interventions includes the same tools as the main set.


Figure 26. Appearance of tools for navigation during otorhinolaryngological interventions additional set

The recorder is made of aluminum, with black and white patterns (markers) applied, which allow it to be detected by security cameras. During recording, the otolaryngologist or neurosurgeon touches patient landmarks with the tip of the recorder and performs a trace registration by sliding the tip of the recorder across the patient's forehead.


Figure 27. Appearance registrar

The maxillary probe is a common instrument for reaching the maxillary sinus area. The maxillary probe includes an ergonomic aluminum attachment with black and white patterns (markers) that allows it to be detected by security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the instrument. Once detected, NaviENT software will display the real-time location of the maxillary probe tip overlaid on the patient's CT scan data.

Figure 28. Appearance maxillary probe

The frontal probe is a common instrument for reaching the frontal sinus area. The front probe includes an ergonomic aluminum attachment with black and white patterns (markers) that allows it to be detected by security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the instrument. After software detection NaviENT will display the real-time location of the frontal probe tip overlaid on the patient's CT scan data.


Figure 29. Appearance of a frontal probe

Bayonet pointer is a bayonet-shaped pointer with a long stainless steel tip that allows the otolaryngologist or surgeon to localize internal anatomical changes. Black and white patterns (markers) allow it to be detected by security cameras. It is important to have a clear line of sight between the security camera unit and the markers on the signpost. After discovery BY NaviENT will display the real-time location of the pointer tip overlaid on the patient's CT scan data.


Figure 30. Appearance of a bayonet pointer

The calibrator is designed to calibrate instrument tips to make them traceable. Black and white patterns (markers) on the calibrator allow it to be tracked using security cameras. To calibrate the instrument tip, the tip must be placed in a small recess on the calibrator, and then both the instrument and calibrator must be presented to the camera. It is important to have a clear line of sight between security cameras and markers. Once the instrument and calibrator are detected, NaviENT software will display a progress bar showing the tip calibration progress.

Figure 31. Calibrator appearance

Dimple to place the instrument tip - recess for placing the tip of the instrument
Calibrator handle – calibrator handle

The patient tracker is made of plastic and is attached to the patient's forehead with double-sided silicone tape. It remains on the patient's forehead throughout the operation. Black and white patterns (markers) allow you to track it using security cameras. It is important to have a clear line of sight between the security camera unit and the black and white patterns on the tracker. Its main purpose is to determine the spatial location of the patient's head.


Figure 32. Patient tracker appearance

Patient tracker marker – patient tracker marker
Double sided tape – double-sided silicone tape
Direction indicator icon – direction indicator icon

The universal marker is designed to be attached using a universal C-shaped clamp to endoscope instruments and allows them to be tracked using security cameras. He has no contact with the patient. Once attached, the user will need to go through a tool tip calibration process. Once calibrated, NaviENT software will display the real-time location of the instrument tip overlaid on the patient's CT scan data.


Figure 33. Appearance of a universal marker

Markers – markers (black and white patterns)

Thumb Screw – screw with a large head

C-Clamp Adapter – universal clamp

The 360 Universal Marker is designed to be attached using a Universal C-clamp to endoscope instruments and allows them to be tracked using security cameras. He has no contact with the patient. Once attached, the user will need to go through a tool tip calibration process. Once calibrated, NaviENT software will display the real-time location of the instrument tip overlaid on the patient's CT scan data.


Figure 34. Appearance of the 360 universal marker

Universal C-clamps include four C-clamps of different sizes. They are designed for attaching universal markers to any surgical instruments.

Figure 35. Appearance of universal clamps

C-Clamp Small – universal C-shaped small clamp

C-Clamp Medium – universal C-shaped medium clamp

C-Clamp Large – universal C-shaped large clamp

C-Clamp X Large – universal clamp C-shaped X large

Universal C-shaped clamps are used for attachment to any surgical instruments with a diameter from 4 mm to 25 mm. Each universal C-shaped clamp has two large-head screws:

- a small, large-head screw for attaching the multi-purpose marker to the clip;
- a long screw with a large head for attaching the clamp to a surgical instrument.

Universal C-shaped clamps are made in 4 different sizes:

- small: for instruments with a diameter of 4-10 mm;
- medium: for instruments with a diameter of 10-15 mm;
- large: for instruments with a diameter of 15-20 mm;
- X large: for instruments with a diameter of 20-25 mm.

The articulated arm is designed to attach the patient tracker to the patient's table. It is made of stainless steel and aluminum and is steam sterilized.

Figure 36. External view of the articulated arm

The dual adapter connects to the third-party Mayfield Skull Stabilization System and has two ports for attaching two articulated arms, one for the cranial reference frame and the other for the biopsy targeting device. It is made of aluminum and is steam sterilized.

Figure 37. Appearance of the double adapter

The cranial reference frame is made of aluminum and is steam sterilized.

Figure 38. Appearance of the cranial reference frame

The biopsy targeting device is used to guide the biopsy needle to the target point. It is made of aluminum and is sterilized with steam.

Figure 39. Appearance of the biopsy targeting device

The 2.2mm reduction tube is made of stainless steel and is steam sterilized.

Figure 40. Appearance of a 2.2 mm reducing tube

The biopsy probe is a short-tipped instrument made of aluminum. Designed to assist in setting up the biopsy targeting device. Supplied sterilized. Sterilization method is steam.

Figure 41. Appearance of the biopsy probe

The straight probe is made of aluminum and is steam sterilized.

Figure 42. External view of the straight probe

To attach the patient tracker, NaviENT equipment is equipped with a silicone tape. This is a self-adhesive, durable tape. It can also be used to attach the Universal Marker and the 360 Universal Marker to a surgical instrument. Before use, the silicone tape is cut into 15-20 cm lengths and steam sterilized individually.

NaviENT equipment is also equipped with an ENT phantom, a phantom stand and a neurosurgical phantom.

Figure 43. Neurosurgical phantom with phantom stand

(Image of ENT phantom Figure 16.1)


Figure 43.1. Case for storing a neurosurgical phantom

Principle of operation

The navigation function in the NaviENT surgical navigation equipment, version Universal neurosurgical and otorhinolaryngological (95-NC-NC) is based on images of the patient obtained before the procedure, in combination with optical measurements of the position of the instrument tips relative to the patient's head. To enable navigation, the patient tracker is attached to a skull holder to provide tracking of the patient's head. Images of the patient are then registered to the position of the patient's head in space by matching the locations of markers marked both on the image and on the patient's head, followed by matching the path traced by the user on the patient's face to a model of the patient's facial surface automatically generated from image data. Subsequently, the NaviENT surgical navigation equipment, version Universal neurosurgical and otorhinolaryngological (95-NC-NC), represents the location of the tip of the required surgical instrument, displaying it on these images using a tracking camera unit and NaviENT software.

EQUIPMENT:

1) Surgical navigation equipment NaviENT, version Universal neurosurgical and otorhinolaryngological (95-NC-NC) assembled:

- 1.1) Trolley for universal neurosurgical and otorhinolaryngological equipment – 1 piece;
- 1.2) Security camera unit – 1 piece;
- 1.3) System unit with pre-installed NaviENT software – 1 piece;
- 1.4) Touch monitor for control optical navigation and planning – 1 piece;
- 1.5) Power cord – 1 piece;
- 1.6) Power adapter – 1 piece;
- 1.7) AccBar camera calibrator – 1 piece;

2) Case for touch monitor – 1 pc;

3) Cover for the security camera unit – 1 piece;

4) A set of instruments for navigation during otorhinolaryngological interventions, up to 10 pieces (if necessary), consisting of:

- patient tracker – up to 10 pcs (if necessary);
- bayonet pointer – up to 10 pcs (if necessary);
- recorder – up to 10 pcs (if necessary);

- front probe – up to 10 pcs (if necessary);
- maxillary probe – up to 10 pcs (if necessary);
- universal marker – up to 10 pcs (if necessary);
- universal marker 360 – up to 10 pcs (if necessary);
- universal C-shaped small clamp – up to 10 pcs (if necessary);
- universal C-shaped medium clamp – up to 10 pcs (if necessary);
- universal C-shaped large clamp – up to 10 pcs (if necessary);
- universal C-shaped clamp X large – up to 10 pcs (if necessary);
- calibrator – up to 10 pcs (if necessary);
- tray for sterilizing instruments for navigation during otorhinolaryngological interventions - up to 10 pcs (if necessary).

5) A set of instruments for navigation during neurosurgical (cranial) interventions, up to 10 pieces (if necessary), consisting of:

- articulated lever – up to 10 pcs (if necessary);
- cranial reference frame – up to 20 pcs (if necessary);
- direct probe – up to 10 pcs (if necessary);
- recorder – up to 10 pcs (if necessary);
- calibrator – up to 10 pcs (if necessary);
- universal marker 360 – up to 10 pcs (if necessary);
- universal C-shaped small clamp – up to 10 pcs (if necessary);
- universal C-shaped medium clamp – up to 10 pcs (if necessary);
- universal C-shaped large clamp – up to 10 pcs (if necessary);
- universal C-shaped clamp X large – up to 10 pcs (if necessary);
- tray for sterilizing instruments for navigation during neurosurgical (cranial) interventions – up to 10 pcs (if necessary).

6) A set of tools for navigation during biopsy, up to 10 pieces (if necessary), consisting of:

- articulated lever – up to 10 pcs (if necessary);
- double adapter – up to 10 pcs (if necessary);
- targeting device for biopsy – up to 10 pcs (if necessary);
- biopsy probe – up to 10 pcs (if necessary);
- universal marker 360 – up to 10 pcs (if necessary);
- reducing tube 2.2 mm – up to 10 pcs (if necessary);
- tray for sterilizing instruments for navigation during biopsy – up to 10 pcs (if necessary).

7) Silicone tape – up to 10 pcs (if necessary);

8) ENT phantom – up to 10 pieces (if necessary) for training;

9) Phantom stand – up to 10 pcs (if necessary) for training;

10) Neurosurgical phantom – up to 10 pieces (if necessary) for training;

11) Operation manual – 1 pc.

4 TECHNICAL SPECIFICATIONS

The operational and technical characteristics of NaviENT equipment must comply with Tables 4-5.

Table 4

Parameter name	Meaning, norm
1	2

Overall dimensions of NaviENT equipment in a compact (folded state of moving parts) state for movement (H×L×W), mm • Otorhinolaryngological (95-NE-NE) • Universal neurosurgical and otorhinolaryngological (95-NC-NC)	1870×440×410 1200×500×400
Overall dimensions of NaviENT equipment in working condition (with moving parts deployed) (H×D×W), mm • Otorhinolaryngological (95-NE-NE) • Universal neurosurgical and otorhinolaryngological (95-NC-NC)	1020×440×410 1900×500×400
Dimensions of NaviENT equipment trolley (L×W), mm • Otorhinolaryngological (95-NE-NE) • Universal neurosurgical and otorhinolaryngological (95-NC-NC)	440×410 500×400
Number of trolley rollers for otorhinolaryngological equipment (95-NE-NE), pcs. • front • rear	2 2
External dimensions of the stopper for rollers of otorhinolaryngological equipment (95-NE-NE) (H×D×W), mm	34×130×58
Internal dimensions of the stopper for rollers of otorhinolaryngological equipment (95-NE-NE) (L×W), mm	66.25×21
Dimensions of trolley rollers of universal neurosurgical and otorhinolaryngological equipment (95-NC-NC), mm Diameter Thickness	125 32
Number of trolley rollers for universal neurosurgical and otorhinolaryngological equipment (95-NC-NC), pcs.	4
Dimensions of the folding trolley arm of otorhinolaryngological equipment (95-NE-NE), mm	635×510
Dimensions of the trolley extension arm of the universal neurosurgical and otorhinolaryngological equipment (95-NC-NC), mm	1000
Dimensions of the lever of the security camera unit of the trolley of the universal neurosurgical and otorhinolaryngological equipment (95-NC-NC), mm Shoulder 1 Shoulder 2	290 370
Time of readiness for work after switching on, min., no more • Otorhinolaryngological (95-NE-NE) • Universal neurosurgical and otorhinolaryngological (95-NC-NC)	5 5
Operating mode	

- Otorhinolaryngological (95-NE-NE) - Universal neurosurgical and otorhinolaryngological (95-NC-NC)	Long lasting Long lasting
Power cord length, mm - Otorhinolaryngological (95-NE-NE) - Universal neurosurgical and otorhinolaryngological (95-NC-NC)	2400 2500
Power cord weight, kg - Otorhinolaryngological (95-NE-NE) - Universal neurosurgical and otorhinolaryngological (95-NC-NC)	0.35 0.56
Rated supply voltage, V - Otorhinolaryngological (95-NE-NE) - Universal neurosurgical and otorhinolaryngological (95-NC-NC)	220-240 100-240
Current consumption, A, no more - Otorhinolaryngological (95-NE-NE) - Universal neurosurgical and otorhinolaryngological (95-NC-NC)	0.6-1.0 2.2
Rated value of supply frequency, Hz - Otorhinolaryngological (95-NE-NE) - Universal neurosurgical and otorhinolaryngological (95-NC-NC)	50-60 50-60
Weight of NaviENT equipment without packaging/with packaging, kg - Otorhinolaryngological (95-NE-NE) - Universal neurosurgical and otorhinolaryngological (95-NC-NC)	30/45 30/33
<i>Characteristics of the security camera unit:</i>	
Overall dimensions of the security camera unit housing (H×D×W), mm	140×280×40
LED Shape	domed
LED size, mm	4.95×8.5
LED IR wavelength, nm	850
Angle of diffusion of IR radiation by LED, degrees.	60
Rated DC voltage of the LED, V	3.1
Number of security cameras, pcs.	3
Security camera lens diameter, mm	20
Rated voltage of security camera, V	12
Security camera power consumption, VA	4
Power adapter model for 95-NE-NE	CUI Inc. model SW124-12-N-P5
Power adapter input specifications: - rated supply voltage, V - rated value of supply frequency, Hz - supply current, A	100-240 50-60 2.2
Power adapter output characteristics:	

- power no more than, W - DC voltage, V	24 12
Length of power cable for security camera unit, mm	120
Case dimensions(L×W), mm	330*245
Case weight, kg	0.3
<i>Management laptop specifications optical navigation and planning "NaviENT" for equipmentOtorhinolaryngological (95-NE-NE):</i>	
Laptop brand	Apple
Laptop model	MacBook Pro
Overall dimensions (H*L×W), mm	180×247.1×358.9
Weight, kg	2.04
Screen diagonal no more than, inch	15.4
Charger with extended cable	MagSafe 2
Extended cable length, mm	1500
<i>Characteristics of the touch monitor for Universal neurosurgical and otorhinolaryngological equipment (95-NC-NC):</i>	
Monitor brand	Acer
Monitor model	T272HUL
Type	color TFT LCD (Touch Screen)
Overall dimensions (L×W), mm	660*410
Display diagonal size, inch	27
Display screen resolution, pixels per inch	(WQHD)2560 x 1440@60 Hz
Nutrition	from power adapter
Case dimensions(L×W), mm	680*470
Case weight, kg	1.2
<i>Characteristics of the system unit for the equipmentuniversal neurosurgical and otorhinolaryngological (95-NC-NC):</i>	
Brand	Intel
Model	NUC8I7BEH
Overall dimensions (L×W), mm	115*115
Monitor connection cable length, mm	325
Rated DC voltage, V	19
Power consumption, W	90
Nutrition	from power adapter
<i>Equipment Power Adapter Specificationsuniversal neurosurgical and otorhinolaryngological (95-NC-NC):</i>	
Brand	Globetek
Model	GTM961600P16019-T3A
Overall dimensions (L×W), mm	240*120
Output characteristics:	
- rated DC voltage, V - direct current, A - power no more than, W	19 8.42 160

<i>Characteristics of the internal DC-DC converter for the Universal neurosurgical and otorhinolaryngological equipment (95-NC-NC):</i>	
Brand	GYVRM
Model	K241205
Input characteristics: • rated DC voltage, V	24
Output characteristics: • rated DC voltage, V • direct current, A • power no more than, W	12 5 60

Note – Tolerances for all characteristics are $\pm 2.5\%$ unless otherwise noted.

Table 5 Characteristics of instruments

Tool name	Overall dimensions, mm	Frequency of application	Type of delivery	Weight, g
1	2			3
Patient tracker	90	reusable	non-sterile	eleven
Bayonet pointer	310	reusable	non-sterile	89.3
Registrar	335	reusable	non-sterile	44.9
Front probe	365	reusable	non-sterile	62.5
Maxillary probe	360	reusable	non-sterile	62.5
Universal marker		reusable	non-sterile	31.4
Universal marker 360	168	reusable	non-sterile	29.4
Universal clamp C-shaped small	42	reusable	non-sterile	18
Universal clamp C-shaped medium	42	reusable	non-sterile	21
Universal clamp C-shaped large	53	reusable	non-sterile	22
Universal clamp C-shaped X large	64	reusable	non-sterile	24
Calibrator	74*50	reusable	non-sterile	84.0
Tray for sterilization of instruments for navigation during otorhinolaryngological interventions	415*335	reusable	non-sterile	2400
Articulated arm	280	reusable	non-sterile	1385
Cranial reference frame	110*115	reusable	non-sterile	126
Straight probe	300	reusable	non-sterile	68
Tray for sterilization of instruments for navigation during neurosurgical (cranial) interventions	415*335	reusable	non-sterile	2400
Double adapter	196	reusable	non-sterile	386
Biopsy targeting device	175	reusable	non-sterile	327

Biopsy probe	218	reusable	non-sterile	79
Reducing tube 2.2 mm	75	reusable	non-sterile	36
Tray for sterilizing instruments for biopsy navigation	415*335	reusable	non-sterile	2400
Silicone tape	89*25	disposable	non-sterile	86
Circle ENT phantom	160	reusable	non-sterile	1268
Phantom stand	400*350	reusable	non-sterile	5100
Circle Phantom neurosurgical	180	reusable	non-sterile	1367
Case for storing an ENT phantom (h*d*w), mm	230*550*315	reusable	non-sterile	6400
Case for storing a neurosurgical phantom (h*d*w), mm	335*610*480	reusable	non-sterile	6400
AccBar chamber calibrator (h*d*w), mm	279*53*50	reusable	non-sterile	145

Note – Tolerances for all characteristics are $\pm 2.5\%$ unless otherwise noted.

NaviENT equipment remains operational at temperatures from $+15^{\circ}\text{C}$ to $+30^{\circ}\text{C}$, relative air humidity from 10 to 90% and atmospheric pressure from 70 to 106 kPa.

NaviENT equipment, packaged in transport containers, is resistant to climatic influences, at temperatures from -25°C to $+45^{\circ}\text{C}$, relative air humidity from 10 to 98% and atmospheric pressure from 70 to 106 kPa.

NaviENT equipment is resistant to impacts caused by rough handling according to IEC 60601-1 and mechanical stress, with the following parameters:

vibration loads – frequency range 10-55 Hz with a displacement amplitude of 0.15 mm;
tools must be drop-resistant from a height of 1 m.

The force required to move equipment on an equipment trolley is no more than 200 N.

All instruments are supplied non-sterile, subject to sterilization, packaged in a sterilization tray (in which the instruments are supplied):

sterilization method - steam;
temperature regime – minimum 132°C ;
full cycle time – 4 minutes;
minimum drying time – 30 minutes.

Silicone tape 15-20 cm long, before use, is suitable for steam sterilization with the parameters described above, and is separately wrapped in paper for sterilization (packaged).

All tools, except the silicone tape, are reusable.

Silicone tape is disposable and should be discarded after use.

The patient tracker cannot be sterilized, only cleaned and disinfected.

The patient tracker can be cleaned using a brush with soapy water, deionized water and a clean soft cloth.

NaviENT equipment can be cleaned using a damp, soft cloth.

The external surfaces of parts of the NaviENT Equipment, such as the keyboard and laptop screen, touch monitor, security camera unit, folding arm, extension arm and security camera unit arm, are resistant to repeated manual disinfection using a soft cloth soaked in 70% isopropyl alcohol.

The patient tracker is resistant to repeated manual disinfection using a nylon-compatible disinfectant such as Metricide, Sporox or Sidex.

Protection against electric shock: ME EQUIPMENT operating from an external Class I electrical supply. Type of applied part: not applicable (not defined as APPLIED PARTS).

NaviENT equipment must meet electromagnetic compatibility requirements (Class A, Group 1).
Ingress protection against water and dust IP X0.
The service life must be 5 years.

5 MATERIALS USED IN MANUFACTURE

In the manufacture of NaviENT equipment the following are used:

Table 6

No.	Name of component, tool		Material	Manufacturer
1	Otorhinolaryngological equipment trolley		Stainless steel 420SS	PokTech RP, China
			ABS plastic PA200	HERBERT Co., Ltd, Thailand
2	Universal neurosurgical and otorhinolaryngological equipment trolley (95-NC-NC)		Stainless steel 420SS	PokTech RP, China
			ABS plastic PA200	HERBERT Co., Ltd, Thailand
3	Calibrator		Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
4	Registrar		Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
5	Universal marker and universal marker 360	Handle	Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
		Tip	Stainless steel 420SS	PokTech RP, China
7	Universal C-shaped clamps	Large head screw	Stainless steel 420SS	PokTech RP, China
		Long with large head screw	Stainless steel 420SS	PokTech RP, China
		Clamp	Stainless steel 420SS	PokTech RP, China
8	Maxillary probe	Handle	Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
		Tip	Stainless steel 420SS	PokTech RP, China
9	Front probe	Handle	Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
		Tip	Stainless steel 420SS	PokTech RP, China
10	Bayonet pointer	Handle	Aluminum AL 6061 T6	Xianfeng Aluminum

				Industries Co., Ltd, China
		Tip	Stainless steel 420SS	PokTech RP, China
eleven	Patient tracker		ABS plastic PA200	HERBERT Co., Ltd, Thailand
12	Self-adhesive silicone tape – transparent		Silicone	Midsun Specialty Products, Germany,
13	Straight probe	Handle	Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
		Tip	Stainless steel 420SS	PokTech RP, China
14	Biopsy probe	Handle	Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
		Tip	Stainless steel 420SS	PokTech RP, China
15	Articulated arm		Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
			Stainless steel 420SS	PokTech RP, China
16	Double adapter		Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
17	Biopsy targeting device		Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
18	Cranial reference frame		Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
19	AccBar Camera Calibrator		Aluminum AL 6061 T6	Xianfeng Aluminum Industries Co., Ltd, China
20	Reducing tube 2.2 mm		Stainless steel 420SS	PokTech RP, China
21	Trays for sterilization		Stainless steel 420SS	PokTech RP, China
22	ENT phantom		ABS plastic PA200	HERBERT Co., Ltd, Thailand

23	Phantom stand	ABS plastic PA200	HERBERT Co., Ltd, Thailand
24	Neurosurgical phantom	ABS plastic PA200	HERBERT Co., Ltd, Thailand
25	Roller stopper	ABS plastic PA200	HERBERT Co., Ltd, Thailand
26	Case for touch monitor	Foamed PVC	HERBERT Co., Ltd, Thailand
27	Case for security camera unit	Foamed PVC	HERBERT Co., Ltd, Thailand

6 USE OF NAVIENT SURGICAL NAVIGATION EQUIPMENT

6.1 Unpacking and installation

Open the package in accordance with the markings on the package and carefully remove the NaviENT equipment and its components:

- Check the NaviENT equipment and its components for any mechanical damage;
- If there is a temperature difference, it is necessary for NaviENT equipment to be at room temperature for at least 8 hours;
- If any problems are found, contact your authorized manufacturer representative immediately.

⚠ Warning

- Save packaging materials for future transportation and storage;
 - Disposal of packaging materials must be carried out in accordance with your local regulations.
-

6.2 Preparing NaviENT surgical navigation equipment for operation

Installation method

Install NaviENT equipment in the intended location of use so that there is a clear line of sight between the security camera unit and the markers on the instruments, as well as the ability to connect to AC power.

Power connection

Connect the power cord to an AC outlet suitable for use with medical equipment.

Inclusion

Before using NaviENT equipment, perform appropriate safety checks.

For surgical navigation equipment NaviENT, Otorhinolaryngological version (95-NE-NE):

- Turn on the laptop by pressing the power button;
- On your desktop, find the NaviENT software application and launch it.

For navigational surgical equipment NaviENT, version Universal neurosurgical and otorhinolaryngological (95-NC-NC):

- Turn on the system unit by pressing the power button;
- Turn on the touchmonitor by pressing the power button;
- On your desktop, find the NaviENT software application and launch it.

Shutdown

To turn off NaviENT equipment, please follow the following procedures:

- Quit working in the NaviENT software;
- Shut down your laptop or system unit as you would on any personal computer;

- Turn off the touchmonitor;
- Disconnect the NaviENT equipment from AC power by unplugging the power cord from the outlet.

7 RULES FOR CARE, CLEANING, DISINFECTION AND STERILIZATION

⚠ Warning

- **Before cleaning NaviENT equipment, turn it off.**

General information on cleaning and disinfection

- When cleaning NaviENT equipment, first turn off the laptop and unplug the power cord. Use a damp, soft, lint-free cloth to clean the outside of your cart and laptop. Avoid getting moisture into any openings. Do not spray liquid directly onto the laptop. Do not use aerosol sprays, solvents, abrasives, or cleaners containing hydrogen peroxide, which may damage the surface of parts.
- Use a damp cloth with 70% isopropyl alcohol to disinfect surfaces that may accidentally come into contact with liquids during surgery, such as the laptop keyboard and display, the security camera unit, and the cart that holds them. When preparing for surgery, ensure that all surfaces that may need to be touched with sterile gloves, such as the security camera unit, laptop trackpad, and cart, are covered with sterile drapes.
- The patient tracker is a reusable item. Between uses, it should be scrubbed with a damp soapy brush and rinsed thoroughly with deionized (DI) water to remove any visible residue and dried with a clean, soft cloth. Inspect the part for visible contamination. If necessary, repeat cleaning. After cleaning the part, it should be disinfected using a nylon-compatible disinfectant such as Metricide, Sporox or Sidex. Patient tracker tapes are single-use items and should be discarded immediately after use. Place new tape on the patient tracker after the patient tracker has been disinfected.
- Universal markers should be attached to the surgical instrument and sterilized together with it using autoclaving.

Cleaning and sterilization before use

1 Before sterilization, the following instruments must be cleaned before use:

- Bayonet pointer
- Registrar
- Front probe
- Maxillary probe
- Universal marker
- Universal marker 360
- Calibrator
- Patient tracker
- Universal clamp C-shaped small
- Universal clamp C-shaped medium
- Large universal C-shaped clamp
- Universal clamp C-shaped X large
- Large head screw
- Long screw with large head
- Articulated arm
- Cranial reference frame
- Biopsy targeting device
- Straight probe
- Biopsy probe

Warnings:

- Clean tools before first use and after each use.

- Before sterilization, remove the Universal Marker from the Universal C-Clip and Large Head Screw(s).
- Remove and properly dispose of the disposable silicone tape used for the Universal Marker and Patient Tracker.

Manual cleaning instructions:

- Rinse the instrument quickly and thoroughly under running tap water to remove all visible dirt. This step should be performed at the point of use to prevent contaminants from drying on the surface of the tool. You can use a soft-bristled brush or a clean cloth to remove dirt.
- Prepare an enzyme-based detergent according to the detergent manufacturer's recommendations, such as Medline One-Part Enzymatic Detergent and pre-soak.
- Immerse the instrument in the enzymatic detergent solution and allow it to soak for at least two (2) minutes.
- Clean the instrument thoroughly with a soft bristle brush and remove any foreign matter.
- Wash crevices, gaps, and blind holes in the instrument with a fresh detergent solution until the solution runs clear.
- Rinse instrument thoroughly in deionized (DI) water for at least one (1) minute to remove any remaining detergent.
- Rinse cracks, gaps and blind holes in the instrument with DI water.
- Dry the instrument with a clean, soft cloth.
- Check the instrument for visible dirt. If necessary, repeat cleaning.

Automatic cleaning instructions:

- Rinse the instrument quickly and thoroughly under running tap water to remove all visible dirt. This step should be performed at the point of use to prevent contaminants from drying on the surface of the tool. You can use a soft-bristled brush or a clean cloth to remove dirt.
- Prepare an enzyme-based detergent according to the detergent manufacturer's recommendations, such as Medline One-Part Enzymatic Detergent and pre-soak.
- Add the cleaning solution to the water tank (based on the size of the tools you will be cleaning and the work level line).
- Place the tools in the basket, slowly submerge it in the bath of cleaner, activate the ultrasonic machine and set the timer for 5-10 minutes and temperature according to the manufacturer's recommendations (do not allow the tools to touch the bottom of the tank).
- Once cleaning is complete, slowly remove them from the washer.
- Rinse instrument thoroughly in deionized (DI) water for at least one (1) minute to remove any remaining detergent.
- Rinse crevices, gaps, and blind holes in the instrument with fresh deionized (DI) water until the solution runs clear.
- Dry the instrument with a clean, soft cloth and place it on a clean bench.
- Check the instrument for visible dirt. If necessary, repeat cleaning.

2 Before use, sterilization must be performed for all instruments, with the exception of the patient tracker.

Warnings:

- Instruments must be properly cleaned before sterilization.
- Before sterilization, remove the plastic backing of the silicone tape.
- Before sterilization, remove the Universal Marker from the Universal C-Clip and Large Head Screw(s).

Instructions:

- Instruments should be individually wrapped in two layers of 1-ply polypropylene film

(Kinguard KC200–510(k) K082554) using sequential envelope folding techniques.

- Cut the silicone tape into 6-inch (15cm) pieces, remove the plastic backing and wrap them individually.
- The sterilization process must be carried out in accordance with the following parameters:
 - Sterilization type: steam
 - Temperature: 132°C
 - Exposure time: 4 minutes
 - Drying time: 30 minutes

3 Before use, the patient tracker must be disinfected.

Warnings:

- Do not sterilize the patient tracker. It will be damaged if heated above 80°C (176°F).

Instructions:

- Follow the cleaning instructions above and thoroughly clean the patient tracker.
- Prepare a nylon-compatible disinfectant (such as Metricide, Sporox or Sidex) according to the manufacturer's recommendations.
- Soak the patient tracker in disinfectant for the time recommended by the manufacturer.
- Remove the patient tracker from the solution and rinse thoroughly with a minimum of 2 liters of sterile water for at least one (1) minute to remove any remaining disinfectant.
- Dry it thoroughly with a soft, sterile, lint-free cloth and remove any remaining sterile water.
- Wrap the Patient Tracker in a sterile, lint-free cloth and place in a sealed bag or container. If the instrument is not to be used immediately, it should be stored in a well-ventilated, dust-free, closed cabinet in a manner that minimizes contamination or damage for 7 days, after which the instrument should be disinfected again before use.

8 MARKING

NaviENT equipment marking

Markings applied to the NaviENT equipment nameplate:

- product name;
- manufacturer's logo;
- product number (version);
- serial number;
- alternating current;
- Rated voltage;
- rated frequency;
- current consumption;
- date of manufacture (year);
- weight;
- symbol "Refer to instructions for use";
- "Attention" symbol;
- "Do not dispose of with household waste" symbol;
- CE marking;
- "Follow operating instructions" symbol;
- name and address of the manufacturer;
- authorized representative in the European Community.

Markings applied to the neurosurgical phantom case:

- product name;
- manufacturer's logo;
- product number (version);
- name and address of the manufacturer;
- "Do not dispose of with household waste" symbol;
- "Atmospheric Pressure Range" symbol;
- "Humidity range" symbol;
- "Temperature range" symbol;
- "Refer to instructions for use" symbol.

Note: separate markings are not applied to the neurosurgical phantom.

Markings applied to the ENT phantom case:

- product name;
- manufacturer's logo;
- product number (version);
- name and address of the manufacturer;
- "Do not dispose of with household waste" symbol;
- "Atmospheric Pressure Range" symbol;
- "Humidity range" symbol;
- "Temperature range" symbol;
- "Refer to instructions for use" symbol.

Note: separate markings are not applied to the ENT phantom.

Note: separate markings are not applied to tools, except for the tools of the additional set, each tool of which is marked with the letter "B".

Note: separate markings are not applied to instruments.

Marking applied to the packaging of NaviENT equipment:

- product name;
- manufacturer's logo;

- product number (version);
- serial number;
- date of manufacture (year);
- weight;
- symbol "Refer to instructions for use";
- "Attention" symbol;
- "Do not dispose of with household waste" symbol;
- CE marking;
- name and address of the manufacturer;
- authorized representative in the European Community;
- Authorized representative on the territory of the Russian Federation;
- Number and date of the registration certificate of the Russian Federation.

Explanation of all marking symbols:

Symbol	Meaning	Description
	Catalog number	Execution option
	Serial number	Product serial number
	Refer to instructions for use	It is necessary to study the instructions for use
	Attention	Notice that precautions must be taken to avoid injury or damage.
	Do not dispose of with household waste	Instructions on the need for disposal at separate collection points for waste recovery and recycling
	Manufacturer	Manufacturer
	Batch code	Indicates the batch code
	Upper temperature limit	Do not store at temperatures above 45 °C
	Not sterile	Indicates non-sterile parts
	Do not reuse	Reuse is not allowed
	fragile	Fragile. Carefully

Symbol	Meaning	Description
	Avoid exposure to sunlight	Avoid exposure to direct sunlight
	Humidity range	Specify the permissible upper and lower limits of relative humidity during transportation and storage
	Atmospheric pressure range	Indication of permissible upper and lower limits of atmospheric pressure during transportation and storage
	Temperature Range	Indication of the maximum and minimum temperatures for storage, transportation and operation
	Operating instructions	Operating instructions must be followed
	Don't push	It is forbidden to push
	Authorized Representative in the European Community	Authorized Representative in the European Community
	date of manufacture	date of manufacture
	Do not stack	Notice regarding the prohibition of stacking of products in packaging
	Keep away from moisture	Instructions to store transport packaging away from rain and in dry conditions
	Turning on the power	Indicates the power button
	Alternating current	Indication on the product plate about the need to connect the product through the appropriate connector only to AC sources
	CE marking	CE marking

9 PACKAGING

The instruments are supplied in a sterilization tray; the sterilization tray along with the instruments is packaged in a plastic bag with a zipper. There is no separate marking on the package.

Phantoms are packaged in plastic cases.

The equipment, along with tools, as well as operational documentation, must be placed in a box, wrapped in polyethylene. The layers of packaged products are separated from each other by a protective material that ensures the safety of this layer; a layer above a layer is not allowed.

The box must be tied with twine or covered with paper tape, paper-based adhesive or polyethylene tape with an adhesive layer, or sealed so that the package cannot be opened without breaking the integrity.

10 SAFETY AND ENVIRONMENTAL REQUIREMENTS

10.1 Safety requirements

Data Security:

- NaviENT equipment loads patient data sets from external data sources using CD/DVD, USB or network storage. Use only CD/DVDs and USB drives from safe sources. Do not connect them to your laptop unless they are from trusted sources as they may contain viruses or malware and may damage your laptop. DVDs and USB drives should only contain patient scan information generated by the scanning devices;
- NaviENT equipment comes with a pre-installed Microsoft firewall and provides basic security when connecting the equipment to the network. Do not connect your laptop to an unsecured wireless network.

Notes:

- NaviENT equipment is equipped with a secure password and user identification number. It is recommended to change the laptop password upon first use (the default password for first login is "Navient") and every 6 months thereafter, and to always lock the laptop display when left unattended;
- It is recommended to turn off wireless access at all times, except for downloading data from local network drives or online support. A laptop is like any other computer on your network and should be protected according to your local IT security policy.

For safety, NaviENT equipment must comply with the requirements of IEC 60601-1.

Protection against electric shock: ME EQUIPMENT operating from an external Class I electrical supply. Type of applied part: not applicable (not defined as APPLIED PARTS).

The security camera unit must additionally comply with the requirements of IEC 62471.

For electromagnetic compatibility, NaviENT equipment must comply with the requirements of IEC 60601-1-2 (Class A, Group 1).

Serviceability in accordance with IEC 60601-1-6 and IEC 62366.

Software safety class in accordance with IEC 62304 – class B.

NaviENT equipment tools must be biologically safe and comply with the requirements of the ISO 10993 group of standards.

10.2 Environmental protection

The materials used in the manufacture of NaviENT equipment must not pose a danger to life, human health and the environment both during operation and after the end of its service life.

The main type of possible hazardous impact on the environment is pollution of atmospheric air in populated areas, soils and waters as a result of:

- accidental leaks (spills) of production materials;
- unorganized disposal of waste on the territory of the manufacturer or outside it;
- arbitrary dumping of waste in places not intended for this purpose.

11 TRANSPORTATION AND STORAGE

11.1 Transportation conditions

Transportation of NaviENT equipment can be carried out by all types of covered vehicles and the rules for the transportation of goods that apply to each type of transport.

Transportation of NaviENT equipment by sea must be carried out in accordance with the “Safety Rules for the Maritime Transportation of General Cargoes”. Type of shipment - containers and small shipments.

NaviENT equipment, during transportation, must be resistant to climatic influences at ambient temperatures from -25°C to $+45^{\circ}\text{C}$, relative air humidity from 10 to 98% and atmospheric pressure from 70 to 106 kPa.

NaviENT equipment, during transportation, must be resistant to mechanical stress, with the following parameters:

- vibration loads – frequency range 10-55 Hz with a displacement amplitude of 0.15 mm.

During transportation, it is prohibited to: push the packaging with NaviENT equipment.

11.2 Storage conditions

NaviENT equipment, when stored, must be resistant to climatic influences at ambient temperatures from $+2^{\circ}\text{C}$ to $+45^{\circ}\text{C}$, relative air humidity from 10 to 98% and atmospheric pressure from 70 to 106 kPa.

Notes on storage conditions:

- Keep away from water;
- Keep away from direct sunlight;
- Do not store near heaters;
- Do not store in dusty places;
- Avoid locations subject to strong shock or vibration.

12 MAINTENANCE AND REPAIR

12.1 NaviENT equipment and associated components do not contain user serviceable parts.

12.2 For repair or replacement of any NaviENT equipment part or component, contact your technical support representative.

12.3 The security camera is specially designed to maintain the high accuracy of the optical calibration performed at the factory over a long period of time. But its accuracy can be compromised by unwanted changes to the equipment for a variety of reasons, such as environmental influences or rough handling. The R-Fine app provides a quick and intuitive method for assessing the accuracy of camera calibration and the ability to refine the calibration if necessary. The R-Fine app will remind the user to perform a calibration every 12 months.

The R-Fine app works by detecting the camera's special AccBar calibrator (included with Navient). The user must place the AccBar camera calibrator in several different areas and orientations, guided by the software using semi-transparent image overlays. Once positioned correctly, the coordinates of the AccBar camera calibrator are used to analyze the accuracy.

12.4 R-Fine Quick Start Guide

When opening the R-Fine app, the user is presented with a translucent image of the AccBar camera calibrator. The user must position the tool as shown on the screen.


- A yellow guide is displayed when the camera cannot detect a stripe.
- A red guide is displayed when a rod is detected but is in the wrong position.
- A green guide is displayed when the rod is installed and oriented correctly.


Try placing each plank on the guide


Sometimes you need to move the bar towards or away from the camera

The bar should be oriented as shown in the figure


In this position the bar is kept pointed towards the camera

Once all samples have been collected, the R-Fine app will determine if the camera is out of alignment. The app will automatically update the calibration file if the improvement exceeds 10%. In this case,

the refinement process should be repeated until the observed improvement is less than 10%. If the improvement is greater than 10% after 3 consecutive attempts, please contact support. This indicates that the camera is not calibrated.

Improvement of more than 10%. Repeat the calibration procedure.

Less than 10% improvement. The revision was successful.

Improvement by more than 10% after 3 attempts in a row. Contact Support.

12.5 Troubleshooting Guide

Various error messages are described and explained here, as well as possible causes and solutions.

Message	Explanation	Solution
Another copy of Navient is running. Please wait until it completes, then restart	You are only allowed to run one instance of the Navient application at a time	Close the running application or wait until it finishes and then try again
The laptop's power cable is disconnected. Make sure the Navient power cord is plugged into an outlet and the power to the laptop and tracking system is turned on.	The laptop power cable is not connected or there is no power.	Make sure the laptop's power cord is connected and the main Navient power cord is plugged into an appropriate outlet.
The camera may not be calibrated. It is highly recommended to run R-Fine every 12 months	To ensure optimal navigation accuracy, the Navient camera must be calibrated every 12 months	Click "Yes" to run R-Fine, or "No" to ignore (not recommended).
Low resolution warning, image spacing is greater than 1.0mm. Navigation accuracy may be reduced	The distance between voxels should be less than 1.0 mm.	Use extreme caution when navigating and perform an accuracy check by touching visible nearby structures before navigating.
Import error, invalid path	You may be trying to import a date from Navient's main data store (where it already resides) or from an inaccessible directory, such as a network drive or portable media that is not currently mapped.	Please correct the path problem and try again.
Failed to import. At least 100 slices required	Navient is capable of processing datasets containing 100 or more slices	Refer to Navient CT Instructions
Import error. There is not enough disk space to import. You want to delete data from storage	Disk is full	To continue working, you must delete existing data sets
Directory exists, overwrite patient data?	The data to be imported already exists in the data warehouse	Select Yes to overwrite the existing data and

Message	Explanation	Solution
		plans file. Select "No" to cancel the import.
Drive error display	There is a problem displaying the specified drive	Please check the disk path and user credentials (if required)
Warning about loading slice stack. More than one stack of slices was found in the folder.	The dataset may be damaged (for example, missing or corrupted files), or multiple scanned files may be stored in the same folder.	The application will open the data and it is recommended that the user ensure the completeness of the patient data by viewing the snippets
Security camera error. Please disconnect and reconnect the camera, then wait a few seconds	MT camera is either disconnected or there is a communication error	Restart the application, and if the application still cannot connect to the camera, restart your computer
Confirm deletion, permanently delete all data of the selected patient(s) from storage?	Removing a dataset from the Navient repository	Selecting "Yes" will delete data from the storage
Not enough disk space	The folder you are exporting data to does not have enough disk space	Free up some space and try again
Directory exists, overwrite patient data?	The exported data already exists in the destination folder	If you select "Yes" the destination folder will be overwritten
Adjust the camera, move Right/Left/Up/Down/Closer/Further	For optimal accuracy, the forehead mark must be within the camera's field of view	In case the markers go out of bounds, the application will warn the user to adjust the camera accordingly
Camera operating outside tested temperature range	The camera may be heating up	You should investigate the cause of the temperature increase, once the problem is resolved, wait until the camera cools down before continuing.
Navigation is inactive	While navigating, if the forehead mark or navigation pointer is not detected (for any reason, such as line of sight interruption), the app will show a sidebar with live camera video.	Check if anything is blocking the FLOW. Make sure there is a clear line of sight between the navigation indicators and the camera body.

12.6 In case of post-warranty and out-of-warranty repair of NaviENT equipment, contact the authorized representative of the manufacturer in the Russian Federation: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address 125445, Moscow, sh. Leningradskoye, 69, building 1, room I floor 1 room 50; phone +7 495 645 47 00; email davletov@im-neuro.ru.

13 SOFTWARE WORKFLOW

Workflow overview

Navient software has been pre-installed on the Navient computer. The shortcut for the Navient application, the Navient Repository and Navient Recordings folder are placed on the desktop.

Note 1: Make sure the main power cord is connected to the power outlet. Once connected a blue LED diode in front of the CameraBox will be on.

Note 2: When launching Navient, a warning message will be displayed if the Navient laptop is not connected to the power. See Appendix 3 for more details.

Navient Repository is a folder where the Navient application stores all imported datasets locally. Section **Ошибка! Источник ссылки не найден.** explains importing patient scan data from CD/DVD to the repository and loading it to the application for guidance.

Navient Recordings is a folder where the Navient application stores videos recordings, screen shots, and application state and log files. Application state file contains information like patient registration, planning and landmarks.

Navient software operation proceeds sequentially through 5 stages, or screens:

1. **Images:** Load the CT or MR images, review and explore images and perform image fusion.
2. **Planning:**
 - a. **Landmarks:** Mark easily identifiable anatomical landmarks on the patient scan.
 - b. **3D modeling:** Generate 3D model of the tumor or sinus.
 - c. **Biopsy:** plan biopsy trajectory or define target points for navigation.
3. **Set Up:** attach the Patient reference, set up the tracking system position and verify the instruments' tip calibration accuracy.
4. **Registration:** Using registrator, point to the marked landmarks on the actual patient's skin, then trace a path on the patient's face to refine the registration.
5. **Navigation:** Observe the position of the tip of various instrument in the Image volume.

Splash screen

When started, Navient displays a “splash” screen that provides the software version number and license information. The splash screen requires the user to select the procedure type. Navient supports trans-cranial and trans-nasal procedures. If the system includes more than one set of instruments, then a dropdown box will be displayed. The user must select the serial number of instruments in order to proceed.

Patient data loading and importing


Once Navient is launched, it opens the data browser window, allowing the user to load, import and delete patient image data. The data browser window can also be initiated from the **Load/Import** button on the bottom of the Navient sidebar.

Navient stores all imported patient scans on a desktop folder called **Navient Repository**. The user can load patient data from the repository or import data from other locations or devices (like CD or DVD and USB drives) to the repository.

Loading data

The Data Browser searches through the Navient Repository folder and shows the search results in two panes. On the top pane it lists the patient information for all patients for which some DICOM-formatted images were found. The browser auto selects the first discovered patient, but the user can click on other rows to change that selection. On the pane below, the discovered image series (slice stack) for the selected patient are shown.

The user can select the desired series and load it with a click on the Load button or with a double click on the corresponding row in the lower pane. The loading automatically triggers the advancement to the IMAGES stage.

- **Import button:** Enables importing a new dataset from devices like CD or DVD and USB drives to the Navient Repository.
- **Load:** Loads the selected patient scan and shows 3D and 2D views.
- **Delete:** Deletes the selected patient data.
- **Refresh:** Refreshes the Data Browser window.
- **Help:** Opens up help window
- **Cancel:** closes the Data Browser window.


Note: Always load the most recent patient data (recommended less than 6 month old) with largest slice number (200 or more) and smallest slice spacing (1mm or less). Navient will show warning messages if the data is more than 2 years old or if it has slice spacing of more than 2mm or if the number of slices are less than 30 slices.

Importing Data

The user can import patient data from CD, DVD or USB drives to the repository.

1. Click on the Import button to open up the folder browser.

2. Select the device that you want to import the data and browse and select the patient data folder.
3. Click **Select Folder** to return to the Data Browser form. Data Browser will display all scans found under the selected folder.
4. Select the scan that you want to import and click on **Import & Load** (see picture below). It will import the patient data to the local **Navient Repository** and then will load the data from the repository.


Note 1: Once data is imported, it will be copied to the **Navient Repository** folder on the desktop. A warning message will be displayed if you try to reimport the same data twice. Once imported, the imported data will be displayed in the Data Browser window.

Note 2: A troubleshooting window will show up if the tracking system is inactive. Follow the instructions displayed in the troubleshooting window (picture below) to resolve the issue.

What is the tracking camera LED status?

Depending on the tracking camera status LED, follow the instructions stated under A, B, or C.

A. If LED is off

Check all 3 connections as displayed in the picture and make sure they are firmly connected. Also try:


1. Unplug connection #1 and plug it to another port on the laptop.
2. Replace the thunderbolt adaptor with a new one.

After connecting the cables, if the status LED is on:

1. Click Retry button to test connection.
2. If the tracker is still inactive, follow instructions B

If the issue still persists, please contact technical support.


B. If LED is constant bright green

First check the tracking cable connections to the laptop, white adaptor and camera and then reboot the laptop. Click the reboot button, and after reboot, re-login with default password "navient" and restart NaviENT.


To avoid this from happening again, please do not remove the white adaptor cable from the laptop.

If the issue still persists, please contact technical support.

C. If LED is blinking

Blinking LED is an indication of a hardware issue, please contact technical support.

Stage 1: IMAGES

Once the patient images are loaded, the application will proceed to the "IMAGES" stage. It displays the 3D patient model as well as 2D axial, coronal and sagittal reformatted slice views and allows the user to rotate, zoom, change contrast and explore the patient scan. Brief information of the patient like ID, Name, Age, Gender and physician name are also displayed on the sidebar. The main goal of this stage is to:

- Verify Patient info:** Make sure the patient info displayed on the side bar panel matches with the actual patient.
- Adjust 2D and 3D images contrasts:** use the toolbar contrast button and adjust the windowing width and center of images. The software will preserve all changes from this stage and will display the exact same contrast during navigation stage.
 - Double clicking on the **Contrast** button resets the contrast to the default value.
 - Click on the **Contrast** button and drag up/down on the 2D or 3D images to change windowing center. Drag left/right to change windowing center.
 - Click on **Presets** button to open the presets drop-down menu. The menu lists available 3D presets options for rendering 3D model. Select an appropriate preset to visualize 3D model as desired.

- F. Correct 3D image orientation:** The patient 3D image must be facing forward as displayed below. Use the Rotate toolstrip to orient the 3D image. This initial orientation will be used as reference for proper presentations in the following stages.




Image Fusion

Navient provides an option to load and merge CT and MR scans of the same patient. the feature is displayed on stage IMAGES. Below are the instructions to load, merge and verify merged images:

1. Load the first modality and confirm the patient info and orient the 3D image.

2. Click on "Add image" button to load the second modality. Once loaded review the patient info and correct the 3D image orientation if needed.


3. Follow the instructions and place 4 landmarks on the CT and MR images. Right click (Tap and hold on touch screen) to place landmarks on the images. It is important to place the CT and MR landmarks as close as possible.

Once all four landmarks are placed, the "Merge" button will be enabled. Click on the button to start image fusion. The process may take a few seconds to complete.

4. Navient will proceed to "Verify" sub stage automatically. Right click (Tap and hold on touch screen) on any anatomy on the images and make sure the software automatically finds the correct corresponding anatomy on the other image.


It is important to review and make sure the two images are fused accurately. Select Axial, Coronal and Sagittal views to review the accuracy of fusion. If the images are not fused accurately, then go back to "Merge" stage and repeat marking the landmarks and the image merging process.

Navient also provides options to overlay the fused images:

- Fused Magenta / Green
- Fused Gray / Gray

The user could also modify the windowing width/level of CT and MR images on 2D and 3D views.


Important Note: The Patient scan must be taken according to instructions provided by Navient Imaging protocol.

FOV (Field of View) coverage area:


- **For trans-cranial procedures**

- Place the patient's head in a foam headrest. If the patient cannot remain still, use a scanner cradle. Separate the patient's head from the scanner cradle with towels or a foam headrest.
- Scan the entire head (Green Area): Include the hard palate, the mouth, tip of the nose, the ears to the top of the head.
- Ensure that the tip of the patient's nose is included in the scan and that it is the most anterior point on the scan.

- **For trans-nasal procedures**

- For the trans-nasal cases, the following anatomical regions must be visible in the patient scan: forehead, dorsum of nose, base of nose and eye laterals.

Scanning protocol:

1. Use square pixels of 0.5mm x 0.5mm – 1.0mm x 1.0mm
2. Use constant slice thickness of 0.5mm – 1.0 mm
3. Axial slices are preferred (Helical scans are also acceptable).
4. Scan contiguous, non-overlapping slices.
5. No motion artifacts or no gantry tilt.
6. Archive the exam from the scanner to a CD or DVD
7. The scan should be uncompressed when saved (non zipped)
8. The scan should be Multi-File DICOM
9. The CD should only contain DICOM data
10. Save scan data into one folder only

Refer to Navient Scanning instructions for a comprehensive Scanning instructions.

Tool strips buttons

All viewports include tool strip buttons that enable basic image manipulations like rotation, zoom, pan, scroll and modifying windowing levels. Select the button on the toolstrip and then drag the mouse on the corresponding viewport.

- **Contrast:** changes the windowing width (Contrast) and level (Intensity) of the images on 2D & 3D views. In order to change windowing width (contrast) drag the mouse left/right on the viewport. In order to change windowing level (intensity) drag Up/Down.
- **Presets:** allows loading and setting various predefined presets to the 3D image. The user can also edit and create new presets.
- **Scroll:** scrolls the slice in/out along the depth dimension of the view. Available only on 2D viewports. In order to scroll, drag up/down.
- **Rotate:** Rotates patient image and is available in both 3D and 2D viewports. In order to rotate the image, select the button and drag on the viewport in the direction you wish to rotate.
- **Rotate 2D:** Rotates patient image around the axis perpendicular to the viewport.
- **Pan:** Pans the images and is available in both 3D and 2D viewports.
- **Zoom:** Zooms in and out the images and is available in both 3D and 2D viewports. In order to zoom in/out the image, select the button and drag up/down on the viewport.
- **Front cut:** Clips the front of the 3D image. In order to change the clip plane, select the button and drag up/down on the 3D viewport. A checkmark next to the front cut button indicates if the feature is active or not.
- **Ruler:** Used for measurement. Select the button, tap and drag on the viewport. Double click on the button to delete all measurements.

14 DISPOSAL INFORMATION

Disposal of NaviENT equipment must be carried out in accordance with SanPiN 2.1.3684-21 for waste class A, including instruments after preliminary disinfection and sterilization, belongs to class A - epidemiological safe waste.

NaviENT equipment must be disposed of in the following cases:

- **end of service life;**
- **confirmation of facts and circumstances that pose a threat to the life and health of health workers and indicate a failure to fulfill the intended purpose.**

All packaging, including transport packaging, must also be recycled. Paper, wood, polyethylene and plastic must be disposed of separately in accordance with the rules for collection, accounting and disposal of packaging.

15 MANUFACTURER WARRANTY

NaviENT equipment is guaranteed for a period of 2 years from the date of sale from the manufacturer's warehouse.

Guaranteed shelf life is 12 months from the date of delivery.

NaviENT equipment is warranted to be free from defects in workmanship and materials. Any defective part will be replaced if the NaviENT equipment has not been tampered with and the NaviENT equipment has been properly operated during this period. Ensure that any malfunction is not due to inadequate care of the NaviENT equipment or failure to follow the operating instructions. In the event of a malfunction, subject to operating requirements in accordance with the operating instructions, the manufacturer carries out routine repairs or replacement of NaviENT equipment or its components.

The manufacturer or the manufacturer's authorized representative in the Russian Federation is not responsible for any personal/property damage caused by erroneous use, operation for purposes other than the intended purpose, and negligence in servicing the product.

Please note that if NaviENT equipment is damaged by the use of unauthorized cleaning or disinfection materials, no warranty will be provided.

Authorized representative of the manufacturer in the Russian Federation: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address 125445, Moscow, sh. Leningradskoye, 69, building 1, room I floor 1 room 50; phone +7 495 645 47 00; email davletov@im-neuro.ru.

WARRANTY CARD

Product category _____

Manufacturer _____

Execution option _____

Serial number/art. _____

Date of sale _____

Organization seller _____

Seller's signature _____

The product has been checked and has no damage.

I have read and agree with the operating rules.

Buyer's signature _____

The warranty period is 2 years from the date of sale from the manufacturer's warehouse.

For warranty questions, please contact: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at 25445, Moscow, sh. Leningradskoye, 69, building 1, room 1 floor 1 room 50; phone +7 495 645 47 00; email davletov@im-neuro.ru.

In case of post-warranty and out-of-warranty repair of the product, contact: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address 25445, Moscow, sh. Leningradskoye, 69, building 1, room 1 floor 1 room 50; phone +7 495 645 47 00; email davletov@im-neuro.ru.

ATTENTION! Make sure that the warranty card is filled out correctly. The serial number of the product/article must exactly match the number on the coupon. All columns of the coupon (except for the columns for the guarantee center) must be completed.

If the coupon is filled out inaccurately or incompletely, the product warranty is considered lost.

The manufacturer's warranty is carried out in accordance with the Civil Code of the Russian Federation on a general basis.

The product should not bear traces of mechanical or thermal influences, as well as traces of exposure to various chemical agents. Operation of the product is permitted only in accordance with the instruction manual for the product.

WARRANTY NOTE

DATE	DESCRIPTION	WARRANTY CENTER

16 INFORMATION ABOUT COMPLAINTS

Complaints should be sent to the authorized representative of the manufacturer in the Russian Federation: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address 125445, Moscow, sh. Leningradskoye, 69, building 1, room I floor I room 50; phone +7 495 645 47 00; email davletov@im-neuro.ru .

17 ELECTROMAGNETIC COMPATIBILITY

1. Electromagnetic Compatibility Specifications

Definition: EMC (electromagnetic compatibility) is the ability of equipment to avoid electromagnetic interference generated by other devices without causing similar electromagnetic interference to other equipment.

Table 1 - Guidance and manufacturer's declaration - electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions		
NaviENT surgical navigation equipment is intended for use in the electromagnetic environment specified below. The purchaser or user of NaviENT Surgical Navigation Equipment should ensure that it is used in the specified electromagnetic environment.		
Electromagnetic Emission Test	Correspondence	Electromagnetic environment - instructions
The group to which ME belongs according to CISPR 11 (GOST CISPR 11-2017)	Group 1	NaviENT surgical navigation equipment uses radio frequency energy only to perform internal functions. RF emissions are low and are not likely to cause interference with nearby electronic equipment
The class to which ME belongs according to CISPR 11 (GOST CISPR 11-2017)	Class A	Surgical navigation equipment NaviENT is suitable for use in all locations, except residential buildings and buildings directly connected to the electrical distribution network that supplies residential buildings
Harmonic components of current according to IEC 61000-3-2 (GOST IEC 61000-3-2-2021)	Class A	
Voltage fluctuations and flicker according to IEC 61000-3-3 (GOST 30804.3.3-2013)	Compliant	

Table 2. Guidance and Manufacturer's Declaration - Immunity

Guidance and Manufacturer's Declaration - Immunity			
NaviENT surgical navigation equipment is intended for use in the electromagnetic environment specified below. The purchaser or user of NaviENT Surgical Navigation Equipment should ensure that it is used in the specified electromagnetic environment.			
Immunity test	Test level according to GOST R IEC 60601-1-2-2014	Compliance level	Electromagnetic environment - instructions
Electrostatic discharges (ESD) according to IEC 61000-4-2 (GOST 30804.4.2-2013)	±6 kV - contact discharge	±2 kV - contact discharge	The floor in the room is made of wood, concrete or ceramic tiles. For floors covered with synthetic material, the relative air humidity is at least 30%. Recommended for use at high relative humidity. If the compliance level is discharged above the specified one, a failure may occur: the security camera unit will turn off. The product needs to be rebooted. If the above method does not help, you must use antistatic wrist straps.
	±8 kV - air discharge	±4 kV - air discharge	
Nanosecond impulse noise according to IEC 61000-4-4 (GOST 30804.4.4-2013)	±2 kV - for power lines	±2 kV - for power lines	Typical network for industrial or medical purposes.
	±1 kV - for input/output lines	-	
High energy microsecond pulse noise according to IEC 61000-4-5 (GOST R 51317.4.5-99)	±1 kV when applying interference in a wire-to-wire manner	±1 kV when applying interference in a wire-to-wire manner	For stable operation, it is necessary to use external surge protection devices, voltage stabilizers or uninterruptible power supplies!
	±2 kV when noise is applied via the wire-to-ground circuit	±2 kV when noise is applied via the wire-to-ground circuit	

Voltage sags, short-term interruptions and voltage changes in input power lines according to IEC 61000-4-11 (GOST 30804.4.11-2013)	<5% Un (voltage dip >95% Un) for 0.5 period	<5% Un (voltage dip >95% Un) for 0.5 period	Typical network for industrial or medical purposes. For continuous operation under voltage surges, it is recommended to use an uninterruptible power supply or battery.
	40% Un (voltage dip 60% Un) for 5 periods	40% Un (voltage dip 60% Un) for 5 periods	
	70% Un (voltage dip 30% Un) for 25 periods	70% Un (voltage dip 30% Un) for 25 periods	
	<5% Un (voltage dip >95% Un) for 5 s	<5% Un (voltage dip >95% Un) for 5 s	
Power frequency magnetic field (50/60 Hz) according to IEC 61000-4-8 (GOST IEC 61000-4-8-2013)	3 A/m	3 A/m	Magnetic fields from the network should not exceed typical levels for typical industrial networks.

NOTE: Un is the voltage level of the electrical network before the test stimulus is applied.

Table 3. Guidance and Manufacturer's Declaration - Immunity

Guidance and Manufacturer's Declaration - Immunity NOT INTENDED TO SUPPORT VITAL FUNCTIONS			
NaviENT surgical navigation equipment is intended for use in the electromagnetic environment specified below. The purchaser or user of NaviENT Surgical Navigation Equipment should ensure that it is used in the specified electromagnetic environment.			
Equipment stability testing	Test level according to GOST R IEC 60601-1-2-2014	Compliance level	Electromagnetic environment - Instructions
			Portable and mobile radio frequency equipment, incl. cables should not be used near the device closer than the recommended distance, calculated by the formula according to the frequency of the transmitter:
Recommended distance			
Conducted interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6 (GOST R 51317.4.6-99)	3B from 150kHz to 80MHz	3B from 150kHz to 80MHz	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$
Radio frequency electromagnetic field according to IEC 61000-4-3 (GOST 30804.4.3-2013)	3V/m from 80MHz to 2.5GHz	3V/m from 80MHz to 2.5GHz	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz
			$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz

Where P is the maximum transmitter output power in watts (W) according to the manufacturer's specifications, and the recommended distance in meters (m).

d - recommended spatial spacing, m;

The field strength of radio waves from fixed radio transmitters, as determined by observations of the electromagnetic environment a) should be lower than the compliance level in each frequency band b)

Interference may occur in the vicinity of equipment marked with:

a) The field strength of radio waves from fixed radio transmitters, such as base stations of radiotelephone networks (cellular/wireless), and land mobile radio stations, amateur radio stations, AM and FM radio broadcast transmitters, television transmitters cannot be determined by calculation with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If measured values at the location of the NaviENT Surgical Navigation Equipment exceed the applicable compliance levels, the operation of the NaviENT Surgical Navigation Equipment should be monitored to verify its normal functioning. If monitoring reveals a deviation from normal operation, additional measures may need to be taken, such as reorienting or relocating the NaviENT Surgical Navigation Equipment.

b) Outside the frequency range from 150 kHz to 80 MHz, the field strength should not exceed (V1) V/m.

Table 4. Recommended separation distances between portable and mobile RF communications and equipment

Recommended separation distances between portable and mobile RF communications and equipment NOT INTENDED TO SUPPORT VITAL FUNCTIONS			
NaviENT surgical navigation equipment is intended for use in an electromagnetic environment in which the levels of radiated interference are controlled. The purchaser or user of NaviENT Surgical Navigation Equipment can avoid the effects of electromagnetic interference by maintaining a minimum separation distance between portable and mobile RF communications equipment (transmitters) and the equipment as recommended below, taking into account the maximum output power of the communications equipment.			
Maximum transmitter output power (W)	Distance depending on transmitter frequency (m)		
	150 kHz + 80 MHz d = $\left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz + 800 MHz d = $\left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz + 2.5 GHz d = $\left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

When determining the recommended spacing values d for transmitters with a rated maximum output power not listed in the table, substitute the rated maximum output power P in watts specified in the transmitter manufacturer's documentation into the given expressions.

Notes: 1 At frequencies of 80 and 800 MHz, higher field strengths apply. 2 The expressions given do not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people. 3 When determining the recommended separation distances d for transmitters with a rated maximum output power not listed in the table, substitute the rated maximum output power P in watts given in the table into the expressions given. transmitter manufacturer documentation

⚠ Warning

Because this NaviENT equipment generates, uses, and radiates radio frequency energy, it may interfere with the operation of other medical and non-medical devices and wireless communications. It is characterized by the creation of electromagnetic interference in the air or cable environment. This product is designed to fully comply with EMC standards and provides reasonable protection against such interference.

⚠ Warning

NaviENT equipment requires special measures to ensure ELECTROMAGNETIC COMPATIBILITY and must be installed and commissioned in accordance with the EMC-related information provided in the instruction manual.

Note:

- **Use of devices such as cellular phones, transceivers, or remote control devices that generate radio waves near this product may cause interference with the operation of this product;**
- **When installing, ensure that this equipment is located as far away from other electrical equipment as possible;**
- **Use only cables provided or designed by our company;**
- **Do not attempt to modify this equipment. Modifications to the product may reduce the EMC performance. Changes include modifications to cables, failure to properly install or position the product, changes to its configuration or components, changes to established equipment operating procedures, etc.**

2. Measures to eliminate EMC problems

If this equipment is found to cause harmful interference, the user is encouraged to correct the problem by one or more of the following measures:

- **Locate other devices as far as possible from this equipment to reduce electromagnetic interference;**
- **To reduce electromagnetic interference, you can adjust the position of the equipment and other devices or the angle between them;**
- **Electromagnetic interference can be reduced by changing the connection point of the NaviENT equipment's power cord;**
- **You can also change the power path to other devices to reduce electromagnetic interference.**

The manufacturer assumes no responsibility for interference caused by using a non-recommended connecting cable or by modifying or altering NaviENT equipment without approval. Unauthorized modification or alteration to NaviENT equipment may void the user's authority to operate it.

List of regulatory documentation used by the manufacturer to confirm compliance:

Standard	Name
ISO 14971	Medical Devices – Application of Risk Management to Medical Devices
ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
IEC 62366	Medical Devices – Application of Serviceability
IEC 62304	Medical Device Software – Software Life Cycle Processes
EN 1041	Information provided by the Manufacturer of Medical Devices
ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing as part of the risk management process
ISO 10993-5	Biological evaluation of medical devices - Part 5: In vitro cytotoxicity tests
ISO 10993-10	Biological evaluation of medical devices - Part 10: Skin irritation and sensitization tests
ISO 15223-1	Medical devices - Symbols to be used on medical device labels, markings and information to be provided - Part 1: General requirements
ISO 17664	Sterilization of medical devices - Information that must be provided by the manufacturer for the processing of refurbished medical devices
ISO 17665-1	Requirements for the development, validation and regular monitoring of the sterilization process for medical devices
IEC 60601-1	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Related standards: Electromagnetic compatibility: Requirements and tests
IEC 62471	Photobiological safety of lamps and lamp systems
NEMA PS 3.1-3.20	Digital Imaging and Communications in Medicine (DICOM) - 2011

The translation of this document from English into Russian was carried out by translator Yulia Aleksandrovna Chigrinova. The translation is correct, accurate and complete.

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©2020 ClaroNav Kolahi Inc. (CKI)