
11006 98952

MANUFACTURER'S STATEMENT

To whom it may concern;

Collagen Matrix, Inc. with a business address of 15 Thornton Road, Oakland, New Jersey, 07436, USA hereby confirms that the attached Instructions for Use (IFU) for SYMBIOS® Collagen Membrane pre-hydrated is a draft copy of the current Instructions for Use (RM0906SYM/Rev B) containing additional specific information for registration in the country of Russia.


Gloria Zulich
Director, Regulatory Affairs

July 18, 2018

Date

Subscribed and sworn to before me this 18th day of July month, 2018 year.

MATEUSZ F. LESZCZAK
NOTARY PUBLIC OF NEW JERSEY
Comm. # 50065913
My Commission Expires 8/10/2022

15 Thornton Road • Oakland, NJ 07436 • USA

CollagenMatrix.com • 201-405-1477

Symbios®
Collagen Membrane pre-hydrated

INSTRUCTIONS FOR USE

Description: Symbios Collagen Membrane pre-hydrated is a pre-hydrated, white, nonfriable collagen membrane matrix derived from bovine pericardium. Symbios Collagen Membrane pre-hydrated is resorbable which eliminates the need for a second surgical procedure that is normally required to remove a non-resorbable membrane. Symbios Collagen Membrane pre-hydrated is sterile, non-pyrogenic and for single use only.

Indications for Use: Symbios Collagen Membrane pre-hydrated is a bioresorbable, implantable collagen material that is intended for use in oral surgical procedures as a resorbable membrane material for use in:

- augmentation around implants placed in immediate extraction sockets, delayed extraction sockets;
- localized ridge augmentation for later implantation;
- alveolar ridge reconstruction for prosthetic treatment;
- filling of bone defects;
- guided bone regeneration in dehiscence defects; and
- guided tissue regeneration procedures in periodontal defects.

Contraindications: Symbios Collagen Membrane pre-hydrated is contraindicated in patients who have:

- acute infections or contaminated wound in the oral cavity;
- known allergy to collagen of animal origin or other bovine-derived products; and/or
- clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which in the physician's judgment will prevent safe implantation or likely healing.

Warning: Clinicians should screen patients prior to surgery to confirm they are not hypersensitive to bovine collagen. Hypersensitivity reactions have been noted with the use of other products containing bovine collagen; therefore, the possibility exists of developing a local sensitivity response to Symbios Collagen Membrane pre-hydrated.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed physician or dentist.

Precautions: As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroidal therapy or currently taking anticoagulants. Patients with clinically significant systemic diseases, indicating a history of anaphylactic reactions, autoimmune diseases, uncontrolled diabetes or severe hypertension have not been implanted with the membrane; therefore, the safety and effectiveness for those patients have not been determined. Nor has it been evaluated in pregnant women, children and/or in patients with conditions involving extremely severe defects with little periodontium or bone.

Symbios Collagen Membrane pre-hydrated cannot be re-sterilized. Open, unused Symbios Collagen Membrane pre-hydrated must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used.

The membrane product is supplied in a pre-hydrated state. If the product is found not properly pre-hydrated upon opening, the product is to be discarded and not rehydrated for clinical use.

Do not use if the product sterilization barrier or its packaging is compromised.

Do not remove the membrane from the package until just prior to surgery.

RM0906SYM- Rev B DRAFT RUSSIA
Corrected 18-07-2018

Adverse Reactions: Possible complications that can occur with any dental surgery include infection, swelling of the intraoral tissue, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing, resorption or ankylosis, with loss of crestal bone height, pain, or complications associated with the use of anesthesia. Minor discomfort may occur for a few days.

Application: Symbios Collagen Membrane pre-hydrated is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The pre-hydrated membrane along with the insert card should be removed from the inner pouch with sterile gloves or instruments.

The bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). Space-making material such as autologous bone, demineralized bone matrix and ceramic materials may be used to fill the defect.

The pre-hydrated membrane is removed from the insert card prior to placement and can be trimmed to the size and shape of the defect using sharp, sterile scissors. The pre-hydrated membrane can be placed as is without further hydration. However, if the clinician prefers, the membrane can be further hydrated in sterile water or sterile saline prior to placement.

Symbios Collagen Membrane pre-hydrated should overlap the walls of the defect by at least 2 mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. The mucoperiosteal flap is sutured over the collagen membrane and the wound should be closed completely to avoid accelerated resorption due to membrane exposure. The membrane is expected to be essentially resorbed in approximately 16 weeks.

Post-operative Procedures: Symbios Collagen Membrane pre-hydrated is completely resorbable and should not be removed. Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent.

The patient should refrain from brushing the treated area for two weeks following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated. The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

Probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.


How Supplied: One (1) membrane per package.

Available sizes: 15 mm x 20 mm, 20 mm x 30 mm, or 30 mm x 40 mm

RM0906SYM- Rev B DRAFT RUSSIA
Corrected 18-07-2018

Additional Country Specific Storage and Product Information: Symbios Collagen Membrane pre-hydrated does not contain any dangerous hazardous substances that would require special handling. Symbios Collagen Membrane pre-hydrated can be disposed of in accordance with all local, state, or federal regulations.

Store at room temperature (15 °C/59 °F to 30 °C/86 °F) in a dry place (relative humidity about 10% to 80%). The shelf life of the device is 12 months from the date of manufacture. The manufacturer guarantees the operability and maintenance of the sterility of the products during the entire shelf life, subject to storage conditions. The product can be transported by any means of transport, such as covered railway wagons, containers, freight vehicles, ship holds, aviation modules in accordance with the rules for these to modes of transport, in compliance with the conditions of transportation without risk of damage to the package. The product does not affect the environment during operation. Minor excursions from the recommended temperature range during storage, shipping, and handling are acceptable. The manufacturer cannot guarantee the product to meet specifications, if stored at temperatures outside of the recommended temperature range for more than 3 days.

Labeling Symbols:

	Manufacturer		Do not use if package is damaged
	Use-by date		Temperature limitation
	Do not re-use		Sterilized using irradiation
	Caution: See instructions for use		Do not re-sterilize
	Batch number	R_x Only	Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist or physician.
	Catalogue number		Authorized representative in the European Community


Manufacturer:
Collagen Matrix, Inc.
15 Thornton Road
Oakland, New Jersey 07436 USA

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Distributor by:

www.dentsplyimplants.com

Authorized representative of the manufacturer in the territory of the Russian Federation: LLC "Dentsply Sirona", Russian Federation, 115432 Moscow, Andropov, Prospekt 18, Building 6 tel.: +7 (495) 725 1087

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

MANUFACTURER'S STATEMENT

To whom it may concern;

Collagen Matrix, Inc. with a business address of 15 Thornton Road, Oakland, New Jersey, 07436, USA hereby confirms that the attached Instructions for Use (IFU) for SYMBIOS® Collagen Membrane SR is a draft copy of the current Instructions for Use (OSCIFU-OU 2014-09 A) containing additional specific information for registration in the country of Russia.

Gloria Zuclich
Director, Regulatory Affairs

June 1, 2018

Date

Sworn to and subscribed
before me this

1st day of June, 2018

MATEUSZ F. LESZCZAK
NOTARY PUBLIC OF NEW JERSEY
Comm. # 50065913
My Commission Expires 8/10/2022

15 Thornton Road • Oakland, NJ 07436 • USA
CollagenMatrix.com • 201-405-1477

SYMBIOS® Collagen Membrane SR

INSTRUCTIONS FOR USE

Description: SYMBIOS® Collagen Membrane SR (slow resorbable) is a white, nonfriable membrane matrix engineered from highly purified type I collagen fibers derived from bovine Achilles tendon. SYMBIOS® Collagen Membrane SR is resorbable which eliminates the need for a second surgical procedure that is normally required to remove a non-resorbable membrane. SYMBIOS® Collagen Membrane SR is supplied sterile, non-pyrogenic and for single use only.

Indication for Use: SYMBIOS® Collagen Membrane SR is a bioresorbable, implantable collagen material that is intended for use in dental surgery procedures as a material for placement in the area of dental implant, bone defect or ridge reconstruction to aid in wound healing post dental surgery.

Contraindications: SYMBIOS® Collagen Membrane SR Collagen is contraindicated in patients who have:

- acute infections or contaminated wound in the oral cavity,
- known allergy to collagen of animal origin or other bovine-derived products,
- clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which in the physician's judgment, will prevent safe implantation or likely healing.

Warning: Clinicians should use care in screening their patients for any known allergies to collagen or bovine-derived products. Hypersensitivity reactions have been noted with the use of other products containing bovine collagen; therefore, the possibility exists of developing a local sensitivity response to SYMBIOS® Collagen Membrane SR.

Precautions: As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroidal therapy or currently taking anticoagulants. Patients with clinically significant systemic diseases, indicating a history of anaphylactic reactions, autoimmune diseases, uncontrolled diabetes or severe hypertension have not been implanted with the membrane; therefore, the safety and effectiveness for those patients have not been determined. Nor has it been evaluated in pregnant women, children and/or in patients with conditions involving extremely severe defects with little periodontium or bone. SYMBIOS® Collagen Membrane SR cannot be resterilized. Open, unused membranes must be discarded. In vivo stability may be adversely affected if re-sterilized. Cross-contamination and infection may occur if re-used. Do not use if the product sterilization barrier or its packaging is compromised.

Safety: The product is manufactured from bovine Achilles tendon, which is classified as tissues with no detected infectivity for Bovine Spongiform Encephalopathy, BSE (World Health Organization Guidelines). The bovine tendon is known to be one of the richest sources of type I collagen that is commercially available. The manufacturing process for the product meets European and International Standards for animal tissue sourcing, handling and inactivation of Spongiform Encephalopathy (SE) pathogens. This process involves a treatment with sodium hydroxide that is a recognized method of inactivation of SE pathogens.

A viral inactivation study for the product's manufacturing process was conducted by an independent laboratory. In this study key manufacturing steps were evaluated for their ability to inactivate the following viral strains: Bovine Viral Diarrhea (enveloped virus) and Porcine Parvoviridae (nonenveloped virus). The study results showed that each of the manufacturing steps evaluated, including the sodium hydroxide treatment, is effective in inactivating these viruses.

Adverse Reactions: Possible complications that can occur with any dental surgery include infection, swelling of the intraoral tissue, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing or dehiscence, root resorption or ankylosis, loss of crestal bone height, pain, or complications associated with the use anesthesia. Minor discomfort may occur for a few days.

Application: SYMBIOS® Collagen Membrane SR is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The membrane should be removed from the inner pouch with sterile gloves or instruments. The periodontal defect or bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). The clinician should perform thorough debridement of the defect. Space-making material such as autologous bone, demineralized bone matrix and ceramic materials may be used to fill the defect. As much tissue as possible should be preserved to allow for primary closure of the wound and correct positioning of the flaps without tension.

OSCIFU-OU 2014-09 REV DRAFT (RUSSIA)

SYMBIOS® Collagen Membrane SR can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution for approximately five minutes prior to the final placement. Then, the device may be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors. The membrane should overlap the walls of the defect by at least 2 mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material. Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a non-cutting needle. Resorbable tacks can also be used to affix the membrane. Primary closure must be achieved and maintained to avoid accelerated resorption due to membrane exposure.

Post-operative Procedures: Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent. The patient should refrain from brushing the treated area for two days following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated. The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.


SYMBIOS® Collagen Membrane SR should be completely resorbed 26 to 38 weeks following surgery. However, probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

How Supplied: One (1) membrane per package; Two (2) packages per box.

Available Sizes: 15 mm x 20 mm
20 mm x 30 mm
30 mm x 40 mm

Additional Country Specific Storage and Product Information: SYMBIOS® Collagen Membrane SR do not contain any dangerous hazardous substances that would require special handling. SYMBIOS® Collagen Membrane SR can be disposed of in accordance with all local, state, or federal regulations.

Store at room temperature (15 °C/59 °F – 30 °C/86 °F) in a dry place (relative humidity about 10% to 80%). The shelf life of the device is 36 months from the date of manufacture. The manufacturer guarantees the operability and maintenance of the sterility of the products during the entire shelf life, subject to storage conditions. The product can be transported by any means of transport, such as covered railway wagons, containers, freight vehicles, ship holds, aviation modules in accordance with the rules for these to modes of transport, in compliance with the conditions of transportation without risk of damage to the package. The product does not affect the environment during operation. Minor excursions from the recommended temperature range during storage, shipping, and handling are acceptable. The manufacturer cannot guarantee the product to meet specifications, if stored at temperatures below -10°C / 14°F or above 49°C / 120°F for more than 10 days.

Labeling Symbols:

Manufacturer

Use-by date

Do not re-use

Caution: See instructions for use

Batch number

Catalogue number

Do not use if package is damaged

Temperature limitation

Sterilized using irradiation

Do not re-sterilize

Rx Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist or physician.

Authorized representative in the European Community

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MDSS GmbH