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## CERTIFICATE

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China Council for the Promotion of International Trade  
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of ZHEJIANG NUROTRON BIOTECHNOLOGY CO., LTD. on the annexed OPERATION MANUAL ENDURO COCHLEAR IMPLANT SYSTEM WITH ACCESSORIES is genuine.

China Council for the Promotion  
of International Trade

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证明书查询网 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

NUROTRON

# Operation Manual Enduro Cochlear Implant System with accessories

LI FANGPING, President  
Zhejiang Nurotron Biotechnology Co., Ltd.

«Zhejiang Nurotron Biotechnology Co., Ltd.»  
No.17 Longtan Road, Yuhang District, HangZhou, Zhejiang, 311121, China

## Table of Contents

|                                                                    |    |
|--------------------------------------------------------------------|----|
| Mode of Operation .....                                            | 6  |
| System Design .....                                                | 6  |
| 2.1 Signal Processing Unit .....                                   | 7  |
| 2.1.1 Ear Hook .....                                               | 7  |
| 2.1.2 Status Indicator Lamp .....                                  | 7  |
| 2.1.3 Microphone .....                                             | 8  |
| 2.1.4 Microphone Cover .....                                       | 8  |
| 2.1.5 Microphone Dust Protection Membrane .....                    | 8  |
| 2.1.6 Radio Frequency Signal Output Interface .....                | 8  |
| 2.1.7 The Volume/Program Selection Switch .....                    | 8  |
| 2.2 Transmitting Coil, Transmitting Wire and Outer Magnet .....    | 9  |
| 2.2.1 Transmitting wire .....                                      | 10 |
| 2.2.2 Transmitting coil .....                                      | 10 |
| 2.2.3 Outer magnet .....                                           | 10 |
| 2.3 BTE rechargeable battery compartment .....                     | 11 |
| 2.3.1 Connector of the rechargeable battery compartment .....      | 11 |
| 2.3.2 Lock of the battery compartment .....                        | 11 |
| 2.3.3 Decorative panel .....                                       | 11 |
| 2.3.4 Specification and Model .....                                | 12 |
| 2.3.5 Way of Charging .....                                        | 12 |
| 2.4 Remote control .....                                           | 12 |
| 2.4.1 Power Key .....                                              | 13 |
| 2.4.2 Volume Key .....                                             | 13 |
| 2.4.3 Touch Screen .....                                           | 13 |
| 2.4.4 Touch Keys .....                                             | 13 |
| 2.4.5 Micro USB Charging Interface .....                           | 13 |
| 2.4.6 Specification and Model .....                                | 14 |
| 2.4.7 Charging Mode .....                                          | 14 |
| 3. Model and Scope of Delivery .....                               | 14 |
| 4. Serial Number .....                                             | 15 |
| 5. Basic product specifications .....                              | 16 |
| 5.1 Basic Product Features .....                                   | 16 |
| 5.2 Operating Environment .....                                    | 17 |
| 5.2.1 Cochlear Implant—Operating Environment .....                 | 18 |
| 5.3 Transportation and Storage Conditions .....                    | 18 |
| 5.3.1 Cochlear Implant—Transportation and Storage Conditions ..... | 18 |
| 5.4 Main Performance Characteristics .....                         | 18 |
| 5.5 Shelf Life and Service Life .....                              | 19 |

|                                                                                                                     |    |
|---------------------------------------------------------------------------------------------------------------------|----|
| Scope .....                                                                                                         | 19 |
| Purpose—Intended Use .....                                                                                          | 19 |
| 2 User and Patient Category .....                                                                                   | 19 |
| 4 Conditions of Use .....                                                                                           | 20 |
| Indications and Contraindications .....                                                                             | 20 |
| 7.1 Contraindications .....                                                                                         | 20 |
| 7.2 Indications for Use .....                                                                                       | 20 |
| 8. Adverse Events .....                                                                                             | 20 |
| Bacterial infectious diseases, inflammatory response .....                                                          | 20 |
| 9. Application .....                                                                                                | 21 |
| 9.1 Assembly and disassembly the BTE rechargeable battery compartment and signal processing unit .....              | 21 |
| 9.2 Change of the BTE hook .....                                                                                    | 22 |
| 9.3 Connecting the Transmitting Coil .....                                                                          | 22 |
| 9.4 Installing the Transmitting Coil .....                                                                          | 23 |
| 9.5 Power On/Off the Speech Processor .....                                                                         | 24 |
| 9.5.1 Power On/Off by Pressing a Key .....                                                                          | 24 |
| 9.5.2 Power On/Off Using the Remote Control .....                                                                   | 24 |
| 9.6 Figure 10. Low battery prompt of the BTE rechargeable battery compartment .....                                 | 25 |
| 9.7 Realization of the remote control features .....                                                                | 26 |
| 9.7.1 Connection with the speech processor .....                                                                    | 26 |
| 9.7.2 Power On/Off the Speech Processor .....                                                                       | 26 |
| 9.7.3 Remote Control Instruction Manual .....                                                                       | 26 |
| 9.7.4 Management of the remote control features .....                                                               | 29 |
| 10. Maintenance, cleaning and repair .....                                                                          | 32 |
| 10.1 Storing a Speech Processor and the Remote Control .....                                                        | 32 |
| 10.2 Requirements for cleaning .....                                                                                | 33 |
| 10.3 Keeping dry .....                                                                                              | 33 |
| 10.4 Disinfection and Sterilization .....                                                                           | 33 |
| 11. Troubleshooting .....                                                                                           | 33 |
| 12. Reference Information .....                                                                                     | 34 |
| 13. Safeguard and Preventive Measures .....                                                                         | 35 |
| 13.1 Replacing Dust Protection Membranes .....                                                                      | 36 |
| 13.2 Precautions for using rechargeable battery compartments .....                                                  | 36 |
| 13.3 Other Preventive Measures .....                                                                                | 36 |
| 13.4 Compatibility Requirements .....                                                                               | 37 |
| 14. Information about the Manufacturer and the Manufacturer's Authorized Representative in the European Union ..... | 37 |
| 15. After-Sales Service Information .....                                                                           | 37 |
| 15.1 Scope of the Service .....                                                                                     | 37 |
| 15.2 Information about the Subdivisions Providing Services .....                                                    | 38 |
| 15.3 Address for Consumers' Requests on the Territory of the Russian Federation .....                               | 38 |

|                                                                                                         |    |
|---------------------------------------------------------------------------------------------------------|----|
| Marking symbols .....                                                                                   | 38 |
| Warranty Liabilities .....                                                                              | 39 |
| Environmental Requirements, Recycling or Disposal Information .....                                     | 40 |
| Reference to List of Standards Used During the Development, Production and Control of the Product ..... | 40 |
| 10. EMC Requirements .....                                                                              | 41 |
| 20.1 Electromagnetic Emission .....                                                                     | 41 |
| 20.2 Electromagnetic Immunity .....                                                                     | 42 |
| 20.3 Isolation Distance .....                                                                           | 45 |

**Warning:**

Implant recipients under the age of 14 or those with intellectual disabilities are not allowed to perform all of the procedures described in this Manual, therefore speech processor must be operated by their guardian/caregiver. A guardian/caregiver is allowed to connect a speech processor, change batteries, perform routine cleaning and maintenance in accordance with this Operation Manual. When the system is running, please do not clean or perform any maintenance of the system.

Implant recipients of 14 years of age or older, literate and having normal intelligence, are allowed to operate the speech processor on their own. The recipients of the implant are its intended users which are allowed to connect the speech processor, change batteries, perform routine cleaning and maintenance according to this Operation Manual. When the system is running, please do not clean or perform any maintenance of the system.

This Operation Manual applies to the Enduro Cochlear Implant System and is suitable for routine patient use.

## Mode of Operation


Cochlear Implant System consists of a cochlear implant speech processor (hereinafter abbreviated as 'speech processor'), a cochlear implant (hereinafter abbreviated as 'implant') and NuroSound\* setting software. The speech processor comprises a signal processing unit, a transmitting wire, a transmitting coil, and battery compartments. It also contains a remote control for signal processing units (hereinafter abbreviated as 'remote control'), a connecting cable and an ear hook. The speech processor is the external part of the Cochlear Implant System, the Cochlear Implant is the implanted part. The remote control can wirelessly control the speech processor, and the speech processor operating environment, and the system battery use environment; the signal processing unit can collect, analyze and process ambient acoustic signals, encode the processed acoustic signals, send them to the implant in the form of radio frequency (RF) signals, and use an implant electrode array to stimulate residual auditory nerves in patients with severe to serious perceptive deafness in order to enable them to perceive sound.

\*Removable media with special software is intended only for a medical specialist audiologist and is not issued to the patient, the user of the Enduro cochlear implantation system with accessories. The software is installed on the specialist's computer once and is used for the initial connection and configuration of the Enduro C accessories cochlear implantation system. Only an audiologist specialist can use removable media with special software during the initial connection of the system and in the future, if necessary. Patients do not use removable media with special software.

Removable media with special software is used during the initial connection of the system and can be used in the future, if necessary, for other patients with the same system.

## 2. System Design

A cochlear implant speech processor is a non-implantable part of a Cochlear Implant System, consisting mainly of a signal processing unit, a transmitting wire, a transmitting coil, a behind-the-ear (BTE) battery compartment, an on-body battery compartment, a remote control etc. Figure 1 shows a behind-the-ear speech processor:


1 Transmission Coil; 2 Transmitting wire; 3 Signal processing unit; 4 BTE rechargeable battery compartment; 5 Remote control

**Figure 1. Behind-the-Ear Type**

## Signal Processing Unit

Figure 2 is a configuration scheme of the signal processing unit. The microphone (7) picks up acoustic signals, and the signal processing unit encodes them to send an output signal to the RF output interface (6). The ear hook (1) can be put on the user's ear, the status indicator lamp (hereinafter abbreviated as 'indicator lamp') (2) displays the operating status of the speech processor, the operating buttons (4 and 5) can be used for volume control, program switching and power on/off command, the dust protection membrane (3) prevents dust and sweat from entering the microphone, and the interface (8) connects the battery compartment and the connecting wire of the compartment.

1 Ear hook; 2 Status indicator lamp; 3 Microphone dust protection membrane; 4 Key “+” for tuning the volume and switching programs; 5 Key “-” for tuning the volume and switching power on/off; 6 RF signal output interface; 7 Microphone; 8 Wired battery compartment interface and settings

**Figure 2. Components of the Signal Processing Unit**

### 2.1.1 Ear Hook

The radian design of the ear hook allows you to have a more comfortable wearing experience and makes the device capable to meet the ear pavilion requirements of most adults or children.

### 2.1.2 Status Indicator Lamp

The status indicator lamp displays the operating status of the speech processor. In Figure 2, the position of Pointer 2 represents the position of the status indicator lamp.

- The green indicator goes on: After the signal processing unit has been connected to a power source, or after using the volume tuning key, or after using the program selection key, or after using the start and stop key, the indicator lamp goes on.

The green indicator lamp goes out after slowly flashing three times: The system initializes upon the start up of the signal processing unit; after reaching stable operation, the indicator lamp goes out.

- The green indicator lamp flashes continuously and slowly: The coil is not right over the implant,

the signal processing unit is searching for the implant. Once the connection is restored, the indicator lamp will flash slowly three times and then go out to switch back to the normal operation status. When the battery capacity is running below the specified value, the alarm prompt starts flashing.

- The green indicator lamp flashes continuously and fast: There are no debug parameters in the signal processing unit, or the debug parameter is incorrect, or there are other possible errors.

## 2.1.3 Microphone

Used primarily to receive signals.

## 2.1.4 Microphone Cover

The microphone cover is primarily used to protect the keys on the signal processing unit. Air vents have been left around the microphone area to make it easier to pick up the acoustic signal.

## 2.1.5 Microphone Dust Protection Membrane

The dust protection membrane of the microphone is located beneath the air vents and is mainly used to prevent dust and moisture from entering the microphone.

## 2.1.6 Radio Frequency Signal Output Interface

The RF signal of the signal processing unit outputs through the interface, and the RF signal is transmitted to the transmitting coil through the transmitting wire connected to the interface.

## 2.1.7 The Volume/Program Selection Switch

The “+” and “-” keys on the signal processing unit (as shown in Figure 3) can be used for volume control, program switching, viewing the running program and switching power on/off, and can feature various functions while different keys being pressed.

Sound volume control function: It is used to tune the sound volume. You can enable this function by briefly pressing the “+” key and the “-” key. After briefly pressing the “+” key, the indicator lamp will flash once, which indicates that the tuning has been successful and the sound volume has been increased; after briefly pressing the “-” key, the indicator lamp will flash once, which indicates that the tuning has been successful and the sound volume has been reduced; sound volume level means Levels 1 to 12. At power-up, the initial sound volume is set to Level 10.


Program switching function: This function enables to switch between different signal environments and processing programs downloaded to the signal processing unit chip from a computer and refers to four programs, including Program 1, Program 2, Program 3 and Program 4. The implant recipient may choose the best program that suits them.

The implant recipient can adjust the running program by pressing and holding the “+” key for 3 seconds, until the indicator lamp flashes once, indicating that the adjustment has been successful; Program 1 is set

original default program.

**the running program:** This function enables to view the running program. Press both the “+” key and “-” key at the same time and release them after the indicator lamp goes on. When the indicator lamp goes on for the first time, it means that the operation has been successful. The number of times the indicator lamp flashes after that corresponds to the number of the running program of the signal processing unit. For example, flashing once indicates that Program 1 is running; flashing twice indicates that Program 2 is running; the rest can be done in the same way.

**The on-off of power supply function.** It is used to power on/off the signal processing unit. You can do this by pressing and holding the “-” key for a while. The indicator lamp will flash to indicate that the operation has been completed successfully.


1 Key Sound Volume “+”; 2 Key Sound Volume “-”

**Figure 3. Keys**

## 2.2 Transmitting Coil, Transmitting Wire and Outer Magnet

The transmitter coil is used primarily for transmitting and receiving RF signals; the transmitting wire bonds the transmitting coil and the signal processing unit; an outer magnet is attracted to a magnet in the implant's receiving coil on the scalp to ensure the proper connection between the implant and the transmitting coil. The design of the transmitting coil and transmitting wire is shown in Figure 4:

1 Connector with 4 slots; 2 Transmitting wire; 3 Transmitting coil; 4 Outer magnet; 5 Male connector

**Figure 4. Scheme of the Transmitting Coil and the Transmitting Wire**

## Transmitting wire

The transmitting wire comes in 5 types of length: 6 cm; 8 cm; 10 cm; 20 cm, recommended for children; 30 cm, recommended for adults; one end of the transmitting wire is connected to the signal processing unit via a four-slot connector; the other end is connected to the transmitting coil via a male connector.

### 2.2.2 Transmitting coil

The transmitting coil must be placed right over the implant. It is attracted by an outer magnet to a magnet in the implant's receiving coil on the scalp to ensure that the transmitting coil is secured in the correct position.

The transmitting coil is bonded to the transmitting wire via a male connector so that it can be replaced by inserting and removing it.

### 2.2.3 Outer magnet

The outer magnet is of the clasp-on type. A clasp-on type pivot magnet is shown in Figure 5. The relative position of the outer magnet and the transmitting coil remains the same after being pivoted and secured. The combined retractive force can be changed by replacing the outer magnet with other models.

The outer magnet has 6 features, including ©-© and ®, marked on the inside of the outer magnet as shown in Figure 5. If the performance of the outer magnet is inadequate, replace the magnet as shown in Figure 5. Rotate the device counterclockwise, to pop the magnet out and insert an outer magnet with suitable features on the transmitting coil.


Figure 5. The Inside of the Outer Magnet


#### Warning

Choking hazard—the device comprises small parts, please supervise babies and toddlers under 3 years of age while they are using the device.

## BTE rechargeable battery compartment

The cochlear implant speech processor is powered by the BTE - a rechargeable battery compartment that contains a rechargeable lithium battery cell. The rechargeable battery compartment contains a 210mAh lithium-ion cell. When the battery is fully charged, the normal lifespan is 6-10 hours. The specific structure is shown in Figure 6:


1 Connector; 2 Lock; 3 Decorative panel

**Figure 6. BTE rechargeable battery compartment**

### 2.3.1 Connector of the rechargeable battery compartment

The battery compartment connector is the interface between the battery and the signal processing unit, having an interface, a silicone rubber O-ring and a metal pin plate, as shown in Figure 7:

**Figure 7 Connector of the rechargeable battery compartment**

### 2.3.2 Lock of the battery compartment

The battery compartment lock is used to connect and lock the battery compartment and the signal processing unit, and to install and remove the battery. When the lock is closed, the battery compartment is secured to the signal processing unit; when the lock is opened, the battery compartment is detached from the signal processing unit and is ready to be removed. As it is shown in Figure 7, the position marked by Pointer 2 corresponds to the locked and closed state. If you need to open the lock, press firmly on the place indicated by Pointer 2.

### 2.3.3 Decorative panel

The decorative panel for BTE battery compartments is located on the front of the rechargeable battery compartment (see item 3 in Figure 7) to give the unit a fancy look.

## Specification and Model

rechargeable battery compartment can be recharged and reused. The nominal voltage is 3.7 V, and its capacity is 210 mAh.

### 3.5 Way of Charging

When a low power prompt appears, please remove the battery compartment and connect it to the charging device to load. Once the battery compartment is fully charged, take it out of the charging device and connect it to the signal processing unit for use. The way of charging is shown in section 8.4 of the Manual.

## 2.4 Remote control

The remote control model is non-contact, and the control is an internal electric drive product with no applied parts. The remote control can perform the following functions:

- Operation control function: It can be paired with a specific signal processing unit. After successful connection, it is able to control the processing unit including volume control, program selection and input mode, mode selection, going to hibernation mode, etc.
- It shows the operating status of the signal processing unit: the volume, the program, charge of the battery compartment and the remote control charge and also displays other related information on the LCD screen.
- It shows the operating condition of the implant, checking the connection between the processing unit and the implant.

A scheme of the general structure of the remote control is shown in Figure 8:


Figure 8. Diagram of the External Structure of the Remote Control

## **Power Key**

Power key is located on the upper right side of the remote control case and has the following functions:

- Press and hold the power key for 5 seconds: In this case, the turning on and off function of the control panel can be performed. During startup, the screen turns on and a welcome interface is displayed to confirm that the remote control has turned on. During shutdown, a shutdown operation is performed in which the screen goes black. A shutdown prompt window is displayed on the screen. After clicking "OK", the screen turns off to identify that the remote control has turned off.
- Press and hold for 1 second: the function of unlocking or manually locking the screen will be enabled; when the remote control is in standby mode and the screen is off, you need to press the power key instantly to display the unlocking interface and unlock the device according to the prompts on the screen; when the remote control screen is on, you need to press the power key instantly so that the screen turns off and the remote control goes to the standby mode.

### **2.4.2 Volume Key**

The volume key is located on the right side of the remote control case; the Volume "+" key is at the top and the Volume "-" key is underneath, they enable the following functions:

- Volume "+" Key: Press the "+" key once so that the volume is increased to Level 1 after the remote control has been successfully connected with the speech processor and the speech processor has gone into its normal operating state.
- Volume "-" Key: Press the key once so that the volume is reduced to Level 1 after the remote control has been successfully connected with the speech processor and the speech processor has gone into its normal operating state.

### **2.4.3 Touch Screen**

A 2.4 inch (60.96 mm) capacitive screen is used.

### **2.4.4 Touch Keys**

Left Function Key: flips the left page/ back to the previous level;

Home Key: displays the main interface of the remote control;

Right Function Key: flips the right page or shows the "Advanced" interface on the display interface of the running program.

### **2.4.5 Micro USB Charging Interface**

The Charging Interface is located on the right side of the remote control case, under the volume control key, it uses the micro USB interface;

ideal conditions, when fully charged once, a standby time of 20 days can be achieved.

### Specification and Model

remote control model is contactless; the water resistance level is IP22; the capacity is 500 mAh.

#### 4.7 Charging Mode

When the remote control is underpowered, please engage the charging wire of the remote control to charge it. While charging the device, a lightning bolt is displayed on the battery symbol at the top right of the screen, and the remote controller can be used normally while charging; when charging is complete, the lightning bolt disappears. When it is impossible to charge the remote controller or the internal battery needs to be replaced, please contact Nurotron for replacement. One may not dismantle or replace it without permission.

## 3. Model and Scope of Delivery

The name of the Cochlear Implant System is Enduro.

The Enduro Cochlear Implant System includes:

I. The Enduro Cochlear Implant System, comprising:

1.1 A cochlear implant, 1 pc., design options: CS-10 implant; CS-10A implant with TS electrode; CS-10A implant with TM electrode; CS-10A implant with TL electrode.

1.2 An Enduro speech processor, comprising: a signal processing unit; a transmitting coil; a transmitting wire, design options: 6 cm, 8 cm, 10 cm, 20 cm, 28 cm (no more than 3 pcs.); BTE battery compartment (no more than 2 pcs.); on-body battery compartment; connecting wire of the on-body rechargeable battery compartment; ear hook (no more than 5 pcs.); remote control.

1.3 A removable drive with special software (if necessary), comprising: a programmer; a micro USB cable; NuroSound software on a USB drive; a configuration cable; a cover.

1.4 Operation Manual

II. Accessories:

1. A charger

2. A charging cable<sup>1)</sup>

3. A configuration cable

4. A signal monitor<sup>2)</sup>

5. A replaceable plate for the signal processing unit (no more than 3 pcs.)

6. A splitter adapter

7. A speech processor holder, design options: large, small (no more than 2 pcs.)

8. A waterproof cover
9. A drying device<sup>3)</sup>
10. A cover to wear the device attached to clothes
11. A BTE rechargeable battery compartment with USB
12. A four-pin splitter adapter
13. A BTE hook
14. Magnet No. 1
15. Magnet No. 2
16. Magnet No. 3
17. Magnet No. 4
18. Magnet No. 5
19. Magnet No. 6
20. Magnet No. 8


<sup>1, 2, 3)</sup> Registration certificate No. RZN 2017/5398 dated 13.08.2019

## 4. Serial Number

The serial number is the only ID symbol for the speech processor manufactured by Nurotron. The location of the mark is shown in Figure 9. The mark should be read as follows:

- a. The naming principles for the marking of the signal processing unit: ED (refers to the Enduro system) 1 (refers to the signal processing unit) 0A (version) 0001 (production serial number starting from 0001) X (security ID number); IP44 (waterproof level);
- b. Method of labeling of the BTE rechargeable battery compartment: ED (refers to the Enduro system) 2 (refers to the BTE battery compartment) 0A (version) DDMYY (date of manufacture); IP44 (waterproof level);
- c. The naming principles for the marking of the transmitting coil: ED (refers to the Enduro system) 3 (refers to the transmitting coil) 0A (version) 0001 (production serial number starting from 0001); IP44 (waterproof level);
- d. The naming principles for the marking of the remote control: ED (refers to the Enduro system) 5 (refers to the remote control) 0A (version) 0001 (production serial number starting from 0001); IP22 (waterproof level);
- e. The naming principles for the marking of the transmitting coil magnet:

There is only one Arabic numeral on the inside of the outer magnet that refers to the magnet strength. The magnet strength is indicated by ©-© and ®. The higher is the number, the stronger is the magnet.

a. Location of the signal processing block marking; b Location of BTE rechargeable battery compartment marking; c. Location of the transmitter coil marking; d. Location of remote control marking; e. Location of the transmitter coil magnet marking.

**Figure 9. Description of the Location of the Main Marking**

## 5. Basic product specifications

### 5.1 Basic Product Features

Table 5-1. Main Enduro system Features

| No. | Name                   | Feature          |                            | Notes |
|-----|------------------------|------------------|----------------------------|-------|
| 1   | Signal processing unit | Outer dimensions | 41.1 mm x 25.3 mm x 9.4 mm |       |
|     |                        | Weight           | No more than 5 g           |       |

|   |                                      |                     |                                 |  |
|---|--------------------------------------|---------------------|---------------------------------|--|
|   | Transmitting coil                    | Coil diameter       | 31.3 mm                         |  |
|   |                                      | Weight              | No more than 9 g                |  |
|   | Transmitting wire                    | Length              | 6 cm, 8 cm, 10 cm, 20 cm, 28 cm |  |
| 5 | BTE rechargeable battery compartment | Outer dimensions    | 19.0 mm x 31.7 mm x 9.83 mm     |  |
|   |                                      | Weight              | No more than 0.48 g             |  |
| 6 | Remote control                       | Outer dimensions    | 94.8 mm x 48.8 mm x 9.9 mm      |  |
| 7 | NuroSound software,                  | Version             | version 1.0.9                   |  |
|   |                                      | Software risk class | Software safety class - B.      |  |

Table 5-2 The main features of the cochlear implant

| No. | Implant                             | CS-10A                  | CS-10A with TS electrode                                                                               | CS-10A with TM electrode                                                                               | CS-10A with TL electrode                                                                               |
|-----|-------------------------------------|-------------------------|--------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|
| 1   | Engineering data                    | Standard electrode (CO) | A thin medium electrode (TS)                                                                           | A thin medium electrode (TM)                                                                           | A thick long electrode (TL)                                                                            |
| 2   | Quantity of contacts                | 24                      | 22                                                                                                     | 24                                                                                                     | 24                                                                                                     |
| 3   | Electrode contact shape             | Saddle bend             | Saddle bend                                                                                            | Saddle bend                                                                                            | Saddle bend                                                                                            |
| 4   | Theoretically exposed electrode, mm | 0.32 mm <sup>2</sup>    | 0.16 mm <sup>2</sup> (apical end)<br>0.22 mm <sup>2</sup> (middle)<br>0.28 mm <sup>2</sup> (basal end) | 0.16 mm <sup>2</sup> (apical end)<br>0.22 mm <sup>2</sup> (middle)<br>0.28 mm <sup>2</sup> (basal end) | 0.16 mm <sup>2</sup> (apical end)<br>0.22 mm <sup>2</sup> (middle)<br>0.28 mm <sup>2</sup> (basal end) |
| 5   | Electrode contact spacing           | 0.85 mm                 | 0.71 mm                                                                                                | 0.85 mm                                                                                                | 1.0 mm                                                                                                 |
| 6   | Length of the array                 | 20.5 mm                 | 16.0 mm                                                                                                | 20.5 mm                                                                                                | 24.0 mm                                                                                                |
| 7   | Apical tip shape                    | Conical shape           | Conical shape                                                                                          | Conical shape                                                                                          | Conical shape                                                                                          |
| 8   | Implantation method                 | Forceps or a fork       | Forceps or a fork                                                                                      | Forceps or a fork                                                                                      | Forceps or a fork                                                                                      |

Table 5-3. Material Safety Data Sheets

| No. | Name             | Data on material                                                                                             |
|-----|------------------|--------------------------------------------------------------------------------------------------------------|
| 1   | Cochlear implant | Platinum-iridium alloy (90% Pt-platinum, 10% Ir-iridium); Titanium; Medical silicone gel and silicone rubber |

## 5.2 Operating Environment

Ambient temperature: +5 °C to +50 °C; Relative humidity: ≤ 85 %.

Atmosphere pressure: 60 kPa ~ 106 kPa

Since the highest temperature in the operating environment is 50 °C, the surface temperature of the device is influenced by the operating temperature, but it will not rise higher than 50.5 °C.

### **Cochlear Implant—Operating Environment**

Ambient temperature from +5 °C to +50 °C, humidity no more than 85 %. After implantation, it withstands temperatures from +32 °C to +42 °C and exposure to blood and other body fluids.

## **5.3 Transportation and Storage Conditions**

Ambient temperature: -20 °C to +50 °C; Relative humidity: ≤ 85 %.

### **5.3.1 Cochlear Implant—Transportation and Storage Conditions**

The sterilized implant must be stored and transported within a humidity range of 15 % to 75 % and a temperature range of 0 °C to +43 °C.

## **5.4 Main Performance Characteristics**

Effective operating distance of the transmitting coil: 2 mm ~ 10 mm;

Effective operating distance of the remote control: 0 m ~ 1.5 m;

The range of response to sound intensity: in the range of 35 dB SPL to ~80 dB SPL, the amplitude of the stimulation pulse and the strength of the incoming sound signal show an increasing response. If the sound intensity exceeds 80 dB SPL, automatic suppression is attained at a stimulation strength of 80 dB SPL;

Signal processing frequency range: 150 Hz ~ 7.8 kHz;

Electrostatic protection: meets the standard EN45502-2-3:2010: The device is able to withstand a 2 kV contact discharge 10 times and a 8 kV air discharge 5 times.

An internal power source and an additional part of the BF type (the working part is the transmitting coil);

Execution mode: Continuous operation;

Waterproof level: remote control: IP22, other parts: IP44.

### **Cochlear Implant:**

Output signal range: 0 μA ~ 1900 μA

Pulse duration: 4 μs ~ 425 μs

Pulse Width Adjustment Scheme: Pulse width can be manually adjusted for individual MAPs through the user interface in the configuration software during the implant recipient's MAP setup.

Single channel stimulation pulse frequency range: 270 Hz – 890 Hz

all stimulation pulse frequency: 20 kHz

transmission bandwidth: 1.067 Mb/s

forward transmission frequency: 16 MHz

backward transmission frequency: 16 MHz

Distance during remote test: Skin graft thickness is 2 mm to 10 mm

## 5.5 Shelf Life and Service Life

The total shelf life of the Enduro Cochlear Implant System is 2 years and is determined by the cochlear implant expiration date. Implant packaging allows maintaining it sterile for 2 years. The service life of the implant is twelve years.

The shelf/service life of the Enduro speech processor is 5 years.

## 6. Scope

Enduro cochlear implant system is a device that digitally processes sound and delivers electrical stimulation to the auditory nerve through a cochlear implant to restore sound perception. The cochlear implant speech processor is a non-implantable part of the Cochlear Implant System, used together with a cochlear implant manufactured by Nurotron; during setup, the speech processor is used in together with the NuroSound configuration software and the programmer. The Enduro cochlear implant system shall be applied in patients aged 12 months and older suffering from severe to significant binaural neurosensory hearing loss. The system can be applied in many places such as home, school, and work unit.

### 6.1 Purpose—Intended Use

The Enduro Cochlear Implant System provides electrical stimulation (through a cochlear implant) to the residual auditory nerve in the cochlea so that the brain of patients with severe to significant neurosensory deafness could perceive sound.

### 6.2 User and Patient Category

The Enduro Cochlear Implant System is designed to be used by physicians and speech audiologists. It shall be applied in patients aged 12 months and older with severe to significant binaural neurosensory hearing loss. Patients under 14 years of age or those with intellectual disabilities are not allowed to perform all operations on their own, but only under the supervision of a guardian/caregiver.

## Conditions of Use

The Enduro Cochlear Implant System is installed and configured by physicians and speech audiologists only in healthcare facilities. The system can be used in many places such as home, school and workplace.

## 7. Indications and Contraindications

### 7.1 Contraindications

Absolute contraindications are:

- 1) A serious inner ear pathology such as Michel deformity or absence of a cochlea;
- 2) Absence of the auditory nerve;
- 3) Deafness caused by a non-cochlear lesion;
- 4) Severe mental disorder;
- 5) Uncured mastoid abscess of the eardrum;

Relative contraindications are:

- 1) Poor physical condition;
- 2) Seizures that can not be controlled.

### 7.2 Indications for Use

- 1) Grade 4 sensorineural hearing loss
- 2) Grade 4 borderline sensorineural hearing loss
- 3) Total hearing loss
- 4) Deafness with residual hearing

## 8. Adverse Events

Bacterial infectious diseases, inflammatory response.

## Application

During the initial setup (power up), the audiologist can assemble your signal processing unit according to the following steps:

- 1) Connect the signal processing unit to the BTE rechargeable battery compartment;
- 2) The transmitting wire shall be connected to the transmitting wire slot on the signal processing unit through a four-way connector for transmitting wires; the transmitting wire shall be connected to the transmitting coil via a male connector.
- 3) Select an outer magnet with the proper strength, and install it on the appropriate location of the transmitting coil;
- 4) Wear the signal processing unit and the transmitting coil for a while and adjust it a convenient position;
- 5) Power on the speech processor and select the appropriate program;
- 6) Adjust the sound volume.

In the course of everyday use, the remote control can be used for the following operation on the speech processor:

- 1) Power on/off the speech processor;
- 2) Recognition of the connection between the speech processor and the remote control;
- 3) Switching programs;
- 4) Volume control;
- 5) Switching the input mode;
- 6) Switching modes;
- 7) Monitoring of the battery charge;
- 8) Detecting the state of the connection.

## 9.1 Assembly and disassembly the BTE rechargeable battery compartment and signal processing unit

Assembly and disassembly of the BTE battery compartment and the signal processing unit are shown in Figure 10:

**Assembly:** Align the signal processing unit connector with the corresponding battery compartment connector and gently push it back to front to make the connectors join.

**Disassembly:** On the top of the decorative plate, press the lock vertically inward and pull the battery compartment back to ensure they are separated from each other.


Figure 10. Assembly and disassembly the BTE rechargeable battery compartment

## 9.2 Change of the BTE hook

As shown in Figure 11, hold the signal processing unit with one hand and the ear hook with the other. Bring the two pieces apart to remove the ear hook; firmly press the backup ear hook to its home position in the opposite direction to that of the removal.



Figure 11. Method for Replacing the BTE Hook

## 9.3 Connecting the Transmitting Coil

Installing the transmitting wire: As shown in Figure 12, hold the signal processing unit with one hand and firmly grip the four-pin transmitting wire connector with two fingers of the other hand. Direct the latter at the right position and insert downward to ensure that there is no gap between the two.

**Figure 12. Connection Between the Transmitting Wire and the Signal Processing Unit**

oving the transmitting wire: As shown in Figure 13, hold the signal processing unit with one hand and gently grip the four-pin transmitting wire connector with two fingers of the other hand. Pull the connector gently upward; do not pull on the wire too much or reverse the four-pin transmitting wire connector.

**Figure 13. Removing the Transmitting Coil and the Signal Processing Unit**

## 9.4 Installing the Transmitting Coil

**Installing the transmitting coil:** As shown in Figure 14, hold the transmitting coil with one hand and grip the transmitting wire pin connector with two fingers of the other hand. Direct the latter at the right position and insert downward to ensure that there is no gap between the two.

**Removing the transmitting coil:** As shown in Figure 15, hold the transmitting coil with one hand and the transmitting wire pin connector with the other hand. Pull the connector out gently.


Figure 14. Installing the Transmitting Coil


Figure 15. Removing the Transmitting Coil

Removing the transmitting coil: As shown in Figure 15, hold the transmitting coil with one hand and the transmitting wire pin connector with the other hand. Pull the connector out gently; do not pull on the wire too much or reverse the four-pin transmitting wire connector.

## 9.5 Power On/Off the Speech Processor

Perform the speech processor power on/off operation using the following two methods:

### 9.5.1 Power On/Off by Pressing a Key

When the signal processing unit is off, you need to press and hold the “-” key on the signal processing unit for 3 seconds, the status indicator lamp will flash 3 times, which indicates that the operation has been completed successfully and the speech processor has turned on;

When the signal processing unit is on, you need to press and hold the “+” key on the signal processing unit for 3 seconds, the status indicator lamp will flash once, which indicates that the operation has been completed successfully and the speech processor has turned off.

### 9.5.2 Power On/Off Using the Remote Control

Under conditions where the remote control has been successfully connected to the speech processor, use the right touch key and select the “BTE Power” (“Battery Power”) module in the remote control and select “Power on” or “Power off”. The progress is shown in green to indicate that the operation has been completed successfully.


Note: if the implant has been connected untimely or unsuccessfully, the speech processor will automatically turn off within approximately 2.5 minutes after turning on the power.

**Figure 10. Low battery prompt of the BTE rechargeable battery**

### partment

The BTE battery compartment can be recharged and reused. When the status indicator lamp on the signal processing unit is flashing continuously and slowly and/or is accompanied by an audible alarm from the signal processing unit, it means that the battery is going into a low power mode. If you notice this, please remove the battery and connect it to the charging holder to load. Once the battery is fully loaded, take it out of the holder and connect it to the signal processing unit for use.

During the loading process, the green battery compartment indicator lamp must be on. When the battery is fully loaded, the indicator lamp will go out. During the loading process, the battery compartment cannot be used. To replace the battery compartment, please press the button for battery insertion and removal; do not pull it out by force. The charging method and the charging status of the BTE battery compartment are shown in Figure 16 and Figure 17, respectively:

1. A charger 2. BTE rechargeable battery compartment  
**Figure 16. Charging method of the BTE rechargeable battery compartment**


Figure 17. Charging method of the BTE rechargeable battery compartment

## 9.7 Realization of the remote control features

### 9.7.1 Connection with the speech processor

- Tap the right touch key to enter the right interface, click and enter the “Settings” - “BTE” (“BTE rechargeable battery compartment”) - “connect” and then click on “Conn1-” to connect to the battery compartment, provided the remote control is turned on and the speech processor is powered;
- If the connection is successful, the left dialog window “Conn1” will display “connected”; otherwise, it will display “unconnected”.

### 9.7.2 Power On/Off the Speech Processor

- Start: Tap the right touch key to enter the right interface, click and enter the “BTE Power” (“Power on of the BTE rechargeable battery compartment”) interface, provided the remote control is powered on and the speech processor is powered: ① If the connection is successful, then click on “Power On” to perform the start-up operation for the speech processor using the remote control.

### 9.7.3 Remote Control Instruction Manual

#### 1) Start of the remote control

During the start the following prompt appears: “The state is powering ... 5/4/3/2/1”. After a successful start, a prompt is displayed stating that “The state is power on.” (“Status - power is on.”); ② if connection failed, click “Power On” and then the BTE interface “Connect” will be displayed) and a prompt will be displayed

ing that the user is performing a power on/off operation for the battery compartment only after the operation is done.

## Power-off of the remote control

During normal use of the speech processor, enter the “BTE Power” interface, click on “Power Off” and perform a power-off operation for the speech processor via the remote control. During the power-off process, a prompt is displayed: “The state is powering down ...”; after a successful power-off, a prompt is displayed that “The state is power off.

## 3) Detecting the connection state

- Provided the remote control is powered on and the speech processor is powered, touch the right touch key to enter the right interface, click and enter the “BTE Power” (“Power on of the BTE rechargeable battery compartment”) interface. In the lower position you can select the function “Detect State.” In some states, if the transmitter coil is lowered, this option can be used to detect abnormalities when wearing the transmitter coils, and the “Implant Fault. Please Check the Implant Status” pop-up window is displayed for the user to check the status After the failure has been corrected, click “Dismiss”.
- If the connection is normal and there are no failures, a prompt is displayed.

## 4) Volume control

### ① Volume key

- During normal speech processor operation, you can use the volume keys on the right side of the remote control to increase or decrease the sound volume of the currently used program.
- When you use the side volume key to adjust the sound volume, the sound volume prompt window is displayed and the current sound volume level is displayed; if the operation is not performed within 3 seconds, the prompt window will automatically disappear.

### ② “Volume” interface

- During a normal speech processor operation, touch the left touch key to enter the left interface, click and enter the “Volume” interface. Click on “+” or “-” to increase or decrease the sound volume.
- The sound volume adjustment range is 1-12. The higher the audio volume level is, the greater the intensity of the audio volume is.

While performing audio volume switching, enter “Volume” interface, read the current audio volume from the speech processor and display the related information; click on “+” or “-” the first time, change the sound volume, and send information into the speech processor; read the reading from the speech processor again (the displayed sound volume is the modified sound volume).

## 5) Switching programs

- In normal speech processor operation, touch the left touch key to enter the left interface, then click

Enter the "Programs" interface to switch the speech processor programs through the control panel.

There are 4 programs for switching. "Program 1" - "Program 4" can proportionally correspond to 4 programs from left to right when downloading programs in NuroSound.

When switching programs, select any program from "Program 1" - "Program 4" to perform the switching.

- When switching programs, enter the "Programs" interface, read the current program from the speech processor and display the related information; click on the plug-in program, send it to the speech processor and display the related information.

## 6) C-tone strategy function

- Under normal speech processor operation, touch the left touch key to enter the left interface, click and enter the "Programs" interface; under the current program, touch the right touch key to enter the "Advance" interface.
- In "C-Tone" column click on the switch key to start a C-tone function.

## 7) Additional noise reduction function

- Under normal speech processor operation, touch the left touch key to enter the left interface, click and enter the "Programs" interface; under the current program, touch the right touch key to enter the "Advance" interface.
- In the "Noise Reduce" column, click on the right switch button to launch the additional noise reduction function.

## 8) Switching the input mode

- Under normal speech processor operation, touch the left touch key to enter the left interface, click and enter the "Input" interface so that the remote control can switch the input mode.
- There are 4 input modes: Mic/audio input (Mic/AUX), mic + external input (Mic + EXT), mic + magnetic induction input (Mic + T-coil) and monitoring. They are switchable. In addition, in "Mic + EXT" and "Mic + T-coil" modes, select the ratio of inputs according to your actual requirements (the default setting is that the ratio of the inputs is 1:3).
- While switching the input mode, enter the "Input" interface, read the input mode of the current program from the speech processor, and display the related information; the first click can change the input mode to send information to the speech processor; after switching programs (for example, Program 1 switches to Program 4), perform the return operation.

This time what the "Input" interface displays is the original input mode of the program as opposed to the input mode switched by the remote control. This indicates that the user operation for the login mode is only temporarily used with the current program, is not stored in the speech processor.

### Switching modes

Under normal speech processor operation, touch the left touch key to enter the left interface, click and enter the "Mode" interface so that the remote control can switch the input modes.

There are 2 modes available for switching (Normal and Rich).

When switching modes, enter the "Mode" interface, read the current mode from the speech processor, and display related information; the first click can change the mode to send information to the speech processor; after switching programs (for example, Program 1 switches to Program 4), perform the return operation and read the speech processor again. What the interface "mode" displays is the original mode of the program as opposed to the mode switched by the leg by the remote control. This suggests that the user operation for the mode is only temporarily used by the current program, not stored in the speech processor.

## 9.7.4 Management of the remote control features

### 1) Power on/off

- Press the power key on the top of the remote control for a long time, display the welcome interface, and then release the power key to perform the startup operation for the remote control; after startup, the content of the display is the main interface by default.
- When the remote control is in startup state and the screen is on, press the power key on the top for a long time until the "Power Off?" dialog window pops up; click on "OK" to turn off the remote control.

### 2) Locking and unlocking the screen

- Briefly press the power key on the top of the remote control to perform a manual screen lock operation.
- Screen auto-lock time is 45 seconds.
- After the remote control screen is locked, briefly press the power key on the top, display unlock interface, and slide on the interface from left to right to unlock the remote control.

### 3) Battery charge indicator

- During normal speech processor operation, touch the right touch key to enter the right interface, click and enter the "Battery" interface, and view the charge of the BTE battery and the remote control battery, as well as the specific battery charge values.
- The charge level can also be displayed on the left (indication of the battery charge level of the BTE battery compartment) and the right (indication of the battery charge level of the remote control) parts of the status bar. However, specific battery charge level values cannot be obtained.

Under normal speech processor operation, touch the right touch key to enter the right interface, click and enter the "Settings" - "Remote" interface to modify the following settings:

## 5) Date and time

Click and enter the “Date & Time” interface, and use “+” and “-” to change the date and time.

Exit after completing the change so that the main interface will immediately display the changed date and time.

## 6) Language

- Click and enter the “Language” interface; click on the appropriate language to switch between languages.

- There are 3 languages available for selection (Chinese, English and Spanish), however Spanish is not included at this time.

- Exit after completing the change so that all interfaces immediately display the changed language.

## 6) Background picture

- Click and enter the “Wallpaper” interface; click on the green “<<” and “>>” to find the wallpaper, define a color for the wallpaper, and click on “Set wallpaper” below the wallpaper with that color to change the wallpaper.

- There are 5 background colors to choose from (dark green, blue, red, yellow and pink).

- Exit after completing the change so that the interface will immediately display the changed wallpaper.

## 7) Reset

- Click on “Reset” so that the remote control immediately performs a factory reset.

## 8) Charge of the remote control

- When the battery level of the remote control is less than 10%, a low battery prompt appears.

There is a micro USB interface on the right side of the remote control, and you can use the corresponding USB cable to charge.

- During the charge, “ ” is displayed in the power display window on the right side of the status bar to indicate that charging is in progress.

- After charging is complete, a prompt window appears, click on “close” and unplug the cable.

## 9) Reboot

- If the remote control shuts down without being able to restart, or is slow during the use, the reboot function is used to force the remote control to return to normal.

- Reboot operation: Press the power key and the “+” key from the side volume keys at the same time, a “shutdown or not” prompt window will pop up. At this time, you do not need to click the shutdown key. Continue pressing the Power key and the “+” key from the side volume keys, and release the key after the remote control screen turns off. Then, the screen turns on which indicates that the reset operation was

ful.

### The status bar and main control page

The screen is divided into two sections, namely the status bar and the main control page, as shown in the following figure:

In the figure: “Status bar” and “Main control page”  
Figure 18 Screen of the main control page

### 11) Memory after changing own parameters

- Change your own parameters of the remote control (date/time, language and wallpaper), and notice if the remote has a memory for changing parameters under different conditions.
- When the remote control is in the unlocked state or turned off, the change of all own parameters is memorized.
- When the remote control is in the reset state, the wallpaper change is not memorized, but the date/time and language changes are memorized.

The detailed information is provided in the following table:

| State of the remote control |                                                            | Memorized or not |          |                    |
|-----------------------------|------------------------------------------------------------|------------------|----------|--------------------|
| Rank                        | Description                                                | Date and time    | Language | Background picture |
| Unlocking                   | Locking after the screen lock                              | ✓                | ✓        | ✓                  |
| Shutdown                    | Restart after shutdown                                     | ✓                | ✓        | ✓                  |
| Reset                       | Automatic shutdown after insufficient charge               |                  |          |                    |
|                             | Reset                                                      |                  |          |                    |
|                             | Reboot                                                     |                  |          |                    |
|                             | Disconnection of the power supply, and the battery removal |                  |          |                    |

## 10. Maintenance, cleaning and repair

To ensure safe and proper operation of the system, it is essential to check its condition daily. It is prohibited to alternate software or components of this system in any way.

Some checks require the use of special tools and/or imply certain risks for the user and therefore require special training. If necessary, contact your physician/audiologist in charge or Nurotron representatives in your region to carry out maintenance and inspections that require special qualifications.

Do not alternate this equipment without permission from Nurotron. Do not replace any parts without authorized technicians, as this may result in the system being unusable. Do not put the battery in an electronic dryer.

Only specialists authorized by Nurotron are allowed to repair individual components of the system during the warranty period. Nurotron shall not accept any responsibility for malfunctions, damage or loss caused by transportation, alteration or repair made by other personnel not authorized by Nurotron.

### 10.1 Storing a Speech Processor and the Remote Control

When the speech processor is not in use, please disconnect the speech processor and the rechargeable battery compartment, and wipe the speech processor, transmission cable, transmission coil and remote control with a clean cloth; please place the speech processor and the remote control in a dedicated electronic dryer or other sealed desiccant boxes.

Wet environment can cause a short circuit in the speech processor and remote control circuitry preventing normal operation. If the speech processor will not be used for a long time, please disconnect the speech processor and the battery compartment and place the speech processor in a sealed plastic bag; if the remote control will not be used for a long time, please turn it off.

Since the remote control has a built-in lithium-ion battery, it must be charged at regular intervals to keep the battery active (the recommended charging frequency is once every 7–10 days).

When the dust protection membrane for the microphones on the key cover plate of the signal processing unit becomes dirty or the sound quality decreases, please contact Nurotron to have it timely replaced.

## Requirements for cleaning

Keep the battery contacts between the speech processor and the battery compartment clean. Wipe the contacts with a cotton ball or a cotton swab dipped in alcohol. Please do not use too much force to avoid damaging the contacts.

When cleaning the surface of the speech processor, wipe it gently with a piece of cloth dampened in a mild detergent solution and avoid penetration of any liquid into the speech processor. Clean the surface frequently as this can effectively prevent dust from entering the speech processor.

When cleaning the surface of the remote control, wipe it gently with a piece of cloth moistened in a mild detergent and avoid penetration of any liquid into the charging port. Wipe the surface frequently, which can effectively prevent dust from entering the remote control.

## 10.3 Keeping dry

If the speech processor is not used at night, it is recommended to wipe it, remove the rechargeable battery compartment and place it in a dryer or other sealed desiccant box for drying (not included in the set).

If you are not going to use the remote control at night, it is recommended that you turn it off, clean the surface and place it in a sealed box with desiccants. Drying frequency:

1) Electronic Dryer: It is recommended to dry the device once every two days (2–3 hours each time). In a damp environment or if the user sweats a lot, the frequency can be increased accordingly.

2) Daily drying: If you are not going to use the device at night, it is recommended that you place it in a sealed box with desiccants.

Warning: If chemical drugs are used as desiccants, please keep them out of the reach of children to prevent children from swallowing these chemical drugs by mistake.

## 10.4 Disinfection and Sterilization

The Enduro Cochlear Implant System is not subject to disinfection or sterilization.

The implant is provided sterile. It undergoes an Ethylene Oxide (EO) sterilization process during its sterile packaging to achieve a Guaranteed Sterility Level (GUS) of  $10^{-6}$  before release.

## 11. Troubleshooting

If the following events occur, your Cochlear Implant System may be malfunctioning:

- No sound is heard when the speech processor is in the startup state;

The sound is interrupted or of poor quality, or it makes the user feel uncomfortable;

Audible alarm sounds;

The status indicator lamp displays a warning prompt;

The remote control indicates that the connection has been unsuccessful;

After entering information via the remote control, there is no response from the speech processor.

In case of troubles, please follow the operating steps; do not apply physical force to any part. After completing each step, please check if you can hear sound.

Restart the remote control and perform the reconnection operation again.

If related troubles still arise after your check, please contact your audiologist or Nurotron.

## 12. Reference Information

| Reference Information | Failure Symptom                                                                   | Check and Solution                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-----------------------|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1                     | The transmitting coil is loose and may fall down.                                 | <ul style="list-style-type: none"> <li>• Use a transmitting coil with a stronger magnetic field.</li> <li>• Cut the connection at the attachment point of the transmitting coil in the appropriate manner.</li> <li>• Contact your audiologist or Nurotron.</li> </ul>                                                                                                                                                                                                                                                                                                                                                    |
| 2                     | No sound is heard                                                                 | <ul style="list-style-type: none"> <li>• Check the position of the transmitting coil, then reboot the system.</li> <li>• Check if the speech processor and power supply are engaged.</li> <li>• Check the position of the volume control.</li> <li>• Replace battery compartments with new ones, reboot the system. Make sure the charge is sufficient.</li> <li>• Check the wire, replace if necessary.</li> <li>• Use a cotton ball soaked in alcohol to clean the metal contacts (metal battery compartment contacts and signal processing unit connector)</li> <li>• Contact your audiologist or Nurotron.</li> </ul> |
| 3                     | The signal processing unit status indicator lamp flashes slowly and continuously. | <ul style="list-style-type: none"> <li>• There is no program in the signal processing unit, or incorrect settings have been made, or there are other possible errors.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 4                     | Status display Fast and continuous flashing of the signal processing unit lamp.   | <ul style="list-style-type: none"> <li>• There is no program in the signal processing unit, or incorrect settings have been made, or there are other possible errors.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| 5                     | Audible alarm sounds.                                                             | <ul style="list-style-type: none"> <li>• The battery compartment charge level is below the prescribed value. Please replace the compartment.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| 6                     | You hear noise                                                                    | <ul style="list-style-type: none"> <li>• Make sure the connection between the transmitting coil and the signal processing unit is correct and that the transmitting location is suitable.</li> <li>• Replace the batteries with new ones, reboot the device.</li> </ul>                                                                                                                                                                                                                                                                                                                                                   |
| 7                     | Sound is fuzzy and weak or of poor quality.                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |

|   |                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|---|-------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   |                                                       | <ul style="list-style-type: none"> <li>• Make sure the volume control is in the correct position.</li> <li>• Check the microphone on the signal processing unit for foreign particles or damage.</li> <li>• Contact your audiologist or Nurotron</li> <li>• Make sure the microphone is not obstructed by other objects.</li> <li>• Replace the battery compartments with new ones, reboot the system. Check whether the poor sound quality was caused by insufficient battery compartment power.</li> <li>• Check whether the microphone cover was dirty or damaged. Replace if necessary.</li> <li>• Contact your audiologist or Nurotron</li> </ul>         |
| 8 | The remote control has been connected unsuccessfully. | <ul style="list-style-type: none"> <li>• Check if the charge level of the remote control is sufficient. If the charge level is insufficient, then please charge it in a timely manner.</li> <li>• Reduce the distance between the remote control and the signal processing unit, and ensure that the length of the connection is less than 1.5 m. Repeat the pairing operation.</li> <li>• Make sure the top of the remote control faces the speech processor and repeat the pairing operation.</li> <li>• Reboot the controller, and repeat the above two operating steps.</li> <li>• Please consult the debugging medical specialist or Nurotron.</li> </ul> |

### 13. Safeguard and Preventive Measures

You can wear the speech processor outdoors in hot or cold weather. In very cold environments, you can wear a hat or wear a scarf over your speech processor. In a very hot environment, please wipe off sweat in time to prevent it from entering the signal processing unit, which will lead to unwanted electrical fault. If you are not wearing the speech processor, please do not place it in hot places (for example, on a windowsill or in a car left in a sunspot).

Before you go bathing or swimming, please remove your speech processor. If a small volume of water or other liquid has accidentally entered into the speech processor, please turn off the processor immediately, wipe it or dry it with an appropriate drying device and do not use it for the next 12 hours. If the processor is still not working after 12 hours, please contact our after-sales personnel.

Please remove your speech processor before applying makeup, powder or hair spray or similar products on your hair, as these substances may damage the signal processing unit.

The signal processing unit must not be exposed to dust and sand. If any dust or sand gets in, please make every effort to shake it out. In serious cases, please contact our after-sales personnel. Do not open the signal processing unit yourself.

## Replacing Dust Protection Membranes

When the dust protection membrane for the microphones on the key cover plate of the signal processing unit becomes dirty or the sound quality decreases, please contact Nurotron to have it timely replaced.

## 13.2 Precautions for using rechargeable battery compartments

These measures, described below, can make the battery compartments demonstrate their optimal performance:

Place the backup battery compartment in a sealed plastic bag to ensure that the metal contact of the battery compartment cannot come into contact with other metal objects such as coins. The above method can effectively prevent abnormal battery discharge, heating or explosion of the battery. Place the battery compartment in a cold and dark environment.

Do not heat. For example, do not place the battery compartment on windowsills or in a car left in a sunspot.

Do not immerse the rechargeable battery compartment in water.

Do not apply excessive force to deform. Do not apply excessive force to assemble the rechargeable battery compartment or drop it on a hard floor.

The rechargeable BTE battery compartment uses a lithium-ion cell. Stay away from open flames. It can explode in an open flame.

In case of leakage, keep your skin or eyes from contact with the escaping liquid. If this liquid comes into contact with your skin or eyes, rinse them immediately with clean water and go to a hospital for treatment. After disposing, please wash your hands immediately.

Do not place rechargeable battery compartments in an electronic dryer.

## 13.3 Other Preventive Measures

When the equipment attains its end-of-life date, it must be disposed of in accordance with the regulations of the local environmental department, or the equipment must be sent to Nurotron company for disposal.

The electromagnetic power emitted by the system of the cochlear implant is very low, therefore, it does not affect the operation of other equipment; in a strong electromagnetic field, it is recommended that you turn off the power of the cochlear implant speech processor.

This medical device is not designed to work in an environment with a high oxygen content.

## 4 Compatibility Requirements

Do not engage the Enduro Cochlear Implant System to other equipment and accessories other than the equipment and accessories provided or supplied by Nurotron. Please charge the batteries with the charger specified by Nurotron.

During the implantation operation, common surgical instruments that meet the requirements should be used. Responsibility for the surgical procedure and the techniques, instruments and techniques used lies with the performing surgeon.

## 14. Information about the Manufacturer and the Manufacturer's Authorized Representative in the European Union

Zhejiang Nurotron Biotechnology Co., Ltd. (Ltd.)

No.17 Longtan Road, Yuhang District, Hangzhou, Zhejiang, 311121, China

Hotline of the National Service: +086-4006-333-571

Website: [www.nurotron.com](http://www.nurotron.com)

MedPass International Limited

Windsor House, Bretforton, Evesham, Worcestershire, WR117JJ, UK

## 15. After-Sales Service Information

### 15.1 Scope of the Service

All buyers of Nurotron products are entitled to the after-sales service provided by Nurotron, including:

- 1) Providing after-sales service specialists to carry out the relevant works.
- 2) Explaining to the buyer general information about Nurotron cochlear products and how to handle them.
- 3) Buyer advice on follow-up professionals and rehabilitation centers; assistance in evaluation of postoperative speech recuperation.
- 4) Timely furnishing customers with updated information on Nurotron cochlear products.
- 5) Filling out a buyer's file, periodic monitoring of the supply of information to them, assistance in

g up requests for information.

Assisting buyers in coordinating with surgical clinics.

Gratuitous repair of Nurotron cochlear products during the warranty period and maintenance during the entire service life of the product.

- 8) Providing a service for the replacement of elements.
- 9) Providing customers at their expense with spare parts and parts for cochlear products of the Nurotron company.

## 15.2 Information about the Subdivisions Providing Services

- 1) After-sales Service Department: Zhejiang Nurotron Biotechnology Co., Ltd Ltd.)
- 2) Address: No.17 Longtan Road, Yuhang District, Hangzhou, Zhejiang, 311121, China

## 15.3 Address for Consumers' Requests on the Territory of the Russian Federation



OOO Van Mourik Medical

109316, Russia, Moscow, Volgogradsky prospect, 47, floor 2, room 234.

+7 (499) 346 08 06

office.vmm@mail.ru

## 16. Marking symbols

|                                                                                     |               |                                                                                     |                                                     |
|-------------------------------------------------------------------------------------|---------------|-------------------------------------------------------------------------------------|-----------------------------------------------------|
|  | Manufacturer  |  | Manufacturing Date                                  |
|  | Serial Number |  | Refer to the instructions for use                   |
|  | CE mark       |  | Authorized representative in the European Community |

|                                                                                     |                             |                                                                                     |                                                                                              |
|-------------------------------------------------------------------------------------|-----------------------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
|    | Fragile, handle carefully   |    | Keep away from sunlight                                                                      |
|    | Keep away from moisture     |    | The working part of type BF                                                                  |
|    | Up                          |    | Working part of type B                                                                       |
|    | Humidity range              |    | Symbol requiring to dispose of the equipment according to the Electronic Equipment Directive |
|    | Use before ....             |    | Temperature limit                                                                            |
|    | Batch Code                  |    | Sterilized with ethylene oxide                                                               |
|   | Do not re-sterilize         |   | Do not use if package is opened or damaged                                                   |
|  | Prohibition of repeated use |  | Quantity                                                                                     |
|  | MRI incompatible device     |  | Degree of protection against penetration of water and solid particles                        |
|  | Electromagnetic field       |                                                                                     |                                                                                              |

## 17. Warranty Liabilities

The Nurotron company, as a manufacturer, is responsible for the quality of the product, subject to the conditions listed below.

The manufacturer shall not be liable for any damage or failures caused by misuse of the product. Correct use implies:

- Compliance with the instructions described in the operating documentation
- Adherence to a proper system maintenance program
- Compliance with international safety standards

The system is warranted against any defects in material and workmanship in accordance with the table

from the date of delivery.

| Product                                              | Warranty Period |
|------------------------------------------------------|-----------------|
| Cochlear implant CS-10A                              | 12 years        |
| Cochlear implant CS-10A with TS electrode            | 12 years        |
| Cochlear implant CS-10A with TM electrode            | 12 years        |
| Cochlear implant CS-10A with TL electrode            | 12 years        |
| Signal processing unit                               | 5 years         |
| Magnet                                               | 1 year          |
| BTE battery compartment                              | 1 year          |
| On-body battery compartment                          | 1 year          |
| Transmitting coil                                    | 1 year          |
| Transmitting wire                                    | 1 year          |
| Connection cable for the on-body battery compartment | 1 year          |
| Remote control                                       | 1 year          |
| Charger                                              | 1 year          |
| Charging cable                                       | 1 year          |
| Signal monitor                                       | 1 year          |
| Splitter adapter                                     | 1 year          |
| BTE battery compartment with USB                     | 1 year          |
| 4-pin splitter adapter                               | 1 year          |

Contact your local Nurotron representative for information about your warranty.

The total guaranteed storage life of the Enduro Cochlear Implant System is 2 years and is determined by the cochlear implant guaranteed storage life.

## 18. Environmental Requirements, Recycling or Disposal Information

During its use, transportation and storage, the device does not have any negative impact on humans and the environment.

It is recommended to keep the packaging material for transportation or repair. If you wish to dispose of the packaging material, please return it to a recognized collection center for recycling.

The device shall be disposed of in accordance with the sanitary rules and regulations of the Russian Federation SanPiN 2.1.3684-21. Hazard class of medical waste for cochlear implant is B, other system and package components belong to hazard class A.

## 19. Reference to List of Standards Used During the Development, Production and Control of the Product

EN ISO 14971:2012; EN ISO 10993-1:2009/AC2010; EN ISO 10993-3:2009; EN ISO 10993-5:2009; EN

10993-6:2010; EN ISO 10993-10:2010; EN ISO 10993-11:2009; IEC 60601-1:2005+CORR.1:2006+CORR.2:2007+A1:20124; EN 60601-1-2:2015; EN 60601-1-11:2010; IEC 60601-1-2:2006/AC:2008; EN 62366:2008; EN ISO 15223-1:2012; EN 1041:2008; ISO 14708-7: 2013; EN 60601-1:2009.

## 20. EMC Requirements

When using this equipment, special precautions related to electromagnetic compatibility (EMC) shall be taken, the equipment shall be installed and operated in accordance with the EMC information provided in the manual.

This equipment may be affected by portable and mobile RF communications equipment.

Except from wires (converters) sold as spare parts for internal components, the use of unspecified accessories and wires (converters) may result in increased emissions and reduced immunity of the equipment or system.

This product or system shall not be operated or stored close to other equipment. If the equipment or system is to be located near or in the same rack with other equipment, the ability of the equipment or system to operate normally in such a configuration shall be checked and verified.

The following wires shall be used to comply with EMC and immunity requirements.

| Wire name         | Length                          |
|-------------------|---------------------------------|
| Transmitting wire | 6 cm, 8 cm, 10 cm, 20 cm, 28 cm |

The basic characteristics shall meet the following conditions: after the immunity test, the maximum DC displacement current due to asymmetry caused by the cochlear implant output waveform shall be less than 0.1  $\mu$ A; the output current range shall be 0 mA to 1.9mA, and the fluctuation range shall not exceed 20 %.

### 20.1 Electromagnetic Emission

Table 20-1. Manufacturer's Guidance and Declaration. Electromagnetic Emissions: Applicable to All Equipment and Systems.

| Manufacturer's Guidance and Declaration: Electromagnetic Emissions. The equipment is expected to be operated in the following specified electromagnetic environment. The buyer or user must ensure that the equipment is operated under the following conditions: |                     |                                        |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------------------------------|
| Radiation test                                                                                                                                                                                                                                                    | Compliance with the | Electromagnetic environment—guidelines |

|                                                                          |               |                                                                                                                                                                                                                                  |
|--------------------------------------------------------------------------|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Radio frequency emission<br>GB 4824 (CISPR 11)                           | Group 1       | This equipment uses radio frequency energy only for its internal functions. Therefore, its RF emissions are very low and there is very little probability of any interference with electronic equipment in the surrounding area. |
| Radio frequency emission<br>GB 4824 (CISPR 11)                           | Class B       | The equipment can be used in all facilities, including low voltage household and utility networks directly connected to residential buildings for domestic use.                                                                  |
| Harmonic radiation<br>GB 17625.1                                         | Class A       |                                                                                                                                                                                                                                  |
| Fluctuations of voltage / Pulsating radiation GB 17625.2 (IEC 61000-3-3) | Extensibility |                                                                                                                                                                                                                                  |

## 20.2 Electromagnetic Immunity

Table 20-2. Manufacturer's Guidance and Statement. Protection Against Electromagnetic Interference. Applicable to All Equipment.

| Manufacturer's Guidance and Declaration. Electromagnetic Immunity. The equipment is expected to be operated in the following specified electromagnetic environment. The buyer or user must ensure that the equipment is operated under the following conditions: |                                                        |                                              |                                                                                                                                                             |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------|----------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunity testing                                                                                                                                                                                                                                                 | IEC60601 investigation level                           | Compliance level                             | Electromagnetic environment - guidelines                                                                                                                    |
| Electrostatic discharge (ESD) GB/T 17626.2 (IEC61000-4-2)                                                                                                                                                                                                        | ± 6 kV for contact charge<br>± 8 kV for for air charge | contact<br>± 6 kV for air<br>± 8 kV          | The floor shall be made of wood, concrete or ceramic tiles. If the floor is covered with composite materials, the relative humidity shall be at least 30 %. |
| • Electrical fast transient (burst) immunity test GB/T 17626.4 (IEC61000-4-4)                                                                                                                                                                                    | ± 2 kV for power cord<br>± 1 kV for I/O line           | ± 2 kV for power cord<br>± 1 kV for I/O line | The mains power supply shall be of adequate quality for use in a typical commercial or hospital environment.                                                |
| Surge GB/T 17626.5 (IEC61000-4-5)                                                                                                                                                                                                                                | ± 1 kV line - ± 2 kV line-ground                       | ± 1 kV line -<br>± 2 kV line-ground          | The mains power supply shall be of adequate quality for use in a typical commercial or hospital environment.                                                |

|                                                                                                                                                                                                              |                                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>drop, short<br/>cut-off and<br/>variation &lt;5 %<br/>0.5 cycle duration<br/>relative to <math>U_T</math>, &gt; 95 %<br/>voltage drop)</p> <p>power supply line<br/>GB/T 17626.11<br/>(IEC61000-4-11)</p> | <p>40 % <math>U_T</math>, last 5 cycles<br/>(relative to <math>U_T</math>, 60 %<br/>voltage drop)</p> <p>70 % <math>U_T</math>, last<br/>25 cycles<br/>(Relative to <math>U_T</math>, 30 %<br/>voltage drop)</p> <p>&lt;5 % <math>U_T</math>, last 5 sec<br/>(Relative to <math>U_T</math>, &gt;<br/>95 % voltage drop)</p> | <p>&lt; 5 % <math>U_T</math>, 0.5<br/>cycle duration<br/>(relative to<br/><math>U_T</math>, &gt; 95 %<br/>voltage drop)</p> <p>40 % <math>U_T</math>, last<br/>5 cycles (relative<br/>to <math>U_T</math>, 60 %<br/>voltage drop)<br/>70 % <math>U_T</math>, last<br/>25 cycles<br/>(relative to <math>U_T</math>,<br/>30 % voltage<br/>drop)<br/>&lt;5 % <math>U_T</math>, last 5<br/>sec<br/>(Relative to<br/><math>U_T</math>, &gt; 95 %<br/>voltage drop)</p> | <p>The mains power supply shall be of<br/>adequate quality for use in a typical<br/>commercial or hospital environment. If<br/>the user needs to ensure uninterrupted<br/>operation of the equipment during a<br/>power failure, it is recommended to<br/>operate the equipment from an<br/>uninterruptible power supply or battery<br/>power.</p> |
| <p>power frequency<br/>magnetic field<br/>(50/60 Hz), GB/T<br/>17626.8 (IEC<br/>61000-4-8)</p>                                                                                                               | <p>3 A/m</p>                                                                                                                                                                                                                                                                                                                | <p>3 A/m</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                      | <p>The power frequency magnetic field shall<br/>have the power frequency magnetic field<br/>level appropriate for typical areas in a<br/>typical commercial or hospital<br/>environment.</p>                                                                                                                                                       |
| <p>Note: <math>U_T</math> is the AC mains voltage before the test voltage is applied.</p>                                                                                                                    |                                                                                                                                                                                                                                                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                                                                                                                                    |

10-3. Manufacturer's Guidance and Declaration. Protection Against Electromagnetic Interference. Supplementary Equipment and Systems

Manufacturer's Guidance and Declaration. Electromagnetic Immunity. The equipment is expected to be operated in the following specified electromagnetic environment. The buyer or user must ensure that the equipment is operated under the following conditions:

| Immunity testing                                                                                                        | IEC60601 investigation level                                                  | Compliance level                | Electromagnetic environment—guidelines                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF conductivity<br><br>GB/T 17262.6<br>(IEC61000-4-6)<br>Radio frequency<br>Radiation<br>GB/T 17262.3<br>(IEC61000-4-3) | 3 V (RMS<br>value)<br>150 kHz –<br>80 MHz<br><br>3 V/m<br>80 MHz –<br>2.5 GHz | 3 V (RMS<br>value)<br><br>3 V/m | During use, the distance between a portable or mobile RF device communication equipment and any part of it, including wires, shall not be shorter than the recommended isolation distance. The distance shall be calculated by the formula corresponding to the frequency of the transmitter.<br>Recommended insulation distance<br>$d = 1.2 \sqrt{150 \text{ kHz}-80 \text{ MHz}}$<br>$d = 1.2 \sqrt{80 \text{ kHz}-800 \text{ MHz}}$<br>$d = 2.3 \sqrt{800 \text{ MHz}-2.5 \text{ GHz}}$ ,<br>where: P — based on the maximum rated output capacity of the transmitter as provided by the transmitter manufacturer, expressed in watts (W);<br>d — the recommended isolation distance, expressed in meters (m).<br>The field strength of a fixed transmitter shall be determined by an electromagnetic field study a, and in each frequency range it shall be below the immunity compliance level b.<br>Interference may occur in the area surrounding the equipment marked with the following symbol.<br> |

Note 1. The formula for the higher frequency range is used when the frequency is 80 MHz and 800 MHz.  
 Note 2. These guidelines may not apply in all circumstances. The propagation of electromagnetic waves is affected by absorption and reflection from buildings, objects and human bodies.

a Field strengths from a fixed transmitter (such as wireless (cellular/radio) phones, mobile radio base stations, amateur radio, AM & FM broadcasting, television and radio) cannot be predicted theoretically with accuracy. To assess the electromagnetic environment of a fixed RF transmitter, it is necessary to consider the study of electromagnetic fields. If the measured field strength in the area where the equipment is located is higher than the applied relevant RF level that has been mentioned above, it is necessary to check whether the equipment is able to operate normally. If abnormal operation is observed, additional measures may be required (for example, changing the position or orientation of the equipment).  
 b. Over the entire frequency range 150 kHz to 80 MHz, field strengths shall be below 3 V/m.

## Isolation Distance

20-4. Recommended Isolation Distances Between Equipment and Portable or Mobile RF Communications Equipment: Applicable to Equipment and Systems Not Involved in Life Support Operations

| Recommended Isolation Distance Between Portable and Mobile RF Communications Equipment. The equipment is intended to be used in an electromagnetic environment where RF radiation-induced disturbance is controlled. In accordance with the maximum rated output capacity of the communications equipment, the buyer or user can prevent electromagnetic interference keeping a minimum isolation distance between portable and mobile RF communications equipment (the transmitter) and the equipment. |                                                                         |                |                 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------|----------------|-----------------|
| The maximum rated output capacity of the transmitter in watts                                                                                                                                                                                                                                                                                                                                                                                                                                           | Isolation distance corresponding to different transmitter frequencies/m |                |                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 150 kHz 80 MHz                                                          | 80 MHz 800 MHz | 800 MHz 2.5 GHz |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 0.12                                                                    | 0.12           | 0.23            |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     | 0.38                                                                    | 0.38           | 0.73            |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 1.2                                                                     | 1.2            | 2.3             |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 3.8                                                                     | 3.8            | 7.3             |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     | 12                                                                      | 12             | 23              |
| To obtain a maximum rated output capacity not listed in the above table, the recommended isolation distance d (expressed in meters (m)) can be determined according to the formula in the corresponding column indicating the frequency of the transmitter. In this case, P is the maximum rated output capacity of the transmitter as provided by the transmitter manufacturer, expressed in watts (W).                                                                                                |                                                                         |                |                 |
| Note 1: The formula for the higher frequency range is used when the frequency is 80 MHz and 800 MHz.                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                         |                |                 |
| Note 2: These guidelines may not apply in all circumstances. The propagation of electromagnetic waves is affected by absorption and reflection from buildings, objects and human bodies.                                                                                                                                                                                                                                                                                                                |                                                                         |                |                 |

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