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## CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade  
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of LIFETECH SCIENTIFIC (SHENZHEN) CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion  
of International Trade

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Authorized

Signature:

Hu Na

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(Date: Apr. 18, 2023)

证明网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

CE  
0344

Approved By

*Xie Yuehui*  
CEO, Xie Yuehui

09.03.2023

Aegisy



Instructions for use  
Aegisy Vena Cava Filter with Delivery System

Lifetech Scientific (Shenzhen) Co., Ltd.

2023

V4.0

|                                                                                     |                                                                                  |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------|
|    | Date of manufacture                                                              |
|    | Manufacturer                                                                     |
|    | Use by                                                                           |
|    | Serial number                                                                    |
|    | Batch code                                                                       |
|    | Sterilized using ethylene oxide                                                  |
|    | Non-pyrogenic                                                                    |
|    | Do not reuse                                                                     |
|    | Consult Instructions For Use                                                     |
|    | Keep away from sunlight                                                          |
|   | Keep dry                                                                         |
|  | Catalogue number                                                                 |
|  | Do not use if package is damaged                                                 |
|  | Do not resterilize                                                               |
|  | Upper limit of temperature range                                                 |
|  | Authorised representative in the European Community                              |
|  | MR conditional                                                                   |
|  | CE mark                                                                          |
|  | Up                                                                               |
|  | Fragile, handle with care                                                        |
|  | Maximum number of identical packages that can be stacked one on top of the other |

## I. Device Description

The Aegisy Vena Cava Filter with Delivery System is designed for percutaneous delivery of a permanent & temporary vena cava filter into the inferior vena cava (IVC).

The filter is laser cut from a nickel-titanium alloy tube (nitinol) and then molded into an asymmetrical lantern-like product. The proximal part resembles a "basket" in the form of six diamonds, and the distal part is shaped like the letter "Y". This design guarantees optimal retention of blood clots. The proximal and distal parts of the filter have one nut each, with which the product is connected to the delivery cable, which ensures controlled extraction. Six clamps, with which the product is attached to the vessel wall, are located at the junction of the distal part and six straight wires. An extraction hook is provided on the nut located at the distal end of the filter (See Figure 1).


The Aegisy Vena Cava Filter with Delivery System is delivered with an Introducer Kit. This Introducer Kit contains a Sheath (6F), a Dilator, a Delivery cable and a Loader. The Introducer Kit is needed to advance the filter through the Sheath to the proper position in the IVC. The ID of Dilator is 1.0 mm, allowing the use of a standard 0.035" (0.89 mm) guidewire. The guidewire is not included in the introducer kit.

Prior to implantation, the filter is connected to the delivery cable and then pulled back into the loading device. First, the dilator is passed over the guidewire through the femoral or jugular vein into the inferior vena cava, after which the guidewire is removed. Next, the filter from the loading device is moved along the dilator with a delivery cable. The dilator is pulled back, after which the filter is opened to the diameter of the inferior vena cava. The anatomical position of the filter is shown in Figure 25. (the proximal end of the filter is perpendicular to the renal vein). In this case, the filter exerts an external radial force on the inner surface of the vessel. Thus, a stable position of the product in the desired location is achieved (see Figure 25). The main design elements of the device are shown in Figure 2. After implantation, the cava filter can retain blood clots, while maintaining the patency of the inferior vena cava. Some dimensional characteristics of the product are presented in table 1.


Table 1

| Product Name and Catalog No.                           | Sheath Length(cm)     |                  | Angiographic Vessel Dilator     |            |
|--------------------------------------------------------|-----------------------|------------------|---------------------------------|------------|
|                                                        |                       |                  | ID(mm)                          | Length(cm) |
|                                                        | 55                    |                  | 1.0                             | 60         |
| Aegisy Vena Cava Filter with Delivery System XJLX326 0 | Recommended Guidewire |                  | Recommended Venous Access Sites |            |
|                                                        | Length (cm)           | Diameter (in)    |                                 |            |
|                                                        | 150                   | 0.035" (0.89 mm) | Femoral or Jugular              |            |

The guidewire used must be complied with ISO11070:2014+AMD1:2018 for sterile single-use intravascular catheter introducers. The guidewire is not included in the kit package.

**Figure 1**

**Figure 2**

## **II. Indications**

Aegisy Vena Cava Filter with Delivery System is indicated for filtrating thrombus after percutaneous placement in the inferior vena cava in the following situations:

- Pulmonary embolism, when anticoagulants are contraindicated
- Failure of anticoagulant therapy in thrombo-embolic diseases
- Emergency treatment following massive pulmonary embolism when anticipated benefits of conventional therapy are reduced
- Chronic, recurrent pulmonary embolism when anticoagulant therapy has failed or is contraindicated.

## **III. Contraindications**

- Patients with risk of septic embolism.
- Patients with uncontrolled infectious disease.
- Patients with an IVC diameter smaller than 20 mm and larger than 28 mm.
- Patients contraindicated for procedures under fluoroscopy.
- Patients with demonstrated hypersensitivity to one of the components of Aegisy
- Patients with deformed IVC.

#### **IV. Recommended Percutaneous Procedure for Aegisy Vena Cava Filter with Delivery System Implantation**

**A pre-placement cavogram is required:**

- To confirm the patency of the vein route from the puncture point to the proper placement position of filter in IVC.
- To evaluate the anatomy of the vena cava and identify any venous anomalies.
- To mark the level of the renal veins.
- To determine the highest level of any thrombus, which may be present.
- To determine the desired level for filter deployment and to mark the intended position with respect to the vertebral bodies.
- To ensure that the diameter of the vena cava (AP projection) at the planned filter deployment site is less than or equal to 28 mm in diameter and larger than or equal to 20 mm.

##### **1. Preparation of the placement accessories**

**Product preparation:** Take Aegisy venous cava filter with delivery system and wire (wire parameters are shown in Table 1)

##### **1.1 Remove all product components from sterile packaging**


**Note:** Contents are supplied sterile. Do not use if product sterile barrier or packaging is broken.

Whether through the femoral or jugular approach, the user must check the screw connection between the delivery cable and the filter before inserting the filter into the vein. Step described in part 1.2 for the femoral approach, part 1.3 for the jugular approach.

##### **1.2. For femoral insertion, check the screw connection between the delivery cable and the filter according to paragraphs 1.2.1-1.2.6 (Figure 3 – Figure 7)**


**Note:** Before checking, tightly connect the loading device and the hemostatic valve by rotating the first clockwise.

**1.2.1. Figure 3:** Pull the filter out of the loading device by pulling out the delivery cable. You can see that the distal part of the filter is connected to the delivery cable.


**Figure 3**

**1.2.2 Figure 4:** Disconnect the filter from the delivery cable by rotating the filter clockwise. Disconnection occurs on the distal nut

**Figure 4**


**1.2.3 Figure 5:** Connect the filter to the delivery cable by rotating the filter counterclockwise. The connection takes place on the distal nut

**Figure 5**


**1.2.4. Figure 6:** Pull the filter and the delivery cable back and forth to check whether the filter is securely attached to the delivery cable

**Note:** If the filter is disconnected from the delivery cable, you must stop using it.

**Figure 6**

**1.2.5. Figure 7:** Loosen the filter connection by a semicircle clockwise to ensure safe disconnection of the filter cable after final placement in the hollow vein.

**Figure 7**

**1.2.6.** After the described steps, insert the filter back into the boot device. Then rinse the product with sterile heparinized saline solution or other suitable isotonic solution.

**1.3.** During the jugular insertion, check the screw connection between the delivery cable and

the filter according to paragraphs 1.2.1, 1.2.2 and 1.3.1-1.3.8 (Figure 3, Figure 4, and Figure 8 – Figure 20).

**Note:** Before checking, tightly fasten the loading device and the hemostatic valve by rotating the loading device clockwise.

1.3.1 After 1.2.1 and 1.2.2, pull the delivery cable out of the loading device, then insert the delivery cable into the loading device, as shown in Figure 8, to the condition shown in Figure 9, so that the screw of the delivery cable can be seen from the side of the hemostatic valve.


Figure 8


Figure 9

1.3.2 Connect the filter to the delivery cable by rotating the filter counterclockwise as shown in Figure 10. The connection takes place on the distal nut




Figure 10

1.3.3 Pull the filter into the boot device, as shown in Figure 11, until the distal part of the filter exits the bootloader, as shown in Figure 12.

**Figure 11**


**Figure 12**

**1.3.4 Disconnect the delivery cable from the filter by rotating it counterclockwise, as shown in Figure 13 and Figure 14. Disconnection occurs on the distal nut**


**Figure 13**


**Figure 14**

**1.3.5 Insert the delivery cable to the filter on the other side of the boot device as shown in Figure 15. Connect the delivery cable tightly to the proximal nut at the proximal end of the filter by rotating clockwise as shown in Figure 16.**


**Figure 15**

**Figure 16**


**1.3.6 Pull the filter and the cable in different sides as shown in figure 17 to check if the filter is firmly attached to the Delivery Cable.**

**Note:** If the filter doff from the Delivery Cable, you should stop using.

**Figure 17**


**1.3.7 Loosen the filter connection by a half circle counter-clockwise rotation as seen from the cranial side of the filter as figure 18, to allow safe disconnection of the cable from the filter after final placement in the vena cava.**


**Figure 18**

**1.3.8 Pull out the delivery cable, as shown in Figure 19, so that the filter completely enters the loading device, as shown in Figure 20.**

**Note:** Make sure the filter is connected to the delivery cable. When the filter is implanted into the vena cava, it will be in the same condition, regardless of femoral or jugular insertion.

**Figure 19**

**Figure 20**

**1.3.9 After the above steps, insert the filter back into the boot device. Then rinse all components with sterile heparinized saline solution or other suitable isotonic solution.**

**1.4 Before and during implantation, a heparinized saline solution is administered to the patient to increase the clotting time by more than 200 seconds. Rinse the wire guide with heparinized saline solution or a suitable isotonic solution through a syringe. Remove the air from the dilator by rinsing it with heparinized saline or a suitable isotonic solution.**

**Note:** This drug reduces the possibility of blood clots forming around the filter during implantation.

## **2. Puncturing the vein**

**Select the appropriate venous access site and perform the puncture using the standard percutaneous puncture technique.**

**Note:** The exact location of the puncture is determined by the patient's anatomy, the position of the blood clot and the doctor's preferences. It is preferable to use femoral access when the venous pathway from the femoral vein to the renal vein is open. If this pathway is not open, it is recommended to use jugular access when the vein pathway for jugular access is open. If the venous pathways for femoral and jugular access are not open, cancel the implantation. The doctor can check whether the venous pathways are open for jugular access and femoral access or not by injecting a contrast agent. Insert the wire guide into the vessel through the puncture needle. Carefully move the wire guide to the desired location, as shown in Figure 21.

**Note:** Figures 21-26 show access through the right femoral vein only as an example, without affecting the procedure of other implantation routes.


Figure 21

### 3. Insertion of the Sheath

**3.1** While holding the wire guide in place, remove the piercing needle and press on the puncture site until the dilator is inserted into the vessel, as shown in Figure 22. Using a rotational motion, push the dilator through the tissue into the femoral or jugular vein and the inferior vena cava. Position the Sheath tip below the renal veins in preparation for an angiographic overview of the IVC (Cavogram).



Figure 22

**3.2** Position the tip of the dilator below the renal vein to prepare for an angiographic survey of the inferior vena cava (cavogram).

**3.3** Determine the diameter of the inferior vena cava at the intended implantation site below the inferior renal vein by injecting a contrast agent. At this time, place the end of the dilator in the right place.

**Note:** the dilator is radiopaque along the entire length

**3.4** While holding the dilator in the desired position, remove the wire guide as shown in **Figure 23**. Aspirate from the sideport extension to remove any potential air.

**Figure 23**


**3.5** Perform aspiration through the dilator opening to remove any possible air.

#### **4. Positioning and release of the filter**


**4.1** Insert the Loader with the filter and the Delivery Cable into the hub with the hemostasis valve until you feel resistance. Then slowly advance the filter into the distal of the Sheath.

**4.2.** Continue to push the filter inside the shell until the marker mark 1 of the delivery cable appears at the end of the Luer connector of the hemostatic valve. This indicates that the filter is located at the distal end of the dilator and completely inside the dilator. (as shown in **Figure 24 (a)**).

**Figure 24 (b):** Mark 1 indicates that the filter is at the distal tip of the Sheath but fully within the Sheath. Mark 2 indicates that the filter is completely deployed out of the Sheath.


**Figure 24 (a)**

**Figure 24 (b)**

**4.3. To deploy the filter, hold the delivery cable in the desired position, then pull the dilator back on the delivery cable until the marker mark 2 of the delivery cable appears at the end of the Luer connector of the hemostatic valve. This indicates that the filter is fully deployed from the dilator. (See Figure 25)**

**Figure 25**

**Note:**

a) if you encounter excessive resistance during the introduction of the filter, stop pressing and make sure that none of the components are damaged. If the dilator or delivery cable is bent, stop implantation, remove the damaged components from the vein and use a new kit to continue implantation.

b) Additional Note: During jugular insertion, you should very carefully check the position of the dilator and filter on the digital subtraction angiography (DSA) image. After confirming the correct position, you can deploy the filter from the dilator as indicated in paragraph 4.3. When the filter is fully deployed, it cannot be pulled back into the dilator, so as not to damage the inferior vena cava. You can adjust the position of the shell and filter according to the DSA image before the filter clips come out of the shell.

**4.2 Ensure that filter is fully deployed.**

**Note:** When femoral access, after deployment of the filter, verify angiographically that the filter is completely and correctly deployed. If the deployment of the filter is incomplete or the position of the filter is not correct or if the filter is tilting, it must be pulled back into the Sheath. After repositioning of the Sheath, the filter can be deployed again. This procedure has to be repeated until the filter is completely deployed in the desired position and shows no sign of tilting.

When jugular access, you should be very careful to check the position of the Sheath and the filter by DSA image. After you have confirmed the position is right, you just can deploy the filter from the Sheath. When the filter completely deploys, don't pulled back the filter into the Sheath to prevent from damage IVC. You can adjust the position of the Sheath and the filter by DSA image before the Secure Anchor is out of the Sheath.

**4.3** Once the Aegisy filter has been deployed, unscrew the Delivery Cable from the filter by counterclockwise rotation of the Delivery Cable and then remove the Delivery Cable from the Sheath.

**4.4** Remove the Sheath, when clinically indicated, by compressing the vessel above the puncture site and slowly withdrawing the Sheath.

**4.5** Now the filter is implanted.

**4.6** Discard the introducer kit and packaging materials.

**V. Recommended Devices and Procedure Required for the Withdrawal of the Aegisy Filter****Indications for removal**

- As a result of various clinical examinations and assessments, DVT (deep vein thrombosis) disappeared or was in a stable condition;
- The risk of PE (pulmonary embolism) is reduced to an acceptable level by continued appropriate basic treatment (treatment or prophylaxis) or changes in clinical conditions.
- Patients are not expected to return to high-risk PE status due to interruption of primary treatment, changes in clinical treatment plans, or changes in clinical conditions.
- Requires filter removal within an appropriate time period for removal.
- The cava filter becomes a potential source of serious disease due to displacement or loss of structural integrity, or may no longer provide embolism protection.

The Aegisy temporary filter can be withdrawn up to and including 14 days after implantation via the femoral vein only.

Recommended devices required to remove the product are presented in Table 2

**Table 2**

| Recommended devices | Dimension |
|---------------------|-----------|
| Angiographic Vessel | 10F       |


|                                       |                                                     |
|---------------------------------------|-----------------------------------------------------|
| <b>Withdrawal Catheter/Introducer</b> |                                                     |
| <b>Guidewire</b>                      | length 150cm, diameter 0.035" (0.89mm)              |
| <b>Snare Kit</b>                      | 15/20/25mm loop diameter, 6F catheter, length 102cm |

Angiographic catheter for vascular diversion and guidewire must comply with the ISO11070+AMD1:2018 sterile disposable sheaths for intravascular catheters, and be registered in the Russian Federation. Compatible loop set registered in the Russian Federation: ONE Snare Endovascular Loop Trapping System: ONE1500 No. PZN 2018/6864 dated 22.02.2018. The above items are not included in the scope of delivery.

#### **Recommended withdrawal procedure:**

##### **1. Placement of the snare**

- Select a puncture site that provides access to the distal end of the filter for removal.
- Puncture the femoral vein using standard percutaneous puncture technique; position the guidewire in the IVC.
- Place a 10F Angiographic Vessel Withdrawal Catheter/Introducer in IVC.
- Install the 6F loop catheter and advance the catheter to the distal nut on the distal end of the filter. (As shown in Figure 26)
- Remove the guidewire.
- Advance the snare through the catheter.


**Figure 26**


##### **2. Capture of the Withdrawal Hook**

- Capture the Withdrawal Hook with the snare loop. (As Figure 27)
- Gently apply pressure to the loop and then move the 6F catheter forward until it contacts the nut. The 6F catheter is now in contact with the distal end of the filter (As Figure 28, Figure 29)

**Warning:** while tightening the loop, press the 6F catheter against the distal end of the filter to prevent movement of the filter.

**Figure 27**

**Figure 28**


Figure 29

### 3. Withdrawal of the filter

- Keep tension on the loop and 6F catheter at the level of the nut. At the same time, advance the 10F angiographic catheter/sheath forward towards the filter until the filter is completely captured. (As Figure 30)
- Remove 10F Angiographic Vessel Withdrawal Catheter with the filter, 6F catheter and snare.
- Now the filter is successfully withdrawn.


Figure 30

## VI. Complications

Procedures requiring percutaneous interventional techniques should be attempted by physicians familiar with the possible complications. Complications may occur at any time

during the procedure.

Possible procedure complications include, but are not limited to the following:

- air embolism
- hematoma at the puncture site
- incorrect positioning of the filter
- incorrect orientation of the filter
- restriction of blood flow
- occlusion of small vessel
- infection
- perforation of the vessel wall
- filter fracture
- thrombus formation
- distal embolization

Possible long-term complications associated with filter implantation include, but are not limited to, the following:

- filter obstruction
- filter perforation of the vena cava wall
- filter migration
- filter fracture
- recurrent pulmonary embolism

Aegisy Vena Cava Filter with Delivery System has undergone fatigue testing simulating 10 years human implant conditions.

#### **VII. Storing and Management**

- Store in a cool, dark, dry place.
- Do not autoclave or resterilize by any other method. Do not expose to temperature above 50°C for any length of time.
- Do not expose to organic solvents.
- Handle with care during transport: keep it dry and cool and do not damage the package

#### **VIII. Packaging and Labeling**

- The Aegisy Vena Cava Filter with Delivery System is packaged in double sealed bag and a colored outer box. The package includes: pass certificate, instruction for use, quality feedback.
- Label information: Product name, sterilization date, use-before date, serial number, bar code, sterilization indication.

#### **IX. Precautions and Warning**

- All components in the Aegisy Filter and Introducer Kit are for single use only. Do not resterilize or reuse. Structural integrity and/or function may be impaired or be lost through cleaning, resterilization, or reuse and may cause adverse patient reactions. Lifetech will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the Aegisy Filter and Introducer Kit.
- As soon as the need for prevention of thromboembolism disappears, the product must

be removed. It is recommended to use a cava filter for 14 days after implantation inclusive. The procedure for removing the cava filter is described in Chapter V. "Recommended Devices and Procedures Required for the Withdrawal of the Aegisy Filter"

- When injecting contrast medium through the Dilator or the Sheath, do not exceed maximum pressure rating of 800 psi.
- After filter implantation any catheterization procedure requiring passage of a device through the filter must be avoided.
- Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are very rare. In addition, the FEA (Finite Element Analysis) and the fatigue test simulate the filter in IVC for 10 years, it can not fracture.
- The Aegisy Vena Cava Filter with Delivery System should only be used by physicians, who are trained in diagnostic and percutaneous interventional techniques. Lifetech will not be responsible for direct and consequential damages or expenses resulting from use by untrained personnel.
- Persons with allergic reactions to nickel titanium (Nitinol) may suffer an allergic response to this implant.
- A patient with an implanted Aegisy Vena Cava Filter with Delivery System can be scanned safely immediately after placement of the device under the following conditions:
  - ✧ Static magnetic field of 3 T or less.
  - ✧ Spatial gradient magnetic field of 720 G/cm or less.
  - ✧ Maximum MRsystem-reported, whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning.

**Note:** MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

- It's preferable to use femoral access when the vein route from the femoral vein to the renal vein is patent. If the vein route from the femoral to the renal vein is not patent, it is recommended to use jugular access when the vein route for jugular access is patent. If the vein routes for femoral access and jugular access are not patent, cancel the implantation. The physician can check whether the vein routes for jugular access and femoral access patent or not by injecting contrast medium.

#### **X. Transportation, storage and operation**

Products are transported by all modes of transport in covered vehicles in accordance with the rules of cargo transportation applicable to this type of transport in accordance with the conditions:

- **Transportation and storage conditions:** The Product should be transported and stored at an air temperature from 0°C to +50°C at a relative humidity of no more than 80% (non-condensing)
- **Operating conditions:** Temperature from +32 to +42 °C with relative humidity up to 100%

Store the product in the manufacturer's packaging in a dry, clean place, away from heat sources and direct sunlight, at a controlled temperature. The shelf life is indicated on the label of each product.

#### **XI. Information about the sterility of the product**

The product is delivered sterile. Sterilization parameters: - ethylene oxide sterilization method in accordance with the requirements of ISO 11135-1 (residual contamination level 10<sup>-6</sup>). Repeated sterilization of the product is prohibited. Repeated use is prohibited.

## **XII. Disposal information**

Dispose of the product in accordance with SanPiN 2.1.3684-21 for class B waste (epidemiologically hazardous waste). After hardware methods of disinfection using physical methods and changing the appearance of waste, which excludes the possibility of their reuse, class B waste can be accumulated, temporarily stored, transported, destroyed and disposed of together with class A waste (epidemiologically safe waste, similar in composition to MSW). The packaging of class B decontaminated medical waste must be labeled indicating that the waste has been decontaminated. In case of non-use of the product, damage to individual packaging and expiration of the storage period, disposal should be carried out in accordance with SanPiN 2.1.3684-21 for class A waste.

## **XIII. Compliance with the standards of the Russian Federation**


1. GOST R ISO 14630-2017 «Inactive surgical implants. General requirements».
2. GOST ISO 11607-1-2018 «Packaging for medical devices subject to terminal sterilization. Part 1. Requirements for materials, sterilization barrier systems and packaging systems».
3. GOST R ISO 15223-1-2020 "Medical devices. Symbols used for labeling on medical devices, on labels and in accompanying documentation. Part 1. Basic requirements".
4. GOST ISO 10555-1-2021 "Sterile intravascular catheters for single use. Part 1. General technical requirements.».
5. GOST R 52770-2016 «Medical devices. Safety requirements. Methods of sanitary-chemical and toxicological tests».
6. GOST 31209-2003 «Containers for blood and its components. Chemical and biological safety requirements and test methods».
7. GOST 4011-72 «Drinking water. Methods for measuring the mass concentration of total iron» p.2.
8. GOST 31870-2012 «Drinking water. Determination of the content of elements by atomic spectrometry» p.4.
9. MUK 4.1.3166-14 «Gas chromatographic determination of hexane, heptane, acetaldehyde, acetone, methyl acetate, ethyl acetate, methanol, isopropanol, acrylonitrile, n-propanol, n-propyl acetate, butyl acetate, isobutanol, n-butanol, benzene, toluene, m-, o- and p-xylenes, isopropylbenzene, styrene, α-methylstyrene in water and aqueous extracts from metals of various compositions».
10. GOST R 55227-2012 «Water. Methods for determining the content of formaldehyde» p.6.
11. MR 2413-81 «Guidelines for the determination of epichlorohydrin in aqueous extracts from polymer materials».
12. MR 1436-76 «Guidelines for the determination of diphenylpropane, as well as some phenols in its presence, during sanitary and chemical studies of products made of polymer materials intended for contact with food».
13. I-880-71 «Instructions for the sanitary and chemical examination of products made of polymer and other synthetic materials intended for contact with food».
14. MVI MN 1924-2003 «Method of gas chromatographic determination of phenol and epichlorohydrin in model media simulating food products».
15. GOST 31214-2016 «Medical devices. Requirements for samples and documentation provided for toxicological, sanitary-chemical tests, tests for sterility and pyrogenicity».
16. GOST ISO 10993-1-2011. «Medical products. Assessment of the biological effect of medical devices. Part 1. Assessment and Research».
17. GOST ISO 10993-2-2009. «Medical products. Assessment of the biological effect of medical devices. Part 2. Requirements for handling animals».


18. GOST ISO 10993-4-2011. «Medical products. Assessment of the biological effect of medical devices. Part 4. Investigation of products that interact with blood».
19. GOST ISO 10993-5-2011. «Medical products. Assessment of the biological effect of medical devices. Part 5. Testing for cytotoxicity: in vitro methods».
20. GOST ISO 10993-6-2011 «Medical devices. Assessment of the biological effect of medical devices. Part 6. Studies of local action after implantation».
21. GOST ISO 10993-7-2016 «Medical devices. Assessment of the biological effect of medical devices. Part 7. Residual content of ethylene oxide after sterilization».
22. GOST ISO 10993-10-2011 «Medical devices. Assessment of the biological effect of medical devices. Part 10. Investigations of irritant and sensitizing effects».
23. GOST ISO 10993-11-2011 «Medical devices. Assessment of the biological effect of medical devices. Part 11. Studies of general toxic action».
24. GOST ISO 10993-12-2015. «Medical products. Assessment of the biological effect of medical devices. Part 12. Sample preparation and control samples».
25. GF RF OFS.1.2.4.0003.15. Sterility.
- 26) OFS.1.2.4.0005.15 Pyrogenicity.

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