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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Shanghai MicroPort
Medical (Group) Co., Ltd. on the annexed INSTRUCTON FOR
USE is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:


Zhou Pingfan

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(Date: Sep. 07, 2021)

证 询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

I affirm
Dr.Zhaohua Chang
Chairman&CEO
Shanghai MicroPort Medical (Group) Co., Ltd

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INSTRUCTIONS FOR USE

Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System

Revision: A.

INSTRUCTIONS FOR USE

Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System

PRODUCT DESCRIPTION

The Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System (Firehawk Liberty™ RTECSS or Firehawk Liberty™) is a device/drug combination product consisting of:

- Drug/Polymer Coating:** a formulation of rapamycin in a biodegradable polymer coating which filled in grooves that located on outside surface of stent struts (Fig. 1.1 and Fig. 1.2).
- A Stent and Delivery System:** L605 CoCr stent pre-mounted to a RX delivery catheter.


Fig. 1.1 The Schematic diagram of Firehawk Liberty™ RTECSS


Fig. 1.2 Construction of Groove loaded drug product

The device characteristics are summarized in Table 1.1 and 1.2.

Table 1.1 Firehawk Liberty™ RTECSS Product Description

Parameter	Characteristic
Stent configuration	shape - cylindrical
Delivery system	Rapid-Exchange (RX) type balloon catheter
Available Stent Lengths (mm)	8, 13, 16, 18, 21, 23, 26, 29, 31, 33, 35, 38
Available Stent Diameters (mm)	2.25; 2.5; 2.75; 3.0; 3.5; 4.0
Stent Material	L-605 cobalt chromium (CoCr) alloy
Drug Coating	An abluminal coating of poly(lactic acid) loaded with Rapamycin, Fig. 1.2
Balloon material	Pebax 3533 & Nylon 12
Radiopaque strip	Platinum and iridium alloy
Hydrophilic coating	PVP
Distances to the Guide wire	250 mm $\pm$ 50 mm

Parameter		Characteristic																
Distal tip of the catheter	outlet																	
	Position of the first mark	950 mm ± 50 mm																
	Position of the second mark	1050 mm ± 50 mm																
Effective catheter length		1400 mm ± 100 mm																
Delivery system design		Single access port to inflation lumen. Guide wire exit notch is located 25 cm from tip. Designed for guide wires ≤ 0.014" (0,36 mm).																
Stent Delivery System Balloon		A semi-compliant balloon with two radiopaque markers to designate the stent																
Balloon inflation pressure		Nominal inflation pressure: 10 atm (1013 kPa) for all sizes; Rated Burst Pressure (RBP): 16 atm (1621 kPa) for 2.25 to 3.5mm Diameters; 14 atm (1419 kPa) for 4.0mm Diameters																
Guiding Catheter Inner Diameter		≥5F (0.056")																
Catheter Shaft Outer Diameter (nominal)		Φ2.25-3.0 mm		Φ3.5-4.0 mm														
	Distal:	2.6 F (0.034"[0.864mm])		2.8F (0.037"[0.940mm])														
	Proximal:	2.3F (0.030"[0.762mm])		2.6 F (0.034"[0.864mm])														
Отношение длины стента к диаметру		При раскрытии стента до минимального и максимального диаметров изменение длины стента должно составлять ≤20%.																
The ratio of stent length to diameter		After the tracking process, the delivery system has no visible damage (for example, kinks), the stents are not loosened or misaligned; The distance between the outer diameter of the balloon and the outer diameter of the twisted stent profile should be ≤ 5% both before and after tracking, as well as with a bend radius of 8 mm.																
Vessel wall compatibility		The outer surface of the stent should fit snugly against the vessel wall, and there should be no deviation or separation of the stent from the vessel wall.																
Visualization		Stents should be visualized by digital subtraction angiography (CSA) or fluoroscopy.																
Resistance to fracture under radially applied load		The table below shows the fracture resistance for stents of different diameters at 12% deformation.																
		<table><tr><th>Nominal diameter (mm)</th><th>12% deformation and fracture resistance (kPa)</th></tr><tr><td>φ2,25</td><td>≥110</td></tr><tr><td>φ2,5</td><td>≥105</td></tr><tr><td>φ2,75</td><td>≥100</td></tr><tr><td>φ3,0</td><td>≥95</td></tr><tr><td>φ3,5</td><td>≥90</td></tr><tr><td>φ4,0</td><td>≥85</td></tr></table>			Nominal diameter (mm)	12% deformation and fracture resistance (kPa)	φ2,25	≥110	φ2,5	≥105	φ2,75	≥100	φ3,0	≥95	φ3,5	≥90	φ4,0	≥85
		Nominal diameter (mm)	12% deformation and fracture resistance (kPa)															
		φ2,25	≥110															
		φ2,5	≥105															
		φ2,75	≥100															
		φ3,0	≥95															
φ3,5	≥90																	
φ4,0	≥85																	
Decreased vessel lumen (Reduction)		≤ 10%																
Durability		The stent must maintain its integrity without any obvious defects or damage for 10 years.																
Bend / Kink		The minimum radius of curvature of a deployed stent should be ≤ 8.0 mm.																
Stent free surface area and stent outer surface area		For stents with a specified diameter, minimum diameter and maximum diameter, the surface area free from the stent should be more than 70% of the total surface area of cylinders of the same diameter and length.																
The appearance of the stent system		The surface of the delivery system must be smooth, free of visible foreign matter; without kinks (especially around the conductor recess). No damage to the tip.																
Stent profile		Stent diameter at nominal pressure (mm)		Profile														
		φ2,25, φ2,5		≤ 0,042 "/1,067 mm														
		φ2,75		≤ 0,043 "/1,092 mm														
		φ3,0		≤ 0,044 "/1,118 mm														
		φ3,5, φ4,0		≤ 0,051 "/1,295 mm														
Stent displacement force		The peak displacement force is at least 2.0 N at an initial distance of 0-5 mm (displacement force is the force required to displace the stent in relation to the balloon without the possibility of recovery).																

Parameter	Characteristic	
Leak	When filling the system with contrast agent, there should be no leaks or damage, such as tears or protrusions on the delivery system. There should be no leaks in the Luer lock, in any area of the connection, or in any other element. Air should not enter through the assembly points of the delivery system during the suction process.	
Bend / Kink	The stent / delivery system should not be kinked during delivery and retrieval.	
Visualization	The radiopaque marker stripes of the balloon catheter should be visualized during access using fluoroscopy.	
Balloon inflation time	Balloon inflation time ≤ 20 seconds	
Balloon deflation time	Balloon deflation time ≤ 20 seconds	
Estimated cylinder wear	After 10 filling and emptying cycles, the delivery system should be free of leaks or ruptures.	
Balloon expansion at both ends of the stent	The diameter of the portion of the balloon outside the stent should not exceed 5% of the outer diameter of the stent at the design burst pressure.	
Hemostasis	At a pressure of 25 kPa and a pushing force of 2 ± 0.2 N, the volume of water flowing out of the hemostatic valve in 20 minutes should not exceed 10 ml.	
Isolation of the drug	Milestone	Selection intensity
	10 min	$\leq 35\%$
	1 hour	22%~58%
	6 hours	41%~82%
	48 hours	$\geq 66\%$
Component size compatibility	Guide wire lumen diameter and travel smoothness: The wire > 0.014 " (0.36mm) can move smoothly without stopping inside the guide wire lumen, and the maximum pushing force does not exceed 0.3 N.	
Factors of stent functioning	- The probability of occurrence of contact, galvanic, pitting and fretting corrosion, the possibility of oxidation: the probability is low, provided by the use of bioinert materials - Tissue perforation hazard: Tissue perforation or rupture can occur if a stent is used that does not match the vessel diameter, or if the precautions regarding stent use are not followed. - Information on biodegradation: low probability, provided by the use of bioinert materials.	

Table 1.2 Firehawk Liberty™ RTECSS Product Matrix and Rapamycin Content

Product Code	Stent Diameter (mm)	Stent Length (mm)	Drug Content (µg)
FR2208	2.25	8	18
FR2508	2.5	8	18
FR2708	2.75	8	21
FR3008	3.0	8	21
FR3508	3.5	8	27
FR4008	4.0	8	27
FR2213	2.25	13	30
FR2513	2.5	13	30
FR2713	2.75	13	35
FR3013	3.0	13	35
FR3513	3.5	13	45
FR4013	4.0	13	45
FR2216	2.25	16	38
FR2516	2.5	16	38
FR2716	2.75	16	41
FR3016	3.0	16	41
FR3516	3.5	16	53
FR4016	4.0	16	53
FR2218	2.25	18	43
FR2518	2.5	18	43
FR2718	2.75	18	47
FR3018	3.0	18	47

Product Code	Stent Diameter (mm)	Stent Length (mm)	Drug Content (µg)
FR3518	3.5	18	60
FR4018	4.0	18	60
FR2221	2.25	21	48
FR2521	2.5	21	48
FR2721	2.75	21	53
FR3021	3.0	21	53
FR3521	3.5	21	68
FR4021	4.0	21	68
FR2223	2.25	23	53
FR2523	2.5	23	53
FR2723	2.75	23	59
FR3023	3.0	23	59
FR3523	3.5	23	76
FR4023	4.0	23	76
FR2226	2.25	26	60
FR2526	2.5	26	60
FR2726	2.75	26	68
FR3026	3.0	26	68
FR3526	3.5	26	87
FR4026	4.0	26	87
FR2229	2.25	29	68
FR2529	2.5	29	68
FR2729	2.75	29	76
FR3029	3.0	29	76
FR3529	3.5	29	98
FR4029	4.0	29	98
FR2531	2.5	31	73
FR2731	2.75	31	79
FR3031	3.0	31	79
FR3531	3.5	31	102
FR4031	4.0	31	102
FR2533	2.5	33	78
FR2733	2.75	33	85
FR3033	3.0	33	85
FR3533	3.5	33	110
FR4033	4.0	33	110
FR2735	2.75	35	91
FR3035	3.0	35	91
ER3535	3.5	35	117
FR4035	4.0	35	117
FR2738	2.75	38	97
FR3038	3.0	38	97
FR3538	3.5	38	125
FR4038	4.0	38	125

Product Code Format: FR xxyy:

xx - stent diameter, mm

yy - stent length, mm.

Type of contact with the body

A stent is a lifelong contact with the internal environment of the body (with blood or lymph).

Delivery system - short-term contact with blood and lymph.

2 PURPOSE AND INDICATIONS FOR USE

2.1. PURPOSE

The Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System (Firehawk Liberty™

RTECSS) is indicated for improving coronary luminal diameter in patients with symptomatic heart disease due to de novo native coronary artery lesions length ≤ 60 mm with reference vessel diameters of ≥ 2.25 to ≤ 4.0 mm.

2. INDICATIONS

The Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System (Firehawk Liberty™ RTECSS) is indicated for improving coronary luminal diameter in patients with symptomatic heart disease due to de novo native coronary artery lesions length ≤ 60 mm with reference vessel diameters of ≥ 2.25 to ≤ 4.0 mm.

3 CONTRAINDICATIONS

The Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System (Firehawk Liberty™ RTECSS) is contraindicated for use in patients:

- Who cannot receive antiplatelet and/or anti-coagulant therapy;
- With lesions that prevent complete angioplasty balloon inflation or proper placement of the stent or stent delivery system;
- With hypersensitivity or contraindication to rapamycin or similar drugs or any other analogue or derivative, cobalt, chromium, nickel, tungsten, or PLA.

4 WARNINGS

- Ensure that the inner package sterile barrier has not been opened or damaged prior to use.
- Judicious patient selection is necessary since the use of the device carries the associated risks of thrombosis, vascular complications, and/or bleeding events.
- This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy.

5 PRECAUTIONS

5.1 General Precautions

- Stent implantation should only be performed by physicians who have received appropriate training.
- Stent placement should be performed at hospitals where emergency coronary artery bypass graft surgery is accessible.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. Long-term outcomes following repeat dilatation of the stent are presently unknown.
- Risks and benefits should be considered in patients with severe contrast agent allergies.
- Do not expose or wipe the product with organic solvents such as alcohol or detergents.
- Care should be taken to control the guiding catheter tip during stent delivery, deployment, and delivery system withdrawal. Before withdrawing the stent delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter advancement into the vessel and subsequent arterial damage.
- Compared to use within the specified Indications for Use, the use of DES in patients and lesions outside of the labeled indications may have an increased risk of adverse events, including stent thrombosis, stent embolization, MI, or death.
- It is very important that the patient comply with post-procedural antiplatelet therapy recommendations. It should be adapted depending on patient profile, clinical indication and intervention characteristics, and in accordance with ESC and ACC/AHA/SCAI guidelines.
- Early discontinuation of prescribed antiplatelet medications could result in a higher risk of stent thrombosis, MI, or death. The risks and benefits of the procedure should be weighed against the possible risks associated with early discontinuation of antiplatelet therapy. Patients who require early discontinuation of antiplatelet therapy should be monitored carefully for cardiac events.
- Orally administered sirolimus combined with cyclosporine is associated with increased serum cholesterol and triglycerides levels.
- When used with cyclosporine medication, The sirolimus mean AUC and mean Cmax may be affected.

Use of Multiple Stents

When multiple stents are required, stent materials should be of similar composition. Placing multiple stents of different materials in contact with each other may increase potential for corrosion. Potential interactions of the Firehawk Liberty™ stent with other drug-eluting or coated stents have not been evaluated.

A patient's exposure to drug and polymer is proportional to the number and total length of implanted stents. Use of more than two Firehawk Liberty™ stents has not received adequate long term clinical evaluation. The maximum amount of stents used needs to be determined by the physician based on the risks versus the clinical benefit for the patient involved.

5.3 Brachytherapy

The safety and effectiveness of the Firehawk Liberty™ stent in patients with prior brachytherapy of the target lesion have not been established. The safety and effectiveness of the use of brachytherapy to treat in-stent restenosis in a Firehawk Liberty™ stent has not been established. Both vascular brachytherapy and the Firehawk Liberty™ stent alter arterial remodeling. The synergy between these two treatments has not been determined.

5.4 Drug Interactions

The effect of drug interactions on the safety or efficacy of the Firehawk Liberty™ stent has not been determined. Several drugs are known to affect the metabolism of Rapamycin, and other drug interactions may be inferred from known metabolic effects. Rapamycin is known to be a substrate for both cytochrome P450 IIIA4 (CYP3A4) and P-glycoprotein.

Consideration should be given to the potential for drug interaction when deciding to implant Firehawk Liberty™ stent if patient is taking a drug that might interact with Rapamycin, or when deciding to initiate therapy if patient has recently received a Firehawk Liberty™ stent.

5.5 Magnetic Resonance Imaging (MRI)

Non-clinical testing has demonstrated the Firehawk Liberty™ stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla
- spatial gradient field of 36 T/m
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg or less at 3 Tesla, Firehawk Liberty™ stent, for 15 minutes of continuous MR scanning.

In non-clinical testing, the single and overlapping (72.4 mm) Firehawk Liberty™ stent respectively produced a temperature rise of 1.2 °C and 3.2 °C at a maximum whole body averaged specific absorption rate (SAR) of ≈ 3.53 W/kg, as assessed by calorimetry for 15 minutes of continuous MR scanning with whole body coil in a 3 Tesla Magnetom Trio, Siemens Medical Solutions (Software: Numaris/4, syngo MR A30) MR Scanner.

Non-clinical testing has not been performed to rule out the possibility of implant migration and RF heating at static gradient magnetic fields stronger than mentioned above.

MR imaging quality may be compromised if the area of interest is in exactly the same area or relatively close to the position of the stent.

5.6 Stent Handling Precautions

- Each stent is for single use only. Do not resterilize or reuse this device. Note the "Use-by" (expiration) date on the product label.
- The foil pouch is not a sterile barrier. The inner header bag (pouch) within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.
- The foil pouch contains 1 packet of oxygen absorber. Please discard this packet after opening the pouch.
- Do not remove the stent from the delivery system as removal may damage the stent and/or lead to stent embolization. These components are intended to perform together as a system.
- Do not induce a vacuum on the delivery system prior to reaching the target lesion. This may cause dislodgement of the stent from the balloon.
- The delivery system should not be used in conjunction with other stents.
- Special care must be taken not to handle or disrupt the stent on the balloon especially during

delivery system removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.

Do not manipulate, touch, or handle the stent with your fingers, which may cause contamination or stent dislodgement from the delivery balloon.

Use only the appropriate balloon inflation media (see Section 9.3.4 – Operator's Instructions, Delivery System Preparation). Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in stent deployment.

5.7 Stent Placement Precautions

5.7.1 Stent Preparation

- Do not pre-inflate the delivery system prior to stent deployment other than as directed. Use the balloon purging technique described in Section 9.3.4 – Operator's Instructions, Delivery System Preparation.
- Use guiding catheters which have lumen sizes that are suitable to accommodate the stent delivery system (see Section 1 – Product Description).

5.7.2 Stent Implantation

- The vessel should be pre-dilated with an appropriate sized balloon. Failure to do so may increase the difficulty of stent placement and cause procedural complications.
- Do not expand the stent if it is not properly positioned in the vessel (see Section 5.8 – Precautions, Stent/System Removal Precautions).
- Implanting a stent may lead to vessel dissection and acute closure require additional intervention (CABG, further dilatation, placement of additional stents, or other).
- Although the safety and effectiveness of treating same coronary artery with more than one Firehawk Liberty™ stents has not been established, if this is performed, place the stent in the distal lesion before the proximal lesion in order to minimize dislodgement risk incurred by traversing through deployed stents.
- Stent placement may compromise side branch patency.
- Do not exceed Rated Burst Pressure (RBP) as indicated on product label. See Table 9.8, Typical Firehawk Liberty™ RTECSS Compliance. Balloon pressures should be monitored during inflation. Applying pressures higher than specified on the product label may result in a ruptured balloon with possible arterial damage and dissection.
- Do not attempt to pull an unexpanded stent back through the guiding catheter, as dislodgment of the stent from the balloon may occur. Remove as a single unit as indicated in Section 5.8 – Precautions, Stent/System Removal Precautions.
- Stent retrieval methods (i.e., using additional wires, snares, and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- Although the stent delivery system balloon is strong enough to expand the stent without rupture, a circumferential balloon tear distal to the stent and prior to complete stent expansion, could cause the balloon to become tethered to the stent, requiring surgical removal. In case of balloon rupture, it should be withdrawn and, if necessary, a new dilatation catheter exchanged over the guide wire to complete the expansion of the stent.
- Ensure the stented area covers the entire lesion/dissection site and that no gaps exist between stents.

5.8 Stent/System Removal Precautions

If any unusual resistance is felt at any time during lesion access before stent implantation, the stent system and the guiding catheter should be removed as a single unit.

When removing the stent system and guiding catheter as a single unit:

- DO NOT retract the stent system into the guiding catheter.
- Position the proximal balloon marker just distal to the guiding catheter tip.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the delivery system to the guiding catheter, and remove the guiding catheter and stent system as a single unit.

Failure to follow these steps and/or applying excessive force to the delivery system can potentially result in loss or damage to the stent and/or delivery system components and/or vasculature.

If it is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

Post-Procedure Precautions

Care should be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guidewire, or balloon catheter to avoid disrupting the stent placement, apposition, geometry.

Post-dilatation: All efforts should be made to assure that the stent is not under dilated. If the deployed stent is not fully apposed to the vessel wall, the stent may be expanded further with a larger diameter balloon. The post-dilatation can be done using a low-profile, high pressure, non-compliant balloon catheter. The balloon should not extend outside of the stented region (see Section 9.7 – Operator's Instructions, Post-Deployment Dilatation of Stent Segments). Do not use the stent delivery balloon for post-dilatation.

- Non-clinical testing has demonstrated that the Firehawk Liberty™ stent is MR Conditional (see Section 5.5 – Precautions, Magnetic Resonance Imaging (MRI)). MR imaging quality may be compromised if the area of interest is in the same area of the position of the stent.
- Antiplatelet therapy should be administered post-procedure (see Section 5.1 – Precautions, General Precautions). Patients who require early discontinuation of antiplatelet therapy (e.g., secondary to active bleeding) should be monitored carefully for cardiac events. At the discretion of the patient's treating physician, the antiplatelet therapy should be restarted as soon as possible.

6 ADVERSE EVENTS

Adverse events may be associated with the implantation of a coronary stent in coronary arteries, but are not limited to the following:

- Allergic reaction to anti-coagulant and/or antiplatelet therapy, contrast medium, or stent materials
- Arrhythmias
- Arteriovenous fistula
- Bleeding
- Cardiac tamponade
- Coronary aneurysm
- Death
- Dissection
- Drug interactions with antiplatelet/anticoagulant/contrast medium
- Emboli, distal (tissue, air, or thrombotic emboli)
- Embolization, stent
- Emergency CABG
- Failure to deliver the stent to the intended site
- Fever
- Heart failure
- hemorrhage
- Hypotension/Hypertension
- Infection, local or systemic
- Myocardial infarction
- Pain, at the access site
- Pseudoaneurysm, femoral
- Restenosis of stented segment
- Stent embolization or migration
- Stent fracture
- Stent thrombosis/occlusion
- Target lesion revascularization
- Target vessels of non-target lesion revascularization
- Total occlusion of coronary artery
- Vessel trauma requiring surgical repair or reintervention

Potential adverse events not captured above, that may be unique to the rapamycin drug coating

- Abnormal liver function tests
- Allergic/immunologic reaction to drug (rapamycin or structurally-related compounds) or the polymer stent coating or its individual components (see Section 1 Product Description)
- Anemia
- Arthralgias
- Diarrhea
- Hypercholesterolemia

- Hypersensitivity, including anaphylactic/anaphylactoid type reactions
- Hypertriglyceridemia
- Hypokalemia
- Infections
- Interstitial lung disease
- Myocardial infarction
- Myocardial ischemia
- Occlusion
- Prolonged angina
- Pseudoaneurysm
- Hematologic dyscrasia (including leukopenia, neutropenia, thrombocytopenia)
- Renal failure
- Restenosis of stented segment (greater than 50% obstruction)
- Stroke
- Vessel spasm
- Vessel perforation

7 INDIVIDUALIZATION OF TREATMENT

The risks and benefits described above should be considered carefully for each patient before use of the Firehawk Liberty™ stent. Patient selection factors to be assessed should include a judgment regarding risk of anti-platelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.

Premorbid conditions that increase the risk of a poor initial result and the risks of emergency referral for bypass surgery (diabetes mellitus, renal failure, and severe obesity) should be reviewed. The relation of baseline and procedural variables to Major Adverse Cardiac Events (MACE) was examined. Multivariable modeling suggested that treatment assignment remained an independent predictor of clinical and angiographic outcomes even after adjusting for other baseline and procedural confounding variables.

8 USE IN SPECIAL POPULATIONS

8.1 Pregnancy

There are no adequate data from the use of sirolimus or Firehawk Liberty™ stents in pregnant women. In animals, sirolimus was embryo- and fetus- toxic at clinically relevant exposures. Teratogenic effects have not been observed. The potential risk for humans is unknown. The Firehawk Liberty™ stent should not be used in pregnant women, unless the clinical condition of the mother requires treatment with the stent.

8.2 Lactation

It is not known whether sirolimus is excreted in human breast milk. Sirolimus is excreted in milk of lactating rats. A risk to the suckling child cannot be excluded. A decision must be made whether to discontinue breast-feeding or to implant the stent taking into account the importance of the stent to the mother.

8.3 Pediatric Use

The safety and effectiveness of the Firehawk Liberty™ stent in pediatric patients below the age of 18 years have not been established.

8.4 Geriatric Use

Clinical studies of the Firehawk Liberty™ stent did not suggest that patients age 65 years and over differed with regard to safety and efficacy compared to younger patients.

Carefully consider whether it is appropriate to use in the above patient populations.

9 OPERATOR'S INSTRUCTIONS

9.1 Inspection Prior to Use

Carefully inspect the sterile package before opening. Do not use if the package has been damaged or opened. Do not use the product after the "Use-by" date.

9.2 Materials Required

Quantity	Material *
	Guiding catheter $\geq 5F$ (1.42 mm, 0.056 inch) inner diameter]
3	20 cc syringe
1000 u/500 cc	Sterile Heparinized Normal Saline (HepNS)
1	Guidewire ≤ 0.014 inch (0.36 mm) diameter]
1	Rotating hemostatic valve with appropriate minimum inner diameter (0.096 inch [2.44 mm])
N/A	Contrast medium diluted 1:1 with heparinized normal saline
1	Inflation device
1	Stopcock (3-way minimum)
1	Torque device
1	Pre-deployment dilatation catheter
1	Appropriate arterial sheath
N/A	Guiding catheter $\geq 5F$ (1.42 mm, 0.056 inch) inner diameter]

* The medical device must have a valid marketing authorization.

9.3 Preparation

9.3.1 Packaging Removal

Note: The foil pouch is not a sterile barrier. The inner pouch within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.

- Peel-open the inner pouch and pass or drop the product into the sterile field using an aseptic technique.
- Carefully remove the stent system from coiled hoop and inspect it for bends, kinks, and other damage. Do not use the product if any damage is noted.
- Carefully remove the protective stent sheath and mandrel together, visually inspect the stent to ensure that it has not been damaged or displaced from its original position (between proximal and distal marker bands) on the balloon.

Note: Do not touch the stent or move the stent with your fingers to avoid contamination of the stent or dislodging it from the balloon of the delivery system.

9.3.2 Rinse Catheter

- Rinse the catheter with sterile heparinized normal saline solution.

9.3.3 Guide Wire Lumen Flush

- Attach the syringe with HepNS to the flushing needle packaged with the catheter.
- Insert the needle into the tip of the catheter and flush the guide-wire lumen with HepNS.

Note: AVOID manipulation of the stent during flushing of the guide-wire lumen, as this may disrupt the placement of the stent on the balloon.

9.3.4 Delivery System Preparation

- 1) Prepare the inflation device or syringe with suitable diluted contrast medium.
- 2) Attach the inflation device or syringe to the stopcock; attach to the balloon inflation port hub.
- 3) Open the stopcock to stent delivery system.
- 4) Leave the inflation device or syringe on neutral.

Note: While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.

9.4 Delivery Procedure

- 1) Prepare the vascular access site according to standard practice.
- 2) Pre-dilate the lesion with a Percutaneous Transluminal Coronary Angioplasty (PTCA) catheter. Limit the longitudinal length of pre-dilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of Firehawk Liberty™ stent.
- 3) Maintain neutral pressure on the inflation device. Open the rotating hemostatic valve as widely as possible.
- 4) Backload the delivery system onto the proximal portion of the guide-wire while maintaining the guide-wire position across the target lesion.
- 5) Advance the stent delivery system over the guide-wire to the target lesion. Use the radiopaque balloon markers to position the stent across the lesion; perform angiography to confirm the position of the stent.
- 6) Sufficiently tighten the rotating hemostatic valve. The stent is now ready to be deployed.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, do not force passage. Remove the entire system as a

single unit. See Section 5.8 – Precautions, Stent/System Removal Precautions for specific delivery system removal instructions.

Deployment Procedure

Before deployment, reconfirm the correct position of the stent relative to the target lesion via the radiopaque balloon markers.

- 1) Attach the inflation device (only partially filled with contrast media) to three-way stopcock and apply negative pressure to purge the balloon of air.
- 2) Turn the stopcock on the catheter to the off position and purge the inflation device of air. Close the side port of the stopcock.
- 3) Under fluoroscopic visualization, inflate the balloon to at least the nominal pressure to deploy the stent. Do not exceed the labeled RBP. Deploy the stent slowly by pressurizing the delivery system until stent is completely expanded. Maintain pressure for 30 seconds. If necessary, the delivery system can be repressurized or further pressurized to assure complete apposition of the stent to the artery wall. Do not exceed the labeled rated burst pressure (RBP) of 16 atm (1621 kPa) for the 2.25 mm – 3.5 mm stents and 14 atm (1419 kPa) for the 4.0 mm stents sizes.
- 4) Stent wall contact should be verified through routine angiography or intravascular ultrasound.
- 5) Fully cover the entire lesion and balloon treated area (including dissections) with the Firehawk Liberty™ stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
- 6) If more than one Firehawk Liberty™ stent is needed to cover the lesion and balloon treated area, adequately overlap stents, taking into account stent foreshortening. Ensure no gaps between stents by positioning the balloon marker bands of the second Firehawk Liberty™ stent inside the deployed stent prior to expansion. See Section 5.8 – Precautions, Stent/System Removal Precautions.
- 7) Deflate the balloon by pulling negative pressure on the inflation device. Allow adequate time, at least 15 seconds, for full balloon deflation. Confirm complete balloon deflation before attempting to remove the delivery system. Ensure the balloon is fully deflated before attempting to move the catheter.
- 8) Confirm that the stent is adequately expanded and in full contact with the artery wall.

9.6 Removal Procedure

- 1) Ensure that the balloon is fully deflated.
- 2) Fully open the rotating hemostatic valve.
- 3) While maintaining the guide wire position and negative pressure on the inflation device, withdraw the delivery system.

Note: Should any resistance be felt at any time during either lesion access or removal of the delivery system post stent implantation, the entire system should be removed as a single unit. See Section 5.8 – Precautions, Stent/System Removal Precautions for specific delivery system removal instructions.

- 4) Tighten the rotating hemostatic valve.
- 5) Repeat angiography to assess the stented area. If post-dilatation is necessary, ensure that the final stent diameter matches the reference vessel diameter. Assure that the stent is not underdilated.

9.7 Post-Deployment Dilatation of Stent Segments

All efforts should be taken to assure that the stent is not underdilated. If the deployed stent size is still inadequate with respect to vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. The stent may be further expanded using a low profile, high pressure, and non-compliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guidewire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region.

CAUTION: Do not dilate the stent beyond the following limits.

Nominal Stent Diameter	Dilatation Limit
2.25 mm to 2.5 mm	2.75 mm
2.75 mm to 3.0 mm	3.5 mm
3.5 mm to 4.0 mm	4.5 mm

Reconfirm stent position and angiographic results. Repeat inflations until optimal stent deployment is achieved.

9.8 In Vitro Compliance Information

Table - Typical Firehawk Liberty™ stent Compliance

Pressure			Stent					
			Inner Diameter (mm)					
atm	kPa		2.25	2.5	2.75	3.0	3.5	4.0**
10	1013	Nominal	220	2.45	2.70	2.95	3.50	3.90
11	1115		225	2.49	2.73	3.00	3.54	3.98
12	1216		230	2.52	2.78	3.05	3.61	4.06
13	1317		233	2.57	2.82	3.09	3.66	4.12
14	1419		236	2.61	2.85	3.12	3.70	4.18
15	1520		239	2.65	2.88	3.16	3.74	4.22
16	1621	RBP*	2.42	2.68	2.91	3.19	3.78	4.27
17	1723		2.44	2.71	2.94	3.22	3.82	4.32
18	1824		2.47	2.73	2.97	3.25	3.86	4.37

*Do not exceed the rated burst pressure (RBP)

**Rated Burst Pressure (RBP): 14 atm (1419 kPa) for 4.0 mm diameters

Note:

- The RBP is based on the results of in vitro testing. At least 99.9 percent of the balloons (with a 95 percent confidence) will not burst at or below their RBP. Use of a pressure monitoring device is recommended to prevent over pressurization.
- The Stent Diameter values listed are actual average stent diameters at the specific balloon inflation pressures obtained during in vitro testing at 37°C and **do not take into account lesion resistance. The stent sizing should be confirmed angiographically..**

10 HOW SUPPLIED

CONTENT

1. Coronary stent system with rapamycin release into the vessel wall Firehawk Liberty™, consisting of:

- Drug-eluting stent
- Delivery System: Rapid Change (RX) Balloon Catheter

2. Balloon extensibility table

3. Instructions for use

STERILE: FOR SINGLE USE ONLY. This product is sterilized with ethylene oxide (EtO) and is nonpyrogenic.

- Do not use if the package is opened or damaged.
- Do not use if labeling is incomplete or illegible.
- Do not resterilize.

If the packaging of the delivered product is opened or damaged, return the product to the manufacturer or to an authorized representative of the manufacturer. Please contact your local Shanghai MicroPort Medical (Group) Co., Ltd. representative for product return information.

11 PERSONNEL REQUIREMENTS

The Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System can be operated by personnel who must be well trained and operate the stent system under the guidance of a suitable healthcare professional familiar with the general benefits and risks of stent system operation..

12 CONDITIONS OF APPLICATION

Scope: medical and treatment-and-prophylactic medical institutions. Use at temperatures between + 32 ° C and + 42 ° C (human body, bloodstream).

13 STERILIZATION

The product is sterilized with ethylene oxide (EtO).
Do not re-sterilize or reuse the device.

14 STORAGE AND TRANSPORTATION CONDITIONS

The product should be stored between 10 °C to 30 °C (50°F to 86°F), in a dry, clean, well-ventilated area free of corrosive gases.

Products should be transported at temperatures between 10 °C and 30 °C (50°F to 86°F). Transport temperature allowed up to 50 °C (122°F) for 2 weeks (products must be protected from pressure, direct sunlight and rainwater or as per contract).

15 DISPOSAL INSTRUCTIONS

Decommissioning of the product is not required.

Disposal of the medical device must be carried out in accordance with the current legislation. When using a medical device, waste is generated that is classified and destroyed (disposed of) as epidemiologically hazardous waste (SanPiN 2.1.3684-21: class B). Products that have become unusable or with an expired shelf life are subject to destruction in accordance with SanPiN 2.1.3684-21 as waste belonging to class "A" (epidemiologically safe waste, close to MSW in composition), in a way that prevents their reuse ...

Dispose of hazardous substances or biocontaminated materials according to your institution's practices. Dispose of all materials in a safe and acceptable manner in accordance with applicable health and safety regulations.

16 SHELF LIFE

Subject to the storage conditions, the shelf life of the product is 2 years.

17 WARRANTY AND LIMITATIONS

Shanghai MicroPort Medical (Group) Co., Ltd. has committed itself to the manufacture of the Firehawk Liberty™ Rapamycin Target Eluting Coronary Stent System. Shanghai MicroPort Medical (Group) Co., Ltd. warrants that all the products shall be free of defects in materials and workmanship upon receipt. Shanghai MicroPort Medical (Group) Co., Ltd. will not be liable for any incidental, special, or consequential loss, damage, or expense resulting, directly or indirectly, from use of its product.

Shanghai MicroPort Medical (Group) Co., Ltd. indicates definitely that this device is intended for one time use only and makes no warranty, representation or condition of any kind, whether expressed or implied (including any warranty of merchantability, suitability or fitness for a particular purpose) respecting the reuse of this product. Moreover, Shanghai MicroPort Medical (Group) Co., Ltd. assumes no responsibility or liability for incidental or consequential damages which may result from such reuse. Contact MicroPort directly with any other questions related to the products.

Firehawk Liberty™ is a trademark of Shanghai MicroPort Medical (Group) Co., Ltd.

If you have any questions regarding the described products, please contact Shanghai MicroPort Medical (Group) Co., Ltd. directly, or an authorized representative of the manufacturer.

18 MANUFACTURER

Shanghai MicroPort Medical (Group) Co., Ltd.

Add: 501 Newton Road ZJ Hi-Tech Park 201203 Shanghai, People's Republic of China

Tel: (86) (21) 38954600

Fax: (86) (21) 50801305

E-mail: customerservice@microport.com

Website: www.microport.com

19 OFFICIAL REPRESENTATIVE IN THE RUSSIAN FEDERATION

Limited Liability Company "ETS-PRO" (LLC "ETS-PRO")

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21 MARKING AND PACKAGING

The Firehawk Liberty™ Rapamycin In-Wall Coronary Stent System is available in five layers and is labeled accordingly:

Layer 1 Coiled tubing: Protects the stent system from kinking and crushing.

Layer 2 The inner bag is a sterile barrier system that prevents microbiological contamination of the product.

Layer 3 Aluminum Foil Bag: Protects the drug-eluting stent from light, oxygen and moisture.

Layer 4 Consumer packaging (cardboard box): Protects packaged products from damage.

Layer 5 Transport Packaging: Protects packaged items during transport.

The explanation of symbols and inscriptions on the marking is given in Section 22 of these instructions.

➤ The label of the inner bag, aluminum foil bag, consumer packaging is additionally equipped with six (6) stickers - adhesive recording labels for attachment to the medical history of the patient receiving the implant (see Figure 2). The stickers contain the following information:

- a) identification of the manufacturer: the logo of the manufacturer "MicroPort";
- b) product information: "Firehawk Liberty" trademark, stent size (diameter x length) mm;
- c) symbol "batch number";
- d) symbol "catalog number"
- f) symbol "expiry date"
- f) product barcode


Figure 2. The layout of the sticker (is an integral part of the marking)

Stickers are part of the labels (perforation is applied along the edge of each sticker on the label - punching holes along the tear line, which makes it easy to separate the sticker from the marking and paste it into the patient's medical history).

➤ The labeling of each stent contains information about:

- the trade mark of the stent "Firehawk Liberty"
- information about the size of the stent in the format, in mm: Ø diameter x length

22 GRAPHIC SYMBOLS AND LABELING

Graphic symbols on the label:

	Do not re-sterilize.		Manufacturer.
	Do not use if package is damaged.		Authorized representative in the European Community.
	Keep away from sunlight.		Date of manufacturer.
	Keep dry.		Use-by date.
	Do not re-use.		Batch code.
	Consult instructions for use.		Catalog number.
	Content: One Unit		Sterilized using ethylene oxide
	MRI conditional.		Non-pyrogen.
	Temperature limit for storage: 10°C~30°C.		CE mark.

Decoding of inscriptions on the marking

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Registration certificate No. _____ dated XX.XX.20XX

I affirm
Dr.Zhaohua Chang
Chairman&CEO
Shanghai MicroPort Medical (Group) Co., Ltd

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