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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Chang Zhaohua, the legal representative of Shanghai MicroPort Medical (Group) Co., Ltd. and the seal of the said company on the annexed INSTRUCTIONS FOR USE are genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Zhou Pingfan

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(Date: Oct. 29, 2019)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Instructions for use

Firehawk™ Rapamycin Target Eluting Coronary Stent System

Document: IFU-03-2019. Rev.D

APPROVED

Shanghai MicroPort Medical (Group) Co., Ltd.
Chen Hua Chang

09.2019

Shanghai MicroPort Medical (Group) Co., Ltd.

NOTICE: The information contained herein is CONFIDENTIAL and PROPRIETARY. The contents of this document are not to be reproduced or made available to third parties without prior consent from MicroPort Medical, and are not to be used in any unauthorized way.

This document has been prepared to provide information for the registration of a medical device in the territory of the Russian Federation. Information is presented in the amount required in accordance with the Rules for the registration of medical devices (Decree of the Government of the Russian Federation dated December 27, 2012 No. 1416), Order of the Ministry of Health of the Russian Federation dated January 19, 2017 No. 11n "On approval of the requirements for the content of technical and operational documentation of the manufacturer (manufacturer)) a medical device.

INSTRUCTIONS FOR USE

Firehawk Rapamycin Target Eluting Coronary Stent System

1. PRODUCT DESCRIPTION

The Firehawk™ Rapamycin Target Eluting Coronary Stent System (Firehawk™ RTECSS or Firehawk™) is a device/drug combination product consisting of

- A Stent and Delivery System: L605 CoCr stent pre-mounted to a RX delivery catheter;
- Drug/Polymer Coating: a formulation of rapamycin in a biodegradable polymer coating which filled in grooves that located on outside surface of stent struts (Fig. 1.1 and Fig. 1.2);


Fig. 1.1 The Schematic diagram of Firehawk™ RTECSS


Рис. 1.2 — Construction of Groove loaded drug product

Stent configuration: shape - cylindrical.

Stent type: Balloon-expandable implant used to treat vascular damage by creating mechanical support to ensure or maintain vessel integrity.

Therapeutic effect: The use of Firehawk™ stent involves minimally invasive interventional therapy during coronary stent implantation to deliver to the lesion site and then expand the stent to open the artery. After the procedure, delivery systems are removed from the body. The stent remains in the vessel for its continued support. Rapamycin is located in the recesses of the stent frame and is slowly released (elutes) to inhibit the proliferation of smooth muscle cells. The coating will degrade and ultimately (9 months after implantation) the stent becomes a holometallic stent. All of these constructs help restore blood flow to the arteries and prevent repeated narrowing or blockage of the treated vessel.

Size range:

The Firehawk™ stent has three sizes: Small size (S) with 6 crowns on a cylindrical ring, designed for stents with nominal diameters of 2.25 mm and 2.5 mm, medium size (M) with 7 crowns on a cylindrical ring, designed for stents with nominal diameters of 2.75 mm and 3.0 mm. Large size (L) with 9 crowns on a cylindrical ring, designed for stents with nominal diameters of 3.5 mm and 4.0 mm. All stent sizes are shown in table 1.1.

Table 1.1 - Dimensions and design of the Firehawk™ stent

Diameter (mm)	Length (mm)											Stent size
	13	16	18	21	23	26	29	31	33	35	38	
Ø2,25	√	√	√	√	√	√	√	—	—	—	—	S
Ø2,5	√	√	√	√	√	√	√	√	√	—	—	
Ø2,75	√	√	√	√	√	√	√	√	√	√	√	M
Ø3,0	√	√	√	√	√	√	√	√	√	√	√	
Ø3,5	√	√	√	√	√	√	√	√	√	√	√	L
Ø4,0	√	√	√	√	√	√	√	√	√	√	√	

"—" - not available

The characteristics of the device are presented in tables 1.2 and 1.3.

Table 1.2 - Description of the Firehawk™ Stent

Characteristic	Value			
Delivery system	Rapid-Exchange (RX) type balloon catheter			
Delivery system design	Single-channel port to the clearance for pumping. The outlet for the conductor is located at a distance of 25 cm from the tip. For conductors with a diameter ≤ 0.014 inches (0.356 mm)			
Available Stent Lengths (mm)	13, 16, 18, 21, 23, 26, 29, 31, 33, 35, 38			
Available Stent Diameters (mm)	2,25; 2,5; 2,75; 3,0; 3,5; 4,0			
Appearance	The surface of the delivery system should be smooth without any visible foreign particles; There should be no twisting (especially around the guide hole). There should be no damage to the tip. Stent carcass elements do not wobble or overlap, the carcass does not tip over, there are no cracks, and there are no visible foreign substances.			
Inner diameter and length of expanded stent at nominal pressure	Nominal Diameter (mm)	The diameter of the expanded stent (mm)	Nominal length (mm)	Length of expanded stent at nominal pressure (mm)
	Ø2,25	2,25±0,22	13	13±3,0
	Ø2,5	2,50±0,25	16	16±3,0
	Ø2,75	2,75±0,27	18	18±3,0
	Ø3,0	3,00±0,30	21	21±3,0
	Ø3,5	3,50±0,35	23	23±3,0
	Ø4,0	4,00±0,40	26	26±3,5
	-	-	29	29±3,5
	-	-	31	31±3,5
	-	-	33	33±3,5
	-	-	35	35±3,5
	-	-	38	38±3,5
Stent Material	L-605 cobalt chromium (CoCr) alloy			
Drug coating	The abluminal coating of poly-DL-lactide (DL-PLA) and rapamycin (CAS 53123-88-9), placed in grooves with a maximum nominal drug content of 125 µg at the largest stent (4.0 × 38 mm). Figure 1.2			
Balloon material	Polymer brand Pebax			
X-ray contrast strip	Pt-Ir Alloy			

Characteristic		Value
The thickness (depth) of the drug coating		The drug (rapamycin in a biodegradable polymer coating) is located in the grooves located on the outer surface of the stent frame. The depth of the grooves is 30 microns.
Distal coating		Hydrophilic
Catheter Effective Working Length		1400 mm
Soft tip length		3,0 mm
Soft Tip Outer Diameter		0.5334 mm
Distances from the distal tip of the catheter	Guide hole position	250 mm (tolerance $\pm$ 50 mm)
	First mark position	950 mm (tolerance $\pm$ 50 mm)
	Second mark position	1050 mm (tolerance $\pm$ 50 mm)
Stent Platform Thickness		For size S, M: 86.4 μ m, for size L: 96.5 μ m
Stent length and diameter ratio		When the stent is stretched from minimum to maximum diameter, the change in stent length should be $\leq$ 20%.
Stent delivery system balloon		A semi-compliant balloon with two radiopaque markers to designate the stent.
Balloon Inflation Pressure		Nominal inflation pressure: 10 atm (1013 kPa) for all sizes; Rated Burst Pressure (RBP): 16 atm (1621 kPa) for 2.25 to 3.5mm Diameters; 14 atm (1419 kPa) for 4.0mm Diameters
Guide Catheter Diameter		The maximum diameter of the compatible catheter guidewire is $\leq$ 0.014 in. (0.36 mm) and the minimum internal diameter of the compatible guide catheter is $\geq$ 1.8 mm (6 F).
Bending / twisting		The minimum radius of curvature of the expanded stent should be $\leq$ 8.0 mm.
Surface area without stent and external stent area		For stents with the declared diameter, minimum diameter and maximum diameter, the surface area without a stent should be more than 70% of the total surface area of the cylindrical part with the same diameter and length.
Magnetic resonance (MR) safety and compatibility		Magnetically induced displacement force: The angle of magnetic deflection should be less than the specific weight of the stent.

Characteristic	Value										
	Moment: The magnetic moment should be less than the specific gravity of the stent $\times$ length of the stent. RF heating: the maximum local temperature rise should not exceed 4 ° C. Artifacts: If the distance of the artifacts L outside the stent lumen is less than 2 mm, the quality of the MR image will be slightly worse. If L is greater than 2 mm, the stent will affect the quality of the MRI.										
Effect / Extension Profile	After the control process, the delivery system has no visible damage (for example, kinks), the stents do not detach and do not move. The distance between the outer diameter of the balloon and the outer diameter of the cross section of the stent should be $\leq$ 5% both before and after the control, and the curvature should be more than 8 mm.										
Adherence to the vessel wall	The outer surface of the stent should fit snugly against the vessel wall, and tilt and separation of the stent from the vessel wall is not allowed.										
Stent profile	Profile of the pressed stent: <table> <tr> <th>Nominal diameter</th><th>Profile</th></tr> <tr> <td>$\varnothing 2,25, \varnothing 2,5$ mm</td><td>$\leq 0,042''$ (1,06 mm)</td></tr> <tr> <td>$\varnothing 2,75$ mm</td><td>$\leq 0,043''$ (1,09 mm)</td></tr> <tr> <td>$\varnothing 3,0$ mm</td><td>$\leq 0,044''$ (1,12 mm)</td></tr> <tr> <td>$\varnothing 3,5, \varnothing 4,0$ mm</td><td>$\leq 0,051''$ (1,30 mm)</td></tr> </table>	Nominal diameter	Profile	$\varnothing 2,25, \varnothing 2,5$ mm	$\leq 0,042''$ (1,06 mm)	$\varnothing 2,75$ mm	$\leq 0,043''$ (1,09 mm)	$\varnothing 3,0$ mm	$\leq 0,044''$ (1,12 mm)	$\varnothing 3,5, \varnothing 4,0$ mm	$\leq 0,051''$ (1,30 mm)
Nominal diameter	Profile										
$\varnothing 2,25, \varnothing 2,5$ mm	$\leq 0,042''$ (1,06 mm)										
$\varnothing 2,75$ mm	$\leq 0,043''$ (1,09 mm)										
$\varnothing 3,0$ mm	$\leq 0,044''$ (1,12 mm)										
$\varnothing 3,5, \varnothing 4,0$ mm	$\leq 0,051''$ (1,30 mm)										
Stent displacement force	The initial bias force must be at least 2.0 N with an initial movement of 0-5 mm (the initial bias force is the force required to create the bias).										
A leak	When filling the system with a contrast agent, there should be no leaks or damage, such as tears or protrusion of the delivery system. There should be no leaks at the Luer tip, in any area of the joint or any other elements. Air must not come from the assembly points of the delivery system during the suction process.										

Characteristic	Value																				
Bending / twisting	The stent / delivery system should not have kinks during delivery and retrieval in a simulated use test.																				
Conductor Guide Lumen Smoothness	The conductor can move smoothly without any obstacles through the clearance of the guide; maximum pushing force shall not exceed 0.3 N.																				
Excessive balloon expansion at both ends of the stent	The diameter of the balloon extending beyond the ends of the stent should be no more than 5% of the outer size of the stent during DDR.																				
Hemostasis	At a pressure of 25 kPa and a pushing force of 2 ± 0.2 N, the volume of water flowing out of the hemostatic valve in 20 minutes should not exceed 10 ml.																				
Balloon filling time	Tank filling time ≤ 20 seconds																				
Balloon Duration	Duration of cylinder emptying ≤ 20 seconds																				
Drug release	<table> <tr> <th>Time point</th><th>Release range</th></tr> <tr> <td>10 min</td><td>$\leq 35\%$</td></tr> <tr> <td>1</td><td>22%~58%</td></tr> <tr> <td>6</td><td>41%~82%</td></tr> <tr> <td>48</td><td>$\geq 66\%$</td></tr> </table>	Time point	Release range	10 min	$\leq 35\%$	1	22%~58%	6	41%~82%	48	$\geq 66\%$										
Time point	Release range																				
10 min	$\leq 35\%$																				
1	22%~58%																				
6	41%~82%																				
48	$\geq 66\%$																				
Joint strength and tube tensile strength	<table> <tr> <th>Validation items</th><th>Minimum Destructive Force (N)</th></tr> <tr> <td>Connection of the distal part of the cylinder</td><td>2</td></tr> <tr> <td>Connection of the proximal part of the balloon</td><td>5</td></tr> <tr> <td>The distal part of the shaft</td><td>5</td></tr> <tr> <td>Adapter tube</td><td>5</td></tr> <tr> <td>Middle part</td><td>3</td></tr> <tr> <td>Conductor outlet</td><td>5</td></tr> <tr> <td>Sheath compound</td><td>5</td></tr> <tr> <td>Connection between Luer tip and strain relief tube</td><td>10</td></tr> <tr> <td>The connection between the proximal part of the rod and the tube to relieve</td><td>10</td></tr> </table>	Validation items	Minimum Destructive Force (N)	Connection of the distal part of the cylinder	2	Connection of the proximal part of the balloon	5	The distal part of the shaft	5	Adapter tube	5	Middle part	3	Conductor outlet	5	Sheath compound	5	Connection between Luer tip and strain relief tube	10	The connection between the proximal part of the rod and the tube to relieve	10
Validation items	Minimum Destructive Force (N)																				
Connection of the distal part of the cylinder	2																				
Connection of the proximal part of the balloon	5																				
The distal part of the shaft	5																				
Adapter tube	5																				
Middle part	3																				
Conductor outlet	5																				
Sheath compound	5																				
Connection between Luer tip and strain relief tube	10																				
The connection between the proximal part of the rod and the tube to relieve	10																				

Characteristic	Value															
	tension															
Visibility	X-ray contrast strips of a balloon catheter should be visible during fluoroscopy.															
Compression Resistance at Radially Applied Load	Compression resistance values for stents of different diameters with 12% strain: <table><tr><th>Nominal Diameter (mm)</th><th>12% deformation and compression resistance HF1 (kPa)</th></tr><tr><td>Ø2,25</td><td>≥ 110</td></tr><tr><td>Ø2,5</td><td>≥ 105</td></tr><tr><td>Ø2,75</td><td>≥ 100</td></tr><tr><td>Ø3,0</td><td>≥ 95</td></tr><tr><td>Ø3,5</td><td>≥ 90</td></tr><tr><td>Ø4,0</td><td>≥ 85</td></tr></table>		Nominal Diameter (mm)	12% deformation and compression resistance HF1 (kPa)	Ø2,25	≥ 110	Ø2,5	≥ 105	Ø2,75	≥ 100	Ø3,0	≥ 95	Ø3,5	≥ 90	Ø4,0	≥ 85
Nominal Diameter (mm)	12% deformation and compression resistance HF1 (kPa)															
Ø2,25	≥ 110															
Ø2,5	≥ 105															
Ø2,75	≥ 100															
Ø3,0	≥ 95															
Ø3,5	≥ 90															
Ø4,0	≥ 85															
Reduction rate	≤ 10%															
Stent shortening	Nominal Length (mm)	Stent shortening														
	10≤L<20	≤ 12%														
	20≤L<30	≤ 10%														
	L≥30	≤ 8%														
Estimated fatigue characteristics of the cylinder	The delivery system must not leak or break after 10 cycles of filling and emptying.															
Luer connector	Outer diameter 4 mm, inner diameter 2 mm, 6% bevel of the conical part of the connection, no thread															
Catheter Shaft Outer Diameter (nominal)	<u>Ø2,25-3,0 mm</u>	<u>Ø3,5-4,0 mm</u>														
	Distal:	2.7F(0.035") 2.9F(0.038")														
	Proximal:	2.3F(0.030") 2.6F(0.034")														
Sterilization packaging	The absence of defects on the outer surface of the inner bag, such as unevenness of the material of the sterile barrier or on it, for example, discoloration, tears, cracks, holes or breaks; continuity of the weld (presence of channels, lack of welding or incomplete welding).															

Table 1.2 — Firehawk RTECSS Product Matrix and Rapamycin Content

Product Code	Stent Diameter (mm)	Stent Length (mm)	Drug Content (µg)
RV2213	2,25	13	30
RV2513	2,5	13	30
RV2713	2,75	13	35
RV3013	3,0	13	35
RV3513	3,5	13	45
RV4013	4,0	13	45
RV2216	2,25	16	38
RV2516	2,5	16	38
RV2716	2,75	16	41
RV3016	3,0	16	41
RV3516	3,5	16	53
RV4016	4,0	16	53
RV2218	2,25	18	43
RV2518	2,5	18	43
RV2718	2,75	18	47
RV3018	3,0	18	47
RV3518	3,5	18	60
RV4018	4,0	18	60
RV2221	2,25	21	48
RV2521	2,5	21	48
RV2721	2,75	21	53
RV3021	3,0	21	53
RV3521	3,5	21	68
RV4021	4,0	21	68
RV2223	2,25	23	53
RV2523	2,5	23	53
RV2723	2,75	23	59
RV3023	3,0	23	59
RV3523	3,5	23	76
RV4023	4,0	23	76
RV2226	2,25	26	60

Product Code	Stent Diameter (mm)	Stent Length (mm)	Drug Content (µg)
RV2526	2,5	26	60
RV2726	2,75	26	68
RV3026	3,0	26	68
RV3526	3,5	26	87
RV4026	4,0	26	87
RV2229	2,25	29	68
RV2529	2,5	29	68
RV2729	2,75	29	76
RV3029	3,0	29	76
RV3529	3,5	29	98
RV4029	4,0	29	98
RV2531	2,5	31	73
RV2731	2,75	31	79
RV3031	3,0	31	79
RV3531	3,5	31	102
RV4031	4,0	31	102
RV2533	2,5	33	78
RV2733	2,75	33	85
RV3033	3,0	33	85
RV3533	3,5	33	110
RV4033	4,0	33	110
RV2735	2,75	35	91
RV3035	3,0	35	91
RV3535	3,5	35	117
RV4035	4,0	35	117
RV2738	2,75	38	97
RV3038	3,0	38	97
RV3538	3,5	38	125
RV4038	4,0	38	125

Stent functioning factors:

- The degree of probability of occurrence of contact, galvanic, pitting and fretting corrosion, the possibility of oxidation: the probability is low, provided by the use of bioinert materials
- Danger of tissue perforation: Perforation or rupture of tissues can occur if a stent is used that is not suitable for the diameter of the vessel, or if precautions regarding the use of the stent are not followed.
- The level of force transmitted by the stent to the surrounding tissue: 0.10-0.15 N / mm
- The period of expected stent ingrowth, penetration into tissues, tissue perforation, stent displacement and migration - about 4 weeks
- Penetration in tissue 50%
- Information on biodegradation: the probability is low, provided by the use of bioinert materials.
- Minimum force to remove the system: at least 2 N.

Type of contact with the body

Stent - lifelong contact with the internal environment of the body (with blood or lymph).

Fast change balloon catheter for stent delivery - short-term contact with blood and lymph.

2. INDICATIONS AND SCOPE

Indications: The Firehawk™ Rapamycin Target Eluting Coronary Stent System (Firehawk™ RTECSS) is indicated for improving coronary luminal diameter in patients with symptomatic heart disease due to de novo native coronary artery lesions length ≤ 60 mm with reference vessel diameters of ≥ 2.25 to ≤ 4.0 mm.

Scope: Medical and medical-preventive medical institutions.

3. USER

The Firehawk™ Stent should be used by qualified medical personnel who are familiar with the general benefits and risks of using the stent system.

4. INDICATIONS AND CONTRAINDICATIONS**INDICATIONS:**

The Firehawk™ Rapamycin Target Eluting Coronary Stent System (Firehawk™ RTECSS) is indicated for improving coronary luminal diameter in patients with symptomatic heart disease due to de novo native coronary artery lesions length ≤ 60 mm with reference vessel diameters of ≥ 2.25 to ≤ 4.0 mm.

CONTRAINDICATIONS:

The Firehawk Rapamycin Target Eluting Coronary Stent System (Firehawk TM RTECSS) is contraindicated for use in patients:

- Who cannot receive antiplatelet and/or anti-coagulant therapy
- With lesions that prevent complete angioplasty balloon inflation or proper placement of the stent or stent delivery system

With hypersensitivity or contraindication to rapamycin or similar drugs or any other analogue or derivative, cobalt, chromium, nickel, tungsten, or PLA.

5. WARNINGS

- Ensure that the inner package sterile barrier has not been opened or damaged prior to use.
- Judicious patient selection is necessary since the use of the device carries the associated risks of thrombosis, vascular complications, and/or bleeding events.
- This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy.

6. PRECAUTIONS**6.1. General Precautions**

- Stent implantation should only be performed by physicians who have received appropriate training.
- Stent placement should be performed at hospitals where emergency coronary artery bypass graft surgery is accessible.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. Long-term outcomes following repeat dilatation of the stent are presently unknown.
- Risks and benefits should be considered in patients with severe contrast agent allergies.
- Do not expose or wipe the product with organic solvents such as alcohol or detergents.
- Care should be taken to control the guiding catheter tip during stent delivery, deployment, and delivery system withdrawal. Before withdrawing the stent delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter advancement into the vessel and subsequent arterial damage.
- Compared to use within the specified Indications for Use, the use of DES in patients and

lesions outside of the labeled indications may have an increased risk of adverse events, including stent thrombosis, stent embolization, MI, or death.

- It is very important that the patient comply with post-procedural antiplatelet therapy recommendations. It is strongly recommended for 12 months in patients who are not at high risk of bleeding per the ACC/AHA/SCAI guidelines.
- Early discontinuation of prescribed antiplatelet medications could result in a higher risk of stent thrombosis, MI, or death. The risks and benefits of the procedure should be weighed against the possible risks associated with early discontinuation of antiplatelet therapy. Patients who require early discontinuation of antiplatelet therapy should be monitored carefully for cardiac events.
- Orally administered sirolimus combined with cyclosporine is associated with increased serum cholesterol and triglycerides levels.
- When used with cyclosporine medication, The sirolimus mean AUC and mean Cmax may be affected

6.2. Use of Multiple Stents

When multiple stents are required, stent materials should be of similar composition. Placing multiple stents of different materials in contact with each other may increase potential for corrosion. Potential interactions of the Firehawk™ stent with other drug-eluting or coated stents have not been evaluated.

A patient's exposure to drug and polymer is proportional to the number and total length of implanted stents. Use of more than two Firehawk stents has not received adequate long term clinical evaluation. The maximum amount of stents used needs to be determined by the physician based on the risks versus the clinical benefit for the patient involved.

6.3. Brachytherapy

The safety and effectiveness of the Firehawk™ stent in patients with prior brachytherapy of the target lesion have not been established. The safety and effectiveness of the use of brachytherapy to treat in-stent restenosis in a Firehawk stent has not been established. Both vascular brachytherapy and the Firehawk™ stent alter arterial remodeling. The synergy between these two treatments has not been determined.

6.4. Drug Interactions

The effect of drug interactions on the safety or efficacy of the Firehawk™ Stent has not been determined. Several drugs are known to affect the metabolism of Rapamycin, and other drug interactions may be inferred from known metabolic effects. Rapamycin is known to be a substrate for both cytochrome P450 IIIA4 (CYP3A4) and P-glycoprotein.

Consideration should be given to the potential for drug interaction when deciding to implant Firehawk™ Stent if patient is taking a drug that might interact with Rapamycin, or when deciding to initiate therapy if patient has recently received a Firehawk™ Stent.

6.5. Magnetic Resonance Imaging (MRI)

Non-clinical testing has demonstrated the Firehawk™ stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla
- Spatial gradient field of 36 T/m
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg or less at 3 Tesla, Firehawk™ stent, for 15 minutes of continuous MR scanning.

In non-clinical testing, the single and overlapping (72.4 mm) Firehawk™ stent respectively produced a temperature rise of 1.2 °C and 3.2 °C at a maximum whole body averaged

specific absorption rate (SAR) of ≈ 3.53 W/kg, as assessed by calorimetry for 15 minutes of continuous MR scanning with whole body coil in a 3 Tesla Magnetom Trio, Siemens Medical Solutions (Software: Numaris/4, syngo MR A30) MR Scanner.

Non-clinical testing has not been performed to rule out the possibility of implant migration and RF heating at static gradient magnetic fields stronger than mentioned above.

MR imaging quality may be compromised if the area of interest is in exactly the same area or relatively close to the position of the stent.

6.6. Stent Handling Precautions

- Each stent is for single use only. Do not resterilize or reuse this device. Note the "Use-by" (expiration) date on the product label.
- The foil pouch is not a sterile barrier. The inner header bag (pouch) within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.
- The foil pouch contains 1 packet of oxygen absorber. Please discard this packet after opening the pouch.
- Do not remove the stent from the delivery system as removal may damage the stent and/or lead to stent embolization. These components are intended to perform together as a system.
- Do not induce a vacuum on the delivery system prior to reaching the target lesion. This may cause dislodgement of the stent from the balloon.
- The delivery system should not be used in conjunction with other stents.
- Special care must be taken not to handle or disrupt the stent on the balloon especially during delivery system removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- Do not manipulate, touch, or handle the stent with your fingers, which may cause contamination or stent dislodgement from the delivery balloon.
- Use only the appropriate balloon inflation media (see Section 9.3.4 – Operator's Instructions, Delivery System Preparation). Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in stent deployment.

6.7. Stent Placement Precautions

6.7.1. Stent Preparation

• Do not pre-inflate the delivery system prior to stent deployment other than as directed. Use the balloon purging technique described in Section 9.3.4 – Operator’s Instructions, Delivery System Preparation.

• Use guiding catheters which have lumen sizes that are suitable to accommodate the stent delivery system (see Section 1.1 – Product Description, Device Component Description).

6.7.2. Stent Implantation

• The vessel should be pre-dilated with an appropriate sized balloon. Failure to do so may increase the difficulty of stent placement and cause procedural complications.

• Do not expand the stent if it is not properly positioned in the vessel (see Section 5.8 – Precautions, Stent/System Removal).

• Implanting a stent may lead to vessel dissection and acute closure require additional intervention (CABG, further dilatation, placement of additional stents, or other).

• Although the safety and effectiveness of treating same coronary artery with more than one Firehawk TM stents has not been established, if this is performed, place the stent in the distal lesion before the proximal lesion in order to minimize dislodgement risk incurred by traversing through deployed stents.

• Stent placement may compromise side branch patency.

• Do not exceed Rated Burst Pressure (RBP) as indicated on product label. See Table 9.1, Typical Firehawk TM RTECSS Compliance. Balloon pressures should be monitored during inflation. Applying pressures higher than specified on the product label may result in a ruptured balloon with possible arterial damage and dissection.

• Do not attempt to pull an unexpanded stent back through the guiding catheter, as dislodgment of the stent from the balloon may occur. Remove as a single unit as indicated in 5.8 Stent/System Removal Precautions.

• Stent retrieval methods (i.e., using additional wires, snares, and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.

• Although the stent delivery system balloon is strong enough to expand the stent without rupture, a circumferential balloon tear distal to the stent and prior to complete stent expansion, could cause the balloon to become tethered to the stent, requiring surgical removal. In case of balloon rupture, it should be withdrawn and, if necessary, a new dilatation catheter exchanged over the guide wire to complete the expansion of the stent.

- Ensure the stented area covers the entire lesion/dissection site and that no gaps exist between stents.

6.8. Stent/System Removal Precautions

If any unusual resistance is felt at any time during lesion access before stent implantation, the stent system and the guiding catheter should be removed as a single unit.

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- DO NOT retract the stent system into the guiding catheter.
- Position the proximal balloon marker just distal to the guiding catheter tip.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the delivery system to the guiding catheter, and remove the guiding catheter and stent system as a single unit.

Failure to follow these steps and/or applying excessive force to the delivery system can potentially result in loss or damage to the stent and/or delivery system components and/or vasculature.

If it is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

6.9. Post-Procedure Precautions

- Care should be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guidewire, or balloon catheter to avoid disrupting the stent placement, apposition, geometry.
- Post-dilatation: All efforts should be made to assure that the stent is not under dilated. If the deployed stent is not fully apposed to the vessel wall, the stent may be expanded further with a larger diameter balloon. The post-dilatation can be done using a low-profile, high pressure, non-compliant balloon catheter. The balloon should not extend outside of the stented region (see Operator's Manual – 9.7 Post-Deployment Dilatation of Stent Segments). Do not use the stent delivery balloon for post-dilatation.
- Non-clinical testing has demonstrated that the Firehawk stent is MR Conditional (see Section

5.5 Magnetic Resonance Imaging (MRI)). MR imaging quality may be compromised if the area of interest is in the same area of the position of the stent.

• Antiplatelet therapy should be administered post-procedure (see Section 5.1 General

Precautions). Patients who require early discontinuation of antiplatelet therapy (e.g., secondary to active bleeding) should be monitored carefully for cardiac events. At the discretion of the patient's treating physician, the antiplatelet therapy should be restarted as soon as possible.

7. ADVERSE EVENTS

Adverse events may be associated with the implantation of a coronary stent in coronary arteries, but are not limited to the following:

- Allergic reaction to anti-coagulant and/or antiplatelet therapy, contrast medium, or stent
- materials
- Arrhythmias
- Arteriovenous fistula
- Bleeding
- Cardiac tamponade
- Coronary aneurysm
- Death
- Dissection
- Drug interactions with antiplatelet/anticoagulant/contrast medium
- Emboli, distal (tissue, air, or thrombotic emboli)
- Embolization, stent
- Emergency CABG
- Failure to deliver the stent to the intended site
- Fever
- Heart failure
- hemorrhage
- Hypotension/Hypertension
- Infection, local or systemic
- Myocardial infarction
- Pain, at the access site
- Pseudoaneurysm, femoral
- Restenosis of stented segment

- Stent embolization or migration
- Stent fracture
- Stent thrombosis/occlusion
- Target lesion revascularization
- C0C03-016 Rev.D Page 10 / 18
- Target vessels of non-target lesion revascularization
- Total occlusion of coronary artery
- Vessel trauma requiring surgical repair or reintervention

Potential adverse events not captured above, that may be unique to the rapamycin drug coating

- Abnormal liver function tests
- Allergic/immunologic reaction to drug (rapamycin or structurally-related compounds) or the
- polymer stent coating or its individual components (see Section 1 Product Description)
- Anemia
- Arthralgias
- Diarrhea
- Hypercholesterolemia
- Hypersensitivity, including anaphylactic/anaphylactoid type reactions
- Hypertriglyceridemia
- Hypokalemia
- Infections
- Interstitial lung disease
- Myocardial infarction
- Myocardial ischemia
- Occlusion
- Prolonged angina
- Pseudoaneurysm
- Hematologic dyscrasia (including leukopenia, neutropenia, thrombocytopenia)
- Renal failure
- Restenosis of stented segment (greater than 50% obstruction)
- Stroke
- Vessel spasm
- Vessel perforation

8. INDIVIDUALIZATION OF TREATMENT

The risks and benefits described above should be considered carefully for each patient before use of the Firehawk™ stent. Patient selection factors to be assessed should include a judgment regarding risk of anti-platelet therapy. Special consideration should be given to those patients with recently active gastritis or peptic ulcer disease.

Premorbid conditions that increase the risk of a poor initial result and the risks of emergency referral for bypass surgery (diabetes mellitus, renal failure, and severe obesity) should be reviewed. The relation of baseline and procedural variables to Major Adverse Cardiac Events (MACE) was examined. Multivariable modeling suggested that treatment assignment remained an independent predictor of clinical and angiographic outcomes even after adjusting for other baseline and procedural confounding variables.

9. USE IN SPECIAL POPULATIONS

9.1. Pregnancy

There are no adequate data from the use of sirolimus or Firehawk™ Stents in pregnant women. In animals, sirolimus was embryo- and fetus- toxic at clinically relevant exposures. Teratogenic effects have not been observed. The potential risk for humans is unknown. The Firehawk stent should not be used in pregnant women, unless the clinical condition of the mother requires treatment with the stent.

9.2. Lactation

It is not known whether sirolimus is excreted in human breast milk. Sirolimus is excreted in milk of lactating rats. A risk to the suckling child cannot be excluded. A decision must be made whether to discontinue breast-feeding or to implant the stent taking into account the importance of the stent to the mother.

9.3. Pediatric Use

The safety and effectiveness of the Firehawk™ stent in pediatric patients below the age of 18 years have not been established.

9.4. Geriatric Use

Clinical studies of the Firehawk™ stent did not suggest that patients age 65 years and over differed with regard to safety and efficacy compared to younger patients.

Carefully consider whether it is appropriate to use in the above patient populations.

10. OPERATOR'S INSTRUCTIONS

10.1. Inspection Prior to Use

Carefully inspect the sterile package before opening. Do not use if the package has been damaged or opened. Do not use the product after the "Use-by" date.

10.2. Required Materials

Ancillary devices required to work with the Firehawk™ stent are presented in table 1.4

Table 1.4. List of required assistive devices

Quantity	Material
1	Guiding catheter [≥ 6 F (1.8 mm, 0.070 inch) inner diameter]
2-3	20 cc syringe
1,000 u/500 cc	Sterile Heparinized Normal Saline (HepNS)
1	Guidewire [≤ 0.014 inch (0.356 mm) diameter]
1	Rotating hemostatic valve with appropriate minimum inner diameter (0.096 inch [2.44 mm])
N/A	Contrast medium diluted 1:1 with heparinized normal saline
1	Inflation device
1	Stopcock (3-way minimum)
1	Torque device
1	Pre-deployment dilatation catheter
1	Appropriate arterial sheath
N/A	Appropriate anticoagulation and antiplatelet drugs

10.3. Preparation

10.3.1 Packaging Removal

Note: The foil pouch is not a sterile barrier. The inner pouch within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner pouch is NOT sterile.

- Peel-open the inner pouch and pass or drop the product into the sterile field using an aseptic technique.
- Carefully remove the stent system from coiled hoop and inspect it for bends, kinks, and other damage. Do not use the product if any damage is noted.
- Carefully remove the protective stent sheath and mandrel together, visually inspect the stent to ensure that it has not been damaged or displaced from its original position (between

- proximal and distal marker bands) on the balloon.

Note: do not manipulate, touch, or handle the stent with your fingers, which may cause contamination or stent dislodgement from the delivery balloon.

10.3.2. Rinse Catheter

- Rinse the catheter with sterile heparinized normal saline solution.

10.3.3. Guide Wire Lumen Flush

- Attach the syringe with HepNS to the flushing needle packaged with the catheter.
- Insert the needle into the tip of the catheter and flush the guide-wire lumen with HepNS.

Note: AVOID manipulation of the stent during flushing of the guide-wire lumen, as this may disrupt the placement of the stent on the balloon.

10.3.4 Delivery System Preparation

- 1) Prepare the inflation device or syringe with suitable diluted contrast medium.
- 2) Attach the inflation device or syringe to the stopcock; attach to the balloon inflation port hub.
- 3) Open the stopcock to stent delivery system.
- 4) Leave the inflation device or syringe on neutral.

Note: While introducing the delivery system into the vessel, do not induce negative pressure on A-C0C03-016 Rev.D Page 13 / 18 the delivery system. This may cause dislodgement of the stent from the balloon.

10.4. Delivery Procedure

- 1) Prepare the vascular access site according to standard practice.
- 2) Pre-dilate the lesion with a Percutaneous Transluminal Coronary Angioplasty (PTCA) catheter. Limit the longitudinal length of pre-dilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of Firehawk Stent.
- 3) Maintain neutral pressure on the inflation device. Open the rotating hemostatic valve as widely as possible.
- 4) Backload the delivery system onto the proximal portion of the guide-wire while maintaining the guide-wire position across the target lesion.

- 5) Advance the stent delivery system over the guide-wire to the target lesion. Use the radiopaque balloon markers to position the stent across the lesion; perform angiography to confirm the position of the stent.
- 6) Sufficiently tighten the rotating hemostatic valve. The stent is now ready to be deployed.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, do not force passage. Remove the entire system as a single unit. See Section 5.8 – Precautions, Stent/System Removal for specific delivery system removal instructions.

10.5. Deployment Procedure

- 1) Before deployment, reconfirm the correct position of the stent relative to the target lesion via the radiopaque balloon markers.
- 2) Attach the inflation device (only partially filled with contrast media) to three-way stopcock and apply negative pressure to purge the balloon of air.
- 3) Turn the stopcock on the catheter to the off position and purge the inflation device of air. Close the side port of the stopcock.
- 4) Under fluoroscopic visualization, inflate the balloon to at least the nominal pressure to deploy the stent. Do not exceed the labeled RBP. Deploy the stent slowly by pressurizing the delivery system until stent is completely expanded. Maintain pressure for 30 seconds. If necessary, the delivery system can be repressurized or further pressurized to assure complete apposition of the stent to the artery wall. Do not exceed the labeled rated burst pressure (RBP) of 16 atm (1621 kPa) for the 2.25 mm – 3.5 mm stents and 14 atm (1419 kPa) for the 4.0 mm stents sizes.
- 5) Stent wall contact should be verified through routine angiography or intravascular ultrasound.
- 6) Fully cover the entire lesion and balloon treated area (including dissections) with the Firehawk Stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
- 7) If more than one Firehawk stent is needed to cover the lesion and balloon treated area, adequately overlap stents, taking into account stent foreshortening. Ensure no gaps between stents by positioning the balloon marker bands of the second Firehawk stent inside the deployed stent prior to expansion. See Section 5.8 Stent/System Removal

Precautions.

- 8) Deflate the balloon by pulling negative pressure on the inflation device. Allow adequate time, at least 15 seconds, for full balloon deflation. Confirm complete balloon deflation before attempting to remove the delivery system. Ensure the balloon is fully deflated before attempting to move the catheter.
- 9) Confirm that the stent is adequately expanded and in full contact with the artery wall.

10.6. Removal Procedure

- 1) Ensure that the balloon is fully deflated.
- 2) Fully open the rotating hemostatic valve.
- 3) While maintaining the guide wire position and negative pressure on the inflation device, withdraw the delivery system.

Note: Should any resistance be felt at any time during either lesion access or removal of the delivery system post stent implantation, the entire system should be removed as a single unit. See Section 5.8 Stent/System Removal Precautions for specific delivery system removal instructions.

- 4) Tighten the rotating hemostatic valve.
- 5) Repeat angiography to assess the stented area. If post-dilatation is necessary, ensure that the final stent diameter matches the reference vessel diameter. Assure that the stent is not underdilated.

10.7 Post-Deployment Dilatation of Stent Segments

All efforts should be taken to assure that the stent is not underdilated. If the deployed stent size is still inadequate with respect to vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. The stent may be further expanded using a low profile, high pressure, and non-compliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guidewire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region.

CAUTION: Do not dilate the stent beyond the following limits.

Nominal Stent Diameter	Dilatation Limit
2.25 mm to 2.5 mm	2.75 mm
2.75 mm to 3.0 mm	3.5 mm

3.5 mm to 4.0 mm	4.5 mm
------------------	--------

Reconfirm stent position and angiographic results. Repeat inflations until optimal stent deployment is achieved.

10.8 In Vitro Compliance Information

Pressure			Stent					
			Inner Diameter (mm)					
atm	kPa		2,25	2,5	2,75	3,0	3,5	4,0**
10	1013	Nominal	2,20	2,45	2,70	2,95	3,50	3,90
11	1115		2,25	2,49	2,73	3,00	3,54	3,98
12	1216		2,30	2,52	2,78	3,05	3,61	4,06
13	1317		2,33	2,57	2,82	3,09	3,66	4,12
14	1419		2,36	2,61	2,85	3,12	3,70	4,18
15	1520		2,39	2,65	2,88	3,16	3,74	4,22
16	1621	RBP *	2,42	2,68	2,91	3,19	3,78	4,27
17	1723		2,44	2,71	2,94	3,22	3,82	4,32
18	1824		2,47	2,73	2,97	3,25	3,86	4,37

*Do not exceed the rated burst pressure (RBP)

**Rated Burst Pressure (RBP): 14 atm (1419 kPa) for 4.0 mm diameters

Note:

- The RBP is based on the results of in vitro testing. At least 99.9 percent of the balloons (with a 95 percent confidence) will not burst at or below their RBP. Use of a pressure monitoring device is recommended to prevent over pressurization.
- The Stent Diameter values listed are actual average stent diameters at the specific balloon inflation pressures obtained during in vitro testing at 37°C and **do not take into account lesion resistance. The stent sizing should be confirmed angiographically.**

11. HOW SUPPLIED

STERILE: FOR SINGLE USE ONLY.

This product is sterilized with ethylene oxide (EtO) and is nonpyrogenic.

- Do not use if the package is opened or damaged.
- Do not use if labeling is incomplete or illegible.
- Do not resterilize.

Return product if package is opened or damaged to MicroPort Returned Goods. Contact your local MicroPort Representative for return information.

CONTENTS

One (1) Firehawk™ Rapamycin Target Eluting Coronary Stent System

One (1) Instructions for Use

One (1) Patient Information Card

PERFORMANCE OPTIONS

Firehawk™ Rapamycin Target Eluting Coronary Stent System, in the following versions:

RV2213, RV2513, RV2713, RV3013, RV3513, RV4013, RV2216, RV2516, RV2716, RV3016, RV3516, RV4016, RV2218, RV2518, RV2718, RV3018, RV3518, RV4018, RV2221, RV2521, RV2721, RV3021, RV3521, RV4021, RV2223, RV2523, RV2723, RV3023, RV3523, RV4023, RV2226, RV2526, RV2726, RV3026, RV3526, RV4026, RV2229, RV2529, RV2729, RV3029, RV3529, RV4029, RV2531, RV2731, RV3031, RV3531, RV4031, RV2533, RV2733, RV3033, RV3533, RV4033, RV2735, RV3035, RV3535, RV4035, RV2738, RV3038, RV3538, RV4038, Consisting of:

- A balloon catheter for stent delivery.
- Stent balloon-expanding cobalt-chrome coronary.
- Instructions for use.

12. PACKAGING

The packaging system for the Firehawk™ stent consists of five elements (from the inside out): a welded rolled tube (a protective sheath that covers the stent / balloon and inserted into the welded rolled tube), inner packaging (inner bag / sterile barrier system), aluminum packaging (bag of aluminum foil / protective packaging), trade packaging (folding box, protective packaging) and transport packaging (corrugated cardboard).

13. STORAGE

Store between 10°C to 30°C (50°F to 86°F); Transportation excursion is permitted to 50°C (122°F) for 2 weeks. Use by the “Use-by” date noted on the package.

The foil pouch contains 1 packet of oxygen absorber. Please discard this packet after opening the pouch.

14. TRANSPORTATION CONDITIONS

The product should be transported at a temperature of 10 ° C to 30 ° C. Keep away from moisture. Keep away from sunlight.

15. WORKING TEMPERATURE

Use at a temperature of + 32 ° C to + 42 ° C (human body, blood flow)

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On the territory of the Russian Federation, the product should be disposed of in accordance with SanPiN 2.1.7.2790-10 for class B wastes (epidemiologically hazardous wastes). After hardware-based disinfection methods using physical methods and changing the appearance of the waste, which excludes the possibility of their reuse, class B waste can be accumulated, temporarily stored, transported, disposed of and disposed of together with class waste (epidemiologically safe waste similar in composition to solid waste). Packaging of Class B medical disinfected waste should be labeled to indicate that the waste has been decontaminated.

17. SHELF LIFE

Under established storage conditions, the expiration date of the Firehawk™ stent is 2 years.

18. WARRANTY AND LIMITATIONS

Shanghai MicroPort Medical (Group) Co., Ltd. has committed itself to the manufacture of the Firehawk™ Rapamycin Target Eluting Coronary Stent System. Shanghai MicroPort Medical (Group) Co., Ltd. warrants that all the products shall be free of defects in materials and workmanship upon receipt. Shanghai MicroPort Medical (Group) Co., Ltd. will not be liable for any incidental, special, or consequential loss, damage, or expense resulting, directly or indirectly, from use of its product.

Shanghai MicroPort Medical (Group) Co., Ltd. indicates definitely that this device is intended for one time use only and makes no warranty, representation or condition of any kind, whether expressed or implied (including any warranty of merchantability, suitability or fitness for a particular purpose) respecting the reuse of this product. Moreover, Shanghai MicroPort Medical (Group) Co., Ltd. assumes no responsibility or liability for incidental or consequential damages which may result from such reuse. Contact MicroPort directly with any other questions related to the products.

Firehawk™ is a trademark of Shanghai MicroPort Medical (Group) Co., Ltd.

In case of accidents, contact your authorized representative (local distributor).

MANUFACTURER:

Shanghai MicroPort Medical (Group) Co., Ltd.

Add: 501 Newton Rd., Z.J Hi-Tech Park, Shanghai 201203, P. R. China

Tel: (86) (21) 38954600

Fax: (86) (21) 50801305

Email: customerservice@microport.com

Website: www.microport.com

Authorized representative of the manufacturer in Russia

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20. SYMBOLS ON THE MARKING OF A MEDICAL PRODUCT

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Registration certificate No. ____ from XX.XX.20XX

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