
**Instructions for use of the medical device:
Catheter for assisted reproductive technologies,
sterile, disposable**

APPROVE

Allwin Medical Devices Inc., USA

Business Development Manager

Digish Mehta

A circular stamp with the text "ALLWIN MEDICAL DEVICES INC." around the perimeter and "USA" in the center. A handwritten signature is written over the stamp.

17/1/2024

ACKNOWLEDGMENT

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.


State of California
County of Orange


On 11/17/2024 before me, Alexandra S. Zambri, Notary Public
(insert name and title of the officer)

personally appeared Digish Menta
who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature  (Seal)

NAME OF MEDICAL DEVICE

Catheter for assisted reproductive technologies, sterile, disposable, in the composition.

INTENDED USE OF THE DEVICE

The catheter for assisted reproductive technologies is intended use for to introduce washed spermatozoa or embryo in the uterine cavity.

INTENDED USER

The product should be used exclusively by professionals specializing (healthcare professional) in the field of assisted reproductive technologies (ART).

MATERIALS NEEDED IN THE WORK, BUT NOT SUPPLIED WITH A MEDICAL DEVICE

Medical sterile syringe, volume 1 ml to 5 ml, medical sterile gloves, vaginal speculum

There are no performance restrictions for these products. The selection and application of these products is carried out by the laboratory regulations.

OPERATING OF USE (CONDITIONS) OF THE MEDICAL DEVICE

Catheters in operation are to be resistant to the effects of climatic factors according to GOST 15150 for the type of execution U6 (from +32 to +42 C).

AREA OF APPLICATION

The catheter for assisted reproductive technologies are used in specialized laboratories performing assisted reproductive technologies.

INDICATIONS FOR USE

- Infertility of the unknown etiology
- Failure of conception due to tubal blockage (for insemination a necessary is tubal patency of at least one tube)
- Cervical factor infertility
- Male factor infertility
- Combined infertility
- Polycystic ovarian syndrome after failure of conception during a year of the conservative treatment
- Endometriosis

CONTRAINDICATIONS FOR USE

- Cervical infection
- Recent pelvic inflammation.
- Sexually transmitted disease
- Pregnancy
- Confirmed or suspected intrauterine device
- Uterine perforation.

If the technique of intrauterine manipulation is violated, complications may occur such as perforation of the uterine wall, bleeding, endometrial damage.

SIDE EFFECTS

Under use of catheter by qualified specialists in compliance with aseptic technique cannot injure endometrium and vaginal mucous membrane, cause bleeding or infection carry.

Sterile single use catheter for assisted reproductive technologies do not have potential side effect.

CAUTION

- To reduce the risk of perforation, if any resistance is felt while inserting the catheter, do not force the catheter against the resistance.
- The visibility of catheter may be compromised with use of lower solution ultrasound equipment.
- The patients should be with a full bladder for performing transabdominal Ultrasound.

WARNING/ PRECAUTION:

- Do not use if package is opened or damaged.
- For single use only.
- Do not re-sterilize.
- Always use washed spermatozoa. The introduction of unwashed spermatozoa into the uterus may result in severe adverse reactions, including anaphylactic shock. Please refer to published medical literature for methods of preparing spermatozoa for intrauterine artificial insemination before performing the procedure.
- To reduce the risk of perforation, if any resistance is felt while inserting the catheter, do not force the it against the resistance.
- Confirm the specification (length) of the catheter prior to use.
- During insertion do not pull back further than the last graduation of the inner catheter, as removing the inner catheter from the outer catheter may result in damage to the endometrial tissue and bleeding

Type of contact with the human body: short-term contact with the mucous membrane of the vagina and uterus during more than 24 hours.

MODIFICATIONS OF EXECUTION END TECHNICAL CHARACTERISTICS

Intrauterine insemination catheter (single-component catheter)

The name of medical device	Order code	Working length without hub (mm)	Catheter distal end outer diameter of the catheter (mm)	Catheter's volume without hub (ml)	Holes on the tip	Packaging
INSEMINA Intra Uterine Insemination Catheter	I-IUI	170.00 mm± 2.00 mm	2.04 mm± 0.05 mm	0.225 ml± 0.015 ml	Two lateral holes on the tip	25 pcs
INSEMINA Intra Uterine Insemination Catheter, Open tip	I-IUI-O	170.00 mm± 2.00 mm	2.04 mm± 0.05 mm	0.225 ml± 0.015 ml	One hole on the tip	25 pcs
CURVED Intra Uterine Insemination Catheter	I-IUI-C	170.00 mm± 2.00 mm	2.04 mm± 0.05 mm	0.225 ml± 0.015 ml	Two lateral holes on the tip	25 pcs
CURVED Intra Uterine Insemination Catheter, Open tip	I-IUI-C-O	170.00 mm± 2.00 mm	2.04 mm± 0.05 mm	0.225 ml± 0.015 ml	One hole on the tip	25 pcs
TRACKABLE Intra Uterine Insemination Catheter with malleable stylet	T-IUI	170.00 mm± 2.00 mm	2.04 mm± 0.05 mm	0.225 ml ± 0.015 ml	Two lateral holes on the tip	25 pcs
TRACKABLE Intra Uterine Insemination Catheter with malleable stylet, Open tip	T-IUI-O	170.00 mm± 2.00 mm	2.04 mm± 0.05 mm	0.225 ml ± 0.015 ml	One hole on the tip	25 pcs

Stylet for TRACKABLE

The name of medical device	Stylet length without connector (mm)
TRACKABLE Intra Uterine Insemination Catheter with malleable stylet	195.00 mm ± 2.00 mm

Intrauterine insemination catheter (two-part catheter)

The name of medical device	Order code	Working length without hub (loading catheter) (mm)	Catheter distal end inner tube outer diameter (loading catheter) (mm)	Internal catheter volume (loading catheter) without hub (ml)	Holes on the tip	Packaging
FLEXI Intra Uterine Insemination Catheter	F-IUI	180.00 mm $\pm$ 2.00 mm	1.60 mm $\pm$ 0.05 mm	0.138 ml $\pm$ 0.015 ml	Two lateral holes on the tip	25 pcs
FLEXI Intra Uterine Insemination Catheter, Open tip	F-IUI-O	180.00 mm $\pm$ 2.00 mm	1.60 mm $\pm$ 0.05 mm	0.138 ml $\pm$ 0.015 ml	One hole on the tip	25 pcs

Embryo Transfer Catheter

The name of medical device	Order code	Working length without hub (loading catheter) (mm)	Catheter distal end inner tube outer diameter (loading catheter) (mm)	Internal catheter volume (loading catheter) without hub (ml)	Packaging
EMBRYO TRANS Embryo Transfer Catheter	EB-ETC23	230.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.101 ml $\pm$ 0.015 ml	10 pcs
EMBRYO TRANS Embryo Transfer Catheter	EB-ETC18	180.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.079 ml $\pm$ 0.015 ml	10 pcs
ECHO TRANS Embryo Transfer Catheter	ET-ETC23	230.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.101 ml $\pm$ 0.015 ml	10 pcs
ECHO TRANS Embryo Transfer Catheter	ET-ETC18	180.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.079 ml $\pm$ 0.015 ml	10 pcs
TWINKLE EMBRYO TRANS Embryo Transfer Catheter	TEB-ETC23	230.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.101 ml $\pm$ 0.015 ml	10 pcs
TWINKLE EMBRYO TRANS Embryo Transfer Catheter	TEB-ETC18	180.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.079 ml $\pm$ 0.015 ml	10 pcs
BULB TRANS ULTRA Embryo Transfer Catheter Bulb Tip	BT-ETC-U	225.00 mm $\pm$ 2.00 mm	0.95 mm $\pm$ 0.05 mm	0.055 ml $\pm$ 0.015 ml	10 pcs
BULB TRANS STAR Embryo Transfer Catheter Bulb Tip	BT-ETC-U-S	225.00 mm $\pm$ 2.00 mm	0.95 mm $\pm$ 0.05 mm	0.055 ml $\pm$ 0.015 ml	10 pcs
PRO EMBRYO TRANS Embryo Transfer Catheter	PEB-ETC	220.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.099 ml $\pm$ 0.015 ml	10 pcs
PRO ECHO TRANS Embryo Transfer Catheter	PET-ETC	220.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.099 ml $\pm$ 0.015 ml	10 pcs
TWINKLE PRO EMBRYO TRANS Embryo Transfer Catheter	TPEB-ETC	220.00 mm $\pm$ 2.00 mm	1.50 mm $\pm$ 0.05 mm	0.099 ml $\pm$ 0.015 ml	10 pcs

Accessories:

The name of the medical device accessories (REF)	Order code	Length	Packaging
Stylet for catheters 23 cm long EMBRYO TRANS / ECHO TRANS / TWINKLE EMBRYO TRANS	STY23	210.00 mm $\pm$ 1.00 mm	10 pcs
Stylet for catheters 18 cm long EMBRYO TRANS / ECHO TRANS / TWINKLE EMBRYO TRANS	STY18	160.00 mm $\pm$ 2.00 mm	10 pcs
Stylet for BULB TRANS ULTRA Bulb Tip / BULB TRANS STAR Bulb Tip	STY-BU	202.00 mm $\pm$ 2.00 mm	10 pcs

Stylet for PRO EMBRYO TRANS /PRO ECHO TRANS /TWINKLE PRO EMBRYO TRANS	STY-PE	213.00 mm ± 1.00 mm	10 pcs
Soft Obturator for PRO EMBRYO TRANS /PRO ECHO TRANS /TWINKLE PRO EMBRYO TRANS	SB-PE	213.00 mm± 1.00 mm	10 pcs

Description of the accessories

The stylet and the obturator are a stainless steel metal rod coated with a polymer with a rounded tip at one end and a connector at the other.

The accessories (stylet, soft obturator) are used for a successful insertion of pre-washed spermatozoa or cultured embryos through the catheter only in individual cases: specific anatomy of the cervical canal (in case of an anatomically tortuous cervical canal of a complex shape) or large amount of mucus in the cervical canal.

In individual cases, in case of anatomically tortuous cervical canal in a patient, the doctor may decide to use a catheter with a stylet for successful passage through the cervical canal. The stylet gives the catheter additional rigidity and the ability to shape the guiding catheter.

In individual cases, if a patient has large amounts of cervical mucus, the doctor may decide to use a catheter with a soft obturator. Cervical mucus can fill the lumen of the guiding catheter and prevent the passage of the loading catheter with the embryo into the uterine cavity. The soft obturator prevents cervical mucus from entering the conductor catheter and thus facilitates the transfer procedure.

These accessories (stylet, soft obturator) are not used independently without the catheter.

Stylet with code STY23 is used (if necessary) with the 23 cm catheters EMBRYO TRANS / ECHO TRANS / TWINKLE EMBRYO TRANS.

Stylet with code STY18 is used (if necessary) with the 18 catheters EMBRYO TRANS / ECHO TRANS / TWINKLE EMBRYO TRANS.

Stylet with code STY-BU is used (if necessary) with the catheters BULB TRANS ULTRA / BULB TRANS STAR.

Stylet with code STY-PE is used (if necessary) with the catheters PRO EMBRYO TRANS / PRO ECHO TRANS / TWINKLE PRO EMBRYO TRANS.

Soft Obturator with code SB-PE is used (if necessary) with the catheters PRO EMBRYO TRANS / PRO ECHO TRANS / TWINKLE PRO EMBRYO TRANS.

Preparation for use and use directions

1. Before to the procedure, check the condition of the cervix prior to the insertion of the catheter, to ascertain the depth of the uterus and any uterine anteversion or retroversion present in individual patients.
2. Place the patient in the dorsal lithotomy position with her feet in stirrups.
3. Position a vaginal speculum within the vagina so the external cervical os can be visualized during catheter insertion.

Caution:

- * use a stylet in case the cervix bends and it seems to be difficult to pass the cervix.
 - * use a soft obturator if there is large amounts of cervical mucus
4. Cleanse the cervix with an appropriate sterilized saline solution.
 5. Affix the adaptor of the catheter to a syringe preloaded with washed spermatozoa or embryo solution
 6. Guide the catheter through the cervix and into the uterine cavity.
 7. Check the catheter reached to the uterus (if necessary use ultrasound guidance).
 8. Slowly press the plunger of the syringe to introduce the spermatozoa or embryo into the uterine cavity.
 9. Following injection, remove the catheter.
 10. In case of embryo transfer check for retention of embryos in the catheter.
 11. Take the vaginal speculum out and keep the patient at rest for further 40-50 minutes.
 12. Discard the syringe, catheter and other materials used as medical waste.

Note: ink marks or depth positioners can be used to facilitate placement catheter.

* For FLEXI catheter steps 5,6, should be:

5. Guide the outer sheath through the cervix into the uterine. Then affix the adaptor of the inner catheter to a syringe preloaded with washed spermatozoa solution.

6. Guide the loading catheter through the inserted conductor catheter to access the uterine cavity.
Use of a medical device with accessories (if necessary):

If necessary and to give the shape to the conductor catheter insert the stylet into the conductor catheter. Ensure the stylet is firmly placed into the catheter hub via the luer slip fit. Shape the conductor catheter and insert the catheter into the uterus through the cervical canal. After insertion, remove the stylet from the conductor catheter. Insert the loading catheter with the spermatozoa or embryos into the conductor catheter. Introduce the spermatozoa or embryos into the uterine cavity by pressing the syringe plunger.

If necessary, so to prevent cervical mucus from entering the conductor catheter, insert a soft obturator into the guide catheter. Ensure the soft obturator is firmly placed into the catheter hub via the luer slip fit. Insert the catheter with the soft obturator into the uterus through the cervical canal. After insertion, remove the soft obturator from the conductor catheter. Insert the loading catheter with the spermatozoa or embryos into the conductor catheter. Introduce the spermatozoa or embryos into the uterine cavity by pressing the syringe plunger.

Note: ink marks or depth positioners can be used to facilitate placement catheter.

Do not reuse, reprocess & resterilize:

Reuse, Reprocessing & resterilization may compromise the structural integrity, can also create risk of contamination & may not give desired result or create complications, infections which may result in injury, illness or death.

STORAGE CONDITIONS, TRANSPORTATIONS AND SHELF LIFE

The shelf life of the medical device are three years after manufacture when stored at 5°C to 30°C°, protected from (sun) light.

Transport conditions: from -30°C ((±3°C) at Uncontrolled RH) to 60°C (±2°C) 15% RH (±5%) up to 24 hours.

PROCEDURE OF UTILIZATION

When using this product, waste classified as class A and B is formed, which should be disposed of according with the requirements of SanPin 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures. Use, expired and unused product is classified as Class B and disposal should happen according to the requirements of SanPin 2.1.3684-21. Packaging of the product after use of the product belongs to waste of class A and is disposed of with household waste.

WARRANTY LIABILITY

The manufacturer guarantees that the product is compliant with the requirements declared by the manufacturer, if the instructions for use, storage are followed. For questions regarding the quality of the product, you must contact the authorized representative: LTD "Gera", Russia, 620026, Ekaterinburg, Karl Marks str., d. 8 of. 416, e-mail: gera.cat96@mail.

REF Catalogue Number **LOT** Batch Code  Date of Manufacture  Use By  Do not re-use  Do not Re-Sterilize **non-toxic**

 Consult instruction for use  Do not use if Packing is damaged  Caution  Keep out of sunlight  Keep Dry **STERILE EO** Sterilized using ethylene oxide

 **apyrogenic**

 **LATEX** Does not contain natural rubber latex

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