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Hauptstraße 28, 15806 Zossen b. Berlin, Germany

www.botiss.com

Attention! This document presents the technical characteristics that are a trade secret of the company. These parameters are not taken into account and are not controlled by the patient when using the product.

This document is intended for professionals with medical education and can be provided upon request.

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Product details.


Specified max. values	
Material Source	PTFE
Colour	Blue
Size	15x20 mm, 20x30 mm, 30x40 mm
Cleanliness	Free of impurities
Corner radius	Rounded, smooth
Thickness	0,08
Sterility (with EO)	sterile
The surface density of the membrane	0,02 g/cm ²
Tensile strength	20-40 N/mm ²
Elongation	250-500%
Cross-linking	No
Resorbable	No
Conformability	The membrane is flexible and conforms to the defect site
Non-pyrogenic	Yes
Biocompatibility	Yes
Shelf Life	3 years
Packaging	Double pouch package

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
- Outer Blister Label:**
- Product Name and Size
 - Temperature Storage Conditions
 - REF number
 - LOT number
 - Sterile EO
 - Do Not Reuse
 - Use By Date
 - CE Mark
 - Do Not use if packaging is damaged
 - Consult instructions for use
 - Do not re-sterilize
 - Information about the manufacturer
 - Keep away from sunlight

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Symbols:

	Do not reuse
	Do not re-sterilize
	Sterilized using ethylene oxide
	Use-by date
	Manufacturer
	Upper Limit of Temperature
	Keep away from sunlight
	Lot Number
	Reference Number

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	Consult instructions for Use.
	Reject product if packaging is damaged or opened

Information in Russian sticker

Russian sticker for final package will be added locally in Russia after approval with information below (information can be changed in accordance with Russian legislation)

Product name (in Russian)

Article:

Lot:

Manufacturer:

Place of manufacture:

Authorized representative, seller:

Sterile. Shelf life: see packaging. For single use only.

Non-pyrogenic

Non-toxic

Do not use if damaged packaging. Refer to instructions for use

ONLY FOR PROFESSIONAL USE IN STOMATOLOGY

RC # _____ date of issue _____

Label of Outer Chipboard box consist the following information:

- Order number
- Shipment number
- Barcode
- Caution about sterile content
- Cautions symbols

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Composition and package size

permamem® non-resorbable PTFE barrier membrane is offered in the following sizes, with rounded edges: 15 x 20 mm / 20 x 30 mm / 30 x 40 mm. permamem® is sterilized by ethylene oxide and made available in double sterile packaging.

Indications

permamem® is a temporarily implantable membrane for use as a space-creating barrier in the treatment of dental defects.

Product description

permamem® is made of 100% polytetrafluorethylene (PTFE), a biologically inert and tissue-compatible material. Irregularities in the texture of the material arise from the technical process and do not limit the quality or functionality of the product. The non-resorbable PTFE barrier membrane is a temporarily implantable material that prevents the integration and passage of bacteria due to the small pores of the material, thus allowing open healing of the membrane. permamem® is intended to prevent ingrowth of undesired soft tissue cells into the bone defect, thereby creating a suitable environment for bone regeneration.

Areas of application

permamem® is intended for use
- as a space-creating barrier in guided bone regeneration (GBR) and in guided tissue regeneration (GTR).
- for the treatment of dental defects.
- to cover bone defects in periodontal, maxillofacial, oral and implant surgeries.

Contraindications

The non-resorbable membrane is not intended for use under strain. permamem® must not be used in the case of inflammation or active infections.

Instructions for use

permamem® should only be used by doctors who are familiar, through appropriate training, with surgical intervention techniques for dental defects.

- Application

After careful opening of the sterile package, permamem® may be removed from the package using sterile instruments and can be cut to the desired shape and size. In order to prevent perforation of the soft tissue it is important not to leave any sharp edges. The membrane should extend 3-4 mm beyond the edges of the defect to ensure membrane stability and protection of the bone substitute material. A minimum distance of 1 mm to the adjacent teeth should be maintained. If required, the membrane may be fixed by sutures, screws or pins.

- Removal

The membrane should be removed after 3-4 weeks. An exposed membrane may be easily removed with tweezers. Patient comfort can be increased using local anaesthesia. If primary closure is obtained during membrane placement, opening of the surgical site will be required to remove the membrane. After removal of permamem®, the primary healing process and the reepithelialization of the regenerating soft tissue is completed within one month.

Use during pregnancy or breastfeeding period

There are no studies on the use of permamem® during pregnancy or while breastfeeding, or on its effect on reproductive ability. Before using permamem®, the treating physician must therefore make an individual consideration of the possible risks to the child against the benefits for the mother.

Use in children or elderly patients

There is no evidence that special precautions need to be taken based on the age of the patients to be treated.

Warnings / precautionary measures

permamem® is a medical device and must therefore be stored out of reach of children. permamem® is intended for single use only. It must not be resterilized or reused due to risk of infection transmission. The bone regeneration takes approximately 6-12 months, depending on the defect and on the bone substitute material used. During this time period, heavy prosthetic load should be avoided.

Storage instructions

permamem® should be stored at temperatures below 30°C. The product must be protected from sunlight.

Shelf life

The expiration date is printed on the box and on the sterile inner packaging. permamem® must not be used past the specified date.

Sterility

permamem® is sterile as long as the packaging remains unopened and undamaged. If the packaging is damaged, permamem® must not be used. The contents of unused packages that have nevertheless been opened or damaged must not be resterilized and must be disposed of, since the functionality and safety of the product can no longer be guaranteed.

Disposal

The product can be disposed together with normal clinic waste. Local guidelines must be followed.

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