
80

Phillip H Jones – Notary Public
Two Chimneys, North Street
Winkfield Windsor SL4 4TE England
Telephone +44 (0) 7860 398485
Email: info@notarypublicassistance.co.uk

I, PHILLIP JONES, Notary Public duly admitted and sworn, carrying on practice at Two Chimneys North Street Winkfield Windsor SL4 4TE England.

CERTIFY THAT I recognise the signature appearing at the foot of the attached Approval to be that of Rachel Swales whose identity is known to me and who has signed the said document in their capacity as Regulatory Affairs Manager of DePuy (Ireland) and DePuy Ireland UC.

In faith and testimony whereof I, the said Notary, have subscribed my name and set and affixed my seal of office at Windsor aforesaid this Thursday, 16th March 2023.

Notary Public

My Commission expires on death

Protocol no. 13011 /23

Regulated by the Faculty Office of the Archbishop of Canterbury

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
Instruction for use Marathon XLPE Cemented Cup, variations

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Further in the text the following names of the medical device can be used: device, cups, Marathon XLPE, The Marathon XLPE Cemented Cup, Marathon, MD, implant.

The medical devices names can be mentioned in the text in abbreviated form for information perception convenience.

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MD classification

Table 1. Potential risk class of MD:

Health Authority, country	Class	Document
Ministry of Health, Russia	IIB	Order of Ministry of Health of Russian Federation №4n from 6 th June 2012

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INTENDED USE AND PRINCIPLE OF OPERATION

The Marathon XLPE Cemented Cup is intended to replace the bearing surface of the acetabulum as part of a total hip replacement; providing increased patient mobility and reduced pain. The implant is used in conjunction with a femoral component and intended for use with bone cement in an appropriately prepared acetabulum.

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INDICATIONS AND CONTRAINDICATIONS FOR MD USE AND POSSIBLE ADVERSE EFFECTS

Indications

Total hip replacement is indicated in the following conditions:

1. A severely painful and/or arthritic joint (typically due to non-inflammatory degenerative joint disease);
2. Failed previous hip surgery;
3. Acute fracture of the femoral head or neck.

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factors listed in WARNINGS AND PRECAUTIONS;

4. Excessive wear or fracture of the bearing components due to: intraoperative damage to the prosthetic components, loose cement, bone fragments, metallic particles, ceramic particles or other factors listed in WARNINGS AND PRECAUTIONS;
5. Early or late infection;
6. Peripheral neuropathies. Subclinical nerve damage may also occur as a result of surgical trauma;
7. Tissue reactions, osteolysis, and/or implant loosening caused by metallic corrosion, allergic reactions, wear or particulate debris (such as loose cement, metallic, polyethylene or ceramic particles);
8. Residual pain in the hip.

Intraoperative

1. Femoral shaft or acetabular perforation, fissure, or fracture, which may require the use of internal fixation;
2. Trochanteric fracture;
3. Damage to blood vessels (e.g. iliac, obturator and femoral artery);
4. Temporary or permanent nerve damage (e.g. femoral, obturator or isolated peroneal nerve);
5. Subluxation or dislocation of the hip joint due to implant size or configuration selection, positioning of components and/or muscle and fibrous tissue laxity;
6. Lengthening or shortening of the affected extremity;
7. Surgical delay.

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Early Postoperative

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism and myocardial infarction;
2. Hematoma and/or delayed wound healing;
3. Subluxation or dislocation;
4. Disassociation of modular components.

Late Postoperative

1. Trochanteric avulsion from excessive muscular tension, weight-bearing, or inadvertent intraoperative weakening of the trochanter;
2. Aggravation of problems in the ipsilateral or contralateral knee and ankle joints due to leg length discrepancy, femoral medialization and/or muscular deficiencies;
3. Acetabular or femoral fracture due to trauma or excessive loading, particularly in the presence of poor bone stock caused by severe osteoporosis, bone defects from previous surgery, intraoperative reaming procedures, or bone resorption;
4. Bone resorption which may contribute to the deterioration of fixation and eventual loosening of the implant;
5. Periarticular calcification or ossification which may lead to a decrease in joint mobility and range of motion;
6. Traumatic arthrosis of the ipsilateral knee secondary to intraoperative positioning, of the extremity during surgery;
7. Subluxation or dislocation;
8. Disassociation of modular components.

INFORMATION ON POTENTIAL USERS OF MD AND

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DIRECTIONS FOR USE

Information for the patient

Patients who are fitted with hip replacements should be warned that the service life of the implant may depend on their weight and activity level.

Directions for use

The Marathon XLPE Cemented Cup instrumentation is designed to accommodate all open surgical approaches including through the back of the hip (posterior approach), the side of the hip (lateral or anterolateral approach), the front of the hip (anterior approach) or through a combination of approaches. The major steps in a total hip arthroplasty procedure include:

- **Preparation of the Cartilage/Bone:** The proximal head and neck of the femur are removed and the channel inside the femur is prepared for the fitting of the femoral stem. The acetabulum is reamed to achieve the optimum bone surface for fixation of the Marathon XLPE Cemented Cup. The focus of reaming is removal of all sclerotic bone, cartilage and soft tissue from the acetabulum in order to facilitate cement interdigitation.

- **Placement/Positioning of the Acetabular Component:** Cup positioning should optimize fixation, range of motion, dislocation resistance and minimize the likelihood of subluxation, impingement and edge loading. This may be assessed during pre-operative planning, acetabular preparation and cup trialing. The target cup inclination (as measured on radiographs) should be at about 45° taking into account local soft tissue and anatomic landmarks. The target cup anteversion (as measured on radiographs) should be 15 – 20° taking into account local soft tissue and anatomic landmarks. Bone cement is used for fixation of the Marathon XLPE Cemented Cups.

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• **Placement/Positioning of Femoral Components:** The femoral stem is inserted into the prepared femoral bone. The stem can be held in the bone with or without cement. The ball of the ball-and-socket hip joint is then inserted on top of the stem.

• **Positioning and Closure:** With the acetabular cup, femoral stem and ball inserted, the hip replacement joint is placed in final position and the incision closed.

Preoperative

THE SURGEON SHOULD DISCUSS ALL PHYSICAL AND MENTAL LIMITATIONS PARTICULAR TO THE PATIENT AND ALL ASPECTS OF THE SURGERY AND THE PROSTHESES WITH THE PATIENT BEFORE SURGERY. The discussion should include the limitations and possible consequences of joint replacement, the necessity to follow the surgeon's instructions postoperatively and the information provided in the patient leaflet.

The preoperative planning and surgical techniques for implantation of these hip replacement components evolved from the surgical experience gained during the development of many hip prostheses. Surgeons should not begin the clinical use of any hip prosthesis before they have thoroughly familiarized themselves with its specific implantation technique. Surgical technique brochures are available from DePuy.

Intraoperative

It is recommended that components at least one size larger and one size

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smaller than were preoperatively determined be available at surgery to accommodate intraoperative selection of the appropriate size(s).

Protective covers should be left on until the components are ready to be implanted. Do not use components if they have been dropped or have impacted a hard surface. Damage to the component may not be visible but could cause early failure of the prosthesis.

When performing revision surgery, it is imperative that any fibrous membrane encapsulating the failed prosthesis be removed in its entirety. Residual fibrous membranes may interfere with implant fixation. Review any pre-existing devices for damage to mating surfaces that would prevent adequate seating or locking of the new device. Improper selection, placement, or positioning of implants may result in unusual stress conditions and a subsequent reduction in the functional life of the implant. Increasing the anteversion of an acetabular component in total hip replacement may result in instability and/or dislocation. Increasing the anteversion of an acetabular component in total hip replacement may result in instability and/or dislocation. Prior to closure, the surgical site should be thoroughly cleansed of bone chips, extraneous bone cement, ectopic bone, etc. Foreign particles at the metal/plastic or ceramic/plastic interface may cause excessive wear. Range of motion should be thoroughly checked for improper mating, instability, or impingement and corrected as appropriate.

Postoperative

Strict adherence by the patient to the surgeon's instructions and warnings is extremely important. Accepted practices should be followed in postoperative care.

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The patient should be released from the hospital with complete written instructions and warnings regarding exercises and therapies and any limitations on their activities. Component malposition or the use of constrained hip or total hip with modular skirted head or tapered sleeve can decrease joint range of motion and can increase the potential for component wear, impingement premature dislocation or revision. In these cases, the surgeon should inform the patient that high range of motion activities should be avoided.

A continuing periodic patient follow-up is recommended because of the unknown functional lifetime of the implant in the specific patient, particularly with respect to the maintenance of implant fixation and bearing surfaces. Radiological investigations should be performed at each follow-up and compared with previous radiological investigations and correlated with the clinical assessment of the patient. If any radiological changes are observed, such as the occurrence of radiolucencies, bone resorption, or any changes in the position of an implant, these changes should be closely monitored to determine whether they are static or progressive and the patient treated appropriately.

Potential users of MD

Orthopaedic surgeons, operating room staff familiar with the orthopaedic procedures and techniques for Total Hip Arthroplasty before using Marathon XLPE Cemented Cups.

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Temperature	from 21 °C to 25°C
Pressure	86 kPa -106 kPa
Operation conditions	
Humidity	100%
Temperature	From 32 °C to 42 °C
Pressure	86 kPa -106 kPa

DESCRIPTION OF MAIN FUNCTIONAL ELEMENTS OF MD

General description of functional elements of MD

Figure 1. Acetabular Cup

The Acetabular Cup is manufactured from cross-linked Ultra High Molecular Weight Polyethylene (UHMWPE) and includes a flange which extends from the perimeter.


The flange is manufactured from UHMWPE and contains 4 holes which are designed to prevent an air pocket forming behind the flange during insertion of the implant.

The flange is designed to pressurise cement during insertion of the implant and to reduce the likelihood of the implant bottoming out in the

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

acetabulum.

The flange is trimmed to size by the surgeon prior to implantation. Concentric lines at 2mm intervals are laser marked on the top surface of the flange to assist the surgeon during trimming.

The rear surface of the Acetabular Cup features a locating hole designed to help define orientation of the Clip-On Marker Wire, and a set of grooves into which the Clip-On Marker Wire is assembled. The set of grooves are also designed to ensure secure fixation of the implant in the cement mantel.


The Marathon XLPE Cemented Cup design includes two types of marker wire.

Figure 2. Embedded Marker Wire

The first marker wire is a length of Stainless Steel wire pre-embedded in the Acetabular Cup and positioned near the perimeter.

Figure 3. Clip-on Marker Wire

The second marker wire is a Titanium Alloy Clip-On Wire supplied with

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the Acetabular Cup

for assembly by the end user. The Clip-On Marker Wire is designed to allow postoperative radiological assessment of position and wear of the implanted Acetabular Cup.

One end of the Clip-on Marker Wire is intended to mate with a locating hole on the acetabular cup. The position of the locating hole is identified by a mark on the top surface of the implant flange. The wire is then rotated to click-fit into the implant grooves, ensuring correct orientation of the Clip-On Marker Wire. Guidance on positioning of the Clip-On Marker Wire is provided in the information card packaged with the implant.

The Marathon XLPE Cemented Cup is available in a range of sizes and is suitable for both left and right side use. "L" and "R" laser etchings on the Acetabular Cup are designed to ensure correct orientation. When the implant is mounted onto the introducer, only one of these two laser etchings should be left exposed, corresponding to the side the implant is orientated for.

The possibility and ways of integrating with other medical devices

The Marathon XLPE cemented cups are compatible with 22, 26, 28, 32 and 36mm diameter standard metal or ceramic heads. The bearing surface has the same diameter as the mating femoral component (head) with a clearance. The bearing geometry (diameter and surface finish) is transferred directly from the predicate devices (Elite and Charnley) and they are identical. Figure 4 below shows the compatibility chart of bearing diameter Vs outer diameter of the cup.

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		Outer diameter of cup, mm						
		38	40	43	45	47	50	53
Bearing diameter, mm	22	x	x	x	x	x	x	x
	26		x	x	x	x	x	x
	28		x	x	x	x	x	x
	32				x	x	x	x
	36						x	x

Figure 4: Compatibility chart of bearing diameter Vs Outer diameter of the cup

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BASIC STANDARDS TO WHICH THE PRODUCTS MEET

The devices have passed verification and validation tests and meet the requirements of the following standards:

Applicable standard	Description
BS EN ISO 14630	Non active Surgical Implants - General Requirements
BS EN ISO 21534	Non-active surgical implants - Joint replacement implants - Particular requirements
ISO 21535	Non Active Surgical Implants- Joint Replacement Implants- Specific Requirements for Hip Joint Replacement Implants
BS EN ISO 14971	Medical Devices - Application of Risk Management to Medical Devices
ISO 13485	Medical Devices – Quality Management Systems – Requirements for regulatory purposes
BS EN 556-1	Sterilisation of medical devices - requirements

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MD SAFETY AND EFFECTIVENESS REQUIREMENTS Warnings and Precautions

CAUTION:

- Implants and trials components from different manufacturers or implant systems should never be used together;
- Hip prosthesis components should never be reimplanted. Even though the implant appears undamaged, the implant may have developed microscopic imperfections which could lead to failure;
- Always use a trial prosthesis for trial purposes. Trials should not be assembled with any components intended for permanent implantation. Trials must have the same configuration size, etc., as the corresponding components to be permanently implanted;
- Single use devices manufactured by DePuy have not been designed to undergo or withstand any form of alteration, such as disassembly, cleaning or re-sterilization, after a single patient use. Reuse can potentially compromise device performance and patient safety;
- Do not alter or modify implants in any way unless as described in the surgical technique;
- Quality bone stock is critical to ensure fixation of the device. Consider bone grafting and other reinforcement techniques, especially when performing revision surgery, to seat and support the device.

CAUTION: The following conditions, singularly or concurrently, tend to impose severe loading on the affected extremity, thereby placing the patient at higher risk for failure of the hip replacement:

1. Obesity or excessive patient weight;
2. Manual labor;

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3. Moderate to high demand sports participation;
4. High levels of patient activity;
5. Likelihood of falls;
6. Alcohol or drug misuse addiction;
7. Other disabilities, as applicable.

CAUTION: The following conditions singularly or concurrently, tend to adversely affect the fixation of hip replacement implants:

1. Marked osteoporosis or poor bone stock;
2. Metabolic disorders or systemic pharmacological treatments leading to progressive deterioration of solid bone support for the implant (e.g., diabetes mellitus, steroid therapies, immunosuppressive therapies, etc.);
3. History of general or local infections;
4. Severe deformities leading to impaired fixation or improper positioning of the implant;
5. Tumors of the supporting bone structures;
6. Allergic reactions to implant materials (e.g., bone cement, metal);
7. Tissue reactions to implant corrosion or implant wear debris;
8. Disabilities of other joints (i.e., knees and ankles).

NOTE: Diabetes, at present, has not been established as a contraindication. However, because of the increased risk for complications such as infection, slow wound healing, etc., the physician should carefully consider the advisability of hip replacement in the severely diabetic patient.

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