**ACKNOWLEDGED**  
as true, reliable, and correct

Bentley InnoMed GmbH  
Lotzenacker 25  
72379 Hechingen  
Germany

[www.bentley.global](http://www.bentley.global)  
Phone: +49 (0)7471 984 995 10  
Fax: +49 (0)7471 984 995 9  
[info@bentley.global](mailto:info@bentley.global)

**EXPLOITATIVE DOCUMENT**

**(BeGraft Coronary Stent Graft System)**

BeGraft Coronary Stent Graft with the Delivery System in the following designs:  
Stent diameters (mm) – 2.50; 2.75; 3.00; 3.50; 4.00; stent lengths – 8; 12; 16; 18; 21; 24  
Stent diameters (mm) – 4.50; 5.00; stent lengths – 16; 18; 21; 24

Bentley InnoMed GmbH  
Lotzenacker 25  
D-72379 Hechingen

*[Signature]*

2017

2017-04-10

**Index**

|                                                                                            |   |
|--------------------------------------------------------------------------------------------|---|
| 1. Name of the Medical Device .....                                                        | 1 |
| 2. General Description .....                                                               | 1 |
| 3. Indications .....                                                                       | 1 |
| 4. Contraindications .....                                                                 | 1 |
| 5. Potential Adverse Effects and Complications .....                                       | 1 |
| 6. Area of Use .....                                                                       | 1 |
| 7. Specification .....                                                                     | 1 |
| 8. Stent Compliance .....                                                                  | 1 |
| 9. Materials of the Main Functional Components .....                                       | 1 |
| 10. Warnings and Precautions .....                                                         | 1 |
| 11. Description and Operation Principles .....                                             | 1 |
| 12. Application Procedure .....                                                            | 1 |
| 13. Storage Conditions .....                                                               | 1 |
| 14. Conditions of Use (Ambient Temperature, Relative Humidity, Atmospheric Pressure) ..... | 1 |
| 15. Transportation Conditions .....                                                        | 1 |
| 16. Packaging and Package Content .....                                                    | 1 |
| 17. Labeling .....                                                                         | 2 |
| 18. Safe Disposal .....                                                                    | 2 |
| 19. Manufacturer Warranty .....                                                            | 2 |
| 20. Shelf Life .....                                                                       | 2 |
| 21. Manufacturer .....                                                                     | 2 |
| 22. Authorized Representative of the Manufacturer .....                                    | 2 |

**1. Name of**

BeGraft Coronar  
Stent diam  
Stent diam

**2. Gener**

The BeGraft C  
• a ball  
which  
• pre-m  
The delivery  
(2.5-4.0 mm)  
of an implant

**3. Indi**

The BeGraft  
coronary arti  
1. Acut  
2. Anci  
Patients bei  
artery bypas

**4. Cor**

- Patients  
or antipl  
lant / ant
- Patients  
access ar
- Patients
- Patients
- Lesions
- Patients
- Lesions  
% residu
- Lesions  
posed st
- Severe s
- Patients
- Diffuse l
- Patients  
terial PT
- Severe a  
(e.g. asp  
/ antiplat

## 1. Name of the Medical Device

BeGraft Coronary Stent Graft with the Delivery System in the following designs:

Stent diameters (mm) – 2.50; 2.75; 3.00; 3.50; 4.00; stent lengths – 8; 12; 16; 18; 21; 24

Stent diameters (mm) – 4.50; 5.00; stent lengths – 16; 18; 21; 24.

## 2. General Description

The BeGraft Coronary Stent Graft System consists of:

- a balloon-expandable cobalt-chromium (L605) stent, covered with an ePTFE graft material, which is
- pre-mounted on the balloon of a rapid exchange (RX) stent delivery system.

The delivery system is compatible with 0.014" guide wires and is available in usable lengths of 143 cm (2.5-4.0 mm) and 153 cm (4.5/5.0 mm). The stent is made by precision laser cutting of a seamless tubing of an implantable cobalt-chrome alloy (L-605).

## 3. Indications

The BeGraft Coronary Stent Graft with the Delivery is indicated for transluminal implantation in coronary arteries or aorto-coronary bypass grafts for the treatment of:

1. Acute perforation or rupture of coronary arteries.
2. Aneurysm of coronary arteries or coronary bypass-vein graft.

Patients being considered for stent graft implantation should be acceptable candidates for coronary artery bypass graft surgery.

## 4. Contraindications

- Patients with hemorrhagic diathesis or other disorders in which the application of anticoagulant or antiplatelet therapy is contraindicated and patients unable or unwilling to tolerate anticoagulant / antiplatelet therapy (e.g. septic ulceration, cerebrovascular complications).
- Patients with proximal atherosclerosis and/or significant vessel tortuosity in whom guide wire access and adequate guiding catheter support is prohibited.
- Patients with diffuse coronary artery disease resulting in low blood flow in the stented area.
- Patients with ostial lesions.
- Lesions with coronary thrombus formation.
- Patients with unprotected left main coronary vessel disease (i. e. > 50 % stenosis).
- Lesions that cannot be successfully predilated (i. e. < 20 % diameter stenosis reduction, and > 50 % residual diameter stenosis).
- Lesions with threatened or abrupt closure produced during attempted predilation prior to proposed stent implantation.
- Severe stenosis of left main coronary artery.
- Patients with totally obstructed coronary arteries.
- Diffuse lesions.
- Patients with known allergies to the used materials: stent material (L605) and/or to the graft material PTFE and/or delivery system (stainless steel) material.
- Severe allergic reaction to the contrast medium and/or the necessary pharmaceutical treatment (e.g. aspirin) or bleeding complications and patients unable or unwilling to tolerate anticoagulant / antiplatelet therapy.

## 5. Potential Adverse Effects and Complications

Potential undesired effects are possible, but not limited to the following:

- Dissection, perforation, injury, rupture, sudden and total occlusion or other damages to the coronary arteries
- Thrombosis of the coronary arteries
- Acute or subacute closure of stent graft
- Pseudoaneurysm formation
- Unstable angina
- Arterial spasm
- Acute myocardial infarction
- Cardiac dysrhythmia, including ventricular fibrillation
- Hypo / hypertension
- Restenosis of the treated artery
- Distal emboli (air, particle, thromboembolic emboli)
- Hemoptysis
- Hemorrhage or hematoma
- Hemorrhagic stroke
- Hematuria
- Bleeding from anticoagulation/antiplatelet medication
- Side effects caused by the concomitant medication
- Allergic reaction to contrast medium
- Allergic reaction to the implanted stent material (L605) and/or PTFE
- Infections, fever
- Arteriovenous fistula
- Emergency surgery of coronary vessel bypass graft, surgical graft removal
- Deformation of stent symmetry on deployment or on re-crossing with another device (IVU, balloon catheter, etc.)
- Stent graft migration
- Intimal tear
- Tissue reaction, necrosis
- Renal failure
- Death.

## 6. Area of Use

Cardiology.

|              |
|--------------|
| Chai         |
| Total the ti |
| Max.         |
| Min.         |
| Outer dista  |
| Outer prox   |
| Crim         |
| (rou         |
| Nom          |
| Rate         |
| Stent        |
| Stent bar    |
| Max          |
| Max          |
| Radi         |
| Flex         |
| Tens         |
| slee         |
| Tens         |
| Tens         |
| (bon         |
| *For         |
| 105:         |

## 7. Specification

| Characteristics                                                                              | Value                                                                                                                  |                                                      |
|----------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------|
|                                                                                              | For stent sizes<br>Ø2.50; 2.75; 3.00; 3.50;<br>4.00                                                                    | For stent sizes<br>Ø 4.50; 5.00                      |
| Total effective length of the catheter (from the tip to the distal end of the strain relief) | 143 cm                                                                                                                 | 153 cm                                               |
| Max. diameter of the guide wire                                                              | 0.014 " (0.36 mm)                                                                                                      | 0.014 " (0.36 mm)                                    |
| Min. size of the guide catheter                                                              | 5 F (inner diameter min 1.40 mm)                                                                                       | 5 F (inner diameter min 1.40 mm)                     |
| Outer diameter of the shaft – distal region                                                  | 2.7 F (0.88 mm)                                                                                                        | 3.2 F (1.05 mm)                                      |
| Outer diameter of the shaft – proximal region                                                | 1.9 F (0.65 mm)                                                                                                        | 1.9 F (0.65 mm)                                      |
| Crimped stent profile (rounded values)                                                       | 1.1 mm (Ø 2.50 mm)<br>1.1 – 1.2 mm (Ø 2.75 mm)<br>1.2 mm (Ø 3.00 mm)<br>1.2 – 1.3 mm (Ø 3.50 mm)<br>1.3 mm (Ø 4.00 mm) | 1.2 – 1.3 mm (Ø 4.50 mm)<br>1.3 – 1.4 mm (Ø 5.00 mm) |
| Nominal inflation pressure (NP)                                                              | 11 bar / 10 <sup>5</sup> Pa                                                                                            | 10 bar / 10 <sup>5</sup> Pa                          |
| Rated burst pressure (RBP)                                                                   | 16 bar / 10 <sup>5</sup> Pa                                                                                            | 14 bar / 10 <sup>5</sup> Pa                          |
| Stent recoil (all sizes)                                                                     | Max. 7 %                                                                                                               | Max. 7 %                                             |
| Stent foreshortening at nominal pressure (11 bar)                                            | Ø 2.50 mm: < 5 %<br>Ø 2.75 – 4.00 mm: < 10 %                                                                           | Ø 4.50 mm: < 10 %<br>Ø 5.00 mm: < 10 %               |
| Max. inflation time                                                                          | 10 s                                                                                                                   | 10 s                                                 |
| Max. deflation time                                                                          | 15 s                                                                                                                   | 25 s                                                 |
| Radial force at 50% compression                                                              | Ø 2.50 mm: min 0.16N/mm stent length<br>Ø 2.75 – 4.00 mm: min 0.13N/mm stent length                                    | Ø 4.50-5.00 mm: min 0.11N/mm stent length            |
| Flex/kink radius:                                                                            |                                                                                                                        |                                                      |
| Distal shaft with crimped stent                                                              | ≥ 9.0 mm                                                                                                               | ≥ 9.0 mm                                             |
| Tensile strength on tip & distal balloon sleeve                                              | ≥ 1.5 N                                                                                                                | N/A*                                                 |
| Tensile strength on proximal balloon sleeve                                                  | ≥ 5 N                                                                                                                  | ≥ 5 N                                                |
| Tensile strength of the side port (bond)                                                     | ≥ 5 N                                                                                                                  | ≥ 5 N                                                |
| Tensile strength of the Luer-Lock sealing (bond)                                             | ≥ 3N                                                                                                                   | ≥ 3N                                                 |

\*For large diameters (4.5 and 5.0 mm) there is no acceptance criterion. Reason: According ISO 10555-1 Appendix B, B.3.1 distal tips less than 3 mm can be excluded from the test piece

## 8. Stent Compliance

The table below shows the ratio pressure / stent outer diameter. The dark grey areas mark the Nominal Pressure (NP – pressure at which the stent outer diameter reaches its nominal diameter), which equal to:

- 10 [bar, atm or 10<sup>5</sup> Pa] for the stent sizes 4.50; 5.00 mm;
- 11 [bar, atm or 10<sup>5</sup> Pa] for the stent sizes 2.50; 2.75; 3.00; 3.50; 4.00 mm.

The light grey areas mark the Rated Burst Pressure (RBP), which is equal to:

- 14 [bar, atm or 10<sup>5</sup> Pa] for the stent sizes 4.50; 5.00 mm;
- 16 [bar, atm or 10<sup>5</sup> Pa] for the stent sizes 2.50; 2.75; 3.00; 3.50; 4.00 mm.

|                | Stent outer diameter (mm) |        |        |        |        |        |        |
|----------------|---------------------------|--------|--------|--------|--------|--------|--------|
| Pressure (bar) | Ø 2.50                    | Ø 2.75 | Ø 3.00 | Ø 3.50 | Ø 4.00 | Ø 4.50 | Ø 5.00 |
| 7              |                           |        |        |        |        | 4.03   | 4.48   |
| 8              | 2.26                      | 2.49   | 2.68   | 3.05   | 3.52   | 4.20   | 4.68   |
| 9              | 2.34                      | 2.58   | 2.80   | 3.20   | 3.69   | 4.37   | 4.85   |
| 10             | 2.42                      | 2.67   | 2.91   | 3.35   | 3.86   | 4.59   | 5.12   |
| 11             | 2.50                      | 2.75   | 3.00   | 3.50   | 4.00   | 4.61   | 5.12   |
| 12             | 2.57                      | 2.82   | 3.07   | 3.62   | 4.13   | 4.70   | 5.23   |
| 13             | 2.63                      | 2.88   | 3.14   | 3.71   | 4.23   | 4.77   | 5.31   |
| 14             | 2.68                      | 2.93   | 3.20   | 3.80   | 4.31   | 4.83   | 5.37   |
| 15             | 2.73                      | 2.98   | 3.26   | 3.86   | 4.39   | 4.88   | 5.45   |
| 16             | 2.78                      | 3.02   | 3.32   | 3.92   | 4.46   | 4.92   | 5.48   |
| 17             | 2.82                      | 3.06   | 3.37   | 3.97   | 4.53   |        |        |
| 18             | 2.87                      | 3.11   | 3.43   | 4.03   | 4.59   |        |        |

Nominal Pressure

RBP Rated Burst Pressure

e Nomin  
l, which

Ø 5.00  
4.48  
4.68  
4.85  
5.12  
5.23  
5.31  
5.37  
6.45  
5.48

## 9. Materials of the Main Functional Components

- Catheter shaft proximal: PTFE coated hypo-tube
- Catheter shaft distal: Polyamide
- Balloon: Modified Polyamide
- Stent Graft: The stent is made of a medical grade cobalt-chrome alloy (L-605). Outer cobalt-chrome stent surface is covered with microporous ePTFE.

| Component |                                 | Component material                                                                                              | Blood contact |
|-----------|---------------------------------|-----------------------------------------------------------------------------------------------------------------|---------------|
| Catheter  | Stent                           | L605 cobalt/chrome                                                                                              | x             |
|           | Graft tube                      | ePTFE AEOS™                                                                                                     | x             |
|           | Soft tip                        | Modified Polyamide<br>Pebax 6333D (yellow)                                                                      | x             |
|           | Balloon                         | Modified Polyamide (Pebax 7233 SA01)                                                                            | x             |
|           | Balloon markers                 | Pt/Ir (90/10)                                                                                                   | -             |
|           | Distal shaft<br>Inner tube      | Outer layer: Pebax 7233D SA01 (transparent)<br>Inner layer: High Density Polyethylene (blue) PETROTHENE LR34000 | x             |
|           | Distal shaft<br>Outer tube      | Rilsan® BESVO A Nylon II 72D (transparent)                                                                      | x             |
|           | Proximal shaft material         | Stainless steel 304 hypo-tube, PTFE coated (black) Duraskin                                                     | x             |
|           | Strain relief of the hub (lucr) | TPE Thermoplastic K TF6STE (transparent)                                                                        | -             |
|           | Hub (lucr)                      | Polyamide (TROGAMID CX7323; transparent)                                                                        | -             |

## 10. Warnings and Precautions

### General

The medical techniques and procedures described here do not entail all the procedures of PTCA acknowledged in medicine.

- The BeGraft Coronary Stent Graft System should only be used by physicians, who are well trained in interventional cardiology and who are experts in the techniques of coronary angioplasty and coronary stent implantation, as well as in management of possible complications.
- Stent graft implantation in patients who are not acceptable candidates for coronary artery bypass graft surgery require careful consideration, including possible hemodynamic support during stenting.
- The use of different stent types in the same vessel is not recommended. This may increase the risk of subacute thrombosis.
- During stent graft implantation, a cardiac surgery team should be on stand-by.
- Use aseptic techniques in order to avoid contamination of the system.
- Use the BeGraft Coronary Stent Graft System before the expiry date noted on the product label.

- Do not expose the system to organic solvents (e.g. alcohol).
- Do not use if the package is open or damaged.
- Ensure that labeled stent outer diameter and length of the stent graft are suitable for the intended procedure. The stent outer diameter must not exceed the diameter of the coronary artery proximally and distally to the aneurysm, rupture or perforation.
- If predilation of the lesion is required, do not overdilate the lesion site. The BeGraft Corona Stent Graft System is not intended for use as primary dilatation device.
- Ensure that the stent graft is positioned between the RO-markers (don't touch) and undamaged after removal of the stent graft protector.
- Ensure that the shaft is not bent nor kinked. Do not attempt to straighten the shaft if it is bent kinked, as it may break. Instead, prepare a new system.
- Administer appropriate anticoagulation and vasodilation therapy before insertion of the BeGraft Coronary Stent Graft System. Continue the anticoagulation treatment after the implantation according to national guide lines and international standards.
- Only use diluted contrast medium.
- Do not overtighten the hemostatic valve, as lumen constriction may occur, possibly impairing inflation/deflation of the balloon.
- Insert the BeGraft Coronary Stent Graft System into the vessel under fluoroscopic guidance only.
- Never attempt to advance the stent graft system against increased resistance.
- If increased resistance during advancement is detected, stop the procedure and only continue when the reason was found. If this is not possible, then carefully pull back the BeGraft Corona Stent Graft System together with the guiding catheter as a unit. Do not attempt to retract the BeGraft Coronary Stent Graft System through the guiding catheter as this may cause the stent dislodge and an embolization could occur.
- Only use an inflation device with a manometer to inflate the balloon and do not use air or a gaseous medium. Do not exceed the Rated Burst Pressure (RBP) mentioned on the label or a compliance chart.
- A balloon burst during implantation may cause migration or deformation of the stent graft or problems to remove the catheter. In such a case, retract the catheter carefully, leaving the guiding wire in place. If the balloon gets caught in the stent graft, the catheter should be removed immediately by cardiac surgery. A new PTCA catheter should be carefully inserted to reposition and completely dilate the stent graft.
- After stent graft delivery, manipulate the catheter only after the balloon has been fully deflated under vacuum.
- When passing an already implanted stent which is located proximally to the lesion, the BeGraft stent can be deformed or dislodged. In such a case, the correct position and expansion of the stent graft must be ensured.
- Recrossing a partially or fully deployed stent graft with adjunct device must be performed with extreme caution to avoid stent graft damage or deformation.
- The use of laser devices, mechanical atherectomy and hyperthermia are not recommended in the stented region.
- The long-term effectiveness of stent graft implantation has not been fully established.

MRI  
Nor  
over  
It ca  
•  
•  
•  
•  
SAR  
•  
The  
F20:  
this  
urati  
in-vi  
were  
The:  
1.5 i  
leng  
after  
Phili  
cons  
RF-  
in le  
MR-  
ly 9  
MR:  
of Be  
Advi  
The  
tion  
and  
Mec  
Hanc  
Ever  
bend  
Do n

### MRI Safety and Compatibility

Non-clinical testing has demonstrated that the **BeGraft Coronary Stent Graft**, in single and in overlapped configuration up to 130 mm in length, is **MR Conditional** as defined in ASTM F2503. It can be scanned safely immediately after implantation under the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Maximum spatial field gradient of 5,700 Gauss/cm (57T/m) or less
- Maximum force product of 96,000,000 G<sup>2</sup>/cm (96T<sup>2</sup>/m)
- Theoretically estimated maximum whole body averaged specific absorption rate (WBA-SAR) of < 2 W/kg for 15 minutes of scanning
- Normal operating mode of the MR-system.

The BeGraft Coronary Stent Graft shows a magnetically induced **displacement force** (ASTM F2052-14) of 8% (4.5° deflection) and torque (ASTM F2213-06) of 1.0% and should not migrate in this MRI environment. **3.0 Tesla RF-induced heating** in non-clinical testing of overlapped configuration up to 130 mm in length (performed according to ASTM F2182-11) produced a maximum in-vitro temperature rise after 15 minutes of 7.0°C when scaled to a WBA-SAR of 2 W/kg. Tests were conducted in a 3.0 T Siemens Magnetom Trio with "Numaris/4 syngo MR B17" software. These calculations do not take into consideration the cooling effects of the blood flow.

**1.5 Tesla RF-induced heating** in non-clinical testing on overlapped configuration up to 130 mm in length (performed according to ASTM F2182-11) produced a maximum in-vitro temperature rise after 15 minutes of 8.9°C when scaled to a WBA-SAR of 2 W/kg. Tests were conducted in a 1.5 T Philips Intera with a "Release 12.6.1.4 2012-05-22" software. These calculations do not take into consideration the cooling effects of the blood flow.

**RF-induced heating general:** the effects of MRI on overlapped Stent Grafts greater than 130 mm in length or Stent Grafts with fractured struts are unknown.

**MR-artefacts:** In non-clinical testing, the image artefact caused by the device extends approximately 9.15 mm from the Stent Graft when imaged with a gradient echo pulse sequence and a 3 Tesla MR system. Therefore, it may be necessary to optimize the MR imaging parameters in the presence of BeGraft Coronary Stent Grafts.

### Advice due to the single-use

The product BeGraft Coronary Stent Graft System is designed for the single-use only. The reconditioning for a re-use is prohibited. The product is exposed to high mechanical stress during its use and possibly affected thereafter. Therefore only a single-use is recommended for a safe application.

### Mechanical damages risk

Handle the device with special care, avoid its fall, strike, and other mechanical impacts.

Every time carefully inspect the whole system for defects prior to use. Inspect the device for cracks, bends, breaks, holes, or other damages.

Do not use the device if it is damaged or has any defect.

## 11. Description and Operation Principles

The BeGraft coronary stent graft with the delivery system consists of two parts: a balloon-expandable cobalt-chromium (L605) stent, covered with ePTFE, which is pre-mounted on the balloon of a rapid exchange (RX) stent delivery system.

### Stent Graft

The balloon-expandable coronary stent with expandable PTFE graft material used here has been designed for the treatment of acute coronary artery perforation, acute coronary artery rupture, coronary artery aneurysm and coronary bypass-vein-graft aneurysm.

### Delivery System

The PTCA-catheter for the stent graft delivery has a shaft system with an integrated balloon on the distal tip. The PTCA-catheter shaft consists of a single-lumen tube, divided into two lumens in the area of the distal tip. The short distal inner lumen of the catheter permits the use of a guide wire (which exits the system proximally through the side-port) to facilitate advancement of the system and through the area of the perforation / rupture / aneurysm. The balloon is inflated with contrast agent via the continuous lumen. The balloon was designed to provide the exact defined diameter at a specific pressure using a semi-compliant balloon material.

For the x-ray visualization of the stent graft position on the balloon in the vessel two radio-opaque platinum/iridium markers are attached according to the stent length. Two additional markers on the proximal catheter shaft help to position the dilatation catheter in relation to the guiding catheter. The marker close to the Luer lock of the dilatation catheter is used for the femoral application, the other one for the brachial application.

## 12. Application Procedure

**Note:** The physician should consult literature on current medical practice of PTCA.

### Accessories

- Femoral or brachial guiding catheter of appropriate size and configuration to match the coronary artery (5-9 F)
- Introducer sheath, compatible with the guiding catheter
- Hemostatic valve
- Sterile heparinized physiological saline
- 60 % contrast medium diluted 50:50 (per volume) with sterile heparinized physiological saline
- 20 cc Luer-Lock syringe (optional)
- Inflation device with manometer
- 3-way stopcock for a minimum pressure of 25 bar
- 0.014" guide wire (Attention: a hydrophilic coated guide wire can get stuck in the balloon catheter when the coating is swelling.)
- Guide wire introducer
- Torquer for the guide wire

### Concomitant Medication

This information is a recommendation only. It is not intended to dictate medical practice. The conditions under which stent grafts have been used to date include:

Administration of aspirin prior to stent graft implantation.

Administration of heparin and nitro-glycerin during the procedure.

Administration of antiplatelet therapy alone or in combination with other anticoagulants up to some months following stent graft implantation.

Most experience with PTCA procedures has been gained with anticoagulant/antiplatelet therapy up to now. Administration of antiplatelet therapy and anticoagulants are depending on the patients' condition and medical history. Prior anticoagulant/antiplatelet administration is an essential adjunct to stent graft implantation. Proper heparin treatment, before and after sheath removal, is also essential for the overall success of the procedure. Before, during and after stent implantation, specific blood values of the patient (e.g. prothrombin time PT partial thromboplastin time PTT, platelet count etc.) must be determined in order to optimize the medication. The selection, duration and dosage of the medication are under the responsibility of the physician treating the patient, as well as the selection of long-term medication after stent graft implantation. Because of an extended endo-thelialization time it may be mandatory to administer a prolonged anticoagulant/antiplatelet therapy.

#### Prior to use

- Please check the sterile packaging of the stent system and all the equipment necessary for this procedure carefully. Do not use, if the packaging is open or damaged.
- Do not use the stent graft system, if it is damaged or kinked.
- Remember that the balloon diameter of the stent graft system at nominal pressure should not be larger than the artery diameter proximal and distal to the stenosis.
- Check the guide wire size according to the product label.

#### Preparation

1. Carefully inspect the entire product for any damage. Prepare the inflation device according to the manufacturer's instructions.
2. Prepare the sterile physiological saline and the diluted contrast medium.
3. Fill a 10 or 20 cc syringe with 4 cc of diluted contrast medium.
4. Carefully remove the stent graft system from the sterile packaging (aseptic technique!). Pull the stent graft system out of the transport coil. Carefully remove the protector with stiletto from the stent graft.

#### Caution:

Ensure that the stent graft is not damaged (e.g. protruding struts) and that it is positioned between the RO-markers (don't touch the implant). Furthermore, the shaft has to be inspected on bends and kinks. For more details please see "Warnings and Precautions" (section 6). Return any damaged product to Bentley InnoMed.

5. Connect the 3-way stopcock (min. 25 bar) with the luer connector to the balloon catheter.
6. Connect the syringe to the stopcock. Position the syringe vertically with the piston pointing downwards.
7. Create a strong vacuum by pulling back the piston. Maintain the vacuum for 15 - 20 seconds and make sure that no more air bubbles can be observed in the diluted contrast medium.
8. Carefully release the vacuum, close the stopcock and remove the syringe.
9. Repeat this procedure until no more air bubbles can be observed.

#### Warning

If it is not possible to create or to maintain a vacuum, check all connections. If you cannot find the cause, do not use this system under any circumstances and return the unit to Bentley InnoMed.

10. Create another vacuum in the system and close the 3-way stopcock.
11. Remove the air from the tubing and connector of the inflation device.

12. Connect the inflation device to the 3-way stopcock of the balloon catheter.
13. Open the stopcock carefully. The inflation device should be in a vertical position in order to avoid air entering the system.
14. Attach a needle to the syringe; insert it into the distal tip of the catheter and flush the guide wire lumen by applying gentle pressure to the syringe.
15. The BeGraft Coronary Stent Graft System is now ready for use.

#### Pre-dilation

If pre-dilation with a suitable PTCA dilatation catheter prior to stenting of the lesion is required, not overdilate the lesion site.

For the pre-dilatation of the vessel you must be aware of atherosclerotic plaques, which may inhibit advancement of the catheter.

#### Implantation of the stent graft

1. Insert a compatible introducer sheath in the femoral or brachial artery using conventional techniques.
2. Insert the guiding catheter into the vessel through the introducer sheath and fix the hemostatic valve to it.
3. Insert the guide wire into the guiding catheter and carefully push it forwards under fluoroscopic control until it is distally to the lesion/defect.
4. After the guide wire has passed the lesion/defect, close the hemostatic valve tightly to prevent displacement of the wire.
5. Advance the BeGraft Coronary Stent Graft System over the guide wire, ensuring that the guide wire exits through the side-port located approximately 25 cm from the catheter tip.
6. Open the hemostatic valve and carefully push the BeGraft Coronary Stent Graft System forwards over the guide wire to the distal tip of the guiding catheter.
7. Continue the procedure under fluoroscopic control. The radiopaque markers facilitate the positioning of the implant in the area of the perforation, rupture or aneurysm.
8. After positioning, inflate the balloon to Nominal Pressure (NP). Keep this pressure for 10 to 30 seconds. Due to the restricted blood flow, the patient can feel heart pain. Adapt the inflation time to the patient's condition.
9. After deployment of the stent graft, deflate the balloon maintaining negative pressure, allow adequate time for full deflation. After full deflation, open the stopcock to allow atmospheric pressure in the balloon.
10. The inflation can be repeated two to three times, but do not exceed the Rated Burst Pressure (RBP). Check the implant for full expansion as long as the balloon catheter has not been retracted (in order) to minimize the risk that part of the balloon got caught in the stent struts.

#### Caution:

- Do not retract the catheter prior to ensuring that the stent graft is implanted correctly.
- Re-dilation with another balloon catheter is principally possible, but with increased risk.
- Do not overtighten the hemostatic valve, since this may restrict the flow of contrast medium into the balloon and thereby increase the inflation/deflation time.
- Do not overstretch the stent graft and do not exceed the labeled Rated Burst Pressure (RBP) of the balloon, as this may cause damage to the graft material and/or the vessel. Permanently check the correct expansion of the stent graft by monitoring with high resolution fluoroscopic equipment.

- Expansion of the stent graft should not be undertaken if the stent graft has not been properly positioned in the area of perforation, rupture or aneurysm. If the position of the stent graft is not optimal, it should be repositioned before inflation.

#### Note:

- It is necessary that guide wire and/or catheter remain in the treated area until the procedure is completed.
- Different contrast media may affect the inflation/deflation time due to different viscosities.
- Multiple angiographic views and measurements should be taken in order to ensure the correct positioning of the implant.

#### Further dilation in the stent graft

- If the primary angiographic result is not optimal, further expansion of the stent graft with the help of a PTCA catheter is possible.
- If the expanded stent graft needs to be recrossed with a guide wire, the distal tip of a flexible guide wire should be formed to a J-shape prior to reinsertion through the stent graft. This technique is recommended so the guidewire will more smoothly pass through the lumen of the stent graft minimizing the risk of getting caught in the stent struts.

#### Balloon Catheter Withdrawal Procedure

- Open the hemostatic valve.
- Hold the guide wire and the hemostatic valve firmly in one hand, while grasping the balloon catheter shaft with the other hand. Maintain the guide wire's position in the vessel by holding the wire stationary and begin pulling the stent graft catheter with the deflated balloon carefully out of the guiding catheter.
- Retract the balloon catheter until the opening for the guide wire lumen is reached. Then proceed in the same way as with an over-the-wire system until the balloon catheter has been totally removed.
- Close the hemostatic valve.

#### Caution:

- Before retracting the stent graft catheter, make sure that the balloon has been completely deflated.
- Observe the patient regularly during the first 30 minutes after stent placement, and perform an angiographic examination for control. If stent graft placement is associated with the onset of thrombus in the stented region, intracoronary administration of a thrombolytic agent is recommended.

### 13. Storage Conditions

- The devices must be stored in a dry, dark and cool place at temperatures from 10°C to 25°C in original packagings. The packaging must not be damaged and should be wrapped into the stretch film to prevent contamination and dusting.
- Relative humidity at the storage place must not be higher than 50%.
- The product must be kept away from direct sunlight.
- Storage must be in closed dry areas, safe from atmospheric precipitations and soil moisture.
- The packed devices must be stored distant from heaters (not less than 1 m).
- Atmospheric pressure 84.0–106.7 kPa (630–800 mm Hg).

### 14. Conditions of Use (Ambient Temperature, Relative Humidity, Atmospheric Pressure)

Stent graft (implant) conditions of use (in accordance with GOST R 50444-92):

temperature: from +32° C – +42 °C.

relative humidity: 100%.

atmospheric pressure: 84.0–106.7 kPa (630–800 mm Hg). Before use, the packaging has to be checked if it has not been damaged or opened. The device must not be used if the inner packaging is open or damaged.

The BeGraft coronary stent graft with the delivery system is intended for single-use only.

The BeGraft coronary stent graft with the delivery system should only be used by physicians, who are well trained in interventional cardiology and who are experts in the techniques of coronary angioplasty and coronary stent implantation, as well as in management of possible complications.

The BeGraft coronary stent graft with the delivery system is beyond repair and maintenance.

### 15. Transportation Conditions

| Temperature                                                                              | Period of time         |
|------------------------------------------------------------------------------------------|------------------------|
| -25° C – +50° C at a humidity level below max. 85%                                       | Not more than 72 hours |
| Note: the device must be kept at the ambient temperature for at least 8 hours before use |                        |

### 16. Packaging and Package Content


The "BeGraft Coronary Stent Graft System" is packaged in the following manner:

- The inner lumen of the stent graft section is protected by a stiletto whereas the tip and outer surface of the stent graft section are covered by a stent protector tube (PTFE natural).
- The catheter with the stent protector is inserted into a dispenser coil.
- The "coiled product" is placed into a header bag (sterile pouch with Tyvek lid) together with compliance card (stent graft / balloon compliance).
- The header bag (containing the product and the compliance card) is then sealed and labelled.
- The EtO indicator sticker is attached to the label.
- The sealed and labelled header bag is placed into a carton box with one IFU (containing 24 images) and one implant information card (packed in a zip lock bag).
- The pull tab of the carton box is closed with a double-sided adhesive tape and the label is affixed.
- The product is then sterilized using ethylene oxide gas (external sterilisation supplier).
- After final release testing and acceptance at Bentley InnoMed, the product is transferred to distribution storage location (ready for shipment to the specific customer).


#### Product box content

|                                                          |          |
|----------------------------------------------------------|----------|
| 1. BeGraft Coronary Stent Graft System (the product)     | 1 piece. |
| 2. Stent Protector (PTFE tubing) with Stiletto (mandrel) | 1 piece. |
| 3. Dispenser Coil                                        | 1 piece. |
| 4. Compliance card                                       | 1 piece. |
| 5. Header bag                                            | 1 piece. |
| 6. Instruction for use (IFU)                             | 1 piece. |
| 7. Implant Information Card                              | 1 piece. |
| 8. Zip Lock Bag                                          | 1 piece. |

1. 2. BeGraft Coronary Stent Graft System (the Product) hub (luer) with strain relief

3. BeGraft Coronary Stent Graft System in Dispenser Coil

4. Compliance Card

**Bentley BeGraft Compliance Chart**


Subtotal from Outer Diameter (mm)

| Pressure (mmHg) | Ø 2.50 | Ø 2.75 | Ø 3.00 | Ø 3.50 | Ø 4.00 | Ø 4.50 | Ø 5.00 |
|-----------------|--------|--------|--------|--------|--------|--------|--------|
| 7               |        |        |        |        |        | 4.00   | 4.48   |
| 8               | 2.78   | 2.68   | 2.88   | 3.05   | 3.53   | 4.20   | 4.68   |
| 9               | 2.78   | 2.58   | 2.80   | 2.70   | 3.00   | 4.37   | 4.85   |
| 10              | 2.62   | 2.63   | 2.81   | 2.88   | 3.86   |        |        |
| 11              |        |        |        |        |        | 4.81   | 5.12   |
| 12              | 2.67   | 2.67   | 3.07   | 3.82   | 4.13   | 4.75   | 5.18   |
| 13              | 2.83   | 2.88   | 3.16   | 3.71   | 4.73   | 4.77   | 5.31   |
| 14              | 2.68   | 2.93   | 3.20   | 3.80   | 4.31   | 4.83   | 5.37   |
| 15              | 2.73   | 2.98   | 3.26   | 3.85   | 4.39   |        |        |
| 16              | 2.78   | 3.03   | 3.32   | 3.93   | 4.48   |        |        |

☐ Nominal Pressure
 ☐ Burst Pressure

© 2016 Bentley Industries Limited All Rights Reserved 430345-007-000

## 5. Header Bag

## 6. Instruction for Use

Fig. 1. The front page of the IfU. Dimensions: (200x200)±2 mm, 80 pages.

## 7. Im An In and U / Lot label: to m show

Fig. 2

## 7. Implant Information Card

An Implant Information Card is supplied with this product. This card must be carefully filled out and the peel-off label of the respective implant, indicating the Manufacturer / Product Name / REF# / Lot#, should be applied to the indicated area on the Implant Information Card (up to 4 peel-off labels can be applied). The patient should carry this Implant Information Card with him and show it to medical personnel in case of further treatments. The view of the Implant Information Card is shown at the next page.


Fig. 2. Implant Information Card: a – front view; b – back view. Dimensions: (217x90)±2 mm.

## 17. Labeling

The device package shows the following required information:

- company name;
- brand name of the device;
- type of the device (generic name), its size and reference number;
- name and registered office address of the manufacturer;
- country of the manufacturer;
- symbols, identifying the information necessary for the proper usage of the the medical device;
- sketch of the device with main dimensional information for the catheter (including recommended guide wire size) and symbols for implant size as well as the information that the catheter type is "rapid exchange";
- quantity of the devices in the package;
- batch code / LOT number;
- use before date
- date (month, year) of manufacturing;
- prohibition against re-use;
- do not use if package is damaged;
- keep away from sunlight;
- keep away from moisture;
- single use only;
- barcode;
- certification information (number of notified body for CE-marking CE0086);
- information of customers' support service;
- number and date of GHOST-R Registration Certificate (provided by a sticker which is put to the product box by Bentley InnoMeds distributor for Russia);
- a sticker: "Mean water permeability of cover material in expanded state: 2.5 ml/min/cm<sup>2</sup>".

The table below shows the label information used for all product sizes

| Feature                                                          | Label Information                                        |
|------------------------------------------------------------------|----------------------------------------------------------|
| Guide catheter max diameter                                      | 0.014" (0,36 mm)                                         |
| Guide catheter min diameter                                      | 5 F                                                      |
| Sterilization method                                             | Sterilized by ethylene oxide                             |
| Instructions for use (IFU)                                       | IFU is put into the box                                  |
| Compliance table (delivered with the product as compliance card) | Ref. stent compliance table (printed on compliance card) |
| Storage time (shelf life)                                        | 3 years                                                  |

The table below shows the label information that depends on the product size.

| Feature                                                                                              | Label Information                  |                                    |
|------------------------------------------------------------------------------------------------------|------------------------------------|------------------------------------|
|                                                                                                      | For stent diameters 2.50 – 4.00 mm | For stent diameters 4.50 – 5.00 mm |
| Catheter effective length (from the catheter tip to the distal end of the transparent strain relief) | 143 cm                             | 153 cm                             |
| Shaft outer diameter, distal                                                                         | 2.7 F (0.88 mm)                    | 3.2 F (1.05 mm)                    |
| Shaft outer diameter, proximal                                                                       | 1.9 F (0.65 mm)                    | 1.9 F (0.65 mm)                    |
| Nominal stent pressure In Vitro                                                                      | 11 atm                             | 10 atm                             |
| Rated burst pressure (RDP)                                                                           | 16 atm                             | 14 atm                             |

Label  
Each  
umen


Symb

### Labels for the patient documentation

Each of the device packages (carton box and sterile bag) has 3 additional labels for the patient documentation with the respective information:


- Reference number
- Date of Manufacture: "respective date"
- Lot number
- Use By: "date"
- Product name "BeGraft Coronary Stent Graft System"
- Company Logo
- Bar code including GTIN, expiry date, and LOT number

### Symbols

- Caution, consult accompanying documents
- Manufacturer/manufactured by
- Date of Manufacture
- Number of Units
- Storage: Store in a dry, dark, and cool place
- Use by (date) (Use Before)
- Do not Use if Package is damaged
- Do not Re-use
- Do not Re-sterilize

- Sterilized by ethylene oxide
- Rapid Exchange (RX)
- Catalogue (Reference) number (REF)
- Batch code/LOT number
- Recommended Guide Wire
- Guide Catheter Min Size
- Pressure (P)
- Nominal Pressure (NP)
- Rated Burst Pressure (RBP)

Expanded Stent Outer Diameter

Expanded Stent Length

Only sterile and non-pyrogenic in unopened packages

MR Conditional

### 18. Safe Disposal

Medical devices must be disposed after usage, after the expiry date is reached, or if the package is damaged and the product cannot be used.

Medical devices must be disposed in accordance with the actual Russian Federation Sanitary Regulations and Standards for medical wastes management, practical guidance for medical products: infection, preconditioning and sterilization procedures.

### 19. Manufacturer Warranty

Bentley InnoMed GmbH guarantees that this product has been manufactured, sterilized and packed with the due care prescribed in the EU Directive for Medical Devices 93/42/EEC. Each product is individually tested and inspected before it is released for packaging. Bentley InnoMed GmbH will repair any device defects which can be proven to have originated during manufacture or packaging are immediately brought to the attention of Bentley InnoMed GmbH or its distributors. This warranty applies in lieu of and ruling out any other verbal or written, express, implicit or other warranty.

As a result of biological differences in individuals, no product is 100% effective under all circumstances. Due to this fact and since we have no control over the conditions, the instruments used in operation and the diagnosis, Bentley InnoMed GmbH and its distributors do not guarantee effectiveness in application; unsatisfactory results also cannot be ruled out. Bentley InnoMed GmbH and its distributors shall not be liable for any, damage or losses of any nature arising directly or indirectly from the use of the device. The product liability expires if any damage caused by transport, storage or improper handling by the operator, or other factors outside the control of manufacturer occurs. Bentley InnoMed GmbH and its distributors shall not be liable for damage arising from re-sterilization or reuse of the product.

### 20. Shelf Life

Shelf life of the medical device – 3 years from the manufacturing date.

### 21. Manufacturer

Bentley InnoMed GmbH, Lotzenacker 25, D - 72379 Hechingen, Germany

[www.bentley.global](http://www.bentley.global)

phone +49 (0)7471 984 995 10,

fax +49 (0)7471 984 995 9.

Email [info@bentley.global](mailto:info@bentley.global)

### 22. Authorized Representative of the Manufacturer

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