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//

Instructions for use of a medical device
Soft tissue filling gel based on hyaluronic acid HYAL SYSTEM ACP

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I, the undersigned **Dr. NICOLAS MARCOLINA Public Notary in Padua**, registered with the Padua District Board of Notaries (*Collegio Notarile del Distretto di Padova*), CERTIFY that:

- Mr. CARLO PIZZOCARO born in Milano (Italy) on 9th August 1969, who I attest is the Chairman of the Board of Directors and Legal Representative of "**FIDIA FARMACEUTICI S.p.A.**" (a Company organised and existing under the laws of Italy), having its registered office in Abano Terme (Padua - Italy), Via Ponte della Fabbrica n. 3/A, registered with the Companies' Registry of Padua with registration number, and fiscal code, 00204260285;-----

whose personal identity, title and signatory power, executed before me in my presence and marked at the foot of the page, the foregoing document.-----

Padua, on the 9th day of February 2024.

A handwritten signature in dark ink, appearing to read 'N. Marcolina', is written over a circular notary seal. The seal is partially visible on the right side of the page, showing some text around its perimeter.

Soft tissue filling gel based on hyaluronic acid HYAL SYSTEM ACP

Instructions for use

Name and composition

Soft tissue filling gel based on hyaluronic acid HYAL SYSTEM ACP consisting of:

1. Pre-filled syringe containing 1 ml of auto-cross-link hyaluronic acid gel, 2.0% - 1 pc.
2. Sterile needle 30G, 13 mm - 1 pc.
3. Sterile needle 30G, 4 mm - 1 pc.
4. Instructions for use - 1 pc.
5. Adhesive sticker indicating name of the product, lot number and expiry date - 2 pcs.

Composition

Each pre-filled syringe of HYAL SYSTEM ACP contains:

Auto-cross-linked hyaluronic acid (ACP) 20 mg

Sodium chloride 9 mg

Water for injections up to 1,1 ml

Intended use

Soft tissue filling gel based on hyaluronic acid HYAL SYSTEM ACP, injected into mid to deep dermis, is intended for a medium lasting correction of minor soft tissue defects such as facial lines of glabellar, periocular, nasolabial, perioral, chin and cheek areas, neck and décolletage areas and in the management of hypertrophic scars (pathological or not) resulting from burns, surgical wounds and skin trauma due to volume replenishment, viscoelastic properties, dense structure, the natural moisturizing function of a hyaluronic acid obtained using auto-cross-link technology (ACP200).

Product description and principle and mechanism of its operation

HYAL SYSTEM ACP is a sterile biodegradable viscoelastic medical device for intradermal injection based on low molecular weight (200 kDa) hyaluronic acid of enzymatic origin (bacterium type: Streptococcus Equi) obtained by internal binding technology ("auto-cross-link", ACP200). HYAL SYSTEM ACP is supplied in single-use pre-filled syringe containing 1 ml of gel.

The crosslinking of the Auto-crosslinked polymer (ACP) occurs between HA molecules only and achieved through a base-catalysed reaction consisting of the activation of HA carboxyl groups by 2-chloro-1-methylpyridinium iodide (CMPI) and subsequent nucleophilic acyl substitution by the hydroxyl groups of HA in organic solvent. CMPI gives a "zero-length" crosslinking, in detail promotes the formation of intra- and intermolecular ester bonds between the carboxyl and the hydroxyl groups of the HA chains and avoided any foreign bond incorporation into the crosslinked matrix of HA (the crosslinking occurs between HA molecules, thus, no residual content of crosslinking agents).

Mechanism of ACP action:

- hyaluronic acid interacts with CD44 receptors involved in the processes of adhesion, proliferation and differentiation of fibroblasts and mesenchymal stem cells (MSC) contained in adipose tissue;
- hydrolysis of ACP attracts macrophages, which secrete growth factors that stimulate the growth of fibroblasts and the formation of a new extracellular matrix;
- hyaluronic acid released from ACP, combining with CD44 receptors, inhibits the formation of collagenase MMP-13 under the action of cytokine IL-1 β .

The 200 kDa hyaluronic acid fraction is released slowly (the bioresorption period is up to 9 days) and, therefore, has a prolonged effect, which allows you to improve and maintain the elasticity and elasticity of the skin for a longer time.

Main physico-chemical and biological characteristics

The medical device corresponds to the parameters specified in the table:

п/п	Parameter	Specifications
1.	Appearance	Pre-filled syringe equipped with plunger rod and backstop. The syringe is filled with a colourless, odourless gel, free from foreign particles.
2.	The origin of hyaluronic acid	Bacterial fermentation (type of bacteria: Streptococcus Equi)
3.	Reticulation ratio	4,5±1,5% m/m
4.	Average molecular weight of the hyaluronic acid fraction as sodium salt	200 kDa
5.	Hyaluronic acid content	16–20 mg (2,0 %)
7.	Osmolality	250–335 mOsm/kg
10.	Dynamic viscosity	10000–25000 мПа·с
11.	pH	5.5-6.7
12.	Extractable volume	1.0–1.2 ml/unit
13.	Sterility	Sterile
14.	Bacterial endotoxins	<0.2 EU/ml
15.	Mass of the filled syringe	5,92 g ± 5%
16.	Luer Lock adapter pull force	Not less than 43 N
17.	Luer Lock adapter torque	10 Ncm or more
18.	Leakage at the syringe tip/needle hub	There must be no leakage.
19.	Cylinder	The piston stop should slide smoothly without stopping.
20.	Needle	30 G, 13 or 4 (+1/-2) mm

Indications

HYAL SYSTEM ACP is indicated for a medium lasting correction of minor soft tissue defects such as facial lines of glabellar, periocular, nasolabial, perioral, chin and cheek areas, neck and décolletage areas and in the management of hypertrophic scars (pathological or not) resulting from burns, surgical wounds and skin trauma.

Application field

Application field: dermatovenereology, cosmetology, plastic surgery.

Potential users

The product is intended for use by qualified medical personnel – doctors.

Injection area

Mid to deep dermis of such as facial lines of glabellar, periocular, nasolabial, perioral, chin and cheek areas, neck and décolletage areas and hypertrophic scars (pathological or not) resulting from burns, surgical wounds and skin trauma.

Dosage and administration

HYAL SYSTEM ACP is intended only for mid to deep dermis injection.

The techniques are reported in the paragraphs below according to the defect to be treated.

The amount injected will depend on the areas to be treated based on the experience and the injection technique of the healthcare professional.

a. Treatment of soft tissue defects

- Inject intradermally the gel, superficially, at mid-level or deeply, according to the aimed results. The mode of administration, site of administration, and principal administration technique (linear or fan) are at the discretion of the healthcare professional. Different injection techniques could be used at healthcare professional's discretion and in different facial areas.

It is generally recommended a linear technique or multiple injections in the middle or deep

dermis. For wide areas of the face, a linear crossed technique is suggested. Fine lines can be treated with superficial injections, using linear and/or serial injection techniques.

- Bring the product to room temperature
- Upon injecting the product, gently massage the treated area to spread it evenly into the tissue
- It is advisable to adopt individual-tailored treatment protocols to fit individual requirements.

The amount of HYAL SYSTEM ACP usually administered in facial wrinkles as nasolabial, glabellar, perioral, chin, cheek and décolletage area range from 0.1 ml up to a maximum of 2 ml.

Clinical data shows the effectiveness of the device, when a treatment cycle consisting of three sessions of HYAL SYSTEM ACP administered two or three weeks apart is performed. It is possible to repeat the treatment cycle if required after 4-6 months.

b. Treatment of hypertrophic scars

- Infiltrate the scar with the product in the intra lesion dermal planes
- Infiltrate the same amount of product in the subdermal level
- Repeat the procedure until complete scar infiltration, using a maximum of 2 ml of product 25 cm² of scar area, up to obtain a transitory ischemia indicated by the transition to a white colour of the lesion

HYAL SYSTEM ACP is usually administered in one session of treatment of hypertrophic scar using a maximum of 2 ml of product for 25 cm² of scar area. Clinical data shows the effectiveness of the device, when a treatment cycle consisting of two sessions of HYAL SYSTEM ACP, administered at 15 days interval is performed.

The expected duration of clinical effect is up to 4-6 months. HYAL SYSTEM ACP is 100% biodegradable. HYAL SYSTEM ACP biodegradation period depends on the individual characteristics of the patient. Degradation products are completely eliminated from human tissues by a natural metabolic pathway.

Directions for use

1. Bring the product to room temperature.
 - **Treatment of soft tissue defects**
2. Remove the protective tip cap of the syringe.
3. Insert one of the needles provided in the box. It is recommended to institute personalized treatments to fulfill patient's needs. It is suggested to use the 4 mm needle to treat areas such as the frontal, peri-ocular, peribuccal and the mandibular supero-inferior line, while it is preferable to use the 13 mm needle to treat other areas such as the zygomatic antero-lateral region and the whole subzygomatic area.
4. Inject intradermally the gel, superficially, at mid level or deeply, according to the aimed results. Different injection techniques could be used at physician's discretion and in different facial areas. It is generally recommended a linear technique or multiple injections in the middle or deep dermis. For wide areas of the face, a linear crossed technique is suggested. Fine lines can be treated with superficial injections, using linear and/or micropomphi injection techniques.
5. After injection, gently massage the treated area to improve product distribution in the tissue.

-Treatment of hypertrophic scars

2. Remove the protective tip cap of the syringe.
3. Insert the 13 mm needles provided in the box.
4. Infiltrate the scar with the product in the intra lesion dermal planes.
5. Infiltrate the same amount of product in the subdermal level.
6. Repeat the procedure until complete scar infiltration, using a maximum of two ml of product for 25 cm² of scar area, up to obtain a transitory ischemia indicated by the transition to a white color of the lesion.

Interactions

HYAL SYSTEM ACP can be used with topical anaesthetics containing lidocaine.

HYAL SYSTEM ACP must not be used concomitantly with disinfectants containing quaternary

Contraindications

HYAL SYSTEM ACP is contraindicated in individuals with known hypersensitivity to any component of this product.

HYAL SYSTEM ACP must not be implanted into blood vessels. Implantation of HYAL SYSTEM ACP into dermal vessels may cause vascular occlusion or embolic phenomena.

HYAL SYSTEM ACP must not be used:

- in anatomical sites with infectious or inflammatory processes
- in individuals with severe allergies manifested by a history of anaphylaxis, or history or presence of multiple severe allergies.

Use in specific populations

The safety and efficacy of HYAL SYSTEM ACP in pregnant women, lactating women who are breast feeding, or in subjects less than 18 years of age have not been established and therefore its use is not recommended in these populations.

Warnings and Precautions

- Healthcare professionals are encouraged to discuss all potential risks of injection with their individuals prior to treatment and ensure that individuals are aware of signs and symptoms of potential complications.
- In order to minimize the risks of potential complications, this product should only be used by healthcare professionals who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Administer in aseptic conditions. As any other percutaneous treatment, application of HYAL SYSTEM ACP may imply a risk of infection.
- A post-administration monitoring time is recommended in order to identify any potential undesirable side effect.
- Use with caution in subjects treated with anticoagulant and/ or antiplatelet agents due to the increased risk of ecchymoses or haematomas in the inoculation site. It is to pay close attention to areas with telangiectasia.
- Do not over correct the skin defects.
- Until any initial swelling and redness have gone away exposure of the treated area to excessive sun, UV lamp exposure, and extreme cold weather should be minimized.
- Sunbathing, sunbeds and intense heat (such as saunas and Turkish baths) should be avoided for 1 day after treatment.
- The product is supplied in blister packed pre-filled single-use syringes. The solution contained in the prefilled syringe is sterile. The blister is a protective packaging and it is not meant to provide a sterile barrier system. The needles provided in packaging are sterile and single use. Do not resterilise after use.
- The product is for single use, that means it is intended to be used once only for a single subject. The assembled syringe must be discarded immediately after use, regardless of whether or not the solution has been completely administered. When multiple injections are performed to the subject, the needle could be damaged; it is advisable to substitute it with a new needle.
- If this product is reprocessed and/or reused, Fidia Farmaceutici cannot guarantee performance, functionality, material structure, or cleanliness or sterility of the product. Reuse could lead to illness, infection and/or serious injury to the individual or user.
- After use, dispose according to applicable national practice.
- Do not use if the package and, in particular, the syringe and needle blisters are damaged.
- Do not utilize the product after the expiration date reported on the packaging
- Keep out of reach of children.
- No clinical data are available on possible interactions with other aesthetic correction treatments (implant of biomaterials, lifting, botulinum toxin injections, laser, chemical peeling) or permanent fillers.
- No clinical data are available on the use of HYAL SYSTEM ACP in patients affected by autoimmune conditions (LE, Scleroderma, etc.), neoplasia or with ongoing immunosuppressor

therapy; depending on the nature of the disease, the healthcare professional shall therefore choose the administration and the associate treatment on an individual case-by-case basis.

Undesirable effects

The treated area may show reactions commonly related to the injection, such as redness, minimal hematomas, oedema, ecchymosis erythema and mild pain, that spontaneously resolve from few hours to few days after application.

Conditions of use, storage, and transportation

Store and transport at temperature between + 2 ° C and + 8 ° C and 60% relative humidity. The short excursions of temperature, such the one that could occur during transportation or storage, do not affect the physico-chemical profile of the device.

Do not freeze.

When applied, the product is resistant to temperatures from 32 to 42°C, exposure to biological fluids and secretions of body tissues with which it comes into contact.

Delivery set

"Soft tissue filling gel based on hyaluronic acid HYAL SYSTEM ACP" is delivered in a cardboard package, containing a blister with 1 (one) pre-filled disposable syringe containing 1 ml of a gel and 2 (two) sterile needles, 30G, 13 mm and 30G, 4 mm, individually packed in plastic capsules, instructions for use and 2 (two) adhesive sticker indicating name of the product, lot number and expiry date.

Sterilization of the medical device

"Soft tissue filling gel based on hyaluronic acid HYAL SYSTEM ACP" is for single use only.

Pre-filled syringes are sterilized with moist heat. Supplied needles are sterilized by radiation method.

Disposal

After use, the product should be disposed of in a dedicated container. It is forbidden to dispose of together with household waste.

The used product should be disposed of as potentially hazardous (contaminated) medical waste (class B). An unused or expired product can be disposed of as normal (uncontaminated) medical waste (class A). Refer to national guidelines and sanitary standards (SanPiN 2.1.3684-21) regarding medical waste treatment.

Manufacturer:

FIDIA Farmaceutici S.p.A.

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35031 Abano Terme (Padova) Italy

Tel.: +39 049 8232111

E-mail: info@fidiapharma.it

Warranty and complaints

"FIDIA PHARMACEUTICI S.P.A." assumes a warranty obligation to the end user in respect of the medical device "Soft tissue filling gel based on hyaluronic acid HYAL SYSTEM ACP" for the absence of defects and compliance with the characteristics declared by the manufacturer within 36 months from the date of manufacture, subject to the conditions specified in the instructions by application.

In the event of undesirable events contact the authorized representative of the manufacturer on the territory of the Russian Federation.


Authorized representative of the manufacturer on the territory of the Russian Federation

Limited liability company "FIDIA FARMA RUSSIA" ("FIDIA RUSSIA" LLC)

123610, Russia, Moscow, intracity territory of municipal district Presnensky, Embankment Krasnopresnenskaya, house 12, room 1/18.

Tel.: +7 (495) 249-54-95, E-mail: claims@fidiapharma.ru

Symbols used on labels

	Consult instructions for use
	Use by
	Batch code
	Sterilized using steam
	Sterilized using irradiation
	Do not use if package is damaged
	Temperature limitation
	Do not reuse
	Name and Address of the manufacturer
	Single sterile barrier system
	Caution

The list of the international documents / standards applied to the medical device

Directive 93/42/EEC, ISO 9001:2015, EN ISO 13485:2016/A11:2021, EN ISO 14971:2019/A11:2021, ISO/TR 24971:2020, EN ISO 14630:2012, Reg 1907/2006, Commission directive EC 2003/94/EC (2003), European Pharmacopoeia (Ph. Eur), current edition, ISO 14644-1:2015, EN 556-1: 2001 / AC:2006, EN ISO 17665-1:2006, ISO/TS 17665-2:2009, EN ISO 11737-1:2018/amd1:2021, ISO 11138-1:2017, EN ISO 11138-3:2017, EN ISO 14161:2009, EN 285:2015 + A1 2021, Ph Eur 2.6.14, Ph Eur 2.6.1, EN ISO 80369-7:2021, ISO 80369-20:2015, ISO 11040-4:2015/Amd1:2020, ISO 11040-5:2012, ISO 11040-8:2016, ISO 11607-1:2020/A11:2022, EN ISO 11607-2:2020/A11:2022, CPMP/QWP/122/02, rev.1 (2007), ICH Q1A(R2), EN ISO 10993-1:2020, EN ISO 10993-2:2022, EN ISO 10993-3:2014, ISO 10993-5:2009, ISO 10993-6:2016, EN ISO 10993-9:2021, EN ISO 10993-10:2023, EN ISO 10993-11:2018, EN ISO 10993-12:2021, EN ISO 10993-16:2017, EN ISO 10993-17:2023, EN ISO 10993-18:2020/Amd1:2022, ISO 10993-23:2021, ICH Q3D, ICH Q3C, EN ISO 14155:2011, MEDDEV 2.7/1 rev.4, EN ISO 15223-1:2016, EN 1041:2008 +A1:2013, MEDDEV 2.7/3 rev.3, MEDDEV 2.12/2 rev.2, IEC 62366: 2015+ Amd1: 2020, IEC/TR 62366-2:2016.

The list of the national standards applied to the medical device

GOST 31214–2016, GOST R 52770–2016, GOST ISO 10993-1-2021, GOST R ISO 10993-2-2009, GOST ISO 10993-5-2011, GOST ISO 10993-6-2021, GOST ISO 10993-10-2011, GOST ISO 10993-11-2021, GPA.1.2.4.0005.15, GOST R 52770–2023, MG 1856–78, MR 1436–76, MGC 4.1.3171–14, MGC 4.1.3166–14, GOST R IEC 62366-1-2021, GOST ISO 11607-1-2018, GOST R ISO 15223-1-2023, GOST R ISO 14630–2017.

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