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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of ATLAS LINK (BEIJING) TECHNOLOGY CO., LTD on the annexed DOCUMENT is genuine.

授权签字
Authorized
Signature: Lyu Cuilian

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Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

OPERATIONAL DOCUMENTATION IN VITRO DIAGNOSTIC MEDICAL DEVICE

pH-Balance FRAUTEST vaginal acidity test

Manufactured by: Atlas Link (Beijing) Technology Co., Ltd. Room 811 Zeyang Plaza, No. 166 Fushi Road, Shijingshan District, 100043 Beijing, PR of China

Signature _____

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

Operational doc

Possible reagents
 * Bacterial vaginosis
 * Mixed infection
 * Sexually transmitted infections

CONTENTS

1. Name of medical device
2. Manufacturer of medical device
3. Authorized representative of manufacturer
4. Information about the intended use of medical device
 - 4.1. Intended use of medical device
 - 4.2. Description of target analyte with indication of the qualitative, semi-quantitative or quantitative type
 - 4.3. Specific disorder, condition or risk factor for the detection, definition or differentiation of which the medical device is designed
 - 4.4. Type of test portion
5. Information about the intended purpose for clinical laboratory diagnostics
6. Functional characteristics and intended use of the medical device
7. Description of the testing procedure using the medical device, description of reagents, calibrators and control material
8. Risks of the use of medical device, contraindications, expected and predictable side effects associated with the use of medical device for its intended purpose
9. The list of materials and special materials that are required for testing (analysis) but are not included in the scope of delivery of medical device
10. Technical characteristics of medical device
11. Description of accessories, medical devices or devices that are not medical, but are intended for use in combination with medical device
12. Information for the identification of medical devices in order to obtain a safe combination and (or) information about the known restrictions on the joint application of medical devices for their intended purpose
13. Information about the presence of medicinal product for human use, the materials of animal and (or) human origin in the medical device
14. Information about the stability characteristics of medical device
 - 14.1. Transportation conditions
 - 14.2. Shelf life and storage conditions
15. Information about the sterile condition of medical device, the method of its sterilization and the sequence of activities in case of the sterile package violation
16. Information about the procedure for processing of medical device for its reuse, including the cleaning, disinfection, packaging and, if necessary, re-sterilization
17. Information for consumers of medical device about the precautions to be taken if necessary, as well as about the precautions and/or measures taken in relation to the potentially infectious material contained in the medical device
18. Information about the intended purpose of medical device only for single application in accordance with its intended use
19. Information about the conditions necessary for the collection, processing and preparation of samples, data on the stability of the analyzed samples, including the conditions and duration of storage, transportation conditions, restrictions with respect to the freezing (thawing) cycles

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

20. Information about the preparation for intended use, and the operation of medical device
21. Description of testing procedure, including the calculations and interpretation of test results and, if necessary, information about the feasibility of confirmatory tests
22. Analytical performance characteristics
 - 22.1 Sensitivity
 - 22.2. The effect of interfering substances
 - 22.3. Repeatability and reproducibility study
23. Clinical efficacy characteristics
24. Biological reference interval for the use of medical device
25. Information about the limitations related to sample that may affect study results
26. Warnings and/or special precautions for the safe disposal of medical device and its accessories

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

The result should be compared with the test from the vagina. Match the color of the solution.

1. Name of medical device

pH-Balance FRAUTEST vaginal acidity test.

2. Manufacturer of medical device.

Atlas Link (Beijing) Technology Co., Ltd., Room 811 Zeyang Plaza, No. 166 Fushi Road, Shijingshan District, 100043 Beijing, PR of China.

3. Authorized representative of manufacturer.

Limited Liability Company BOLEAR. Address: Russia, 15201, Moscow, Kashirsky proezd, 13, room XIV. Tel./fax +7 (495) 787-92-68, e-mail: info@bolear.ru, www.bolear.ru.

4. Information about the intended use of medical device.

4.1. Intended use of medical device


The test is intended for the quantitative determination of the vaginal fluid acidity (pH of the vaginal medium).

4.2. Description of target analyte with indication of the qualitative, semi-quantitative or quantitative type.

The acidic vaginal pH of 3.8 to 4.5 is a basic requirement for the optimal functioning of the body's own system to protect the vagina. This system effectively avoids contamination by pathogens and the occurrence of vaginal infections. The most important and natural defense against the vaginal problems is a healthy vaginal microflora. There are many factors that regulate the composition and dynamics of the vaginal microbial ecosystem and support the natural protective functions of the vagina. The lactobacilli are the one of the important factors of this.

Under the influence of lactobacilli, the glycogen secreted by the upper wall of the vagina decomposes with the formation of glucose sugar, and then the lactic acid. This results in the acidic environment in the vagina with pH value usually between 3.8 and 4.5. The acidic environment in the vagina forms an effective protection of the vaginal mucosa against the pathogenic microorganisms.

Low pH level is optimal for the growth of lactobacilli, and high pH level is optimal for the growth of pathogenic bacteria

A deviations of pH to the alkaline side (more than 4.5) indicates an imbalance in the vaginal microflora and may be one of the indicators of the risk of vaginal infection development (colpitis/vaginitis) or the most common disease – bacterial vaginosis.


Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

Possible reasons for vaginal pH changes are:

- * Bacterial vaginosis;
- * Mixed infections;
- * Sexually transmitted diseases;
- * Premature rupture of the foetal membranes;
- * Estrogen deficiency and antibiotic treatment;
- * Intestinal dysbiosis;
- * Some types of contraception: intrauterine devices, products containing antiseptics;
- * Menopause;
- * Improper or excessive intimate grooming.

If the pH is normal, but the characteristic symptoms are present (itching, burning, discharge of white color and curdled consistency), these manifestations may be the symptoms of candidal vaginitis (milk thrush).

Vaginal pH at different stages of a woman's life

4.3. Specific disorder, condition or risk factor for the detection, definition or differentiation of which the medical device is designed.

- Perform the test no earlier than 12 hours after sexual activity.
- Perform the test at least 12 hours after using vaginal medicines (vaginal suppositories, creams, gels, etc.).
- Perform the test only 3-4 days after the end of the critical days.
- Perform the test 15 minutes after urinating because leftover urine can lead to false test results.
- Do not wash or shower immediately prior the testing.
- Use each test only once.

4.4. Type of test portion.

The test is intended solely for the analysis of vaginal secretion.

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

15. Information
sterilization and
The test is man

5. Information about the intended purpose for clinical laboratory diagnostics.

The test is a disposable medical device for self-diagnosis in the medical institutions and at home.

Target use profile:

- 1) type of user: woman, no special education;
- 2) education: special education is not required;
- 3) knowledge: reading and understanding of Europeanized Arabic numerals;
- 4) linguistic understanding: reading the instructions in Russian;
- 5) experience: no experience is required.

6. Functional characteristics and intended use of the medical device.

The test is intended for the quantitative determination of the vaginal fluid acidity (pH of the vaginal medium). The pH-Balance FRAUTEST vaginal acidity test is an auxiliary device for determining the pH value, and not the presence of any infection.

7. Description of the testing procedure using the medical device, description of reagents, calibrators and control materials.


Test application description.

Preparation for administration.


1. Wash your hands before testing.
2. Open the packaging and remove the test from the foil bag just before testing. Touch the test only in the area of the handle. Avoid contact of the pH indicator zone with liquids.

Administration.

1. Grasp the round handle with your thumb and forefinger as shown in Fig. 1, while the window with the pH indicator on the applicator should be on the side of the index finger. If you use the test this way, then the execution will be correct.

2. Sit with your knees apart or one foot on a toilet seat or other firm surface.
3. Gently part the labia and gently insert the applicator into the vagina with your free hand, just as you would use a tampon (Fig. 2), to a depth of about 2.5 cm, but no deeper than the area of the round handle. The insertion of the applicator should be completely painless.


4. Press the inserted part of the applicator against the back of the vagina with light pressure (the wall that faces the rectum) and hold for 10 seconds so that the pH indicator is well moistened with vaginal secretion. Adequate wetting with secretion is very important for the evaluation of the test.

5. Gently remove the test from the vagina and evaluate the result.

Evaluation of the test results.

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

The result should be assessed immediately after removing the test from the vagina, while the measurement area is still wet. Match the color in the test area with the color scale on the foil package and select the most similar area (Figure 3). The darker the color in the test area, the higher the pH.

After evaluating the result, you can dispose the used test along with your normal household waste.

The test results can be entered on a special plate for further monitoring of the pH level.

Testing	1	2	3	4	5
pH level					

- ✓ If the pH indicator does not turn into a uniform color after testing, this indicates that the test was not performed properly. A possible reason is that the pH indicator has not sufficiently contacted with the vaginal secretion. It is recommended to repeat the testing with a new test.
- ✓ If the pH level is out of the normal range once, then repeat the measurement the next day.
- ✓ If the pH level is again outside the normal range, it is recommended to see a doctor to determine the cause of the increased pH value. Remember that many factors can temporarily change the pH level of the vagina without causing any disease.

If the following symptoms occur and your vaginal pH is too high, you should see your doctor:

- Bleeding, including spotting (not during menstruation);
- Suspected discharge of amniotic fluid;
- Burning sensation when urinating;
- Itching or burning in the vagina;
- Increased discharge and the appearance of a "fishy", unpleasant odor.

8. Risks of the use of medical device, contraindications, expected and predictable side effects associated with the use of medical device for its intended purpose.

- Perform the test no earlier than 12 hours after sexual activity.
- Perform the test at least 12 hours after using vaginal medicines (vaginal suppositories, creams, gels, etc.).
- Perform the test only 3-4 days after the end of the critical days.
- Perform the test 15 minutes after urinating because leftover urine can lead to false test results.
- Do not wash or shower immediately prior the testing.
- Use each test only once.
- The test only detects changes in pH levels, not the presence of any infection.
- An acidic pH value is not 100% protection against infections. If you notice painful symptoms despite the normal pH, see your doctor.
- Do not use the test after the expiration date (see the date on the package).
- Do not start any treatment before discussing the test result with your doctor.
- If the test is not used properly, it can lead to rupture of the hymen in women who are not yet sexually active.

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

22. Analytical
22.1. Sensitivity
Objective.
Verify the

9. The list of materials and special materials that are required for testing (analysis) are not included in the scope of delivery of medical device.

Not required.

10. Technical characteristics of medical device.

The pH-Balance FRAUTEST vaginal acidity test is a vaginal applicator device (Fig. 4) with a pH measuring strip integrated into it. The vaginal device consists of a round handle and an insertion tube. On the one hand, at the tip of the insertion tube, there is a window in which the indicator zone of the pH strip is located (the zone for measuring the pH level).

The round handle allows you to touch safely the vaginal devices. The vaginal device is inserted about 2.5 cm (1 inch) into the vagina and the pH measurement zone is gently pressed against the back of the vagina. This ensures that the pH measurement zone is in contact with the vaginal secretion. The vaginal device is then removed from the vagina and the pH is read.


Figure 4. pH-Balance FRAUTEST vaginal acidity test.

11. Description of accessories, medical devices or devices that are not medical, but are intended for use in combination with medical device.

Not required.

12. Information for the identification of medical devices in order to obtain a safe combination and (or) information about the known restrictions on the joint application of medical devices for their intended purpose.

The test is not intended to be used in combination with other medical devices.

13. Information about the presence of medicinal product for human use, the materials of animal and (or) human origin in the medical device.

The test does not contain the medicinal products, and materials of animal and/or human origin.

14. Information about the stability characteristics of medical device.

14.1. Transportation conditions.

The transportation is carried out by any type of transport according to the requirements of transportation applicable to this type of transport, subject to temperature conditions from +2°C to +30°C in a dry place, protected from direct sunlight and freezing.

14.2. Shelf life and storage conditions.

Shelf life is 24 months from the date of manufacture, subject to the observance of the storage and transportation conditions. The test should not be used after the expiration date.

Store in sealed containers at temperatures between + 2 ° C and + 30 ° C in a dry place protected from direct sunlight. Do not freeze!

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

15. Information about the sterile condition of medical device, the method of its sterilization and the sequence of activities in case of the sterile package violation.

The test is manufactured in the non-sterile form.

16. Information about the procedure for processing of medical device for its reuse, including the cleaning, disinfection, packaging and, if necessary, re-sterilization.

The test is a disposable medical device and disposed of after use.

17. Information for consumers of medical device about the precautions to be taken if necessary, as well as about the precautions and/or measures taken in relation to the potentially infectious material contained in the medical device.

- Perform the test no earlier than 12 hours after sexual activity.
- Perform the test at least 12 hours after using vaginal medicines (vaginal suppositories, creams, gels, etc.).
- Perform the test only 3-4 days after the end of the critical days.
- Perform the test 15 minutes after urinating because leftover urine can lead to false test results.
- Do not wash or shower immediately prior the testing.
- Use each test only once.
- The test only detects changes in pH levels, not the presence of any infection.
- An acidic pH value is not 100% protection against infections. If you notice painful symptoms despite the normal pH, see your doctor.
- Do not use the test after the expiration date (see the date on the package).
- Do not start any treatment before discussing the test result with your doctor.
- If the test is not used properly, it can lead to rupture of the hymen in women who are not yet sexually active.

18. Information about the intended purpose of medical device only for single application in accordance with its intended use.

The test is intended for single use.

19. Information about the conditions necessary for the collection, processing and preparation of samples, data on the stability of the analyzed samples, including the conditions and duration of storage, transportation conditions, restrictions with respect to the freezing (thawing) cycles.

The test is intended solely for the analysis of vaginal secretion.

20. Information about the preparation for intended use, and the operation of medical device.

1. Wash your hands before testing.
2. Open the packaging and remove the test from the foil bag just before testing. Touch the test only in the area of the handle. Avoid contact of the pH indicator zone with liquids.

21. Description of testing procedure, including the calculations and interpretation of test results and, if necessary, information about the feasibility of confirmatory tests.


Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

Operational c

Table 2.


Administration.

1. Grasp the round handle with your thumb and forefinger as shown in Fig. 1, while the window with the pH indicator on the applicator should be on the side of the index finger. If you use the test this way, then the execution will be correct.

2. Sit with your knees apart or one foot on a toilet seat or other firm surface.

3. Gently part the labia and gently insert the applicator into the vagina with your free hand, just as you would use a tampon (Fig. 2), to a depth of about 2.5 cm, but no deeper than the area of the round handle. The insertion of the applicator should be completely painless.


4. Press the inserted part of the applicator against the back of the vagina with light pressure (the wall that faces the rectum) and hold for 10 seconds so that the pH indicator is well moistened with vaginal secretion. Adequate wetting with secretion is very important for the evaluation of the test.

5. Gently remove the test from the vagina and evaluate the result.

Evaluation of the test results.

The result should be assessed immediately after removing the test from the vagina, while the measurement area is still wet. Match the color in the test area with the color scale on the foil package and select the most similar area (Figure 3). The darker the color in the test area, the higher the pH.

After evaluating the result, you can dispose the used test along with your normal household waste.

The test results can be entered on a special plate for further monitoring of the pH level.

Testing	1	2	3	4	5
pH level					

✓ If the pH indicator does not turn into a uniform color after testing, this indicates that the test was not performed properly. A possible reason is that the pH indicator has not sufficiently contacted with the vaginal secretion. It is recommended to repeat the testing with a new test.

✓ If the pH level is out of the normal range once, then repeat the measurement the next day.

✓ If the pH level is again outside the normal range, it is recommended to see a doctor to determine the cause of the increased pH value. Remember that many factors can temporarily change the pH level of the vagina without causing any disease.

If the following symptoms occur and your vaginal pH is too high, you should see your doctor:

- Bleeding, including spotting (not during menstruation);
- Suspected discharge of amniotic fluid;
- Burning sensation when urinating;
- Itching or burning in the vagina;
- Increased discharge and the appearance of a "fishy", unpleasant odor.

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

22. Analytical performance characteristics.

22.1. Sensitivity.

Objective.

Verify the ability of the test to provide the same sensitivity as BV-pH precision panel (pH 4.0 or pH 4.4).

Regulatory documents and compliance with their requirements.

Guidelines for the use of in vitro diagnostic medical device prepared by the Food and Drug Administration (FDA).

Guidelines prepared by the State Food and Drug Administration

This study complies with all the applicable laws and regulations.

Materials.

- Precision panel BV-pH, pH levels 4.0, 4.4, 4.7, 5.5, 7.0.

Note: The accuracy was determined by adding a solution prepared from citric acid. All the samples were confirmed with a pH meter. 10 replications for each sample in accordance with the package insert were tested. The results are provided below.

- pH-Balance FRAUTEST vaginal acidity test.

Lot number: 20180926, 20181101, 20181209.

- Bacterial vaginosis detection card (Amsel method) manufactured by Shenzhen SkiariBio-Tech Co., Ltd.

Lot number: 20180629, expiration date: May 8, 2020.

Study design.

Check the samples for pH 4.0, 4.4, 4.7, 5.5, 7.0 with 3 different lots, repeat the tests 30 times for each sample. The recorded results for the sensitivity test should be the same as for the standard panel and commercial test.

Evaluation criteria.

Negative (-)

If the convex pH measurement zone on the applicator turns pistachio green upon exposure to the sample. The result is negative.

Positive (+)

If the convex pH measurement zone on the applicator turns yellow, bluish green or dark green upon exposure to the sample. The result is positive.

Analytical sensitivity criteria:

pH level	color	result
pH 4.0	Yellow	-(30/30)
pH 4.4	Yellow-green	-(30/30)
pH 4.7	Dark yellow-green	+(30/30)
pH 5.5	Blue-green	+(30/30)
pH 7.0	Blue	+(30/30)

Results.

The results of testing are provided in Tables 1-2.

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

pH 4.7
pH 5.5
pH 7.0

Table 1. pH-Balance FRAUTEST vaginal acidity test

pH level	lot: 20180926		lot: 20181101		lot: 20181209	
	color	result	color	result	color	result
pH 4.0	Yellow	-(30/30)	Yellow	-(30/30)	Yellow	-(30/30)
pH 4.4	Yellow-green	-(30/30)	Yellow-green	-(30/30)	Yellow-green	-(30/30)
pH 4.7	Dark yellow-green	+(30/30)	Dark yellow-green	+(30/30)	Dark yellow-green	+(30/30)
pH 5.5	Blue-green	+(30/30)	Blue-green	+(30/30)	Blue-green	+(30/30)
pH 7.0	Blue	+(30/30)	Blue	+(30/30)	Blue	+(30/30)

Conclusion.

The pH-Balance FRAUTEST vaginal acidity test is able to identify all the positive samples and show the same sensitivity as the standard panel including the BV-pH precision panel (pH 4.0 or pH 4.4) and the commercial test. No false positive results were observed.

22.2. The effect of interfering substances.

Regulatory documents and compliance with their requirements.

Guidelines for the use of in vitro diagnostic medical device prepared by the Food and Drug Administration (FDA).

Guidelines prepared by the State Food and Drug Administration

This study complies with all the applicable laws and regulations.

Materials.

- pH-Balance FRAUTEST vaginal acidity test

Lot number: 20180926, 20181101.

- Samples with pH 4.4, 4.7 and 5.5.
- Human protein Mucin-1/MUC-1.

Study design.

The test device was evaluated using the samples from two different lots. The test devices were analyzed using the samples with pH 4.4, 4.7 and 5.5 and human protein Mucin-1/MUC-1 1 µg/ml. Testing was repeated five times for each sample. The tests were carried out according to the instruction leaflet. The results of studies are provided in Table 2.

Evaluation criteria.

Negative (-)

If the convex area of the applicator is yellowish-green when measuring the pH of the sample.

Positive (+)

The result is positive if the convex area of the applicator when measuring the pH of the sample is yellow, bluish-green or darker in color.

Results.

The results of testing are provided in Table 2.

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

Table 2.

	Sample	Lot	
		20201028	20201029
Control group	pH: 4.4	5/5 pH 4.4	5/5 pH 4.4
	pH: 4.7	5/5 pH 4.7	5/5 pH 4.7
	pH: 5.5	5/5 pH 5.5	5/5 pH 5.5
Human protein Mucin-1/MUC-1 1 µg/ml	pH: 4.4	5/5 pH 4.4	5/5 pH 4.4
	pH: 4.7	5/5 pH 4.7	5/5 pH 4.7
	pH: 5.5	5/5 pH 5.5	5/5 pH 5.5

Conclusion.

The interfering substances do not affect the test results.

22.3. Repeatability and reproducibility studies.

Objective.

To confirm that the repeatability and reproducibility of the test meet the specified requirements.

Regulatory documents and compliance with their requirements.

Guidelines of the Food and Drug Administration (FDA) for in vitro diagnostic medical devices.

Guidelines of the National Medical Products Administration (NMPA).

This study complied with all the applicable laws and regulations.

Materials.

- Precision panel for determination of pH in bacterial vaginosis (BV-pH); the pH value is 4.0, 4.4, 4.7, 5.5, 7.0.

Note: The accuracy was determined using the solutions with a known amount of substance prepared using the citric acid and disodium phosphate. All the values of the samples were confirmed with a pH meter. 10 replications for each sample in accordance with the package insert were tested. The results of studies are provided in Tables 4-1, 4-2, 4-3, 4-5, 4-6, 4-7, 4-8, 4-9

- pH-Balance FRAUTEST vaginal acidity test, lot number: 20180926, 20181101, 20181209.

Design of studies.

Testing of three different lots of tests with a BV-pH precision panel by three different operators in three different experimental locations. For each panel number, repeat 10 tests for each lot. 3 operators for reproducibility, 1 operator on 5 different days for repeatability. The test conditions for repeatability and reproducibility are illustrated in Tables A-1, A-2 and A-3.

Table A-1. Operator 1 at location 1.

Lots	Lot 1	Lot 2	Lot 3
Days			
Day 1	10 replications of each panel number	10 replications of each panel number	10 replications of each panel number
Day 2	10 replications of each panel number	10 replications of each panel number	10 replications of each panel number
Day 3	10 replications of each panel number	10 replications of each panel number	10 replications of each panel number
Day 4	10 replications of each panel number	10 replications of each panel number	10 replications of each panel number
Day 5	10 replications of each panel number	10 replications of each panel number	10 replications of each panel number

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

pH 7.0
pH 4.0
pH 4.4

Table A-2. Operator 2 at location 2.

Lots	Lot 1	Lot 2	Lot 3
Days			
Day 1	10 replications of each panel number	10 replications of each panel number	10 replications of each panel number

Table A-3. Operator 3 at location 3.

Lots	Lot 1	Lot 2	Lot 3
Days			
Day 1	10 replications of each panel number	10 replications of each panel number	10 replications of each panel number

Evaluation criteria.

Negative (-)

If the convex area of the applicator is yellowish-green when measuring the pH of the sample.

Positive (+)

The result is positive if the convex area of the applicator when measuring the pH of the sample is yellow, bluish-green or darker in color.

Study criteria:

pH level	color	result
pH 4.0	Yellow	-(30/30)
pH 4.4	Yellow-green	-(30/30)
pH 4.7	Dark yellow-green	+(30/30)
pH 5.5	Blue-green	+(30/30)
pH 7.0	Blue	+(30/30)

When testing the same sample, the results should not differ significantly for three different lots, performed by different operators at different experimental locations. The accuracy is equal to or greater than 95%.

Table 4-1. Test results. Operator 1, location 1. Lot: 20180926

		Test 1	Test 2	Test 3	Test 4	Test 5
Day 1	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
Day 2		Test 1	Test 2	Test 3	Test 4	Test 5
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green

Name	pH-Balance FRAUTEST
Number	OD-RU-B-pH
Version	1

	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 1	Test 2	Test 3	Test 4	Test 5
Day 3	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
Day 4	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 1	Test 2	Test 3	Test 4	Test 5
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
Day 5	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 1	Test 2	Test 3	Test 4	Test 5

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-pH
	Version	1

Operational doc

pH 4.0
pH 4.4
pH 4.7
pH 5.5

Table 4-2. Test results. Operator 1, location 1. Lot: 20181101

		Test 1	Test 2	Test 3	Test 4	Test 5
Day 1	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
Day 2		Test 1	Test 2	Test 3	Test 4	Test 5
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
Day 3		Test 1	Test 2	Test 3	Test 4	Test 5
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
Day 4		Test 1	Test 2	Test 3	Test 4	Test 5
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green

	pH 7.0	Blue	Blue	Blue	Blue	Blue
		<i>Test 1</i>	<i>Test 2</i>	<i>Test 3</i>	<i>Test 4</i>	<i>Test 5</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
Day 5		<i>Test 6</i>	<i>Test 7</i>	<i>Test 8</i>	<i>Test 9</i>	<i>Test 10</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue

Table 4-3. Test results. Operator 1, location 1. Lot: 20181209

		<i>Test 1</i>	<i>Test 2</i>	<i>Test 3</i>	<i>Test 4</i>	<i>Test 5</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
Day 1		<i>Test 6</i>	<i>Test 7</i>	<i>Test 8</i>	<i>Test 9</i>	<i>Test 10</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		<i>Test 1</i>	<i>Test 2</i>	<i>Test 3</i>	<i>Test 4</i>	<i>Test 5</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
Day 2		<i>Test 6</i>	<i>Test 7</i>	<i>Test 8</i>	<i>Test 9</i>	<i>Test 10</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		<i>Test 1</i>	<i>Test 2</i>	<i>Test 3</i>	<i>Test 4</i>	<i>Test 5</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
Day 3		<i>Test 6</i>	<i>Test 7</i>	<i>Test 8</i>	<i>Test 9</i>	<i>Test 10</i>
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-ph
	Version	1

Operational do
- All the test
the manufacture
- Visual inte
performed by

	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
Day 4		Test 1	Test 2	Test 3	Test 4	Test 5
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
Day 5	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 1	Test 2	Test 3	Test 4	Test 5
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue
		Test 6	Test 7	Test 8	Test 9	Test 10
	pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
	pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
	pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
	pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
	pH 7.0	Blue	Blue	Blue	Blue	Blue

Table 4-4. Test results. Operator 2, location 2. Lot: 20180926

	Test 1	Test 2	Test 3	Test 4	Test 5
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue
	Test 6	Test 7	Test 8	Test 9	Test 10
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue

Table 4-5. Test results. Operator 2, location 2. Lot: 20181101

	Test 1	Test 2	Test 3	Test 4	Test 5
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue
	Test 6	Test 7	Test 8	Test 9	Test 10

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue

Table 4-6. Test results. Operator 2, location 2. Lot: 20181209

	Test 1	Test 2	Test 3	Test 4	Test 5
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue
	Test 6	Test 7	Test 8	Test 9	Test 10
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue

Table 4-7. Test results. Operator 3, location 3. Lot: 20180926

	Test 1	Test 2	Test 3	Test 4	Test 5
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue
	Test 6	Test 7	Test 8	Test 9	Test 10
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue

Table 4-8. Test results. Operator 3, location 3. Lot: 20181101

	Test 1	Test 2	Test 3	Test 4	Test 5
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue
	Test 6	Test 7	Test 8	Test 9	Test 10
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRA
	Number	OD-RU-B-ph
	Version	1

26. Warning accessories.
The used tests
The unused tests

Table 4-9. Test results. Operator 3, location 3. Lot: 20181209

	Test 1	Test 2	Test 3	Test 4	Test 5
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue
	Test 6	Test 7	Test 8	Test 9	Test 10
pH 4.0	Yellow	Yellow	Yellow	Yellow	Yellow
pH 4.4	Yellow-green	Yellow-green	Yellow-green	Yellow-green	Yellow-green
pH 4.7	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green	Dark yellow-green
pH 5.5	Blue-green	Blue-green	Blue-green	Blue-green	Blue-green
pH 7.0	Blue	Blue	Blue	Blue	Blue

Conclusion.

The tests showed the correct results for all the control samples. No differences in the results of different tests for each lot, different lots of samples, different operators and different test locations were observed. This indicates that both the repeatability and reproducibility among different lots and operators have reached 100% which shows that the test meets the specified requirements.

23. Clinical efficacy characteristics.

Objective.

To confirm the diagnostic sensitivity and specificity of the test.

Regulatory documents and compliance with their requirements.

Guidelines for the use of in vitro diagnostic medical device prepared by the Food and Drug Administration (FDA).

Guidelines prepared by the State Food and Drug Administration

This study complies with all the applicable laws and regulations.

Materials.

- The potential samples were collected at the Third Clinical Hospital of Hebei Medical University. pH-Balance FRAUTEST vaginal acidity test, lot No. 20181209, date of manufacture: December 9, 2018

- Bacterial vaginosis determination card (Amsel method) manufactured by Shenzhen Shieraybayo-Tech Co., Ltd.

Lot No.: 20180629, expiration date: May 8, 2020

Study design.

I. Evaluation of samples collected from patients at the Third Clinical Hospital of Hebei Medical University.

- Subjects were female patients who were hospitalized due to the discovery of a potential gynecological infection;

- In total, about 105 samples were collected.

II. Experts and clinical laboratories.

Third Clinical Hospital of Hebei Medical University.

III. Requirements to samples.

All the samples were confirmed.

The samples had to be selected at random.

IV. Testing performance.

Operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

- All the testing was performed by clinical specialists in each clinical laboratory according to the manufacturer's instructions using the approved samples.

- Visual interpretation of the results of the pH-Balance FRAUTEST vaginal acidity test was performed by an independent clinical specialist.

- The diagnostic center was obliged to sum up the results and send them to Atlas Link

Evaluation criteria.

The pH of normal vaginal discharge is between 3.8 and 4.5.

Results.

To carry out the correlation, 105 samples were taken. 51 positive samples and 54 negative samples were confirmed using the bacterial vaginosis determination card (Amsel method) manufactured by Shenzhen Shieraybayo-Tech Co., Ltd. These samples were randomly selected and tested using the device for determining the vaginal acidity manufactured by ATLAS LINK. The test results are provided in Table 5.

Table 5. Test results.

Results of test performed using the comparison method								
		PH: 4.0	PH: 4.4	PH: 4.7	PH: 5.0	PH: 5.3	PH: 5.5	PH: 7.0
Results of test performed using the test	PH: 4.0	17	2					
	PH: 4.4	1	33	1				
	PH: 4.7		1	12	2			
	PH: 5.0			1	13			
	PH: 5.3					9		
	PH: 5.5						4	
	PH: 7.0							3
Total		18	36	14	15	9	6	3
Matching of the same results		90.5%						
Matching within +/- 1 percent		100%						
Matching of positive results		98.0%						
Matching of negative results		98.1%						

Results.

The clinical evaluation was performed to confirm the clinical efficacy of the instrument for the use of pH-Balance FRAUTEST vaginal acidity test. In total, 105 samples were tested as the control reagents using the bacterial vaginosis determination card (Amsel method) manufactured by Shenzhen Shieraybayo-Tek Co., Ltd. the diagnostic sensitivity of which is 98.0%; diagnostic specificity is 98.1%; overall compliance is 98.05%.

24. Biological reference interval for the use of medical device.

The biological reference interval is pH 4.0-7.0

25. Information about the limitations related to sample that may affect study results.

- Perform the test no earlier than 12 hours after sexual activity.
- Perform the test at least 12 hours after using vaginal medicines (vaginal suppositories, creams, gels, etc.).
- Perform the test only 3-4 days after the end of the critical days.
- Perform the test 15 minutes after urinating because leftover urine can lead to false test results.
- Do not wash or shower immediately prior the testing.

Balance FRA
OD-RU-B-pH
1

operational documentation for in vitro diagnostic medical device	Name	pH-Balance FRAUTEST
	Number	OD-RU-B-pH
	Version	1

26. Warnings and/or special precautions for the safe disposal of medical device and its accessories.

The used tests should be disposed as household waste.
The unused test due to the expiration date is disposed as household waste.

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