The State of Maryland

Office of the Secretary of State

This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille does not certify the content of the document for which it is issued.

Apostille

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1. Country: United States of America
This public document
2. has been signed by **Scott A. Poyer**
3. acting in the capacity of **Clerk of the Circuit Court**
4. bears the seal/stamp of the **Clerk of Circuit Court of Anne Arundel County**

Certified

5. at Annapolis, Maryland
6. the **6th day of December, 2018**
7. by The Secretary of State of Maryland

8. No. 470347

9. Seal

10. Signature

Secretary of State

State of Maryland, Anne Arundel County, Sct.

I, SCOTT A. POYER, Clerk of the Circuit Court for Anne Arundel County, Maryland, a court of record, do hereby certify that **Linda Davis** was a commissioned/appointed and qualified Notary Public commencing on the 1st day of May, 2017.

In Testimony Whereof, I have hereunto set my hand and affixed the seal of the court December 6, 2018

Scott A. Poyer

Scott A. Poyer, Clerk
Circuit Court for Anne Arundel County, Maryland

Signed

Subscribed and sworn to before me this

6 day December, 20 18.

[Signature]
Notary Public

LINDA DAVIS
Notary Public-Maryland
Anne Arundel County
My Commission Expires
May 14, 2021

Affidavit

State of Maryland)
) SS:
County of Anne Arundel)

Before me, the undersigned notary public, this day, personally,
appeared to me known or proven, who being duly sworn according to law,
deposes the following:

Yoshiko Yano

I, _____, hereby certify that the attached
document is

- ☒ an original document
☐ a true copy of an original document

to the best of my knowledge and belief.

Signed _____

Subscribed and sworn to before me this

6 day December, 20 18.

Notary Public

LINDA DAVIS
Notary Public-Maryland
Anne Arundel County
My Commission Expires
May 14, 2021

Bayer Medical Care Inc
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
(412) 767-2400
www.radiology.bayer.com

To whom it may concern,

December 3, 2018

To be submitted to the Competent Health Authorities in Federal Republic of Russia

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Best regards,

William Hohn

Regulatory Affairs International Strategist – Europe
December 3, 2018

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©2017 Bayer. All Rights Reserved. Reproduction of this manual is strictly
prohibited without express written consent of Bayer Medical Care Inc.

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Bayer, the Bayer Cross, FluiDots, MEDRAD, Stellant, and MEDRAD Stellant may be
registered trademarks of Bayer in the US and other countries.

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UNITED STATES

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

English	5
Deutsch	5
Français	6
Italiano	6
Ελληνικά	7
Svenska	7
Nederlands	8
Dansk	8
Norsk	9
Português	9
Español	10
日本語	10
Magyar	11
Русский	11
Polski	12
Türkçe	13
繁體中文	14

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©2010-2014, 2017 Bayer. All Rights Reserved.
Reproduction of this manual is strictly prohibited without express written consent of Bayer Medical Care Inc.

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Bayer, the Bayer Cross, FluiDots and MEDRAD Stellant may be registered trademarks of Bayer in the US and other countries.

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ОО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Read the information contained in this instruction manual. This information will assist you in using the device in a safe manner.

Safety Notice: This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

Warnings: The contents of this package are intended for the delivery of contrast media or saline. Do not use for single-use on one patient only with Stellant Injectors.

Precautions: These devices are not intended for parent use, drug infusion, chemotherapy, or use for which the device is not indicated.

Caution: Rx Only

Warnings:

Contamination can cause death or serious injury. Do not connect a patient to the injector until the air has been cleared from the syringe plunger. Carefully read the instructions for use of the FluidDots indicators (where applicable) to reduce the chance of air embolism.

Contamination may result in biological contamination, product degradation and/or performance issues. Properly discard items after use, or if there is any possibility of contamination may have occurred.

Sterility will be compromised, and patient may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.

Contamination can occur if syringes are used to store contrast media. Use loaded syringes only. Discard unused loaded syringes.

Warnings labeled for single use, please note: This product is intended for single use only. Do not reuse, reprocess or reuse. The disposable syringe has been designed and validated for single use only. Re-use of the single use disposable syringe poses risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential patient injury could result if syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged.

Caution if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged contents are used. Visually inspect contents and package before each use.

Caution for operator injury may result from contrast leaks or tubing ruptures. Ensure that the fluid is not open; do not exceed pressures identified on front of the package. Use of greater pressures or obstructions in the fluid path may result in leaks or ruptures.

Caution for patient injury could result if syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged.

Caution in handling and inserting spike into contrast bottle. The spike is sharp and may cause personal injury.

Contamination may occur if the end of the spike is not properly sealed.

Automatic prime will not function properly and an injection may result if SJS-LP-60-T-J is used without an adapter. Do not use SJS-LP-60-T-J if an adapter is not installed on the injector.

Warnings:

Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten. This will help minimize leaks, disconnection, and component damage.


Refer to the Injector Operation Manual for further instructions.

Installing a Syringe
Insert the syringe until it snaps into place. The piston automatically engages the plunger, and advances it to the front of the syringe (most models). To reduce the volume and quantity of air bubbles drawn into the syringe during loading, a Fluid Displacement Device from Bayer is recommended.

Operator vigilance and care, coupled with a set procedure, is essential to minimizing the possibility of air embolism. Point the injector head up during loading. Point the injector head down during an injection.

Caution: To help avoid an air injection, syringes from Bayer are equipped with FluidDots indicators. FluidDots indicators could be observed as part of an arming procedure. When the FluidDots are viewed through an empty

syringe, the dots appear as small narrow ellipses. When viewed through a full syringe, the dots become larger, almost round. See below:

Empty Syringe Filled Syringe

To minimize air embolization risks, ensure that one operator is designated the responsibility of filling the syringe(s). Do not change operators during the procedure. If an operator change must occur, ensure that the new operator verifies that the fluid path is purged of air.

Loading and Priming a Syringe Using a Spike or a Quick Fill Tube (QFT)

The syringe can be loaded manually or automatically on the MEDRAD Stellant injector. See the MEDRAD Stellant Operation Manual for more detailed instructions.

MEDRAD Stellant Manual Loading:

1. Install a new syringe.
2. Install the spike or QFT onto the end of the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
3. Insert the spike or QFT into the fluid source and fill the syringe with fluid, using the load buttons.
4. Expel air.
5. Remove spike or QFT. Connect the disposable tubing set.
6. Follow the instructions in the section "Connector Tube Installation."

MEDRAD Stellant Integral Autoload with a spike or QFT:

1. Set and lock protocol.
2. Install a new syringe.
3. Install the spike or QFT onto the end of the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
4. Insert the spike or QFT into the fluid source and press Autoload.
5. Press Fill.
6. Expel air.
7. Remove spike or QFT. Connect the disposable tubing set.
8. Follow the instructions in the section "Connector Tube Installation."
9. Press Start/Hold.

NOTE: Patency check through aspiration using the injector manual knob is not possible when using a check valve on the end of the connector tube. If aspiration is important, remove the check valve from the connector tube and directly connect the connector tube to the catheter.

Removing a Syringe(s)

1. Disconnect the disposable tubing set from the vascular entry device. The disposable tubing set does not need to be disconnected from the syringe.
2. Rotate the syringe approximately 1/4 turn counter-clockwise and gently pull the syringe out of the injector head, discarding the syringe with disposable tubing set.

Note: Once the syringe is removed from the injector, the piston will automatically retract (most models).

Note: In order to remove the syringe, the last piston motion must be in the forward direction, which is typical. If you cannot remove the syringe, rotate the injector head manual knob approximately one revolution in the forward direction, then repeat Step 2.

Connector Tube Installation

1. Remove the connector tube from the package, removing the dust covers on the luer fittings.
2. Ensure all air is purged from the syringe.
3. Attach the connector tube to the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
4. Ensure that the connector luer fitting is secured to the tip of the syringe and verify that the tubing is not kinked or obstructed.

- If you are using a connector tube with T-connector, attach the straight portion of the T-connector to the contrast (Syringe A) and the extension to the saline (Syringe B). If the T-connector is connected to Syringe B, priming will not completely fill the connector tubing.
5. Connect the Prime Tube to the patient end of the tubing set, then prime the tubing by pressing Prime.
 6. Ensure all air is purged.
 7. Rotate the injector head downward.
 8. Remove the Prime Tube, connect to the patient and press Check For Air.

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Certified

5. at Annapolis, Maryland
6. the **28th day of January, 2019**
7. by The Secretary of State of Maryland
8. **No. 473914**

9. Seal

10. Signature

John C. W. W. W.

Secretary of State

State of Maryland, Anne Arundel County, Set

I, **SCOTT A. POYER**, Clerk of the Circuit Court for Anne Arundel County, Maryland, a court of record, do hereby certify that **Linda Davis** was a commissioned, appointed and qualified Notary Public commencing on the 1st day of May, 2017.

In Testimony Whereof, I have hereunto set my hand and affixed the seal of the court January 28, 2019

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
Best regards,

William Hohn Regulatory Affairs International Strategist – Europe

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