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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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
THIS IS TO CERTIFY THAT: the seal of XIAMEN COMPOWER MEDICAL TECH. CO., LTD. on the annexed DOCUMENT is genuine.

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Authorized

Signature:

Wei Shanshan

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(Date: Apr. 03, 2024)

证明书查询网址: site for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED
Xiamen Compower Medical Tech. Co., Ltd.

Organization
Unit 301, No. 16 Xianghong Road, Xiāng'an Torch Industrial Zone,
361101 Xiamen, People's Republic of China

Address

Director

Occupation


Frank Lin

Surname, Name

18/03/2024

Date DD.MM.YYYY

**Signature
Stamp**

Frank Lin

DRAFT

Non-reusable artificial manual breathing unit AMBU YUKIgroup, of different
versions

Operational manual

Contents

1	Medical device name	3
2	Information about the medical device manufacturer	4
3	Indented use and principle of operation of the medical device, stated by the manufacturer	4
4	Conditions of use	4
5	List of indications for use, contraindications, possible side effects	4
6	Description, parameters and characteristics of a medical device	5
7	Instructions for use	12
8	Installation and method of application	13
9	Maintenance and current repair	15
10	Scope of supply	15
11	Package	16
12	Marking	17
13	Operation conditions	19
14	Transportation	19
15	Storage	19
16	Warranty obligations	20
17	Procedure for the disposal and destruction of a medical device	20
18	Contacts	20

1 Medical device name

Non-reusable artificial manual breathing unit AMBU YUKIgroup, of different versions (hereinafter referred to as the device).

Versions:

1. **Non-reusable artificial manual breathing unit for adults consists of:**
 - breathing bag for adults.
 - PVC breathing mask for adults, sizes: No.4, 5.
 - patient valve.
 - PEEP valve (if necessary).
 - oxygen tube (if necessary).
 - reservoir bag 2500 (if necessary).
 - Godel airway, sizes: 2, 3, 4 (if necessary).
 - mouth expander (if necessary).
 - consumer packing (plastic box with transparent lid and carrying handle).
 - operational manual.

Short name: Non-reusable AMBU bag for adults

2. **Non-reusable artificial manual breathing unit for children consists of:**
 - breathing bag for children.
 - PVC breathing mask for children, sizes No. 1, 2, 3.
 - patient valve.
 - PEEP valve (if necessary).
 - oxygen tube (if necessary).
 - reservoir bag 1600 (if necessary).
 - Godel airway, sizes: 000, 00, 0, 1, 2, 3 (if necessary).
 - mouth expander (if necessary).
 - consumer packing (plastic box with transparent lid and carrying handle).
 - operational manual.

Short name: Non-reusable AMBU bag for children.

3. **Non-reusable artificial manual breathing unit for newborns consists of:**
 - breathing bag for newborns.
 - PVC breathing mask for newborns, sizes: No. 00 (42 mm), No. 0, 1.
 - patient valve.
 - PEEP valve (if necessary).
 - oxygen tube (if necessary).

- reservoir bag 1600 (if necessary).
- model airway, sizes: 000, 00, 0 (if necessary).
- mouth expander (if necessary).
- consumer packing (plastic box with transparent lid and carrying handle).
- operational manual.

Short name: Non-reusable AMBU bag for newborns

2 Information about the medical device manufacturer

Xiamen Compower Medical Tech. Co., Ltd.

Unit 301, No. 16 Xianghong Road, Xiang'an Torch Industrial Zone, Xiamen, 361101, China

Tel.: +86 592 5799319, Fax: +86 592 5796765, e-mail: sales-compower@cpmedi.com, web-site: <http://www.cpmedi.com/en/>

3 Indented use and principle of operation of the medical device, stated by the manufacturer

The device is used for forced artificial lung ventilation (ALV) during cardiopulmonary resuscitation (CPR) in newborns, children and adults with unstable breathing.

The product is a non-invasive device and is used exclusively for temporary ventilation due to the absence of any monitoring of the delivered volume of air, monitoring of pressure in the respiratory tract and the fatigue of the process.

Principle of operation is based on the forced supply of air to the lungs through the facial mask by manually compressing the ventilation bag. and the device consists of a large elastic chamber, manually ventilated, a gas tank, tubes and a connector for attaching to a facial mask or endotracheal tube.

4 Conditions of use

The device is used by specialist doctors, in emergency medical services in ambulances, intensive care units, during internal movements of the patient and incidents with mass casualties.

5 List of indications for use, contraindications, possible side effects

5.1 Indications for use, contraindications for use, possible side effects

5.1.1 Indications for use:

- artificial myoplegia in a complex anesthetic aid and in intensive therapy of convulsive, hyperventilation and hyperthermic syndromes,

- apnea or hypoventilation due to disorders of the central regulation of respiration (clinical death, cerebrovascular accident, brain damage, etc.).
- apnea or hypoventilation due to a violation of the mechanical apparatus of respiration (respiratory muscle damage, pneumothorax, including surgery, etc.).
- acute respiratory failure due to restrictive and obstructive lung lesions, accompanied by increased work of the respiratory muscles that absorb a lot of oxygen.

5.1.2 Contraindications to use:

- tense pneumothorax, or mediastinal emphysema without decompression and drainage.
- moderate or more serious hemoptysis.
- severe case of lung cyst or bullous lung disease.
- hypovolemic shock without replenishing the volume of circulating blood.
- acute myocardial infarction.
- suffocating respiratory failure caused by inhaling a foreign body.

5.1.3 Possible side effects:

- regurgitation and aspiration of gastric contents.

6 Description, parameters and characteristics of a medical device

6.1 Device description

The device can include:

- breathing bag for newborns: – 1 pcs.
- breathing bag for children: – 1 pcs.
- breathing bag for adults: – 1 pcs.
- PVC breathing mask, sized: No. 00 (42 mm), 0, 1, 2, 3, 4, 5 – from 2 to 3 pcs.
- patient valve – 1 pcs.
- PEEP valve (if necessary) – 1 pcs.
- oxygen tube (if necessary) – 1 pcs.
- reservoir bag: 1600 or 2500 (if necessary) – 1 pcs.
- Godel airway, sizes: 000, 00, 0, 1, 2, 3, 4 (if necessary) – from 3 to 6 pcs.

- **mouth expander (if necessary) – 1 pcs.**
- **consumer packing (plastic box with transparent lid and carrying handle).**
- **operational manual.**

The appearance of the device is shown in Figure 1, a schematic image is in Figure 2.


Figure 1 – The appearance of the device


Figure 2 – A schematic image of the device

The product is intended for one-time use in one patient, it is not sterile, is not subject to repeated sterilization.

The breathing bag is a self-expanding PVC bag with an outlet valve and a patient valve, on the nozzle of which a facial mask is attached. The structured surface of the respiratory bag prevents slipping in the doctor's hands. A 360-degree swivel adapter facilitates the provision of medical care, the adapter is equipped with a petal valve of the patient and a pressure limiting valve in order to prevent lung barotrauma.

PEEP valve is used for monitoring of pressure in the patient lungs.

The reservoir bag is designed to save oxygen during ventilation..

The PVC breathing mask is equipped with an inflatable cuff with a nipple valve and a Luer adapter.

Oxygen tube 2 m long with a star-shaped inner lumen, resistant to sticking during kinks, is a tube with connectors at the ends for connection to an oxygen source (balloon) and to the intake valve of the breathing bag.

The Godel airway helps to ensure the passage of air into the respiratory tract. The airway should keep the patient's airways open so that the tongue (partially or completely) does not cover the epiglottis cartilage, which could make it difficult for the patient to breathe.

The Godel airway has a bite block with a color designation of the dimensions specified in Table 1, which is designed for easy identification, prevention of biting and blocking of the airway. It is made of two compatible materials with different properties – polyethylene and polyoxymetelene, connected during casting, with different properties: the inner solid to stiffen the structure and the outer plastic to prevent tooth damage.

Table 1 – Color coding of the dimensions of the Godel airway

Length, mm	Dimension	Color marking
40	000	pink
50	00	blue
60	0	black
70	1	white
80	2	green
90	3	yellow
100	4	red

The mouth expander is designed to safely open the patient's mouth to ensure the patency of the respiratory tract in life-threatening conditions.

The PEEP valve provides positive pressure at the end of exhalation – it is an auxiliary means of lung ventilation to maintain the vital capacity of the lungs. At the end of exhalation, PEEP creates a pressure of 2 to 20 cm H₂O to counteract passive emptying of the lungs and maintain airway pressure above atmospheric. Thus, the introduction of PEEP can open destroyed or unstable alveoli and increase lung capacity and gas exchange surface area.

If it is necessary to perform an ALV with the use of an additional volume of oxygen, the device is connected to an oxygen cylinder using an oxygen tube. One end of the tube is connected to the oxygen source, the other end of the tube is connected to the connector of the intake valve cover, and the reservoir bag is attached to the intake valve cover.

An artificial lung ventilation device, if it is connected to a medical oxygen source, can supply additional oxygen using a reservoir bag. When used, oxygen enters the reservoir during compression of the breathing bag.

6.2 Parameters and characteristics

Technical characteristics of the device are specified in Table 2.



Table 2 – Technical characteristics

Parameter name	Sizes of breathing bags		
	For newborns	For children	For adults
1 Respiratory volume of the breathing bag, ml	from 20 to 150	from 150 to 450	from 600 to 1000
2 Volume of the breathing bag, ml	280±100	550±200	1500±200
3 The volume of the dead space of the respiratory bag, ml, no more	7	20	65
4 Dimensions of the breathing bag (length × diameter), mm	180×75, permissible deviation ±20	180×100, permissible deviation ±20	251×140, permissible deviation ±20
5 Dimensions of the assembled device (length × width × height), mm	645×195×75, permissible deviation ±20	645×195×100, permissible deviation ±20	730×225×130, permissible deviation ±20
6 Mass of the assembled device, g, not more	900		
7 The time of complete straightening of the breathing bag, sec., no more	1	1	2
8 Outer diameter of the pressure limiting valves, mm	25.0±0.5		
9 Inner diameter of the pressure limiting valves, mm	18.4±0.2		
10 Length of the pressure limiting valves, mm	25.0±0.5		
11 Flow capacity of the pressure limiting valve, cm H ₂ O	30-45	35 – 50	40 – 60
12 Inhalation resistance, no more	5 cm H ₂ O (at 50 l/min)		
13 Exhalation resistance, no more	5 cm H ₂ O (at 50 l/min)		
14 The volume of the reservoir bag, ml, no more	1600	1600	2500
15 Dimensions of the reservoir bag (length × width × height), mm	275×195×0.16 permissible deviation ±10%	270×195×0.16 permissible deviation ±10%	295×225×0.16 permissible deviation ±10%
16 Patient connector	according to ISO 5356-1: Ø22/15 mm		
17 Outer diameter of the oxygen tube, mm	6.0 ± 0.2		

19 Inner diameter of the oxygen tube, mm	4.0 ± 0.2		
20 Length of the oxygen tube, mm	2000±150		
20 The volume level of noise generated by the product, dBA, no more than	30		
21 Dimensions of the mouth expander (length × diameter), mm	82.4 × 59.2 permissible deviation ± 0.5		
22 Dimensions of PEEP valve (length × diameter), mm	150 × 37.6 permissible deviation ± 2.0		
23 The flow capacity of the PEEP valve, ml/min	500 ± 10		
24 Oxygen concentration:			
24.1 Oxygen supply rate, l/min, no more	15	15	15
24.2 Respiratory rate, respiratory movements/min, no more	20	25	12
24.3 Inhalation volume, ml, no more	60	150	600
24.4 Oxygen concentration with a reservoir bag, %, no more*	85		
24.5 Oxygen concentration without reservoir bag, %, no more*	85		
25 Patient body weight, kg	not more 10	from 10 to 30	not less 40

*Test methods according to ISO 10651-4, the oxygen source must provide at least 85% oxygen concentration in the oxygen-air mixture at a flow of 2-8 l/min.

Technical characteristics of the breathing masks are specified in Table 3.

- 1 – Patient connector
 2 – Luer-Lock Adapter
 A, B, C – dimensions are specified in Table 3.

Figure 3 - Image of silicone breathing masks

Table 3 – Technical characteristics of PVC breathing masks

Mask type	Connector dimensions, mm	Size A, mm	Size B, mm	Size C, mm	Mass, g
1 For newborns (size No. 00)	Ø15	38±2	42±2	42±2	11±0.5
2 For newborns (size No. 0)	Ø15	40±2	48±2	48±2	12.3±0.5
3 Nursing infant (size No.1)	Ø15	54±2	62±2	82±2	26±0.5
4 Standard for children (size No.2)	Ø22	64±2	74±2	96±2	45±0.5
5 For children (size No. 3)	Ø22	65±2	68±2	108.5±2	43±0.5
6 For adults M (size No. 4)	Ø22	75±2	77±2	125±2	59±0.5
7 For adults L (size No.5)	Ø22	80±2	84±2	137±2	73±0.5

Technical characteristics of Godel airway are specified in Table 4.

Table 4 – Technical characteristics of Godel airway

Dimensions	Inner diameter, mm	Outer diameter, mm	Length, mm	Mass, g, not more	Length of the reinforcing liner, mm
000	5.0±0.2	10±0.8	40±2	2.57±0.4	38.0±2.5
00	5.5±0.2	11±0.8	50±2	3.9±0.4	48.0±2.5
0	6.0±0.3	11±0.8	60±3	5.8±0.4	58.0±3.0
1	6.0±0.3	14±0.8	70±3	9±0.4	68.0±3.0
2	6.5±0.3	14±0.8	80±4	10.7±0.4	78.0±4.0
3	8.5±0.3	15±0.8	90±4	13.2±0.4	88.0±4.0
4	8.0±0.3	17.5±0.8	100±5	15.7±0.4	98.0±5.0

Materials, applied for device manufacturing are specified in Table 5.

Table 5 – Materials

No.	Name of the component part	Material	Brand	Manufacturer
1	breathing mask	PVC	MRI 2310-15	Changzhou Healthcare, China
2	Patient valve	PC	L1250Y	TEIJIN CHEMICALS LTD, China
3	PEEP valve	PC	L1250Y	TEIJIN CHEMICALS LTD, China
4	breathing bag	PVC	MI-3250LT-18	Changzhou Healthcare, China
5	output valve	PC	L1250Y	TEIJIN CHEMICALS LTD, China
6	oxygen tube	PVC	MT-1182-16	Changzhou Healthcare, China
7	reservoir bag	EVA	8038	Fujian Refining & Petrochemical Co., China
8	Godel airway	Polyethylene	17703-010	Formosa Plastics, China
9	Mouth expander	Polyethylene	17703-010	Formosa Plastics, China

Note:

PVC – polyvinyl chloride.

PC – polycarbonate.

EVA – ethylene vinyl acetate.

6.3 Methods and conditions of sterilization, terms of preservation of sterility

The product is non-reusable, not sterile and is not subject to repeated sterilization.

7 Instructions for use

7.1 General issues

7.1.1 Before operating with the product, read the Operating Manual.

7.1.2 The product is intended for use by personnel trained in the technique of pulmonary resuscitation.

7.2 Safety requirements

7.2.1 Before using the product, it is necessary to clean the patient's respiratory tract.

7.2.2 After unpacking and assembling, it is necessary to check the operability of the product.

7.2.3 It is necessary to check the operability of the valve by observing, during resuscitation, the respiratory movements of the chest, the color of the skin of the face and the mucous membrane of the lips.

7.2.4 Do not use the product in the presence of toxic or explosive substances in the air! If the device sucks gas from such an atmosphere and allows it to be inhaled by the patient, then using the device in a polluted environment can be dangerous if measures are not taken to prevent such suction or inhalation, for example, using a filter.

7.2.5 At high oxygen concentrations – it is forbidden to use an open flame and smoking near the product – this can lead to fire!

7.2.6 It is forbidden to use oil or lubricants near oxygen devices – this may lead to a fire!

7.2.7 The product is disposable and intended for use in one patient. It is not subject to repeated sterilization. Repeated use may cause patients to develop a cross infection.

8 Installation and method of application

8.1 Device assembly

8.1.1 Check the integrity of the package. Carefully remove the product.

8.1.2 Fully straighten the reservoir bag before use.

8.1.3 Squeezing the breathing bag, visually check the operability of the patient's valve.

8.1.4 Connect the patient's valve to the breathing mask adapter.

8.1.5 In case of resuscitation with oxygen supply, fully deploy one connector of the the oxygen to the oxygen source and the other to the oxygen port on the intake valve of the breathing bag. Connect the tank bag.

8.1.6 Set the oxygen supply level to no more than 15 l/min or as directed by a doctor.

8.2 Health check

Close the pressure limiting valve, block the patient's valve connection with your thumb. Bring a breathing bag filled with compressed air to your ear. Squeeze the

bag, the bag should be hard to the touch, the sound of leakage should not be heard. Open the pressure limiting valve. Remove the thumb from the patient's valve, squeeze the breathing bag, the sound of air passing through the patient's valve should be heard.

8.3 Method of application

8.3.1 Open the patient's mouth, clear the respiratory tract of foreign objects and liquids. It is recommended to use resuscitation suction. Tilt the patient's head back and extend the lower jaw until the neck is straightened as completely as possible and the airways open.

In cases where it is not possible to open the patient's mouth, use a mouth expander: pulling the corner of the mouth and cheek to the bottom, insert the narrow edge of the mouth expander between the molars and screw the distance between the upper and lower jaw clockwise.

8.3.2 For more adequate ventilation, it is recommended to install an airway.

8.3.2.1 Insert the airway into the oral cavity with the curved side to the tongue.

8.3.2.2 When the airway reaches the back wall of the pharynx, turn it 180° – it will press the root of the tongue and the epiglottis, creating free airway patency. The criteria for the correct position of the product are free spontaneous breathing and unhindered ventilation of the lungs.

8.3.2.3 Make sure that the airway does not move the tongue backwards and thus does not block the airway.

8.3.3 Firmly press the breathing mask to the patient's face, closing the mouth and nose, tilt the head back as much as possible, holding the lower jaw with the same hand with which you hold the breathing mask. Gently squeeze the breathing bag and observe the expansion of the chest.

8.3.4 Loosen the pressure on the breathing bag and let the chest get in. Repeat 12-20 times a minute, for children – 30 times a minute.

If you feel an increase in resistance, check the patency of the respiratory tract or the correct position of the patient's head. If using the product, it is not possible to achieve adequate ventilation, immediately proceed to the expiratory technique (mouth-to-mouth, mouth-to-nose).

8.3.5 If the patient begins vomiting during mask ventilation, it is necessary to immediately clear the airways, then compress the respiratory bag several times in a row until normal rhythmic ventilation resumes.

8.3.6 For the adult breathing bag of the size range, it is necessary to install a pressure limiting valve for 60 cm H₂O.

8.3.7 For a bag of respiratory neonatal and pediatric size range, it is necessary to install a pressure limiting valve for 40 cm H₂O.

When the valve is open, a hissing sound may be heard.

In some patients, high pressure during ventilation can cause lung damage or stomach bloating.

8.3.8 If it is necessary to increase the pressure, it is necessary to press and turn the button during compression of the breathing bag – the valve is disconnected.

8.3.9 Adjust the oxygen supply level according to the Specification or for medical reasons.

8.3.10 After the end of use, remove the product in accordance with medical recommendations and protocols.

8.3.11 The used product must be disposed of.

9 Maintenance and current repair

The product is not subject to maintenance and repair.

10 Scope of supply

Scope of supply of the product is specified in Table 6.

Table 6 – Scope of supply

Name	Quantity, pcs.
1 Non-reusable artificial manual breathing unit for adults, consisting of:	1
– breathing bag for adults.	1
– PVC breathing bag for adults, sizes: No. 4, 5.	2
– patient valve.	1
– PEEP valve (if necessary).	1
– oxygen tube (if necessary).	1
– reservoir bag 2500 (if necessary).	1
– Godel airway, sizes: 2, 3, 4 (if necessary).	3
– mouth expander (if necessary).	1

- consumer packing (plastic box with transparent lid and carrying handle)	1
- operational manual	1
2 Non-reusable artificial manual breathing unit for children, consisting of:	1
- breathing bag for children.	1
- PVC breathing bag for children, sizes: No. 1, 2, 3.	3
- patient valve.	1
- PEEP valve (if necessary).	1
- oxygen tube (if necessary).	1
- reservoir bag 1600 (if necessary).	1
- Godel airway, sizes: 000, 00, 0, 1, 2, 3 (if necessary).	** 6
- mouth expander (if necessary).	1
- consumer packing (plastic box with transparent lid and carrying handle)	1
- operational manual	1
3 Non-reusable artificial manual breathing unit for newborns, consisting of:	1
- breathing bag for newborns.	1
- PVC breathing bag for children, sizes: No. 00 (42 mm), 0, 1.	2
- patient valve.	1
- PEEP valve (if necessary).	1
- oxygen tube (if necessary).	1
- reservoir bag 1600 (if necessary).	1
- Godel airway, sizes: 000, 00, 0 (if necessary).	3
- mouth expander (if necessary).	1
- consumer packing (plastic box with transparent lid and carrying handle)	1
- operational manual	1

11 Package

Each component part of the product, except for the breathing bag, must be packed in an individual polyethylene bag. The edges of the package are secured with various types of fasteners: adhesive flap, zip lock, adhesive tape or sealed.

The composite products are packed in a plastic box with a transparent lid and a carrying handle, along with an instruction manual. Overall dimensions of the plastic

length × width × height), – 270× 160× 140 mm, permissible deviation is ± 20 mm. Box weight is 0.9 ± 0.1 kg.

The plastic box is packed in a polyethylene bag and placed 12 pieces each in a transport package, which is a cardboard box.

12 Marking

12.1 The pressure limiting valve is marked with the indication of the maintained pressure, on the breathing bag – the size of the bag and information about sterilization, on the breathing mask – the maxi size, on the reservoir bag – the name and volume.

12.2 The following markings must be applied to the plastic box for storage and carrying:

- name and/or trademark of the manufacturer and its address (indicated by the corresponding symbol).
- name of the manufacturer's authorized organization and contact details.
- device name.
- device dimensions.
- date of manufacture (indicated by the corresponding symbol).
- expiration date (indicated by the corresponding symbol).
- batch code (indicated by the corresponding symbol).
- catalog number.
- sterility information (indicated by the corresponding symbol).
- information about the phthalate content (indicated by the corresponding symbol).
- warning about the need to refer to the operating manual (indicated by the corresponding symbol).
- warning about not using the product if the packaging is damaged (indicated by the corresponding symbol).
- warning about the prohibition of repeated use (indicated by the corresponding symbol).
- designation of storage conditions.
- registration certificate number and date.

12.3 Marking of transportation package includes

- name and/or trademark of the manufacturer.

- device name.
- year and month of packaging.
- designation of storage and transportation conditions.

Designation of applied symbols is specified in Table 7.

Table 7 – Designation of applied symbols

Symbol image	Symbol designation
	Batch code
	Catalog number
	Manufacturing date
	Manufacturer
	Expiration date ...
	Refer to User manual
	Do not use if the package is damaged
	Prohibition of repeated use
	Does not contain phthalates
	Non-sterile
	Keep away from moisture
	Avoid exposure to sunlight

	Temperature range
---	-------------------

13 Operating conditions

The device is used indoors or on open surfaces under the following climatic conditions:

air temperature from minus 18°C to plus 50 °C.

relative air humidity from 40 to 95%.


atmospheric pressure 84.0-106.7 kPa (630-800 mm Hg).

14 Transportation

The packaged device is transported by all types of closed vehicles, in accordance with the rules of cargo transportation applicable to this type of transport.


The device is transported under the following climatic conditions:

air temperature from minus 40 to plus 60 °C.


relative air humidity from 40 to 95%.


atmospheric pressure 84.0-106.7 kPa (630-800 mm Hg).

15 Storage

The device is stored indoors under the following climatic conditions:

air temperature from minus 40 to plus 60 °C

relative air humidity from 40 to 95%.

atmospheric pressure 84.0-106.7 kPa (630-800 mm Hg).

It is necessary to store the device away from heating devices, in a place protected from the sun.

16 Warranty obligations

The manufacturer guarantees that the device complies with the quality indicators and characteristics stated in the documentation, subject to the rules of transportation and storage throughout the shelf life.

The manufacturer is not responsible for incidental or consequential damages, damages or expenses directly or indirectly related to the use of this device; he is not responsible for products that are re-processed or sterilized, reused, exposed to climatic and mechanical factors during transportation and storage, other than those specified in the manufacturer's documentation, and does not give any guarantees, express or implied, including, but not limited to, guarantees of commercial quality and suitability for intended use in with respect to such products.

The expiration date of the product is 3 years from the date of manufacture.

The shelf life of the product is 3 years from the date of production.

17 Procedure for the disposal and destruction of a medical device

After using the product for its intended purpose, it must be disposed of as potentially hazardous waste of Class B.

Upon expiration of the expiration date and/or violation of the integrity of the sealed packaging, the medical device must be disposed of as Class A medical waste.

Disposal is carried out in accordance with the rules in force in the country of application at the time of disposal.

18 Contacts

In case of questions related to the use of the product, as well as in case of claims against the manufacturer, the consumer can contact the specified address:

Authorized organization for quality issues:

The limited liability company (OOO) ARMEDIKA, Address: 620075, Russia,
Ekaterinburg, ul. Shevchenko, 20, prem. 7, Tel.: +7 (343)-286-42-73, E-mail:
armika@yukigroup.ru, Web-site: yukigroup.ru.

Annex A

Information about regulatory documentation

List of applied international standards

- EN ISO 13485:2016 Medical products. Quality management systems. System requirements for regulatory purposes.
- EEC Directive 93/42/EEC (certificate No. G2063744 0018, valid until "26" 05 2024).
- EN 1041:2013 Information prepared by the manufacturer accompanying medical devices.
- EN ISO 14971:2019 Medical devices. Application of risk management to medical devices.
- EN ISO 10651-4:2002 Artificial lung ventilation medical devices. Part 4. Specific safety requirements, considering the basic functional characteristics of artificial lung ventilation devices for manual actuation.
- EN ISO 10651-5:2006 Artificial lung ventilation medical devices. Part 5. Safety requirements, considering the main functional characteristics of artificial lung ventilation devices for reviving.
- EN ISO 5364:2016 Equipment for anesthesia and artificial respiration. Oropharyngeal air ducts.
- ISO 10993-1:2018 Evaluation of biological medical devices. Part 1. Assessment and testing within the framework of the risk management process.
- ISO 10993-5:2009 Evaluation of biological medical devices. Part 5. Cytotoxicity tests in vitro.
- ISO 10993-10:2010 Evaluation of biological medical devices. Part 10. Tests for irritation and allergic skin reaction.
- Certificate of conformity ISO 13485:2016, registration no. Q6 063744 0017, date of issue "01" 10 2019.

List of standards, applied on the territory of the Russian Federation

- GOST R IEC 62366-1-2021 Medical products. Part 1. Application of usability engineering to medical devices.
- GOST R ISO 10651-4-2015 Lung ventilators. Part 4. Particular requirements for operator-powered resuscitators

- GOST R ISO 10651-5-2015 Lung ventilators for medical use. Part 5. Particular requirements for basic safety and essential performance for gas-powered emergency resuscitators
- GOST 31057-2012 Resuscitators intended for use with humans. General technical requirements and test methods
- GOST R ISO 15223-1-2020 Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements
- GOST ISO 10993-1-2011 Medical devices. Assessment of the biological effect of medical devices. Part. 1. Evaluation and research
- GOST R ISO 10993-2-2009 Medical devices. Biological evaluation of medical devices. Part 2. Animal welfare requirements
- GOST ISO 10993-5-2011 Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity
- GOST ISO 10993-10-2011 Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and delayed-type hypersensitivity

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