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Nantong Egens Biotechnology Co. Ltd

Address: Building 15, Building 12 (west), No.1692 Xinghu Avenue
Nantong Economy & Technology Developments Zone
Jiangsu Province China
Tel.: 86-0513-853-27458, Fax: 86-0513-85327120

Instruction for use of medical device in vitro (for registration in Russia)

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Manufacturer:

Nantong EGENS Biotechnology Co., Ltd
Building 15, Building 12(west) , No.1692 Xinghu Avenue
Nantong Economy & Technology Development Zone
Jiangsu Province China
Tel.: 86-0513-85327458, Fax: 86-0513-85327120
e-mail: egens@egens-bio.cn

Authorized representative in Russia:

FC GRAND CAPITAL Ltd.
Russian Federation, 125284, Moscow, Leningradsky avenue, h. 31A, bld. 1,
floor/room 14/1
Tel.: +7 (495)258-34-28
e-mail: info@grand-capital.pro

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- Pregnancy test-strip – 1 ps;
- Instruction for use – 1 ps.

2) Pregnancy test Bona Mente! С добрыми намерениями! in variants of execution in composition of:


- Pregnancy test-strips– 2 pics;
- Instruction for use – 1 ps.

1.2. Intended use: Pregnancy Test is a rapid chromatographic immunoassay for the qualitative detection of human Chorionic Gonadotropin (hCG) in urine samples to aid in the early detection of pregnancy by lay users.

1.3. Scope of application: Clinical laboratory diagnostics. The medical device is used independently by patients at home. The medical device is intended for self-use by persons without medical training at home. The pregnancy test is intended for one-time use.

1.4. Principle: Pregnancy Test for is a self-performing immunoassay designed for the qualitative determination of human chorionic gonadotropin (HCG) in urine for early detection of pregnancy. Human Chorionic gonadotropin (HCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. During normal pregnancy, HCG can be detected in urine as early as 7 days following conception, doubling every 1.3 to 2 days. At the time of the last missed menstrual period, urine HCG levels are about 100mIU/ml with peak levels of 100,000 to 200,000mIU/ml seen at the end of the first trimester. The presence of HCG soon after conception and its subsequent increase in concentration during early gestational growth make it an ideal marker for the early detection of pregnancy.

The method employs a unique combination of monoclonal dye conjugate and monoclonal-solid phase antibodies to selectively identify the HCG in the test samples with an extremely high degree of sensitivity. At 5 minutes, levels of HCG as 25mIU/ml can be detected.

A membrane with an absorbent pad over a strip of fiber glass paper is impregnated with a colloidal conjugate of gold particles and monoclonal antibodies to HCG. Other absorbent pads at the end of the assay absorb excess sample fluid. The urine sample is introduced into a chromatographic membrane (A). As it contacts the membrane, the sample dissolves the lyophilized conjugate (B). In a reactive sample, the HCG antigen will attach to the antibodies in the colloidal solution. As the conjugate moves forward on the membrane, anti-HCG monoclonal antibodies affixed on the test zone ("T") (C) will bind the HCG-gold conjugate complex, forming a pink line ("T"). Any sample will cause a pink colored line to appear in the control zone ("C") (D). This line is formed by the binding of polyclonal antibodies (Anti-mouse IgG) affixed onto the control zone to the sample-colloidal gold conjugate. Presence of this line indicates that the test has been carried out correctly.

Validity: 3 years.

2. Specification and materials of the product

2.1. Parts list

- Test Strip
- Desiccant
- Foil pouch
- Carton box
- Packing tape (which is provided after registration of medical device in Russian Federation)

3. Contents

Each box contains: 1 (one) or 2 (two) pouches and the operating instruction.

Each pouch contains: the Pregnancy Test Strip and the Desiccant.

4. Precautions

- FOR IN-VITRO DIAGNOSTIC USE ONLY (IVD), NOT for internal use.
- Read directions for use carefully before performing this test. Perform each step exactly as described in IFU.
- Do not use beyond the expiration date marked on the foil pouch.
- Each one step HCG pregnancy Test can only be used once; do not re-use the one step HCG pregnancy Test.
- Do not use the Test if the foil pouch is damaged.
- Once the foil pouch has been opened, the Test should be used immediately.
- Discard each one step HCG pregnancy Test after use. Treat urine samples and used test devices as if they were potentially infectious. Avoid contact with skin.
- KEEP OUT OF REACH OF CHILDREN.


5. TEST PROCEDURE

SAMPLE COLLECTION AND HANDLING

To perform this test, collect a urine sample into a clean and dry container. Urine samples collected first thing in the morning are ideal because they generally contain higher concentration of HCG. However, urine samples collected at any time of the day may also be used.

If it is impossible to perform a pregnancy test immediately after collecting urine, storage of biomaterial is no more than 8 hours at room temperature 20-25°C.


TO CARRY OUT THE TEST

Absorbent cotton Test line Control line

DIRECTIONS FOR USE

- Bring the sealed pouch to room temperature.
- Remove the strip from the foil pouch.
- Immerse the strip into the urine with the arrow end pointing toward the urine cup. Do not cover the urine over the MAX (maximum) line.

- Take the strip out of the cup after a minimum of 15 seconds and lay it flatly on a non-absorptive, clean surface.
- Read the result in 5 minutes.
- Discard the test device after a single use.

Depending on the concentration of HCG in the urine sample, positive results may be seen within 40 seconds. **Wait the full five minutes to confirm the result.** Do not read the one step HCG pregnancy Test after 5 minutes as this may give an incorrect test result.

This is a single use test. Please safely dispose of the strip, regard this as infectious material and dispose of the test appropriately after use.

HOW TO READ THE RESULT

Positive: If two pink/purple colored lines appear on the one step HCG pregnancy Test, one in the control area and one in the test area, then this means HCG has been detected in the urine and there is a strong possibility of pregnancy. The lines seen may be faint and one line may be lighter than the other; however this is still a positive result. The lines do not have to match in strength or color.

Negative: If one colored pink/purple line appears in the control area on the one step HCG pregnancy Test and there is no line in the test area, then this indicates a negative result. This means no HCG has been detected in the urine. This indicates you are either not pregnant or you have tested too early. If you are not sure repeat the test in 48 hours.

Invalid: If no lines appear anywhere on the control zone, the test is invalid and repeat the test using another strip.

Invalid Test Results can be caused by:

- The absorbent tip not being wetted thoroughly.
- Reading the test result too early or too late. A 5 minutes reaction time is required.

For further information, see LIMITATIONS OF THE PROCEDURE.

6. QUALITY CONTROL

Built in Quality Control Features:

The appearance of a line in the control area indicates that the one step HCG pregnancy Test is performing properly and serves as a procedural control. It is recommended that a positive hCG control (containing 25-250mIU/mL hCG) and a negative hCG control (containing 0mIU/mL hCG) are evaluated to verify proper test performance when a new shipment of tests is received.

7. Performance characteristics

Sensitivity

Pregnancy Test (during one-stage) will display positive results with specimens containing HCG at levels of 25mIU/ml or greater.


Accuracy

Comparison studies on the Pregnancy Test with a legally marketed device were performed in-house and in a clinical reference laboratory. Positive and negative results were compared and the correlation was > 99%.

Specificity

The following compounds exhibited no interference when dissolved in urine, which had HCG levels of 0 and 25 mIU/ml.

hLH.....	500mIU/ml
hFSH.....	1000mIU/ml
hTSH.....	1000μIU/ml

There is no interference observed for the substances in the urine at below concentration.

Paracetamol	20 mg/dl	Glucose	2000 mg/dl
Ascorbic acid	20mg/dL	Bilirubin	2 mg/dl
Caffeine	20mg/dL	Atropine	20mg/dL
Phenothiazine or its derivatives*	20mg/dL	Hemoglobin	500 mg/dl
EDTA	80mg/dL	Thiophene or its derivatives**	20mg/dL
Acetylsalicylic acid	20mg/dL	Ketone	20mg/dL
Ethanol	1,0 %	Tetracycline	20mg/dL
Albumin	2000 mg/dl		

*Compounds belonging to the phenothiazine group, according to the chemical classification.

** Compounds containing a thiophene ring.

No interference observed when the urine samples at pH 4 ~ 9; Density, 1.000 ~ 1.035g/mL.

Sensitivity and specificity

The pregnancy test strip allows to determine the presence of HCG in concentrations from 25 mIU/ ml.

During clinical trials, clinical efficacy was evaluated, and the results obtained on the Test product were proven to be consistent with the results obtained on the Comparison Product:

- Diagnostic sensitivity – 100 %,
- Diagnostic specificity – 100 % .

Cross-Reactivity

There's no cross-reactivity observed with hLH (500mIU/mL), hFSH (1000mIU/mL), hTSH (1000µIU/mL).

Intra-Assay/ Inter-Assay

Assay were carried out to determine assay reproducibility using replicates of 10 tests for each of three lots by using the specimen from the HCG QC panel, and observe the color of the T line. The color of the T line is the same.

Hook study

The product still shows positive results at HCG concentration of 2000IU/ml.


Studies with Lay person and clinical

100 lay persons with variant ages, educations participant the study, 3 sites involved the lay person study, the test performed by lay persons after reading the IFU. The reference assay is Acon Biotech (Hangzhou) Co., Ltd. And the conformity is 100%.

Stability studies

The stability studies included the Accelerated Stability, Real-time stability, after first opening stability, and transportation stability.

Accelerated Stability Summary

Accelerated stability study was conducted. Three pilot production lots will be stored at 50°C to demonstrate product stability. A series of stability tests were performed at 1, 20, 30, 40, 50, 60, 70, 80, 90, 100, 110, 120, 125 days. Three Calibrators (0mIU/ml HCG, 5mIU/ml HCG, another HCG sample is prepared according to the product sensitivity) were tested in triplicate according to the instructions specified in the package insert.

Results show that the product is stable at 50°C for 125 days. The data was plotted on an Arrhenius plot and the shelf life was determined to be at least 36 months from the date of manufacture.

Real-time stability summary

Real Time stability study was conducted. Tests were evaluated using samples from three different batches. The real-time stability study is set up at the room temperatures (4°C, and 30°C).

A series of stability tests were performed at 1, 4, 7, 10, 13, 16, 19, 22, 25, 27, 29, 31, 33, 35, 36, 37 months. At each predetermined time-point, the appropriate number of devices are pulled out and allowed to return to room temperature in our laboratory. Three Calibrators (0mIU/ml HCG, 5mIU/ml HCG, another HCG sample is prepared according to the product sensitivity) were tested in triplicate according to the instructions specified in the package insert.

Results show that the product is stable at 4°C, and 30°C for 37 months. The data support the indication of product shelf life was determined to be at least 36 months from the date of manufacture.

After first opening stability summary

After first opening stability study was conducted. Tests were evaluated using samples from three different batches. The After first opening stability study is set up at 4 conditions (Humidity: 50%~70%, 6h after Opening, Humidity: 50%~70%, 24h after Opening, Humidity: ≥70%, 6h after Opening, and Humidity: ≥70%, 24h after Opening)

A series of stability tests were performed at these 4 conditions. At each predetermined point, three Calibrators (0mIU/ml HCG, 5mIU/ml HCG, another HCG sample is prepared according to the product sensitivity) were tested in triplicate according to the instructions specified in the package insert.

Results show that the product is stable at Humidity: 50%~70%, 6h after Opening.

Transportation stability summary

Transportation stability study was conducted. Tests were evaluated using samples from three different batches. The Transportation stability study is set up at the extreme temperatures (-20°C, and 40°C) for 35 days.

A series of stability tests were performed at 1, 6, 14, 21, 24, 28, 31, 25 days since the study started. At each predetermined time-point, the appropriate number of devices are pulled out and allowed to return to room temperature in our laboratory. Three Calibrators (0mIU/ml HCG, 5mIU/ml HCG, another HCG sample is prepared according to the product sensitivity) were tested in triplicate according to the instructions specified in the package insert.

Also, a simulate transport stability study according to ASTM-D4169-16 was performed. Results show that the product is stable at cold, hot and vibration transport environment.

Limitation of the procedure

Performing the test too early can yield a false negative due to low HCG levels. In this case, another urine specimen should be obtained at least 48 hours later and tested. HCG levels may remain detectable for several weeks after normal delivery, delivery by caesarean section, or spontaneous abortion.

If urine sample is too dilute (ie: low specific gravity) it may not contain a representative level of HCG. If pregnancy is still suspected, wait 48 hours.

Collect a urine sample first thing in the morning and retest it with another strip test. As is true with any diagnostic procedure, the user should evaluate data obtained by the use of this kit with consideration to other clinical information and consult their doctor for the final diagnosis of pregnancy before making any decision of medical relevance.

There are some situations which will effect the test result: some drugs; some diseases, e.g. myoma of uterus; the menopausal population and some hormones, such as the raising concentration of LH, FSH and TSH; etc. Consult the doctors for final confirmation.

8. Transportation and storage

The temperature at transportation should be between -20°C to $+40^{\circ}\text{C}$ with a range of $\pm 2^{\circ}\text{C}$ and Relative Humidity 80% with a range of $\pm 5\%$. Devices packaged by manufacturer may be transported by all modes of transport in covered vehicles in accordance with the rules of transportation of goods for each type of transport.

Store in the original manufacturer's packaging for the labeled shelf life at $+4^{\circ}\text{C}$ to $+30^{\circ}\text{C}$, in an enclosed dry space to protect from rain, heat, direct sunlight.

Do not freeze.

Environment protection requirements

With proper usage, storage and disposal Pregnancy Test does not adversely impact the environment.



9. Disposal instructions

Pregnancy Test should be used in accordance with the instruction for use. The product should be disposed after a single use by a user.

In household use the product should be disposed as household waste. Do not put it down the drain.

In all other cases the Pregnancy Test should be disposed in accordance with administrative and/or local government policy or regulation.

Labeling symbols

	For in vitro diagnostic use only		See instruction for use
	Store between 4°C – 30°C		Do not reuse
	Batch number		Expiry date
	Keep away from sunlight		Tests per kit
	Keep dry		Manufacturer
	Date of manufacturing		

10. The form of claims (complaints)

Information about the claim:

The person/organization that filed the claim		Date of filing the claim	
--	--	--------------------------	--

Sales region		Phone number of the person who filed the claim	
Name of the medical device		Specifications	
Production batch number		LOT number of the medical device	
Number of claims		Total number of transferred products (LOTs)	
Return of the medical device	☐ Yes, time _____ ☐ No	Recipient (seller)	

Method of filing a claim: ☐ Verbally; ☐ By phone; ☐ In writing form/by fax (to attach); ☐ By email (to attach); ☐ Other _____;

Type of claim: ☐ Appearance; ☐ Size; ☐ Functional characteristics; ☐ Reliability; ☐ Reliability; ☐ Labelling; ☐ Other _____;

Is the medical device within the expiration date?: ☐ Yes ☐ No _____;

Was the medical device in use?: ☐ Yes ☐ No _____.

Has the reported adverse event resulted in serious injury or death? ☐ Yes ☐ No

Is the reported undesirable phenomenon related to a malfunction of the equipment, and can a malfunction of such equipment lead to serious injury or death? ☐ Yes ☐ No

Detailed information about the claim:

Signature/Date: _____

Description of the assessment:

- 1) Has the reported adverse event resulted in serious injury or death? ☐ Yes, ☐ No
- 2) Is the reported undesirable phenomenon related to a malfunction of the device, and also can a malfunction of such a device lead to serious injury or death? ☐ Yes, ☐ No
- 3) Can the nature of the reported adverse event lead to serious injury or death? ☐ Yes, ☐ No

Note: when answering "yes" to any of the above questions, it is necessary to notify the competent authorities of the country or region in which the reported undesirable phenomenon occurred;

Evaluation result:

1. Should an undesirable phenomenon be reported to the State in which it occurred

- ☐ Yes, report to ☐ FDA;
☐ relevant authorities in the EU;
☐ the relevant state organization;
☐ the relevant authority of another country or region.
☐ No, the incident does not belong to the undesirable phenomena that need to be reported.

2. Has there been an investigation?

- ☐ Yes.
☐ No (specify the reason) _____

Assessment conducted/date:

Approved/Date:

Investigation (including product characteristics, labeling and packaging):

Archived sample of the medical device:

Information about production control:

Return of goods:

Other:

Conclusion of the investigation:

Signature/Date:

Measures taken:

Are corrective and preventive actions required::

- ☐ Yes
☐ No (specify the reason):

Approved/Date:

Response (with appendices, if available):

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NOTARIAL CERTIFICATE

(2024) STNTZ. Zi, No. 4387

Applicant: Nantong Egens Biotechnology Co., Ltd.,
Unified Social Credit Code: 91320600718557959C,
address: Block A, Building 15, No. 1692, Xinghu Avenue,
Nantong Development Zone.

Legal Representative: Ou Weijun, male, born on
September 16, 1967, Citizen's ID Card No.
320106196709162011.

Authorized Agent: Yu Dan, female, born on July 19,
1979, Citizen's ID Card No. 320602197907191549.

Issue under notarization: seal

This is to certify that the seal of Nantong Egens
Biotechnology Co., Ltd., affixed on the foregoing English
document, is found to be authentic.

Notary: Qin Xiaowen (signature)

Nantong Nantong Notary Public Office,

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NOTARIAL CERTIFICATE

(2024) STNTZ. Zi, No. 4388

Applicant: Nantong Egens Biotechnology Co., Ltd.,
Unified Social Credit Code: 91320600718557959C,
address: Block A, Building 15, No. 1692, Xinghu Avenue,
Nantong Development Zone.

Legal Representative: Ou Weijun, male, born on
September 16, 1967, Citizen's ID Card No.
320106196709162011.

Authorized Agent: Yu Dan, female, born on July 19,
1979, Citizen's ID Card No. 320602197907191549.

Issue under notarization: English translation conforms
to the original document in Chinese

This is to certify that the English translation of the
foregoing (2024) STNTZ. Zi, No. 4387 Notarial Certificate
conforms to the original document in Chinese.

Notary: Qin Xiaowen (signature)
Nantong Nantong Notary Public Office,
Jiangsu Province

The People's Republic of China (seal)
April 2, 2024

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Nantong Notary Public Office, Jiangsu Province, Китай

Nantong Egens Biotechnology Co., Ltd.

Адрес: Building 15, Building 12 (west), No.1692 Xinghu Avenue,
Nantong Economy & Technology Developments Zone,
Jiangsu Province, Китай
Тел.: 86-0513-853-27458, Факс: 86-0513-85327120

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Nantong EGENS Biotechnology Co., Ltd.

Building 15, Building 12 (west), No.1692 Xinghu Avenue,
Nantong Economy & Technology Development Zone,
Jiangsu Province, Китай

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