BEYOĞLU BEŞİNCİ NOTERLİĞİ
Meclisi Mebusan Cad. No: 33
Salıpazarı-Beyoğlu/İSTANBUL
Tel: (0212) 251 69 84-85 Fax: 249 24 00
beyoglu5noter@gmail.com

№01143

İlgili Makama

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	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

INSTRUCTIONS FOR USE

INVADER PTCA Balloon Dilatation Catheter

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

1. General information

1.1 Name

«INVADER™ PTCA Balloon Dilatation Catheter

In versions:

1. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,25x6 mm;
2. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,25x10 mm;
3. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,25x12 mm;
4. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,25x15 mm;
5. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,25x 20 mm;
6. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,50x6 mm;
7. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,50x10 mm;
8. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,50x12 mm;
9. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,50x15 mm;
10. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 1,50x20 mm;
11. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,00x6 mm;
12. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,00x10 mm;
13. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,00x12 mm;
14. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,00x15 mm;
15. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,00x20 mm;
16. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,00x25 mm;
17. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,25x6 mm;
18. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,25x10 mm;
19. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,25x12 mm;
20. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,25x15 mm;
21. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,25x20 mm;

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

22. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,25x25 mm;
23. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x6 mm;
24. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x10 mm;
25. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x12 mm;
26. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x15 mm;
27. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x20 mm;
28. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x25 mm;
29. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x30 mm;
30. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,50x34 mm;
31. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x6 mm;
32. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x10 mm;
33. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x12 mm;
34. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x15 mm;
35. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x20 mm;
36. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x25 mm;
37. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x30 mm;
38. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 2,75x34 mm;
39. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x6 mm;
40. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x10 mm;
41. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x12 mm;
42. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x15 mm;
43. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x20 mm;
44. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x25 mm;
45. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x30 mm;

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

46. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,00x34 mm;
47. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,25x10 mm;
48. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,25x12 mm;
49. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,25x15 mm;
50. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,25x20 mm;
51. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,25x25 mm;
52. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,25x30 mm;
53. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,25x34 mm;
54. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,50x10 mm;
55. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,50x12 mm;
56. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,50x15 mm;
57. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,50x20 mm;
58. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,50x25 mm;
59. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,50x30 mm;
60. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,50x34 mm;
61. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,75x10 mm;
62. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,75x12 mm;
63. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,75x15 mm;
64. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,75x20 mm;
65. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,75x25 mm;
66. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,75x30 mm;
67. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 3,75x34 mm;
68. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,00x10 mm;
69. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,00x12 mm;

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

70. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,00x15 mm;
71. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,00x20 mm;
72. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,00x25 mm;
73. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,00x30 mm;
74. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,00x34 mm;
75. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,50x10 mm;
76. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,50x12 mm;
77. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,50x15 mm;
78. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,50x20 mm;
79. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,50x25 mm;
80. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,50x30 mm;
81. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 4,50x34 mm;
82. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 5,00x10 mm;
83. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 5,00x12 mm;
84. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 5,00x15 mm;
85. INVADER PTCA Balloon Dilatation Catheter with working length of 144 cm and balloon (diameter x length): 5,00x20 mm».

1.2 Device description

INVADER PTCA Balloon Dilatation Catheter (hereinafter referred to as the catheter, PTCA catheter, device) is a sterile, single-use, monorail, semi-compliant rapid exchange (Rx) dilatation catheter for dilatation (predilatation) of target vessels at lesion sites during balloon angioplasty and / or coronary artery stenting operations.

The balloon catheter consists of the proximal tube and distal tube with the balloon as a working element.

The proximal tube is a single-lumen metal tube with outside polytetrafluoroethylene (PTFE) coating. There are two marker bands on the proximal tube enabling the operator to visually control the depth of balloon insertion relative to the puncture site (point of entry into the artery) for the brachial (92 cm) or femoral (102 cm) artery, respectively. One end of the proximal tube contains a hub (made in the form of a Luer port) connecting to the device for inflating/deflating

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

the catheter balloon (indeflator) and pressure control in the balloon. The second end of the proximal tube is rigidly connected to the distal tube of the balloon catheter.

The distal tube of the catheter has two lumens, one for inflation / deflation of the catheter working element (balloon), and one for the guidewire, which serves as a guide for putting the catheter into the target vessel. The lumen for the guidewire begins at the soft atraumatic tip and ends after 27 cm with an exit to the distal tube outer surface. The distal part of the catheter, except for the balloon and soft atraumatic tip, has a hydrophilic coating. The distal part of the catheter contains two radiopaque marker bands labeling the beginning and end of the balloon working length, which allows the operator to accurately position the balloon at the target vessel lesion during the catheter use.

The balloon achieves its nominal dimensions (diameter) when filled the contrast fluid through the Luer port at the proximal catheter end. Each device having a nominal balloon diameter and length is accompanied with a table of the balloon diameter compliance depending on the pressure created in the balloon by a pressure delivery device. The nominal pressure (NP) is a pressure at which the balloon achieves its nominal diameter. The rated burst pressure (RBP) is a pressure at which 99.95% of the balloons will remain intact. The range of pressure values from nominal (NP) to rated burst pressure (RBP) is called the operating pressure range of the balloon.

The catheter is compatible with 0.356 ± 0.025 mm (0.014" inch) interventional guidewires and guiding catheters with an outside diameter of at least 5F.

The effective catheter length is 144 cm. The range of catheters includes balloons from 1.25 mm to 5.0 mm in diameter and from 6 mm to 34 mm in length (see Table 1).

Table 1. Catheter balloon size matrix

		Balloon Length (mm)							
		6	10	12	15	20	25	30	34
Balloon Diameter (mm)	1.25	✓	✓	✓	✓	✓	—	—	—
	1.50	✓	✓	✓	✓	✓	—	—	—
	2.00	✓	✓	✓	✓	✓	✓	—	—
	2.25	✓	✓	✓	✓	✓	✓	—	—
	2.50	✓	✓	✓	✓	✓	✓	✓	✓
	2.75	✓	✓	✓	✓	✓	✓	✓	✓
	3.00	✓	✓	✓	✓	✓	✓	✓	✓
	3.25	—	✓	✓	✓	✓	✓	✓	✓
	3.50	—	✓	✓	✓	✓	✓	✓	✓

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

	3.75	—	✓	✓	✓	✓	✓	✓	✓
	4.00	—	✓	✓	✓	✓	✓	✓	✓
	4.50	—	✓	✓	✓	✓	✓	✓	✓
	5.0	—	✓	✓	✓	✓	—	—	—

Note: The specified balloon diameter is outer diameter, and the balloon length is operating length, both parameters apply to an inflated balloon at nominal pressure.

The catheter is a sterile single-use medical device (no maintenance required).

The catheters are sterilized with ethylene oxide gas.

The catheters do not contain materials of animal or human origin.

The catheters do not contain drugs.

The catheter localization in the patient's body may be done by means of fluoroscopy.

Radiation, electromagnetic and magnetic fields do not affect the device.

1.3 Information about the manufacturer and production site

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1.4 Device classification

According to Directive MDD 93/42/EEC of June 14, 1993. Annex IX, Rule 6, the balloon dilatation catheter is an invasive, sterile medical device used in direct contact with the central circulatory system. Therefore, it is classified as a medical device of Class III risk. The device is CE marked.

1.5 Type of body contact

The catheter is a device having a short-term contact with the circulating blood.

2. Intended use, operating places and potential consumers

Intended use: For use in percutaneous transluminal coronary angioplasty (PTCA) to dilate a stenotic coronary artery by controlled balloon inflation at its distal end.

Operating places: Medical treatment and preventive care institutions.

Potential consumers: Physicians in medical treatment and preventive care institutions.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

3. Operating principle

The device is a percutaneous transluminal coronary angioplasty (PTCA) catheter using hydraulic pressure to dilate stenotic lesions and facilitate expansion of the stent after delivery. After flushing the lumen of the guide and filling the inflation/deflation lumen of the balloon with contrast agent solution, the catheter is inserted over the 0.356 ± 0.025 mm (0.014" inch) PTCA guide through the guiding catheter and placed across the lesion. Radiopaque markers help to correctly position the balloon across the lesion. In accordance with standard PTCA procedures, the balloon is inflated with a 50% contrast medium and 50% heparinized saline solution using an inflator. If necessary, depending on the fluoroscopy data, inflation of the balloon is repeated. Care must be taken not to exceed the balloon rated burst pressure.

After the procedure is completed, the balloon is deflated and the PTCA balloon catheter is removed.

	INVADER PTCA Balloon Dilation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

4. Indications for use

- The catheters are intended for dilating stenotic atherosclerotic lesions in coronary arteries or vascular bypass grafts to increase myocardial perfusion;
- Patients should be eligible for coronary bypass surgery;
- The catheters are indicated for patients with single-artery non-calcified atherosclerotic lesions that can be dilated using a PTCA catheter.
- This procedure can also be indicated for certain patients who have multi-artery disease, and for patients who have undergone aorta-coronary bypass surgery but still have:
 - recurrent symptoms;
 - progressive coronary artery disease;
 - stenosis or obstruction in bypass grafts.

5. Contraindications

The catheter use is contraindicated for patients:

- not eligible for coronary bypass surgery
- with fully obstructed coronary arteries
- with diffuse lesions
- with severe stenosis of the left main coronary artery.

6. Adverse events

- Dissection of the coronary artery
- Tearing, perforation of, or damage to the coronary artery
- Complete obstruction of the coronary artery or bypass graft
- Thrombosis of the coronary artery
- Unstable angina
- Acute myocardial infarction
- Restenosis of the dilated artery
- Spasm of the coronary artery
- Arrhythmias, including ventricular fibrillation
- Hemorrhage and hematoma
- Drug reactions, allergic reaction to contrast agents used during the procedure;
- Hypertension or hypotension
- Infection
- Arteriovenous fistula
- Embolism
- Death
- Urgent coronary artery bypass graft surgery
- Renal failure

7. Personnel qualification

The catheters should only be used by experienced physicians trained in PTCA operations.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

8. Warnings / Precautions

- The catheters are designed for single use only, reuse is prohibited.
- Resterilization is prohibited.
- Do not use a catheter that has been exposed to the sunshine or stored in non-recommended conditions.
- Do not use a catheter after the expiry date printed on its packaging.
- Use diluted contrast medium only.
- Do not use air or any other gas to inflate the balloon.
- Check the packaging for any damage. Do not use the catheter if the integrity of the sterilization packaging is violated.
- When using the catheter, avoid its coming into contact with non-sterile surfaces.
- Do not exceed the rated burst pressure indicated in the Instructions for Use when inflating the balloon.
- The inflated balloon diameter at the points limiting the stenotic area from the proximal and distal side should not exceed the determined value.
- The catheters should only be used by experienced physicians trained in PTCA operations.
- Appropriate anticoagulation and vasodilatation therapy should be given prior to catheterization.
- The PTCA operation should only be performed in medical treatment and preventive care institutions having necessary equipment for emergency coronary bypass surgery in case of severe complications.
- Do not tighten the hemostatic adapter in the Y-connector, as this can compress the shaft, thus impeding the inflation and deflation of the balloon.
- All procedures to be done after the catheter introducing into the body should be performed under fluoroscopic control. Do not move the catheter until the balloon is completely deflated. If any resistance is encountered during the procedure, stop and try to identify the cause and then gradually advance the balloon catheter avoiding kinking of its proximal part. If you fail to identify the cause, remove the entire system.
- Do not use Ethiodol or Lipiodol contrast media.
- Do not treat the catheter with organic solvents (such as alcohol, etc.).

9. Delivery set

The device delivery set includes:

- Catheter in the sterilization packaging - 1 pc.;
- Instructions for Use - 1 pc.

10. Instructions for Use

10.1 Selecting the right catheter size

- When choosing a catheter, special attention should be paid to the balloon diameter. The diameter of the catheter balloon inflated at nominal pressure should not exceed the coronary artery diameter of the at the points adjacent to the stenotic area on the proximal and distal side. The catheter shall be chosen is by a qualified physician with specialized medical education based on the indications, type and features of the forthcoming manipulations with the device.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

10.2 Preparation

- Prepare an inflation solution of a 1:1 mixture of contrast material and sterile saline solution.
- Fill a 20 cc syringe (not included in the delivery set) with approximately 10 cc of the saline solution.
- Place a needle (not included in the delivery set) on the syringe, and then carefully insert the needle through the distal tip of the catheter and flush the guidewire lumen.
- Attach the stopcock (not included in the delivery set) to the catheter hub.
- Attach the syringe to the stopcock and then turn the syringe vertically so that its plunger is at the top. Pull the plunger and purge the balloon of air.
- Continue suction with the syringe for 15-20 seconds, making sure that no air bubbles appear in the diluted contrast material.
- Carefully stop suction.
- Repeat the process if required.
- Exert negative pressure and then switch off the stopcock.
- Make sure that no air enters the system, and attach the inflation device to the catheter. Switch on the stopcock, inflate the balloon to nominal pressure, and then check the integrity of the catheter.
- Exert negative pressure and then switch off the stopcock.

10.3 Insertion of the catheter

- Using the hemostatic valve (not included in the delivery set), attach the Y-connector to the guiding catheter already inserted in the artery.
- Pass the guidewire (not included in the delivery set) through the guiding catheter (not included in the delivery set), advance and then place it as required.
- Gently tighten the hemostatic valve of the Y-connector on to the guidewire so as to prevent unintended movement of the guidewire.
- Insert the rear tip of the guidewire in the distal tip of the catheter and attach it to the catheter. The rear tip of the guidewire will come out at approximately 27 cm proximal to the catheter's distal tip.

WARNING: Kissing balloon: Bench tests have shown that two catheters with 4 mm diameter can be inserted in a 6F (or larger) guide catheter.

- Loosen the hemostatic valve of the Y-connector and advance the catheter towards the distal part of the guiding catheter.

WARNING: Always advance the PTCA catheter when it is deflated and placed over the guidewire.

- Do not tighten the hemostatic valve of the Y-connector too much, as this may prevent the passage of the contrast material through the balloon inflation lumen, thus resulting in the prolongation of the balloon deflation/inflation time.
- Using standard angioplasty techniques, continue the procedure under fluoroscopy. Radiopaque markers can help to better position the balloon in the stenosis.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

WARNING: It is recommended that the guidewire and/or PTCA catheter remain within the lesion until the dilatation is completed. The balloon deflation/inflation time may vary depending on the viscosity of the contrast medium.

10.4 Removal of the catheter

- Loosen the hemostatic valve of the Y-connector.
- Hold the hemostatic valve and guidewire with one hand, and the catheter shaft with the other hand.
- Whilst keeping the guidewire immobile in order to maintain its position in the coronary artery, remove the PTCA catheter from the guiding catheter.

WARNING: During the removal, check the position of the guidewire by means of fluoroscopy.

- Switch on the hemostatic valve of the Y-connector.
- If necessary, prepare another PTCA catheter and repeat the dilatation process.

11. Device specifications

The device specifications are presented in Table 2.

Table 3. Device specifications

#	Specification	Value
Device specifications		
01	Effective length of the catheter (L)* * detailed dimensioning specifications of the device are given in clause 12	1440 +0/-150 mm
02	Designation of the cylinder length (for indication in the documentation and on the label) * The dimensions of the main parts of the catheter balloon are shown in Table A2	6, 10, 12, 15, 20, 25, 30, 34 mm $\pm$ 10 %
03	Designation of the cylinder diameter (for indication in the documentation and on the label) * The effective diameter of the cylinder corresponding to the pressure created inside is shown in Table 4.	1.25, 1.50, 2.00, 2.25, 2.50, 2.75, 3.00, 3.25, 3.50, 3.75, 4.00, 4.50, 5.00 mm $\pm$ 6 %
04	Effective length of the balloon at nominal pressure (NP)	6 mm: $5 \leq L \leq 7$ mm 10 mm: $9 \leq L \leq 11$ mm 12 mm: $11 \leq L \leq 13$ mm 15 mm: $14 \leq L \leq 16$ mm 20 mm: $19 \leq L \leq 21$ mm 25 mm: $24 \leq L \leq 26$ mm 30 mm: $29 \leq L \leq 31$ mm 34 mm: $33 \leq L \leq 35$ mm

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

#	Specification		Value
05	Outer diameter of the distal tube		2,7F* / 0,9 mm $\pm$ 0,25 mm
06	Outer diameter of the proximal tube		1,78F* / 0,6 mm $\pm$ 0,25 mm
07	Outer diameter of the distal tip		0,017" / 0,43 mm $\pm$ 0,25 mm
08	Atraumatic tip shape		must be rounded, conical
09	Minimum internal diameter of a compatible conductive catheter		5F* / 1,422 mm $\pm$ 1%
10	Compatible guidewire diameter		0.356 $\pm$ 0.025 mm (0.014" inch)
11	Distance from the middle of the first radiopaque balloon marker to the middle of the second radiopaque balloon marker	Balloon length, mm	Distance between the marker bands, mm
		6 mm	6 $\pm$ 0,1 mm
		10 mm	10 $\pm$ 0,1 mm
		12 mm	12 $\pm$ 0,1 mm
		15 mm	15 $\pm$ 0,1 mm
		20 mm	20 $\pm$ 0,1 mm
		25 mm	25 $\pm$ 0,1 mm
		30 mm	30 $\pm$ 0,1 mm
		34 mm	34 $\pm$ 0,1 mm
12	Nominal pressure of the balloon (NP)		8 atm (811 kPa)
13	Rated burst pressure (RBP)		16 atm (1621 kPa)
14	Inflation time of the balloon		$\leq$ 15 s
15	Deflation time of the balloon		$\leq$ 15 s
16	Burst strength of the balloon		The balloon must withstand a minimum of 10 rated burst pressure cycles
17	Tensile force of tested parts of the catheter: catheter tip detachment		$\geq$ 3.0 N
18	Tensile force of tested parts of the catheter: break of the connection between the proximal and distal tube of the catheter**		$\geq$ 5.0 N
19	Tensile force of tested parts of the catheter: break of the connection of the proximal tube and hub**		$\geq$ 15.0 N
20	Appearance of the catheter		The catheter should be free of kinks, cuts and contaminants. The thermal welding sites should be even and have smooth transitions. The balloon should have no damage. The Luer port must be free of burrs, air bubbles, damage, drips, cracks, etc.
21	Stylet	Appearance	Should be free of burrs, damage, cracks, etc.
		Diameter	0,36 $\pm$ 0,07 mm
		Length	295 $\pm$ 3 mm
22	Clip	Appearance	Should be free of burrs, air bubbles, damage, drips, cracks, etc.
		Number of holes	3 pieces.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

#	Specification		Value
		Diameter of holes	$3,8 \pm 0,1$ mm
		Length	$19 \pm 1,5$ mm
		Width	$8,5 \pm 0,15$ mm
		Height	$5,3 \pm 0,15$ mm
23	Strain relief	Appearance	Must be free of burrs, air bubbles, damage, drips, cracks, etc.
		Outer diameter	$4,5 \pm 0,05$ mm
		Inner diameter	$3,2 \pm 0,07$ mm
		Length	145 ± 5 cm

* Value according to the French scale (Charrière scale)
** Connection points of the catheter parts are shown in Figure 3

Effective balloon diameters corresponding to the pressure created inside are presented in Table 3.

Table 3. Balloon compliance table

Pressure (atm)	Pressure (kPa)	Diameter of the balloon (mm)												
		1.25	1.50	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.50	5.0
5	507	1,17	1,37	1,88	2,03	2,40	2,63	2,86	3,02	3,34	3,47	3,84	4,29	4,64
6	608	1,20	1,40	1,92	2,09	2,46	2,68	2,93	3,09	3,41	3,57	3,91	4,38	4,79
7	709	1,22	1,45	1,97	2,15	2,51	2,73	2,98	3,16	3,46	3,66	3,96	4,45	4,91
8*	811	1,24	1,48	2,01	2,21	2,56	2,76	3,02	3,22	3,51	3,74	4,01	4,51	5,02
9	912	1,27	1,52	2,04	2,26	2,60	2,80	3,06	3,26	3,55	3,80	4,05	4,55	5,11
10	1013	1,30	1,55	2,07	2,30	2,63	2,83	3,09	3,30	3,58	3,85	4,08	4,59	5,18
11	1115	1,33	1,58	2,10	2,33	2,66	2,85	3,12	3,33	3,60	3,90	4,11	4,62	5,24
12	1216	1,36	1,61	2,12	2,37	2,68	2,87	3,14	3,36	3,63	3,93	4,14	4,65	5,30
13	1317	1,40	1,64	2,15	2,40	2,71	2,90	3,16	3,40	3,65	3,96	4,17	4,68	5,34
14	1419	1,43	1,67	2,18	2,43	2,74	2,92	3,19	3,42	3,68	3,99	4,19	4,70	5,39
15	1520	1,47	1,71	2,21	2,47	2,77	2,94	3,22	3,46	3,70	4,02	4,22	4,74	5,42
16**	1621	1,51	1,75	2,24	2,50	2,80	2,97	3,24	3,49	3,73	4,05	4,24	4,77	5,47
17	1723	1,55	1,79	2,28	2,54	2,84	3,00	3,27	3,52	3,76	4,08	4,27	4,80	5,51
18	1824	1,59	1,83	2,31	2,57	2,87	3,03	3,29	3,56	3,78	4,11	4,30	4,83	5,54
19	1925	1,65	1,88	2,35	2,61	2,90	3,06	3,32	3,59	3,81	4,13	4,33	4,86	5,57
20	2027	1,73	1,92	2,39	2,64	2,94	3,09	3,35	3,62	3,84	4,15	4,36	4,89	5,61

*Nominal pressure of the balloon (NP)

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

**** Rated burst pressure of the balloon (RBP)**

Corrosion resistance

Metal parts of catheters are not corrodible. After testing catheters according to the method specified in Annex A of ISO 10555-1 standard, metal parts of the catheter shall be free from corrosion.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

12. Schematic representation of the device


Figure 1. Schematic representation of the catheter

Table A1. Main parts of the catheter with their sizes

#	Ref	Description	Value
1	A	Effective length of the catheter	1440 +0/-150 mm
2	B	Distance from the atraumatic tip to the Rx port	270 ± 5 mm
3	C	Distance from the atraumatic tip to the first marker band of the proximal tube	920 ± 60 mm
4	D	Distance from the atraumatic tip to the second marker band of the proximal tube	1020 ± 60 mm
5	E	Outer diameter of the distal tip	0,017" / 0,43 mm ± 0,25 mm
6	F	Outer diameter of the distal tube	2,7F / 0,9 mm ± 0,25 mm
7	G	Outer diameter of the proximal tube	1,78F / 0,6 mm ± 0,25 mm
8	H	Catheter balloon	* Detailed dimensioning specifications of the balloon are shown on Figure 2
9	I	Rx port	For information only
10	J	Luer port	Inner diameter is 3,9 ± 0,5 mm

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

			Outer diameter is $7,8 \pm 0,5$ mm Cone 6%
11	K	Length of the proximal tube marker band	10 ± 1 mm


Figure 2. Schematic representation of the catheter balloon

Table A2. Main parts of the catheter balloon with dimensions

№	Letter symbol	Description	Value	
1	A	Atraumatic tip length	$2,75 \pm 0,25$ mm	
2	B	Distance from the end of the balloon cones to the middle of radiopaque markers	Diameter of the balloon, mm	Value, mm
			1.25	$1,25 \pm 0,1$
			1.50	$1,50 \pm 0,1$
			2.00	$2,00 \pm 0,1$

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

№	Letter symbol	Description	Value	
			2.25	$2.25 \pm 0,1$
			2.50	$2.50 \pm 0,1$
			2.75	$2.75 \pm 0,1$
			3.00	$3.00 \pm 0,1$
			3.25	$3.25 \pm 0,1$
			3.50	$3.50 \pm 0,1$
			3.75	$3.75 \pm 0,1$
			4.00	$4.00 \pm 0,1$
			4.50	$4.50 \pm 0,1$
			5.00	$5.00 \pm 0,1$
3	C	Distance from the middle of the first marker to the middle of the second radiopaque balloon marker	* dimensioning specifications see in item 11 of Table 3	
4	D	Length of the radiopaque marker band	$1,016 \pm 0,102$ mm	


Figure 3. Connection points of the catheter parts

A - junction of the atraumatic tip, B - junction of the proximal part of the balloon, C - junction of the catheter proximal and distal tube, D - junction of the proximal tube and connecting module






	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

13. Labeling

Labeling is performed by printing the basic data about the device on the labels attached to the device sterilization, packaging, and transportation boxes.

13.1 Labeling of the sterilization packaging

The labeling of the sterilization packaging should contain:

Description	Symbol
Company Logo	-
Product Name (including the Russian version)	-
Schematic Representation of the Product with Dimensions in mm	-
Balloon Compliance Table	-
Designation of the cylinder diameter (mm)	
Cylinder length designation (mm)	
Catalogue Number	
Batch Code	
Date of Manufacture	
Use by	
Minimum Guiding Catheter Size	
Rated Burst Pressure	
Compatible Guidewire Diameter	
Sterilization with ethylene oxide	
Contents	
Storage temperature range +10°C - +40°C	
Avoid exposure to sunlight.	
Avoid exposure to moisture	
Caution.	
Prescription only-device restricted to use by or on order of a physician.	

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

Description	Symbol
Consult instructions for use.	
Do not reuse.	
Do not use if package is damaged.	
Do not resterilize	
Barcode	-
Identification labels	-
Manufacturer	
Authorised Representative in the European Community	
CE Mark	
Company Website	-
Date and number of the Registration Certificate (specified on an additional label)	-
Information about the Authorized Representative in Russia (specified on an additional label)	-
Label abbreviations: atm - atmosphere mm - millimeter NP - Nominal Pressure RBO - Rated Burst Pressure "- inch	

There must be two identification stickers on the sterilization packaging label. Identification stickers should bear the following information: brand name of the catheter, balloon diameter (mm), balloon length (mm), Catalogue Number, Batch Code, and Use by symbols, as well as barcode. After percutaneous transluminal coronary angioplasty, the physician shall put one identification sticker in the patient's medical record, and another one in the operation log.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024



Figure 4. Labeling of the sterilization packaging

13.2 Labeling of the packaging box

Labeling of the packaging box must contain:

fvgtDescription	Symbol
Company Logo	-
Product Name (including the Russian version)	-
Schematic Representation of the Product with Dimensions in mm	-
Balloon Compliance Table	-
Designation of the cylinder diameter (mm)	
Cylinder length designation (mm)	

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

fvgtDescription	Symbol
Catalogue Number	
Batch Code	
Date of Manufacture	
Use by	
Minimum Guiding Catheter Size	
Rated Burst Pressure	
Compatible Guidewire Diameter	
Sterilization with ethylene oxide	
Contents	
Storage temperature range +10°C - +40°C	
Avoid exposure to sunlight	
Avoid exposure to moisture	
Caution.	
Prescription only-device restricted to use by or on order of a physician.	
Consult instructions for use.	
Do not reuse.	
Do not use if package is damaged.	
Do not resterilize	
Barcode	-
Identification labels	-
Manufacturer	
Authorised Representative in the European Community	
CE Mark	
Company Website	-

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

fvgtDescription	Symbol
Date and number of the Registration Certificate (specified on an additional label)	-
Information about the Authorized Representative in Russia (specified on an additional label)	-
Label abbreviations: atm - atmosphere mm - millimeter NP - Nominal Pressure RBO - Rated Burst Pressure "- inch	

There must be two identification stickers on the packaging box label identical to those on the sterilization packaging, which stickers are duplicate. Those identification stickers should bear the following information: brand name of the catheter, balloon diameter (mm), balloon length (mm), Catalogue Number, Batch Code, and Use by symbols, as well as barcode. After percutaneous transluminal coronary angioplasty, the physician shall put one identification sticker in the patient's medical record, and another one in the operation log.

The labeling of the packaging box for devices delivery to the Russian Federation should be provided with an additional sticker.

Layout of the additional sticker to be put on the packaging box for devices delivery to the Russian Federation:

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Content translation:

INVADER PTCA Balloon Dilatation Catheter, 144 cm working length with balloon (diameter x length): (size indicated on the packaging) Article: indicated on the packaging

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

13.3 Labeling of the transportation packaging

The handling signs "Top", "Keep away from moisture", "Fragile, careful" should be applied to the transport package.

Figure 10. Labeling of the transportation box

14.4 Luer port labeling

The Luer port must bear the following information: brand name of the catheter, balloon diameter (mm), balloon length (mm).

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

14. National and international standards

The list of standards used in the medical device development and production is given in Table 5.

Table 4. List of applicable standards

No.	Standard Number	Standard Title
01	EN 556-1:2001 / AC:2006	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
02	EN 1041:2008 A1:2013	Information supplied by the manufacturer of medical devices
03	EN 1422:1997+A1:2009	Sterilizers for medical purposes - Ethylene oxide sterilizers - Requirements and test methods
04	EN 62366:2008	Medical devices - Application of usability engineering to medical devices
05	EN 20594-1:1993	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Part 1: General requirements (ISO 594-1:1986)
06	EN ISO 10555-1:2009	Sterile, single-use intravascular catheters - Part 1: General requirements (ISO 10555-1:1995, including Amd 1:1999 and Amd 2:2004)
07	EN ISO 10993-1:2009 / AC:2010	Biological evaluation of medical devices - Part 1: Evaluation and testing
08	EN ISO 10993-4:2009	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood
09	EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
10	EN ISO 10993-7:2008/ AC:2009	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals
11	EN ISO 10993-11: 2009	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
12	EN ISO 10993-12:2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
13	EN ISO 11135-1:2007	Sterilization of health care devices - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024


No.	Standard Number	Standard Title
		devices
14	EN ISO 11607 - 1:2009	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
15	EN ISO 11607 - 2 :2006	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
16	EN ISO 13485 : 2016	Medical devices - Quality management systems - Requirements for regulatory purposes
17	EN ISO 14971: 2012	Medical devices - Application of risk management to medical devices
18	EN ISO 25539-2:2009	Cardiovascular implants - Endovascular devices - Part 2: Vascular stents
19	ISO 80369-7:2016	Small-bore connectors for liquids and gases in healthcare applications — Part 7: Connectors for intravascular or hypodermic applications
20	ISO 6009: 2016	Hypodermic needles for single-use - Colour coding for identification
21	ISO 10555-1-2013	Intravascular catheters - Sterile and single-use catheters - Part 1: General requirements
22	ISO 10555-4-2013	Intravascular catheters - Sterile and single-use catheters - Part 4: Balloon dilatation catheters
23	ISO 10993-10:2010	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
24	ISO 11135: 2014	Sterilization of health care devices - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices
25	ISO 13485 : 2016	Medical devices - Quality management systems - Requirements for regulatory purposes
26	ISO 14644-1:2015	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness by particle concentration
27	ISO 14644-2:2015	Cleanrooms and associated controlled environments -

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

No.	Standard Number	Standard Title
		Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration
28	ISO 14644-3:2005	Cleanrooms and associated controlled environments - Part 3: Test methods
29	ISO 14644-4:2001	Cleanrooms and associated controlled environments - Part 4: Design, construction and start-up
30	ISO 14698-1:2003	Cleanrooms and associated controlled environments - Biocontamination control - Part 1: General principles and methods
31	ISO 15223-1: 2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
32	BS EN ISO 10993-10: 2013	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
33	BS EN ISO 25539-2:2012	Cardiovascular implants - Endovascular devices - Part 2: Vascular stents
34	CAN/CSA-ISO 13485:2003	Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes
35	ASTM D 618- 13	Standard Practice for Conditioning Plastics for Testing
36	ASTM D 638-14	Standard Test Method for Tensile Properties of Plastics
37	ASTM D 642 -15	Standard Test Method for Determining Compressive Resistance of Shipping Containers, Components, and Unit Loads
38	ASTM D 999 - 08	Standard Test Methods for Vibration Testing of Shipping Containers
39	ASTM D 3078-02 (Reapproved 2008)	Standard Test Method for Determination of Leaks in Flexible Packaging by Bubble Emission
40	ASTM D 4169 – 2016	Standard Practice for Performance Testing of Shipping Containers and Systems (non-monetary information)
41	ASTM D 4332-14	Standard Practice for Conditioning Containers, Packagings, or Packaging Components for Testing

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

No.	Standard Number	Standard Title
42	ASTM D 4728 – 17	Standard Test Method for Random Vibration Testing of Shipping Containers
43	ASTM D 5276 – 98 – Reapproved 2009	Standard Test Method for Drop Test of Loaded Containers by Free Fall (mandatory and non-monetary information)
44	ASTM D 5487-2016	Standard Test Method for Simulated Drop of Loaded Containers by Shock Machines
45	ASTM D 6344 – 04 – Reapproved 2017	Standard Test Method for Concentrated Impacts to Transport Packagings
46	ASTM F 88- F 88 M-15	Standard test method for seal strength of flexible barrier materials ¹
47	ASTM F 640-12	Standard Test Methods for Determining Radiopacity for Medical Use (non-monetary information)
48	ASTM F 1886- F 1886M -09	Standard Test Method for Determining Integrity of Seals for Flexible Packaging by Visual Inspection
49	ASTM F 1980 – 2016	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
50	ASTM F 1984-99 (Reapproved 2013)	Standard Practice for Testing for Whole Complement Activation in Serum by Solid Materials
51	ASTM F 2065-00 (Reapproved 2010)	Standard Practice for Testing for Alternative Pathway Complement Activation in Serum by Solid Materials
52	ASTM F 2096-11	Standard Test Method for Detecting Gross Leaks in Packaging by Internal Pressurization (Bubble Test)
53	ASTM F 2097-16	Standard Guide for Design and Evaluation of Primary Flexible Packaging for Medical Device
54	ANSI ASQ Z1.4:2013	Sampling Procedures and Tables for Inspection by Attributes
55	ANSI ASQ Z1.9 - 2008	Sampling Procedures and Tables for Inspection by Variables for Percent Nonconforming

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

15. Transportation, storage and use conditions

15.1 Transportation conditions

Transportation conditions:

Ambient temperature: - 5°C to + 40°C.

Relative air humidity: no more than 80%.

Atmospheric pressure: 84 to 106.7 kPa

Keep away from moisture and direct sunlight.

15.2 Storage conditions

Store in a cool, dry and dark place in the original packaging.

Storage conditions:

Ambient temperature: + 10°C to + 40°C.

Relative air humidity: no more than 80%.

Atmospheric pressure: 84 to 106.7 kPa

Keep away from moisture and direct sunlight.

15.3 Use conditions

The devices are manufactured for use in a biological environment in the temperature range from 32 °C to 42 °C.

16. Sterilization, cleaning and disinfection

The device is supplied sterile. Sterilization parameters: Gas sterilization method (ethylene oxide).

The device is non-pyrogenic.

The device is non-toxic.

The device is a sterile single-use medical device (cleaning and disinfection is not applicable).

Re-sterilization of the device is prohibited.

Reuse is prohibited.

22. Disposal information

After use, all components of the device must be disposed of in accordance with the guidelines governing disposal of such devices and applicable waste control regulations.

Expired devices, damaged devices or their parts shall be destroyed as Class A medical waste (epidemiologically safe waste similar in composition to municipal solid waste) in accordance with SanPiN 2.1.3684.

Used devices shall be destroyed as Class B medical waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.3684.

	INVADER PTCA Balloon Dilatation Catheter	Rev. 03
	Instructions for Use	Date: January 22, 2024

23. Guaranteed Shelf Life

Guaranteed shelf life of the devices is 3 years from the date of sterilization.

24. Warranty

The manufacturer guarantees that all of its devices at the time of delivery are free from defects due to the material quality and manufacturing process.

The manufacturer emphasizes that this device is for single use only (no maintenance is applicable), and makes no warranties, representations or conditions, either express or implied (including any warranties of merchantability or usability for a particular purpose) relating to the reuse of this device. Moreover, the manufacturer assumes no obligations or liability for incidental or consequential damage that may result from reuse of the device.

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As to questions on the medical device quality within the Russian Federation, please contact the authorized representative of the manufacturer in the Russian Federation.

Send claims to:

Authorized representative of the manufacturer in the Russian Federation:

Limited Liability Company Research and Production Association "Deost" (LLC NPO "Deost")

Legal address of the authorized representative of the manufacturer in the Russian Federation:

Russia, 142290, Moscow Region, Pushchino, Prospect Nauki, building 3, suite 4, office 22, floor 4

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