Serial No.5260-7-23.

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Hereby certify that I was presented with:

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☐ a copy that is not the original and is not a true copy of the original

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☐ computerized information (please specify) from the database of

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I confirm, with my signature and seal, that the attached document marked with the letter/ numberA. is an exact copy of the document or information presented to me.

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«Approved»

Alpha-Bio Tec Ltd.

Regulatory Affairs Project Manager

Shani Aroch

02 July 2023
Data


Signature

Stamp

Shani Aroch
RA Team Leader, PRRC

INSTRUCTION FOR USE FOR MEDICAL DEVICE

«Stomatological instruments for dental implantology with accessories»

Name of medical device

Dental instruments for dental implantology with accessories

1. Surgical drill, variations: 65007 CD2-8; 65008 CD2-10; 65009 CD2-11.5; 65010 CD2-13; 65011 CD2-16; 65012 CD2.8-8; 65013 CD2.8-10; 65014 CD2.8-11.5; 65015 CD2.8-13; 65016 CD2.8-16; 65017 CD3.2-8; 65018 CD3.2-10; 65019 CD3.2-11.5; 65020 CD3.2-13; 65021 CD3.2-16; 65022 CD3.6-8; 65023 CD3.6-10; 65024 CD3.6-11.5; 65025 CD3.6-13; 65026 CD3.6-16; 65027 CD4.1-8; 65028 CD4.1-10; 65029 CD4.1-11.5; 65030 CD4.1-13; 65031 CD4.1-16; 65032 CD4.5-8; 65033 CD4.5-10; 65034 CD4.5-11.5; 65035 CD4.5-13; 65036 CD4.5-16; 65070 CD2.4-8; 65071 CD2.4-10; 65072 CD2.4-11.5; 65073 CD2.4-13; 65074 CD2.4-16.
2. Tissue punch, variations: 65003 TPS; 65004 TPL.
3. Crestal drill, variations: 65005 CDS; 65006 CDL.
4. Pin, variations: 65047 LP; 65048 CPS; 65049 CPL.
5. Milling drill: 65050 MCD1.5.
6. Implant amount, variations: 65037 IMS; 65038 IML; 65055 IMC; 65064 IMSS; 65065 IMSL.
7. Driver, variations: 65061 IMCD; 65062 IMSD; 65063 IMLD; 4057 HTD1.5; 4058 HTD1.5S; 4059 HHS1.5; 4060 HHL1.5; 4168 HT1.5; 4140 G-ITDL2.5; 4141 G-ITDM2.5; 4142 G-ITDS2.5; 4143 GITL2.5/1.25; 4144 GITM2.5/1.25; 4145 GITS2.5/1.25; 4146 MITD2.5-IH; 4147 MITD2.1-CHC.
8. Implant mount screw 65039 IMHS.
9. Implant mount extension 65042 IMX.
10. Handpiece insertion adapter 65044 HIA.
11. Implant amount extractor 65045 IME.
12. Implant mount screw - CHC: 65056 IMCS.
13. Sleeve L/S adapter driver: 65057 SAD.
14. Sleeve for adapter 65058 SLSA.
15. Sleeve, variations: 66012 SLS; 66013 SLL; 66014 SLSE.

Accessories:

1. Box: 65002 DWT, 4699 OBSA-CS.

Product description

Important information!

All instruments are supplied nonsterile and must be cleaned and disinfected before use after packaging removing.

These guidelines apply to all stainless steel or titanium products. Products are referred to implant drivers, transfers and similar instruments, trephines, burs, healing formers, titanium alloy abutments and all other surgical instruments made from stainless steel and titanium alloys.

Alpha-Bio Tec dental drills from the group of instruments for guided surgery include drills with external irrigation in various diameters and lengths. Drills are made of stainless surgical steel. Using them a dental surgeon has better control the process of creating a channel for implantation of the required diameter and depth. The Alpha-Bio Tec Dental Drills in the Guided Surgery Toolbox are available with and without laser marking. The diameter and length of each drill is laser engraved on the shank. Each drill has a built-in stop that prevents the user from drilling deeper than required.

Attention!

1. Depth marks on laser marked drills refer to the length of the drill tip (measured from the tip of the drill to the middle of the indicator line). The length of the drill tip depends on its diameter and must be considered when preparing for an osteotomy.
2. In drills without laser marking, the specified drill length does not include the length of the triangular cutting edge. When determining the depth of drilling, it is necessary to add 0.3-0.8 mm to it, especially when manipulating near anatomical structures.

Intended for use

The product is intended for use in surgical and restorative dental practice.

Indications

Alpha-Bio Tec Guided Surgery drills and other instruments are intended for preparation of the cavity for Alpha-Bio Tec implants placement.

Alpha-Bio Tec Guided Surgery drills must be used with the help of a hand piece. Recommended working speed is from several rpm to 800 rpm (high speed use over 1000 rpm is discouraged).

The use of sharp drills, sufficient irrigation, in-and out drilling motion, short cutting cycles, waiting for the bone to cool, and use of pilot drills in successively increasing sizes are essential.

Contraindications

- Patients with pre-existing known allergy to stainless steel material.

Potential adverse events that may occur include the following:

- infection
- bleeding
- irritation.

Precautions

- Always wear gloves when handling contaminated instruments.
- Eye protection must be worn to prevent particulate matter.
- A surgical mask must be worn to protect against inhalation of the generated dust or aerosol.
- Alpha-Bio Tec Guided Surgery Drills may only be used by licensed dentists who have the proper specialization, skill, and training to work with dental implants and prosthetic elements.

Cautions:

Alpha-Bio Tec GSTK Drills and parts are sold non-sterile and must be cleaned, disinfected and sterilized prior to first clinical use and between subsequent uses in accordance with the Instructions for Reprocessing/Sterilization section of the User Instructions.

- All packaging must be removed prior to processing. The drills should be inspected for damage or wear before each use and any defected drill should be discarded.
- Alpha-Bio Tec instruments including drills may be damaged by alkaline and acidic detergents. Always use a neutral pH detergent solution for cleaning.
- Cleaning, disinfecting and sterilization must only be performed with specific materials, specifically; detergents containing the following substance should be avoided: Oxalic acid and high concentrate Chloride.
- Ensure that the drill is fully seated and gripped in the handpiece collet before use.
- Maintain hand piece in good working order and correctly lubricated.
- Do not exceed the maximum speeds indicated in this User Instructions.
- Avoid excessive drilling speed and/or drilling duration in order to avoid overheating and associated complications.
- Move the drill continuously when in use to avoid localized heating.
- Insert the drills in the disinfectant solution immediately after surgery to prevent encrustation of blood, secretions, etc.
- Do not store wet or damp instruments including drills.
- Do not sterilize, clean or disinfect instruments made with different metals in the same

cycle.

- Steam Sterilization in autoclave is recommended to prevent dulling.
- Alpha-Bio Tec instruments including drills should not be exposed to temperatures more than 140°C.
- The operating surgeon / practitioner in charge must carefully read the cleaning & sterilization instructions prior to use.
- Alpha-Bio Tec is not liable for complications, other negative effects or damages that might occur for reasons such as unsuitable choice of material or handling thereof, unsuitable use or handling of the instruments, asepsis and so on. The operating surgeon is responsible for any such complications or other consequences.

Information about potential consumers of the medical device

Intended user profile:

A dentist who completed a postgraduate course in the field of implantology.

Application area:

Dentistry.

These products should only be used by certified dental surgeons/dental surgeons.

Procedures with the use of these products must be carried out within the framework of a dental office provided for surgical operations and located in professional dentistry.

Method of using of medical device

Principle of operation

Computer Guided Surgery

Alpha-Bio Tec supports dental professionals with a range of innovative solutions for implant applications and surgical template production for computer-guided surgery.

☐ Maximum safety, whereby clinical decisions are made in advance, prior to the surgery, deploying software tools and cutting-edge visual technologies, resulting in less invasive, flapless surgery and guided implant insertion via a template.

Planning, inserting implants and surgical template preparation, dedicated for computer guided implantation (Guided Surgery) in cases of partially edentulous or edentulous, using a safe, fast and easy method.

☐ Maximum precision and matching between planning and execution, ensuring implant placement predictability, prior to performing the surgery, and in many cases, enabling the preparation of temporary rehabilitation surgery based on planning only.

☐ Executing a less invasive, flapless surgery and guided insertion of the implants through the template. Drilling is performed using the designated drills from the GSTK and continues until contacting the template stopper. The GSTK drills are available in various lengths. In cases where various length drills are required, drill until the GSTK drills' built-in stoppers come into contact with the sleeves in the surgical template. Refer to the specifications in the relevant protocol for each implant.

In the case of guided implant placement through the template, the implant should be placed through an appropriate sleeve. A designated guide is screwed to the implant.

After connecting the designated implant guide, the implant is inserted using a ratchet, until the stopper of the designated guide comes in contact with the sleeve in the template. After implant placement, the designated guide can be removed by loosening the screw. After the screw is loosened, use the extractor provided with the GSTK and turn the extractor until the guide disconnects from the implant. After disconnecting all of the guides, the fixation pins can be removed (if used) and the surgical template can also be removed.

Technical characteristics

Product code	Product name	Product description	Specification	Raw materials	Detailed description
1. Surgical drill, variations:					
65007 	CD2-8	Surgical drill D2.0mm – L8mm	Mass: 1.82g Height: 33.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 58.5mm Length of working part: 8mm Diameter: 2mm Shrank type: type 1 – diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65008 	CD2-10	Surgical drill D2.0mm – L10mm	Mass: 1.85gr Height: 35.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 64.5mm Length of working part: 10mm Diameter: 2mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65009 	CD2-11.5	Surgical Drill D2.0mm – L11.5mm	Mass: 1.87gr Height: 37mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 69mm Length of working part: 11.5mm Diameter: 2mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65010 	CD2-13	Surgical Drill D2.0mm – L13mm	Mass: 1.87gr Height: 38.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 73.5mm Length of working part: 13mm Diameter: 2mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65011	CD2-16	Surgical Drill D2.0mm – L16mm	Mass: 1.93gr Height: 41.5mm Hardness: 48-54 HRC Roughness: 32 μinch	Stainless steel 410S	Surgical drill with extended body which precisely fit the

			General length:82.5mm Length of working part:16mm Diameter: 2mm Shrank type: type 1- diameter 2.35 mm		guide sleeve
65012 	CD2.8-8	Surgical Drill D2.8mm – L8mm	Mass: 1.95gr Height: 33.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length:58.5mm Length of working part:8mm Diameter: 2.8mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65013 	CD2.8-10	Surgical Drill D2.8mm – L10mm	Mass: 2.01gr Height: 35.5mm Hardness: 48-54 HRC Roughness: 32 μinch Hardness: 48-54 HRC Roughness: 32 μinch General length:64.5mm Length of working part:10mm Diameter: 2.8mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65014 	CD2.8-11.5	Surgical Drill D2.8mm – L11,5mm	Mass: 2.04gr Height: 37mm Hardness: 48-54 HRC Roughness: 32 μinch Hardness: 48-54 HRC Roughness: 32 μinch General length:69mm Length of working part:11.5mm Diameter: 2.8mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65015 	CD2.8-13	Surgical Drill D2.8mm – L13mm	Mass: 2.09gr Height: 38.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length:73.5mm Length of working part:13mm Diameter: 2.8mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65016 	CD2.8-16	Surgical Drill D2.8mm – L16mm	Mass: 2.16gr Height: 41.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length:82.5mm Length of working part:16mm Diameter: 2.8mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65017	CD3.2-8	Surgical Drill	Mass: 2.01gr	Stainless steel	Surgical drill

		D3.2mm – L8mm	Height: 33.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 58.5mm Length of working part: 8mm Diameter: 3.2mm Shrank type: type 1- diameter 2.35 mm	410S	with extended body which precisely fit the guide sleeve
65018 	CD3.2-10	Surgical Drill D3.2mm – L10mm	Mass: 2.09gr Height: 35.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 64.5mm Length of working part: 10mm Diameter: 3.2mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65019 	CD3.2-11.5	Surgical Drill D3.2mm – L11,5mm	Mass: 2.14gr Height: 37mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 69.5mm Length of working part: 11.5mm Diameter: 3.2mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65020 	CD3.2-13	Surgical Drill D3.2mm – L13mm	Mass: 2.19gr Height: 38.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 73.5mm Length of working part: 13mm Diameter: 3.2mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65021 	CD3.2-16	Surgical Drill D3.2mm – L16mm	Mass: 2.3gr Height: 41.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 82.5mm Length of working part: 16mm Diameter: 3.2mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65022 	CD3.6-8	Surgical Drill D3.6mm – L8mm	Mass: 2.12gr Height: 33.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 58.5mm Length of working part: 8mm Diameter: 3.65mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65023	CD3.6-10	Surgical Drill D3.6mm – L10mm	Mass: 2.21gr Height: 35.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 64.5mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve


				Length of working part: 10mm Diameter: 3.65mm Shrank type: type 1- diameter 2.35 mm		
65024 	CD3.6-11.5	Surgical Drill D3.6mm – L11.5mm		Mass: 2.27gr Height: 37mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 69mm Length of working part: 11.5mm Diameter: 3.65mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65025 	CD3.6-13	Surgical Drill D3.6mm – L13mm		Mass: 2.33gr Height: 38.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 73.5mm Length of working part: 13mm Diameter: 3.65mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65026 	CD3.6-16	Surgical Drill D3.6mm – L16mm		Mass: 2.46gr Height: 41.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 82.5mm Length of working part: 16mm Diameter: 3.65mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65027 	CD4.1-8	Surgical Drill D4.1mm – L8mm		Mass: 2.18gr Height: 33.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 58.5mm Length of working part: 8mm Diameter: 4.1mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65028 	CD4.1-10	Surgical Drill D4.1mm – L10mm		Mass: 2.28gr Height: 35.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 64.5mm Length of working part: 10mm Diameter: 4.1mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65029 	CD4.1-11.5	Surgical Drill D4.1mm – L11.5mm		Mass: 2.37gr Height: 37mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 69mm Length of working part: 11.5mm	Нержавеющая сталь 410S	Хирургическое сверло с удлиненной рабочей частью, которая точно





			Diameter: 4.1mm Shrank type: type 1- diameter 2.35 mm		соответствует шаблону направляющей штулки
65030 	CD4.1-13	Surgical Drill D4.1mm – L13mm	Mass: 2.42gr Height: 38.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 82.5mm Length of working part: 16mm Diameter: 4.1mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65031 	CD4.1-16	Surgical Drill D4.1mm – L16mm	Mass: 2.56gr Height: 41.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 82.5mm Length of working part: 16mm Diameter: 4.1mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65032 	CD4.5-8	Surgical Drill D4.5mm – L8mm	Mass: 3.17gr Height: 33.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 58.5mm Length of working part: 8mm Diameter: 4.5mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65033 	CD4.5-10	Surgical Drill D4.5mm – L10mm	Mass: 3.32gr Height: 35.5mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 64.5mm Length of working part: 10mm Diameter: 4.5mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65034 	CD4.5-11.5	Surgical Drill D4.5mm – L11.5mm	Mass: 3.42gr Height: 37mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 69mm Length of working part: 11.5mm Diameter: 4.5mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65035	CD4.5-13	Surgical Drill D4.5mm – L13mm	Mass: 3.50gr Height: 38.5mm Hardness: 48-54 HRC	Stainless steel 410S	Surgical drill with extended body which




			Roughness: 32 μ inch General length: 73.5mm Length of working part: 13mm Diameter: 4.5mm Shank type: type 1- diameter 2.35 mm		precisely fit the guide sleeve
65036 	CD4.5-16	Surgical Drill D4.5mm – L16mm	Mass: 3.70gr Height: 41.5mm Hardness: 48-54 HRC Roughness: 32 μ inch General length: 82.5mm Length of working part: 16mm Diameter: 4.5mm Shank type: type 1- diameter 2.35 mm	Stainless steel 410S	Surgical drill with extended body which precisely fit the guide sleeve
65070 	CD2.4-8	Surgical Drill D2.4mm – L8mm	Mass: 1.88gr Height: 33.5mm Hardness: 53-57 HRC Roughness: 32 μ inch General length: 58.5 Length of working part: 8 Diameter: 2.4mm Shank type: type 1- diameter 2.35 mm	Stainless steel 440C	Surgical drill with extended body which precisely fit the guide sleeve
65071 	CD2.4-10	Surgical Drill D2.4mm – L10mm	Mass: 1.92gr Height: 35.5mm Hardness: 53-57 HRC Roughness: 32 μ inch General length: 64.5mm Length of working part: 10mm Diameter: 2.4mm Shank type: type 1- diameter 2.35 mm	Stainless steel 440C	Surgical drill with extended body which precisely fit the guide sleeve
65072 	CD2.4-11.5	Surgical Drill D2.4mm – L11.5mm	Mass: 1.96gr Height: 37mm Hardness: 53-57 HRC Roughness: 32 μ inch General length: 69mm Length of working part: 11.5mm Diameter: 2.4mm Shank type: type 1- diameter 2.35 mm	Stainless steel 440C	Surgical drill with extended body which precisely fit the guide sleeve
65073 	CD2.4-13	Surgical Drill D2.4mm – L13mm	Mass: 1.98gr Height: 38.5mm Hardness: 53-57 HRC Roughness: 32 μ inch General length: 73.5mm Length of working part: 13mm Diameter: 2.4mm Shank type: type 1- diameter 2.35 mm	Stainless steel 440C	Surgical drill with extended body which precisely fit the guide sleeve
65074	CD2.4-16	Surgical Drill D2.4mm – L16mm	Mass: 2.03gr Height: 41.5mm Hardness: 53-57 HRC Roughness: 32 μ inch General length: 82.5mm Length of working part: 16mm	Stainless steel 440C	Surgical drill with extended body which precisely fit the guide sleeve

			Diameter: 2.4mm Shrank type: type 1- diameter 2.35 mm		
2. Tissue punch, variations:					
65003 	TPS	Tissue Punch Small	Mass: 1.2gr Height: 30.3mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 30.3mm Length of working part: 30.3mm Diameter: 4.10mm	Stainless steel 410S	Punch drill intended for removing a small portion of the gingiva, exposing the bone before drilling (small diameter punched hole)
65004 	TPL	Tissue Punch Large	Mass: 1.4gr Height: 30.3mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 30.3mm Length of working part: 30.3mm Diameter: 5.45mm	Stainless steel 410S	Punch drill intended for removing a small portion of the gingiva, exposing the bone before drilling (small diameter punched hole)
3. Crestal drill, variations:					
65005 	CDS	Crestal Drill Small	Mass: 1.7gr Height: 29mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 29mm Length of working part: 3.5mm Diameter: 4.10mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	The primary drill intended for drilling only through the crestal bone (3 mm drilling depth) - small
65006 	CDL	Crestal Drill Large	Mass: 2.33gr Height: 29mm Hardness: 48-54 HRC Roughness: 32 μinch General length: 29mm Length of working part: 3.5mm Diameter: 5.45mm Shrank type: type 1- diameter 2.35 mm	Stainless steel 410S	The primary drill intended for drilling only through the crestal bone (3 mm drilling depth) - small
4. Pin, variations:					
65047	LP	Lateral Pin	Mass: 0.67gr Height: 28.5mm Hardness: 32-36 HRC Roughness: 32 μinch General length: 28.5mm Length of working part: 24.5mm Diameter: 4.9mm	Ti-6Al-4V ELI	Metal pin intended for laterally fixing the guide to the jaw bone



					
65048 	CPS	Crestal Pin Small	Mass: 1.91gr Height: 21.8mm Hardness: 48-54 HRC Roughness: 32 µinch General length: 21.8mm Length of working part: 5.5mm Diameter: 4.1mm	Stainless steel 410S	Metal pin intended for vertically fixing the guide to the jaw bone (small body diameter)
65049 	CPL	Crestal Pin Large	Mass: 2.70gr Height: 21.8mm Hardness: 48-54 HRC Roughness: 32 µinch General length: 21.8mm Length of working part: 5.5mm Diameter: 5.45mm	Stainless steel 410S	Metal pin intended for vertically fixing the guide to the jaw bone (small body diameter)
5. Milling drill					
65050 	MCD1.5	Milling drill, D1.5mm	Mass: 0.66gr Height: 35mm Hardness: 48-54 HRC Roughness: 32 µinch General length: 63.5mm Length of working part: 8mm Diameter: 1.5mm Shank type: type 1- diameter 2.35 mm	Stainless steel 410S	1.5 mm diameter drill preparing a 1.5 mm hole to fit the following Lateral pin insertion (65047)
6. Implant mount, variations:					
65037 	IMS	Implant Mount Small	Mass: 0.67gr Height: 15.5mm Hardness: 32-36 HRC Roughness: 32 µinch General length: 15.5mm Length of working part: 15.5mm Diameter: 4.8mm	Ti-6AL-4V ELI	Small IH mounter connected to the implant by 65039 screw for implant insertion through the guide
65038 	IML	Implant Mount Large	Mass: 1.03 gr Height: 15.5mm Hardness: 32-36 HRC Roughness: 32 µinch General length: 15.5mm Length of working part: 15.5 mm Diameter: 4.8mm	Ti-6AL-4V ELI	Large IH mounter connected to the implant by 65039 screw for implant insertion through the guide
65055	IMC	Implant Mount CHC	Mass: 0.71 gr Height: 15.3mm Hardness: 32-36 HRC	Ti-6AL-4V ELI	CHC mounter connected to the implant by

			Roughness: 32 µinch General length: 15.3mm Length of working part: 15.3mm Diameter: 4.8mm		65056 screw for implant insertion through the guide
65064 	IMSS	Implant Mount CS Small	Mass: 0.64 gr Height: 15.85mm Hardness: 32-36 HRC Roughness: 32 µinch General length: 15.85mm Length of working part: 15.85mm Diameter: 4.8mm	Ti-6AL-4V ELI	Small CS mounter connected to the implant by the 65039 screw for implant insertion through the guide
65065 	IMSL	Implant Mount CS large	Mass: 1.03gr Height: 15.85mm Hardness: 32-36 HRC Roughness: 32 µinch General length: 15.85mm Length of working part: 15.85mm Diameter: 4.8mm	Ti-6AL-4V ELI	Large CS mounter connected to the implant by the 65039 screw for implant insertion through the guide
7. Driver, variations:					
65061 	IMCD	Implant Mount CHC Driver	Mass: 2.35gr Height: 20.17mm Diameter: 6.2mm Hardness: 53-57 HRC Roughness: 32 µinch General length: 20.17mm Length of working part: 20.17mm Hexagons width: 2.25mm Sharp edges R: 1.25mm	Stainless steel 440C	Grip driver for CHC implants
65062 	IMSD	Implant Mount Int. Hex Driver Small	Mass: 2.39gr Height: 19.9mm Diameter: 6.2mm Hardness: 53-57 HRC Roughness: 32 µinch General length: 19.9mm Length of working part: 19.9mm Hexagons width: 2.74mm Sharp edges R: 1.43mm	Stainless steel 440C	Grip driver for IH implants (Small diameter body)
65063 	IMLD	Implant Mount Int. Hex Driver Large	Mass: 3.06gr Height: 19.9mm Diameter: 6.2mm Hardness: 53-57 HRC Roughness: 32 µinch General length: 19.9mm Length of working part: 19.9mm Hexagons width: 2.74mm Sharp edges R: 1.43mm	Stainless steel 440C	Grip driver for IH implants (Large diameter body)
4057 	HTD1.5	Short Hex Driver 1.5 mm	Mass: 1.95 g General length: 23mm Length of working part: 14.5mm Hexagon width: 1.5 mm Diameter: 7.2mm Sharp edges R: 0.03x45° Hardness: 50-57 HRC Roughness: 32 µinch (0.8 microns)	Hardening stainless steel SANDVIK 1RK91	Fits hexagonal 6.35 mm or square 4 mm ratchet or surgical screwdriver
4058	HTD1.5S	Pro Short Hex Driver 1.5 mm	Mass: 1.72 g General length: 15.5mm	Hardening stainless steel	Fits hexagonal 6.35 mm or

			Length of working part: 7mm Hexagon width: 1.5mm Diameter: 7.2mm Sharp edges R: 0.03x45° Hardness: 50-57 HRC Roughness: 32 µinch (0.8 microns)	SANDVIK 1RK91	square 4 mm ratchet or surgical screwdriver
4059 	HHS1.5	Pro Hand Hex Screw Driver 1.5 mm	Mass: 2.86 g General length: 18.1mm Length of working part: 7mm Hexagon width: 1.5 mm Diameter: 7.85mm Sharp edges R: 0.03x45° Hardness: 50-57 HRC Roughness: 32 µinch (0.8 microns)	Hardening stainless steel SANDVIK 1RK91 (body) + Ti-6Al-4V (handle)	For manual use
4060 	HHL1.5	Pro Hand Hex Long Screw Driver 1.5 mm	Mass: 4.27 g General length: 28.1mm Length of working part: 13 mm Hexagon width: 1.5 mm Diameter: 7.85mm Sharp edges R: 0.03x45° Hardness: 50-57 HRC Roughness: 32 µinch (0.8 microns)	Hardening stainless steel SANDVIK 1RK91 (body) + Ti-6Al-4V (handle)	For manual use
4168 	HT1.5	Motor Mount Hex Driver 1.5 mm	Mass: 0.60gr General length: 21 mm Length of working part: 7mm Hexagons width: 1.5 mm Diameter: 2.35 0/-0.02 mm Sharp edges R: 0.03x45° Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	To be used with a contra angle motor
4140 long 	G- ITDL2.5	Internal Hex Insertion Drivers Long	Mass: 2.88gr General length: 27mm Length of working part: 18.5mm Hexagons width: 2.43 mm Diameter: 7.2mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	Compatible with hexagonal 6.35 mm or square 4 mm ratchet or surgical screwdriver
4141 stand ard 	G- ITDM2.5	Internal Hex Insertion Drivers Standard	Mass: 2.23gr General length: 18mm Length of working part: 9.5mm Hexagons width: 2.43 mm Diameter: 7.2mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	Compatible with hexagonal 6.35 mm or square 4 mm ratchet or surgical screwdriver
4142 short 	G- ITDS2.5	Internal Hex Insertion Drivers Short	Mass: 1.74gr General length: 14.5mm Length of working part: 6.0mm Hexagons width: 2.43mm Diameter: 7.2mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	Compatible with hexagonal 6.35 mm or square 4 mm ratchet or surgical screwdriver
4143, long	GITL2.5/1 .25	Internal Hex Contra-Angle Driver Long motor mount 2.5/1.25 mm	Mass: 0.83gr General length: 25.0mm Length of working part: 25.0mm Hexagons width: 2.43 mm Diameter: 2.74mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	Implant driver for IH platform which is used to insert the implant into the osteotomy and compatible with

					motor mounts.
4144 	GITM2.5/ 1.25	Internal Hex Contra-Angle Driver Standard	Mass: 0.62gr General length: 20.0 mm Length of working part: 20.0mm Hexagons width: 2.43mm Diameter: 2.74mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	Implant driver for IH platform which is used to insert the implant into the osteotomy and compatible with motor mounts.
4145 short 	GITS2.5/ 1.25	Internal Hex Contra-Angle Driver Short motor mount 2.5/1.25 mm	Mass: 0.59gr General length: 18.0mm Length of working part: 18.0mm Hexagons width: 2.43 mm Diameter: 2.74mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	Implant driver for IH platform which is used to insert the implant into the osteotomy and compatible with motor mounts.
4146 	MITD2.5- IH	Manual Driver	Mass: 4.07gr General length: 25.5mm Length of working part: 11.5mm Hexagons width: 2.43mm Diameter: 7.854mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	For manual use
4147 	MITD2.1- CHC	Manual Driver	Mass: 4.07gr General length: 25.5mm Length of working part: 11.5mm Hexagons width: 2.09 mm Diameter: 7.854mm Hardness: 53-57 HRC Roughness: 32 µinch	Stainless steel 440C	For manual use
8. Implant mount screw					
65039 	IMHS	Implant Mount Screw	Mass: 0.39 gr Height: 22.1mm Hardness: 32-36 HRC Roughness: 32 µinch General length: 22.1mm Length of working part: 22.1mm Diameter: 2.75mm	Ti-6AL-4V ELI	Screw for mounter to implant connection
9. Implant mount extention					
65042 	IMX	Implant mount extention	Mass: 1.40 gr Height: 6mm General length: 12mm Length of working part: 4mm Diameter: 5mm	hardening stainless steel SANDVIK 1RK91	Extension part intended to prolong the mounter connection in cases of tight area which could not be reached by our standard tooling

10. Handpiece insertion adapter					
65044	HIA	Handpiece Insertion Adaptor	Mass: 1.42gr Height: 6mm General length: 23.1mm Length of working part: 14.1mm Diameter: 6mm	hardening stainless steel SANDVIK 1RK91	Mounter insertion tool which fit motor-mount handpiece connection
					
11. Implant mount extractor					
65045	IME	Implant Mount Extractor	Mass: 1.00 gr Height: 32mm Hardness: 32-36 HRC Roughness: 32 µinch General length: 32mm Length of working part: 32mm Diameter: 6mm	Ti-6Al-4V ELI	Metal extractor intended to aid in removing the mounters
					
12. Implant mount screw - CHC					
65056	IMCS	Implant Mount Screw CHC	Mass: 0.29gr Height: 21.5mm Hardness: 32-36 HRC Roughness: 32 µinch General length: 21.5mm Length of working part: 21.5mm Diameter: 2.75mm	Ti-6Al-4V ELI	Screw for mounter to implant connection
					
13. Sleeve L/S adaptor driver					
65057	SAD	Sleeve L/S Adaptor Driver	Mass: 2.37gr Height: 20.5mm Hardness: 53-57 HRC Roughness: 32 µinch General length: 20.5mm Length of working part: 20.5mm Diameter: 6mm	Stainless steel 440C	Insertion tool intended for insertion of the small sleeve adaptor (65058) to the large sleeve
					
14. Sleeve for adaptor					
65058	SLSA	Small Sleeve for Adaptor	Mass: 0.31gr Height: 4mm Hardness: 53-57 HRC Roughness: 32 µinch General length: 4mm Length of working part: 4mm Diameter: 5.5mm	Stainless steel 440C	Small metal sleeve submerge and fix in the guide
					
15. Sleeve, variations:					
66012	SLS	Sleeve small	Mass: 0.03gr General length: 4mm Length of working part: 4mm Diameter: 5mm	Stainless steel 630 (17-4PH)	Used to guide and position the drills and implants in the planned position and depth.
					
66013	SLL	Sleeve large	Mass: 0.36gr Height: 4mm	Stainless steel 630 (17-4PH)	Used to guide and position the

			General length: 4mm Length of working part: 4mm Diameter: 7.3mm		drills and implants in the planned position and depth.
66014 	SLSE	Securing sleeve	Mass: 0.03gr Height: 3.8±0.1mm General length: 7.6±0.05mm Length of working part: 3.5mm Diameter: 3.2mm	Stainless steel 630 (17-4PH)	Used to secure the surgical guide in place with a lateral pin
Accessories:					
1. Box					
65002 	DWT	DW GSTK Box	Mass: 670gr Dimensions: Length: 269mm Height: 148.3mm Height of cover: 56.2mm	- Radel (R-5000): Cover; - silicon (SID 501): Strips; - Radel polyphenylsulfone (PPSU) (also for 4613): Base.	GSTK Box
4699 	OBSA-CS	CS Advanced Surgical Organizer Box	Mass: 480gr Dimensions: Length: 185.06mm Width: 133.49mm Height of cover: 26.5mm	- Radel (R-5000): Cover; - silicon (SID 501): Strips; - Radel polyphenylsulfone (PPSU) - base	Organizer Box

Delivery set of medical devices

Dental instruments, packed in an individual (primary) package, are delivered to the consumer by the piece.

Instructions for the use of products are provided to the consumer upon delivery by an authorized representative upon request.

Information about cleaning, disinfection, sterilization

Cleaning, disinfection and sterilization before the first clinical application

1. Immerse components in an approved disinfectant. Refer to the manufacturer instructions for use, do not use disinfectants containing chlorine, ammonia and aldehydes.

2. Wash the used instruments at least three times in distilled water and dry them with a clean lint-free towel. Alternately, dry the instruments by compressed air.

3. Place the clean instruments in the Surgical Organizer Box.

4. Sterilize the instruments using a steam autoclave according to autoclave manufacturer's recommendations (at a temperature of 132°C / 207°F for 7 minutes followed by a 30-minute dry cycle vapor pressure of up to 2 Bar).

Distilled water should be used to avoid surface stains. Make sure before use that the elements, inside the autoclave, are not rusted. Chemclave is NOT recommended.

5. Do not exceed the autoclave manufacturer's stated maximum load.

Cleaning, disinfection and sterilization after working

1. Place instruments in the disinfectant solution immediately after surgery to prevent encrustation of blood, secretions, etc.

2. Meticulously eliminate all the post-operative residual tissue or bone from used instruments including drills, by immersing and rinsing in cold water (<40°C).

3. Clean the used instruments thoroughly first by immersing them in mild, pH-neutral enzymatic solution for up to 20 minutes or until rinsing can be performed.

Avoid contact with phenol alcohol, chlorine, acid or quaternary ammonia.

Rinse the instruments under a hard stream of cold water (<40°C). Avoid water with high concentration of chlorine.

4. Flush the drill lumen with hypodermic needle.

5. Scrub the used instrument with a soft nylon brush. Make sure to meticulously eliminate all the post-operative residues (blood, bone etc.) from the irrigation channel.

6. Immerse components in an approved disinfectant. Refer to the manufacturer instructions for use, do not use disinfectants containing chlorine, ammonia and aldehydes.

7. Wash the used instruments at least three times in distilled water and dry them with a clean lint-free towel. Alternately, dry the instruments by compressed air.

8. Place the clean instruments in the Surgical Organizer Box.

9. Sterilize the instruments using a steam autoclave according to autoclave manufacturer's recommendations (at a temperature of 132°C / 207°F for 7 minutes followed by a 30 minute dry cycle vapor pressure of up to 2 Bar).

Distilled water should be used to avoid surface stains. Make sure before use that the elements, inside the autoclave, are not rusted. Chemclave is NOT recommended.

10. Do not exceed the autoclave manufacturer's stated maximum load.

Visual inspection:

After cleaning, disinfection, sterilization and drying, inspect the product for unacceptable damage such as corrosion, discoloration, pitting, cracked seals.

Reprocessing limitations

Repeated reprocessing in line with these instructions has only minimal effect on Alpha-Bio Tec GSTK drills. Drill end of life is normally determined by wear and attrition during use rather than by the effects of reprocessing. Drills should be visually inspected for signs of wear prior to reuse.

Maintenance instructions

- A visual inspection is carried out to exclude contamination.
- A visual inspection is carried out for signs of damage and wear. The cutting edges must be free of defects, and the product itself must have the correct shape.
- Drills Tools showing signs of damage and corrosion must be disposed of.
- Discoloration or blunting of the incisal edge may occur after steam sterilization.
- Resterilization reduces cutting efficiency over time.
- Drills are recommended to be replaced after 6-10 uses.

Labelling

Labelling of medical device, which is supplied in packaging, contains the following information:

- Product name;
- LOT;

- Catalogue number;
- caution;
- instruction for use;
- trademark and manufacturer's address;
- «non sterile»;
- Manufacturing date;
- «for practitioners only»;
- quantity.

Table of applicable symbols

	Caution
	Catalogue number
LOT	Series code
	CE marking
	«for practitioners only»
	Non sterile
	Manufacturer
	Manufacturing date
	Consult instruction for use
	Instructions for use are available at the website: www.e-labeling.eu/ABT
	Instructions for use are available by phone request

Transport, operation and storage conditions

During operation instruments are resistant to temperatures from +32 to + 42C and relative humidity of 80%. Products are resistant to the effects of biological fluids and secretions of body tissues with which they come into contract during operation.

The products may be used in public or private dental clinics and only by surgeons who have been trained and have sufficient experience in this field.

Transportation is carried out by transport in accordance with local legislation.

Transport conditions: temperature from -50 °C to + 50 °C and relative humidity 75% at 15 °C.

During storage, the drill must be kept in a dry, clean environment at a temperature of 15÷30°C.

Exposure to conditions that affect the marking of the drill or its expiration date, change its surface or geometric properties (for example, excessive shaking, pressure, heat and UV radiation) should be avoided.

All necessary measures should be taken to prevent re-contamination before subsequent use.

List of applied standards

Number	Name
MEDDEV 2.7.1 Rev. 4	Clinical evaluation: Guide for manufacturers and notified bodies
EN ISO 13485	Medical devices — Quality management systems — Requirements for regulatory purposes
EN ISO 15223-1	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
EN ISO 14971	Medical devices - Application of risk management to medical devices
EN ISO 7405	Dentistry. Evaluation of biocompatibility of medical devices used in dentistry
EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 17664-1	Processing of health care products — Information to be provided by the medical device manufacturer for the processing of medical devices — Part 1: Critical and semi-critical medical devices
ISO 1797	Dentistry — Shanks for rotary and oscillating instruments
EN ISO 13504	Dentistry — General requirements for instruments and related accessories used in dental implant placement and treatment
ISO 17665-1	Sterilization of health care products — Moist heat — Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
EN 62366 - 1	Medical devices - Part 1: Application of usability engineering to medical devices
EN 1639	Dentistry - Medical devices for dentistry - Instruments
ASTM F136	Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401)
ASTM D4169	Standard Protocol for Performance Testing of Shipping Containers and systems

Warranty

Alpha-Bio Tec warrants that all products will be free of defects in materials and/or workmanship. This warranty applies to the original purchaser only. There are no warranties, express or implied, except this warranty, which is given in lieu of any other warranties, express or implied, including any implied warranty of fitness for a particular purpose.

The operating surgeon / practitioner in charge must carefully read the user instructions prior using the drills.

Alpha-Bio Tec is not liable for complications, other negative effects or damages that might occur for reasons such as incorrect indications or surgical technique, unsuitable choice of material or handling

thereof, unsuitable use or handling of the instruments, asepsis and so on. The operating surgeon is responsible for any such complications or other consequences.

Procedure and conditions for medical device disposable

Potentially contaminated or no longer usable medical devices must be safely disposed of as medical (clinical) waste in accordance with local public health guidelines, national laws and policies, and national laws and policies.

Segregation, recycling, or disposal of packaging material must be carried out in accordance with the local laws of the country and government legislation on packaging and packaging waste, as applicable.

Disposal should be carried out in compliance with the requirements for class B epidemiologically hazardous waste, in accordance with sanitary rules and SanPiN 2.1.3684-21.

Information about manufacturer

Manufacturer

Alpha-Bio Tec Ltd.,

7 Hatnufa St., Kiryat Arye, POB 3936, Petach Tikva, 4951025, Israel

Tel.: +972-3-9291055

Authorized representative in Russian Federation:

Nobel Biocare Russia LLC

Address: 109004 Moscow, Stanislavskogo str., 21, building 2.

Tel.: +7 495 974 77 55, Fax: +7 495 974 77 66

E-mail: cs-team.russia@nobelbiocare.com

In case of appearance of undesirable events, which have signs of an adverse events, it's necessary to send a message to authorized representative on the territory of the Russian Federation.

Date of issue or last revision of the instructions for use May 2023

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