**IMAXEON**

for. x 20

November 1, 2021

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JD

November 1, 2021

To be submitted to the Competent Health Authorities of the Russian Federation.

We, Imaxeon Pty Ltd., Australia (Unit 1, 38-46 South Street, Rydalmere, New South Wales 2116, Australia) hereby state that the attached Instruction for use for MD "CENT-PL 250 cm Patient line for MEDRAD Centargo injection system" (Линия 250 см для одного пациента CENT-PL для системы инъекционной MEDRAD Centargo) is true.

Best regards,

S. Jafari

Sajjad Jafari

Regulatory Affairs Manager  
Unit 2, 38-46 South Street  
Rydalmere NSW 2116 Australia  
Mobil: +61 439078849  
E-mail: [sajjad.jafari@bayer.com](mailto:sajjad.jafari@bayer.com)

I hereby certify this to be a true and correct copy  
of the original document.

John Howard Fisher - Notary Public,  
Parramatta NSW.

JD  
6/11/2021**IMAXEON**

Imaxeon Pty Ltd  
ABN 93 093 950 906

Rydalmere Metro Centre Unit 1, 38-46 South St Rydalmere NSW 2116 Australia  
PO Box 6033 Silverwater NSW 1811 Australia

t: 612 8845 4999  
f: 612 8845 4936

e: [info@imaxeon.com](mailto:info@imaxeon.com)  
w: [www.imaxeon.com](http://www.imaxeon.com)

RESTRICTED

# APOSTILLE

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- |                                            |                                   |
|--------------------------------------------|-----------------------------------|
| 1. Country                                 | Australia                         |
| 2. This public document has been signed by | John Howard Fisher                |
| 3. acting in the capacity of               | Notary Public                     |
| 4. bears the seal/stamp of                 | John Howard Fisher, Notary Public |
| 5. at Sydney Passport Office               | Certified                         |
| 7. by Luke O'Neill                         | 6. the 18th day of November, 2021 |
| 8. No. 1 PPT-6A-142718                     |                                   |
| 9. Seal/Stamp                              |                                   |
| 10. Signature                              |                                   |

*This Apostille only certifies the authenticity of the signature (where applicable) and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. This Apostille can be verified at <https://orao.dfat.gov.au/pages/verifyapostille.aspx>*

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## ENGLISH

**Important Safety Notice**

This device is intended to be used by individuals with adequate training and experience in diagnostic studies in x-ray applications.

PLEASE FOLLOW YOUR INSTITUTION'S POLICY FOR MEDIA ADMINISTRATION AND ASEPTIC TECHNIQUE.

**Indications for Use**

The MEDRAD® Centargo CT Injection System is indicated for the specific purpose of injecting intravenous contrast media and flushing solutions into humans for diagnostic studies in x-ray applications.

**Contraindications**

None known.

**NOTE:** Refer to the injector system operation manual for additional operating instructions and instructions for configuring the injector system.

**NOTE:** Refer to contrast media manufacturer's package insert for specific indications and use time.

**NOTE:** If cleaning of disposables is desired, use only 70% isopropyl alcohol.

Exposure to other cleaning agents or disinfectants may lead to component failure.

**NOTE:** For products labeled as sterile fluid path, sterile barrier system refers to dust caps.

Report any serious incident that has occurred in relation to this device to Bayer (radiology.bayer.com/contact) and to your local European competent authority (or, where applicable, to the appropriate regulatory authority of the country in which the incident has occurred).

**Operating temperature:** 16°C to 28°C

**Operating humidity:** 20% to 90%

**Warning**

Follow aseptic technique at all times. Improper handling can cause infection.

**For Patient Line and Replacement Spike, inspect contents and package prior to use.** Do not use if sterile package is opened or damaged. Patient or operator injury may result.

**For Day Set, inspect components prior to use. Do not use if any dust caps are missing or displaced.** Do not remove dust caps until ready to make connections. Early removal of dust caps can result in contamination.


**Do not reinstall, re-sterilize, or reprocess Day Set.** Potential device failure includes significant component deterioration and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be reinstalled, re-sterilized, or reprocessed.

**The Patient Line is intended for single use only. Do not re-sterilize, reprocess, or reuse.** Potential device failure includes significant component deterioration and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be re-sterilized, reprocessed, or reused.

**Do not use the Day Set for longer than 24 hours.** Over-use poses risks to the patient, including contamination and injury due to device failure.

**Do not connect the Patient Line to the patient until all trapped air has been cleared.** Do not modify or attempt to circumvent the operation of the air sensors. Air embolization can cause death or serious injury.


**Check the labeling of any disposables for their maximum pressures. If none are provided, do not use. Ensure that the programmed pressure limit does not exceed the maximum labeled pressures, but is not too low so as to compromise the quality of the study.** Use of greater pressures can result in fluid leaks or tubing ruptures and patient or operator injury.

**Install New Day Set:**

1. Press button to unlock injector door.
2. Insert new Day Set straight into injector.


**NOTE:** Maximum use time of Day Set is 24 hours.

3. Push down on lock levers. Lock levers click when Day Set is locked into place.
4. Push small-diameter tube into outlet air sensor and close outlet air sensor door.


**Warning**

Use care in handling and inserting spikes into saline and contrast. The spikes are sharp and may cause injury.


Do not remove dust caps until ready to make connections. Early removal of dust caps can result in contamination.

**Install Spike Adapters onto Injector:**

1. Snap on spike adapters from Day Set onto assigned inlet air detectors.
- NOTE:** Ensure spike adapters are aligned with corresponding fluid reservoirs.
- NOTE:** Adapters click when locked into place.
2. Ensure outlet air sensor door is closed.
3. Close injector door.

**Load Contrast:**


1. Remove spike dust cap (if present).
2. Spike contrast onto middle or far right spike.
- NOTE:** Spike adapter may be removed from injector to ease spiking.
3. Secure contrast.
- NOTE:** For bottles, use clamp to secure. For bags, use clamp hook to hang.
4. Select **Injector** on display, and select one contrast fluid icon.
5. Select contrast type on display, and press **FILL**.
6. For second contrast, repeat steps 1-5.

**Load Saline:**


1. Remove spike dust cap (if present).
2. Spike saline on far left spike.
- NOTE:** Spike adapter may be removed from injector to ease spiking.
3. Secure saline.
- NOTE:** For bottles, use clamp to secure. For bags, use clamp hook to hang.
4. Select **Injector** on display, and select saline fluid icon.
5. Press **FILL** on display.

**Warning**

Do not reuse or reconnect the Patient Line. Cross-contamination can cause infection.

**Install Patient Line:**


1. Remove orange Day Set cap (if present).
2. Insert Patient Line until it clicks.
- NOTE:** System primes Patient Line automatically. If lights are solid blue, Patient Line is primed and ready. If lights are red, Patient Line is not primed. Refer to the operation manual for troubleshooting.

**Prepare and Inject Patient:**

1. Check Patient Line for air.
- NOTE:** If additional prime fluid is needed, press and hold **Advance Saline** button on injector. Saline will be pushed through Patient Line.
2. Disconnect patient end of Patient Line from injector.
3. Connect Patient Line to patient.
4. Perform injection.

**Complete Patient Exam:**

1. After injection is complete, press **End Exam** on display.
2. Disconnect used Patient Line from injector and patient.
3. Discard used Patient Line, per facility guidelines.
4. Insert new Patient Line until it clicks.

**Eject Day Set:**


1. (Optional) Empty Day Set. Refer to the operation manual for instruction.
2. Press **Eject Day Set** on display.
3. Press button to unlock injector door.
4. Lift lock levers.
5. Remove small-diameter tube from outlet air sensor, and close outlet air sensor door.
6. Discard contrast and saline containers.
7. Remove and discard used Day Set, per facility guidelines.
8. Ensure outlet air sensor door is closed, then close injector door.

**Within 24 Hours:****Replace Fluids As Needed:**

- ◆ For additional contrast, repeat steps 1-5 in "Load Contrast."
- ◆ For additional saline, repeat steps 1-5 in "Load Saline."

**Replace Spike If Required:**

1. Press both levers on left and right side of spike adapter.
2. While pressing, remove used spike.
3. Insert new spike into spike adapter.
- NOTE:** When installing a new spike, the circular orange vent must align with indent on back of adapter.

**Empty Prime Container As Needed:**

1. Press button to unlock injector door.
2. Remove prime container and dispose of fluids per facility guidelines.
3. Replace prime container, and close injector door.

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© 2018-2021 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.
Bayer, the Bayer Cross, MEDRAD, Centargo, MEDRAD Centargo, and Imaxeon are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.