This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille does not certify the content of the document for which it is issued.

Apostille

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1. Country: United States of America
This public document
2. has been signed by **Barbara H. Meiklejohn**
3. acting in the capacity of **Clerk of the Circuit Court for Montgomery County**
4. bears the seal/stamp of the **Circuit Court for Montgomery County**

Certified

5. at Annapolis, Maryland
6. the **26th day of July, 2019**
7. by The Secretary of State of Maryland
8. No. 492028
9. Seal
- 
- The seal of the Secretary of State of Maryland is partially visible at the bottom of the page. It is a circular gold-colored seal with the words "OFFICE OF THE SECRETARY OF STATE" and "MARYLAND" inscribed around the perimeter. The seal is positioned to the right of the word "Seal" in the line number 9.

John C. Wickersham

Secretary of State

Affidavit

State of Maryland)
) SS:
County of Anne Arundel)

Before me, the undersigned notary public, this day, personally, appeared to me known or proven, who being duly sworn according to law, deposes the following:

I, Vera Chawla, hereby certify that the attached document is


- ☒ an original document
☐ a true copy of an original document

to the best of my knowledge and belief.

V. Chawla
Signed

Subscribed and sworn to before me this

24 day of July, 2019.

Notary Public 

To whom it may concern,

July 21, 2019

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
(412) 767-2400
www.radiology.bayer.com

To be submitted to the Competent Health Authorities in Federal Republic of Russia

We, Bayer Medical Care Inc., USA (1 Bayer Drive Indianola Pennsylvania 15051-0780 USA) hereby state that attached document (Addendum to Instruction for use) (Multi-Patient Disposable Set MEDRAD Avanta consisting of) is true.

Best regards,

William Hohn

Regulatory Affairs International Strategist – Europe

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Apostille

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1. Country: United States of America
This public document
2. has been signed by **Robert P. Duckworth**
3. acting in the capacity of **Clerk of the Circuit Court for Anne Arundel County**
4. bears the seal/stamp of the **Circuit Court for Anne Arundel County**

Certified

5. at Annapolis, Maryland
6. the **1st day of October, 2018**
7. by The Secretary of State of Maryland
8. No. 465286
9. Seal

10. Signature

Secretary of State

Affidavit

State of Maryland)
)SS:
County of Anne Arundel)

Before me, the undersigned notary public, this day, personally,
appeared to me known or proven, who being duly sworn according to law,
deposes the following:

Yoshiko Yano

I, _____, hereby certify that the attached
document is

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
27 Sept, 2018

To whom it may concern:

To be submitted to the Competent Health Authorities in Federal Republic of Russia
We, Bayer Medical Care Inc., USA (1 Bayer Drive Indianola Pennsylvania 15051
USA) hereby state that attached document (Instruction for use) (Multi-Patient
Disposable Set MEDRAD Avanta consisting of) is true.

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UNITED STATES

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

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© 2010-2014, 2017 Bayer. All Rights Reserved.
Reproduction of this manual is strictly prohibited without
express written consent of Bayer Medical Care Inc.

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Bayer, the Bayer Cross, MEDRAD Avanta, MEDRAD, Avanta
and FluiDots may be registered trademarks of Bayer in the
US and other countries.

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Introduction: Read the information contained in this brochure. Understanding the information will assist you in using the sterile Disposable Sets properly.

Important Safety Notice: The MEDRAD Avanta Sterile Disposable Sets are intended to be used by medical professionals with adequate training and experience in the operation of the MEDRAD Avanta Fluid Management Injection System and angiographic procedures and techniques.

Intended Use: The MEDRAD Avanta Sterile Disposable Sets are specifically intended for use in the x-ray angiography environment. They are designed to administer intravascular radio-opaque contrast compounds, and common flushing agents at various volumes and flow rates into humans for use in diagnostic and interventional angiographic procedures performed in cardiology, radiology, and vascular surgery.

Contraindications: For use only with the MEDRAD Avanta Fluid Management Injection System. The MEDRAD Avanta Sterile Disposable Sets are not intended to be used for chemotherapy and are not intended to administer fluids other than intravascular contrast agents and common flushing solutions.

Restricted Sale: Rx Only.

Warnings

- ⚠ For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.
- ⚠ Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector, or attempt an injection until all air has been removed from the syringe and fluid path. Carefully read the instructions for loading and the use of FluidDots indicators (where applicable) to verify that fluid is present and to reduce the chance of air embolization during an injection.
- ⚠ Do not attempt to aspirate fluid from the waste port stopcock when open to the PIV. Doing so may increase the risk of inducing air into the fluid path.
- ⚠ Do not attempt to aspirate fluid from the hemodynamic port. Doing so may increase the risk of inducing air into the fluid path.
- ⚠ Before connecting the Single-Patient Sterile Disposable Set (SPAT) to the patient, ensure all stopcocks and open ports are closed to air and all air has been removed from the fluid path before injection. Improper manipulation of the waste port stopcock as well as any fluid port left open to air may increase the risk of introducing air into the fluid path.
- ⚠ Biological contamination can result from reusing disposable items or failure to follow sterile technique during setup or use. Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.
- ⚠ Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.
- ⚠ Patient or operator injury may result from mishandling of the sharp spike. Use care in handling and inserting the spike into the contrast bottle and the saline bag/bottle.

- ⚠ **Storage of filled syringes can promote bacterial growth.** Syringes from Bayer are designed to be filled just prior to the procedure. Discard unused filled syringes.
- ⚠ **Patient injury could result if syringe is not properly engaged.** Improper engagement may cause air embolization or under-volume delivery.
- ⚠ **Patient or operator injury may result from contrast media leaks or ruptures to the syringe or connector tubing.** Use of greater pressure or occlusions in the fluid path may result in leaks or ruptures to the syringe or connector tubing. Ensure that the fluid path is open; do not exceed appropriate pressure.
- ⚠ **Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe.** Do not remove the plunger to fill the syringe.
- ⚠ **Inadvertent aspiration can occur if the Multi-Patient disposable contrast tubing set is kinked or blocked during a syringe fill process.** Verify multi-patient contrast tubing is free of kinks before attempting to fill the syringe.
- ⚠ **Patient injury can result if the patient is connected during any contrast or saline purge.** Ensure the patient is disconnected before any purge.
- ⚠ **Blood can contaminate the multi-patient portion of the disposable set.** Ensure that blood does not appear past the pressure isolation valve in the single patient disposable set.
- ⚠ **Ensure vent cap on spike(s) are open if using bottles.** If closed, a vacuum may be created and air may be introduced.
- ⚠ **Do not exceed the five patients per day usage requirement for the Multi-patient Disposable Set (MPAT).** This is the usage approved by appropriate regulatory bodies.
- ⚠ **Use of non-Bayer supplied disposables, including administration sets and additions to administration sets, such as but not limited to bleed back control devices and pressure transducers may cause patient injury if not properly connected or flushed.** These devices must be compatible with your system. Refer to manufacturer's instructions for proper use of these devices.
- ⚠ **Use of extension tubing connected to the Single Patient Disposable Set may cause patient injury and will degrade system performance.** Do not connect extension tubing to the Single Patient Disposable Set.

Cautions


- ⚠ **Component damage may occur if disposables are not installed properly.** Assure all connections are secure. This will help minimize leaks, disconnections, air introduction, and component damage. Do not use tools to over tighten connections or to assist in removal of disposables.
- ⚠ **Do not use methyl alcohol on disposable sets.** Use of methyl alcohol will degrade the integrity of the disposable sets.
- ⚠ **All contrast and saline media should be used in accordance with the manufacturer's indication for use.**

NOTE: To perform the following procedures, the operator must select the Case Setup screens on the DCU. Refer to the injector operation manual for additional information.

Multi-Patient Sterile Disposable Set (MPAT) Installation (usable for up to five patients). The Multi-Patient Disposables Set consists of the following items. (Also see MPAT Setup illustration):

Vented Contrast Spike	2nd HPCT w/ male luer
Contrast Drip Chamber	Saline Spike
Low Pressure Contrast Tube (LPCT)	Saline Drip Chamber

High Pressure Contrast Tube (HPCT)	Low Pressure Saline Tube (LPST)
Contrast Flow Control Valve (CFCV)	Post Drip Chamber Saline Tube
150 ml Multi-Patient Syringe	

General Setup Instructions

Be sure that the patient is disconnected and the injector is on.

1. Install syringe. (See Install Multi-Patient Syringe below).
2. Install saline tubing. (See Install Saline Tubing below).
3. Install contrast tubing. (See Install Contrast Tubing below).
4. Rotate head to upright position, fill syringe and remove dust caps.
5. Purge multi-patient disposables. Catch the discharged fluid in a suitable container and tap disposable connection points and pressure jacket to remove air.
6. Rotate head downward.
7. Attach MEDRAD Avanta Sterile Hand Controller, and maintain sterility, if applicable.

NOTE: Refer to Operation Manual for MEDRAD Avanta Fluid Management Injection System for further information.

After step 7 in the General Setup Instructions is complete, go to SPAT Installation.

Install Multi-Patient Syringe

NOTE: Ensure patient is not connected.

Select Multi-Patient Syringe from Case Setup.

1. Retract piston.
2. Insert syringe into pressure jacket and raise drop front arms.
3. Engage the piston to the plunger.

Install Saline Tubing

NOTE: Ensure patient is not connected.

Select Multi-Patient Disposables from Case Setup, then select Install Saline Tubing.

1. Install drip chamber.
2. Route tubing through pump and close door.
3. Route tubing around peg.
4. Insert tubing into gross air detector and close door.
5. Spike saline bag.
6. Fill saline drip chamber completely and snap into position before use.

NOTE: To differentiate the saline and contrast lines of the MPAT, the contrast line will be identified as the line connected to the stopcock that is placed in the CFCV.

Install Contrast Tubing

Select Install Contrast Tubing from Case Setup.

1. Install drip chamber.

2. Connect tubing to syringe.
3. Install stopcock.
4. Insert tubing into gross air detector and close door.
5. Spike contrast bottle.
6. Fill contrast drip chamber completely and snap into position before use.

Go to step 4 in General Setup Instructions.

Single-Patient Sterile Disposable Set (SPAT) Installation (must be replaced for each patient use). The Single-Patient Disposable Set consists of the following items (Also see SPAT Setup illustration):

Rotating luer (A)	High Pressure Saline Tube w/ male luer connector (D)
Stopcock for waste port (B)	High Pressure Contrast Tube w/ female luer connector (E)
Pressure Isolation Valve (PIV) Assembly / Hemodynamic pressure monitoring port (C)	

Single-Patient Sterile Disposable Set (AVA 500 SPAT Angio) is indicated for use when a user attachable transducer is unnecessary and does not contain item (C) pictured below.

Select Single-Patient Disposables from Case Setup.

Be sure that the patient is disconnected.

1. Remove all dustcaps.
2. Connect the SPAT saline tube to MPAT saline tube.
3. Connect the SPAT contrast tube to MPAT contrast tube.

To purge air out of SPAT:

NOTE: Ensure patient is not connected.

When using AVA 500 SPAT L continue with step 4. When using AVA 500 SPAT ANGIO proceed to step 5.

4. Attach transducer to hemodynamic port, verify the transducer stopcock is open, and the waste port stopcock is off to waste port.
 - a. Attach 20cc saline syringe to the transducer port.
 - b. Hold SPAT parallel to the patient with the PIV and 20cc syringe connection facing upward at a 45° angle.
 - c. Forcefully flush through the transducer and PIV assembly until all air is expelled from both components. Flick or tap PIV to facilitate air removal (Repeat as necessary).
 - d. Close transducer port stopcock. Remove 20cc syringe.
5. Initiate contrast and saline purge, tap SPAT to MPAT saline connection point and wasteport stopcock to remove air while purging is in progress.
6. Repeat steps above as needed, until all air is purged.
7. Carefully inspect the entire SPAT and MPAT to verify all air is removed before connecting to patient.

NOTE: Disconnect the SPAT from the MPAT immediately following each procedure.

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©2017 Bayer. All Rights Reserved. Reproduction of this manual is strictly prohibited without express written consent of Bayer Medical Care Inc.
Bayer, the Bayer Cross, MEDRAD Avanta, MEDRAD, Avanta, and FluiDots may be registered trademarks of Bayer in the US and other countries.