To whom it may concern,

Date: 07.06.2019


To be submitted to the Competent Health Authorities in Federal Republic of Russia

We, Imaxeon Pty Ltd, Australia (Unit 1, 38-46 South Street, Rydalmere, New South Wales 2116 Australia) hereby state that attached document (Instructions for use) (Наборы одноразовые для системы инъекционной MEDRAD Salient) is true.

Note: Instructions for use (IFU) includes not only the sets but also other Imaxeon products that is "Extension tubing" are specified in the IFU. "Extension tubing" will be registered separately.


Best regards,

A handwritten signature in black ink, appearing to read 'Anhua Hu'.

Anhua Hu  
Regulatory Affairs Manager  
Imaxeon Pty Ltd

SIGNED IN MY PRESENCE

A handwritten signature in black ink, appearing to read 'John Howard Fisher'.

JOHN HOWARD FISHER  
NOTARY PUBLIC  
PARRAMATTA NSW  
AUSTRALIA

26 JUNE, 2019

## APOSTILLE

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- |                              |                                                                                |
|------------------------------|--------------------------------------------------------------------------------|
| 1. Country                   | Australia                                                                      |
| This public document         |                                                                                |
| 2. has been signed by        | John Howard Fisher                                                             |
| 3. acting in the capacity of | Notary Public                                                                  |
| 4. bears the seal/stamp of   | John Howard Fisher, Notary Public                                              |
| Certified                    |                                                                                |
| 5. at Sydney Passport Office | 6. the 15th day of July, 2019                                                  |
| 7. by Milan Stekovic         | Department of Foreign Affairs and Trade<br>Sydney Passport Office<br>Australia |
| 8. No. UPPT-F3-94278         |                                                                                |
| 9. Seal/Stamp                | 10. Signature                                                                  |

*This Apostille only certifies the authenticity of the signature (where applicable) and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. This Apostille can be verified at <https://orao.dfat.gov.au/pages/verifyapostille.aspx>*

**MEDRAD® Salient**

**MEDRAD® Vistron Plus**

Contrast Injection System

## Disposable sets and Connector Tubing INSTRUCTIONS FOR USE



| Catalog No. | Description                                                                           |
|-------------|---------------------------------------------------------------------------------------|
| ZY6320      | Salient Syringe 190 ml and quick fill tube                                            |
| ZY6321      | Salient Syringe 190 ml and spike                                                      |
| ZY6322      | Salient Syringe 190 ml and quick fill tube                                            |
| ZY6323      | Salient Syringe 190 ml and spike                                                      |
| ZY6324      | Salient Syringe 190 ml and spike<br>150 cm 300 psi coiled tube set – single           |
| ZY6325      | Salient Syringe 190 ml and quick fill tube<br>150 cm 300 psi coiled tube set – single |
| ZY5151      | 150 cm 300 psi coiled tube set – single                                               |
| ZY5152      | 150 cm 300 psi coiled tube set – dual                                                 |

©Copyright MN090002-10 Rev Date 18 Jan 2019  
Imaxeon Pty Ltd. All Rights Reserved.  
Imaxeon® and FluiDots® are registered trademarks.

| Language         | Page |
|------------------|------|
| English          | 3    |
| Deutsch          | 3    |
| Français         | 4    |
| Italiano         | 4    |
| Ελληνικά         | 5    |
| Svenska          | 5    |
| Nederlands       | 6    |
| Dansk            | 6    |
| Norsk            | 7    |
| Português-Brasil | 7-8  |
| Español          | 8    |
| 日本語              | 9    |
| Slovenščina      | 9    |
| Български        | 10   |
| Srpski           | 10   |
| Polski           | 11   |
| Türkçe           | 11   |
| Magyar           | 12   |
| Hrvatski         | 12   |
| 简体中文             | 13   |
| Русский          | 14   |
| Português        | 15   |

Imaxeon Pty Ltd  
Unit 1, 38-46 South Street  
Rydalmere NSW 2116  
Australia  
T: +61-2-8845-4999  
F: +61-2-8845-4936  
W: www.imaxeon.com

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**Introduction:**

Read the information contained in this section. Understanding the information will assist you in operating the device in a safe manner. **Important Safety Notice:** This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

**Intended Use:** The contents of this package are intended to be used in the delivery of contrast media or saline. They are indicated for single-use on one patient only with Salient and Vistron Plus Injectors. The transfer set is intended for use for one container of media only, and must be discarded when the container becomes empty, or after six hours, whichever comes first. Refer to the media manufacturer's instructions for use for additional indications.

**Contraindications:** These devices are not intended for multiple patient use, drug infusion, chemotherapy, or any other use for which the device is not indicated.

**Restricted Sale:** Rx Only

**Warnings**

- Air embolization can cause death or serious injury to the patient. Do not connect a patient to the injector until all trapped air has been cleared from the syringe and fluid path. Carefully read the instructions for loading and the use of FluidDots® indicators to reduce the chance of air embolism.
- Biological contamination can result from reusing disposable items or failure to follow aseptic technique. Properly discard disposable items after use, or if there is any possibility that contamination may have occurred.
- Syringe sterility will be compromised, and patient infection may result, if the plunger is removed from the syringe. Do not remove the plunger to fill the syringe.
- Bacterial contamination can occur if syringes are used to store contrast media. Use loaded syringes immediately. Discard unused loaded syringes.
- Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.
- Patient or operator injury may result from contrast media leaks or tubing ruptures. Ensure that the fluid path is open; do not exceed pressures identified on the front of the package. Use of greater pressures or occlusions in the fluid path may result in leaks or ruptures.
- Patient injury could result if syringe is not properly engaged. Do not load or inject unless the syringe is properly engaged.

**Cautions**

- Component damage or leaks may occur if not installed properly. Ensure all connections are secure; do not overtighten. This will help minimize leaks, disconnection, and component damage.

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
- Refer to the Injector Operation Manual for further instructions. **NOTE:** To minimize bubbles when loading a syringe through a transfer set, do not exceed a reverse load rate of 7 ml/sec.

**Installing a Syringe**

Insert the syringe until it snaps into place. The piston automatically engages the plunger, and advances it to the front of the syringe. Operator vigilance and care, coupled with a set procedure, is essential to minimizing the possibility of an air embolism. Point the injector head up during loading. Point the injector head down during an injection.

To help avoid an air injection, IMAXEON syringes are equipped with FluidDots® indicators. FluidDots indicators should be observed as part of an arming procedure. When the FluidDots are viewed through an empty syringe, the dots appear as small narrow ellipses. When viewed through a full syringe, the dots become larger, almost round.

To minimize air embolization risks, ensure that one operator is designated the responsibility of filling the syringe(s). Do not change operators during the procedure. If an operator change must occur, ensure that the new operator verifies that the fluid path is purged of air.

**Loading and Priming a Syringe Using a Quick Fill Tube (QFT)**

The syringe can be loaded manually or automatically on the Salient injector. See the Salient Operation Manual for more detailed instructions.

**Salient Manual Loading:**

1. Install a new syringe.
2. Install the QFT onto the end of the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
3. Insert the QFT into the fluid source and fill the syringe with fluid, using the load buttons.
4. Expel air.
5. Remove QFT. Connect the disposable tubing set.
6. Follow the instructions in the section "Connector Tube Installation."

**Salient Integral Autoload with a QFT:**

1. Install a new syringe.
2. Select the syringe on the fill screen and press the "Auto" Fill tab
3. Install a QFT onto the end of the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
4. Insert the QFT into the fluid source and press the start button on the "Auto" Fill tab.
5. Expel air.

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
6. Remove QFT. Connect the disposable tubing set.
  7. Follow the instructions in the section "Connector Tube Installation."
  8. Press Next button to set the injection parameters
- Note:** Patency check through aspiration using the injector is not possible when using a check valve on the end of the connector tube. If aspiration is important, remove the check valve from the connector tube and directly connect the connector tube to the catheter.

**Filling with Spike**

Prepare the bottle or bag of contrast and/or flushing solution. Push the bottle or bag onto the spike until the seal is punctured.

**Removing a Syringe(s)**

1. Disconnect the disposable tubing set from the vascular entry device. The disposable tubing set does not need to be disconnected from the syringe.
2. Rotate the syringe approximately 1/4 turn counter-clockwise and gently pull the syringe out of the injector head, discarding the syringe with disposable tubing set.

**Note:** Once the syringe is removed from the injector, the piston will automatically retract.

**Note:** In order to remove the syringe, the last piston motion must be in the forward direction, which is typical. If you cannot remove the syringe, press the "Back" button until the fill screen is displayed, press the retract button, then repeat Step 2.

**Connector Tube Installation**

1. Ensure all air is purged from the syringe.
2. Remove the connector tube from the package. Remove the dust covers on the luer fittings. Attach the connector tube to the syringe, 1/4 turn to 1/2 turn maximum. Do not install with excessive force.
3. Ensure that the connector luer fitting is secured to the tip of the syringe and verify that the tubing is not kinked or obstructed.

If you are using a connector tube with T-connector, attach the straight portion of the T-connector to the contrast (Syringe A) and the extension to the saline (Syringe B). If the T-connector is connected to Syringe B, priming will not completely fill the connector tubing.

4. Ensure all air is purged.
5. Rotate the injector head downward.
6. Connect to the patient and press ARM.
7. The injector will ask for CHECK FOR AIR.

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