Serial No.5048-03-23.

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Hereby certify that I was presented with:

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- ☒ the original document
- ☐ a true copy of the original document
- ☐ a copy that is not the original and is not a true copy of the original
- ☐ computerized information (please specify) from the database of
- Drawn up in the languageENGLISH in which I am fluent.

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I confirm, with my signature and seal, that the attached document marked with the letter/ numberA. is an exact copy of the document or information presented to me.

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Alpha-Bio Tec Ltd.

Regulatory Affairs Project Manager

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Instruction for use	
Dental implants Alpha-Bio	

1. MEDICAL DEVICE NAME

Dental Implants Alpha-Bio:

I. Dental implants Dual Fit Implant (DFI) with internal hex connection (IH):

1. DFI IH diameter 3,75 mm, length 8,0 mm.
2. DFI IH diameter 3,75 mm, length 10,0 mm.
3. DFI IH diameter 3,75 mm, length 11,5 mm.
4. DFI IH diameter 3,75 mm, length 13,0 mm.
5. DFI IH diameter 4,2 mm, length 8,0 mm.
6. DFI IH diameter 4,2 mm, length 10,0 mm.
7. DFI IH diameter 4,2 mm, length 11,5 mm.
8. DFI IH diameter 4,2 mm, length 13,0 mm.
9. DFI IH diameter 3,3 mm, length 8,0 mm.
10. DFI IH diameter 3,3 mm, length 10,0 mm.
11. DFI IH diameter 3,3 mm, length 11,5 mm.
12. DFI IH diameter 3,3 mm, length 13,0 mm.
13. DFI IH diameter 5,0 mm, length 8,0 mm.
14. DFI IH diameter 5,0 mm, length 10,0 mm.
15. DFI IH diameter 5,0 mm, length 11,5 mm.
16. DFI IH diameter 5,0 mm, length 13,0 mm.

II. Dental implants Dual Fit Implant (DFI) with conical standard connection (CS):

1. DFI CS diameter 3,75 mm, length 8,0 mm.
2. DFI CS diameter 3,75 mm, length 10,0 mm.
3. DFI CS diameter 3,75 mm, length 11,5 mm.
4. DFI CS diameter 3,75 mm, length 13,0 mm.
5. DFI CS diameter 4,2 mm, length 8,0 mm.
6. DFI CS diameter 4,2 mm, length 10,0 mm.
7. DFI CS diameter 4,2 mm, length 11,5 mm.
8. DFI CS diameter 4,2 mm, length 13,0 mm.
9. DFI CS diameter 5,0 mm, length 8,0 mm.
10. DFI CS diameter 5,0 mm, length 10,0 mm.
11. DFI CS diameter 5,0 mm, length 11,5 mm.
12. DFI CS diameter 5,0 mm, length 13,0 mm.

Further in the text, abbreviated names of registered medical devices will be used: Implant, DFI implant.

2. DESCRIPTION

Dental Implants Alpha-Bio are made of biocompatible Titanium alloy (Ti-6Al-4V ELI). The implants are treated with NanoTec surface treatment. Cover screw is delivered with implant package. The Cover screw is used to seal the implant concealed under the soft tissue for the osteointegration period. The cover screws are made of Titanium alloy (Ti-6Al-4V ELI). Cover screws are color coded in accordance with the implant's connection.

Instruction for use	
Dental implants Alpha-Bio	

Table 1. Color code for cover screw

IMPLANT TYPE	IMPLANT CONNECTION	COVER SCREW
DFI	Internal Hex (IH)	Blue
DFI	Conical Standard (CS)	Green

Implants package identification labels indicates implant type, diameter and length are coded by color

3. INTENDED USE

Designed to replace the tooth/teeth root and for supporting fixed or removable dentures in order to improve and to restore chewing function.

4. APPLICATION AREA

Dentistry. Intended for professional use only.

5. INDICATIONS

Endosseous implants are indicated for the reconstruction of deficient dental apparatus (tooth/teeth) by acting as a supporting bone fixed device for dental prosthesis in the lower or upper jaw. The dental prosthesis may be fixed or a removable, single or multiple unit. Implant can be used in immediate, one or two stage procedures, assuming an adequate primary stability is achieved.

The cover screw is indicated for use with implants, with its corresponding implant connection.

6. CONTRAINDICATIONS

It is contraindicated to use Alpha-Bio dental implants and cover screws in:

- Patients who are medically unfit for an oral surgery.
- Patients in whom a physical limitation does not allow a safe support of functional/parafunctional loads (limited implants size, numbers and/or positions)
- Patients who are allergic or hypersensitive to the following materials: titanium alloy Ti-6Al-4V ELI (titanium, aluminum, vanadium).

7. WARNINGS

Over tightening of the implant can result in:

- Clinical damage to the implantation site in the short and/or long run.

Instruction for use	
Dental implants Alpha-Bio	

- Breakage and/or damage to the implant.

It is recommended to follow carefully the drilling protocol as described according to implants and drill type to avoid over or under drilling.

Failure to recognize the proper length and size of osteotomies may lead to wrong drilling depth, resulting in permanent injury to nerves or other vital structures. Drilling beyond the depth intended for lower jaw surgery may potentially result in permanent numbness to the lower lip and chin or lead to a life-threatening hemorrhage in the floor of the mouth.

Besides the mandatory precautions for any surgery such as asepsis, during drilling in the jawbone, one must avoid damage to nerves and vessels by referring to anatomical knowledge and preoperative radiographs.

Because of the small size of the cover screw, care must be taken that it will not be swallowed or aspirated by the patient. It is appropriate to use specific supporting tools to prevent aspiration of loose parts (e.g. a throat shield). It is recommended to use the prosthetic “grip” driver when placing the screw in the implant.¹

8. CAUTION

General:

One hundred percent implant success cannot be guaranteed. Especially, nonobservance of the indicated limitations of use and working steps may result in failure. Treatment by means of implants may lead to loss of bone, biologic or mechanical failures including fatigue fracture of implants.

Close cooperation between surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment.

It is strongly recommended that Alpha-Bio Tec. surgical instruments and prosthetic components are used only with Alpha-Bio Tec. implants, as combining components that are not dimensioned for correct mating can lead to mechanical and/or instrumental failure, damage to tissue or unsatisfactory aesthetic results.

It is strongly recommended that clinicians, new as well as experienced users of implants, prosthetics, and associated software, always go through special training before undertaking a new treatment method. Alpha-Bio Tec. offers a wide range of courses for various levels of knowledge and experience. For more information, please visit www.alpha-bio.net.

When using a new device/treatment method for the first time, working with a colleague who is experienced with the new device/treatment method may help avoid possible complications. Alpha-Bio Tec. has a global network of mentors available for this purpose. It is especially important to achieve proper stress distribution through adaptation and fitting of the crown or bridge by adjusting occlusion to the opposing jaw. In addition, avoid excessive transverse loading forces, particularly in immediate loading cases.

¹ RC № FSZ 2011/10261(ref. 4052 / 4053)

Instruction for use	
Dental implants Alpha-Bio	

Before Surgery:

Careful psychological, evaluation followed by clinical and radiological examination must be performed to the patient prior to surgery to determine the suitability of the patient for treatment. Special attention must be given to patients who have local or systemic factors that could interfere with the healing process of either bone or soft tissue or the osseointegration process (e.g., cigarette smoking, poor oral hygiene, uncontrolled diabetes, orofacial radiotherapy, steroid therapy, infections in the neighboring bone or other situations that may jeopardize success as described in the literature special caution is advised in patients who receive bisphosphonate therapy).

In general, implant placement and prosthetic design must accommodate individual patient conditions. In case of bruxism, other parafunctional habits or unfavorable jaw relationships, reappraisal of the treatment option may be considered.

With respect to pediatric patients, routine implant treatment is not recommended till complete growth of facial, TMJ, jaws and cranium.

Pre-operative hard tissue or soft tissue deficits may yield a compromised esthetic result or unfavorable implant angulations.

All components, instruments and tooling used during procedure must be maintained in good condition and care must be taken that instrumentation does not damage implants or other components.

At Surgery:

Care and maintenance of instruments are crucial for a successful treatment. Sterilized instruments not only safeguard your patients and staff against infection but are also essential for the outcome of the total treatment.

Because of the small sizes of the devices, care must be taken that they are not swallowed or aspirated by the patient.

The number of implants and angulation are introduced in professional literature and the user should follow these recommendations.

After implant installation, the surgeon's evaluation of bone quality and primary stability will determine when implants may be loaded. Lack of adequate quantity and/ or quality of remaining bone, infection and generalized diseases may be potential causes for failure of osseointegration both immediately after surgery, or after osseointegration is initially achieved.

Bending movements may create forces that are known to be most unfavorable, as they can potentially jeopardize the long-term stability of an implant-supported restoration. In order to decrease bending movements, the distribution of forces should be optimized by cross-arch stabilization, minimizing distal cantilevers, having a balanced occlusion as well as decreased cuspal inclination of the prosthetic teeth and should be very well balanced.

While modifying the restoration, use copious irrigation and appropriate protection equipment. Avoid inhalation of dust created while performing the above-mentioned procedure.

After Surgery:

To help secure the long term-treatment outcome, it is advised to provide comprehensive patient

Instruction for use	
Dental implants Alpha-Bio	

follow up after implant treatment and to inform about appropriate oral hygiene.

9. POTENTIAL USERS

Potential consumers of the medical device are male and female patients with complete or partial edentulism with a properly developed and documented growth phase of the jawbone.

The medical devices are intended for use by a qualified dentist who has completed a professional course in implant surgery.

Procedure conditions:

The product is intended for use in a medical institution (in dental clinics), in a professional, equipped dental office.

10. SIDE EFFECTS

Placement of dental implants and restorative components is an invasive surgical procedure that may be associated with the following residual risks: pain, swelling, infection, bleeding, hematoma and inflammation.

Drilling into the jawbone or misplacing an implant may (in rare cases) result in fenestration/fracture of bone, damage to neighboring structures, nerve damage causing temporary or permanent loss of sensation of the lower lip and/or tongue, damage to blood vessels, perforation of nasal cavity and/ or maxillary sinus and local and systemic infection depending on drilling and implant location.

During the placement of an implant and use of its devices, the pharyngeal (gag) reflex may be triggered in sensitive patients.

Dental implants are part of a multi-component system that replaces tooth structures. As a result, side effects similar to those occurring in natural teeth might develop overtime around replaced teeth such as but not limited to mucositis, calculus, peri- implant it is, fistulas, ulcers, soft tissue hyperplasia, soft/hard tissue recession/loss. Some patients may experience mucosal discoloration. During submerged healing period, bone may grow over the cover screw. In some cases, the cover screw may become prematurely exposed.

During soft tissue healing period, gingiva may grow over the healing abutment.

11. NOTICE REGARDING SERIOUS INCIDENTS

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the authorized representative of the manufacturer in your country and to your national authority.

The contact information to report a serious incident is as follows:

Nobel Biocare Russia LLC, Moscow, st. Stanislavsky, 21, building 2, 109004;

Contact phone: +7 (495) 974-77-55.

Website: www.nobelbiocare.ru

Instruction for use	
Dental implants Alpha-Bio	

12. SURGICAL PROCEDURE AND HANDLING PROCEDURE


Drilling:

The implant site is prepared in a sequential procedure, using drills of increasing diameters with depth indication laser marking that provide a reading of the desired drilling depth. Drills should be replaced when their cutting efficiency is reduced. All preparation of bone tissue must be carried out under ample irrigation, with saline solution and using an intermittent drilling technique. The drilling protocol should be adapted to the implant diameter and bone quality at the site.


The instructed drilling steps are suggested for dense/cortical bone. In low density bone, fewer steps may be performed. Drill sequence is recommended according to bone type and implant type.

13. PACKAGE. SET OF DELIVERY.


Product packaging is a modern and easy-to-use package designed for maximum comfort, absolute sterility and superiority ergonomics.

Picture 1. Primary MD package

Instruction for use	
Dental implants Alpha-Bio	

Picture 2. Implant removing from package

1. Open the package - tear the cardboard, pull the tyvek, remove the inner holder.
2. Open the cap – this can be accomplished with an easy on hand operation.
3. Insert the preferred tool/driver (connected either to a ratchet or a contra-angle, or manually) in a continuous clockwise turn until hex is found.
4. Implant is now connected and prepared to be removed from the package. Do it carefully while verifying that there is no gap between the driver and the implant's hex. Implant may be inserted directly to site using the preferred driver (ratchet driver, motor mount driver, or surgical screwdriver).
5. Implant installation should continue until desired position is achieved, using ABT's recommended torque maximum of 35Ncm.
6. For Immediate loading, use the appropriate torque decided by the surgeon. The implant should withstand occlusal forces. Depending on chosen surgical protocol, place a cover screw or an abutment

Set of delivery

The implant together with the cover screw, packed in an individual (primary) package, is delivered to the consumer in a group cardboard box. Instructions for the use of products are provided to the consumer upon delivery by an authorized representative upon request.

14. MATERIALS

Dental implants and cover screws: titanium alloy Ti-6Al-4V ELI (titanium, aluminum, vanadium)

15. STERILITY AND REUSABILITY INFORMATION

Dental implants and cover screws are delivered sterile for single use only. Do not use after the labeled expiration date.

Instruction for use	
Dental implants Alpha-Bio	

The sterilization method is radiation.

Caution: Do not use device if the packaging has been damaged or previously opened.

Caution: Dental implants and cover screws are single use products and should not be reprocessed. Reprocessing could cause loss of mechanical, chemical and / or biological characteristics. Reuse could cause infection.

Caution: Do not use device if the packaging has been damaged or previously opened.

16. PERFORMANCE REQUIREMENTS AND LIMITATIONS

To achieve the desired performance of Alpha-Bio Tec.'s dental implants, cover screws and healing abutments the products must be used according to the Instructions for use and/or in accordance with other Alpha-Bio Tec. components and the Instructions for Use of each product.

To confirm the compatibility of other Alpha-Bio Tec. products with the dental implants, cover screws and healing abutments check the color coding, dimensions, connection type and/or any other applicable marking on the product/product package.

17. IMPORTANT – DISCLAIMER OF LIABILITY

This product is part of an overall concept and may only be used in conjunction with the associated original products according to the instructions and recommendation of Alpha-Bio Tec. (ABT). Non-recommended use of products made by third parties in conjunction with Alpha-Bio Tec. products will void any warranty or other obligation, express or implied, of Alpha-Bio Tec. The user of Alpha-Bio Tec. products has the duty to determine whether or not any product is suitable for the particular patient and circumstances. Alpha-Bio Tec. disclaims any liability, express or implied, and shall have no responsibility for any direct, indirect, punitive or other damages, arising out of or in connection with any errors in professional judgment or practice in the use of Alpha-Bio Tec products. The user is also obliged to study the latest developments in regard to this Alpha-Bio Tec. product and its applications regularly. In cases of doubt, the user has to contact Alpha-Bio Tec. Since the utilization of this product is under the control of the user, they are his/her responsibility.

Alpha-Bio Tec. does not assume any liability whatsoever for damage arising thereof.

18. OPERATING, TRANSPORT AND STORAGE CONDITIONS. SHELF LIFE. PATIENT CONTACT.

The device must be stored and transported in dry conditions in the original packaging at room temperature and not exposed to direct sunlight. Incorrect storage and transportation may influence device characteristics leading to failure.

The shelf life (sterility) of implants is 5 years.

Transportation of products should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport.

Patient contact

Instruction for use	
Dental implants Alpha-Bio	

Product of long-term contact (more than 30 days). Implantable product in contact with bone and soft tissues.

19. THE PROCEDURE AND CONDITIONS FOR THE DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.

Safely discard potentially contaminated or no longer usable medical devices as healthcare (clinical) waste in accordance with local healthcare guidelines, country and government legislation or policy.

Separation, re-cycling or disposal of packaging material shall follow local country and government legislation on packaging and packaging waste, where applicable.

Implant is delivered in a sterile, single use package (including cover screw) and intended for single use.

Medical devices in the Russian Federation must be disposed of in local public/private organizations in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of the territories of urban and rural settlements, for water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures.


If the integrity of the sterile packaging is damaged, as well as when the expiration date has expired, the product is disposed of in accordance with the rules for handling epidemiologically safe medical waste, close in composition to municipal solid waste (class A).

The product that has come into contact with the patient's body is disposed of in accordance with the rules for handling epidemiologically hazardous medical waste (class B).

The medical product does not adversely affect the human body and the environment during operation, storage and transportation in accordance with the instructions for use.

20. SYMBOLS DESCRIPTION

Table 2. Symbols on the package description

SYMBOL	DESCRIPTION
	Sterile. Sterilized using irradiation
	Do not reuse
	Caution!

Instruction for use	
Dental implants Alpha-Bio	
	Consult instructions for use
	Do not use if package is damaged
	Use by date
UDI	<i>unique device identification</i>
LOT	Batch code
REF	Catalogue number
	Manufacturing date
	Manufacturer
Rx only	Federal (USA) law restricts the sale of this device to, or on the order of, a licensed physician or dentist
Made in	Made in

21. TECHNICAL REQUIREMENTS OF MEDICAL DEVICE

The main technical characteristics of the medical device "Alpha Bio Dental Implants" are presented in Table 3.

The deviation of values not shown on the drawings is: ± 0.05 mm (for values from 0.5 mm to 6 mm) and ± 0.1 mm (for values from 6 mm to 30 mm) according to DIN ISO 2768.

Instruction for use	
Dental implants Alpha-Bio	

Table 3. Main technical characteristics and materials for the manufacture of the medical product "Dental Implants Alpha Bio

Device name	Device technical characteristics
I. Dental implant Dual Fit Implant (DFI) with internal hex connection (IH):	
1. DFI IH diameter 3,75 mm, length 8,0 mm.	Mass:0.23±0.03 gr Diameter:3.85±0.02mm General length:8±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
2. DFI IH diameter 3,75 mm, length 10,0 mm.	Mass:0.30±0.03 gr Diameter:3.85±0.02mm General length:10±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
3. DFI IH diameter 3,75 mm, length 11,5 mm.	Mass:0.37±0.03 gr Diameter:3.85±0.02mm General length:11.5±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
4. DFI IH diameter 3,75 mm, length 13,0 mm.	Mass:0.43±0.03 gr Diameter:3.85±0.02mm General length:13±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
5. DFI IH diameter 4,2 mm, length 8,0 mm.	Mass:0.30±0.03 gr Diameter:4.2±0.02mm General length:8±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
6. DFI IH diameter 4,2 mm, length 10,0 mm.	Mass:0.39±0.03 gr Diameter:4.2±0.02mm General length:10±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
7. DFI IH diameter 4,2 mm, length 11,5 mm.	Mass:0.45±0.03 gr

Instruction for use	
Dental implants Alpha-Bio	

	Diameter: 4.2 ± 0.02 mm General length: 11.5 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
8. DFI IH diameter 4,2 mm, length 13,0 mm.	Mass: 0.54 ± 0.03 gr Diameter: 4.2 ± 0.02 mm General length: 13 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
9. DFI IH diameter 3,3 mm, length 8,0 mm.	Mass: 0.20 ± 0.03 gr Diameter: 3.7 ± 0.02 mm General length: 8 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
10. DFI IH diameter 3,3 mm, length 10,0 mm.	Mass: 0.27 ± 0.03 gr Diameter: 3.7 ± 0.02 mm General length: 10 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
11. DFI IH diameter 3,3 mm, length 11,5 mm.	Mass: 0.32 ± 0.03 gr Diameter: 3.7 ± 0.02 mm General length: 11.5 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
12. DFI IH diameter 3,3 mm, length 13,0 mm.	Mass: 0.38 ± 0.03 gr Diameter: 3.7 ± 0.02 mm General length: 13 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
13. DFI IH diameter 5,0 mm, length 8,0 mm.	Mass: 0.43 ± 0.03 gr Diameter: 4.95 ± 0.02 mm General length: 8 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm

Instruction for use	
Dental implants Alpha-Bio	

14. DFI IH diameter 5,0 mm, length 10,0 mm.	Mass:0.53±0.03 gr Diameter:4.95±0.02mm General length:10±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
15. DFI IH diameter 5,0 mm, length 11,5 mm.	Mass:0.63±0.03 gr Diameter:4.95±0.02mm General length:11.5±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7 mm, Cover screw (IH) length 4,9 mm
16. DFI IH diameter 5,0 mm, length 13,0 mm.	Mass:0.73±0.03 gr Diameter:4.95±0.02mm General length:13±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (IH) diameter 3,7mm, Cover screw (IH) length 4,9mm
II. Dental implant Dual Fit Implant (DFI) with conical standard connection (CS):	
1. DFI CS diameter 3,75 mm, length 8,0 mm.	Mass:0.22±0.03 gr Diameter:3.85±0.02mm General length:8±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
2. DFI CS diameter 3,75 mm, length 10,0 mm.	Mass:0.29±0.03 gr Diameter:3.85±0.02mm General length:10±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
3. DFI CS diameter 3,75 mm, length 11,5 mm.	Mass:0.35±0.03 gr Diameter:3.85±0.02mm General length:11.5±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
4. DFI CS diameter 3,75 mm, length 13,0 mm.	Mass:0.42±0.03 gr Diameter:3.85±0.02mm General length:13±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm

Instruction for use	
Dental implants Alpha-Bio	

5. DFI CS diameter 4,2 mm, length 8,0 mm	Mass:0.29±0.03 gr Diameter:4.2±0.02mm General length:8±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
6. DFI CS diameter 4,2 mm, length 10,0 mm.	Mass:0.37±0.03 gr Diameter:4.2±0.02mm General length:10±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
7. DFI CS diameter 4,2 mm, length 11,5 mm.	Mass:0.44±0.03 gr Diameter:4.2±0.02mm General length:11.5±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
8. DFI CS diameter 4,2 mm, length 13,0 mm.	Mass:0.53±0.03 gr Diameter:4.2±0.02mm General length:13±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
9. DFI CS diameter 5,0 mm, length 8,0 mm.	Mass:0.42±0.03 gr Diameter:4.95±0.02mm General length:8±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
10. DFI CS diameter 5,0 mm, length 10,0 mm.	Mass:0.52±0.03 gr Diameter: 4.95±0.02mm General length:10±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15mm, Cover screw (CS) length 5,5mm
11. DFI CS diameter 5,0 mm, length 11,5 mm.	Mass:0.62±0.03 gr Diameter:4.95±0.02mm General length:11.5±0.1mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) length 5,5mm

Instruction for use	
Dental implants Alpha-Bio	

12. DFI CS diameter 5,0 mm, length 13,0 mm.	Mass: 0.71 ± 0.03 gr Diameter: 4.95 ± 0.02 mm General length: 13 ± 0.1 mm Hardness: no less than 30 HRC Roughness: no more than 0.8 micro-meters Cover screw (CS) diameter 3,15 mm, Cover screw (CS) length 5,5 mm
---	--

Insertion torques should not exceed 50Ncm and shall be minimum 35Ncm in immediate loading procedures.

22. WARRANTY OBLIGATIONS

To achieve the desired performance of Alpha-Bio Tec.'s dental implants, cover screws and healing abutments the products must be used according to the Instructions for Use for each Alpha-Bio Tec product.

To confirm the compatibility of other Alpha-Bio Tec. products with the dental implants, cover screws and healing abutments check the color coding, dimensions, connection type and/or any other applicable marking on the product/product package.

For questions about the quality of a medical product, manufactured by Alpha-Bio Tec Ltd. please contact:

Authorized representative on the territory of the Russian Federation:

LLC «Nobel Biocare Russia»,

Russian Federation, 109004 Moscow, Stanislavskogo str., 21, bld. 2,

Telephone: +7 (495) 974-77-55 fax +7 (495) 974-77-66

Web-site: <https://www.nobelbiocare.com>

23. REQUIREMENTS FOR THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES

Not provided. Products are sterile single use.

24. LIST OF APPLIED STANDARDS

EN ISO 13485:2016/AC:2018	Medical devices — Quality management systems — Requirements for regulatory purposes BS EN ISO 13485:2016 + A11:2021 / ISO 13485:2016
FDA QSR	21 CFR 820 Quality System regulation for medical device federal register
EN 1641:2009	Dentistry - Medical devices for dentistry - Materials
EN 1642:2011	Dentistry - Medical devices for dentistry - Dental implants
EN ISO 14971:2019/A11:2021	Medical devices - Application of risk management to medical devices
ASTM F2503-20	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment
EN 62366-1:2015 + AC:2015 A1:2020	Medical devices - Part 1: Application of usability engineering to medical devices (BS EN 62366-1:2015+A1:2020 / IEC 62366:2015)
MDD 93/42/EEC with	Medical Device Directive

Instruction for use	
Dental implants Alpha-Bio	

amendment Directive 2007/47/E	
EN ISO 14155: 2020	Clinical investigation of medical devices for human subjects - Good clinical practice (BS EN ISO 14155:2020 / ISO 14155:2020)
MedDev 2.7/1 rev 4 June 2016	Clinical evaluation: Guide for manufacturers and notified bodies
MedDev 2.12-1 rev 8 Jan 2013	Guidelines on a medical devices vigilance system
MedDev 2/12 rev 2:2012 Januar 2012	Guidelines on postmarketing clinical testings
MedDev 2.5/5 rev 3:1998	Translation procedure
ISO 10451:2010	Dentistry — Contents of technical file for dental implant systems
ASTM F136-13 2013	Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401)
EN ISO 10993-1 2020	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process (BS EN ISO 10993-1:2020 /ISO 10993-1:2018)
EN ISO 7405:2018	Dentistry. Evaluation of biocompatibility of medical devices used in dentistry
EN ISO 14698-1 2003	Cleanrooms and associated controlled environments - Biocontamination control - Part 1: General principles and methods (ISO 14698-1:2003)
EN 11137-1:2015 + A2:2019	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices (BS EN ISO 11137-1:2015+A2:2019 / ISO 11137-1:2006+Amd 1:2013)
ISO 11137-2:2013 + Amd1:2022	Sterilization of health care products — Radiation — Part 2: Establishing the sterilization dose
EN ISO 11737-1:2018+ A1:2021	Sterilization of health care products - Microbiological methods Part 1: Determination of a population of microorganisms on products – (BS EN ISO 11737-1:2018+A1:2021 /ISO 11737-1:2018)
EN ISO 11737-2: 2020	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process (BS EN ISO 11737-2:2020 /ISO 11737-2:2019).
EN ISO 14664-1:2015	Cleanrooms and associated controlled environments — Part 1: Classification of air cleanliness by particle concentration
EN ISO 14664-2:2015	Cleanrooms and associated controlled environments. Part 2. Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration
ISO/TS 13004:2022	Sterilization of health care products-Radiation-Substantiation of a selected sterilization dose—Method Vdmax.
EN ISO 14937:2009	Sterilization of health care products — General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices
EN ISO 11607-1: 2020	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

Instruction for use	
Dental implants Alpha-Bio	

	(BS EN ISO 11607-1:2020 / ISO 11607-1:2019).
EN ISO 11607-2: 2020	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes (BS EN ISO 11607-1:2020 / ISO 11607-1:2019)
ASTM F 1929-15	Standard Test Method for Detecting Seal Leaks in Porous Medical Packaging by Dye Penetration
ASTM F1980-21	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices.
EN 868-5:2018	Packaging for terminally sterilized medical devices. Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods.
ASTM F2096 (2019)	Standard Test Method for Detecting Gross Leaks in Packaging by Internal Pressurization (Bubble Test).
ASTM D4169 2022	Standard Protocol for Performance Testing of Shipping Containers and Systems
ASTM F 1929-12	Standard Test Method For Detecting Seal Leaks In Porous Medical Packaging By Dye Penetration
EN ISO 20417-1: 2021	Medical devices — Information to be supplied by the manufacturer
EN ISO 15223-1:2021	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements
EN 556-1:2001/AC:2006	Sterilization of medical devices. Requirements for medical devices to be designated "STERILE". Requirements for terminally sterilized medical devices.
ISO 14801:2016	Dentistry — Implants — Dynamic fatigue test for endosseous dental implants.
ISO 2859-1:1999	Sampling procedures for inspection by attributes. Part 1: Sampling schemes indexed by acceptance quality limit for lot-by-lot inspection

25. THE NAME AND LEGAL ADDRESS OF THE ORGANIZATION-DEVELOPER, PRODUCER OF MEDICAL DEVICES AND AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER:

Developer:

Alpha-Bio Tec Ltd.

7 Hatnufa St., Kiryat Arye, POB 3936, Petach Tikva, 4951025, Israel

Manufacturer:

Alpha-Bio Tec Ltd.

7 Hatnufa St., Kiryat Arye, POB 3936, Petach Tikva, 4951025, Israel

Manufacturing site:

Alpha Bio Tec Ltd.,

Instruction for use	
Dental implants Alpha-Bio	

21 HaReches Blvd. POB 264 Modi'in 7171102, Israel

Authorized representative of the manufacturer in the territory of the Russian Federation:

Nobel Biocare Russia LLC, Moscow, st. Stanislavsky, 21, building 2, 109004;

Tel.: +7(495)974 77 55, Fax: +7(495)974 77 66

Date of issue or last revision of the instructions for use

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