**Milestone Srl**


Headquarter
Via Fatebenefratelli, 1/5
24010 Sorisole (BG) - Italy
Tel +39 035 4128264
Fax +39 035 575498

Valbrembo Plant
Corso Europa Unità, 6
24030 Valbrembo (BG) - Italy
Tel +39 035 0401600

marketing@milestonemedsrl.com
www.milestonemed.com

C.F. - P. IVA - Reg. Imp. 01879330163
R.E.A. Bergamo 244134
Cap. Soc. € 1.000.000,00 int. vers.
UNI EN ISO 9001 - UNI EN ISO 13485 CERTIFIED

APPROVED
Milestone S.r.l., Italy

(Diego CORTESI, General Director)

2024-03-01

MILESTONE S.r.l.
Via Fatebenefratelli, 1/5
24010 SORISOLE (BG)
P. IVA 01879330163

For the Russian Federation

Attention to

Roszdraznadzor

In For the Russian Federation

Dear Sirs,

attached to this letter, there is the Operational Documentation of the Manufacturer for the FormSAFE in-vitro diagnostic medical device, professional use only.

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We certify that the signer is a duly authorized signatory on behalf of the
Company for any documents and acts used in trade with foreign countries.

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Ass.te Amm.vo Baiocchi

MILESTONE
H E L P I N G
P A T I E N T S

**FormSAFE container with fixative for collection,
transportation and storage of biological materials**

User Manual - MM127

Thank you for choosing one of our products and welcome to the ever more crowded club of Milestone laboratory instrumentation users.

We are sure that you will be completely satisfied with this new product entering your laboratory.

We invite you to read carefully this user manual and to keep it in reach for convenient and fast consulting.

For any clarification or request for assistance, please contact our representative in your country:

Authorized representative in the Russian Federation:

BioLine LLC

Legal address: Russia, 197022 Saint Petersburg, Municipal district Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, structure 3H/5H, room 208

Postal address: Russia, 197022, Saint Petersburg, 23 Professora Popova Street, Let. E

E-mail: main@bioline.ru

Tel: (812) 320-49-49

Service Hotline: 8-800-555-49-40.

or contact

Legal manufacturer:

Milestone S.r.l., Italy

Via Fatebenefratelli, 1/5, IT-24010 Sorisole (BG), Italy

Tel. +39.035.412 8264

Fax +39.035.575498

Web-site: www.milestonemed.com

E-mail: marketing@milestonemedsrl.com

Manufacturing site:

DIAPATH S.P.A., Italy

Via Savoldini, 71

IT-24057 Martinengo (BG)

Please read the user manual before using the product.

FormSAFE

Table of Contents

1. Name of the medical device	5
2. Description.....	5
2.1. <i>Composition of reagent in FormSAFE:</i>	6
2.2. <i>Delivery set of the product</i>	7
3. FormSAFE intended use	7
3.1. <i>FormSAFE intended use and indications</i>	7
3.2. <i>Contraindications and side effects</i>	8
3.3. <i>Type of specimen analyzed</i>	8
4. Information on potential consumers (personnel qualification)	8
5. Technical specifications of FormSAFE	8
6. Warnings and precautions for use.....	9
7. Principle of using.....	10
8. Information on the conditions required for the collection, processing and preparation of specimens, data on the stability of the analyzed specimens, including storage conditions and duration, restrictions on freezing (thawing) cycles	10
9. Installation, assembly, configuration, preparation for use	11
10. Room requirements	11
11. Operate with FormSAFE	11
12. List of materials that are required to work with FormSAFE container, but are not included in the medical device delivery set	12
13. Information needed to identify medical products in order to obtain a safe combination and (or) information on known restrictions on joint use of medical devices.....	12
14. Information on the presence in the medical product of drugs for medical use, materials of animal and (or) human origin	12
15. Storage and transportation conditions.....	12
16. Packaging.....	13
17. Labelling	15
18. Symbols used	18
19. Waste disposal.....	19
20. Responsibility of the Manufacturer.....	19
21. Repair and service information	19
22. Cleaning and disinfection	20


23. Warranty and shelf life.....	20
24. Information on the initial release and revision of the user manual.....	20
25. List of standards applied.....	20

FormSAFE**1. Name of the medical device**

Name of the medical device: FormSAFE container with fixative for collection, transportation and storage of biological materials (hereafter referred to as the FormSAFE container(s), FormSAFE container(s) with fixative, container(s), FormSAFE, medical device(s) or product).

2. Description


Milestone s.r.l. has introduced the concept of a two-layered liquid container to reduce the risk of formalin exposure to the staff.

- 60 ml container
- 30 ml 10 % Neutral Buffered Formalin, containing 4% formaldehyde (w/v)
- 6 ml Floating shield fluid

FormSAFE container provides an excellent alternative to traditional prefilled containers when personnel are concerned with formalin exposure during collection, transportation and storage.

A second fluid component maintains a position above the filled formalin and acts as a protective floating shield which prevents fixative fumes from escaping the container.

This floating shield fluid is immiscible with formalin, having a lower specific gravity, and covers the top surface to form a protective film. This proprietary fluid is nontoxic and inert to biospecimens: if the container is agitated, the floating layer quickly reforms.

Conventional pre-filled
container (toxic fumes)FormSAFE
(no toxic fumes)

The floating shield fluid protects the operator from the exposure to toxic fumes.

The FormSAFE container are designed to be as a primary fixative: the presence of a "biphasic" floating protection shield prevents the escape of harmful fumes of formalin preventing exposure to the user.

Medical product for in vitro diagnostics	<i>IVD medical device according to Regulation 2017/746/EU.</i>
Potential risk class	<i>I</i>
Type of the medical device in accordance with the Russian nomenclature classification	<i>259900</i>
Catalog number	<i>51425 - FormSAFE, box including N. 50 of 60 ml pre-filled containers 51424 – single FormSAFE pre-filled container (no orderable)</i>

The device is intended for clinical laboratory diagnostics.

The FormSAFE container with fixative is for only one-time use.

2.1. Composition of reagent in FormSAFE:

Name	CAS-No.	Concentration
Composition of fixative (30 ml)		
<i>Water</i>	<i>7732-18-5</i>	<i>>98%</i>
<i>Formaldehyde</i>	<i>50-00-0</i>	<i>4%</i>
<i>Methanol</i>	<i>67-56-1</i>	<i><2.5%</i>
Composition of floating shield fluid (6 ml)		
<i>Alkanes (Isoparaffin)</i>	<i>90622-58-5</i>	<i>>90%</i>

Refer to the Safety Data Sheet for more specification.

2.2. Delivery set of the product

FormSAFE container with fixative for collection, transportation and storage of biological materials

Composed of:

1. FormSAFE container with fixative for collection, transportation and storage of biological materials - 50 pcs/pack.
2. User Manual in paper - 1 pc.

3. FormSAFE intended use

3.1. FormSAFE intended use and indications

Intended use:

The FormSAFE container with fixative for collection, transportation and storage of biological materials is a disposable auxiliary medical device for in vitro diagnostics intended for collection, transportation and storage of human biological materials (biopsy and autopsy materials, including surgical specimens) for purpose of its further processing and preparation of a histological specimen for conducting microscopic analysis and various types of clinical laboratory studies of collected specimens of human biological materials.

The FormSAFE container is intended only for professional use in clinical diagnostic laboratories of healthcare institutions for in vitro diagnostics.

Functional purpose	Auxiliary product for in vitro diagnostics
Description of a target analyte	Not applicable
Specific pathology, condition or risk factor for the detection, identification or differentiation of which the IVD device is intended	It is used for pathological anatomical studies to detect pathological changes in human tissues.
Population, demographic aspects of the medical device use	The FormSAFE container is not intended for a specific group of patients. It is allowed to use the product for fixation and storage of biological tissues of all populations, patient groups of different age.
Information on analytical sensitivity (limit of detection), analytical specificity, diagnostic sensitivity and diagnostic specificity	Not applicable. FormSAFE container is not an analytical or diagnostic kit of reagents, but it is auxiliary product in diagnostics.
Clinical performance characteristics: diagnostic sensitivity and diagnostic specificity	
Biological reference interval of use	
Information on interfering substances or sample limitations that may interfere with study result	The result of the study may be affected by untimely placement of the specimen in a container, failure to comply with storage and transportation conditions. Biological material should be stored only in the fixative.
Description of measurement procedures, metrological traceability of calibrators and	Not applicable. The delivery set of FormSAFE container

control materials values	does not include any calibrators and control materials.
--------------------------	---

Indications for use: 10% Neutral Buffered Formalin used in the FormSAFE is intended as a general purpose fixative for all tissues to preserve cells and tissue components in a “life-like state” and to allow for the preparation for histopathology analysis.

3.2. Contraindications and side effects

Contraindications (limits) for use: it is not allowed to use expired products, a damaged container with loss of tightness, if a color of fixative is changed.

Side effects: not applicable.

3.3. Type of specimen analyzed

Types of analyzed specimens for fixation: all human tissues: biopsy and autopsy material, including surgical material.

4. Information on potential consumers (personnel qualification)

The medical device is intended for use in medical and healthcare facilities in laboratories of the following profile: histological / pathoanatomical / pathomorphological / morphological; for clinical and laboratory diagnostics.

The medical device is intended for use by professional users (doctor, medical laboratory technicians of clinical diagnostic laboratory, pathoanatomical department of healthcare institution).

The FormSAFE container refers to medical devices that do not have direct or indirect contact with the human body.

The medical device, FormSAFE container, is not sterile and cannot be sterilized.

All employees must carefully read the User Manual before using FormSAFE and complete training.

5. Technical specifications of FormSAFE

The plastic container is made from polyethylene terephthalate (PET), and a masterbatch (PET based) is added if necessary. The white cap, which contains a multi-layer liner inside, is made of polypropylene (PP).

No.	Technical specifications	Values
1	Height of container, mm	$61.0 \pm 5\%$
2	Diameter of container, mm	$44.0 \pm 5\%$
3	Height of cap, mm	$11.4 \pm 5\%$

4	Diameter of cap, mm	40.2 ± 5%
5	Internal diameter of container cap, mm	38.0 ± 5%
6	Container cap thread	1 turn thread
7	Torque effect of container cap	min. 120 N*cm/ max. 380 N*cm
8	Container wall thickness, mm	0.75 ± 5%
9	Color of cap	white
10	Nominal volume of container	60.0 ml ± 3.0 ml
11	Volume of fixative in container	36.0 ml ± 3.5 ml
12	Appearance of fixative	Two-layered transparent liquid, colorless, without floating particles and/or white sediment
13	Odor	Characteristic
14	pH	7.0 - 7.4
15	Formalin volume	30.0 ml ± 2.0 ml
16	Isoparaffin volume	6.0 ml ± 1.5 ml
17	Weight of container with fixative, g	50.0 ± 10%

6. Warnings and precautions for use

The product should be handled in accordance with current regulations.

The device is intended for use by specialized technical personnel.

The product is ready for use.

Read carefully the information related to the Hazard Indications and to the Safety Advice Phrases, as per the Safety Data Sheet and the label. Always consult the Safety Data Sheet, where you will find the information pertaining to the risks arising from the product, the precautionary measures to be adopted during use, the first aid and intervention measures in case of accidental release.

1. Not for internal use.
2. Do not drink or inhale solution from container.
3. If swallowed seek medical advice immediately.
4. Avoid that liquid does not make direct contact with skin.
5. When skin become exposed to solution, rinse thoroughly with clear water with soap immediately and consult a doctor if necessary.
6. When you are working with FormSAFE use resistant gloves, goggles and laboratory coat/ apron.
7. After using FormSAFE containers, they should be treated as biohazard material.
8. Precautions must be observed when contacting biomaterial. During work, gloves and special laboratory clothing should be used, since biomaterial samples should be considered as potentially infectious, capable of retaining and transmitting HIV, hepatitis B virus or any other infectious agent for a long time.
9. Do not use the device if its integrity is damaged or the cap and container are deformed.
10. The device cannot be reused.
11. The device must be used before its expiration date.

7. Principle of using

Tissue specimens should be placed in FormSAFE container immediately after removal. Prompt fixation will prevent tissue necrosis, autolysis, putrefaction, and all the degenerative processes. FormSAFE container is large enough to accommodate specimens up to approximately 6cm³ and to allow a ratio up to 5 times the volume of 10% Neutral Buffered Formalin to the volume of the specimen.

If you use specimens already in cassettes in order to ensure a proper fixation insert only one cassette in the FormSAFE container.

Ensure the specimen is completely covered with the reagent.

Details on Principle of Formalin Fixation:

The term “fixation” means the inactivation of the autolytic enzymes and bacteria that cause the putrefaction of tissues and the provision of structural stability to the chemical components of the protoplasm. Fixatives can denature proteins and render the other cellular components insoluble. 10% Buffered Formalin solution is one of the main non-coagulant fixatives; it hardens the protein gels without separating the water from the proteins, and fixes protoplasm without causing the formation of a microscopic spongy mesh. It has a very low oxidation potential, penetrates tissues quite fast and hardens them gradually. With regard to proteins, it binds with NH₂ groups and leaves the hydrophilic groups intact. It can form a methylene bridge between two protein chains, binding to the amine groups of the side-chains. It neutralises basic groups, increasing the acidity of proteins; that is why tissues fixed in formalin display greater affinity with basic rather than acidic dyes. It does not precipitate the DNA, it conserves the majority of lipids and renders phospholipids insoluble to lipid solvents; it does not fix soluble carbohydrates and dissolves glycogen and urea.

8. Information on the conditions required for the collection, processing and preparation of specimens, data on the stability of the analyzed specimens, including storage conditions and duration, restrictions on freezing (thawing) cycles

Tissue specimens should be placed in a FormSAFE container immediately after removal to prevent all the degenerative processes. The volume of the specimen should be no more than 6 cm³.

Preparation of the patient and collection of biopsy material is carried out in accordance with the current local rules for conducting biopsy studies.

Autopsy material is removed in accordance with the current local rules for conducting pathology studies.

When removing the material, it is necessary to prevent its possible mechanical damage, such as compression with tweezers, crushing when cutting with a blunt tool, etc.

The material should not show signs of degenerative processes (putrefaction).

Biomaterial should be stored only in the fixative solution.

Before further use, the container with a specimen should be stored at a temperature between +15 °C and +25 °C, a storage period is no more 12 months. Protect from freezing!

9. Installation, assembly, configuration, preparation for use

Not applicable.

The product is ready for use.

10. Room requirements

The use of FormSAFE container should be carried out in rooms equipped with forced-air and exhaust ventilation and ventilated rooms.

11. Operate with FormSAFE

The product is ready for use.

Do not shake immediately before use.

If the container is shaken before use, wait at least 1 minute before opening the container, in such a way that the protective floating shield reforms and then it prevents the fumes of the formalin to come out.

1. Identify the FormSAFE container with patient ID number using the proper labelling system at your medical facility. It is also available a space on the container label for indicating the patient ID number.

Space for patient ID number.

2. Unscrew the cap and place it by inside down on absorbent paper, in order that the excess of the floating shield (if any) is absorbed by the paper.

3. Put the tissue specimen immediately after removal into the FormSAFE container so that it is completely covered with the fixative and close the cap tightly to the full end of the thread to

avoid spillage during transportation.

4. FormSAFE container can now be transferred to the histopathology lab for analysis.

Keep the container in a vertical position during the transportation.

5. If the FormSAFE container has been shaken during the transportation wait some seconds to allow the “two-phases” stratification and restore the floating shield fluid before opening.

12. List of materials that are required to work with FormSAFE container, but are not included in the medical device delivery kit

FormSAFE container with fixative is an independent medical product, additional materials are not required to work with them.

13. Information needed to identify medical products in order to obtain a safe combination and (or) information on known restrictions on joint use of medical devices

FormSAFE container with fixative is not used in combination with other medical devices.

14. Information on the presence in the medical product of drugs for medical use, materials of animal and (or) human origin

FormSAFE container with fixative does not contain drugs for medical use, materials of animal and (or) human origin.

15. Storage and transportation conditions

FormSAFE container with fixative should be stored in closed ventilated rooms at temperature between +15°C and +25°C, away from heat sources, avoid exposure to sunlight.

Transportation of products is carried out by all types of covered transport in accordance with the rules for the carriage of goods in force for this type of transport, regulating the transportation of dangerous goods. Transport temperature from +2°C to +25°C, avoid exposure to sunlight.

In relation to transportation at temperatures below 15° C, provided that it is allowed, it is possible to see a slight turbidity of the reagent. This does not affect its integrity and functioning in the fixation phase. Before working with the containers received from transportation (when they arrive from the delivery) be sure to let the containers reach the working temperature of the environment (we recommend a minimum of one day).

FormSAFE container with a specimen should be transported at temperature between +2 °C and +25 °C and stored at temperature between +15 °C and +25 °C, a storage period is no more than 12 months. Transportation by all means of transport is allowed.

16. Packaging


FormSAFE container with nominal volume of 60 ml is filled with a fixative solution (30+6 ml).

FormSAFE container consists of a plastic container made of polyethylene terephthalate (PET) and a white plastic cap made of polypropylene (PP) and contained a multi-layer liner inside.

Part of container	Material
Container	polyethylene terephthalate (PET)
Cap	polypropylene (PP)
Cap liner	multi-layer: facing layer - aluminum foil; inner layers - polyethylene terephthalate (PET), polyolefin layer and sealing layer.

The FormSAFE containers are packed in 2 cardboard holders (pallets) of 25 pcs each, placed in a cardboard box (consumer packaging) of 50 pcs (2x25 pcs).

Dimensions of consumer packaging: 324 mm x 322 mm x 133 mm ($\pm 5\%$).


Weight of consumer packaging with FormSAFE containers: 2.90 kg ($\pm 10\%$).

Each consumer packaging contains a User Manual in paper.

The maximum number of consumer packages in a stack is 5.

17. Labelling

Container label

Label on the consumer packaging – a box (50 containers)

The symbols of **UDI**, Data Matrix (2D) code, lot number **LOT** and expiration date  are the variable information.

Explanations of the symbols are provided in section 18.

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

18. Symbols used

Symbol	Symbol name
	Top. Indicate the correct vertical position of the packaging during transportation and storage.
	Refer to the User Manual. An instruction accompanied by this symbol provides a quite important information and require more attention.
	An instruction accompanied by this symbol provides a cautionary statement: WARNING.
	Biohazard.
	This product complies to European community directives.
	IVD medical device.
	Catalog number
	Lot number
	Best before:
	It is prohibited to re-use. Indicates a medical device for one-time use.
	The delivery set is enough for <50> tests. Quantity 50 pcs.
	Storage temperature range in warehouse
	Storage temperature range only for transportation
	Legal manufacturer
	Exclamation sign. Harmful substances.
	Dangerous for human health.
	Unique device identification
	Stacking limit

19. Waste disposal

Handle with care and store in a cool dry space using a tightly closed container. Vapours may collect in empty container. Treat empty containers as hazardous.

Waste material should be disposed of in an approved incinerator or in a designated landfill site, in compliance with all federal, provincial and local government regulations.

For further information, please refer to the SDS (Safety Data Sheet).

For the Russian Federation, it is recommended to dispose of the used FormSAFE container in accordance with the rules implemented by "SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste classified as class B (epidemiologically hazardous waste).

It is recommended to dispose of the unused expired FormSAFE container in accordance with the rules implemented by "SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for maintenance of the territories of urban and rural settlements, for water bodies, potable water and potable water supply, air, soils, living quarters, operation of industrial, public premises, organization and implementation of sanitary and anti-epidemic (preventive) measures" for the medical waste classified as class Г (toxicological hazardous waste).

The FormSAFE container can have potential health effects and cause environmental pollution if not correctly disposed.

20. Responsibility of the Manufacturer

Milestone S.r.l. is the legal manufacturer that is responsible for the quality, efficiency and safety of the registered medical product for use, normal operation and possible negative results of using it, except for cases of negligence, or incorrect use through the fault of user, or force majeure.

21. Repair and service information

The medical product (FormSAFE container with fixative) is single use, and it is not a subject of technical service (maintenance) and repair.

For additional services or support, please contact Customer Support Service or your local distributor.

Customer's Support Service

Tel.: +39 035 572735

E-mail: customersupport@milestonemedsrl.comAuthorized representative in the Russian Federation:**BioLine LLC**Legal address: Russia, 197022 Saint Petersburg, Municipal district Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, structure 3H/5H, room 208Postal address: Russia, 197022, Saint Petersburg, 23 Professora Popova Street, Let. EE-mail: main@bioline.ru

Tel.: +7(812) 320-49-49

Fax: +7(812) 320-49-40

Hotline: 8-800-555-49-40

In case of any undesired reactions or events during the use of the product, please immediately contact with the Customer's Support or your local distributor.

22. Cleaning and disinfection

Not applicable.

The medical product (FormSAFE container with fixative) is single use, and it is not a subject of cleaning and disinfection.

23. Warranty and shelf life

Milestone S.r.l. guarantees that the product has been subjected to a validation procedure based on our in-house testing standards, and that it complies with the technical specifications and/or agreed characteristics warranted in the label. The manufacturer guarantees high quality of product during all term until the expiration date under condition of storage and recommendations observance.

The shelf life of the FormSAFE container is 12 months when stored as recommended. The expiration date is specified on the label.

The FormSAFE container cannot be used after the expiration date.

24. Information on the initial release and revision of the user manual

Initial issue. Rev. 01 dd. 05.2023

25. List of standards applied

The FormSAFE container with fixative complies with the following standards:

- Regulation 2017/746/EC on in vitro diagnostic medical devices.
- Regulation (EC) No 1272/2008 – classification, labelling and packaging of substances and mixtures.
- Regulation (EC) No 1907/2006 - REACH (Registration, Evaluation and Authorisation of Chemicals).

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