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APOSTILLE

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1. COUNTRY: UNITED STATES OF AMERICA
2. THIS PUBLIC DOCUMENT HAS BEEN SIGNED BY:
JANICE M LOPEZ
3. ACTING IN THE CAPACITY OF:
NOTARY PUBLIC OF NEW JERSEY
4. BEARS THE SEAL/STAMP OF:
JANICE M LOPEZ, NOTARY

CERTIFIED

5. AT TRENTON, NEW JERSEY
6. THE 30TH DAY OF DECEMBER 2021
7. BY: Elizabeth Maher Muoio
State Treasurer
8. NO: A780273
9. SEAL/STAMP:
10. SIGNATURE

Elizabeth Maher Muoio
State Treasurer

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UP Regulators Affairs International

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Instruction for use
CereLink® ICP
Monitor

JANICE M. LOPEZ
NOTARY PUBLIC OF NEW JERSEY
Commission # 50048048
My Commission Expires 10/19/2026

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Explosion Hazard: Do not use in the presence of flammable materials (e.g., anesthetics, solvents, cleaning agents and endogenous gases).

Electrical Shock Hazard:

- Use only Codman approved power supplies manufactured for this device. Use of another power supply may not provide electrical isolation from supply mains and protection against electrical hazards.
- Do not remove side, front or rear panels. Contact manufacturer's representative for service and repair.

The ICP Monitor contains a lithium-ion battery inside (Replacement battery, REF 82-6824), which available for separate purchase also.

Warnings for battery:

- Do not open the battery compartment or attempt to replace the battery while monitoring a patient or while the device is connected to the AC Power Supply.
- Use in extreme conditions (e.g., extreme temperature, high humidity, deep discharge) may decrease battery performance.
- Do not puncture the battery.
- Do not dispose of in fire.

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- Do not short circuit battery contacts as battery may get hot, leak, ignite or explode.
- Use only Codman approved batteries listed in Recommended Accessories and Reordering Information section. Use of another battery may present a risk of fire or explosion.
- The device may become warm during normal operation and surface temperatures may reach up to 50°C. Avoid prolonged skin contact with the device (less than 5 minutes) to reduce heat-related concerns. The device complies with user-accessible surface temperature limits defined by the International Standards for Safety (IEC 60601-1).
- Never sterilize the Battery.

For Replacement Power Supply:

Electrical shock hazard:

- The Replacement Power Supply is intended for indoor use only.
- Do not disassemble the Replacement Power Supply in any way. Unauthorized modifications to it can cause it to malfunction or fail in use.
- To prevent injury to the patient and/or user, do not connect the Replacement Power Supply to anything other than the CereLink® ICP Monitor.
- Do not immerse the Replacement Power Supply in any type of liquid.
- Check the Replacement Power Supply for proper function prior to use. Replace if necessary.
- The Replacement Power Supply must be disposed of in accordance with hospital procedures and regulations.
- Ensure the CereLink® Power Supply is grounded:
 - Use only the CereLink Power Supply, which provides three or more electrical pins/contacts depending on the geographical region.

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- Connect the CereLink to only “hospital grade” or “hospital only” receptacles.
- Do not attempt to bypass the grounding contact on the CereLink Power Supply by using an adapter.

NOTE: Failure to properly connect the CereLink® ICP Monitor to a grounded power supply can affect the accuracy and performance of the pressure readings.

For ICP Extension Cable:

- The Extension Cable is MR Unsafe. Do not bring the Extension Cable into the MR environment.
- Modification of the Extension Cable is not permitted. Unauthorized modifications to the Extension Cable can cause a malfunction resulting in serious patient injury.
- Visually inspect the Extension Cable before and after each use for the following:
 - Cable insulation is intact and no cracks are visible.
 - Connectors are undamaged; connector pins are not bent.
 - Connectors must be dry.
- If the Extension Cable does not comply with these conditions or has any other problems that might affect function, do not use the Extension Cable.
- The connectors for the Extension Cable must be properly aligned with the receptacles on the Monitor before pushing these parts together. DO NOT twist the connector. Twisting the connector in this manner will damage the pins and can lead to product malfunction.
- Do not immerse the Extension Cable in any liquid or allow liquid to enter the connectors. Do not use the Extension Cable if connectors become wet; erratic results may occur.

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- Calculation (a.k.a., Pressure Time Dose or PTD) where $PTD = \text{Time} \times \text{ICP measurement}$.
- Historic mean ICP histogram with time above threshold calculation
- Continuous operation on AC power
- Internal rechargeable battery for at least 2 hours of continuous operation
- User programmable visual and audible mean ICP alarms that activate if the ICP valve exceeds a user-specified limit
- Real-time data streaming of mean ICP and/or ICP waveform through the USB connection
- Data storage and export via USB:
 - 24 hours of ICP waveform
 - 14 days of mean ICP
 - 10 screen-shots
 - Data markers (points in the data identified by the user as special events)
- One-touch transducer zeroing function
- Free-standing positioning or attachable to an IV pole
- 7" (177,8 mm) color touch screen

Description of monitor:

The CereLink® ICP Monitor (REF 82-6820) packaged with a ICP Extension Cable (REF 82-6845), rechargeable battery (REF 82-6824, installed in the monitor and also available for additional purchase) and replacement power supply (REF 82-6822). Also the package could include Adapter (REF EXPORTCAB), if necessary. The devices are provided to the user non-sterile. All the devices are available for separate order also.

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
Figure 1: Front and Side of ICP Monitor

A Power Button

- Press for 1 second to turn device on.
- Blue light – Device powered on.
- Press for 3 second to turn device off.
- No light – Device powered off.

B External Power Source / Battery Charging Indicator

- Green light – Device connected to external power supply.
- No light – Device not connected to external power source.

C Touch Screen

D System Failure Indicator

- Red light – System issue. Please refer to Troubleshooting section for details.

E Rubber Feet

F Mute Alarm Button:

- Push to mute alarm for 3 minutes
- No light – No alarm
- Yellow light flashing – Alarming
- Yellow light steady on – Alarming and Muted

G ICP Sensor Extension Cable Connection (yellow)

H Patient Bedside Monitor Interface Cable Connection (blue)

I High Speed USB Ports

J Alarm Beacon

- No light – No alarm
- Yellow light flashing – Alarm or Alarm Muted
- Yellow light steady on – System Failure

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
Figure 2: Back of ICP Monitor

- | | |
|---|--|
| <p>(A) System Information Label
See Labeling Symbols Glossary section.</p> <p>(B) Battery Access
See Battery Replacement section.</p> <p>(C) IV Pole Clamp</p> <p>(D) Power Input</p> | <p>(E) Power Cable Retention Mechanism
Refer to Figure 3.</p> <p>(F) Recessed Hand Hold</p> <p>(G) Cooling Vents</p> |
|---|--|


Figure 3: Power Cable Retention Mechanism Usage

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MECHANICAL	
Box Dimensions ($\pm 5\%$)	166 mm H x 222 mm W x 60 mm D (118 mm D with clamping)
Box Weight ($\pm 5\%$)	1.8 kg
Table Top	YES
Connectors	Connectors Keyed and color coded ICP Input: yellow Bedside Monitor: blue
Push Button Depress Force	between 100g and 700g

POWER	
Power Supply	Medical grade Input: 100 - 240V AC, 50 - 60Hz Output: 12V DC, 2.5A, 30W • Earth connection is only functional, providing 2 means of protection from line and neutral to the functional earth
Battery Type	3.7 V (5Ah) Li-Ion Polymer, rechargeable
Battery Life	2 hours (fully charged, new battery)
Battery Charge time	3 hours (for fully discharged battery)
Low battery alert	Approximately 15 minutes of battery life remaining

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INPUT PRESSURE CHANNEL

Operating Range	-50 to 150mmHg
Sensor Excitation Voltage	5VDC +/- 0.01V
Pressure Signal Sensitivity	5uV/Vex/mmHg
Input Impedance	10 MΩ nominal
Pressure Signal Band Width	0 to 25Hz (-3dB)
Zero Range	+/- 50 mmHg
Input Channel Accuracy	+/- 0.5mmHg over the range -50 to 50mmHg and +/-1% over the range 50 to 150mmHg (exclusive of sensor)
Input Channel Resolution	0.1mmHg

DISPLAY

Type	TFT 1024 (H) x (RGB) x 600 (V) (±5%) LCD with projected capacitive touchscreen
Size (±10.00 mm)	7.0" (= 177,0 mm)
Viewing Angle	+/- 30°

VISUAL INDICATOR

Power indicator Color	Blue (400nm to 700nm)
Power indicator Intensity	light in the intensity of 0.2 to 0.4 lumens at 0.1m +/- 0.01m.

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AC/Battery Status Indicator Color	Green (490nm to 540nm)
AC/Battery Status Indicator Intensity	intensities of light at 0.1m +/- 0.01m: Green (0.2 to 0.9 lux)
Alarm Indicator Color	Yellow (565nm to 595nm)
Alarm Indicator Intensity	intensities of light at 0.1m +/- 0.01m: Yellow (0.2 to 0.4 lux)
Alarm Mute Visual Indicator Color	Yellow (565nm to 595nm)
Alarm Mute Visual Indicator Intensity	light in the intensity of 0.2 to 0.4 lumens at 0.1m +/- 0.01m.
Fault Visual Indicator Color	Red (620nm to 670nm)
Fault Visual Indicator Intensity	intensity of 0.4 to 1.0 lumens at 0.1m +/- 0.01m.

ALARMS

Alarm Parameter	Mean ICP
High Alarm Range	-50 to 150 mmHg
Low Alarm Range	-50 to 150 mmHg
High/Low Alarm Step Resolution	1 mmHg
High Alarm Default Value	20 mmHg
Low Alarm Default Value	0 mmHg
Audible Alert Indicator Volume Range	in the range of 30 to 60 dBA, measured at a distance of 1m +/- 0.01m, from any orientation.

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USER INTERFACE

Display Mean ICP	YES
Mean ICP Resolution	1 mmHg
Display ICP Waveform	YES
ICP Waveform Amplitude	-50 to 50 mmHg; -10 to 20 mmHg; 0 to 30mmHg; 0 to 50 mmHg; 0 to 100 mmHg; 0 to 150 mmHg.
Display mean ICP Trend	YES
ICP Trend amplitude	-50 to 50mmHg; -10 to 20mmHg; 0 to 30mmHg; 0 to 50mmHg; 0 to 100mmHg; 0 to 150mmHg.
ICP Trend window	3hr, 12hr, 24hr, 48hr, 6 days
ICP Trend time scroll	14 days
Display High/Low Alarms	YES
Histogram screen	Two percentage thresholds. Calculate and display a percentage histogram of the amount of time mean ICP values have been within the predetermine ranges with respect to the time scale selected. The available ranges are: <0, 0 to 5, 5 to 10, 10 to 15, 15 to 20, 20 to 25, 25 to 30, 30 to 35, 35 to 40, 40 to 45, 45 to 50 and >50 mmHg.

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Histogram window	3hr, 12hr, 24hr, 48hr, 6days
Screen Shots	YES
Time and Date	YES
Multilanguage	YES English, Simplified Chinese, French, German, Italian, Portuguese, Russian, Spanish, Japanese, Swedish, Finnish, Danish, & Norwegian.

STORAGE

Mean ICP	14 days, 1 sample every 10 seconds
Mean ICP Alarms	14 days. 1 sample every 10 seconds
Waveform	24 hours full resolution
Screen	50 Screen Shots (10 for each main screen type)

SAFETY

Input type	Type BF equipment, floating input protected against defibrillator discharge
Ingress Protection Rating	IP22
Fire Retardant Rating	UL 94 V-0

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Port	USB 2.0
Streaming	USB to RS232 Serial port setup: Speed – 115200 Data – 8 bit Parity - none Stop Bits - 1 bit Flow Control - none Mean ICP and Waveform
Data Download	USB, .csv file, 24 hours waveform, 14 days mean ICP trend, 50 Screen Shots, Alarms, and Data Markers.

Replacement Power supply (REF 82-6822) used to power the medical device. The power supply cable is protected by a braided shield that connects to earth ground via a filtering choke. The package contain additional one cord with US blade and 3 blades depends of country:

Diameter	5,4 mm $\pm$ 3%
Cable Length	4000 mm $\pm$ 150 mm
Conductors	One (1) conductor surrounded by braided shield connected to neutral
Specifications:	Class II double insulated 120 VAC floating power supply (i.e., 2-pin connection plug)
Type	Wall Plug-In type
Input Voltage Range:	100-240V

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CereLink® ICP Extension Cable (REF 82-6845) used to connect monitor to ICP sensors. This cable is reusable and is supplied nonsterile.

Weight ($\pm$ 5%)	120 g
Length ($\pm$ 25mm)	3000 mm
maximum force to mate the ICP sensor Connector	$\leq$ 30.0 N
Demating Force (The force required to disconnect the cable from the ICP sensor)	From $\geq$ 4.5N to $\leq$ 35N

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Replacement battery (REF 82-6822) is a lithium-ion battery designed for use with the CereLink® ICP Monitor (REF 82-6820) and used to power the medical device.


Type	Rechargeable Li-Ion Polymer Battery
Typical Capacity	5000mAh
Used cell type:	Li-Polymer 726090
Battery configuration:	1S1P with PCM, NTC, Connector, Lead wires
Heavy metal content:	Mercury free, Cadmium free, Lead free
Internal resistance: (1kHz, fully charged)	approx. 140mΩ
Normal Volt	3.7V
Limit Charge Volt	4.2V ±0.03V (4.1V recommended due to safety and cycle life => approx. 10% capacity loss)
Standard Discharge Current:	1.000mA (0.2C-rate) @ 3.0V
Max. Discharge Current:	5.000mA (1.0 C-rate) @ 3.0V
Discharge Environment	From -10°C to 60°C
Standard Charge Current:	1.000mA (0.2C-rate) @ 4.2V
Max. Charge Current	2.500mA (0.5 C-rate) @ 4.2V
Charge Environment	0°C to +40°C
Charge Time (CCCV Mode):	3 h (± 10%)

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USB-to-RS232 Adapter cable (REF EXPORTCAB) used to digitally export ICP mean trend data from an ICP monitor to a remote device.

1 – connector and cable


2 – Connection of Adapter to Monitor through USB port

During streaming, the monitor reports the average ICP value in mmHg, determined on the basis of current ICP dynamics data, with an interval of one minute. The USB to RS232 adapter cable transfers information at 9600 baud, 8 data bits, 1 stop bit, even pair.

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


The overview of serial connector provide below:

RS-232 Pinout – DB9 Male DTE		
PIN	Signal Name	RS-232 Signals
1	DCD (Data Carrier Detect)	Input
2	RD (Receive Data)	Input
3	TD (Transmit Data)	Output
4	DTR (DTE Ready)	Output
5	SG (Signal Ground)	Ground

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
Connectivity type	USB to RS232
USB port type	USB 2.0
Length of cable ($\pm 5\%$)	950 mm
Connector dimensions in mm ($\pm 0,5$ mm)	L 52 x W 32 x H 16,4 mm
Weight, gm with cable ($\pm 5\%$)	104.3 gm
USB type	B (F)
Adapter type	DB-9 (M)
Speed of streaming of information	115200
Data	8 bit
Parity	none
Stop bits	1 bit
Flow control	none

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When connected to a ICP Sensors, the ICP Monitor provides a numeric display of the mean ICP, ICP waveform and mean ICP trend.

Operator could use only special ICP Sensors, Russian Registration Certificate № P3H 2021/14639.

The technical specification of sensors provided below:

SENSING TECHNOLOGY – ICP Sensors	
Pressure Sensor Technology	Strain gauge
Configurations	Basic Kit, Metal Bolt Kit, Plastic Bolt Kit, Ventricular Catheter Kit, Ventricular Catheter Kit with Tuohy-Borst adapter
Zeroing	Automatic, stored in sensor memory
Zeroing medium	Sterile water or sterile saline
Reconnection to other CereLink® ICP Monitors	YES
MRI	1.5T and 3T
Pressure Accuracy	a) ± 2 mmHg over the range -20 mmHg to 20 mmHg b) $\pm 10\%$ of reading over the range -50 to -20 and 20 to 150mmHg

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CereLink® ICP monitor could be connected to external patient monitor with compatible Interface cable. ICP Monitor could be connect with external patient monitors manufacturer be following manufacturers and which have the following technical specifications.

Manufacture by:

- DRAEGER\SIEMENS (connection through Integra Interface cable REF 82-6880)
- PHILIPS (connection through Integra Interface cable REF 82-6881)
- GE (connection through Integra Interface cable REF 82-6882)
- Spacelabs (connection through Integra Interface cable REF 82-6883)
- Nihon-Kohden (connection through Integra Interface cable REF 82-6887)
- Fukuda (connection through Integra Interface cable REF 82- 82-6889 and 82-6888)

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- GE Datex-Ohmeda (connection through Integra Interface cable REF 82-6884)

EXTERNAL PATIENT MONITOR INTERFACE	
Output Signal Sensitivity	5uV/V _{ex} /mmHg
Input Excitation Voltage Range	2.5 to 6.5 VDC (incl. symmetric supply voltages)
Output Accuracy	+/- 0.5 mmHg over the range -50 to 50 mmHg and +/-1% over the range 50 to 150 mmHg (exclusive of sensor)
Output Calibration Signals	0, 20 and 100 mmHg
Input Impedance	350 to 500 Ω
Source Impedance	350 Ω
Output Pressure Signal Band Width	0 to 25Hz (-3dB)
Zero Range	± 50 mmHg
Expected External Patient Monitor Isolation	Type BF or Type CF. The external monitor must provide this type of isolation to ensure patient safety

Refer to the external Patient Monitor Instructions for Use before connecting CereLink® ICP Monitor. Make sure the external monitor pressure input is compatible with CereLink® ICP Monitor output.

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The ICP Monitor can stream data through the USB port using a USB to serial communication. Computers used for streaming data must be connected to a grounded power source, as describe below:

- Computer connected to CereLink® requires an earth ground connection which provides three or more electrical pins/contacts depending on the geographical region.
- Connect only to "hospital grade" or "hospital only" receptacles.
- Do not attempt to bypass the grounding contact on the computer power supply by using an adapter.

Earth ground symbol. Connection plug must have three or more electrical pins/contacts.

The serial port settings required to connect a PC to the ICP Monitor through a terminal application or other customized serial program are as follows.

- Speed: 115200
- Data: 8 bit
- Parity: none
- Stop: 1 bit
- Flow Control: none

The communication protocol is unidirectional. The ICP Monitor will only output data and no responses shall be sent from the receiving device (i.e., PC Host should be in a "listen only" mode). The ICP Monitor automatically detects in which USB port the USB/Serial connection is present. Either USB port may be used based on user preference. The data stream is always active except when the Data Download menu is open. Steaming will resume as soon as the Data Download menu is closed. Data packets are transmitted at various frequencies depending on the data type.

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System Status Packet Description:

Description	String	Size (bytes)	Details
Start of Packet	*	1	-
Product ID	CERELINK	8	-
Device SN	CLKxxxxxxxx	10	-
Protocol Version	V01	3	-
Packet Type	ICPNC ICPNZ or FAULT	5	ICPNC = ICP sensor Not Connected ICPNZ = ICP sensor Not Zeroed or FAULT = System is in a Fault State
Date	YYYY/MM/DD	10	-
Time	HH:mm:ss	8	-
FAULT Error Code	XXXX	4	Only available for FAULT Packet Type

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Message Counter	00001	5	Increment with each message Reset when 99999 is reached
Packet Checksum	0xFFFF (initial value)	4	CRC16 $x^{16} + x^{12} + x^5 + 1$
End of Packet	\r\n	2	-

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Package includes a ICP Monitor with a pre-installed battery inside (1 pcs) and a replaceable power supply unit (from 1 to 5 pcs, depending on the customer's order). If necessary, the delivery set is also completed with additional extension cables, additional batteries, an additional USB-to-RS232 Adapter cable. The quantity of each of the products, if necessary, varies from 1 to 5 pieces, depending on the order of the customer.
Each product comes with an IFU.

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1) Unpack the ICP Monitor.

WARNING: Use only manufacturer approved accessories with the ICP Monitor. Mounting and installation must be carried out by trained personnel.

Remove the contents from the shipping box and verify the items. After unpacking the contents, inspect the items for any signs of damage. If any damages are discovered, notify the carrier and the supplier. Retain all shipping contents and boxes for examination.

2) Prepare the Power supply.

Remove the Power Cord from the package and attach the region-specific adapter plug or cord to the power supply.

3) Plug the ICP Monitor into AC Power.

Connect the power supply to the ICP Monitor by inserting the connector end into the Power Input. Insert the grounded electrical plug (with three or more electrical pins/contacts depending on the geographical region) of the Power Cord into a grounded wall outlet

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and ensure the Power Indicator light is green. Use the Power Cable retention mechanism to prevent the Power Cord from detaching from the ICP Monitor.

Note: Use only provided CereLink® ICP Power Supply.

- Use only the CereLink® power connection, which provides three or more electrical pins/contacts depending on the geographical region.
- Connect only to “hospital grade” or “hospital only” receptacles.
- Do not attempt to bypass the grounding contact on the CereLink® Power Supply by using an adapter.
- Connection plug must have three or more electrical pins/contacts.

4) Charge the Battery to Full Capacity.

Keep the ICP Monitor plugged in for at least 5 hours.

5) Turn on the ICP Monitor.

Press the Power Button on the front of the ICP Monitor for at least 1 second. The Power Button will illuminate to indicate that the device is on. After approximately 20 seconds, while performing internal self-checks, the ICP Monitor will display the “Codman CereLink” logo. Confirm that the battery is fully charged.

6. Specify the Language.

From the Home Screen, select Settings, then select Set Language. Choose the desired language.

7. Set the Time and Date.

From the Home Screen, select Settings, then select Set Time/Date. Adjust the time and date as needed.

8. Set the Alarm Volume and Screen Brightness.

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3. Connecting a New ICP Sensor.


A: Connect the ICP Extension Cable to the yellow connector on the ICP Monitor. The ICP Extension Cable should be autoclaved according to the associated instructions for use. Follow your institution's procedures and maintain sterile technique.

NOTE: Cable Connectors are push/pull. DO NOT twist.

- Extension Cable connector tab and arrows should face forward, in same direction as the monitor screen.
- Push yellow connector into yellow port until you hear a click.
- To remove cable, grab gray textured part of connector and pull away from monitor.

B: While maintaining sterile technique, connect the ICP Sensor to the ICP Extension Cable. The connection is a tight fit with no audible or tactile click.

C: Once the sensor is connected, the monitor will play a sound and will display the zero symbol.

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
ICP Sensors must be “zeroed” before implantation. The zeroing process is performed to calibrate the Sensor to atmospheric pressure prior to measuring ICP. Follow the zeroing steps included in the Instructions for Use packaged with the ICP Sensor and Extension Cable.

CAUTION: Do NOT submerge the tip of the ICP Sensor vertically in a deep pool or cup of sterile water or sterile saline. Doing so will impose a hydrostatic pressure on the ICP Sensor diaphragm that is higher than atmospheric pressure, resulting in an inaccurate zero reference.

A: Prepare to zero the device by laying the tip of the ICP Sensor flat in a shallow pool of sterile water or sterile saline.

B: Access the sensor zeroing menu by selecting the Home button from any of the main screens and then the Zero Sensor option.

C: While keeping the tip of the ICP Sensor flat and still in the sterile water or sterile saline, press the zero button on the ICP Monitor screen.

NOTE: The device is equipped with an algorithm that will detect when the sensor is ready to be zeroed.

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

Zero Sensor

1. Keep sensor tip wet during entire zeroing process (see sensor instructions)

2. Hold sensor still. Press to zero.


Ready to Zero

Sensor not Zeroed

- Zero button will NOT be available if sensor is moving.
- Zero button will NOT be available if sensor is already implanted.

Zero Sensor

1. Keep sensor tip wet during entire zeroing process (see sensor instructions)

2. Hold sensor still. Press to zero.


Not Ready
Sensor implanted or not still

Sensor not Zeroed

D: When sensor zeroing has been successful, the date and time will be displayed on the screen:

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E: To confirm the Mean ICP value is zero, return to the ICP Waveform Screen. If necessary, repeat the zeroing steps.

F: While maintaining sterile technique, disconnect the ICP Sensor from the ICP Extension Cable before implantation and after zeroing. Follow the implantation steps listed in the Instructions for Use packaged with the ICP Sensor.

5. Reconnecting an ICP Sensor after Implantation.

After the ICP Sensor has been zeroed and surgically implanted, reconnect sensor to the ICP Extension Cable on the side of the ICP Monitor.

NOTE: The sensor can be reconnected to any CereLink® ICP Monitor. It is not required to use the same monitor used for the zeroing process. Historic trend data will not be transferred if the sensor is connected to a new CereLink® device.

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Existing Sensor Data.

Data from two different sensors cannot be stored on the same monitor. Prior sensor data **MUST** be deleted before using a new sensor. When a new ICP Sensor is detected, the device will:

1. Confirm deletion of data.
- OR
2. Recommend disconnecting sensor to prevent data loss.

NOTE: If the sensor is disconnected, monitoring will **NOT** resume until the previous sensor is reconnected, or the stored data is downloaded through the USB

Taking a "Screen Shot".

To save an image displayed on the ICP Monitor screen, press the Screen Shot button. Up to 50 images (10 for each main screen) can be stored at the same time. After 10 images, of the same screen type, the oldest image in memory will be overwritten. Images can be downloaded to a USB memory drive via the Data Download option on the Home menu.

Using Data Markers.

Use the Data Marker button on right side of the screen to identify and save in memory special events in real time. For example, use the Data Marker before and after patient transport, administration of medication, change in dosage, or any other

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interventions or procedures. Every time the Data Marker button is selected, a unique number (1-99) will be displayed on the screen to identify the event marker. Log this number on the patient history file along with the event information. The Mean ICP Trend Screen will display all event markers on a pop-up screen at the time the data was tagged. The pop-up screen will read "Data Marker #X, Set in trend screen". After 99 stored events, the data marker counter will reset to one. The Data Marker ID is linked to the mean ICP trend data and can be downloaded via USB.

Data Download.

The device can store the following data:

- The last 24 hours of the ICP waveform (100 samples per second)
- Last 14 days of mean ICP values (1 sample every 10 seconds)
- Last 14 days of physiological alarm conditions
- Last 14 days of event markers
- 50 screen shots

Data will be overwritten once the memory capacity has been reached, starting from the oldest data point. A warning will be displayed prior to overwriting data except for ICP waveform. The ICP waveform will be replaced every 24 hours. Make sure to download waveform information as frequently as necessary. Stored data will not be lost if the sensor is disconnected or the monitor is powered down.

To Download Data:

1. Insert a digital memory device into the USB port on the right side of the monitor.

USB device should meet the following specifications:

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When one device is female and one device is male, use a straight-through cable.

3. Loopback Test.

Loopback pins 2 and 3. Using Hyper Terminal or similar program, connect to the COM port. Set the desired baud rate. Ensure that Hyper Terminal local echo is OFF. Transmit data. If the same character string is returned, the test is good.

Transporting the Patient – Battery Power.

CAUTION: The ICP Monitor is protected against ingress of dripping fluids when tilted up to 15 degrees from its upright position. However, the fluid ingress protection may be compromised when not positioned upright, for example if transporting while on the patient bed. Always ensure that the cooling vents are unblocked and kept away from fluid sources.

The ICP Monitor includes a rechargeable lithium-ion battery, that supplies power to the ICP Monitor for at least 3 hours when the battery is fully charged. It is recommended to leave the ICP Monitor connected to a power supply during routine use. The battery should only be used for short periods of time when a power supply is not available, e.g., during patient transport.

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To attach the Sensor to the Extension Cable and then connect the Extension Cable to the Monitor, follow these steps:

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WARNING: The ICP Monitor provides a blue output connector that allows connectivity to an external Patient Monitor.

Use the appropriate Interface Cable manufactured by Integra to connect to the patient monitor.

NOTE: The patient monitor calibration can be performed without a sensor connected. In this case, after calibration, the CereLink® ICP Monitor will transmit an ICP value of -50mmHg to the Patient Monitor.

Connect Interface Cable.

NOTE: Connectors are push/pull. DO NOT twist.

- Connector tab and arrows should face forward, in same direction as the screen.
- Push blue connector into blue port until you hear a click.
- To remove cable, grab textured part of connector and pull away from monitor.

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
Select Patient Monitor Setup.

Once the External Patient Monitor is recognized by the ICP Monitor, the "Patient Monitor Setup" option will be available in the Home Screen.

NOTE: You can access the Patient Monitor Setup menu by selecting the Home button from any of the main screens.

Start Calibration Process.

- a. Start the calibration process by selecting the "0" button on the CereLink® ICP Monitor, as provided below. The monitor will send a 0mmHg signal to the External Monitor.

- b. Zero the pressure input on the External Patient Monitor.
- c. Wait until the External Patient Monitor displays "0mmHg".

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




Return to the Codman CereLink® ICP Monitor and press the arrow to continue the synchronization process, as provided below:

Use the available calibration signals (0, 20 and 100mmHg) on the ICP Monitor to confirm the correct calibration and scaling of the External Bedside Monitor:

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
Refer to the Instructions for Use of the External Bedside Monitor to optimize the calibration.

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ICP Waveform Screen (default screen):

(A) Mean ICP Status and Numeric Displays:


25 ICP Sensor connected, zeroed, and acquiring data in millimeters of mercury (mmHg). ICP value will flash if alarming.

→0← ICP value will flash if alarming. ICP Sensor connected, but not zeroed.

-- ICP Sensor not connected.

(B) Alarm Threshold Indicators and Alarm Status:

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HIGH/LOW alarm thresholds and current mean ICP indicator (triangle). Mean ICP Indicator will flash if alarming.

Alarm ON. Alarm will sound if the mean ICP value exceeds the alarm threshold.

Alarm OFF. Alarm will NOT sound if the mean ICP value exceeds the alarm threshold.

C

Dashed Lines in Graph:

Alarm threshold value.

D

ICP Waveform

E

ICP Waveform Scale Adjustment:

Press the button to cycle through the available ICP amplitude scales (0 to 30, 0 to 50, 0 to 100, 0 to 150, -50 to 50, -10 to 20).

F

Current Time and Date

G

Time and Date the Sensor was Previously Zeroed

H

Sensor ID: Available when a sensor is connected or when there is data stored in the device. Provides a unique number to identify the data.

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I Patient Bedside Monitor Indicator:
 Bedside Monitor NOT detected.
 Bedside Monitor Detected. Symbol flashes when monitor is not calibrated.

J External Power Status Indicator:
 AC power NOT detected.
 AC power present and charging the battery.

K Battery Status (Level) Indicator:
 Battery level full.
 Battery level low, 40% remaining.
 Battery level low, 20% remaining. A warning will be displayed when 15 minutes of battery life remains.
 Battery level critical. Charge immediately. If power is disconnected the device may shut down.

L Screen Toggle Button:
 Access to trend, boxplot and histogram screens.

M Screen Shot Button:
 Screen Images will be stored in memory and can be downloaded to a USB stick.

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Mark special events on the trend screen. Refer to Using Data Markers section for more information.


Home Button:

Refer to Home Screen section for more information.

Mean ICP Trend Screen:

The ICP Trend Screen displays the historic mean ICP values over a period of time. Additionally, the Pressure Time Dose (PTD) is also presented on this screen. PTD is a measure to quantify the burden and time spent above an ICP threshold. PTD represents the area under the curve above a user-defined threshold within the displayed time frame (refer to Analysis Threshold Setup Screen section for information on how to set the two available thresholds). For example, in the graphic below, the two thresholds selected are “22 mmHg” and “15 mmHg”.

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
A

Trend/Boxplot button:

Toggle between trend and boxplot view.

B

PTD Thresholds:

Dashed lines represent PTD thresholds.

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E

Data Markers:

Icon marks the point in time special events occurred. Refer to Using Data Markers section for more information.

F

Trend Time Label:

Time axis corresponding to the time the trend data was acquired.

G

Timescale unit of measure:

Amount of time between two grid line divisions.

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 Waveform Trend Arrow: Right:
 Press arrow to navigate trend data.

 Timescale button:
 Press the button to cycle through the available time scales (3 hr, 12 hr, 24 hr, 48 hr, 6 days).

 Trend Date Label:
 Time axis corresponding to the date the trend data was acquired.

 ICP Trend Stats:
Section where the PTD calculations are displayed:
1. PTD = Pressure Time Dose (mmHg * hours)
2. (xx) = Pressure Threshold at which PTD is calculated (mmHg).

Mean ICP Boxplot Screen:

The ICP Boxplot Screen shows the same data from the mean ICP trend screen, but as a graphical representation of the statistical data. The box plot graphic uses descriptive statistics to display the distribution of the mean ICP data based on quartiles. The quartiles split the data into four equal groups based on five key values: minimum (0%), first quartile (25%), median (50%), third quartile (75%), and maximum (100%). The minimum and maximum represent the true extreme values in the data set. The median represents the middle number of the data set, the first quartile is the middle number between the minimum and the median and the third quartile is the middle number between the maximum and the

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

median. Therefore, the box includes 50% of the data points around the median point. Additionally, the mean value within each time division is represented as green diamonds.

The functionality of the Boxplot resembles that from the Mean ICP Trend screen.

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The ICP Histogram screen shows the distribution of mean ICP values over a period of time. Additionally, the Time Above Threshold is calculated. Time above threshold is the total time the mean ICP has been above one of two user-defined thresholds within the selected timeframe (refer to Analysis Threshold Setup Screen section for information on how to set the thresholds).

NOTE: The same thresholds used for the PTD calculations in the Mean ICP Trend Screen are used for the Time Above Threshold calculation.


- A** ICP Time Above Threshold:
1. T = Time above threshold (hours, minutes and %).
 2. (xx) = Pressure Threshold at which T is calculated (mmHg).

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Press the button to cycle through the available time scales (3 hr, 12 hr, 24 hr, 48 hr, 6 days).

Ⓒ

T Thresholds:

Dashed Lines indicate T Thresholds

Ⓓ

Histogram:

Distribution of mean ICP within the selected time window.

Home screen:

The Home Screen provides access to several frequently used monitor functions and is accessed by pressing the Home Button on the User Interface Screen.

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
Zero Sensor:

Pre-implantation Zeroing of the ICP Sensor. Only available if a sensor is detected.

Alarm Setup:

Set mean ICP alarm limits.

Patient Monitor Setup:

Synchronize ICP Monitor with external patient monitor. Only available if an external Patient Monitor is detected.

Analysis Setup:

Set Pressure Time Dose (PTD) and Time Above Threshold limits.

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
Set the current time and date.

Language:

Set the user interface display language.

Technician Access:

Information for Codman Service and Repair personnel.

Home:

Exit the Settings menu and return to Home menu.

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
This screen is accessed by selecting HOME -> ALARM SETUP.

Use this screen to select the HIGH and LOW alarm for the mean ICP value.


The HIGH and LOW alarm can be turned ON/OFF with the switch on the left (green indicates alarm is ON).

Alarm values must be within -50 to 150mmHg and the HIGH alarm value must always be greater than the LOW value.

Analysis Threshold Setup Screen:

This screen is accessed by selecting HOME -> ANALYSIS SETUP.

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

Use this screen to select the HIGH and LOW levels used for the Pressure Time Dose (PTD) and Time Above Threshold calculations.

The HIGH and LOW levels can be turned ON/OFF with the switch on the left (green indicates that the alarm is ON).

The threshold values must be within -50 to 150 mmHg and the High value must always be greater than the Low value.

Data Download Screen:

This screen is accessed by selecting HOME -> DATA DOWNLOAD.

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

Use this screen to select which data to download to an external USB memory drive.

Note that more than one item can be selected for download. Refer to Data Download section for additional information.

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The ICP Monitor can download data to an external USB flash drive with the following specifications:

- USB 2.0
- 8 to 32 GB capacity
- FAT32 format

Data available for download include:

- Last 1 hour or 24 hours of ICP Waveform data (100 samples per second)
- Trend data, including last 24 hours or all 14 days of:
- Mean ICP values (1 sample every 10 seconds)
- Physiological alarm conditions
- Event markers
- Last or all available screen shots (50 maximum)

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Data is stored in the following format:

- Waveform data - Tab delimited text file with 3 columns:

Date	Time	ICP [mmHg]
YYYY/MM/DD	HH:mm:ss.xx	XX.X

- Trend data - Tab delimited text file with 5 columns:

Date	Time	ICP [mmHg]	Alarm[H/L]	Marker[1-99]
YYYY/MM/DD	HH:mm:ss	XX.X	X	X

- Images – Bitmap (BMP) file

The data will be stored in the USB flash drive in files following a specific file naming convention:

- Waveform data:

CereLink-wave-SNxxxxxxx-YYYYMMDD-HH:mm:ss(#)

- Trend data:

CereLink-trend-SNxxxxxxx-YYYYMMDD-HH:mm:ss

- Screen shots:

CereLink-IMG-zzzz-SNxxxxxxx-YYYYMMDD-HH:mm:ss

- Where:

YYYY = Year

MM = Month

DD = Day

HH = Hour (24 hour format)

mm = minutes

ss = seconds

(#) = file number

xxxxxxx = Sensor serial number (uses "x" if no sensor is connected)

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File Name Examples:

- Monitoring started 12/01/2017 at 3:45 PM
- Sensor Serial number = 1234567
- Download initiated to USB flash drive at 12/6/2017 3:54 PM
- Patient has been monitored for 5 days and 9 minutes.

Scenario 1: Last 24 hours of trend data:

- CereLink-trend-SN1234567-20171206-155400

Scenario 2: All available trend data:

- CereLink-trend-SN1234567-20171206-155400

Note: Scenarios 1 and 2 will create the same file name, but with different amount of data

Scenario 3: Last 1 hour waveform:

- CereLink-wave-SN1234567-20171206-155400(1)

Scenario 4: Last 24 hours waveform:

- CereLink-wave-SN1234567-20171206-155400(1)
- CereLink-wave-SN1234567-20171206-155400(2)
- CereLink-wave-SN1234567-20171206-155400(3)
- : : : : : :
- CereLink-wave-SN1234567-20171206-155400(23)
- CereLink-wave-SN1234567-20171206-155400(24)

Note: Scenarios 3 and 4 will create the same file name for the first hour data, containing the same information

Scenario 5: 2 images in memory, one taken before sensor was connected:

- CereLink-IMG-wave-SN1234567-20171205-114600
- CereLink-IMG-wave-SNxxxxxxx-20171201-150100

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Alarms.

The ICP Monitor activates three types of audio/visual alarms to indicate problems or situations that require immediate attention:

1. Physiological Alarms.

CAUTION: The Physiological alarm can be disabled by manually turning it off or by selecting extreme limits.

Alerts the user that the patient's mean ICP value has exceeded the specified alarm threshold set by the user. Physiological alarms can be muted temporarily for 3 minutes by pressing the mute button on the front of the device. These alarms will be stored along with the historic trend data and can be retrieved via USB.

NOTE: Alarm functionality can be tested any time a sensor is connected by setting the alarm limit below the ICP value.

2. System Alerts.

Alerts the user of situations that require user action, including:

- a. New ICP Sensor connected. User must acknowledge if old data can be erased. Otherwise, sensor must be disconnected.
- b. ICP Sensor disconnected from the Extension Cable.
- c. ICP Sensor failure.
- d. Screen shot limit reached. Oldest image will be overwritten.
- e. Mean ICP trend data limit reached. Oldest data will be overwritten.
- f. Battery charge at 10%.

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Alerts the user that there is a system malfunction and in general, the device becomes unusable. Contact authorized representative for further instructions.

Alarm Priority.

Alarm Type	Alarm Priority	Action	Audio Indicator	Visual Indicator
System Failure Alarms	Medium Will override all other alarms.	Patient ICP will not continue to be monitored. Troubleshoot system as it may be unusable.	Burst of 5 pulses with a maximum volume of 60dBA (during active ICP monitoring).	Fault indicator will illuminate and yellow screen will be displayed identifying the error code. Alarm beacon will illuminate yellow (during active ICP monitoring).
Physiological Alarms	Medium Will override System Alerts.	Patient ICP has exceeded the alarm limit. Prompt attention to patient is necessary to prevent a potential harm to occur.	Burst of 5 pulses with a maximum volume of 60dBA.	Alarm Beacon will flash yellow ICP value will flash. Active alarm: Mute button will flash yellow as sound is played. Audio Paused: Mute

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				button will illuminate yellow (not flashing).
System Alerts	Low	ICP monitoring may stop for technical reasons (e.g. sensor disconnected, low battery). Prompt attention to device is necessary to prevent a potential harm to occur.	Burst of 3 pulses with a maximum volume of 60dBA.	Popup message will appear on the screen indicating the alert type and waiting for user action.

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System Problems.

Problem	Cause	Recommendation
Battery not charging or not holding charge.	Defective battery.	• Replace battery. • Refer to Replacing the Battery section.
Monitor freezes and does not turn off.	Overall system software or hardware error.	• Press and hold the power button for several seconds to force shutdown of the ICP Monitor. • If problem persists, contact Authorized representative
Wrong date and time.	Internal back up battery defective.	• Internal battery needs replacement. • Contact Authorized representative

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Annual Preventive Maintenance and Testing is recommended to inspect for any possible wear that will degrade CereLink® ICP Monitor System's essential performance or operation, including Electrical Safety Testing per 62353:2014.

Perform the periodic maintenance procedures according to the schedule provided below:

Task	Frequency
Wiping exterior surfaces	Daily, when in use (by operator)
Sterilization of Extension cable	Before each use
Preventative maintenance	Annually (by manufacturer representative)
Replace the battery	Every 3 years or when battery performance has deteriorated (by operator)

Preventive Maintenance.

While this device does not require periodic repair or calibration, general care, handling and maintenance will extend the life of the device. Such practices the user must check include:

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- The integrity of equipment, plastic covers, LCD and interconnecting cables is sound.
- External equipment markings remain legible.
- The equipment is free of dirt or debris likely to impair safety or performance (follow the instructions for wiping the device).
- The controls and display function in accordance with the Instruction Manual.
- The rechargeable battery is charge normally.

After this check (or after the product is out of order), if any damage / unfitness for use is found, the product should be handed over to an authorized representative of the manufacturer to assess the possibility of repair. Disposable items cannot be repaired.

The manufacturer will provide schematics, parts list, and other technical information to authorized representative. For service, repair or replacements, contact your local Integra representative. Equipment under warranty should be returned to Integra LifeSciences Corporation. It will be repaired or replaced without charge to the purchaser.

Repair.

- The Extension Cable is not a repairable product; contact your local representative to order a replacement.

Replacing the Battery.

WARNING: The ICP Monitor contains a lithium ion battery:

- Do not open the battery compartment or attempt to replace the battery while monitoring a patient or while the device is connected to the AC power Supply.
- Do not puncture the battery.
- Do not dispose of in fire.

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- Do not short circuit battery contacts as battery may ignite, explode, leak or get hot.
- Use only Integra approved batteries listed in IFU. Use of another battery may present a risk of fire or explosion.

1. Remove Screws:

Remove four (4) screws from back of ICP Monitor using cross head screwdriver.

2. Remove Battery Door:

CAUTION: Do not open the battery compartment or attempt to replace the battery while the device is connected to the AC Power Supply.

Remove battery door to expose battery and connector cable.

3. Disconnect Battery:

A. Remove the old battery by gently pulling the battery from the monitors connector slots.

B. Dispose of the old battery according to local regulations.

4. Replace Battery:

Use of another battery may present a risk of fire or explosion. Plug white connector into connector slot in back of monitor.

5. Reassemble Battery Door:

Replace four (4) screws from back of ICP Monitor using cross head screwdriver.

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Wipe the exterior of devices as necessary with a soft, damp cloth or sponge. Use one of the following to remove debris, stains or adhesives that stick to the surface:

- Combination wipes (Quaternary Ammonium / Isopropyl Alcohol) or mild hospital cleaning detergent/ bactericide.
- 70% Isopropyl alcohol

Sterilization before usage:

In a properly functioning, calibrated autoclave steam sterilizer, effective sterilization can be achieved with the below parameters:

Temperature	Duration	Drying Time
132° C	4 minutes	30 minutes

The Extension Cable has been tested to withstand 100 autoclave cycles.

Cable performance beyond 100 autoclave cycles has not been evaluated. Do not subject cable to additional weight when placed in autoclave since this may compress and deform the connectors.

Note:

- Remove all bedclips (if present) before sterilization
- Do not immerse the Extension Cable in any liquid or allow liquid
- to enter the connectors. Do not use the Extension Cable if connectors become wet; erratic results may occur.

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ENVIRONMENTAL	
Operating Temperature Range	From 5° C to 40° C
Operating Humidity Range	20% to 80% non-condensing

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Shipping Temperature Range	From -10° C to 35°C, +8°C recommended
Shipping Humidity Range	15% to 85% non- condensing
Shipping Atmospheric Pressure Range	500hPa to 1060hPa

ENVIRONMENTAL

Storage Temperature Range	From 0° C to 50°C
Storage Humidity Range	15% to 85% non- condensing
Storage Atmospheric Pressure Range	500hPa to 1060hPa

Note: store battery cool and dry at 50 - 80% DOD state.

For ICP Extension Cable:

ENVIRONMENTAL

Shipping Temperature Range	From 0° C to 50°C
Shipping Humidity Range	15% to 85% non- condensing
Shipping Atmospheric Pressure Range	500hPa to 1060hPa

ENVIRONMENTAL

Storage Temperature Range	From 0° C to 50°C
Storage Humidity Range	15% to 85% non- condensing
Storage Atmospheric Pressure Range	500hPa to 1060hPa

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Although the Extension Cable is not itself subject to electrostatic discharge (ESD), it is connected to a component (Sensor) that is subject to ESD and could be damaged during a discharge.

CAUTION: Exposure to electrostatic discharge (ESD) energy could damage the Sensor. High levels of ESD could damage the electronic components and cause the Sensor to be rendered inaccurate or inoperable. Take all precautions to reduce the buildup of electrostatic charge during the use of this product.

- Provide patient grounding; e.g., grounding straps on gurneys.
- Avoid the use of materials that could generate ESD during patient movement and transport; e.g., nylon transfer boards with bedding.
- Before touching the patient, caretakers should discharge ESD buildup by touching a grounded metal surface, such as a bedrail.

It is recommended that all hospital personnel in contact with these devices

receive an explanation of the ESD symbol and training in ESD precautionary procedures. Training should include, at a minimum, an introduction to electrostatic discharge, when and why it occurs, precautionary measure, and the damage that can be done to electronic components if touched by a user who is electrostatically charged.

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this device may affect its life expectancy.

The ICP Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ICP Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ICP Monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $D = (3.5/\sqrt{P})(\sqrt{P})$	80 MHz to 800 MHz $D = (3.5/\sqrt{P})(\sqrt{P})$	800 MHz to 2.5 GHz $D = (7/\sqrt{P})(\sqrt{P})$
0.01	.012	.012	0.23
0.1	.037	.037	.074
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

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The CereLink® ICP Monitor, battery, power supply, cable, adapter is reusable, non-sterile devices. Therefore, there is no expiry date and shelf-life is not applicable for these devices.

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DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply

to, and manufacturer shall not be responsible for, any loss arising in connection with the purchase or use of any product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as, in INTEGRA's judgment, to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by manufacturer. THIS LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON INTEGRA'S PART, AND INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY REPRESENTATIVE OR OTHER PERSON TO ASSUME FOR IT ANY OTHER LIABILITY IN CONNECTION WITH INTEGRA'S PRODUCTS. INTEGRA DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESSED OR IMPLIED WARRANTY TO PATIENTS.

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Medical device consists of electrical components. Consider recycling at the end of service life. The product can be returned to

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