Certified copy

**Notar
Hartmut Pfeifer
St. Wendel**

I hereby certify that the following is a true
copy of the Copy I saw today.

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Extract N. 3 "2.2 Instructions For Use (IFU)"

Reference to:

STED (Summary of Technical Documentation) for

FX-Dialysers sterilised by steam – FX classix
Russia, Revision December 2022

Product name	Product code
FX 50 classix	F00002385
FX 60 classix	F00002386
FX 80 classix	F00002387
FX 100 classix	F00002388

The extract replaces the corresponding predecessor documents in section
"2.2 Instructions For Use (IFU)" of the original dossier

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1. Name of the Medical Device

FX classix dialyser, models:

- 1) Dialyser FX 50 classix with Instructions for Use
- 2) Dialyser FX 60 classix with Instructions for Use
- 3) Dialyser FX 80 classix with Instructions for Use
- 4) Dialyser FX 100 classix with Instructions for Use

2. Information on the Manufacturer of the Medical Device

2.1 Manufacturer

Fresenius Medical Care AG & Co. KGaA

Else-Kröner-Straße 1, 61352 Bad Homburg, Germany

Phone: +49 6172 609-0

Fax: +49 6172 609-2422

Email: <https://www.freseniusmedicalcare.com/en/contact/>

2.2 Authorized Representative of the Manufacturer in the Russian Federation
Joint Stock Company "Fresenius SP" (JSC "Fresenius SP")

Legal address: Ugreshskaya street 2, building 6, room 09

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+49 6172 609-0

www.freseniusmedicalcare.com

3. Purpose of the medical device, indicating the potential consumer

FX classix dialysers are designed for single use in chronic haemodialysis or haemo(dia)filtration and uses by trained and certified health care professionals.

Indication

The FX dialysers are designed for single use in chronic haemodialysis, haemofiltration or haemo(dia)filtration for the renal replacement therapy in case of patients with stage 5 chronic kidney diseases (CKD5).

4. Functional Characteristics and Purpose of the Medical Device

FX classix dialysers are designed for single use in chronic haemodialysis or haemo(dia)filtration respectively.

Table 1: FX classix dialysers: Summary data of in vitro performance tests

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In vitro: acc. to ISO 8637-1; in vitro data are likely to differ from in vivo data due to the patient's blood composition and clinical settings.

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5. Risk of Application of the Medical Device

Certain side effects may occur during dialysis and may result from factors specific to the patient, operating parameters, equipment, priming procedure, dialysis solution, dialyser, anticoagulation, medication etc. Therefore the selection of a dialyser and the selection and monitoring of treatment parameters based on individual patient characteristics, therapy tolerance and clinical requirements as well as compliance with the water and dialysis fluid standards are essential to minimise side effects. The FX classix dialyser is designed for high performance dialysis. In patients not treated with this dialyser before and incident patients starting HD or HDF therapy, treatment intensity shall be gradually increased to permit adequate adaptation. Hypersensitivity or hypersensitivity-like reactions have been observed during dialysis mainly in the first weeks of treatment with the dialyser. Symptomatology can vary and may include: dyspnoea, chest congestion, bronchospasm, respiratory arrest, hypotension, tachycardia, urticaria, erythema, flushing, angioedema, ocular hyperaemia, pruritus, abdominal pain, nausea, convulsions and unconsciousness. Carefully monitor patients who have not previously been treated with the dialyser, or who have shown possible hypersensitivity symptoms during previous treatments, or who have a history of allergy including asthma. If severe hypersensitivity or hypersensitivity-like reactions occur, the dialysis must be discontinued and the blood from the extracorporeal system must not be returned to the patient. Initiate appropriate emergency medical treatment.

The dialyser is steam sterilised and thus contains no sterilisation residues. For patients with known ethylene oxide hypersensitivity all components of the extracorporeal circuit should be ethylene oxide-free.

Contraindications:

Special contraindications are unknown. Generally, contraindications for haemodialysis are applicable. Patients with known hypersensitivity to any of the dialyser's material must not be treated with the dialyser.

6. Technical Characteristics of the Medical Device

Table 2: Technical characteristics of FX classix dialysers

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* For additional information see the section under "WARNINGS" in the Instructions for Use.

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8. Information on the Presence of Medicinal Products, Substances of Animal or Human Origin in the Medical Device

Not applicable.

9. Information on the Order of Installation, Setting, Tuning, Calibration and Other Actions Necessary to Put the Medical Device into Operation

It is recommended to introduce an anticoagulant to the extracorporeal circuit. Anticoagulant requirements may vary with the patient's condition, application site, dialyser characteristics and treatment modality. Nature, amount and method of application of an anticoagulant must be prescribed by the responsible physician.

HAEMOFILTRATION

Priming

Affix dialyser in a vertical position. Aseptically connect the arterial and venous blood lines to the lower and upper dialyser ports respectively. Attach the dialysate tubing connectors with the inlet to the upper and the outlet to the lower dialysate port. Fill dialyser according to the instructions provided with the machine (there is no need to turn the dialyser). We recommend bag-free priming with ONLINEplus (priming volume: 500 mL). Alternatively, the blood compartment can be filled and de-aerated in recirculation mode with isotonic saline solution (e.g. 500 mL bag). Ensure complete de-aeration of the blood compartment.

Patient Connection

Connect arterial line to the patient's circulation. Allow blood to flow in the blood line and dialyser (pump speed approx. 100 mL/min). Connect venous blood line to the patient. Adjust treatment parameters to desired setting. The weight loss should be monitored.

Termination

Perform reinfusion of the blood at the end of the treatment according to the instructions provided with the machine. We recommend the bag-free ONLINEplus procedure. Alternatively, isotonic saline solution (e.g. 500 mL bag) can be used for reinfusion. The blood should be reinfused completely.

10. Room Requirements Where the Medical Device Is to Be Installed and Requirements on the User for Installation of the Medical Device

The requirements for installation (setup) of the device for carrying out a dialysis define the room requirements for the installation of the medical device. FX classix Dialysers are intended to be used only by healthcare professionals.

11. Information for Check of Correctness of Installation (Setup) of the Medical Device and its Readiness for Safe Work of Operation

- a) Contents and frequency of maintenance, including cleaning and disinfection of the medical device

Not applicable.

- b) List provided by the manufacturer of the medical device, information, identification keys, passwords of access, programs necessary for installation, adjustment, operation and maintenance of the medical device

Not applicable.

- c) List of expendables (components, reagents) and procedure of their application and replacement

Not applicable.

- d) Necessity of calibration for ensuring appropriate and safe work of the medical device during the term of its service

Not applicable.

- e) Methods of risk minimization connected with installation, calibration or maintenance of the medical device

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WARNINGS

Due to the high hydraulic permeability of the membrane, the dialysers must only be used with dialysis machines which enable precise volumetric fluid control. Refer to dialysis machine manufacturer's technical and safety instructions for the use of high permeability membranes.

During postdilution-HDF and postdilution-HF with manual setting of the filtration flow rate the maximum filtration flow (substitution rate + ultrafiltration rate) must not exceed 25 % of the effective blood flow rate. This is to keep the blood water content within a certain uncritical range and to avoid the risk of haemoconcentration and coagulation in the extracorporeal circuit. The TMP has to be monitored.

The dialyser is intended for single-use only. Re-use may be hazardous to both the patient and operator. Cleansing solutions and disinfectants may damage materials employed for the housing, potting and membrane. Safety of use can no longer be guaranteed and the manufacturer assumes no liability.

Blood flow rates below the recommended blood flow range can be set. The coagulation should be monitored and the decrease of the clearance values might be steeper than expected.

Blood flow rates above the recommended blood flow range can be set. This might result in a decrease of the sieving coefficient of middle molecular substances (e.g. β_2 - m, myoglobin) and an increase of the sieving coefficient of albumin when performing HDF or HF. The increase of the middle molecule clearance might be lower than expected from raised blood flow and filtration flow rates.

WARRANTY

Products with manufacturing defects will be replaced if the defect is reported with details of the lot number.

The manufacturer will not be liable for any misuse, improper handling, non-compliance with instructions for use and cautionary notes and for any damage incurred subsequent to the manufacturer's delivery of the dialyser.

- f) Information on installation, adjustment (setup), control, calibration and other actions necessary to put the medical device into operation and for its correct application

For more information please refer to section 11 subsection e) "Methods for risk minimisation in the context of installation, calibration or maintenance of the medical device".

- g) Information on the list of the main operation characteristics (application) of the medical device, on transport conditions and storage

The user has to transport and store the FX classix dialysers steady against the influence of climatic factors, i.e. from +5 °C - +30 °C. The shelf-life is three years, calculated from the date of manufacture (see also use-by date on the product and box label).

For more information, please refer to product or box label for:

REF	Catalogue number	Use-by date
Date of manufacture	Do not re-use	
Manufacturer	Temperature limit	
Sterilized using steam. Blood pathway sterile.	Use only on machines with exact UF-control.	
Caution	Units	
Non-pyrogenic	Do not use if the package is damaged	
Keep dry	Fragile, handle with care	
LOT	Batch code	

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h) List of national standards applied by the manufacturer of the medical device

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* Additional information according to paragraph 6.2.1, 6.2.2, 6.2.3 of GOST ISO 8637-2012 can be provided to the user upon request.

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12. Sterility of the Medical Device

The dialyser is steam sterilised and thus contains no sterilisation residues.

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14. Combined Use of the Medical Device with Other Devices

FX classix dialysers are designed for use during kidney replacement therapy in combination with haemodialysis machines. Use only devices with exact ultrafiltration control (UF control).

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1) Recommended lines:

– Consumables for the device "Artificial kidney": bloodlines and dialysis accessories, Reg. cert. № FSZ 2010/07339;

2) Recommended machines:

– Apparatus "Artificial kidney" type 4008, models 4008 B, 4008 E, 4008 H, 4008 S, 4008 ADS, Reg. cert. no. RZN 2014/1643;

– Apparatus "Artificial kidney" type 5008 with accessories, Reg. cert. No. RZN 2014/1645;

– Apparatus "Artificial kidney" type 5008S with accessories, Reg. cert. No. FSZ 2008/01823.

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Not applicable.

16. Safety

- a) Malfunctions of the medical device, failure in its work or deviations in functioning which can influence the safety of the medical device, including those that can be identified by external signs

Use only if unit package is intact, sealing caps are in place and the dialyser is undamaged. Dialysers must not be used after expiry date (see label). Each dialyser is checked for integrity prior to leaving the factory. If a blood leak should arise, the dialyser must be exchanged.

- b) Impacts on functioning of the medical device of the external factors connected with application of the medical device in a combination with other medical devices and (or) the equipment, or such predictable factors as external electromagnetic fields, electrostatic categories, radiation (electromagnetic, ionizing, other), the atmospheric pressure and its differences, humidity and environmental temperature

FX classix dialysers produced by Fresenius Medical Care are intended for use by healthcare professionals (doctors and nurses), conditions of the environment and the requirements are strictly established by clinical experts according to specific instructions of the medical institution.

- c) Risk of the electromagnetic interference created by a medical device for other medical devices, the equipment and means of communication during the carrying out and assessment of results of diagnostics, treatment or at its application to destination (for example, the electromagnetic radiation of a medical device exerting impact on other equipment)

Not applicable.

17. Warnings and/or Precautions Taken by the User When Using a Medical Device Containing a Medicinal Product, Material of Animal or Human Origin, Materials That Are Carcinogenic, Mutagenic or Toxic, Possible Isolation or Leaching of Which Leads to Sensibilization, Allergic Reaction or Negatively Affects Reproductive Function

The dialyser is steam sterilised and thus contains no sterilisation residues. For patients with known ethylene oxide hypersensitivity all components of the extracorporeal circuit should be ethylene oxide-free. For more information on warnings and precautionary measures please refer to section 11 subsection e) "Methods for risk minimisation in the context of installation, calibration or maintenance of the medical device".

18. Warnings and Precautions Taken by the User when Disposing the Medical Device, Accessories and Consumables Used Together with it (if any), Including Information about the Infectious, Microbial, Environmental or Physical Hazards of the Medical Device

After the use FX classix Dialysers can be potentially dangerous to patients, healthcare professionals and the environment. Information on the correct disposal of the waste materials can be found in "Procedure and conditions for disposal or destruction of medical device".

19. Consultation of a Healthcare Professional

If you observe any side effects, please report them immediately (see also "Risk of Application of the Medical Device"). The healthcare professional will initiate appropriate medical treatment.

20. Information on this Document

Initial release of the text: 2022-11-17

Date of revision of the text: 2022-11-17

21. Procedure and conditions for disposal or destruction of medical device

Disposal or destruction of medical devices is carried out in accordance with the classification, rules for collection, use, neutralisation, placement, storage, transportation, accounting and disposal of medical waste established by the authorised federal executive body.

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Warning and / or special precautions for the safe disposal of a medical device and its accessories:

- a) Infectious or microbial risks, including the possibility of contamination of consumables with infectious agents of human origin:

Single-use medical devices that have come into contact with human blood or other biological fluids and can potentially be infectious. Used disposable products are subject to collection in containers of different colors, depending on the class of medical waste to which they belong.

- b) Environmental risks associated with potentially hazardous materials and substances:

Not applicable.

- c) Physical risks, including the possibility of explosion or fire:

Not applicable.

Disposal of used dialyzers in contact with the internal environment of the body is carried out in accordance with the procedure provided for by SanPiN 2.1.3684-21 for Class B waste. If dialyzers were used to assist patients infected (or suspected to be infected) with epidemiologically dangerous infections, then the disposal of dialyzers should be carried out in accordance with the procedure provided for by SanPiN 2.1.3684-21 for Class V waste. Disposal of an expired product that has not been in contact with the patient (who have not had contact with the biological fluids of patients, infectious patients) is carried out in accordance with the procedure provided for by SanPiN 2.1.3684-21 for Class A waste. Disposal of packaging is carried out in accordance with the procedure provided for Class A waste. The requirements of Sanitary Rules for the treatment of solid municipal waste are applied to the treatment of medical waste of Class A.

Transportation and disposal of disposable medical devices, as well as other types of medical waste, must be carried out by a licenced company.

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