
Certified copy

**Notar
Hartmut Pfeifer
St. Wendel**

I hereby certify that the following is a true
copy of the Copy I saw today.

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**FRESENIUS
MEDICAL CARE**

**Addition to
the Instruction for use
for multiFiltrate PRO tubings
(July-2022)
Date of Issue:
May-2023**

Addition to the Instructions for Use for

for multiFiltratePRO tubings (July-2022)

*(these additions apply to STED and extract from STED No. 2.
Only for Russian Federation)*

St. Wendel, 22-May-2023

**FRESENIUS
MEDICAL CARE**
FRESENIUS MEDICAL CARE AG & Co. KGaA
Else-Kröner-Str. 1
61352 Bad Homburg, Germany

M. Wenke-Pfeifer

I. V. Marina Wenke-Pfeifer
Senior Manager Regulatory Affairs Medical Devices
Fresenius Medical Care AG & Co. KGaA
Else-Kröner-Str. 1
61352 Bad Homburg, Germany

C. Betz

I. V. Dr. Christian Betz
Expert Regulatory Affairs Medical Devices
Fresenius Medical Care AG & Co. KGaA
Else-Kröner-Str. 1
61352 Bad Homburg, Germany

Addition to the Instruction for use for for multiFiltratePRO tubings (July-2022)

*(these additions apply to STED and extract from STED No. 2.
Only for Russian Federation)*

St. Wendel, 05.06.2023

I.V. Marina Wenke-Pfeifer
Senior Manager Regulatory Affairs Medical Devices

Fresenius Medical Care AG & Co. KGaA
Else-Kröner-Str. 1, 61352 Bad Homburg, GERMANY

I.V. Christian Betz
Expert Regulatory Affairs Medical Devices

Fresenius Medical Care AG & Co. KGaA
Else-Kröner-Str. 1, 61352 Bad Homburg, GERMANY

1. Name of the medical device

multiFiltratePRO tubings are tubing systems for the extracorporeal blood purification. The variants of the multiFiltratePRO tubings are as follows:

multiFiltratePRO tubings are tubing systems for the extracorporeal blood purification. The variants of the multiFiltratePRO tubings are as follows:

1. multiFiltratePRO Ci-Ca HD (both heparin and citrate-calcium anticoagulation; for HD therapy), composed by:
--

Page 1 / 15

1	arterial line with citrate line and heparin line - 1 pc
2	venous line with calcium line - 1 pc
3	filtrate line - 1 pc
4	dialysate line - 1 pc
5	tray to facilitate handling - 1 pc
6	IFU - 1 pc
7	recirculation connector - 1 pc
2. multiFiltratePRO Ci-Ca HDF (both heparin and citrate-calcium anticoagulation; for HDF therapy), composed by:	
1	arterial line with citrate line and heparin line - 1 pc
2	venous line with calcium line - 1 pc
3	filtrate line - 1 pc
4	dialysate line - 1 pc
5	substitute line - 1 pc
6	tray to facilitate handling - 1 pc
7	IFU - 1 pc
8	redirculation connector - 1 pc
3. multiFiltratePRO HDF (heparin anticoagulation; for HDF therapy), composed by:	
1	arterial line with heparin line - 1 pc
2	venous line - 1 pc
3	filtrate line - 1 pc
4	dialysate line complete with substitute line - 1 pc
5	tray to facilitate handling - 1 pc
6	IFU - 1 pc
7	recirculation connector - 1 pc
4. multiFiltratePRO TPE (Therapeutic Plasma Exchange therapy), composed by:	
1	arterial line with heparin line - 1 pc
2	venous line - 1 pc
3	filtrate line - 1 pc
4	plasma line (plasma line is a mix of the dialysate and substitute line) - 1 pc
5	tray to facilitate handling - 1 pc
6	IFU - 1 pc
7	recirculation connector - 1 pc
Abbreviation:	
IFU: Instruction for Use; Pc: piece(s)	

2. Purpose

The multiFiltratePRO tubings are designed to be used with the active medical device multiFiltratePRO for extracorporeal blood purification treatment. The main treatment of the multiFiltratePRO is the Continuous Renal Replacement Therapy (CRRT), specifically used in ICUs, enabling the removal of toxins by convection, diffusion or a combination of both.

- multiFiltratePRO Ci-Ca HD extracorporeal blood purification tubing system with heparin and regional citrate anticoagulation for continuous veno-venous hemodialysis (CVVHD).
- multiFiltratePRO Ci-Ca HDF extracorporeal blood purification tubing system with heparin and regional citrate anticoagulation for continuous veno-venous hemodiafiltration (CVVHDF).

- multiFiltratePRO HDF extracorporeal blood purification tubing system with heparin anticoagulation for continuous veno-venous hemodialysis (CVVHD); for continuous veno-venous hemofiltration (CVVH) with pre- or post-dilution; for continuous veno-venous hemofiltration (CVVH) with pre-post dilution; for continuous venous-venous hemodiafiltration (CVVHDF).
- multiFiltratePRO TPE extracorporeal blood purification tubing system with heparin anticoagulation for therapeutic plasma exchange (TPE).

3. Scope of delivery

multiFiltratePRO tubings are tubing systems for the extracorporeal blood purification. The variants of the multiFiltratePRO tubings are as follows:

multiFiltratePRO HDF (heparin anticoagulation; for HDF therapy), composed by:	
1	arterial line with heparin line - 1 pc
2	venous line - 1 pc
3	filtrate line - 1 pc
4	dialysate line complete with substitute line - 1 pc
5	tray to facilitate handling - 1 pc
6	IFU - 1 pc
7	recirculation connector - 1 pc
multiFiltratePRO CI-Ca HDF (both heparin and citrate-calcium anticoagulation; for HDF therapy), composed by:	
1	arterial line with citrate line and heparin line - 1 pc
2	venous line with calcium line - 1 pc
3	filtrate line - 1 pc
4	dialysate line - 1 pc
5	substitute line - 1 pc
6	tray to facilitate handling - 1 pc
7	IFU - 1 pc
8	recirculation connector - 1 pc
multiFiltratePRO CI-Ca HD (both heparin and citrate-calcium anticoagulation; for HD therapy), composed by:	
1	arterial line with citrate line and heparin line - 1 pc
2	venous line with calcium line - 1 pc
3	filtrate line - 1 pc
4	dialysate line - 1 pc
5	tray to facilitate handling - 1 pc
6	IFU - 1 pc
7	recirculation connector - 1 pc
multiFiltratePRO TPE (Therapeutic Plasma Exchange therapy), composed by:	
1	arterial line with heparin line - 1 pc
2	venous line - 1 pc
3	filtrate line - 1 pc
4	plasma line (plasma line is a mix of the dialysate and substitute line) - 1 pc
5	tray to facilitate handling - 1 pc
6	IFU - 1 pc
7	recirculation connector - 1 pc
Abbreviation:	

4. List of national (in Russian Federation) standards applied by the manufacturer:

Decree of the Government of the Russian Federation of December 27, 2012 № 1416 "On approval of the rules of state registration of medical devices.

Ministry of Health of the Russian Federation № 885n of 30.08.2021

"Ministry of Health of the Russian Federation. Order. On Approval of the Procedure for Conformity Assessment of Medical Devices in the Form of Technical Tests, Toxicological Investigations, and Clinical Tests for the Purpose of State Registration of Medical Devices".

Federal Law No. 323-FZ of 21.11.2011 "On the Fundamentals of Health Protection in the Russian Federation".

GOST R ISO 14155-2014 Clinical trials. Good Clinical Practice

GOST ISO 8638-2012

Extracorporeal blood circuit for haemodialysers, haemodiafilters and haemoconcentrators. Technical requirements and test methods

GOST ISO 11607-1-2018

Packaging for terminally sterilized medical devices. Part 1. Requirements for materials, sterile barrier systems and packaging systems

GOST R IEC 62366-1-2021

Medical devices. Part 1. Application of usability engineering to medical devices

GOST R ISO 15223-1-2020

Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements

GOST 31214-2016

Medical devices. Requirements for samples and documentation presented for toxicological tests, sanitary and chemical analyses, tests for sterility and pyrogenicity

GOST ISO 10993-12-2015

Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials

GOST R 52770-2016

Medical devices. Safety requirements. Methods of sanitation-chemical and toxicological tests.

GOST ISO 10993-4-2020

Medical devices. Biological evaluation of medical devices. Part 4. Selection of tests for interactions with blood.

GOST ISO 10993-5-2011

Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity.

GOST ISO 10993-10-2011

Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and delayed-type hypersensitivity.

GOST ISO 10993-11-2021

Medical devices. Biological evaluation of medical devices. Part 11. Tests for systemic toxicity.

Russia State Pharmacopoeia 13 (XIII), GPM.1.2.4.0005.15

Pyrogenicity

5. Disposal

Based on the characteristics of the morphological composition, the hazard class of this medical device waste is defined as epidemiologically hazardous waste class B (SanPiN 2.1.3684-21).

Any use of materials is allowed only if there is a hyperlink to the Sanitary and Epidemiological Requirements for the Handling of Medical Waste).

The used product may constitute a biologically contaminated product and must be disposed of in accordance with current local legislation.

The packaging is to be disposed of as domestic waste.

Unused product and packaging are not harmful to the environment and are disposed of in accordance with local legislation as solid domestic waste.

6. Indications and contraindications**Medical indication**

- Acute kidney injury requiring continuous renal replacement therapy
- Volume overload requiring continuous renal replacement therapy
- Intoxications requiring continuous renal replacement therapy
- Conditions requiring the exchange of blood plasma via TPE
- Conditions requiring CRRT combined with hemoperfusion to eliminate further pathogens from blood
- Conditions requiring CRRT combined with ECCO2R for additional CO2 Removal

Product related contraindications

Hypersensitivity to any of the material of the tubings as listed under composition in the IFU or to ethylene oxide residues.

In severe cases treatment must be discontinued, blood must not be returned to the patient, and appropriate emergency medical care must be initiated. In severe cases of hypersensitivity, terminate the treatment, do not return the blood to the patient, and initiate appropriate emergency medical treatment.

Therapy related contraindication

- The contraindications of all the tubing systems, filters, dialysers and solutions intended for use in a therapy must be observed.
- Inadequate vascular access to the venous blood circuit, e.g., with a double-lumen central venous catheter, so that the blood flow necessary for treatment cannot be achieved.

Systemic anticoagulation is often contraindicated for patients with bleeding or a high risk of bleeding. Here, measures should be taken within the treatment regimen which permit a renal replacement treatment with little or no anticoagulation. A regional citrate anticoagulation can be applied where appropriate.

An established disordered citrate metabolism is an absolute contraindication.

7. Intended patient population

The devices have been specified by the manufacturer for the purpose of treating patients with a body weight of 40 kg and more, irrespective of their age, under consideration of the specified technical data of the device and the single-use items used (e.g. delivery rates, fill volumes).

8. Intended user group and intended environment

The disposable must only be used by individuals with the appropriate training, knowledge, and experience on the proper operation and handling and for whom proof of instruction can be shown (e.g. physicians, nurses/paramedics, anesthesiology/resuscitation, transfusion medicine). Operate in rooms suitable for dialysis located in professional healthcare facilities. Normative and local regulations must be observed.

9. Description of the base tray for ease of use

The tray ensures an easy handling of the product from the user's point of view. On the tray, all the parts of the multiFiltratePRO tubing set are fixed, some in a definitive way, some other in a temporary way. It is designed to be connected with and fixed to the multiFiltratePRO Machine.

Photo: Basic tray for ease of use



Photo: example of tray on the dialysis machine

Characteristics of the tray

	Dimension
Width	495 mm $\pm$ 5%
Height	480 mm $\pm$ 5%
Weight	120 g $\pm$ 5%

10. Description of the recirculation connector

Pre-packaged recirculation connector, used for recirculation during the preparation phase after filling/priming to connect the arterial and venous lines included in the Tubing system for extracorporeal blood purification multiFiltratePRO.

Figure: drawing of the recirculation connector



Figure: Photo of the recirculation connector

Characteristics of the recirculation connector

	Dimension
Length with cap	43,5 mm $\pm$ 5%
Length without cap	42 mm $\pm$ 5%
Internal diameter	4,3 mm $\pm$ 5%
External diameter	6,8 mm $\pm$ 5%
Force of the plug opening	Not more than 5 N
Thread pitch	2,1 mm $\pm$ 5%
Weight	1.23g $\pm$ 5%

11. Operating conditions

Operating characteristics are given according to the information in the Instructions for Use of the multiFiltratePRO extracorporeal detoxification machine with accessories, RU No. FSZ 2010/08757 of 11.01.2023.

Return line pressure measurement	The hydrophobic filter in the return line is determined to be fully wetted when the return pressure sensor (blue) fails to detect any pressure fluctuations over at least 4 consecutive blood pump cycles. Pressure monitoring is deactivated for blood flow rates of less than 50 ml/min.
Access pressure	Display range: -300 to +300 mmHg Resolution: 5 mmHg Accuracy: 10 mmHg No blood detected: Access pressure alarm window size: -300 to +300 mmHg Blood detected: Access pressure alarm window size: +40 to +200 mmHg Default value adjustable in User setup, Factory setting: +200 mmHg If the access pressure drops below the lower limit, the access line clamp will remain open to allow the pressure in the system to disperse. In the event of any further pressure alarm, the clamp will close.
Return pressure (safety system against external blood leaks)	Display range: -100 to +500 mmHg Resolution: 5 mmHg Accuracy: 10 mmHg No blood detected: Return pressure alarm window size: -100 to +500 mmHg Blood detected: Return pressure alarm window size: +40 to +200 mmHg Default value adjustable in User setup, Factory setting: +100 mmHg Alarm window position can be set over a range from +10 to +500 mmHg (switchover to between -100 and +500 mmHg possible in the event of an alarm, if extending low return pressure alarm is activated in Service setup) Factory setting: extending lower return pressure limit is deactivated.
Pre-filter pressure	Display range: -50 to +750 mmHg Resolution: 5 mmHg Accuracy: 10 mmHg No blood detected: Pre-filter pressure alarm window size: -50 to +750 mmHg

TMP (CRRT)

Blood detected:

Pre-filter pressure alarm window size: +40 to +200 mmHg

Default value adjustable in User setup,

Factory setting: +200 mmHg

Display range: -300 to +500 mmHg

Lower alarm limit: -60 mmHg

Upper alarm limit: +520 mmHg

Accuracy: 20 mmHg

Displayed only in Pressure / alarm history tab of treatment screen.

The TMP is calculated and displayed according to the following formula:

$$\text{TMP} = ((P_{\text{ven}} + P_{\text{pre F}}) / 2) - P_{\text{fil}} + \text{Offset}$$

TMP = transmembrane pressure

P_{ven} = return pressure

$P_{\text{pre F}}$ = pre-filter pressure

P_{fil} = filtrate pressure

Offset = 20 mmHg (correction value to compensate for hydrostatic pressure differences)

TMP (TPE)

Display range: -60 to +270 mmHg

Pressure alarm windows

Lower alarm limit: -60 mmHg

Upper alarm limit: +50 mmHg to maximum upper alarm limit

Maximum upper alarm limit can be defined in User setup between +50 and +100 mmHg

Accuracy: 20 mmHg

The TMP is calculated and displayed according to the following formula:

$$\text{TMP} = ((P_{\text{ven}} + P_{\text{pre F}}) / 2) - P_{\text{fil}} + \text{Offset}$$

TMP = transmembrane pressure

P_{ven} = return pressure

$P_{\text{pre F}}$ = pre-filter pressure

P_{fil} = filtrate pressure

Offset = 20 mmHg (correction value to compensate for hydrostatic pressure differences)

Blood pump

Spring-loaded rollers, fully occluding, pressure-limited to 2 bar for standard line with pump segment 6.4 x 1.8 (when using the prescribed tubing systems).

Warning

Risk for the patient as a result of insufficient detoxification

If the access pressure before the blood pump reaches extreme negative values, the blood flow can be reduced, which will impair the effectiveness of the treatment.

> Take suitable steps to avoid an extreme negative access pressure.

Delivery rate:

CRRT: 10 to 500 ml/min

CRRT with citrate anticoagulation: 10 to 200 ml/min

	TPE: 10 to 300 ml/min Resolution: 10 ml/min Flow accuracy over Pressure range ≥ -300 mmHg ≤ 10 % Standard line with pump segment 6.4 x 1.8 mm System accuracy of the delivered blood volume: ± 10 % considered over the entire treatment duration and valid in typical treatment situations. Blood pump stop alarm: time-based standstill monitoring as a safety system against blood loss through clotting. Alarm delay when blood pump stops: 1 minute (during treatment) 3 minutes (while patient is being connected or disconnected) Alarm repeat if blood pump standstill continues: every 60 s
Fill level detector	Method: Capacitive measurement Switching point 13 mm, ± 4 mm from upper edge
Optical detector	Method: Infrared transmission Distinguishes between: – No blood detected (NaCl solution or air in tubing) – Blood detected (blood in tubing)
Air bubble detector	Method: Ultrasound transmission measurement through tubing Detects: – Air bubbles – Blood foam (air/blood mixture) – Microbubbles Air alarm in the following cases: – Blood foam – Microbubbles – Blood flow rate < 100 ml/min: Air bubble: volume ≥ 20 μl – Blood flow rate ≥ 100 ml/min: 10 air bubbles with a bubble volume of < 50 μl each, or 1 air bubble with a bubble volume of ≥ 50 μl, or blood foam The above data is based on a worst-case assumption, at a blood flow rate of 0 to 500 ml/min using the prescribed tubing systems. Full sensitivity at the maximum blood flow is achieved with patients weighing upwards of 45 kg. In order to ensure a similar sensitivity with patients weighing less than 45 kg in a worst-case scenario (level in bubble catcher has dropped), select a lower maximum blood pump rate according to the following table.

General limit value: 0.03 (ml/min) per kg		
Patient weight	Max. infused air for lowest possible hazard	Limited max. blood flow (condition: wetted)
40 kg	1.2 ml/min	458 ml/min
From 45 kg	≥ 1.35 ml/min	≥ 500 ml/min

Heparin pump

Pump type: syringe pump
 Delivery rate: 0.5 to 25 ml/h
 Resolution: 0.1 ml/min
 Accuracy: $\pm 5\%$ for delivery rate 1 to 25 ml/h measured over 2 hours with up to 1.2 bar counterpressure. At delivery rates < 1.0 ml/h, the tolerance could exceed the specified $\pm 5\%$.
 Bolus administration: 0.1 to 5 ml in increments of 0.1 ml (preset maximum bolus volume is 5 ml. This parameter can be set to a lower volume in the System Parameters).
 Bolus rate: 30 ml/min

Audible tone

Sound pressure level settings of the audible alarm:
 Volume range: 50 to 80 dB ± 5 dB
 Factory setting: ≥ 65 dB
 High-priority alarm: 60 to 80 dB ± 5 dB
 Medium priority alarm: 60 to 80 dB ± 5 dB

Ci-Ca drip counter

Measuring range: 0 to 5 drips per second (independently for citrate and calcium)
 Measuring method: optical
 To permit drips to be detected accurately, the fluid level must be within or below the markings.

Citrate pump

Pump type: roller pump
 Delivery accuracy: $\pm 10\%$
 Delivery rate: 10 to 600 ml/h, depending on citrate/blood ratio.
 Dose can be set.
 Concentration of citrate per litre of delivered blood: 2 to 6 mmol/l in 0.1 mmol/l increments
 Default value: 4.0 mmol/l

Calcium pump

Pump type: roller pump
 Delivery accuracy: $\pm 10\%$, at delivery rates < 6 ml/h the deviation can be $\pm 20\%$
 Delivery rate: off, 1 to 100 ml/h, depending on calcium/filtrate ratio.

Dose can be set.
 Concentration of calcium per litre filtrate: 0 to 3 mmol/l in 0.1 mmol/l increments
 Default value: 1.7 mmol/l

The Ci-Ca pumps run at a higher delivery rate (400 ml/h) while the Ci-Ca tube segments are being inserted/removed and the tubing system is being filled.

Ci-Ca fill level detector

Function: for detecting and differentiating between a full or empty Ci-Ca drip chamber (independently for citrate and calcium).
 Measuring method: optical

Cassette detector

To permit a filled drip chamber to be detected accurately, the fluid level must be within or above the markings.

Differentiates between cassettes with and without a Ci-Ca system using a colour sensor and colour codes on the cassettes.

Cassette without Ci-Ca: blue marking

Ci-Ca cassette: yellow marking

12. Operating time of tubing systems

The maximum application time is 72 hours or 1,250 liters. The disposable must be replaced after the maximum usage time or respectively after a specific blood volume pumped as indicated on the primary packaging.

13. Flow rate of fluids

Flow rate	Depends on the method of procedure	
	Blood Flow	0/10 to 500 ml/min $\pm$ 10%
	Substitute flow	0/10 to 80 ml/min
	Dialysate Flow	0/10 to 80 ml/min
	Citrate Flow	0/10 to 600 ml/h
	Calcium flow	0/1 to 100 ml/h
	Ultrafiltration rate (UF rate)	0/10 to 180 ml/min
	Filtrate Flow	0/10 to 180 ml/min
	Net UF rate	0/10 to 990 ml/h

14. Priming Volumes

The filling volumes (priming volumes) of each finished product are as follows:

	multiFiltratePRO Ci-Ca HD	multiFiltratePRO Ci-Ca HDF	multiFiltratePRO HDF	multiFiltratePRO TPE
Priming volume (ml)	116 $\pm$ 5%	116 $\pm$ 5%	116 $\pm$ 5%	117 $\pm$ 5%

15. Shelf life

The product has a shelf life of 3 years.

16. Pressure limit value

	Range
Blood pressure (inlet)	- 300 ÷ 300 mm Hg.
Arterial pressure (out)	- 50 ÷ 750 mm Hg.
Venous pressure	- 100 ÷ 500 mm Hg.

17. Technical Specifications

Limit value of blood flow	500 ml/min
----------------------------------	------------

Connections to auxiliary devices	Fracture force of connection not less than 100 N Elements for connection with auxiliary devices are terminated with a threaded connection with an inner cone of 6%
Connection to the injection unit	The injection port must remain leak-tight after being pierced six times with a 0.8 mm inner diameter needle (No. 21)
Connection to a vascular access device	The connection of the trunk to the vascular access device must be provided by the "Luer-lock" type coupling with 6% taper.

18. Weight and dimensions specifications

The weight/dimension specifications are not necessary for the operation of the appliance, and therefore are not required to be mentioned in the Instructions for Use. This information can be provided to the user on request.

19. Connection diagrams

Schemes (drawings), explaining the connection of the lines to the device to enable their safe use and reduce the risks for the consumer are shown in the Operating Instructions for the "multiFiltratePRO extracorporeal detoxification device with accessories" (RU No. FSZ 2010/08757 of 11.01.2023) item. 4 "Operation" and item 7 "Functional description", and are not required in the Instructions for Use for the multiFiltratePRO extracorporeal blood detoxification tubing system.

20. The tubing systems are compatible with the following devices in use:

The tubing systems are compatible with the following devices in use:

- 1 multiFiltratePRO extracorporeal purification apparatus with accessories RU No. FSZ 2010/08757 of 11.01.2023.
- 2 Gas exchange module multiFiltrate ECCO2R (undergoing the registration procedure as part of the registration of a new software version of multiFiltratePRO extracorporeal purification apparatus).
- 3 multiEffluent Bag 10 L (under registration);
- 4 Ultraflux dialyzer (hemofilter) RU № FSZ 2010/07931
- 5 Solution for hemofiltration Multilak, 2 mmol/l Potassium, RU № LP-000706 from 28.09.2011
- 6 Moisture for hemofiltration Multilak, 4 mmol/l Potassium, RU № LP-000706 from 28.09.2011
- 7 Liquid for hemofiltration Multibik, 2 mmol/l Potassium, RU № LP-000706 from 28.09.2011
- 8 Solution for hemofiltration Multibik, 4 mmol/l Potassium, RU № LP-000706 from 28.09.2011
- 9 Dialysis solution Multiplus K+2 mmol/L, 5000 mL, RU № FSZ 2012/12411 from 25.07.2012
- 10 Dialysis solution (dialysate) for hemodialysis with citrate-calcium anti-coagulation (RU № FSZ 2011/09209 from 25.02.2011):
 - Ci-Ca Dialysate K 2, 5000 ml. - Ci-Ca Dialysate K2 5000 ml - Dialysate solution (dialysate) for hemodialysis with citrate-calcium anti-coagulation, potassium content 2 mmol/l
 - Ci-Ca Dialysate K 4, 5000 ml. - Ci-Ca Dialysate K4 5000 ml - Dialysate solution (dialysate) for hemodialysis with citrate-calcium anti-coagulation, potassium content 4 mmol/l
- 11 Dialyzer (Hemofilter) Ultraflux (Ультрафлюкс). Versions (RU No. FSZ 2010/07931 of 22.09.2010):
 - AV 400 S

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12 PlasmaFlux P1 dry, plasmaFlux P2 dry. RU № RZN 2015/2389 of 03.12.2020.

13 Double-lumen catheter proVencare in kits (RU № RZN 2017/6656 of 29.12.2017):

- Double-lumen catheter in kits (available sizes, diameter, F/length, mm: 6,5/75; 6,5/100; 6,5/125; 8/100; 8/120; 8/150; 11/150; 11/200; 11/250; 12/150; 12/200)
- High-flow dual-lumen catheters in kit (available sizes, diameter, F/length, mm: 11,5/150; 11,5/200; 11,5/250; 13/150; 13/175; 13/200; 13/250).

14 Any three-component syringe for pump or injector that is registered in the territory of the Russian Federation. We most commonly use the following medical product:

Single-use syringes with and without needles: 1 ml, 2 ml, 2.5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 30 ml, 50 ml, 60 ml, 100 ml, 150 ml. (RU No. FSZ 2010/08748 of 31.12.2010).

21. Information about sterilization date on the label

The sterilization process of medical devices is carried out no later than within 30 calendar days from the date of their production.

The date of sterilization of medical devices is traceable based on the production date in label. Thus, the requirements for packaging, labeling and operational documentation of clause 6.1.1 of GOST ISO 8638-2012 "Extracorporeal blood circuit for haemodialysers, haemodiafilters and haemoconcentrators. Technical requirements and test methods" regarding the date of sterilization are fulfilled.

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