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Amep

4

INSTRUCTION FOR USE

Version 3

Date: March 28, 2023

1. Name

NITHYA sterile intradermal implant based on equine collagen Type I, models:

- 1) NITHYA implant 3x70 mg including:
 - glass vial with 70 mg of collagen powder - 3ea
 - instruction for use - 1 ea
- 2) NITHYA implant 200 mg.
 - glass vial with 200 mg of collagen powder - 1ea
 - instruction for use - 1 ea.

2. Composition of medical device and accessories

The medical device is composed at 100% of collagen powder, sterile; pyrogen-free. The limit content of bacterial endotoxins not more than 0.5 U/ml.

The CAS number of collagen powder is: 9007-34-5.

Type I collagen is characterized by its complete and complex molecular form: a protein in the typical triple-helical structure.

3. Intended use

The medical device is indicated for injection into the mid-deep dermis. The medical device helps restore optimal connective tissue conditions after traumatic events and post-traumatic lesions. It is indicated as an adjuvant in skin biorevitalization (dermo-regeneration) and to promote the regeneration of the connective tissue of the dermis by creating optimal conditions for the physiological neoformation of collagen. Also indicated in the treatment of facial and body defects related to ageing over time (chrono and photo-ageing) and for the remodeling of genetic defects such as, for example, skin hypotrophy and hypotonia, reduced skin elasticity.

4. Packing set

- Each pack of NITHYA sterile intradermal implant based on equine collagen Type I, model NITHYA implant 3x70 mg contains n.3 of glass vials each containing 70 mg of equine Type I collagen powder, sterile. The packaging provides 1 leaflet, too.
- Each pack of NITHYA Powder of equine collagen Type I, sterile model NITHYA implant 200 mg contains n.1 of glass vial containing 200 mg of equine Type I collagen powder, sterile. The packaging provides 1 leaflet, too.

The needles and syringes are not provided in package, but 30-32 G ½" - 4-13 mm needles are recommended. For NITHYA sterile intradermal implant based on equine collagen Type I, model 3x70 mg 5 syringes of 1 mL are recommended; for NITHYA sterile intradermal implant based on equine collagen Type I, model 200 mg - 5 to 10 syringes (depending on the treated area) of 1 mL syringes are recommended.

The choice of syringes and needles is at discretion of physician.

Packaging materials

No.	Item	Material	Composition	Grade	Size	Manufacturer
1.	Vial capacity 10 ml	Glass	SiO ₂ 71-73% Al ₂ O ₃ 2-4% CaO 10-12% Na ₂ O 0-3% K ₂ O 12-14% P ₂ O ₅ 0-1% MgO 0-1%	Yellow glass type 1	height 53 mm ± 1 mm; the diameter of body 25 mm ± 1 mm; outer diameter of the neck 20 mm ± 1.0 mm; inner diameter of the neck 13.6 mm ± 0.5 mm;	Shvets BAS Via G. Turin, 6 28158 Milan - I
2.	Drug	Rubber	Chlorobutyl	PA1936	cap diameter 30 mm ± 1 mm; base diameter 13.0 mm ± 0.5 mm; height 9.0 mm ± 0.5 mm;	PHARMA TRADE S Zone Ind. Loc. Poin. 66054 V/ Italy
3.	Flip-off cap	Aluminium	Aluminium alloy	6011	outer diameter (assembly) 23.6 mm ± 0.5 mm; height (assembly) 6.0 mm ± 1.0 mm; diameter 20 mm ± 0.5 mm;	PHARMA TRADE S Zone Ind. Loc. Poin. 66054 V/ Italy
	cap	Plastic	Polypropylene	HP456H	outer diameter 23.0 mm ± 0.5 mm; inner diameter 20.0 mm ± 0.5 mm; height 1.0 mm ± 0.5 mm;	

5. Technical characteristics

- Weight of the powder:
 - for NITHYA sterile intradermal implant based on equine collagen Type I, model NITHYA implant 3x70 mg: 70 mg ± 5 mg for each vial;
 - for NITHYA sterile intradermal implant based on equine collagen Type I, model NITHYA implant 300 mg: 280 mg ± 10 mg for each vial.
- Total weight of finished vial:
 - for NITHYA sterile intradermal implant based on equine collagen Type I,

model 3x70 mg: 25g ± 10%;

-for NITHYA sterile intradermal implant based on equine collagen Type I, model 300 mg: 25g ± 10%;

Size of collagen powder particles: 0.020 µm to 2000.00 µm

6. Risk class

Medical device of Class III potential risk.

7. Indications

- Skin and dermis rejuvenation;
- Dermal regeneration after trauma or posttraumatic lesions;
- Treatment of age-related skin imperfections affecting face, neck, décolleté and body;
- Treatment of tissue laxities to correct body and face imperfections such as skin hypertrophy and hypotonia, or reduced skin elasticity;
- Reshaping hereditary body and face conditions such as skin hypertrophy or hypotonia, or reduced skin elasticity;
- Restoring and/or correcting signs of facial fat loss (Spondylosis) in people with human immunodeficiency virus;
- Restoring and/or correcting facial acne scars.

8. Contraindications

- Do not inject the product in the periorbital area (eyelids, crow's feet) or in the glabellar area (forehead).
- Do not inject into blood vessels (intravascular use).

Do not use:

- in patients with vascular collagen diseases or autoimmune connective tissue diseases (scleroderma);
- in patients with epilepsy not controlled by medication;
- in patients prone to developing hypertrophic scars;
- in patients with known hypersensitivity to collagen;
- in patients with porphyria;
- during pregnancy or breast-feeding;
- in children;
- on skin areas affected by inflammation and/or infection (acne, herpes, etc.);
- in immediate association with a laser treatment, deep chemical peel or dermabrasion. In case of superficial peel, injection is not recommended if a significant inflammatory reaction is present.
- coagulation disorders

9. Possible side effects

The product is hypoallergenic. No side effects associated with the product have been reported.

Patients should be informed of potential side effects secondary to administration of this medical device by means of injection; these might appear immediately or some time after administration.

They include (but are not limited to):

- slight stinging sensation, which disappears in a few minutes;
- inflammatory reactions (such as redness, swelling, edema, erythema, etc.), sometimes in conjunction with itching, or pain to touch, or both. These reactions can last for a week;
- hematoma;
- hardening or nodules at the injection site;
- coloration or discoloration of the treated area;
- cases of necrosis of the glabellar region, abrasion, granuloma, and immediate or delayed hypersensitivity after collagen injections have been described in the literature. It is therefore imperative to take these potential risks into account;
- inflammatory reactions lasting longer than a week or other side effects must be immediately reported by the patients to their practitioner, who will provide adequate treatment to eliminate such reactions;

Any other undesirable side effect associated with injection of the NITHYA must be reported to the distributor and / or manufacturer.

10. Mode of application

The product should only be used by medical personnel.

For use of NITHYA sterile intradermal implant based on equine collagen Type I, model NITHYA implant 3x70 mg dilute the collagen powder contained in the vial in 5 ml of saline solution or double distilled water in order to obtain a suspension. For use of NITHYA sterile intradermal implant based on equine collagen Type I, model NITHYA implant 206 mg vials are to be diluted in 10 ml saline solution or double distilled water (w.f.) in order to obtain a suspension.

The suspension should be prepared immediately before injection and is intended for single-use; any remaining solution must be discarded. The product should be injected with a fine needle (30-32G) into the deep dermis or dermo-epidermal junction. To ensure perfect sealing and avoid product leakage during administration, make sure that the needle is securely fastened to the syringe. Evenly distribute the product in the area to be treated and then apply a light massage. Effect of 1 procedure lasts not less than 1 month.

The treatment should be performed 3 times at 15-days intervals. The course of treatment can be repeated at least 3 months after the last application. Do not overcorrect. The amount of suspension to be administered, as well as the frequency of administration and technique are at the discretion of the doctor, depending on the area to be treated.

Operating temperature range of MD is from +32°C to +2°C

11. Precautions of use

• Before starting the treatment, the patient should be informed of the device's indications and contraindications, as well as of its incompatibilities and possible adverse effects. As with any other intracutaneous treatment, the product may pose a

• Before injecting, thoroughly disinfect the area to be treated. The usual precautions for intradermal infiltrations should be observed.

• Inject slowly. Failure to comply with these precautions may lead to needle detachment and/or product leakage in the later back area.

• The amount of injected product will depend on the area to be treated.

• After injection, it is important to massage the treated area in evenly spread the product.

Do not use in case of known individual hypersensitivity to the product.

• Do not inject in skin areas with active infection or inflammation.

• Injection of the product in the area below the eyes should only be performed by practitioners specifically trained in this technique and with solid knowledge of the physiology of this area.

• The product is not recommended during pregnancy and breast-feeding, although specific contraindications are not known.

• The physician's technical ability is essential for the treatment's success, therefore this medical device should only be administered by physicians who have received specific training in the technique of filler injections.

• Altering the product or using it in any way that is not in accordance with the conditions of use described in this document may affect its sterility and performance, which will not be guaranteed any longer.

• Make sure that the package and the glass vial are intact before use. The product is for single use; do not reuse any remaining solution after its first opening. Quality and sterility can be guaranteed only for vials/packages sealed at the origin.

• Dispose of this product responsibly after use.

• Check the product's expiry date and do not use it after it has expired.

• Keep this medicine out of the sight and reach of children

12. Methods of cleaning

Before injection, proceed to thoroughly disinfection the area to be treated.

13. Sterility

Gamma irradiation.

14. Storage and transportation conditions and Terms of use

To store and transport NITHYA at a temperature below 25 °C, in dry environment and in its original packaging. Protect from sunlight, heat and handle with care. Humidity not more 60%.

No particular caution is required during transportation. It is carried out by all types of covered transport. Products are placed in cardboard boxes.

List of applicable standards

The medical device " NITHYA sterile intradermal implant based on equine collagen Type I obtained the EC Certificate in compliance with the following national and international applicable standards:

- 93/42/EC Directive concerning the medical device, recognised in Italy by D.Lgs n. 46 of 24.02.1997;

2007/47/EC Directive 0405.05.2007 amending the 93/42/EC Directive, recognised in Italy by D.Lgs n. 27 of 25.01.2010;

- UNI EN ISO 14971/2013 - Application of risk management to medical devices;
- ISO 9001/2015 - Quality risk management;
- UNI EN ISO 22442/2009 - Medical devices utilizing animal tissues and their derivatives;
- UNI EN ISO 1841/2009 - Information supplied by the manufacturer of medical devices;
- UNI EN ISO 15223-1/2012 - Symbols to be used with medical device labels, labelling and information to be supplied;
- UNI EN ISO 13485/2016 - Quality management systems - Requirements for regulatory purposes;
- UNI EN ISO 10993-5:2009 and UNI EN ISO 10993-10:2010 - Biological evaluation of medical devices - Biocompatibility tests;
- UNI EN ISO 11137-1, 2, 3 - Sterilization of Medical Devices.

15. Shelf Life

5 years since the date of manufacture.

16. Labelling of sterile consumer packaging

REF/PA	
End-use device	
Type I collagen	
Opromed S.p.A. / Corso Venezia 5 20121 Milan (Italy) / Tel: 39 02 8855698	Manufacturer
☒	Single use
☒	Do not use if packaging is damaged
☒	Reading information contained within the box
☒	Storage temperature: below 25 °C
STERILE	Sterile Product: sterilized by irradiation
CE mark	CE Mark
LOT	Lot number
EXP	Shelf life (month/year)

Vital Information:	
RTMHA	-
Collagen powder	-
70 mg / 300 mg	-
EURORESEARCH S.r.l. Corso Venezia 5 20121 Milano (Italy)	-
CE 0870	CE Mark
lot.	-
Exp.	Expiration date (month/year)

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17. Disposal of medical device
Products that have expired and packaging materials should be destroyed as medical waste of Class A (epidemiologically non hazardous waste, class in composition to the municipal solid waste) in accordance with standard practice in the medical institutions, and in the manner prescribed by local, federal and state rules and regulations.

18. Warranty
The manufacturer, Euroresearch S.r.l. (from Italy), ensures the compliance of medical device with the requirements of the Technical File, if it follows conditions of use, storage and transportation. The guaranteed shelf life is 5 years from the date of manufacturing.
The product in the manufacturer's packaging should be stored in the warehouses of the supplier at a temperature below 25°C in a dry place protected from light, humidity no more than 60%. Manufacturer assumes absolute responsibility in case of defective goods and replaces them.

19. Manufacturer and manufacturing site
Manufacturer: Euroresearch S.r.l., Corso Venezia 5, 20121 Milano, Italy
Manufacturing site: HighTech Contract Manufacturing Medical, Nieuwenhof 30-413, Nijmegen, 6546 BA, The Netherlands

20. Authorized representative in Russia
LLC Heavy Polite - Russia, 125280, Moscow, Novy Arbat 34, bld 1, placement 1, floor 4, office 7. Tel. +7 495 918 77 88; info@easmed.com

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