Instructions for use PD-Paed Plus Paediatric PD system

GB

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RUS

**FRESENIUS
MEDICAL CARE**

Fresenius Medical Care AG & Co. KGaA
61346 Bad Homburg, Germany
Phone: +49 6172-609-0

-- for registration purposes --

Bad Homburg, 05.02.2018
Place and date

**FRESENIUS
MEDICAL CARE**
Fresenius Medical Care AG & Co. KGaA
61346 Bad Homburg - Germany

i. A. Dorothea Hoffmann
Associate Medical Device Regulatory

Fresenius Medical Care AG & Co. KGaA
61346 Bad Homburg - Germany

i. V. Annette Manck
Senior Expert Medical Device Regulatory Affairs

Fresenius Medical Care AG & Co. KGaA
61346 Bad Homburg - Germany

Instructions for use

PD-Paed Plus

Paediatric PD system

GB

General Information

Refer to product or box label for reference number of PD-Paed Plus.

Indications

The PD-Paed Plus is a system to perform peritoneal dialysis for the treatment of chronic and acute renal failure in premature babies and infants up to 10 kg of body weight and/or a peritoneal fill volume up to 200 ml. It is intended for clinical use only. The PD-Paed Plus is intended to be used only with stay•safe® bicaVera®, stay•safe® balance or stay•safe® CAPD/DPCA peritoneal dialysis fluid from Fresenius Medical Care.

Contraindications

Special contraindications are unknown. In general, contraindications for peritoneal dialysis are applicable.

Side effects

Specific side effects are not known. In general, side effects for peritoneal dialysis are applicable.

Potential adverse reactions of the treatment mode

Infections and infestations

very common: peritonitis indicated by a cloudy effluent. Later abdominal pain, fever, and malaise (generally feeling unwell) may develop, in very rare cases, generalised blood poisoning (sepsis). The patient should seek medical advice immediately. The bag with the cloudy effluent should be closed with a sterile cap and assessed for microbiological contamination and white blood cell count.

Skin exit site and tunnel infections indicated by redness, oedema, exudations, crusts and pain at the catheter exit site. In case of skin exit site and tunnel infections the attending physician should be consulted as soon as possible.

Respiratory, thoracic and mediastinal disorders not known: dyspnoea caused by the elevated diaphragm.

Gastrointestinal disorders very common: hernia,

common: abdominal distension and sensation of fullness uncommon: diarrhoea, constipation.

Injury, poisoning and procedural complications.

common: in- and outflow disturbances of the dialysis solution, shoulder pain.

Materials

Tubes: medical grade soft polyvinylchloride (PVC).

Connectors and other components: polyvinylchloride (PVC), polycarbonate (PC), acrylonitrile butadiene styrene (ABS), polypropylene (PP), polybutylenterephthalate (PBT), polytetrafluorethylene (PTFE), phenylalanine (PHE), silicon, synthetic rubber (polyisoprene), high density polyethylene (HDPE), acrylic, polystyrene (PS). The PD-Paed Plus is DEHP-free and does not contain any natural latex.

Directions for use

The PD-Paed Plus is a closed system with a combined inflow ① and outflow

② part and a drainage bag ③.

Cautions

- The PD-Paed Plus is intended for maximal 24 hours use and for single use only. Usage for periods longer than 24 hours and reuse may be hazardous to the patient. Cleansing solutions and disinfectants may damage materials employed for the PD-Paed Plus. Safety and performance of use can no longer be guaranteed and the manufacturer assumes no liability.
- The PD-Paed Plus must only be used in conjunction with the peritoneal dialysis fluids listed above by healthcare professionals with the appropriate training, knowledge and experience and for whom proof of instruction can be shown. In addition, refer to the corresponding training manual.
- Misuse or improper handling may lead to severe complications!
- Do not use if sterile package is damaged or caps are not in place.
- Unpack and handle system according to aseptic technique.
- Ensure that all caps and connections are secure. Inspect the peritoneal dialysis system and the accessories for leaks during the priming phase and treatment, taking corrective measures (e.g. tightening Luer-Lock/stay•safe® connections) or exchanging the system if necessary.
- Do not use after expiry date (see label).
- Do not use needles with a diameter larger than 20 gauge to puncture injection ports.
- The plastics used may be incompatible with drugs or disinfectants. Please check compatibility before use.

Application

Assemble the PD-Paed Plus as follows:

- Attach drainage bag to yellow connector of drainage system at outflow part.

Install the PD-Paed Plus on the infusion pole:

organizer.

- Place inflow burette in upper burette holder. Place outflow burette in lower burette holder.
- Connect the peritoneal dialysis fluid bags to the PD-Paed Plus. For the handling of the PD fluid bags please refer to the respective product insertion leaflets.
- Avoid kinking of tubes.

Warranty

The manufacturer warrants that the PD-Paed Plus has been manufactured in accordance with its specifications and in compliance with the GMP regulatory guidelines.

Products with manufacturing defects will be replaced if the defect is reported with the lot number.

The manufacturer shall not be liable for any misuse, improper handling, non compliance with instructions for use and cautionary notes and for any damage incurred subsequent to the manufacturer's delivery of the PD-Paed Plus.

Expiry date

The expiry date of the medical device is three years from the manufacturing date.

Storage condition

Store in the temperature range from 5 to 30 degrees C avoiding water, direct sunlight, and heat and moisture. Atmospheric pressure, low humidity (normal environmental humidity).

Transportation

The packed devices can be transported by all types of covered vehicles according to the local laws of transportation of medical devices and according to the indications stated on the product label.

Operating conditions

Since the product PD-Paed Plus produced by Fresenius Medical Care is intended to be used by the clinical personnel (physicians and nurses), the environmental conditions and requirements are strictly stated by clinical specialists according to hospital specific instructions.

Disposal

The hazard class of medical waste disposal (device not used) on the basis of morphological characteristics of the composition is CER 07 0 13. After use the product must be disposed according to local rules for infective waste.

The sizes of the product components

#	Description	Length (mm)	Width (mm)	Diameter (mm)
1.	Drainage bag 5lt	387	300	NA
2.	Burette 230 ml	177	48	49,6
3.	Measure chamber	176	48	49,6
4.	Drip chamber with insert for valve	80	17	17,8
5.	Patient connector PC assembled recessed cone	66,5	30	13 10,3
6.	Housing in roller for tube 3.1 x 4.3	55	11	NA
	Roller	NA	NA	13,85 10,55
7.	Stay Safe connector PD Paed	34,3	10.4	10.4 8.8 7.84 6.7 6.4 4
8.	PVC Tube 4,3x6,8	40	NA	4,3 6,8
		100		
		1500		
		80		
		300		
		300		
		150		
		50		
		1300		
9.	PVC Tube 3,5x5,5	50	NA	3,5
		50		5,5
10.	PVC Tube 3,1x4,3	100	NA	3,1 4,3

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No. 149 of the deed register for 2018

Upon the notary's question as to any prior involvement within the meaning of Sec. 3 Para. 1 Sent. 1 No. 7 German Notarial Recording Act (Beurkundungsgesetz – BeurkG) the below-named persons denied any such prior involvement.

I hereby certify, that the above are the true signatures, respectively acknowledged in my presence, of

Mrs. Dorothea Hoffmann, born 4 March 1966,
and Mrs. Annette Manck, born 18 April 1973
business address of both: Else-Kröner-Straße 1, 61352 Bad Homburg v.d.Höhe,
Germany
- both personally known to me -

acting on behalf of the partnership limited by shares

Fresenius Medical Care AG & Co. KGaA,
Else-Kröner-Straße 1, 61352 Bad Homburg v.d.Höhe, Germany.

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