
Certification for Information in Dossier

The submitted dossier for the Russian application contains the following:

- Instructions for use for Advantage hip prostheses:

Product	IFU #
Advantage Reload Cementless Acetabular Cups	F192
Advantage 3P Cementless Acetabular Cups	F192
Advantage 3P Plasma TiHA Cementless Acetabular Cups	F192
Advantage Cemented Acetabular Cups	F192
Advantage Revision Cementless Acetabular Cups	F192
Advantage Inlays	F192
Advantage E1 Inserts	F192
Advantage Screws	F126
Advantage Pin	F192

I certify that the documents which are submitted are true and accurate to the best of my knowledge. In addition, Biomet France takes responsibility for any legal obligations related to this submission.

Yannick LE GUAY
Regulatory Affairs Manager
Biomet France
November 2, 2017

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1. Product name and description

This box contains a medical device from the AVANTAGE ACETABULUM range intended to replace the natural bone acetabulum during total hip arthroplasty.

This range is made up of the following components:

- AVANTAGE RELOAD Cup P0460Pxx
- AVANTAGE CEMENTED Cup P04630xx
- AVANTAGE 3P Cup P04610xx
- AVANTAGE REVISION Cup P0464Hxx
- AVANTAGE 3P Ti+HA Cup P0461PXX
- AVANTAGE Insert P05600xx / P05610xx
- AVANTAGE E1 Insert P0560Exx / P0561Exx
- AVANTAGE Pin P0461070

2. Materials used for the implant

Name	Material	Standards
AVANTAGE RELOAD Cup	Stainless steel	ISO 5832-9: 07
AVANTAGE 3P Ti+HA Cup	Hydroxyapatite Ceramic	ASTM F1185-03
AVANTAGE CEMENTED Cup	Rough pure titanium	ASTM 1580-07
AVANTAGE 3P Cup	Stainless steel	ISO 5832-9: 07
AVANTAGE REVISION Cup	Hydroxyapatite Ceramic	ASTM F1185-03
AVANTAGE Insert	ARCOM UHMWP	ASTM F648-10
AVANTAGE E1 Insert	E1 HXLPE (highly cross-linked polyethylene enhanced with vitamin E)	ASTM F648-10
AVANTAGE Pin	Stainless steel	ISO 5832-1: 08

3. Choice of implants

The use of radiographic templates during preoperative planning and instrumentation specific to the surgical technique will aid in the choice of the most appropriate sizes for the patient's morphology.

4. Indications

AVANTAGE range acetabula can be used for the following indications of total hip arthroplasty:

- Primitive or secondary coxarthrosis;
- Inflammation of the hip; rheumatoid arthritis, etc.;

- Sequelae from previous hip operations, osteotomies, etc.;
- Congenital hip dysplasia;
- Femoral neck fracture;
- Femoral head necrosis.
- Revision of the bearing with a RECAP/MAGNUM cup or an AVANTAGE cup (see below).

Revision of the bearing – AVANTAGE cup:

Only AVANTAGE E1 P0560Exx and P0561Eyy (xx=44 to 66, yy=50 to 66) inserts may be used for a revision of the bearing (revision of the head and insert) when an AVANTAGE cup has already been fitted.

Revision of the bearing – RECAP/MAGNUM cup:

The AVANTAGE E1 P0560Ezz and P0561Ezz (zz=44 to 66) inserts may be used for revisions of the bearing when a RECAP/MAGNUM cup has already been fitted following the failure of the metal/metal association.

Note : AVANTAGE E1 inserts – sizes 44/46/48 Ø 28 mm - are only intended to be used with a RECAP/MAGNUM cup that has already been fitted and not for an AVANTAGE cup revision.

AVANTAGE range acetabula are double mobility and are therefore particularly recommended for patients presenting the following profiles:

- Patients of any age with a high risk of dislocation.
- Pathology exposing a patient to a risk of dislocation:
 - Neck fracture;
 - Tumor surgery;
 - Neuromuscular disorder (Parkinson's, buttock muscle paralysis, etc.).
- Hip prosthesis revision especially for recurring dislocations.

AVANTAGE REVISION range cup can be used for the following indications of total hip arthroplasty revision:

- Aseptic loosening.
- Septic loosening in connection with a replacement in one or two steps depending on the surgeon's therapeutic choice for treating the infection.
- Acetabular bone defect with or without damage to the bony ring and/or destruction of the area (Paprosky Classification: 2B – 2C – 3A).
- Periprosthetic fracture of the coxo-femoral joint.
- Unexplained persistent pain.
- Recurring dislocation.

The choice between the different versions of AVANTAGE ACETABULUM depends on the acetabulum bone defect

which gives the following grades using the Paprosky classification:

AVANTAGE ACETABULUM	Paprosky Classification
AVANTAGE RELOAD Cup + AVANTAGE or AVANTAGE E1 Insert	1 – 2A
AVANTAGE 3P Cup and 3P Ti+HA + AVANTAGE or AVANTAGE E1 Insert + AVANTAGE Pins	1 – 2A – 2C
CEMENTED AVANTAGE Cup + AVANTAGE or AVANTAGE E1 Insert	1 – 2A
AVANTAGE REVISION + AVANTAGE or AVANTAGE E1 insert	2B – 2C – 3A

5. Expected performance

- Restoration of joint function and pain reduction
- Osteo-integration for hydroxyapatite-coated medical devices.

6. Contraindications

- Local or systemic infections.
- Severe muscular, neurological or vascular deficiencies of the extremity involved.
- Bone destruction or poor bone quality that is likely to affect implant stability (Paget's disease, osteoporosis).
- Any concomitant disease that is likely to affect implant function.
- Allergy to any of the components of the implant.
- Patients incapable of following surgeon's recommendations.
- Patients weighing over 110 kg for the AVANTAGE cup / AVANTAGE insert association.
- Patients weighing over 100 kg for the RECAP/MAGNUM cup/ AVANTAGE E1 insert association.
- The AVANTAGE Insert and AVANTAGE E1 Insert are not cleared for use with the Magnum Tri-Spike Acetabular Cup.

7. Factors likely to compromise implant success and device performance

- Significant deformations, congenital dislocation.
- Local bone tumors.
- Ligament instability or serious and untreatable muscle contractures.
- Patient history of infections or repeated falls.

CE 0459

F192

Notice updated: 11/2014

CE marking: 1998

- Substance dependence.
- Poor positioning of insert in cup.
- Trauma possibly causing an implant to break
- Poor orientation of acetabulum.

8. Adverse reactions

- Wearing of the acetabular components' polyethylene joint surfaces has been observed following total hip replacement. Higher wear rates may be caused by particles of cement, metal or other debris that may lead to abrasion of joint surfaces. Significant wear may shorten the prosthesis' useful life and lead to premature replacement of the prosthesis' worn components.

For insert revisions (the cup remaining in place), an inspection of the cup (no deep scratches, etc.) and thorough cleaning are required to remove all debris (bone, metal, tissue, cement, etc.) which could lead to excessive and/or premature wearing of the insert. The surgeon will evaluate whether the cup is showing signs of wear requiring its complete removal or the replacement of the insert alone.

For insert revisions, it is recommended that the prosthetic femoral head should be revised to prevent any excessive and/or premature wearing of the components.

- In all cases of joint replacement, asymptomatic localized progressive bone resorption (osteolysis) may be noted around prosthetic components as a result of foreign body reactions triggered by particles. These particles are generated by the interaction between the various components, as well as between the components and bone, mainly through mechanisms of wear, adhesion and fatigue. Other particles may also be produced by the wearing of another body. Osteolysis may lead to successive complications requiring the removal and replacement of prosthetic components.
- Although rare, cases of metal intolerance following joint replacement have been observed. Implantation of foreign material in tissues may result in histological reactions involving the formation of macrophages and fibroblasts.
- Dislocation or subluxation of prosthetic components due to improper positioning and/or migration of components can occur. Muscle and fibrous tissue laxity can also contribute to these conditions. An unsuitable choice, incorrect positioning, poor alignment or incorrect fixation of the cup may result in excessive stresses on the insert, which can reduce the life expectancy of the prosthesis.

- Prosthetic components can come loose or migrate following trauma.
- When used in association with the RECAP/MAGNUM cup in patients weighing over 100 kg, a blocking of the mobility of the insert in the cup (large joint) may occur.

9. Shelf life and sterility

Implants are supplied sterile and packaged individually in double wrapping. Sterilization is carried out by exposure to a minimum dose of 25 kGy of gamma irradiation. The expiration date is printed on the label. Do not use implants after the expiration date.

SINGLE USE - DO NOT RESTERILIZE IMPLANTS

The performance of the device cannot be guaranteed if implants are reused after a re-sterilisation process (as plastic parts can become distorted or change structurally during the sterilisation process, and surfaces can become damaged; moreover, the efficacy of the decontamination method and integrity of the coating cannot be guaranteed).

10. Specific instructions for use

- Packaging must not show signs that could indicate a defect in the sterility and/or integrity of the medical device.
- Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants that have been, even momentarily, placed in a different patient. Furthermore, reusing an implant, could cause patient contamination
- When handling implants in the operating room, all necessary precautions must be taken to avoid damaging the implants (scratches, etc.).
- The implantation must be performed by an orthopedic surgeon experienced in lower limb surgery, according to the manufacturer's recommended surgical technique.
- For bearing revisions, following the removal of the prosthetic head, check for the absence of damage to the neck and/or cone of the femoral stem opposite.
- Do not use inserts made from ArCom Polyethylene (P05610xx and P05600xx) for bearing revisions already in place.
- Use the instrumentation designed by the manufacturer for inserting this implant.

- Do not use components of this prosthetic system with components from another manufacturer other than the one indicate in § 11
- Comply with size relationships between various related prosthetic parts.
- The surgeon must consider the design and size of the related stem to prevent a cam effect.
- Never allow antibiotics or acid local antiseptics to come into contact with the hydroxyapatite coating.
- Do not cement a hydroxyapatite-coated cup.
- The possible effects of a magnetic resonance imaging environment on this device have not been determined. This product has not been subjected to heating or migration tests in a magnetic resonance imaging environment.
- Dirty implants should be treated as biological waste
- Comply with size relationships between various related prosthetic parts: insert and cup sizes must be identical.

11. Compatible medical devices

- Specific instrumentation:

Instrumentation name	Reference
AVANTAGE UPGRADE INSTRUMENT KIT	110026613
AVANTAGE instrumentation	A0900046
Additional AVANTAGE Revision instrumentation	A0900067
MINI-INVASIVE instrumentation for acetabulum	A0900164
MINI-INVASIVE instrumentation for AVANTAGE acetabulum	A0900165
Instrumentation for removing heads and stems	31-601498

- Femoral stems from the Biomet range: AURA/ARPEGE, AURA RR, AURA II, ART, DB-10, EXCEPTION, UP-TION, PERTROCH+, MAINSTREAM MÜLLER, GTS, TAPERLOC, BIMETRIC, STANMORE, F40, PPF, MALLORY HEAD, CONELOCK, Oxford Revision, Reach.
- RECAP/MAGNUM acetabula: 1578xx / 1579xx / 1308xxxx in association with an AVANTAGE E1 insert.
- Biomet Femoral heads, 22.2 and 28mm diameter, made of stainless steel, alumina ceramic or cobalt chromium with a 12/14 taper and alumina ceramics (BioloX Delta) or cobalt-chromium with a Type 1 taper
- A CEMENTED AVANTAGE cup can be used with an AVANTAGE acetabular reinforcement cross in cases of revision of the acetabular part of a total hip prosthesis:

AVANTAGE acetabular reinforcement cross
P0464xxx / P0465xxx

CE 0459

F192
Notice updated: 11/2014
CE marking: 1998

- AVANTAGE screws:
AVANTAGE Ø 4.5mm screws should be used with AVANTAGE 3P & REVISION cups and 3P Ti+HA cups:

Name	Reference
AVANTAGE Ø 4.5mm screws	P06060xx

- Radiographic templates:
AVANTAGE RELOAD Acetabulum templates: ref. CA057
AVANTAGE 3P Acetabulum templates: ref. CA056
AVANTAGE CEMENTED Acetabulum templates: ref. CA080 and ref. CA086.
AVANTAGE REVISION template : ref. CA079.
AVANTAGE 3P Ti+HA Acetabulum templates :

Preoperative template designation	Reference
AVANTAGE 3P Plasma Ti HA cup template 100%	XRAY000356
AVANTAGE 3P Plasma Ti HA cup template 110%	XRAY000357
AVANTAGE 3P Plasma Ti HA cup template 115%	XRAY000358
AVANTAGE 3P Plasma Ti HA cup template 120%	XRAY000359

- PMMA-type surgical bone cement used with AVANTAGE CEMENTED cups.
- REGENEREX cages:
The AVANTAGE CEMENTED cup can be used with a REGENEREX cage, when revising the acetabular component of a total hip prosthesis.
The table below shows the corresponding size in the REGENEREX range for each size of AVANTAGE cemented cup:

REGENEREX CAGE		AVANTAGE CEMENTED CUP	
Size	Reference	Size	Reference
T60	Pt-1250-60	T44	P0463044
T62	Pt-1252-62	T46	P0463046
T64	Pt-1254-64	T48	P0463048
T66	Pt-1256-66	T50	P0463050
T68	Pt-1258-68	T52	P0463052
T70	Pt-1260-70	T54	P0463054
T72	Pt-1262-72	T56	P0463056

T74	Pt-1264-74	T58	P0463058
T76	Pt-1266-76	T60	P0463060

12. Storage

Products should be stored in their original packaging.

13. Additional instructions for the surgeon

The surgeon should be aware that the anatomic, chemical and physical characteristics of living tissue limit the possibilities in the form and choice of material. Consequently, it should be noted that, like bone that no longer carries out its functions, prosthetics have their limits.

The proper implantation of the hip prosthesis is essential for successful arthroplasty and long-lasting results. This is why the detailed recommendations for surgical technique must be followed.

14. Patient information

Patients should be informed that hip arthroplasty is a major operation, and their hopes should not be raised too high. Patients should be warned not to overwork their prosthesis by performing unreasonable activities that could cause shocks, and to monitor their weight.

All trademarks herein are the property of Biomet, Inc. or one of its subsidiaries unless otherwise indicated.
Biologix trademark is the property of CERAMTEC.

Manufacturer:
Biomet France
Plateau de Lautagne
26000 Valence
Tel.: +33 (0) 4 75 75 91 00
Fax: +33 (0) 4 75 75 91 01

Warning, see instruction leaflet

Do not reuse

Batch code

Catalogue reference

Use by

Sterilized by irradiation

CE 0459

F192
Notice updated: 11/2014
CE marking: 1998

CANCELLOUS BONE FIXING SCREW AVANTAGE SCREW

1. Description of the product

This box contains a BIOMET France fixing screw.
BIOMET France fixing screws are available as follows :

- The AVANTAGE screw: P06060xx
- The cancellous bone fixing screw P06011xx.

This product is for single use only.

2. Materials used for the implant

Item	Material	Standard
AVANTAGE screw	Stainless steel	ISO 5832-1 2007
Cancellous bone fixing screw	Titanium alloy	ISO 5832-3 1996

3. Choice of implants

Radiographic templates for pre-operative planning together with the ancillary for the operating technique make it possible to select the most appropriate sizes for the morphology of the patient.

4. Indications

BIOMET France fixing screws are designed to be implanted in combination with the implants described in Point 11.

The AVANTAGE cortical bone screw is used for fixing the acetabular cups, armatures and support rings for the hip and for fixing the osteosynthesis plates.

The cancellous bone fixing screw is designed to complete the anchorage of the acetabular cup for the hip or the tibial implant for the knee.

5. Expected performance

For the performance of this medical device, refer to the performance of the implant with which it will be combined.

6. Contraindications

- Local or systemic infections.
- Serious muscular, neurological or vascular deficiencies of the limb concerned.
- Bone destruction or poor quality bone likely to affect the stability of the implant (Paget's disease, osteoporosis, etc.)
- Any concomitant disorder likely to affect the function of the implant.

- Allergy to one of the implant components

7. Factors likely to compromise the success of implantation and the performance of the device

- Significant deformities, congenital dislocation.
- Local bone tumours.
- Serious, untreatable ligament instability or muscle contracture.
- Anamnesis indicating repeated infections and falls.
- Drug and alcohol dependency

8. Undesirable effects

- In all cases of joint replacement, it is possible to note progressive, asymptomatic, localised bone resorption (osteolysis) around the prosthetic components following the reaction to a foreign body caused by particles. These are generated by the interaction between the various components as well as between the components and the bone and are mainly due to the mechanisms of wear, adhesion and fatigue. Other particles can also be produced by wear in another body. Osteolysis can cause successive complications requiring the removal and replacement of components of the prosthesis.
- Even though these are rare, cases of intolerance to the metal following joint replacement have been observed. The implantation of foreign materials in the tissues may cause histological reactions involving the formation of macrophages and fibroblasts.
- Dislocation or sub-luxation of the prosthetic components may result in poor positioning and/or displacement of the components. Relaxation of the muscles and fibrous tissues may also contribute to these phenomena.
- The prosthetic components may loosen or move following trauma.

9. Shelf life and sterility

The implants are delivered sterile and packed as single units in double packaging. Sterilisation is performed using irradiation (gamma ray) at a minimal dose of 25 kGy. The expiry date is printed on the label. The implants must not be used after the date shown.

DO NOT RE-STERILISE THE IMPLANTS

Use of the device would no longer guarantee the achievement of performance if it has been re-sterilised (distortion or structural change in the plastic parts during the sterilisation, damage to the surface condition,

questions concerning the efficacy of the decontamination method and questions concerning the integrity of the coating).

10. Specific instructions for use

- The packaging must not show traces that would call into question the sterility and/or integrity of the medical device.
- Never use damaged or explanted implants. If reused, the implant's long-term performances in terms of restoring functions and mobility cannot be guaranteed. In addition, the use of an explant could infect the patient.
- Every precaution must be taken when handling implants in the operating theatre to avoid damaging them (scratches, etc.).
- The implantation must be performed by an orthopaedic surgeon experienced in lower limb surgery.
- Use the special ancillary designed by the manufacturer for insertion of the associated implant.
- Do not use items from this prosthetic system with items from other systems and/or other manufacturers, other than those listed in Point 11.
- The possible effects of a magnetic resonance imaging environment on this device have not been determined. This product has not been subjected to heating or migration tests in a magnetic resonance imaging environment.
- Any soiled implant must be treated as biological waste.

11. Compatible medical devices

The medical devices associated with the AVANTAGE screws are Biomet France devices compatible with 4.5 mm diameter screws such as :

- Biomet France acetabular cups and the associated ancillary
- Biomet France osteosynthesis plates and the associated ancillary
- Biomet France armatures and the associated ancillary
- Biomet France support rings and the associated ancillary

The medical devices associated with the cancellous bone fixing screws are Biomet France devices compatible with 6.5 mm diameter screws such as :

- Biomet France acetabular cups and the associated ancillary
- Biomet France tibial implant and the associated ancillary

CE 0459

F126
Update of the instructions: 08/2015
First put on market: 1992

The products must be stored in their original packaging.

13. Additional information for surgeons

The surgeon must be aware that the particular anatomical, chemical and physical characteristics of living tissue limit the possibilities in terms of shape and the choice of materials. Consequently it must be specified that, just as bones cease to fulfil their function, prostheses have their limits.

Correct implantation of the hip prosthesis is essential as the success of the arthroplasty and the longevity of the results depend upon it. For this reason, it is important to comply with the recommendations detailed in the operating technique.

It is essential to perform regular radiographic examinations to prevent subsequent serious complications.

14. Patient information

All patients should be informed that hip arthroplasty remains a major operation, for which it would be stupid to give them false hope. Patients should be warned not to overtax their prosthesis by undertaking unreasonable activities that could cause knocks and to monitor their weight.

All trademarks herein are the property of Biomet, Inc. or its subsidiaries unless otherwise indicated.

 Caution, see the instructions for use

 Do not re-use

 Batch code

REF Catalogue reference

 Use by

 Sterilisation using irradiation

Manufacturer :
Biomet France
Plaisau de Lautagne
26000 Valence
France
Tel.: +33 (0) 4 75 75 91 00
Fax: +33 (0) 4 75 75 91 01

F126
Update of the instructions: 08/2015
First put on market: 1992

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