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Mr. **Matthias BÜRGER**, Nationality: Germany, in Winterthur, Switzerland,
personally known to us,

who is entered in the Register of Commerce of the Kanton of Zürich as person with the right to sign
per procuration individually for the

Zimmer GmbH, company with limited liability with registered head office in Winterthur,
Switzerland.

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The inspection of the commercial register has taken place directly before the official certification by
internet inquiry.

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Certification Information in Dossier

The submitted dossier of the Russia application **NCB implants for osteosynthesis - non-sterile** contains the following:

- Indication for use № D011500204
- Indication for use № D011500201

I certify that the documents which are submitted are true and accurate to the best of my knowledge. In addition, Zimmer GmbH takes responsibility for any legal obligations related to this submission.

24 July 2015 

Matthias Bürger
Quality Assurance & Regulatory
Compliance VP
Zimmer GmbH

Official Certification

Seen for authentication of the foregoing signature of

Mr. **Matthias BÜRGER**, born 21st July 2015, Nationality: Germany, residing at Maierriedweg 16, 8408 Winterthur, Switzerland, personally known to us.

The signature was acknowledged before us by an authorized third party, Mr. J. Fetibegovic.

Winterthur, 24th July 2015

BK no. 1856

Fee CHF 20.00

NOTARIAT OBERWINTERTHUR-WINTERTHUR

Martin Schmuki, Notary Public

Barcode	D011500201	Ed. 2014/10	
Non-sterile Metal Internal-Fixation Implants			
 Zimmer GmbH Sulzerallee 8 8404 Winterthur, Switzerland Telephone +41/ (0)52 262 60 70 Fax +41/ (0)52 262 01 39 www.zimmer.com		Representative in the USA: Zimmer, Inc. 1800 West Center Street Warsaw, Indiana, 46580, USA	
   			

ENGLISH

The CE mark is valid only if it is also printed on the product label.

Important information for the operating surgeon

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications that may arise from use of the device in circumstances outside of Zimmer's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.

DESCRIPTION

Zimmer implants for internal fixation are available in stainless steel and/or titanium versions, depending on the product.

- Stainless steel implants are made from *Protasul*[®]-S (XCrNiMo – ISO 5832-1 / ASTM F 138/139) or *Protasul*-S30 (XCrNiMnMoNbN – ISO 5832-9 / ASTM F 1586).
- For titanium implants, pure titanium (ISO 5832-2 / ASTM F67), titanium alloy *Protasul*-64 (Ti6Al4V - ISO 5832-3 / ASTM F136) or titanium alloy *Protasul*-100 (Ti6Al7Nb - ISO 5832-11 / ASTM F1295-92) are used.

For detailed information refer to the corresponding product information and/or surgical technique.

INDICATIONS

The indications are product dependent and are listed in the corresponding description of the surgical technique and/or in the product information.

CONTRAINDICATIONS

- All concomitant diseases that may impair the fixation of the implant and/or the success of the intervention.
- Lack of bone substance or poor bone quality impairing stable fixation of the implant.
- Severe muscular, neural or vascular diseases that endanger the success of the intervention.
- Allergy to the implanted material.
- Acute or chronic, local or systemic infections.
- Further product dependent contraindications are mentioned in the corresponding description of the surgical technique and/or in the product information.

⚠ WARNINGS

- Do not reuse. This device is single use only. Reuse may adversely affect device performance.
- Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- Do not use any component if damage is found or caused during setup or insertion.
- Implants and implant parts must only be combined with components belonging to the same system. No liability is accepted for products of third parties.
- Use only instruments specifically designed for use with these devices to help ensure accurate surgical implantation.
- The load-bearing capacity of the implant can be compromised by notching, scratching, or striking.
- Do not mate titanium alloy components with stainless steel. (only applicable for the USA)
- Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions reducing the service life of the implants.
- Do not use these products for other than labeled indications (off-label use).

PRECAUTIONS

- Continued surveillance for new or recurrent infection should be continued as long as the device is in place.
- These implants are designed for temporary usage and are generally removed once the bone can resume its normal function. The physician in charge is responsible for the decision to remove the implants. The risk of adverse effects such as secondary infections, allergies, material fatigue (breakage), implant failure and/or impaired blood circulation, increases with duration of time in situ, non-union, patient weight and activity level.
- Implants must not be machined or altered in any way, unless this is expressly envisaged in the design and in the surgical technique.
- **MAGNETIC RESONANCE (MR) SAFETY AND COMPATIBILITY:** In the event of exposure to foreseeable environmental conditions such as magnetic fields the user and/or patient should be informed of the following:
 - This device has not been evaluated for safety and compatibility in the MR environment.
 - This device has not been tested for heating or migration in the MR environment.
 - There is a potential for heating and migration in the MR environment.
 - There is the potential for metal implants to create MR imaging artifacts in the vicinity of the implant.

ADVERSE EFFECTS

The following adverse effects have been reported:

- Displacement and loosening of implant or implant components
- Deep wound infections
- Vascular complications
- Pseudarthrosis
- Nerve damage
- Corrosion of metal implants
- Wear
- Fatigue fracture
- Inflammatory reactions and osteolysis
- Bone fractures
- Metal sensitivity

STERILIZATION INSTRUCTIONS

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for items supplied non-sterile, for reprocessing reusable devices, or for sterile items that were opened but unused.
 - See Zimmer booklet 97-5000-170-00 (available from your distributor or at zimmer.com) for detailed sterilization guidelines.
- Do not reuse instruments or devices labeled for single use only. Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated

with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.

- All-metal implants without an articulating surface may be resterilized multiple times. This is subject to the exceptions listed below.

DO NOT RESTERILIZE

- Single use only components that have been contaminated with body fluids or debris or previously implanted.
- Components with a packaging expiration date that has been exceeded.
- Do not use the original packaging or one of its components for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing the pouch or wrap.
- Sets of devices may be loaded into general or dedicated trays. Use a medical grade sterilization pouch or wrap (double wrap method - ANSI/AAMI ST79). Sterilization should be performed immediately prior to use.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment chemicals or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.

Recommended Sterilization/Resterilization Specifications

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

Single All-Metal Implants

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Gravity Displacement	121°C (250°F)	30 minutes	15 minutes
Gravity Displacement	132°C (270°F)	15 minutes	
UK Pre-vacuum/Pulsating Vacuum ¹	134°C (273°F)	3 minutes	
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	

¹ This cycle is not for use in the USA.

All-Metal Implants In Sterilization Cases and Trays

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
UK Pre-vacuum/Pulsating Vacuum ¹	134°C (273°F)	3 minutes	30 minutes
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	
Pre-vacuum/Pulsating Vacuum ^{1,2}	134°C (273°F)	18 minutes	

¹ This cycle is not for use in the USA.

² This cycle is not to be used for the inactivation of prions.

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

STORAGE AND HANDLING

- Protective caps or other protective devices must not be removed until immediately before use.
- Assembled implants and implant parts must be thoroughly cleaned and dried before they are stored and placed in cases.
- Non-sterile provided implants are to be stored in the appropriate packaging that protects the implants from damage and must be examined for possible damage before use.

- Implants, implant parts and instruments that can no longer be used may be returned to the manufacturer for proper disposal free of charge (not applicable for the USA).

PATIENT COUNSELING INFORMATION

Complications and/or failure of temporary internal fixation devices are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity or trauma can result in loosening, wear, and/or fracture of the implant. The patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced or removed. The implant is not as strong, reliable, or durable as natural, healthy bone, and will fail under normal weight bearing or load bearing in the absence of complete bone healing.

All names, trademarks, service marks and logos referred to within this package insert are the property of Zimmer, Inc. and/or its respective subsidiaries.

Caution: Federal law restricts this device to sale by or on order of a physician.

not sterile

Symbol for «Contents packed without sterilization»

zimmer

D011500204 Ed. 07/13

NCB® Polyaxial Locking Plate System

Zimmer GmbH
Sulzerallee 8
8404 Winterthur, Switzerland
Telephone +41/ (0)52 262 60 70
Fax +41/ (0)52 262 01 39
www.zimmer.com

Representative in the USA:
Zimmer, Inc.
1800 West Center Street
Warsaw, Indiana, 46580, USA

System: SOP-100 / LB9 / D0115002FREIS 13.01.2014 15:49:37

TR SV ES PT IT EL DE FR NL DA EN

ENGLISH

The CE mark is valid only if it is also printed on the product label.

NCB® Polyaxial Locking Plate System

Important information for the Operating Surgeon

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g. product literature, written surgical technique). Zimmer is not liable for complications arising from the use of the device outside of its indicated uses, surgical technique or judgment, product selection, and similar matters outside the control of Zimmer.

DESCRIPTION

The NCB Polyaxial Locking Plate System is a plate solution for the treatment of complex fractures of the distal femur, proximal humerus and proximal tibia. An additional extension tuberculum minus plate can be assembled to the humeral plate with cerclage wire technique for fixation of AP tuberculum minus fractures.

NCB holes allow for polyaxial screw placement (30°) with subsequent screw locking. Before locking the screws can act as lag screws and be used for fracture reduction. In the locked mode, the plate acts as an internal fixator without contact between the plate and the bone surface reducing the risk of periosteal blood supply impairment.

Femoral, humeral and tibial plates, cancellous, cannulated and locking screws, blind screw inserts and spacers are made from titanium alloy *Protasul*™-64 (Ti6Al4V-ISO 5832-3/ASTM F136). The tuberculum minus plate is made from pure titanium *Protasul*-Ti (ISO 5832-2/ASTM F67), cortical screws from titanium alloy *Protasul*-100 (Ti6Al7Nb ISO 5832-11/ASTM F1295) and the cerclage wire from 316L stainless steel *Protasul*-S (ISO 5832-1/ASTM F138-F139).

INDICATIONS

The NCB Polyaxial Locking Plate System is indicated for temporary internal fixation and stabilization of fractures and osteotomies of long bones.

System: SOP-100 / LB9 / D0115002FREIS 13.01.2014 15:49:37 EN

CONTRAINDICATIONS

- All concomitant diseases that may impair the fixation of the implant and/or the success of the intervention.
- Lack of bone substance or poor bone quality which makes stable seating of the implant impossible.
- Acute or chronic, local or systemic infections.
- Allergy to the implanted material.
- Severe muscular, neural or vascular diseases that endanger the extremities involved.

WARNINGS

- Do not reuse. This device is single patient use only. Reuse may adversely affect device performance.
- Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- Do not use any component if damage is found or caused during setup or insertion.
- Implants and implant parts must only be combined with components belonging to the same system. No liability is accepted for products of third parties that are used by the purchaser or user.
- Use only instruments specifically designed for use with these devices to help ensure accurate surgical implantation.
- The load-bearing capacity of the implant can be compromised by notching, scratching, or striking.
- Do not bend the NCB plates.
- Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions reducing the service life of the implants.
- **To maintain the integrity of the screw / NCB plate construct, two or more screws are required to be placed into each load bearing fracture fragment. If the fracture is in the distal or proximal cancellous region of the bone, use of all screws is recommended.**
- Do not use this product for other than labeled indications (off-label use).

PRECAUTIONS

- Continued surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- These implants are designed for temporary usage and are generally removed once the bone resumes its normal function. The physician in charge is responsible for the decision to remove the implants. The risk of adverse effects such as secondary infections, allergies, material fatigue (breakage), implant failure and/or impaired blood circulation, increases with the time the implant is implanted, delayed union, non-union, patient weight and activity level.
- Implants must not be machined or altered in any way, unless this is expressly envisaged in the design and in the surgical technique.
- This Implant has not been evaluated for safety and compatibility in the MR environment. The Implant has not been tested for heating or migration in the MR environment.

ADVERSE EFFECTS

The following adverse effects have been reported:

- Displacement and loosening of implant
- Deep wound infections
- Vascular complications
- Pseudoarthrosis
- Subclinical nerve damage
- Corrosion of metal implants
- Wear
- Fatigue fracture
- Inflammatory reactions and osteolysis
- Bone fractures
- Metal sensitivity
- Avascular necrosis

STERILIZATION

- These devices are provided either sterile or non-sterile.
- Gamma irradiation is indicated by the "Sterile-R" symbol on the labeling. These devices remain sterile as long as the package integrity has not been violated.
- Inspect each package prior to use and do not use the component if any seal or cavity is damaged, breached, or if the expiration date has been exceeded.
- Once opened, the component must be used, discarded, or resterilized.

- If the packaging is damaged or the sterility expiration date has been reached, the implants must be returned to the manufacturer. (Not applicable for the USA.)

STERILIZATION INSTRUCTIONS

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for items supplied non-sterile, for reprocessing reusable devices, or for sterile items that were opened but unused.
 - See Zimmer booklet 97-5000-170-00 (available from your distributor or at zimmer.com) for detailed sterilization guidelines
- Do not reuse instruments or devices labeled for single use only. Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- All-metal implants without an articulating surface may be resterilized multiple times. This is subject to the exceptions listed below.
- **DO NOT RSTERILIZE:**
 - Single use only components that have been contaminated with body fluids or debris or previously implanted.
 - Components with a packaging expiration date that has been exceeded.
 - Components containing UHMWPE or HDPE.
- Do not use the original plastic cavities or lids for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing the pouch or wrap.
- Sets of devices may be loaded into general or dedicated trays. Use a medical grade sterilization pouch or wrap (double wrap method—ANSI/AAMI ST79). Sterilization should be performed immediately prior to use.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment chemicals or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.

Recommended Sterilization/Resterilization Specifications

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads

Single All-Metal Implants

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Gravity Displacement	121°C (250°F)	30 minutes	15 minutes
Gravity Displacement	132°C (270°F)	15 minutes	
UK Pre-vacuum/Pulsating Vacuum	134°C (273°F)	3 minutes	
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	

This cycle is not for use in the United States.

All-Metal Implants In Sterilization Cases and Trays

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
UK Pre-vacuum/Pulsating Vacuum	134°C (273°F)	3 minutes	30 minutes
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	
Pre-vacuum/Pulsating Vacuum ^{1,2}	134°C (273°F)	18 minutes	

This cycle is not for use in the United States.

² This cycle is not to be used for the inactivation of prions.

Please contact Zimmer at the following number if you have additional questions. In the USA, c 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-26 6131.

STORAGE AND HANDLING

- Implants that are provided sterile must be stored unopened in their original packaging.

System: SOP-100 / LB9 / D0115002FREIS 13.01.2014 15:49:37

- Protective caps or other protective devices must not be removed until immediately before use.
- Assembled implants and implant parts must be thoroughly cleaned and dried before they are stored and placed in cases.
- Implants, implant parts and instruments that can no longer be used may be returned to the manufacturer for proper disposal free of charge. (Not applicable for the USA.)

PATIENT COUNSELING INFORMATION

Complications and/or failure of temporary internal fixation devices are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients and/or with patients who fail to follow through with the required rehabilitation program. Physical activity can result in loosening, wear and/or fracture of the implant. The prospective implant patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient's life, or any particular length of time. Because temporary internal fixation devices are not as strong, reliable, or durable as natural, healthy tissues/bones, all such devices may need to be replaced at some point.

All names, trademarks, service marks and logos referred to within this package insert are the property of Zimmer, Inc. and/or their respective subsidiaries.

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©System. SOP-100 / LB9 / D015002FEN 13.01.2014 15:49:37