ZIMMER BIOMET

Zimmer GmbH

Sulzerallee 8,

CH-8404 Winterthur,

Switzerland


www.zimmer.com

Certification for Information in Dossier

The submitted dossier for the Russian application **NCB implants for osteosynthesis - non-sterile plates** contains the following:

- Instructions for use D011500204
- Instructions for use D011500290
- Instructions for use 87-6203-801-01

I certify that the documents which are submitted are true and accurate to the best of my knowledge. In addition, Zimmer GmbH takes responsibility for any legal obligations related to this submission.

Matthias Bürger
Quality Assurance & Regulatory
Compliance VP
Zimmer GmbH

Official Certification

Seen for authentication of the foregoing signature of

Mr. **Matthias BÜRGER**, Nationality: Germany, in Winterthur, Switzerland,
personally known to us,

who is entered in the Register of Commerce of the Kanton of Zurich as person with the right to sign
per procuration individually for the

Zimmer GmbH, company with limited liability with registered head office in Winterthur ZH,
Switzerland.

The signature was acknowledged before us by an authorized third party, Ms. V. Räber.

The inspection of the commercial register has taken place directly before the official certification by
internet inquiry.

Winterthur, 23rd May 2017
BK no. 1272
Fee CHF 30.00

NOTARIAT OBERWINTERTHUR-WINTERTHUR

~~Hanspeter Lee, Deputy Notary Public~~

Martin Schmuki
Notar

Barcode with Ed.	D011500204	Ed. 2016/09	
NCB® Polyaxial Locking Plate System			
 Zimmer GmbH Sulzerallee 8 8404 Winterthur, Switzerland Telephone +41/ (0)58 854 80 00 Fax +41/ (0)52 244 86 70 www.zimmer.com			
Representative in the USA: Zimmer, Inc. 1800 West Center Street Warsaw, Indiana, 46580, USA			
   			

ENGLISH

The CE mark is valid only if it is also printed on the product label.

NCB® Polyaxial Locking Plate System

Important information for the operating surgeon

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g. product literature, surgical technique). Zimmer is not liable for complications that may arise from the use of the device in circumstances outside Zimmer's control including, but not limited to, deviations from the device's indicated uses or surgical technique.

PRODUCT COMPATIBILITY

Plates, screws and accessories must only be combined with components listed in the product-specific surgical technique.

DESCRIPTION

The NCB Polyaxial Locking Plate System is a plate solution for the treatment of complex fractures of the distal femur, proximal humerus and proximal tibia. NCB holes allow for polyaxial screw placement (30°) with subsequent screw locking. Before locking, the screws can act as lag screws and be used for fracture reduction. In the locked mode, the plates act as an internal fixator without contact between the plate and the bone surface reducing the risk of periosteal blood supply impairment.

An additional extension tuberculum minus plate can be assembled to the humeral plate with cerclage wire technique for fixation of AP tuberculum minus fractures. NCB femoral, humeral and tibial plates, cortical, cancellous and cannulated screws, locking cap, blind screw inserts and spacers are made from titanium alloy Protasul®-64 (Ti6Al4V - ISO 5832-3/ASTM F136). The tuberculum minus plate is made from pure titanium Protasul-Ti (ISO 5832-2/ASTM F67), cortical screws from titanium alloy Protasul-100 (Ti6Al7Nb ISO 5832-11/ASTM F1295) and the cerclage wire from 1.4441 stainless steel Protasul-S (ISO 5832-1/ASTM F138-F139).

INDICATIONS

The NCB Polyaxial Locking Plate System is indicated for temporary internal fixation and stabilization of fractures and osteotomies of long bones.

CONTRAINDICATIONS

- All concomitant diseases that may impair the fixation of the implant and/or the success of the intervention.
- Lack of bone substance or poor bone quality which makes stable seating of the implant impossible.
- Acute or chronic, local or systemic infections.
- Allergy to the implanted material.
- Severe muscular, neural or vascular diseases that endanger the extremities involved.

WARNINGS

- Do not reuse. Implants are for single use only. Reuse may adversely affect device performance.
- Re-use of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- Do not use any component if damage is found or caused during setup or insertion.
- Implants and implant parts must only be combined with components belonging to the same system. No liability is accepted for products of third parties that are used by the purchaser or user.
- Use only instruments specifically designed for use with these devices to help ensure accurate surgical implantation.
- The load-bearing capacity of the implant can be compromised by notching, scratching, or striking.
- Do not bend the NCB plates.
- Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions reducing the service life of the implants.
- **To maintain the integrity of the screw / NCB plate construct, two or more screws are required to be placed into each load bearing fracture fragment. If the fracture is in the distal or proximal cancellous region of the bone, use of all screw holes is recommended.**
- Do not use this product for other than labeled indications.

PRECAUTIONS

- Surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- These implants are designed for temporary usage and are generally removed once the bone can resume its normal function. The physician in charge is responsible for the decision to remove the implants. The risk of adverse effects such as secondary infections, allergies, material fatigue (breakage), implant failure and/or impaired blood circulation, increases with the time the implant is implanted, delayed union, non-union, patient weight and activity level.
- Implants must not be machined or altered in any way, unless this is expressly envisaged in the design and in the surgical technique.
- **MAGNETIC RESONANCE (MR) SAFETY AND COMPATIBILITY:** In the event of exposure to foreseeable environmental conditions such as magnetic fields the user and/or patient should be informed of the following:
 - This device has not been evaluated for safety and compatibility in the MR environment.
 - This device has not been tested for heating or migration in the MR environment.
 - There is a potential for heating and migration in the MR environment.
 - There is the potential for metal implant to create MR imaging artifacts in the vicinity of the implant.

ADVERSE EFFECTS

The following adverse effects have been reported:

- Avascular necrosis
- Bone fractures
- Corrosion of metal implants
- Displacement and loosening of implant
- Impingement
- Implant breakage
- Inflammatory reactions and osteolysis
- Intra-articular screw perforation
- Irritating screw tip
- Joint instability
- Joint stiffness
- Loss of reduction
- Mal-, delayed- and non-union

- Metal sensitivity
- Pain
- Post traumatic osteoarthritis
- Superficial and deep wound infections
- Tendon rupture
- Vascular and neurological complications
- Wear

STERILITY

These devices are provided either sterile or non-sterile.

For non-sterile-packed devices:

Devices provided non-sterile are indicated by the  symbol on the label. It is highly recommended that non-sterile implants are thoroughly cleaned. If the packaging is damaged then ensure that the product has not also been damaged. All implants must be sterilized prior to use. Please follow the cleaning and sterilization instructions described below.

For sterile-packed devices:

Gamma irradiation is indicated by the  symbol on the label. These devices remain sterile as long as the package integrity has not been violated and the expiration date has not been exceeded. Inspect each package prior to use and do not use the component if any seal or cavity is damaged or breached or if the expiration date has been exceeded. Once opened, the component must be used, discarded, or resterilized. If the packaging is damaged or the sterility expiration date has been reached, the implants must be returned to the manufacturer (not applicable for the USA).

RESTERILIZATION INSTRUCTIONS

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for items supplied non-sterile, for reprocessing reusable devices, or for sterile items that were opened but unused.
 - See Zimmer booklet 97-5000-170-00 (available in English from your distributor or at www.zimmer.com) for detailed sterilization guidelines.
- Do not reuse instruments or devices labeled for single use only.
- All-metal implants without an articulating surface may be resterilized multiple times. This is subject to the exceptions listed below.
- **DO NOT RESTERILIZE:**
 - Single use only components that have been contaminated with body fluids or debris or previously implanted.
 - Components with a packaging expiration date that has been exceeded.
- Sterilize only completely cleaned and dried products.
- Do not use the original packaging or one of its components for resterilization. Single devices should be packed in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing the pouch or wrap.
- Sets of devices may be loaded into general or dedicated trays. Use a medical grade sterilization pouch or wrap (double wrap method - ANSI/AAMI ST79). Sterilization should be performed immediately prior to use.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment, chemicals, or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.

Recommended Sterilization/Resterilization Specifications

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

Single, All-Metal Implants

Steam Sterilization

Type	Temperature	Exposure Time	Minimum Dry Time
UK Pre-vacuum/Pulsating Vacuum ¹	134 °C (273 °F)	3 minutes	15 minutes
Pre-vacuum/Pulsating Vacuum	132 °C (270 °F)	4 minutes	

¹ This cycle is not for use in the USA.

All-Metal Implants In Sterilization Cases and Trays

Steam Sterilization

Type	Temperature	Exposure Time	Minimum Dry Time
UK Pre-vacuum/Pulsating Vacuum ¹	134 °C (273 °F)	3 minutes	30 minutes
Pre-vacuum/Pulsating Vacuum	132 °C (270 °F)	4 minutes	
Pre-vacuum/Pulsating Vacuum ^{1,2}	134 °C (273 °F)	18 minutes	

¹ This cycle is not for use in the USA.

² This cycle is not to be used for the inactivation of prions.

For further information on cleaning, refer to Zimmer booklet 97-5000-170-00 – “Orthopaedic Reusable Devices”, available in English from your distributor or at www.zimmer.com.

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

STORAGE AND HANDLING

- It is recommended to store the product in a dry place at room temperature.
- Implants that are provided sterile must be stored unopened in their original packaging.
- Before sterile implants are removed from their packaging, the protective packaging must be examined for possible damage to the sterile barrier. If an expiration date for sterility of the product is indicated, this must be observed. If the sterile barrier system is damaged or the sterility expiration date has been reached, the implants must not be used.
- Implants are extremely sensitive to damage. Therefore, careful handling is strongly recommended. Implants, implant parts and instruments that can no longer be used may be returned to the manufacturer for proper disposal free of charge (not applicable for the USA).

PATIENT COUNSELING INFORMATION

Complications and/or failure of temporary internal fixation devices are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients and/or with patients who fail to follow through with the required rehabilitation program. Physical activity or trauma can result in loosening, wear and/or fracture of the implant. The prospective implant patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant is not as strong, reliable or durable as natural, healthy bone and will fail under normal weight bearing or load bearing in the absence of complete bone healing.

All trademarks and logos referred to within this package insert are the property of Zimmer, Inc. or its affiliates.

not sterile

Symbol for “Contents packed without sterilization»

Caution: Federal law restricts this device to sale by or on the order of a physician.

D011500290 Ed. 07/14

NCB® Straight Narrow Shaft Plate

Zimmer GmbH
Sulzerallee 8
8404 Winterthur, Switzerland
Telephone +41/ (0)52 262 60 70
Fax +41/ (0)52 262 01 39
www.zimmer.com

Representative in the USA:
Zimmer, Inc.
1800 West Center Street
Warsaw, Indiana, 46580, USA

- TR
- SV
- ES
- SK
- SR
- RO
- PT
- PL
- NO
- LT
- IT
- HU
- EL
- DE
- FR
- ET
- NL
- DA
- CS
- HR
- BG
- EN

ENGLISH

EN

The CE mark is valid only if it is also printed on the product label.

NCB® Straight Narrow Shaft Plate

Important information to the operating surgeon

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). Zimmer is not liable for complications that may arise from use of the device in circumstances outside of Zimmer's control including, but not limited to, product selection and deviations from the device's indicated uses or surgical technique.

PRODUCT COMPATIBILITY

Plates, screws and accessories can only be combined with components listed in the product-specific Surgical Technique.

DESCRIPTION

The *NCB* (Non-Contact-Bridging) Straight Narrow Shaft Plate is a line of polyaxial locking plates for the treatment of humeral and tibial shaft fractures, including periprosthetic fractures. The *NCB* System technology allows for polyaxial screw placement (30° cone) with screw locking achieved with the use of locking caps that are threaded into the plate holes. In the locked mode the *NCB* Straight Narrow Shaft Plate acts as an internal fixator without contact between the plate and the bone surface, reducing the risk of periosteal blood supply impairment. This Non-Contact-Bridging concept can be specifically controlled through the use of 1, 2, or 3 mm spacers, which are threaded into the plate holes prior to plate insertion. Plates, screws, blind screw inserts, spacers and locking caps are made of titanium alloy.

INDICATIONS

The *NCB* Straight Narrow Shaft Plate is indicated for temporary internal fixation and stabilization of humeral and tibial shaft fractures and osteotomies, including:

- Periprosthetic fractures
- Comminuted fractures
- Fractures in osteopenic bone
- Nonunions
- Malunions

CONTRAINDICATIONS

- All concomitant diseases that may impair the fixation of the implant and/or the success of the intervention.
- Lack of bone substance or poor bone quality which makes stable seating of the implant impossible.
- Acute or chronic, local or systemic infections.
- Allergy to the implanted materials.
- Severe muscular, neural or vascular diseases that endanger the extremities involved.
- Loose prosthesis, which requires immediate revision.

WARNINGS

- Do not reuse. This device is single use only. Reuse may adversely affect device performance.
- Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- Failure to immobilize a delayed or nonunion of bone will result in excessive and repeated stresses which can cause eventual bending or breakage of the device.

- Do not mix dissimilar metals and alloys that can accelerate corrosion and enhance fracture of implants. It is important that mechanical fit and metallurgical compatibility be considered in selecting mating components.
- Do not use this product for other than labeled indications (off-label use).
- Plates, screws and accessories must only be combined with components listed in the surgical technique. No liability is accepted for products of third parties that are used by the purchaser or the user. It is the responsibility of the user to make sure that the systems are compatible.
- Contact between screws, plates, cables, drilling tools and a previously implanted prosthesis is to be avoided as it may result in the generation of metal particle debris and damage to the existing prosthesis.
- Use only instruments specifically designed for use with these devices to help ensure accurate surgical implantation.
- The load-bearing capacity of the implant can be compromised by notching, scratching or striking the device.
- Do not use any technique for implantation other than that described by Zimmer.
- Do not use any component if damage is found or caused during set up or insertion. Damage could subject the plate to increased stresses, which may lead to device failure.
- Do not use screws in this system for attachment or fixation of the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.
- Drilling and inserting screws into the outer edge of the cortical bone may lead to bone cortex breakage.
- If only unicortical screws are used, additional cables are required.

PRECAUTIONS

- Continued surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- These implants are designed for temporary usage and are generally removed once the bone can resume its normal function. The physician in charge is responsible for the decision to remove the implants. The risk of adverse effects such as secondary infections, allergies, material fatigue (breakage), implant failure and/or impaired blood circulation, explantation related complications increases with the time the implant is implanted, delayed union, non-union, patient weight and activity level.
- Implants must not be machined or altered in any way, unless this is expressly envisaged in the design and in the surgical technique.
- To prevent thread stripping and allow for adequate screw purchase, pass screws as centrally through the bone as possible. In addition, irrigation may be used to help prevent thermal necrosis as the screws are inserted.
- It must be ensured that already implanted intramedullary systems, their fixation and/or anchorage elements are not affected or damaged by screws or drill bits.
- Be aware that bending the plate may decrease its fatigue strength. Furthermore, the locking mechanism of the NCB Plate hole may be damaged and, therefore, may no longer function. Do not use a hole that has been altered by plate contouring for locking.
- **MAGNETIC RESONANCE (MR) SAFETY AND COMPATIBILITY:** IN THE EVENT OF EXPOSURE TO FORESEEABLE ENVIRONMENTAL CONDITIONS SUCH AS MAGNETIC FIELDS THE USER AND/OR PATIENT SHOULD BE INFORMED OF THE FOLLOWING:
 - This device has not been evaluated for safety and compatibility in the MR environment.
 - This device has not been tested for heating or migration in the MR environment.
 - There is a potential for heating and migration in the MR environment.
 - There is the potential for metal implants to create MR imaging artifacts in the vicinity of the implant.

ADVERSE EFFECTS

- Displacement and loosening of implant
- Deep wound infections
- Vascular complications
- Pseudoarthrosis

- Nerve damage
- Corrosion of metal implants
- Wear
- Fatigue fracture
- Inflammatory reactions and osteolysis
- Bone fractures
- Metal sensitivity
- Avascular necrosis

STERILITY

- These devices are provided either sterile or non-sterile.
- Gamma irradiation is indicated by the **[STERILE]** symbol on the labeling. These devices remain sterile as long as the package integrity has not been violated and before the expiration date is reached.
- Inspect each package prior to use and do not use the component if any seal or cavity is damaged, breached, or if the expiration date has been exceeded.
- Once opened, the component must be used, discarded, or resterilized.
- If the packaging is damaged or the sterility expiration date has been reached, the implants must be returned to the manufacturer (not applicable for the USA).

STERILIZATION INSTRUCTIONS

- These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines. They should be used for items supplied non-sterile, for reprocessing reusable devices, or for sterile items that were opened but unused.
 - See Zimmer booklet 97-5000-170-00 (available from your distributor or at zimmer.com) for detailed sterilization guidelines.
- Do not reuse instruments or devices labeled for single use only. All- metal implants without an articulating surface may be resterilized multiple times. This is subject to the exceptions listed below.
- DO NOT RESTERILIZE:
 - Single use only components that have been contaminated with body fluids or debris or previously implanted.
 - Components with a packaging expiration date that has been exceeded.
- Do not use the original packaging (e.g. plastic cavities or lids) for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing.
- Sets of devices may be loaded into general or dedicated trays. Use a medical grade sterilization pouch or wrap (double wrap method - ANSI/AAMI ST79). Sterilization should be performed immediately prior to use.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment chemicals or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.

Recommended Sterilization/Resterilization Specifications

Follow the sterilizer manufacturer's instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

Single, All-Metal Implants

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Gravity Displacement	250°F (121°C)	30 minutes	15 minutes
Gravity Displacement	270°F (132°C)	15 minutes	
UK Pre-vacuum/Pulsating Vacuum ¹	273°F (134°C)	3 minutes	
Pre-vacuum/Pulsating Vacuum	270°F (132°C)	4 minutes	

¹ This cycle is not for use in the USA.

All-Metal Implants In Sterilization Cases and Trays

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
UK Pre- vacuum/Pulsating Vacuum ¹	273°F (134°C)	3 minutes	30 minutes
Pre-vacuum/Pulsating Vacuum	270°F (132°C)	4 minutes	
Pre-vacuum/Pulsating Vacuum ^{1,2}	273°F (134°C)	18 minutes	

¹ This cycle is not for use in the USA.

² This cycle is not to be used for the inactivation of prions.

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

STORAGE AND HANDLING

- Implants/instruments must be stored unopened in their original packaging and must be examined for possible damage before use.
- Before sterile implants/instruments are removed from their packaging, the protective packaging must be examined for possible damage to the sterile barrier. If an expiration date for sterility of the product is indicated, this must be observed. If the sterile barrier system is damaged or the sterility expiration date has been reached, the implants/instruments must not be used.
- Implants/instruments are extremely sensitive to damage. Even small scratches or marks left by impacts on the surfaces will cause excessive wear and may give rise to complications. Extremely careful handling is therefore strongly recommended.

PATIENT COUNSELING INFORMATION

Complications and/or failure of implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity or trauma can result in loosening, wear, and/or fracture of the implant. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities, and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant is not as strong, reliable, or durable as natural, healthy bone, and will fail under normal weight bearing or load bearing in the absence of complete bone healing.

All names, trademarks, service marks and logos referred to within this package insert are the property of Zimmer, Inc. and/or their respective subsidiaries.

Rx only
physician.

Caution: Federal law restricts this device to sale by or on the order of a

NOT STERILE Symbol for «Contents packed without sterilization»

Barcode	87-6203-801-01	Ed. 2014/09	
NCB® PERIPROSTHETIC FEMUR POLYAXIAL LOCKING PLATE SYSTEM			
 Zimmer GmbH Sulzerallee 8 8404 Winterthur, Switzerland Telephone +41/ (0)52 262 60 70 Fax +41/ (0)52 262 01 39 www.zimmer.com	 Zimmer, Inc. 1800 West Center Street Warsaw, Indiana, 46580, USA	EC REP Zimmer, U.K. Ltd. 9 Lancaster Place South Marston Park Swindon, SN3 4FP, UK	
 0086   			

ENGLISH

The CE mark is valid only if it is also printed on the product label.

NCB® PERIPROSTHETIC FEMUR POLYAXIAL LOCKING PLATE SYSTEM

Before using a product placed on the market by Zimmer, the operating surgeon should study carefully the following recommendations, warnings and instructions, as well as the available product-specific information (e.g., product literature, written surgical technique). The manufacturer, the importer and the suppliers of Zimmer products are not liable for complications or other effects that may occur for reasons such as incorrect indications or surgical technique, use of the device outside of its indicated uses, unsuitable choice of material or handling thereof, unsuitable use or handling of the instruments, asepsis and so on. The operating surgeon is responsible for any such complications or other consequences.

PRODUCT COMPATIBILITY

Plates, screws and accessories can only be combined with components listed in the product-specific Surgical Technique.

DESCRIPTION

The NCB (Non-Contact Bridging) Periprosthetic Femur Polyaxial Locking Plate System is a line of polyaxial locking plates for the treatment of femur fractures, particularly of periprosthetic femur fractures. It consists primarily of Proximal Femur Plates, Distal Femur Plates, and Curved Shaft Plates. The NCB System technology allows for polyaxial screw placement (30° cone) with screw locking achieved with the use of locking caps that are threaded into the plate holes. In the locked mode the NCB Periprosthetic Plate acts as an internal fixator without contact between the plate and the bone surface, reducing the risk of periosteal blood supply impairment. This Non-Contact Bridging concept can be specifically controlled through the use of 1, 2, or 3 mm spacers, which are threaded into the plate holes prior to plate insertion. Plates, screws, blind screw inserts, spacers and locking caps are made of titanium alloy.

INDICATIONS/INTENDED USE

The NCB Periprosthetic Femur Polyaxial Locking Plate System is indicated for temporary internal fixation and stabilization of fractures and osteotomies of long bones, including:

- Periprosthetic fractures
- Comminuted fractures
- Supracondylar fractures
- Fractures in osteopenic bone
- Nonunions
- Malunions

CONTRAINDICATIONS

- All concomitant diseases that may impair the fixation of the implant and/or the success of the intervention.
- Lack of bone substance or poor bone quality which makes stable seating of the implant impossible.
- Acute or chronic, local or systemic infections.
- Allergy to the implanted materials.
- Severe muscular, neural or vascular diseases that endanger the extremities involved.
- Loose prosthesis, which requires immediate revision.

WARNINGS

- Do not reuse. This device is single use only.
- Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- Do not use this product for other than labeled indications (off-label use).
- Plates, screws and accessories must only be combined with components belonging to the same system. No liability is accepted for products of third parties that are used by the purchaser or the user. It is the responsibility of the user to make sure that the systems are compatible.
- Contact between screws, plates, cables, drilling tools and a previously implanted prosthesis is to be avoided as it may result in the generation of metal particle debris and damage to the existing prosthesis.
- Use only instruments specifically designed for use with these devices to help ensure accurate surgical implantation.
- The load-bearing capacity of the implant can be compromised by notching, scratching or striking the device.
- Do not use any technique for implantation other than that described by Zimmer.
- Do not use any component if damage is found or caused during set up or insertion. Damage could subject the plate to increased stresses, which may lead to device failure.
- Do not use screws in this system for attachment or fixation of the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.
- Drilling and inserting screws into the outer edge of the cortical bone may lead to bone cortex breakage.
- If only unicortical screws are used, additional cables are required.
- Do not insert three screws in one diagonal three hole pattern.

PRECAUTIONS

- Continued surveillance for new or recurrent sources of infection should be continued as long as the device is in place.
- These implants are designed for temporary usage and are generally removed once the bone can resume its normal function. The physician in charge is responsible for the decision to remove the implants. The risk of adverse effects such as secondary infections, allergies, material fatigue (breakage), implant failure and/or impaired blood circulation increases with the time the implant is implanted, delayed union, non-union, patient weight and activity level.
- Implants must not be machined or altered in any way, unless this is expressly envisaged in the design and in the surgical technique.
- To prevent thread stripping and allow for adequate screw purchase, pass screws as centrally through the bone as possible. In addition, irrigation may be used to help prevent thermal necrosis as the screws are inserted.

- It must be ensured that already implanted intramedullary systems, their fixation and/or anchorage elements are not affected or damaged by screws or drill bits.
- The *NCB* Periprosthetic Femur Polyaxial Locking Plates are anatomically shaped. Be aware that bending the plate may decrease its fatigue strength. Furthermore, the locking mechanism of the *NCB* Plate hole may be damaged and, therefore, may no longer function. Do not use a hole that has been altered by plate contouring for locking.
- **MAGNETIC REONANCE (MR) SAFETY AND COMPATIBILITY:** In the event of exposure to foreseeable environmental conditions such as magnetic fields the user and/or patient should be informed of the following:
 - This device has not been evaluated for safety and compatibility in the MR environment.
 - This device has not been tested for heating or migration in the MR environment.
 - There is a potential for heating and migration in the MR environment.
 - There is the potential for metal implants to create MR imaging artifacts in the vicinity of the implant.

ADVERSE EFFECTS

- Displacement and loosening of implant
- Deep wound infections
- Vascular complications
- Pseudoarthrosis
- Nerve damage
- Corrosion of metal implants
- Wear
- Fatigue fracture
- Inflammatory reactions and osteolysis
- Bone fractures
- Metal sensitivity
- Avascular necrosis

STERILITY

- These devices are provided either sterile or non-sterile.
- Gamma irradiation is indicated by the "Sterile-R"  symbol on the labeling. These devices remain sterile as long as the package integrity has not been violated and before the expiration date is reached.
- Inspect each package prior to use and do not use the component if any seal or cavity is damaged, breached, or if the expiration date has been exceeded.
- Once opened, the component must be used, discarded, or resterilized.
- If the packaging is damaged or the sterility expiration date has been reached, the implants must be returned to the manufacturer. (Not applicable for the USA.)

STERILIZATION INSTRUCTIONS

These sterilization instructions are consistent with ANSI/AAMI/ISO standards and guidelines.

- They should be used for items supplied non-sterile, for reprocessing reusable devices, or for sterile items that were opened but unused.
 - See Zimmer booklet 97-5000-170-00 (available from your distributor or at zimmer.com) for detailed sterilization guidelines.
- Do not reuse instruments or devices labeled for single use only. Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- All-metal implants without an articulating surface may be resterilized multiple times. This is subject to the exceptions listed below.
- **DO NOT RESTERILIZE:**
 - Single use only components that have been contaminated with body fluids or debris or previously implanted.
 - Components with a packaging expiration date that has been exceeded.
- Do not use the original packaging or one of its components for resterilization. Single devices should be packaged in a medical grade sterilization pouch or wrap which conforms to the

recommended specifications for steam sterilization provided in the table below. Ensure that the pouch or wrap is large enough to contain the devices without stressing the seals or tearing the pouch or wrap.

- Sets of devices may be loaded into general or dedicated trays. Use a medical grade sterilization pouch or wrap (double wrap method - ANSI/AAMI ST79). Sterilization should be performed immediately prior to use.
- Items made from titanium and titanium alloys can form oxide layers from steam boiler treatment chemicals or detergent residues. While these oxides are biocompatible, they can obliterate etchings and stampings.

Recommended Sterilization/Resterilization Specifications

Follow the sterilizer manufacturer’s instructions for loading patterns and selection of sterilization parameters. Drying times vary according to load size and should be increased for larger loads.

Single All-Metal Implants

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
Gravity Displacement	121°C (250°F)	30 minutes	15 minutes
Gravity Displacement	132°C (270°F)	15 minutes	
UK Pre-vacuum/Pulsating Vacuum ¹	134°C (273°F)	3 minutes	
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	

¹ This cycle is not for use in the USA.

All-Metal Implants In Sterilization Cases and Trays

Steam Sterilization

Type	Minimum Temperature	Minimum Exposure Time	Minimum Dry Time
UK Pre-vacuum/Pulsating Vacuum ¹	134°C (273°F)	3 minutes	30 minutes
Pre-vacuum/Pulsating Vacuum	132°C (270°F)	4 minutes	
Pre-vacuum/Pulsating Vacuum ^{1,2}	134°C (273°F)	18 minutes	

¹ This cycle is not for use in the USA.

² This cycle is not to be used for the inactivation of prions.

Please contact Zimmer at the following number if you have additional questions. In the USA, call 1-800-348-2759. For calls outside the USA, call the local international access code +1-574-267-6131.

PATIENT COUNSELING INFORMATION

Complications and/or failure of prosthetic implants are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients who fail to follow through with the required rehabilitation program. Physical activity can result in loosening, wear, and/or fracture of the implant. The prospective implant patient must be counseled about the capabilities of the implant and the impact it will have on his or her lifestyle. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail or need to be replaced. The implant may not last the rest of the patient’s life, or any particular length of time. Because prosthetic implants are not as strong, reliable, or durable as natural, healthy tissues/bones, all such

devices may need to be replaced at some point.

All names, trademarks, service marks and logos referred to within this package insert are the property of Zimmer, Inc. and/or their respective subsidiaries.

The manufacturer for each *NCB* Periprosthetic Femur Polyaxial Locking Plate System component is listed in the table below:

Component	Manufacturer
Plates	Zimmer GmbH
Unicortical and Deep Thread Screws	Zimmer, Inc.
All other Screws	Zimmer GmbH
Locking Caps and Spacers	Zimmer GmbH

not sterile

Symbol for «Contents packed without sterilization

R_xonly

Caution: Federal law restricts this device to sale by or on order of a physician

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© «Зиммер ГмбХ»