7 Monalee Manor, Knocknacarra, Galway

T: 086 8152932

E: jtxseymour@gmail.com

TO ALL TO WHOM THESE PRESENTS SHALL COME


I, James T Seymour Notary Public duly authorised, admitted and sworn for the County and City of Galway, Ireland and practising at 7 Monalee Manor, Knocknacarra, Galway, County of Galway, Ireland

DO HEREBY CERTIFY that the paper writing hereto annexed and upon which I have signed my name and affixed my seal of office is a true copy of a Declaration Letter issued by Medtronic Inc

AND I HAVE CERTIFIED as herein at the request of Medtronic Ireland having its place of business at Parkmore Business Park West in the City of Galway, Republic of Ireland, which is the registered business name of Medtronic Vascular Galway Unlimited Company having its registered office at 10 Earlsfort Terrace, Dublin 2, Ireland which is a subsidiary of Medtronic plc (formerly known as Medtronic Inc.) and for the benefit of the proper authorities of Russia and of all to whom it shall come.

IN FAITH AND TESTIMONY whereof I have hereunto set my hand and affixed my Seal of Office.

Dated at 7 Monalee Manor, Knocknacarra, Galway
In the County of Galway
In Ireland
This 9th day of February 2024

James T. Seymour
Notary Public

7 Monalee Manor, Knocknacarra
Galway
Ireland
Commissioned for Life.

Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879

Declaration letter

To whom it concerns:

It is hereby certified that the document attached is the Operating documentation (Instruction for Use and Appendix to the Instruction for Use) for the Symplicity Spyral™ Multi-Electrode Renal Denervation Catheter.

Yours sincerely,

**Zameer
Judge**

Digitally signed
by Zameer Judge
Date: 2024.02.09
11:29:20 Z

Signed on behalf of Medtronic Inc.

Name: Zameer Judge

Title: Sr. Regulatory Affairs Specialist

Medtronic

Symlicity Spyral™

Multi-Electrode Renal Denervation Catheter
Катетер за бубречна денервација с множествено електрода
Viceelektroodový katéter pro renální denervaci
Katheter med flere elektroder til renal denervering
Multi-Elektroden-Renaldenervierungskatheter
Καθετήρας νεφρικής απενεύρωσης πολλαπλών ηλεκτροδίων
Catèter multi-electrode de denervación renal
Miltielektroodiline neeruarterite denervatsioonikateeter
Munuksien denervatiossa käytettävä monielektroodinen katetri
Catheter de dénervation rénale multi-électrodes
Katheter za renalnu denervaciju s više elektroda
Többelektrodsos renális denervációs katéter
Catetere multi-electrodo per denervazione renale
알티전극 신장 신경 차단 카테터
Katheter mit mehreren Elektroden zur Nierendennervation
Daudzelektrodu nieru denervācijas katērs
Katheter za renalnu denervaciju s više elektroda
Multi-elektrodekateter voor renale denervatie
Multi-elektrodekateter for renal denervering
Wielopokłowy cewnik do denervacji nerek
Cateter de desnervação renal multi-eletródo
Cateter de desnervação renal multipolos
Kateter multi-elektrod de denervare renală
Многосекторный катетер для абляции почечных артерий
Katéter na renálnu denerváciu s viacerými elektródami
Kateter za renalnu denervaciju z več elektrodami
Katheter for renal denervering med flere elektroder
Çok Elektrodlu Renal Denervasyon Kateteri
Биполярноэлектродный катетер для нервової денервації
多电极血管交感神经技术导管

Instructions for Use • Инструкции за употреба • Návod k použití • Bruksanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso •
Käyttöohjeet • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • 사용 지침 • Neudodijne instrukcije •
Lietojuma pamācība • Упатства за употреба • Gebrauchsanweisung • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização •
Instruções de utilização • Instrucţiuni de utilizare • Инструкции по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu •
Bruksanvisning • Kallanin Tallennet • Инструкції з експлуатації • 使用說明

CE
2797

$\overline{C \in}$

STERILE	R
---------	---

②

②

LOT

Quantity / Количество / Množství / Antal / Menge / Ποσότητα / Cantidad / Kogus / Määrä / Quantité / Količina / Mennyiség / Quantità / 수량 / Kiekis /
Daudzums / Колличина / Aantal / Antall / Ilość / Quantidade / Quantidade / Cantitate / Количество / Množstvo / Količina / Količina / Mängd / Miktar / Кількість /
数量

Manuel Control / Руководство / Ручно одлагање / Manuel stop / Manuelle Steuerung / Xapodavetou stágyou / Control manual / Ykätsä julkimise /
Manuelinen ohjaus / Controllo manual / Ruchno upravljajnje / Manuelis kassaförord / Controllo manuale / 手動 / Runkine valdymas / Manuelis valdība /
Ручное управление / Handmatige bediening / Manueli bediening / Stawowe ręczne / Controlle manual / Controllo manuale / Comandă manuală / Ручное
управление / Manuelis valdymas / Ruchno upravljajnje / Ruchno kontrola / Manueli kontrolli / Manuel Control / Ручное управление / 手動制御

Symplcity Spyrat™

Multi-Electrode Renal Denervation Catheter

1. Product Description

The Symplcity Spyrat™ multi-electrode renal denervation catheter is designed to be used with the Symplcity G3™ renal denervation radio frequency (RF) generator. The catheter connects to the generator using the integrated cable attached to the catheter handle. The catheter requires the use of a 0.36-mm (0.014-in) guidewire for delivery, preferably without hydrophilic coating. For a straighter electrode array during delivery, Medtronic recommends using an extra support guidewire such as the Medtronic Thunder™ guidewire. In addition, a standard dispersive electrode (also known as a neutral electrode, return electrode pad, or grounding pad) must be placed on the patient and connected to the generator for the therapy to be delivered. Other ancillary devices include the generator, remote control, and DVI-D cable, which all come in the package, as well as the optional foot switch and cart which are both available separately. The catheter has an effective length of 117 cm and is compatible with a 6-Fr guide catheter, and is designed for treating in vessels with diameters ranging from 3 to 8 mm. As shown in Table 1, the catheter features 4 gold radiopaque electrodes at the spiral (helical) distal end. In the straight configuration, the distance between the electrodes is 6.5 mm. The electrodes are deployed into a spiral (helical) shape by partially retracting the guidewire proximal to the spiral section of the catheter. The treatment length (the distance between electrodes 1 and 4) of the catheter changes with the vessel diameter (Table 1). A radiopaque tip marker is located 1 mm proximal to the catheter tip and assists in the positioning of the catheter using fluoroscopic guidance. The catheter also features a straightening tool that facilitates safe insertion of the guidewire into the catheter (Figure 1). This tool is located near the handle and slides along the catheter shaft to straighten the distal end.

Table 1.


Symplcity Spyrat™ multi-electrode renal denervation catheter electrode configuration		Treatment length: Distance between electrodes 1 and 4 as a function of deployed diameter.	
		Vessel Diameter (mm)	Treatment Length (mm)
	1. Electrode 1 (Distal)	3	21
	2. Electrode 2	4	20
	3. Electrode 3	5	20
	4. Electrode 4 (Proximal)	6	19
		7	18
		8	17

Figure 1. Straightening tool used over the distal portion of the Symplcity Spyrat™ catheter

The generator is represented in Figure 2. The front panel touch screen shows information such as temperature, impedance, ablation time, and messages. The front panel also features an RF activation button. Electrode identifiers on the generator screen correspond to each electrode on the catheter (refer to Table 1 and Figure 2). The generator touch screen and remote control allow the user to navigate different options, such as the selection/deselection of channels, viewing previous ablation data sets, or selecting the left or right kidney. The generator uses an automated algorithm to control the power and time settings for use with the catheter. Refer to the G3 generator user manual for further information.

Figure 2. Representative image of the Symplcity G3™ renal denervation RF generator displaying the impedance monitoring screen

2. Intended Use

The Symplcity G3™ Renal Denervation RF Generator when used with the Symplcity Spyrat™ multi-electrode renal denervation catheter is intended to deliver low-level radio frequency (RF) energy through the wall of the renal artery to denervate the human kidney.

3. Indications for Use

The Symplcity Spyrat™ multi-electrode renal denervation catheter is indicated for the treatment of uncontrolled hypertension.

4. Conditions for Use

- The catheter is intended for single use only.
- The catheter is intended for use only with the Symplcity G3™ renal denervation RF generator.
- The catheter has a limited shelf life. The product must be used on or before the Use By date provided on the label.
- Prior to use, the product should be stored in a cool, dry place. The product should not be exposed to organic solvents, ionizing radiation, or ultraviolet light. Carefully inspect the sterile package for damage before opening. Do not use if the package has been damaged or opened.

5. How Supplied

- The catheter is contained in a dual tray configuration. The inner tray retains the catheter while the outer tray with a sealed Tyvek® lid provides a sterile barrier.
- The generator, remote control, DVI-D cable, optional foot switch, and optional cart are nonsterile and reusable. These items are supplied separately from the catheter.
- Items that are required to complete the treatment, though not supplied, include a 0.36-mm (0.014-in) guidewire, a standard adult-size dispersive electrode, and a sterile bag to cover the remote control if used in the sterile field, as well as any other standard items used to aid percutaneous transluminal catheterization in renal arteries. In addition, the following accessories are needed to gain access to the target vessels: a 6-Fr guide catheter, an introducer sheath, a stopcock sidearm, and a Tuohy-Borst adapter.

6. Risks

Intervention risks: death, cardiopulmonary arrest, heart rhythm disturbances including bradycardia, formation of blood clots and/or embolism (which may result in ischemic events such as myocardial infarction, pulmonary embolism, stroke, kidney damage, or peripheral ischemia), retroperitoneal hematoma, hematoma, bruising, bleeding, arterial damage, arterial spasm, arterial stenosis, arterial dissection or perforation, renal artery aneurysm, kidney perforation, pseudoaneurysm, AV fistula, pain, skin burns, and thermal injury to the vasculature or other structures from energy application.

Contrast agents, narcotics, anxiolytics, other pain medications and anti-vasospasm agents are also commonly used during the procedure or after the procedure; use of these agents is associated with known risks.

Other treatment risks include: proteinuria, hematuria, electrolyte disturbances, worsening renal function, hypotension, hypertension, orthostatic hypotension, hypotension causing end organ hypoperfusion, nausea, and vomiting.

Biological hazards: risks of infection, toxicity, abnormal hematology profile, allergy, hemorrhage, and pyrogenicity.

Environmental hazards: consistent with standard hospital protocols for proper use and disposal of biological wastes.

Radiation hazards: consistent with normal use of x-ray during interventional procedures.

7. Warnings and Precautions

7.1. Related to the Use of Radio Frequency in Catheterization Laboratories

- Radio frequency surgery uses high-frequency output. Do not perform procedures if flammable or explosive media are present, such as flammable anesthetics or skin preparation agents.

- Interference produced by the operation of high-frequency surgical equipment may adversely influence the operation of other electronic medical equipment such as monitors and imaging systems.
- Radio frequency surgery may produce a hazardous electrical output. This equipment is for use only by qualified medical personnel trained in the use of this equipment.

7.2. Related to Interventional Techniques

- A thorough understanding of the technical principles, clinical applications, and risks associated with vascular access techniques and percutaneous transluminal catheterization in renal arteries is necessary before using this device. Physicians should be familiar with techniques used to mitigate potential procedural problems that could be encountered while treating renal arteries such as arterial dissection or perforation, or kidney perforation. Ensure that accessories and products that are typically used in such situations are available.
- Ensure that the guide catheter is flushed with heparinized saline between each treatment.
- Prior to use, do not flush the catheter lumen or the catheter while in the hoop. Do not wipe the spiral section of the catheter.
- Be careful to avoid excessive exposure of the patient to contrast agents.
- Avoid using ionized contrast agent when performing renal denervation.

7.3. Related to Patient

- The catheter has not been evaluated in patients who are pregnant, nursing, plan to become pregnant, and in patients with Type I diabetes mellitus, prior renal angioplasty, indwelling renal stent, aortic grafts, or abnormal renal anatomy.
- Avoid use of the catheter in individuals in whom a reduction in blood pressure would be considered hazardous (such as those with hemodynamically significant valvular heart disease).
- Implantable pacemakers and implantable cardioverter/defibrillators (ICDs) may be adversely affected by RF ablation. Consider deactivating ICDs during ablation, having temporary external sources of pacing and defibrillation available during ablation, and performing a complete analysis of the implanted device's function after ablation.
- Avoid treating in arteries with a diameter less than 3 mm or greater than 8 mm.
- Avoid treating in arteries with significant disease or with flow-limiting obstructions.
- Safe use of monopolar radio frequency surgery demands proper coupling of the dispersive electrode to the patient. Follow all of the manufacturer's directions for the placement of the dispersive electrode and for proper insulation between the patient and any metal surfaces.
- Failure to achieve good skin contact with the entire adhesive surface of the dispersive electrode may result in a burn or high impedance measurements.
- The patient should not come into contact with metal parts that are grounded or have an appreciable capacitance to ground (such as operating table supports, etc.). The use of antistatic sheeting is recommended for this purpose.
- The patient's heart rate may drop during the ablation procedure. The administration of medication such as atropine might be needed.
- The patient is likely to experience pain when radio frequency energy is delivered. Proper pain medication should be administered at least 5 minutes prior to ablating renal nerves.
- Nitroglycerine should be administered to prevent vessel spasm during ablation, if not contraindicated.
- An activated clotting time (ACT) of at least 250 seconds should be maintained during the procedure; proper use of systemic anticoagulation medication is required.

7.4. Related to Ablation Catheter and Generator

- The generator should be powered on and allowed to complete a system self-test prior to placing the catheter within the vasculature.
- The catheter is intended for single patient use only. Do not resterilize or reuse. Reuse, reprocessing or resterilization may compromise device integrity and functionality and may create the risk of transmission of infectious diseases from one patient to another, which may result in injury, illness, or death of the patient.
- Do not advance the catheter against resistance.
- Avoid advancing the guidewire too distal to reduce risk of damaging kidney. Similarly, Medtronic recommends using guidewires that are not hydrophilic coated to prevent unintentional damage to the kidneys or renal arteries.
- The catheter should be deployed while using fluoroscopy. Avoid torquing the catheter beyond 180 degrees to prevent guidewire entanglement.
- Avoid treating renal artery segments that may be in close proximity to structures such as the ureter or kidney.

7.5. Related to RF Treatment

- Remove any guidewires that are not contained within the catheter (such as a buddy wire) from the treatment site before activating the RF delivery.
- In the event of multiple ablations, the most distal ablation at the subsequent treatment site should be located approximately 5 mm proximal to the most proximal ablation performed during the preceding ablation.
- Do not ablate if the electrodes are in contact with each other per fluoroscopic observation.
- Avoid occluding renal blood flow during RF delivery.
- During RF delivery, do not move the catheter or guidewire and do not inject saline or contrast agent.
- Increased vessel reactivity, such as spasm, may be encountered when treating in areas with reduced blood flow.
- In the event that the generator stops delivering energy due to high temperature, take an image of the vessel to ensure there is no spasm or occlusion prior to repositioning the catheter in a different section of the artery.
- Do not touch a catheter electrode and the dispersive electrode at the same time during energy delivery as this may result in superficial skin burns.
- Do not allow a catheter electrode or dispersive electrode to come into contact with a metal instrument or surface during energy delivery as this may result in superficial skin burns.

Please consult the generator user manual for additional warning and precautions.

8. References

- The physician should consult recent literature on current medical practices for renal artery interventions.
- The physician should consult the generator user manual for more details on how to operate the generator.

9. Instructions for Use

Closely follow these Instructions for Use and consult the generator user manual for additional instructions for use.

9.1. Equipment and Procedure Preparation

1. Install the generator on a cart or table. For proper equipment ventilation, position the generator more than 30 cm (12 in) away from a wall and do not cover the generator while in use.
2. If the use of a remote control and/or foot switch is desired, connect the remote control and/or foot switch into the respective receptacles on the rear panel of the generator. If desired, the information displayed on the touch screen can also be projected on a cathlab monitor by connecting the DVI-D cable between the rear panel of the generator and the cathlab monitor.
3. Plug the power cable into the back panel of the generator and turn it on by pressing the on/off switch also located on the back panel. Make sure that no catheter is connected to the generator while the generator is being turned on.
4. Check for any system indicator messages or warnings (such as fault or status lights). Following a system self-test, the system is in the STANDBY state and no measurements are possible. After a successful self-test, the front panel will display a screen prompting the user to connect a catheter to the generator.
5. Gather the accessories needed for the procedure, such as dispersive electrode, 5-Fr guide catheter, introducer sheath, 0.36-mm (0.014-in) guidewire, stopcock sidearm, Tuohy-Borst adapter, as well as any other standard items used to aid percutaneous transluminal catheterization in renal arteries.
6. Gather the medications needed for the procedure, such as pain medications, atropine, nitroglycerine, and heparin.

9.2. Patient Preparation

1. Prepare the patient using standard techniques for electrosurgery and catheterization. Ensure the patient's entire body, including extremities, is insulated from contact with grounded metal parts. Closely follow instructions provided by the manufacturer of the dispersive electrode. The dispersive electrode should be placed on the thigh or other nonbony area of the body and should be outside of the angiographic field of view. Shave the patient's skin if necessary for good contact between the dispersive electrode and the skin. Failure to achieve good skin contact by the entire adhesive surface of the dispersive electrode may result in a burn or high impedance measurements. Do not apply the dispersive electrode where fluid may pool.

Note: A disposable dispersive electrode specified for adult use that complies with IEC 60601-2-2 is required.

2. Connect the dispersive electrode to the generator using the receptacle located on the side panel.
3. Ensure that the patient has intravenous (IV) access for drug administration during the procedure. Prior to starting the procedure, administer appropriate systemic anticoagulation (such as heparin) to the patient. An activated clotting time (ACT) of at least 250 seconds should be maintained during RF delivery.
4. Administer pain medication at least 5 minutes prior to ablation. Check vital signs throughout the procedure.
5. Prepare the patient for catheter placement using standard interventional techniques. Connect the Tuohy-Borst adapter to the guide catheter and the stopcock sidearm to the Tuohy-Borst adapter.

Note: If the remote control is being used, insert it into a sterile bag and place it within the sterile field using standard aseptic techniques.

6. Under fluoroscopy, inject contrast in both renal arteries to assess anatomy.
7. Determine whether the arteries are suitable for treatment.


9.3. Catheter Insertion in Renal Artery

1. Using aseptic technique, carefully remove the seal on the outer tray and place the inner tray containing the catheter into the sterile field.
2. Once the tray containing the catheter is in the sterile field, carefully remove the lid by pulling on the lid's pull tab to gain access to the catheter and integrated cable.
3. Remove the coiled cable from the tray and place on a stable sterile surface. Grip the catheter handle with one hand and the hoop with the other hand. Carefully remove the handle and hoop from the tray and place on the stable sterile surface next to the coiled cable.
4. Remove the twist-tie clip from the coiled portion of the cable and pass the integrated cable out of the sterile field for an assistant to connect the cable to the appropriate receptacle on the side panel of the generator. The cable should be secured to the table or drape using a towel clamp, hemostats, or equivalent to help prevent movement of the catheter and handle.
5. An assistant outside the sterile field must perform patient selection on the touch screen (new patient or same patient).
6. Advance a 0.35-mm (0.014-in) guidewire into the target vessel.
 - It is recommended to use only guidewires with a flexible distal tip that are not hydrophilic coated to avoid kidney perforation.
7. Remove the catheter from the hoop, ensure that the straightening tool stays with the handle when pulling the catheter out of the hoop. Inspect the catheter for damage.
 - If the catheter is damaged, do not use.
 - Do not advance the catheter into the hoop after full or partial removal from the hoop. If advanced, fully remove the catheter from the hoop and inspect for damage. If damaged, replace the catheter.
 - Prior to use, do not flush the catheter lumen or the catheter while in the hoop. Do not wipe the spiral section of the catheter.
8. Slide the straightening tool over the spiral portion of the catheter as illustrated in Figure 1, making sure that approximately 5 mm of the catheter tip still protrudes from the distal end of the straightening tool.
 - If excessive resistance is felt while advancing the straightening tool over the spiral section of the catheter, stop, retract the straightening tool, and assess for damage.
 - If the electrodes or the distal end of the catheter are damaged, replace the catheter.
9. Squeeze the distal flare of the tool to secure the catheter. Carefully insert the proximal end of the guidewire through the tip of the catheter. Continue to pass the guidewire through the catheter until the guidewire exits through the rapid exchange port. This exit port is located 30 cm proximal to the distal tip of the catheter.
 - If the guidewire does not exit from the rapid exchange port, remove the guidewire from the catheter and reinsert the guidewire while assessing for device breaches.
 - If the catheter is breached or damaged, replace the catheter and guidewire.
10. Once the guidewire has exited the rapid-exchange port, return the straightening tool by the handle to prevent interference with the guidewire.
11. Use nitroglycerine prior to advancing the catheter in the artery to reduce risk of arterial spasm, if not contraindicated.
12. Advance the catheter over the guidewire through the guide catheter.
 - If using a 55-cm guide catheter, the catheter tip will exit the guide catheter when the shaft marker enters the rotating hemostatic valve.
13. When all four electrodes exit the guide catheter, the impedance monitoring screen (Figure 4) will then be displayed.

Note: If the display does not continue to the impedance monitoring screen, follow these steps:



 - a. Check the catheter position and ensure that all 4 electrodes are outside of the guide catheter.
 - b. Verify appropriate dispersive electrode connection and contact with patient.
 - c. If the previous steps do not result in the display of the impedance monitoring screen, try moving the dispersive electrode to the patient's flank. If needed, replace the dispersive electrode.

9.4. Achieving adequate wall contact

1. Guidewire inserted beyond the distal tip (spiral not deployed). Electrode 1 is located in the bifurcation.
2. Guidewire retracted proximal to the proximal most electrode (spiral deployed). Electrode 1 is now in the main artery.

Figure 3. Device placement within the renal artery

1. Adequate wall contact. All 4 electrode impedance values are stable as shown by an overall linear impedance tracing at all electrodes.
2. Inadequate wall contact. Cyclic, large amplitude tracing is observed on electrode 2 in particular, and on electrode 1. Catheter adjustments are necessary to achieve adequate wall contact.

Figure 4. Making adequate contact with the artery

1. Under fluoroscopic guidance, advance the catheter until the distal electrode is located in the renal artery (Figure 3).
2. Under fluoroscopic guidance, deploy the spiral by retracting the guidewire into the device until the guidewire tip is proximal to electrode 4 (Figure 3, image 2). Make sure the guidewire does not completely exit the rapid exchange port.
3. Adequate wall contact is assessed by the physician and is achieved when the following two conditions are met:
 - a. Deployment of the distal end appears adequate when observed angiographically.
 - b. Impedance values at each electrode are stable through at least one respiratory cycle (Figure 4, image 1).
 - If wall contact does not appear to be adequate per the two criteria above, it is recommended to slightly adjust electrode positions. To do so, slightly torque the catheter clockwise and/or slightly move the catheter forward. These small maneuvers should improve electrode apposition against the vessel wall.
 - If these small adjustments do not improve wall contact, reinsert the guidewire in the distal end of the catheter and change the device location in the artery.
4. If an electrode is not located within the renal artery, or if any electrode deploys in an unsuitable location (such as the ostium of a small vessel or an adrenal gland feeder), deselect (turn off) these electrodes by pressing the electrode number button on the remote control or on the generator touch screen. By deselecting these individual electrodes, RF energy will not be delivered to these electrodes when RF is activated.

Note: Deselection must happen when all electrodes are outside the guide catheter and are displaying impedance values.
5. If desired, for annotation purposes, the left or right kidney can be selected for the treatment by pressing the icons on the generator touch screen or by depressing the kidney button on the remote control. Pressing the button on the remote control will alternate between the left and right kidney selection.

9.5. Performing ablation procedure

1. Once electrodes are well apposed angiographically and impedance values and tracings are stable, RF energy can be delivered to the treatment site. This is done by pressing any of the following: the RF button on the generator front panel, the RF button on the remote control, or the foot switch. The generator delivers power for a target duration of 60 seconds using an automated algorithm and will cease power delivery upon completion of the treatment after 60 seconds. The timer begins counting up and the LED indicator remains blue while RF therapy is

- being delivered. At any point during the procedure, delivery of RF energy can be stopped by pressing the RF button on the generator front panel, pressing the RF button on the remote control, or depressing the foot switch.
- Note:** If the ablation does not initiate due to high-impedance values, first check the catheter position, then check the contact of the dispersive electrode, and finally try moving the dispersive electrode to the patient's flank.
2. If the generator stops delivering RF energy to one or more electrode(s) prior to reaching the 60-second treatment duration, an additional RF application may be performed from the electrode(s) that did not complete treatment, at the same location. First image the artery to ensure that it is safe to perform an ablation. Using the touchscreen, deselect electrodes that completed a 60-second cycle. If needed, perform a slight adjustment to the catheter to ensure proper wall contact, then initiate ablation again.
- Note:** The generator may automatically stop delivering RF energy if certain conditions are detected. A system indicator message or code will appear on the display (see the generator user manual). In the case of a hardware fault condition, the generator will activate a red LED indicator light, emit an audio alert, and display a fault code, if applicable (see the generator user manual for more information about indicator messages and codes).
3. If multiple treatments are to be performed in one artery, move the catheter proximally by pulling it back while taking care to avoid diseased or calcified areas of the vessel. A slight clockwise rotation while pulling back can be applied to ease the motion. All treatments should be located at least 5 mm proximal to any prior treatment location.
4. Once the treatment is completed on one side, take an image of the artery, then advance the guidewire carefully out the tip of the catheter to straighten the spiral distal end.
5. Retract the straightened catheter into the guide catheter. If treating another vessel, reposition the guide catheter within the next vessel. Repeat the procedure for positioning the catheter and delivering treatments.
- If excessive resistance is felt between the guide catheter and electrodes while retracting, consider adjusting the guide catheter position in the vessel to align the catheter coaxially with the guide catheter tip.
 - Ensure that the guide catheter is flushed with heparinized saline periodically, or, at a minimum, between each treatment. Whenever flushing the guide catheter, wait at least 3 seconds to allow the temperature and impedance measurements to stabilize before initiating the next treatment.

9.6. Postprocedure

1. Upon completion of all treatments, straighten the distal end by advancing the guidewire, and then withdraw the straightened catheter completely from the guide catheter.
2. Retract the guidewire and guide catheter from the sheath at the same time.
3. Remove the introducer sheath from the artery and use standard of care procedures to achieve hemostasis at the puncture site.
4. Dispose of the devices in accordance with local hospital, administrative, and/or other government policies.

10. Disclaimer of Warranty

ALTHOUGH THE SYMPLECTICITY SPYRAL™ MULTI-ELECTRODE RENAL DENERVATION CATHETER, HEREINAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC AND ITS AFFILIATES, HEREINAFTER REFERRED TO AS "MEDTRONIC," HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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2020JUL14 at 9:43 a.m. Doc number: M957336A001 Rev. AA [Ruled]
Printhead j - 8 5 x 11.0 inches

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2020/JUL/14 at 9:43 a.m. Doc number: M957336A001 Rev. AA [Autoc]
Printsize: j - 8.5 x 11.0 inches

master:20180508, editDocMpd:20180827, editDocKaz:20171221, editDocStyle_graph:20180316, editDocStyle:20180308, AnnotatDoc_global:20180215, master_pageset:20180222, pagespecs:20180128,
body_outset:20171025, ImagePageLayout:20171023, table:20180806, characterAut:20180903, 3of9:20170211, 3of9_avg:20170211

2020/JUL/14 at 9:43 a.m. Doc number: M967336A001 Rev. AA [ifudoc]
Printspec j - 8.5 x 11.0 inches

master:20180508, alldocMod:20180927, alldocMod2:20171221, alldocStyle_graph:20180316, alldocStyle:20180308, templates_global:20180215, master_pageset:20180222, pagespecs:20180129,
body_extnt:20171025, languageLookup:20171025, table:20180606, characterJoin:20180903, 3of9:20110211, 3of9_svg:20110211

2020JUL14 at 9:13 a.m. Doc Number: M957358AU01 Rev: AA [Hidden]
Printsize | - 8.5 x 11.0 inches

header-20180508_miscdata02-20180527_miscdata02-20171231_abbrevSys_green-20180316_abbrevSys-20180208_tomatoes_gfshad-20180213_header_request-20180222_header-20180129-
body_saint-20171025_honeygate.soup-20171025_mba-20180504_character-20180503_3d6-20110211_3d6_img-20110211

master:20180508, alldocMod:20180927, alldocMod:20171025, languageLookUp:20171025, lnhb:20180608, characterJoin:20180903, 3of9:20110211, 3of9.svg:20110211
body, extent:20171025, lnhb:20180608, characterJoin:20180903, 3of9:20110211, 3of9.svg:20110211
master:20180508, alldocMod:20180927, alldocMod:20171025, languageLookUp:20171025, lnhb:20180608, characterJoin:20180903, 3of9:20110211, 3of9.svg:20110211
body, extent:20171025, lnhb:20180608, characterJoin:20180903, 3of9:20110211, 3of9.svg:20110211

2020/JUL/14 at 9:43 a.m. Doc number: M967336A001 Rev. AA [rudo]
Printspec] - 8.5 x 11.0 inches

2020.06.14 at 9:43 a.m. Doc number: 16867306A001 Rev: AA [ru:doc]
Printed on: j - 8.5 x 11.0 inches

Printer: 20180508_ahs-cs-hd-20180327_ahs-cs-hd-20171221_ahs-cs-hd_20180316_ahs-cs-hd_20180308_ahs-cs-hd_20180215_main_page-set-20180323_page-set-20180128_body_main-201711025_hungary-set-20171025_mah-20180304_czech-rep-hd-20180303_mah-20180211_3rd_page-20110211

2020/JUL/14 at 9:43 a.m. Doc number: M987336/AJ001 Rev AA [Rusdoc]
PrintSpec: 7.85 x 11.0 inches

master: 20180508_#docMod2 20180927_#docMod2 20171221_#docStyle_#trunc 20180316_#docStyle 20180308_#template_global 20180215_master_#pageSet 20180227_#pageSpecs 20180129_#body_#extend 20171025_#languageSet.docx 20171025_#table 20180408_#characterSet 20180903_#obj 20110211_#obj_#svg 20110211

master:20180508, alldocMod:20180927, alldocMod2:20171221, alldocStyle_graph:20180316, alldocStyle:20180306, templates_global:20180215, master_pageset:20180222, pagespecs:20180129, body_extent:20171025, language_lookup:20171025, table:20180606, characterJoin:20180903, 3dF:20110211, 3dF_svg:20110211

2020/JUL/14 at 9:43 a.m. Doc number: M95736A001 Rev. AA [rtdoc]
Printspec] - 8.5 x 11.0 inches

2020JUL14 at 0:43 a.m. Doc number: 1857335AUD1 Raw AA [pdfs]
Printemps j - 6.5 x 11.0 inches

reader_20180508, attached_20180527, attached_20171221, attachStyle_green 20180716, attachStyle_20180208, attached_style_20180215, reader_style_20180222, images_20180128,
book_style_20171025, images_20171025, table_20180008, chessboard_20180007, book_20110211, book_20110211

2020/JUL/14 at 9:43 a.m. Doc number: M957336A001 Rev. AA [ftudoc]
Printspec] - 8.5 x 11.0 inches

master:20180508, alldocMod:20180927, alldocMod2:20171221, alldocStyle_graph:20180316, alldocStyle:20180308, templates_global:20180215, master_pageset:20180222, pagespacs:20180129,
body_extnt:20171025, linguisticLookUp:20171025, lnbk:20180606, characterJoin:20180903, 3of9:20110211, 3of9_svg:20110211

2020/JUL/14 at 9:43 a.m. Doc number: M967336A001 Rev. AA [fluidoc]
Printspec | - 8.5 x 11.0 inches

master:20180508, alldocMod:20180927, alldocMod2:20171121, alldocStyle_graph:20180316, alldocStyle:20180308, templates_global:20180215, master_pageset:20180222, pagespecs:20180129,
body_extint:20171025, languageLookup:20171025, table:20180606, characterJoin:20180903, 3of9:20110211, 3of9_avg:20110211

2020/JUL/14 at 0:43 a.m. Doc number: M967336AD01 Rev. AA [RuKuc]
Printspec: j - 8.5 x 11.0 inches

header:20180508, altDocMod:20180827, altDocMod2:20171221, altDocStyle_graph:20180216, altDocStyle:20180308, template_global:20180215, master_pageant:20180222, pagespecs:20180128,
body_v:altent:20171025, language:lookup:20171026, table:20180806, checkstyle:con:20180803, 3aB:20110211, 3aB_avg:20110211

2020JUL/14 at 9:43 a.m. Doc number: 8667336A001 Rev AA (kdsod)
Printspec I - 8.5 x 11.0 inches

master_20180508, attached_20180527, attached_20171121, addc_style_green_20180204, addc_style_20180204, template_global_20180715, master_jagged_20180222, jaggedness_20180128,
body_color_201711025, jaggedness_color_201711025, table_20180608, crumpled_20180603, 3d_20180711, 3d_20180711, 3d_20180711

2020JUL14 at 9:43 a.m. Doc number: M967336A001 Rev. AA [Mudoc]
Printsize: J - 8.5 x 11.6 inches

master:20180508, aRDocMod:20180927, aRDocMod2:20171221, aRDocStyle_graph:20180318, aRDocStyle:20180308, template_global:20180215, master_pageend:20180223, pagespecs:20180128
body_extent:20171025, languageLocale:20171025, isbn:20180608, characterSet:20180802, 3d6:20110214, 3d8_avg:20110211

2026/JUL/14 at 8:43 a.m. Doc number 18967236A001 Raw_AA [hidden]
Printspace) - 8.9 x 11.0 inches

master_20180308, allDocId=20180307, allDocId=20171221, allDocId= print_20180316, allDocId= 20180308, language_gibet=20180215, master_language=20180222, language=20180728,
book_start=20171023, language=20171025, book=20180308, connection=20180303, book=20171011, book_end=20171011

2020/JUL/14 at 8:43 a.m. Doc number: M957356A/001 Rev. AA [fludoc]
PrintSpec] - 8.5 x 11.0 inches

master:20180508_alldocMod:20180927_alldocMod:20171025_alldocMod:20171025_alldocStyle_graph:20180316_alldocStyle:20180308_template_global:20180215_master_pageset:20180222_pagespecs:20180129_body_extnt:20171025_languageLookup:20171025_table:20180606_characterDom:20180903_3of9:20110211_3of9_avg:20110211

2020/JUL/14 at 9:43 a.m. Doc number: M96736A001 Rev. AA [fluidoc]
PrintSpec | - 8.5 x 11.0 inches

master:20180508, alldocMod:20180927, alldocMod2:20171221, alldocStyle_graph:20180316, alldocStyle:20180308, template_global:20180215, master_page:20180222, pagespec:20180129,
body_ext:20171025, language:20171025, table:20180606, character:20180903, 3d:20110211, 3d:20110211

1

[illegible]

2020/11/14 at 8:43 a.m. Doc Number: M057336A001 Rev. A4 (Rusoc)
Printed on: - 8.5 x 11.0 inches

2020.JUL/14 at 9:43 a.m. Doc number: M96738A001 Rev. AA [fudoc]
Printspec j - 8.5 x 11.0 inches

master:20180508_allDocMod:20180927_allDocMod:20171221_allDocStyle_graph:20180316_allDocStyle:20180308_templates_global:20180215_master_page:20180222_pagespec:20180129_body_extent:20171025_languageLookup:20171025_title:20180606_character:20180903_309:20110211_309_avg:20110211

2020JUL14 at 9:43 a.m. Doc number: 508735AU01 Rev. AA [Russia]
Printer: j - 8.5 x 11.0 inches

master:20180308, addtojob:20180327, addtojob:20171221, addtojob:20180316, addtojob:20180308, template:global:20180315, master:global:20180327, pagesize:20180328,
body:master:20171025, language:English:20171025, date:20180409, checksum:20180903, info:20180111, job_img:20180111

[illegible]

2020.04.14 at 9:43 a.m. Doc number: 6097336A07 Rev. AA (Public)
Printed on: 4.5 x 11.0 inches

2020/JUL/14 at 9:43 a.m. Doc number: M967336A001 Rev. AA (Auto)
Printed on: 8.5 x 11.0 inches

master:20180808, allDocMod:20180827, allDocMod2:20171221, allDocStyle_graph:20180316, allDocStyle:20180306, templates_global:20180213, master_page:20180222, pagespecs:20180129,
body_style:20171025, languageLayout:20171025, table:20180806, characterJoin:20180903, 3d:20110211, 3d_svg:20110211

20204JUL1/14 at 8:43 a.m. Doc number: M967538A001 Rev. AA [Russian]
PrintDate: j - 8.5 x 11.0 inches

reader: 20180508 addDocId: 20171231, addDocStyle: 20180306, template: jrcid: 20180215, reader_jrcid: 20180222, pageDate: 20180728
body: content: 20171026, language: 20171026, label: 20180008, characterSet: 20180008, 2019 20110211, 2019_20110211

2020JUL/14 at 9:43 a.m. Doc number: M957356A001 Rev. AA [fluted]
Printspec | - 8.5 x 11.0 inches

master:20180508_allDocMod:20180927_allDocMod:20171221_allDocStyle_graph:20180316_allDocStyle:20180308_templates_global:20180215_master_pageset:20180222_pagespecs:20180129_body_extnt:20171025_languageLookup:20171025_table:20180606_characterJoin:20180903_3of9:201710211_3of9_svg:201710211

2020/JUL/14 at 9:43 a.m. Doc number: MB67336A001 Rev. AA [fildoc]
Printspec.j - 8.5 x 11.0 inches

master_20180508, alldocMod_20180927, alldocMod_20171121, alldocStyle_graph_20180316, alldocStyle_20180308 templates_global_20180215, master_pageset_20180222, pagespecs_20180129
body_extnt_20171025, languageLookup_20171025, table_20180606, characterJoin_20180903, 3of9_20171021, 3of9_svg_20171021

[illegible]

2020/JUL/14 at 9:43 a.m. Doc Number 1697396A01 Ray AA (AutoDoc)
 Postage: 1 - \$5 x 11.0 inches

Medtronic

Manufacturer:
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com

Medtronic Ireland
Parkmore Business Park West
Galway, Ireland

Authorized Representative in the European Community
Medtronic Ireland

Parkmore Business Park West
Galway, Ireland
Tel. +353-91708800
Fax +353-91707324

Australia Sponsor
Australia
Medtronic Australia Pty Ltd
87 Winslow Road, North Ryde,
NSW 2113 Australia

"MB07336A001"

© 2020 Medtronic
MB07336A001 Rev. AA
2020-07-14

Medtronic

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879

Appendix to the Instructions for Use M967336A001
Simplicity Spyral Multi-Electrode Renal Denervation Catheter

Appendix to the Instructions for Use Simplicity Spyral Multi-Electrode Renal Denervation Catheter

PRESCRIPTION OF MEDICAL DEVICE

Intended to deliver low-level radio frequency (RF) energy through the wall of the renal artery to denervate the human kidney.

INFORMATION ABOUT STICKER

Information on Russian Sticker:

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Warning and informing symbols and inscriptions (if applicable)

Sample of the sticker:

Name of medical device	
Manufacturer:	Code UPN
Country of origin:	Code CFN
Manufacturer's authorized representative::	
Regestretion Certificate № _____ dd _____	
Date of manuf./Lot/UBD: See on pack	See IFU See packaging symbols

SHELF LIFE

3 years

**REQUIREMENTS FOR THE PROTECTION OF THE ENVIRONMENT IN THE
APPLICATION OF MEDICAL DEVICES. PROCEDURE FOR
UTILIZATION/LIQUIDATION OF MD**

After use, this product may be a potential biohazard (hazard class B). Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional and national laws and regulations.

**THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS,
WHICH MATCHES MEDICAL DEVICE**

EN ISO 13485, EN ISO 14971, EN IEC 60601-1, EN IEC 60601-2-2, EN IEC 60601-1-2, EN IEC 62366-1, ISO 10555-1, ISO 2233, ISO 11607-1, ISO 11607-2, IEC 60417, EN ISO 15223-1, ISO 7010, EN 556-1, EN ISO 11137-1, EN ISO 11137-2, EN ISO 11737-1, EN ISO 11737-2, EN ISO 10993-1.

CUSTOMER SERVICE

Contact for all questions, including complaints, repair, and service of medical devices the
Authorized Representative in Russia:

Medtronic, LLC,

123112, Moscow, in.ter.c. Presnensky Municipal District, Presnenskaya nab., 10, floor 9, office III,
room 41, Russian Federation,

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