INSTRUCTIONS FOR USE FOR BELOTERO® SOFT LIDOCAINE MERZ AESTHETICS

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Geneva, the 11th September 2023

1. Description

BELOTERO Soft Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer.

BELOTERO Soft Lidocaine contains 0.3 % of lidocaine hydrochloride.

2. Presentation

BELOTERO Soft Lidocaine is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition:

Cross-linked sodium hyaluronate: 20 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: 1 ml

4. Intended Use/Indications

4.1. Intended Use

BELOTERO Soft Lidocaine is an injectable biodegradable implant intended for enhancement of skin tissue to fill skin depression.

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4.2. Indications

BELOTERO Soft Lidocaine is indicated for injection into the superficial to mid-dermis for treatment of perioral fine lines.

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BELOTERO Soft Lidocaine is designed to be injected into the superficial to mid-dermis by authorized practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection in order to minimize the risk of potential complications.

BELOTERO Soft Lidocaine can be injected with the blanching technique.

Inject BELOTERO Soft Lidocaine slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles.

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BELOTERO Soft Lidocaine must be injected under appropriate aseptic conditions into healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.

To ensure optimal use of BELOTERO Soft Lidocaine, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and / or leakage of the gel at the Luer-lock connection during injection.

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6. Contra-indications

BELOTERO Soft Lidocaine is contra-indicated:

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Do not inject BELOTERO Soft Lidocaine into blood vessels.

Do not inject BELOTERO Soft Lidocaine into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.

Do not inject BELOTERO Soft Lidocaine into an area previously treated with a permanent dermal filler.

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In the absence of available clinical data on tolerance of the injection of BELOTERO Soft Lidocaine in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Soft Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject BELOTERO Soft Lidocaine in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

BELOTERO Soft Lidocaine can be used in combination with other Belotero® products during the same session but in different facial areas. Instructions for use of each product should be followed.

No clinical data is available on the injection of BELOTERO Soft Lidocaine into patient with a Fitzpatrick skin type V/VI.

BELOTERO Soft Lidocaine can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as not only benefits but also adverse events can be cumulative and causality of adverse events could

become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of BELOTERO Soft Lidocaine into an area already treated with other aesthetic products or procedures.

BELOTERO Soft Lidocaine must not be used in association with other aesthetic techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Soft Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Soft Lidocaine with the above-mentioned treatments.

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Injection of BELOTERO Soft Lidocaine into patients with a history of previous herpetic eruption may be associated with herpes reactivation and HHV related diseases (e.g. pityriasis rosea).

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Soft Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

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Do not transfer BELOTERO Soft Lidocaine into another container and do not add other substances to the product.

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8. Warnings

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Soft Lidocaine does not come into contact with such substances.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication, vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.
- As with all dermal filler procedures, BELOTERO Soft Lidocaine should not be used in vascular rich areas. Use in these areas, such as glabella, nose, temples and periorbital regions, has resulted in cases of vascular embolization and symptoms consistent with ocular vessel occlusion.

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• Adverse events:

In occasional cases, one or more of the following may occur in conjunction with the use of products of the Belotero portfolio either immediately or as a delayed reaction: acne cystic, milia, skin dryness (rough facial skin, exfoliation), injection site erosion, inflammation, shivering, fatigue, lymphatic system disorder, rash, burning sensation, injection site sore/warmth/fever, pruritus/itching, urticaria, hematoma, telangiectasia, ecchymosis, edema (including lymph edema), headache/cephalgia, tumefaction, tension, swelling (including persistent swelling), hyper- or hypo-pigmentation, angioedema, induration, blister, vesicle, papule, lump/ bump (visible and/or palpable material) or nodule (including inflammatory nodules), mass, granuloma (including inflammatory signs and foreign body

reactions), necrosis, ischemia, vascular occlusion, embolization, infarction, Tyndall effect (including translucent chords), hypersensitivity, allergic reactions (including asthma attack, Quincke's edema, anaphylactic shock or throat tightening) to one of the product's components (e.g. hyaluronic acid, BDDE, lidocaine hydrochloride), oral and dental disorders, nervous system impairment, impairment of the otorhinolaryngological system (e.g. nasal congestion, oropharyngeal pain, dysguesia, rhinorrhea, epistaxis, sinusitis, transient hearing loss), mastication pain, parotid gland enlargement, muscle twitching, muscle injury/disorder, nausea, vomiting, circulatory collapse, presyncope, peripheral venous disease, hot flush, anxiety caused by trypanophobia, patient dissatisfaction and disappointment (due to lack of or reduced performance, decreased firmness/response, undesirable aesthetic effect), injection site discharge, device migration, product distribution issue (e.g. product accumulation), injection site indentation, superficial vein prominence, overcorrection or cranial nerve disorder (e.g. cranial nerve paralysis, facial paralysis, trigeminal neuralgia).

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10. Assembly of needle to syringe

For optimal use of BELOTERO Soft Lidocaine it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. **Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor** between the thumb and the index fingers.

2. Grasp the protective cap with the other hand and unscrew it.
3. **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

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
Terms of warranty:

Anteis SA, 18, Chemin des Aulx/1228 Plan-les-Ouates, Geneva, Switzerland, as the manufacturer of Belotero Soft Lidocaine guarantees the compliance of the product with all parameters and properties specified in this document, subject to the conditions of transportation, storage, and use in accordance with the instructions for use for two years.

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


Manufacturer of BELOTERO Soft Lidocaine

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INSTRUCTIONS FOR USE FOR BELOTERO® BALANCE LIDOCAINE MERZ AESTHETICS

30.08.2023

Belotero line, viscoelastic implant based on cross-linked sodium hyaluronate

for correction of depressions of tissue, in the following versions:

1) Belotero SOFT Lidocaine, consisting of:

- implantable gel, in a 1.0 ml syringe - 1 pc.;
- sterile needle 30G ½ (0,3 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use - 1 pc.

2) Belotero BALANCE:

Belotero BALANCE Lidocaine;

Belotero LIPS Contour,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) and 30G ½ (0,3 x 13) - 1 pc. each;
- label - 2 pcs.;
- Instructions for use - 1 pc.

3) Belotero INTENSE:

Belotero INTENSE Lidocaine;

Belotero LIPS Shape,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use - 1 pc.

Seen exclusively for the legalization of
the above signature of Mrs Anne-Sophie
Dreano.

*The undersigned notary assumes no
responsibility for the validity and
content of the present document.*

Geneva, the 11th September 2023

1. Description

BELOTERO Balance Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer.

BELOTERO Balance Lidocaine contains 0.3 % of lidocaine hydrochloride.

2. Presentation

BELOTERO Balance Lidocaine is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition:

Cross-linked sodium hyaluronate: 22.5 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: 1 ml

4. Intended Use / Indications

4.1. Intended Use

BELOTERO Balance Lidocaine is an injectable biodegradable implant intended for enhancement of skin tissue to fill skin depression.

The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

4.2. Indications

BELOTERO Balance Lidocaine is indicated for injection into the superficial to mid-dermis for treatment of naso-labial folds, marionette lines, perioral lines and moderate oral commissures. BELOTERO Balance Lidocaine is indicated for lip enhancement.

5. Posology and administration method

BELOTERO Balance Lidocaine is designed to be injected into the superficial to mid-dermis by authorized practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy (explicitly about the distribution pattern of labial arteries) at and around the site of injection in order to minimize the risk of potential complications.

BELOTERO Balance Lidocaine can be injected with the blanching technique. Submucosal or subcutaneous injection is recommended for lip enhancement. In the tear trough area, inject deep on the bone below the orbicularis oculi muscle.

Inject BELOTERO Balance Lidocaine slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles.

The risk of an intravascular injection can be reduced by different strategies, including aspiration prior to injection, utilizing lower volumes and serial injections in high-risk areas, treating one side at a time, pinching/tenting the skin to provide more space superficial to the branches of the main arteries, and manual occlusion of the origin of the supratrochlear vessels with the non-dominant finger. Blunt cannulas may reduce, but not eliminate the risk.

General recommended injection techniques are for example: linear or serial threading, fanning, cross-hatching or serial (micro)puncture. The quantity of product to be injected depends on the area to be corrected.

BELOTERO Balance Lidocaine must be injected under appropriate aseptic conditions into healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.

To ensure optimal use of BELOTERO Balance Lidocaine, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and / or leakage of the gel at the Luer-lock connection during injection.

If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.

The practitioner should note that the use of 30G½" needle requires greater force to inject the gel than 27G½" needle.

The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.

The graduations on the syringe label are only intended for orientation for the user.

Gently massage the treated area after the injection to distribute the product uniformly.

Before treatment, the patient's suitability for the treatment and the patient's need for pain relief (topical anaesthetics, ice packs, distraction techniques, local anaesthetic injections, or nerve blocks depending on the injection site and size of needle used), should be assessed.

6. Contra-indications

BELOTERO Balance Lidocaine is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride, BDDE or to amide-type local anaesthetics,
- In pregnant and breast-feeding women,
- In patients under 18 years old,
- In patients presenting a general infection,
- In patients presenting an active auto-immune disease.

Do not inject BELOTERO Balance Lidocaine into blood vessels.

Do not inject BELOTERO Balance Lidocaine into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.

Do not inject BELOTERO Balance Lidocaine into an area previously treated with a permanent dermal filler.

Do not inject BELOTERO Balance Lidocaine into the glabellar or nose region.

7. Precautions for use

Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.

In the absence of available clinical data on tolerance of the injection of BELOTERO Balance Lidocaine in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Balance Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject BELOTERO Balance Lidocaine in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

BELOTERO Balance Lidocaine injected in the forehead (excluding glabellar) area may be associated with an increased risk for intravascular complications and the consequences of local vascular occlusion, embolization, vision impairment, blindness ischemia, necrosis or infarction.

BELOTERO Balance Lidocaine injected in the periorbital region may be associated with an increased frequency and severity of side effects and adverse events.

BELOTERO Balance Lidocaine can be used in combination with other Belotero® products during the same session but in different facial areas. Instructions for use of each product should be followed.

Limited clinical data is available on the injection of BELOTERO Balance Lidocaine into patients with a Fitzpatrick skin type V/VI.

BELOTERO Balance Lidocaine can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as not only benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of BELOTERO Balance Lidocaine into an area already treated with other aesthetic products or procedures.

BELOTERO Balance Lidocaine must not be used in association with other aesthetic techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Balance Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Balance Lidocaine with the above-mentioned treatments.

Patients using anti-coagulation, anti-platelet, or thrombolytic medications (e.g. warfarin), anti-inflammatory drugs (oral/injectable corticosteroids or non-steroidal anti-inflammatory drugs (NSAIDs; e.g. aspirin, ibuprofen)), or other substances known to increase coagulation time (vitamins or herbal supplements, e.g. Vitamin E, garlic, Ginkgo biloba and St. John's Wort), from 10 days pre- to 3 days post-injection may have increased reactions of hematomas, nodules or bleeding at the injection site.

Injection of BELOTERO Balance Lidocaine into patients with a history of previous herpetic eruption may be associated with herpes reactivation and HHV related diseases (e.g. pityriasis rosea).

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Balance Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.

It should be noted that the presence of lidocaine may cause local redness, hypersensitivity or transient loco-regional numbness.

For normal healthy adults it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in signs of the central nervous system or cardiovascular toxicity.

When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or

agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.

Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.

Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the inner packaging has been opened or damaged.

Do not transfer BELOTERO Balance Lidocaine into another container and do not add other substances to the product.

Only the gel is sterile, but not the outside of the syringe.

Discard the syringe, the remaining product and the needles in the appropriate container after use.

Do not re-sterilize and do not reuse due to the associated risks including infection.

The patient must avoid applying makeup (including skin care products) for at least 12 hours after treatment as well as avoid saunas, peeling, Turkish baths and prolonged exposure to the sun, UV rays, extreme heat and cold for 2 weeks after the treatment. Patients should also avoid putting pressure on and/ or handling the treated area and should avoid strenuous physical activity following treatment.

The patient must avoid drinking alcohol for 24 hours before and after treatment. Alcohol may cause the blood vessels to dilate and cause more bruising.

8. Warnings

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Balance Lidocaine does not come into contact with such substances.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication, vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

9. Side effects and adverse events

Patients must be informed by the practitioner about possible side effects and adverse events before treatment.

• Side effects:

Injection site reactions may occur following injection into the skin but disappear spontaneously within a few days. This includes swelling, nodule or lump/bump, bruising/purpura, hematoma, ecchymosis, induration, erythema/redness, tenderness, pain, discoloration and pruritus/itching, tingling, paraesthesia, numbness, hypoaesthesia, scabbing, needle mark and discomfort or irritation. These injection site reactions are generally of mild or

moderate intensity. A transient bleeding may also occur at the injection site and usually stops spontaneously as soon as the injection is finished.

• Adverse events:

In occasional cases, one or more of the following may occur in conjunction with the use of products of the Belotero portfolio either immediately or as a delayed reaction: acne cystic, milia, skin dryness (rough facial skin, exfoliation), injection site erosion, inflammation, shivering, fatigue, lymphatic system disorder, rash, burning sensation, injection site sore/warmth/fever, pruritus/itching, urticaria, hematoma, telangiectasia, ecchymosis, edema (including lymph edema), headache/cephalgia, tumefaction, tension, swelling (including persistent swelling), hyper- or hypo-pigmentation, angioedema, induration, blister, vesicle, papule, lump/ bump (visible and/or palpable material) or nodule (including inflammatory nodules), mass, granuloma (including inflammatory signs and foreign body reactions), necrosis, ischemia, vascular occlusion, embolization, infarction, Tyndall effect (including translucent chords), hypersensitivity, allergic reactions (including asthma attack, Quincke's edema, anaphylactic shock or throat tightening) to one of the product's components (e.g. hyaluronic acid, BDDE, lidocaine hydrochloride), oral and dental disorders, nervous system impairment, impairment of the otorhinolaryngological system (e.g. nasal congestion, oropharyngeal pain, dysguesia, rhinorrhea, epistaxis, sinusitis, transient hearing loss), mastication pain, parotid gland enlargement, muscle twitching, muscle injury/disorder, nausea, vomiting, circulatory collapse, presyncope, peripheral venous disease, hot flush, anxiety caused by trypanophobia, patient dissatisfaction and disappointment (due to lack of or reduced performance, decreased firmness/response, undesirable aesthetic effect), injection site discharge, device migration, product distribution issue (e.g. product accumulation), injection site indentation, superficial vein prominence, overcorrection or cranial nerve disorder (e.g. cranial nerve paralysis, facial paralysis, trigeminal neuralgia).

Rare cases of the following adverse events have been reported with hyaluronic acid products such as infection (e.g. erysipelas, phlegmon, cellulitis, including open or draining wounds and (dental) abscess, impetigo, pustules), chronic infection (including biofilm formation), scarring, persistent skin discoloration, sensory dysfunction, non-thrombotic lung embolism as well as sarcoid granuloma formation in subjects with hepatitis C and interferon treatment, cerebral injuries (e.g. intracranial penetration, subarachnoid hemorrhage), strabismus, ophthalmoplegia, iris adhesions, cataract, conjunctival hemorrhage, eyelid ptosis and lacrimation.

The risk of granuloma, ischemia, necrosis and vascular occlusion is higher with deep injections and high volumes.

Isolated cases of visual impairment or blindness following unintentional intra-arterial injection have been reported in literature.

Patients should be instructed to report any side effects which last for more than one week and any adverse event as soon as it occurs to his/her practitioner, especially if patient has changes in his/her vision, signs of a stroke (including sudden difficulty speaking, numbness or weakness in his/her face, arms, or legs, difficulty walking, face drooping, severe headache, dizziness, or confusion), white appearance of the skin, or unusual pain during or shortly after treatment. The practitioner may then refer the patient to the appropriate treatment.

Patients with lighter skin types are more likely to develop injection-related adverse events. However, patients with skin of color are more likely to develop post-inflammatory hyperpigmentation and /or hypertrophic scar/keloid formation following injection procedures. Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. itching, swelling, erythema, inflammation.

- There is no known interaction with other local or loco-regional anaesthetics.

10. Assembly of needle to syringe

For optimal use of BELOTERO Balance Lidocaine it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
2. Grasp the protective cap with the other hand and unscrew it.
3. **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

11. Transportation, storage, use


Store between 2°C and 25°C. Protect from light and freezing.

Avoid mechanical shocks.

The manufacturer has no requirements to humidity and atmospheric pressure

12. References

Updated documentation may be available from ANTEIS S.A., Switzerland.

Backstop in the right position during injection

Disposal

In the territory of the Russian Federation used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with SanPiN 2.1.3684-21 „Sanitary and epidemiological requirements to the maintenance of the territory of urban and rural settlements, to waterbodies, to drinking water and water supply of citizens, to atmosphere air, to soils, to living areas, to exploitation of manufacturing and communal areas, to organization and realization of sanitary and anti-epidemic (preventive) events” as belonging to hazard class B.

Requirements for environmental protection

The device should be disposed of in accordance with any applicable local/ national regulations. With proper disposal, there is no harm to the environment.

Terms of warranty:

Anteis SA, 18, Chemin des Aulx/1228 Plan-les-Ouates, Geneva, Switzerland, as the manufacturer of Belotero Balance Lidocaine guarantees the compliance of the product with all parameters and properties specified in this document, subject to the conditions of transportation, storage, and use in accordance with the instructions for use for two years.

Warning fragments

Caution

Consult instructions for use

Do not use if package is damaged.

Single use product. Do not re-use


Open the blister by pulling the blister lid following the arrow.

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.

Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Temperature limit of storage: 2 °C – 25 °C

Date of manufacture

Use-by date


Batch code

CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number.

Manufacturer

Authorized Representative in the European Community

Manufacturer of BELOTERO Balance Lidocaine

Anteis SA

18 Chemin des Aulx
1228 Plan-les-Ouates
Switzerland

European Community Representative:

Merz Aesthetics GmbH
Eckenheimer Landstrasse 100
60318 Frankfurt
Germany

Distributed by:

Authorized Representative in RF

Merz Pharma LLC

3rd street Yamskogo Polya 2 bld.26, 2nd floor, premises VII, room 29;

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30.08.2023

INSTRUCTIONS FOR USE FOR BELOTERO® LIPS CONTOUR MERZ AESTHETICS

Belotero line, viscoelastic implant based on cross-linked sodium hyaluronate
for correction of depressions of tissue, in the following versions:

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Geneva, the 11th September 2023

1. Description

BELOTERO Lips Contour is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer.

BELOTERO Lips Contour contains 0.3 % of lidocaine hydrochloride.

2. Presentation

BELOTERO Lips Contour is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition:

Cross-linked sodium hyaluronate: 22.5 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: 0.6 ml

4. Intended Use/Indications

4.1. Intended Use

BELOTERO Lips Contour is an injectable biodegradable implant intended for enhancement of skin tissue to fill skin depression.

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4.2. Indications

BELOTERO Lips Contour is indicated for injection into the superficial to mid-dermis for treatment of lip contour, perioral lines and moderate oral commissures. BELOTERO Lips Contour is indicated for lip enhancement.

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BELOTERO Lips Contour is designed to be injected into the superficial to mid-dermis by authorized practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy (explicitly about the distribution pattern of labial arteries) at and around the site of injection in order to minimize the risk of potential complications.

BELOTERO Lips Contour can be injected with the blanching technique. Submucosal or subcutaneous injection is recommended for lip enhancement. In the tear trough area, inject deep on the bone below the orbicularis oculi muscle.

Inject BELOTERO Lips Contour slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles.

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BELOTERO Lips Contour must be injected under appropriate aseptic conditions into healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.

To ensure optimal use of BELOTERO Lips Contour, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and / or leakage of the gel at the Luer-lock connection during injection.

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6. Contra-indications

BELOTERO Lips Contour is contra-indicated:

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Do not inject BELOTERO Lips Contour into blood vessels.

Do not inject BELOTERO Lips Contour into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.

Do not inject BELOTERO Lips Contour into an area previously treated with a permanent dermal filler.

Do not inject BELOTERO Lips Contour into the glabellar or nose region.

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In the absence of available clinical data on tolerance of the injection of BELOTERO Lips Contour in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Lips Contour on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject BELOTERO Lips Contour in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

BELOTERO Lips Contour injected in the forehead (excluding glabellar) area may be associated with an increased risk for intravascular complications and the consequences of local vascular occlusion, embolization, vision impairment, blindness ischemia, necrosis or infarction.

BELOTERO Lips Contour injected in the periorbital region may be associated with an increased frequency and severity of side effects and adverse events.

BELOTERO Lips Contour can be used in combination with other Belotero® products during the same session but in different facial areas. Instructions for use of each product should be followed.

Limited clinical data is available on the injection of BELOTERO Lips Contour into patients with a Fitzpatrick skin type V/VI.

BELOTERO Lips Contour can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as not only benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of BELOTERO Lips Contour into an area already treated with other aesthetic products or procedures.

BELOTERO Lips Contour must not be used in association with other aesthetic techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Lips Contour must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Lips Contour with the above-mentioned treatments.

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Injection of BELOTERO Lips Contour into patients with a history of previous herpetic eruption may be associated with herpes reactivation and HHV related diseases (e.g. pityriasis rosea).

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Lips Contour on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

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When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.

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Do not transfer BELOTERO Lips Contour into another container and do not add other substances to the product.

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The patient must avoid applying makeup (including skin care products) for at least 12 hours after treatment as well as avoid saunas, peeling, Turkish baths and prolonged exposure to the sun, UV rays, extreme heat and cold for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area and should avoid strenuous physical activity following treatment.

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8. Warnings

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Lips Contour does not come into contact with such substances.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication, vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

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• Side effects:

Injection site reactions may occur following injection into the skin but disappear spontaneously within a few days. This includes swelling, nodule or lump/bump, bruising/purpura, hematoma, ecchymosis, induration, erythema/redness, tenderness, pain, discoloration and pruritus/itching, tingling, paraesthesia, numbness, hypoaesthesia, scabbing, needle mark and discomfort or irritation. These injection site reactions are generally of mild or moderate intensity. A transient bleeding may also occur at the injection site and usually stops spontaneously as soon as the injection is finished.

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Patients with lighter skin types are more likely to develop injection-related adverse events. However, patients with skin of color are more likely to develop post-inflammatory hyperpigmentation and /or hypertrophic scar/keloid formation following injection procedures.

Patients with specific ethnic characteristics, e.g. Asian population, should be informed of a higher risk of tissue reactions, e.g. itching, swelling, erythema, inflammation.

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10. Assembly of needle to syringe

For optimal use of BELOTERO Lips Contour it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

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

Disposal

In the territory of the Russian Federation used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with SanPiN 2.1.3684-21 „Sanitary and epidemiological requirements to the maintenance of the territory of urban and rural settlements, to waterbodies, to drinking water and water supply of citizens, to atmosphere air, to soils, to living areas, to exploitation of manufacturing and communal areas, to organization and realization of sanitary and anti-epidemic (preventive) events" as belonging to hazard class B.

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The device should be disposed of in accordance with any applicable local/national regulations. With proper disposal, there is no harm to the environment.

Terms of warranty:

Anteis SA, 18, Chemin des Aulx/1228 Plan-les-Ouates, Geneva, Switzerland, as the manufacturer of Belotero Lips Contour guarantees the compliance of the product with all parameters and properties specified in this document, subject to the conditions of transportation, storage, and use in accordance with the instructions for use for two years.

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Manufacturer of BELOTERO Lips Contour

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Geneva, the 11th September 2023

1. Description

BELOTERO Intense Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer.

BELOTERO Intense Lidocaine contains 0.3 % of lidocaine hydrochloride.

2. Presentation

BELOTERO Intense Lidocaine is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition:

Cross-linked sodium hyaluronate: 25.5 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: 1 ml

4. Intended Use / Indications

4.1. Intended Use

BELOTERO Intense Lidocaine is an injectable biodegradable implant intended for enhancement of skin tissue to fill skin depression.

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4.2. Indications

BELOTERO Intense Lidocaine is indicated for injection into the deep dermis for treatment of naso-labial folds and marionette lines. BELOTERO Intense Lidocaine is also indicated for lip enhancement.

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BELOTERO Intense Lidocaine is designed to be injected into the deep dermis by authorized practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection in order to minimize the risk of potential complications. Submucosal or subcutaneous injection is recommended for lip enhancement.

Inject BELOTERO Intense Lidocaine slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles.

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BELOTERO Intense Lidocaine must be injected under appropriate aseptic conditions into healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.

To ensure optimal use of BELOTERO Intense Lidocaine, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and / or leakage of the gel at the Luer-lock connection during injection.

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The use of the enclosed 27G½" needle is recommended, as a smaller needle diameter would require a greater force to inject the implant.

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6. Contra-indications

BELOTERO Intense Lidocaine is contra-indicated:

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Do not inject BELOTERO Intense Lidocaine into blood vessels.

Do not inject BELOTERO Intense Lidocaine into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.

Do not inject BELOTERO Intense Lidocaine into an area previously treated with a permanent dermal filler.

Do not inject BELOTERO Intense Lidocaine into the glabellar or nose region.

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In the absence of available clinical data on tolerance of the injection of BELOTERO Intense Lidocaine in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Intense Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject BELOTERO Intense Lidocaine in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

BELOTERO Intense Lidocaine injected in the temples area may be associated with an increased risk for intravascular complications and the consequences of local vascular occlusion, embolization, vision impairment, blindness, ischemia, necrosis or infarction.

BELOTERO Intense Lidocaine can be used in combination with other Belotero® products during the same session but in different facial areas. Instructions for use of each product

No clinical data is available on the injection of BELOTERO Intense Lidocaine into patient with a Fitzpatrick skin type V/VI.

BELOTERO Intense Lidocaine can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as not only benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of BELOTERO Intense Lidocaine into an area already treated with other aesthetic products or procedures.

BELOTERO Intense Lidocaine must not be used in association with other aesthetic techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Intense Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Intense Lidocaine with the above-mentioned treatments.

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Injection of BELOTERO Intense Lidocaine into patients with a history of previous herpetic eruption may be associated with herpes reactivation and HHV related diseases (e.g. pityriasis rosea).

In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Intense Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

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Do not transfer BELOTERO Intense Lidocaine into another container and do not add other substances to the product.

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8. Warnings

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Intense Lidocaine does not come into contact with such substances.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication, vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

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• Adverse events:

In occasional cases, one or more of the following may occur in conjunction with the use of products of the Belotero portfolio either immediately or as a delayed reaction: acne cystic, milia, skin dryness (rough facial skin, exfoliation), injection site erosion, inflammation,

shivering, fatigue, lymphatic system disorder, rash, burning sensation, injection site sore/warmth/fever, pruritus/itching, urticaria, hematoma, telangiectasia, ecchymosis, edema (including lymph edema), headache/cephalgia, tumefaction, tension, swelling (including persistent swelling), hyper- or hypo-pigmentation, angioedema, induration, blister, vesicle, papule, lump/ bump (visible and/or palpable material) or nodule (including inflammatory nodules), mass, granuloma (including inflammatory signs and foreign body reactions), necrosis, ischemia, vascular occlusion, embolization, infarction, Tyndall effect (including translucent chords), hypersensitivity, allergic reactions (including asthma attack, Quincke's edema, anaphylactic shock or throat tightening) to one of the product's components (e.g. hyaluronic acid, BDDE, lidocaine hydrochloride), oral and dental disorders, nervous system impairment, impairment of the otorhinolaryngological system (e.g. nasal congestion, oropharyngeal pain, dysguesia, rhinorrhea, epistaxis, sinusitis, transient hearing loss), mastication pain, parotid gland enlargement, muscle twitching, muscle injury/disorder, nausea, vomiting, circulatory collapse, presyncope, peripheral venous disease, hot flush, anxiety caused by trypanophobia, patient dissatisfaction and disappointment (due to lack of or reduced performance, decreased firmness/response, undesirable aesthetic effect), injection site discharge, device migration, product distribution issue (e.g. product accumulation), injection site indentation, superficial vein prominence, overcorrection or cranial nerve disorder (e.g. cranial nerve paralysis, facial paralysis, trigeminal neuralgia).

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10. Assembly of needle to syringe

For optimal use of BELOTERO Intense Lidocaine it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. Firmly hold the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
2. Grasp the protective cap with the other hand and unscrew it.
3. Push & Twist the needle on the syringe until a resistance is felt. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

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
Disposal

In the territory of the Russian Federation used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with SanPIN 2.1.3684-21 „Sanitary and epidemiological requirements to the maintenance of the territory of urban and

rural settlements, to waterbodies, to drinking water and water supply of citizens, to atmosphere air, to soils, to living areas, to exploitation of manufacturing and communal areas, to organization and realization of sanitary and anti-epidemic (preventive) events" as belonging to hazard class B.

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Terms of warranty:

Anteis SA, 18, Chemin des Aulx/1228 Plan-les-Ouates, Geneva, Switzerland, as the manufacturer of Belotero Intense Lidocaine guarantees the compliance of the product with all parameters and properties specified in this document, subject to the conditions of transportation, storage, and use in accordance with the instructions for use for two years.

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


Manufacturer of BELOTERO Intense Lidocaine

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(Signature)

30.08.2023

INSTRUCTIONS FOR USE FOR BELOTERO® LIPS SHAPE MERZ AESTHETICS

Belotero line, viscoelastic implant based on cross-linked sodium hyaluronate for correction of depressions of tissue, in the following versions:

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Seen exclusively for the legalization of the above signature of Mrs Anne-Sophie Dreano.

The undersigned notary assumes no responsibility for the validity and content of the present document.

Geneva, the 11th September 2023

1. Description

BELOTERO Lips Shape is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer.

BELOTERO Lips Shape contains 0.3 % of lidocaine hydrochloride.

2. Presentation

BELOTERO Lips Shape is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

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4. Intended Use/Indications

4.1. Intended Use

BELOTERO Lips Shape is an injectable biodegradable implant intended for enhancement of skin tissue to fill skin depression.

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4.2. Indications

BELOTERO Lips Shape is indicated for submucosal or subcutaneous injection into the lips. BELOTERO Lips Shape is indicated for injection into the deep dermis for treatment of perioral lines and severe oral commissures.

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BELOTERO Lips Shape is designed to be injected into the deep dermis by authorized practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy (explicitly about the distribution pattern of labial arteries) at and around the site of injection in order to minimize the risk of potential complications. Submucosal or subcutaneous injection is recommended for lip enhancement.

Inject BELOTERO Lips Shape slowly and not too fast to apply the least amount of pressure necessary, according to the appropriate injection technique using the provided needles.

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BELOTERO Lips Shape must be injected under appropriate aseptic conditions into healthy, non-inflamed skin. Before injection, thoroughly disinfect the area to be treated.

To ensure optimal use of BELOTERO Lips Shape, it is recommended to assemble the needle according to the diagrams below. Improper assembly may lead to a separation of the needle and syringe and / or leakage of the gel at the Luer-lock connection during injection. If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.

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6. Contra-indications

BELOTERO Lips Shape is contra-indicated:

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Do not inject BELOTERO Lips Shape into blood vessels.

Do not inject BELOTERO Lips Shape into skin areas presenting active cutaneous inflammation or infection due to e.g. immunological, allergic, bacterial, fungal or viral causes.

Do not inject BELOTERO Lips Shape into an area previously treated with a permanent dermal filler.

Do not inject BELOTERO Lips Shape into the glabellar or nose region.

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In the absence of available clinical data on tolerance of the injection of BELOTERO Lips Shape in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO Lips Shape on a case-by-case basis depending on the nature of the disease as well as the associated treatment as it may worsen the existing patient health condition. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

It is recommended not to inject BELOTERO Lips Shape in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or keloids.

BELOTERO Lips Shape injected in the temples area may be associated with an increased risk for intravascular complications and the consequences of local vascular occlusion, embolization, vision impairment, blindness, ischemia, necrosis or infarction.

BELOTERO Lips Shape can be used in combination with other Belotero® products during the same session but in different facial areas. Instructions for use of each product should be followed.

No clinical data is available on the injection of BELOTERO Lips Shape into patient with a Fitzpatrick skin type V/VI.

BELOTERO Lips Shape can be used in combination treatments such as with botulinum toxin and/or calcium hydroxylapatite (Radiesse®) only if injected in different facial areas. Practitioners should be experienced and patients appropriately selected as not only benefits but also adverse events can be cumulative and causality of adverse events could become difficult to determine. Instructions for use, depth of injection and appropriate recommendation of each product should be followed. No clinical data are available on the injection of BELOTERO Lips Shape into an area already treated with other aesthetic products or procedures.

BELOTERO Lips Shape must not be used in association with other aesthetic techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO Lips Shape must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO Lips Shape with the above-mentioned treatments.

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Injection of BELOTERO Lips Shape into patients with a history of previous herpetic eruption may be associated with herpes reactivation and HHV related diseases (e.g. pityriasis rosea). In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO Lips Shape on a case-by-case basis depending of the nature of the disease as well as the associated treatment.

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Do not transfer BELOTERO Lips Shape into another container and do not add other substances to the product.

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8. Warnings

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO Lips Shape does not come into contact with such substances.
- Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vascular complication, vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures. Practitioners should immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

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• Adverse events:

In occasional cases, one or more of the following may occur in conjunction with the use of products of the Belotero portfolio either immediately or as a delayed reaction: acne cystic, milia, skin dryness (rough facial skin, exfoliation), injection site erosion, inflammation, shivering, fatigue, lymphatic system disorder, rash, burning sensation, injection site sore/warmth/fever, pruritus/itching, urticaria, hematoma, telangiectasia, ecchymosis, edema (including lymph edema), headache/cephalgia, tumefaction, tension, swelling

(including persistent swelling), hyper- or hypo-pigmentation, angioedema, induration, blister, vesicle, papule, lump/ bump (visible and/or palpable material) or nodule (including inflammatory nodules), mass, granuloma (including inflammatory signs and foreign body reactions), necrosis, ischemia, vascular occlusion, embolization, infarction, Tyndall effect (including translucent chords), hypersensitivity, allergic reactions (including asthma attack, Quincke's edema, anaphylactic shock or throat tightening) to one of the product's components (e.g. hyaluronic acid, BDDE, lidocaine hydrochloride), oral and dental disorders, nervous system impairment, impairment of the otorhinolaryngological system (e.g. nasal congestion, oropharyngeal pain, dysphagia, rhinorrhea, epistaxis, sinusitis, transient hearing loss), mastication pain, parotid gland enlargement, muscle twitching, muscle injury/disorder, nausea, vomiting, circulatory collapse, presyncope, peripheral venous disease, hot flush, anxiety caused by trypanophobia, patient dissatisfaction and disappointment (due to lack of or reduced performance, decreased firmness/response, undesirable aesthetic effect), injection site discharge, device migration, product distribution issue (e.g. product accumulation), injection site indentation, superficial vein prominence, overcorrection or cranial nerve disorder (e.g. cranial nerve paralysis, facial paralysis, trigeminal neuralgia).

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10. Assembly of needle to syringe

For optimal use of BELOTERO Lips Shape it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

1. **Firmly hold** the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.
2. Grasp the protective cap with the other hand and unscrew it.
3. **Push & Twist** the needle on the syringe until a resistance is felt. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
4. Keep holding the Luer-lock and remove the sheath from the needle.

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
Avoid mechanical shocks.

The manufacturer has no requirements to humidity and atmospheric pressure.

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Disposal

In the territory of the Russian Federation used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with SanPiN 2.1.3684-21 „Sanitary and epidemiological requirements to the maintenance of the territory of urban and

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Terms of warranty:

Anteis SA, 18, Chemin des Aulx/1228 Plan-les-Ouates, Geneva, Switzerland, as the manufacturer of Belotero Lips Shape guarantees the compliance of the product with all parameters and properties specified in this document, subject to the conditions of transportation, storage, and use in accordance with the instructions for use for two years.

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
Manufacturer of BELOTERO Lips Shape

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Authorized Representative in RF

Merz Pharma LLC

3rd street Yamskogo Polya 2 bld.26, 2nd floor, premises VII, room 29;
125124, Moscow

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