Anteis S.A.
18 Chemin des Aulx
CH-1228 Plan-les-Ouates, Geneva
Switzerland

To whom it may concern:

Medical Device:

Belotero line viscoelastic implant based on cross-linked hyaluronate for correction of depressions of tissue

We, Anteis S.A., a company organized and existing under laws of Switzerland, situated at Chemin des Aulx 18, CH-1228 Plan-les-Ouates, Geneva, Switzerland, in person of Anne-Sophie DREANO, Regulatory Affairs Expert, hereby declare that the attached documents, Instructions for Use for Belotero line viscoelastic implant based on cross-linked hyaluronate for correction of depressions of tissue in the following versions:

- 1) Belotero SOFT Lidocaine, consisting of:
 - implantable gel, in a 1.0 ml syringe - 1 pc.;
 - sterile needle 30G ½ (0.3 x 13) - 2 pcs.;
 - label - 2 pcs.;
 - Instructions for use – 1 pc.
- 2) Belotero BALANCE Lidocaine:
Belotero BALANCE Lidocaine;
Belotero LIPS Contour, consisting of:
 - implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
 - sterile needle 27G ½ (0.4 x 13) and 30G ½ (0.3 x 13) - 1 pc. each;
 - label - 2 pcs.;
 - Instructions for use – 1 pc.
- 3) Belotero INTENSE Lidocaine:
Belotero INTENSE Lidocaine;
Belotero LIPS Shape, consisting of:
 - implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
 - sterile needle 27G ½ (0.4 x 13) - 2 pcs.;
 - label - 2 pcs.;
 - Instructions for use – 1 pc.,

are original and confidential.

Enclosed IFUs:

Product	File name	Number of pages
Belotero Soft Lidocaine	IfU_Belotero Soft Lidocaine_v4_RU_190107	8 pages
Belotero Balance Lidocaine	IfU_Belotero Balance Lidocaine_v3_RU_190107	8 pages
Belotero Lips Contour	IfU_Belotero Lips Contour_v2_RU_190107	8 pages
Belotero Intense Lidocaine	IfU_Belotero Intense Lidocaine_v3_RU_190107	9 pages

ante

Instructions

BELOTERO

Belotero Lips Shape

IfU_Belotero Lips Shape_v2_RU_190107

8 pages

Signature:

Date: 07 January 2019

DREANO Anne-Sophie
Anteis Regulatory Affairs Expert

Seen exclusively for the legalization of the
beside signature of Mrs Anne-Sophie Dreano.

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MERZ AESTHETICS

INSTRUCTIONS FOR USE FOR BELOTERO® SOFT LIDOCAINE

Belotero line viscoelastic implant based on cross-linked sodium hyaluronate for correction of depressions of tissue, in ~~the~~ following versions:

- 1) Belotero SOFT Lidocaine, consisting of:
 - implantable gel, in a 1.0 ml syringe - 1 pc.;
 - sterile needle 30G ½ (0.3 x 13) - 2 pcs.;
 - label - 2 pcs.;
 - Instructions for use – 1 pc.

- 2) Belotero BALANCE Lidocaine:
 Belotero BALANCE Lidocaine;
 Belotero LIPS Contour,
 consisting of:
 - Implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
 - sterile needle 27G ½ (0.4 x 13) and 30G ½ (0.3 x 13) - 1 pc. each;
 - label - 2 pcs.;
 - Instructions for use – 1 pc.

- 3) Belotero INTENSE Lidocaine:
 Belotero INTENSE Lidocaine;
 Belotero LIPS Shape,
 consisting of:
 - implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
 - sterile needle 27G ½ (0.4 x 13) - 2 pcs.;
 - label - 2 pcs.;
 - Instructions for use – 1 pc.

Intended Use

Belotero line in the following versions: Belotero SOFT Lidocaine, Belotero BALANCE Lidocaine, Belotero LIPS Contour, Belotero INTENSE Lidocaine, Belotero LIPS Shape Is an injectable biodegradable implant indicated for filling of depressions due to facial atrophic scars, filling of facial wrinkles and folds as well as for lip enhancement. The presence of lidocaine aims to reduce local pain associated **with** the injection of the gel and to improve patient comfort.

For use by medical practitioner.

For successful treatment, it is essential that the practitioner has received a special training on the injection techniques for soft tissue augmentation and for increasing the volume of skin tissues. A good knowledge of the anatomy and physiology of the site to be treated is required.

Scope of application

Dermatology

1. Description

BELOTERO® Soft Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. BELOTERO® Soft Lidocaine contains 0.3% of lidocaine hydrochloride.

2. Presentation

BELOTERO® Soft Lidocaine is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, **two** traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate: 20.0 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: gel volume

The volume of the gel in each syringe is stated on the external box.

4. Indications

- BELOTERO® Soft Lidocaine is an injectable resorbable implant indicated to fill fine facial lines and wrinkles as well as for lip enhancement. It is also suitable for the correction of facial atrophic scars.
- The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

5. Posology and administration method

- BELOTERO® Soft Lidocaine is designed to be injected into the superficial to mid-dermis by a legally approved practitioner.
- BELOTERO® Soft Lidocaine can be used for all skin types.
- For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for filling wrinkles and folds. A good knowledge of the anatomy and physiology of the site to be treated is required.
- The treatment must be carried out under appropriate aseptic conditions.
- BELOTERO® Soft Lidocaine must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.
- It is recommended to use one of the supplied needles.
- To ensure optimal use of BELOTERO® Soft Lidocaine, it is recommended to assemble the needle according to the diagrams below.
- BELOTERO® Soft Lidocaine can be used in combination with other BELOTERO® products during the same session. Instructions for use of each product should be followed.

- If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.
- BELOTERO® Soft Lidocaine should be injected slowly. The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.
- Gently massage the treated area after the injection to distribute the product uniformly.

6. Contra-indications

BELOTERO® Soft Lidocaine is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride or to amide-type local anesthetics,
- In pregnant and breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection.
- Do not inject BELOTERO® Soft Lidocaine into blood vessels.
- Do not inject BELOTERO® Soft Lidocaine into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).
- Do not inject BELOTERO® Soft Lidocaine into an area previously treated with a permanent dermal filler.

7. Precautions for use

- Before treatment the patient must be informed about the device, its contra-indications and possible side effects.
- In the absence of available clinical data on tolerance on the injection of BELOTERO® Soft Lidocaine in patients with antecedents or with an active autoimmune disease or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO® Soft Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.
- It is recommended not to inject BELOTERO® Soft Lidocaine in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or cheloids.
- BELOTERO® Soft Lidocaine can be injected in the glabellar area only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area could potentially lead to local vascular occlusion, ischemia and necrosis.
- BELOTERO® Soft Lidocaine can be injected in the periorbital region only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects.

ante

- BELOTERO® **Soft** Lidocaine must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO® Soft Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO® Soft Lidocaine with the above-mentioned treatments.
- No clinical data is available on the injection of BELOTERO® Soft Lidocaine into an area already treated with another filling product.
- Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.
- Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.
- In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO® Soft Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.
- Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.
- It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.
- For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.
- When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.
- Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.
- Do not transfer BELOTERO® Soft Lidocaine into another container and do not add other substances to the product.
- Only the gel is sterile, but not the outside of the syringe.
- BELOTERO® Soft Lidocaine must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.
- Discard the syringe and the remaining product after use.
- Do not re-sterilize and do not reuse due to the associated risks including infection.

- The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

8. Incompatibilities

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO® Soft Lidocaine does not come into contact with such substances.
- There is no known interaction with other local or loco-regional anesthetics.

9. Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases, one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema, or pain sometimes accompanied by itching in the treated area. These reactions may last for a week,
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate and lidocaine hydrochloride,
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

10. Assembly of needle to syringe

For optimal use of BELOTERO® Soft Lidocaine, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

- (1) **Firmly hold** the glass cylinder of the **syringe and the Luer-lock adaptor** between the thumb and the index fingers.
- (2) Grasp the protective cap with the other hand and unscrew it.
- (3) **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
- (4) Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

12. Disposal


Used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with medical procedures adopted in medical institutions and local, state or federal regulations.

13. Requirements for environmental protection


The device should be disposed of in accordance with any applicable local/national regulations. With proper disposal, there is no harm to the environment.

14. References

Updated documentation may be available from ANTEIS SA, Switzerland.

Backstop in the right position during injection

Caution

Consult instructions for use

Do not use if package is damaged

Single use product. Do not re-use

Open the blister by pulling the Tyvek lid following the arrow

Sterile. Sterilized by moist heat.
Only the gel is sterile, but not the outside of the syringe.

Stenle. Sterilized by irradiation.
Only the needle itself is sterile, but not the outside of the needle

Temperature limit of storage: 2 °C – 25 °C

Date of manufacture

Use-by date

Batch code

CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number.

Manufacturer

Manufacturer of BELOTERO® Soft Lidocaine:

ANTEIS SA

18 Chemin des Aulx
CH-1228 Plan-les-Ouates
Geneva, Switzerland

Distributed by:

Authorized Representative in RF

LLC «Merz Pharma»: 125040, Moscow, str. 3-ya Yamskogo polya, 2, bld. 26, 2nd floor,
premises VII, room 29; tel. +7 (495) 653 8555; www.merz.ru

Date of the instructions for use: 2015-04-13

..

2

BELOTERO® BALANCE LIDOCAINE

MERZ AESTHETICS

INSTRUCTIONS FOR USE FOR BELOTERO® BALANCE LIDOCAINE

Belotero line, viscoelastic Implant based on cross-linked sodium hyaluronate for correction of depressions of tissue, in the following versions:

1) Belotero SOFT Lidocaine, consisting of:

- implantable gel, in a 1.0 ml syringe - 1 pc.;
- sterile needle 30G ½ (0,3 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

2) Belotero BALANCE Lidocaine:

Belotero BALANCE Lidocaine;

Belotero LIPS Contour,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) and 30G ½ (0,3 x 13) - 1 pc. each;
- label - 2 pcs.;
- Instructions for use – 1 pc.

3) Belotero INTENSE Lidocaine:

Belotero INTENSE Lidocaine;

Belotero LIPS Shape,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

Intended Use

Belotero line in the following versions: Belotero SOFT Lidocaine, Belotero BALANCE Lidocaine, Belotero LIPS Contour, Belotero INTENSE Lidocaine, Belotero LIPS Shape is an injectable biodegradable implant indicated for filling of depressions due to facial atrophic scars, filling of facial wrinkles and folds as well as for lip enhancement. The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

For use by medical practitioner.

For successful treatment, it is essential that the practitioner has received a specific training on the injection techniques for soft tissue augmentation and for increasing the volume of skin tissues. A good knowledge of the anatomy and physiology of the site to be treated is required.

1. Description

BELOTERO® Balance Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. BELOTERO® Balance Lidocaine contains 0.3% of lidocaine hydrochloride.

2. Presentation

BELOTERO® Balance Lidocaine is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate: 22.5 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: gel volume

The volume of the gel in each syringe is stated on the external box.

4. Indications

- BELOTERO® Balance Lidocaine is an injectable resorbable implant indicated to fill moderate facial wrinkles and folds as well as for lip enhancement. It is also suitable for correction of facial atrophic scars.
- The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

5. Posology and administration method

- BELOTERO® Balance Lidocaine is designed to be injected into the superficial to mid-dermis by a legally approved practitioner.
- BELOTERO® Balance Lidocaine can be used for all skin types.
- For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for filling wrinkles and folds. A good knowledge of the anatomy and physiology of the site to be treated is required.
- The treatment must be carried out under appropriate aseptic conditions.
- BELOTERO® Balance Lidocaine must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.
- It is recommended to use one of the supplied needles.
- To ensure optimal use of BELOTERO® Balance Lidocaine, it is recommended to assemble the needle according to the diagrams below.

- The practitioner should note that the use of 30G½ needles requires greater force to inject the gel than 27G½ needle.
- BELOTERO® Balance Lidocaine can be used in combination with other BELOTERO® products during the same session. Instructions for use of each product should be followed.
- If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.
- BELOTERO® Balance Lidocaine should be injected slowly. The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.
- Gently massage the treated area after the injection to distribute the product uniformly.

6. Contra-indications

BELOTERO® Balance Lidocaine is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride or to amide-type local anesthetics,
- In pregnant and breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection.
- Do not inject BELOTERO® Balance Lidocaine into blood vessels.
- Do not inject BELOTERO® Balance Lidocaine into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).
- Do not inject BELOTERO® Balance Lidocaine into an area previously treated with a permanent dermal filler.

7. Precautions for use

- Before treatment the patient must be informed about the device, its contra-indications and possible side effects.
- In the absence of available clinical data on tolerance and efficacy of the injection of BELOTERO® Balance Lidocaine in patients with antecedents or with an active autoimmune disease, or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO® Balance Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.
- It is recommended not to inject BELOTERO® Balance Lidocaine in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or cheloids.
- BELOTERO® Balance Lidocaine can be injected in the glabellar area only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area could potentially lead to local vascular occlusion, ischemia and necrosis.

- BELOTERO® Balance Lidocaine can be injected in **the** periorbital region only by trained experienced practitioners **with** a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects.
- BELOTERO® Balance Lidocaine must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO® Balance Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO® Balance Lidocaine with the above-mentioned treatments.
- No clinical data is available on the injection of BELOTERO® Balance Lidocaine into an area already treated with another filling product.
- Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.
- Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.
- In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO® Balance Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.
- Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.
- It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.
- For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.
- When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.
- Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.
- Do not transfer BELOTERO® Balance Lidocaine into another container and do not add other substances to the product.
- Only the gel is sterile, but not the outside of the syringe.
- BELOTERO® Balance Lidocaine must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.
- Discard the syringe and the remaining product after use.

- Do not re-sterilize and do not reuse due to the associated risks including infection.
- The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

8. Incompatibilities

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO® Balance Lidocaine does not come into contact with such substances.
- There is no known interaction with other local or loco-regional anesthetics.

9. Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain sometimes accompanied by itching in the treated area. These reactions may last for a week.
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate and lidocaine hydrochloride.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

10. Assembly of needle to syringe

For optimal use of BELOTERO® Balance Lidocaine, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

(1) **Firmly hold** the glass cylinder of the **syringe and the Luer-lock adaptor** between the thumb and the index fingers.

(2) Grasp the protective cap with the other hand and unscrew it.

(3) **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.

(4) Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

12. Disposal



Used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance **with** medical procedures adopted in medical institutions and local, state or federal regulations.

13. Requirements for environmental protection





The device should be disposed of in accordance **with** any applicable local/national regulations. With proper disposal, there is no harm to the environment.

14. References

Updated documentation may be available from ANTEIS SA, Switzerland.

Backstop in the right position during injection

Caution

Consult instructions for use

Do not use if package is damaged

Single use product. Do not re-use


Open the blister by pulling the Tyvek lid following the arrow

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.

Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Temperature limit of storage: 2 °C – 25 °C

Date of manufacture

Use-by date

Batch code

CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number.

Manufacturer

Manufacturer of BELOTERO® Balance Lidocaine:

ANTEIS SA

18 Chemin des Aulx
CH-1228 Plan-les-Ouates
Geneva, Switzerland

Distributed by:

Authorized Representative in RF

LLC «Merz Pharma»: 125040, Moscow, str. 3-ya Yamskogo polya, 2, bld. 26, 2nd floor,
premises VII, room 29; tel. +7 (495) 653 8555; www.merz.ru

Date of the instructions for use: 2015-04-13

2

BELOTERO® LIPS CONTOUR

MERZ AESTHETICS

INSTRUCTIONS FOR USE FOR BELOTERO® LIPS CONTOUR

Belotero line viscoelastic implant based on cross-linked sodium hyaluronate for correction of depressions of tissue, in the following versions:

1) Belotero SOFT Lidocaine, consisting of:

- implantable gel, in a 1.0 ml syringe - 1 pc.;
- sterile needle 30G ½ (0,3 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

2) Belotero BALANCE Lidocaine:

Belotero BALANCE Lidocaine;

Belotero LIPS Contour,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) and 30G ½ (0,3 x 13) - 1 pc. each;
- label - 2 pcs.;
- Instructions for use – 1 pc.

3) Belotero INTENSE Lidocaine:

Belotero INTENSE Lidocaine;

Belotero LIPS Shape,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

Intended Use

Belotero line in the following versions: Belotero SOFT Lidocaine, Belotero BALANCE Lidocaine, Belotero LIPS Contour, Belotero INTENSE Lidocaine, Belotero LIPS Shape is an injectable biodegradable implant indicated for filling of depressions due to facial atrophic scars, filling of facial wrinkles and folds as well as for lip enhancement. The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

For use by medical practitioner.

For successful treatment, it is essential that the practitioner has received a special training on the injection techniques for soft tissue augmentation and for increasing the volume of skin tissues. A good knowledge of the anatomy and physiology of the site to be treated is required.

1. Description

Belotero® Lips Contour is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. Belotero® Lips Contour contains 0.3% of lidocaine hydrochloride.

2. Presentation

Belotero® Lips Contour is presented in a single use pre-filled glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate 22.5 mg/ml

Lidocaine hydrochloride: 3.0 mg/ml

Phosphate buffer pH 7 q.s.: gel volume

The volume of the gel in each syringe is as stated on the external box.

4. Indications

- Belotero® Lips Contour is an injectable resorbable Implant indicated to correct facial atrophic scars and for lip enhancement (e.g. lip contour, perioral fine lines, moderate oral commissures).
- The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

5. Posology and administration method

- Belotero® Lips Contour is designed to be injected into the superficial to mid-dermis and the mucous membrane of the lips by a legally approved practitioner.
- Belotero® Lips Contour can be used for all skin types.
- For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for filling wrinkles and folds and for lip enhancement. A good knowledge of the anatomy and physiology of the site to be treated is required.
- The treatment must be carried out under appropriate aseptic conditions.
- Belotero® Lips Contour must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.
- It is recommended to use one of the supplied needles.
- To ensure optimal use of Belotero® Lips Contour, it is recommended to assemble the needle according to the diagrams below.

Instructions for use –Belotero® Lips Contour - Russia

- The practitioner should note that the use of 30G½ needles requires greater force to inject the gel than 27G½ needle.
- Belotero® Lips Contour can be used in combination with other Belotero® products during the same session. Instructions for use of each product should be followed.
- If the needle becomes obstructed and the injection pressure is too high, stop the injection and change the needle.
- Belotero® Lips Contour should be injected slowly. The quantity of the gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.
- Gently massage the treated area after the injection to distribute the product uniformly.

6. Contra-indications

Belotero® Lips Contour is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride or to amide-type local anesthetics,
- In pregnant and breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection.
- Do not inject Belotero® Lips Contour into blood vessels.
- Do not inject Belotero® Lips Contour into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).
- Do not inject Belotero® Lips Contour into an area previously treated with a permanent dermal filler.

7. Precautions for use

- Before treatment the patient must be informed about the device, its contra-indications and possible side effects.
- In the absence of available clinical data on tolerance and efficacy of the injection of Belotero® Lips Contour in patients with antecedents or with an active autoimmune disease, or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject Belotero® Lips Contour on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.
- It is recommended not to inject Belotero® Lips Contour in patients with a history of streptococcal diseases and in patients pre-disposed to hypertrophic scars or cheloids.
- Belotero® Lips Contour can be injected in the glabellar area only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area could potentially lead to local vascular occlusion, ischemia and necrosis.

ante

Instructions

Discard

- Belotero® Lips Contour can be injected in the periorbital region only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects.
- Belotero® Lips Contour must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, Belotero® Lips Contour must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of Belotero® Lips Contour with the above-mentioned treatments.
- No clinical data is available on the injection of Belotero® Lips Contour into an area already treated with another filling product.
- Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.
- Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.
- In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject Belotero® Lips Contour on a case-by-case basis depending of the nature of the disease as well as the associated treatment.
- Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.
- It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.
- For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.
- When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.
- Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.
- Do not transfer Belotero® Lips Contour into another container and do not add other substances to the product.
- Only the gel is sterile, but not the outside of the syringe.
- Belotero® Lips Contour must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.

- Discard the syringe and the remaining product after use.
- Do not re-sterilize and do not reuse due to the associated risks including infection.
- The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

8. Incompatibilities

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that Belotero® Lips Contour does not come into contact with such substances.
- There is no known interaction with other local or loco-regional anesthetics.

9. Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain sometimes accompanied by itching in the treated area. These reactions may last for a week.
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate and lidocaine hydrochloride.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

10. Assembly of needle to syringe

For optimal use of Belotero® Lips Contour, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

(1) **Firmly hold** the glass cylinder of the **syringe and the Luer-lock adaptor** between the thumb and the index fingers.

(2) Grasp the protective cap with the other hand and unscrew it.

(3) **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.

(4) Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

12. Disposal



Used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with medical procedures adopted in medical institutions and local, state or federal regulations.

13. Requirements for environmental protection


The device should be disposed of in accordance with any applicable local/national regulations. With proper disposal, there is no harm to the environment.

14. References

Updated documentation may be available from ANTEIS SA, Switzerland.

Backstop in the right position during injection

Caution

Consult instructions for use

Do not use if package is damaged

Single use product. Do not re-use

Open the blister by pulling the Tyvek lid following the arrow

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.

Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Temperature limit of storage: 2 °C – 25 °C

Date of manufacture

Use-by date

Batch code

CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number.

Manufacturer

Manufacturer of Belotero® Lips Contour:

ANTEIS SA

18 Chemin des Aulx
CH-1228 Plan-les-Ouates
Geneva, Switzerland

Distributed by:

Authorized Representative in RF

LLC «Merz Pharma»: 125040, Moscow, str. 3-ya Yamskogo polya, 2, bld. 26, 2nd floor,
premises VII, room 29; tel. +7 (495) 653 8555; www.merz.ru

Date of the instructions for use: 2015-04-13

A large, stylized blue ink signature that spans across the middle of the page.

BELOTERO® INTENSE LIDOCAINE

MERZ AESTHETICS

INSTRUCTIONS FOR USE FOR BELOTERO® INTENSE LIDOCAINE

Belotero line viscoelastic implant based on cross-linked sodium hyaluronate for correction of depressions of tissue, in the following versions:

1) Belotero SOFT Lidocaine, consisting of:

- implantable gel, in a 1.0 ml syringe - 1 pc.;
- sterile needle 30G ½ (0,3 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

2) Belotero BALANCE Lidocaine:

Belotero BALANCE Lidocaine;

Belotero LIPS Contour,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) and 30G ½ (0,3 x 13) - 1 pc. each;
- label - 2 pcs.;
- Instructions for use – 1 pc.

3) Belotero INTENSE Lidocaine:

Belotero INTENSE Lidocaine;

Belotero LIPS Shape,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

Intended Use

Belotero line in the following versions: Belotero SOFT Lidocaine, Belotero BALANCE Lidocaine, Belotero LIPS Contour, Belotero INTENSE Lidocaine, Belotero LIPS Shape is an injectable biodegradable implant indicated for filling of depressions due to facial atrophic scars, filling of facial wrinkles and folds as well as for lip enhancement. The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

For use by medical practitioner.

For successful treatment, it is essential that the practitioner has received a special training on the injection techniques for soft tissue augmentation and for increasing the volume of skin tissues. A good knowledge of the anatomy and physiology of the site to be treated is required.

1. Description

BELOTERO® Intense Lidocaine is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. BELOTERO® Intense Lidocaine contains 0.3% of lidocaine hydrochloride.

2. Presentation

BELOTERO® Intense Lidocaine is presented in a single use pre-filled sterile glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate:	25.5 mg/ml
Lidocaine hydrochloride:	3.0 mg/ml
Phosphate buffer pH 7 q.s.:	gel volume

The volume of the gel in each syringe is stated on the external box.

4. Indications

- BELOTERO® Intense Lidocaine is an injectable resorbable implant indicated to fill deep wrinkles and folds, as well as to restore and enhance soft tissue volume (e.g. contours of the face, lip volume etc.). It is also suitable for correction of facial atrophic scars.
- The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

5. Posology and administration method

- BELOTERO® Intense Lidocaine is designed to be injected into the deep dermis by a legally approved practitioner.
- BELOTERO® Intense Lidocaine can be used for all skin types.
- For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for soft tissue augmentation. A good knowledge of the anatomy and physiology of the site to be treated is required.
- The treatment must be carried out under appropriate aseptic conditions.

- BELOTERO® Intense Lidocaine must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.
- It is recommended to use one of the supplied needles.
- To ensure optimal use of BELOTERO® Intense Lidocaine, it is recommended to assemble the needle according to the diagrams below.
- The use of the enclosed 27G½ needle is recommended, as a smaller needle diameter would require a greater force to inject the implant.
- BELOTERO® Intense Lidocaine can be used in combination with other BELOTERO® products during the same session. Instructions for use of each product should be followed.
- If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle.
- BELOTERO® Intense Lidocaine should be injected slowly. The quantity of gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.
- Gently massage the treated area after the injection to distribute the product uniformly.

6. Contra-indications

BELOTERO® Intense Lidocaine is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride and amide-type local anesthetics,
- In pregnant or breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection.
- Do not inject BELOTERO® Intense Lidocaine into blood vessels.
- Do not inject BELOTERO® Intense Lidocaine into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).
- Do not inject BELOTERO® Intense Lidocaine into an area previously treated with a permanent dermal filler.
- Do not inject BELOTERO® Intense Lidocaine in the glabellar region.

7. Precautions for use

- Before treatment the patient must be informed about the device, its contra-indications and possible side effects.
- In the absence of available clinical data on tolerance and efficacy of the injection of BELOTERO® Intense Lidocaine in patients with antecedents or with an active autoimmune disease or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject BELOTERO® Intense Lidocaine on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to

these patients and to not inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.

- It is recommended not to inject BELOTERO® Intense Lidocaine in patients with a history of streptococcal diseases or in patients pre-disposed to hypertrophic scars or cheloids.
- BELOTERO® Intense Lidocaine can be injected in the periorbital region only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects.
- BELOTERO® Intense Lidocaine can be injected in the nose only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may lead to local vascular complications such as ischemia or necrosis.
- BELOTERO® Intense Lidocaine must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, BELOTERO® Intense Lidocaine must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of BELOTERO® Intense Lidocaine with the above-mentioned treatments.
- No clinical data is available on the injection of BELOTERO® Intense Lidocaine into an area already treated with another filling product.
- Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.
- Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.
- In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject BELOTERO® Intense Lidocaine on a case-by-case basis depending of the nature of the disease as well as the associated treatment.
- Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.
- It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.
- For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.
- When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.

- Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.
- Do not transfer BELOTERO® Intense Lidocaine into another container and do not add other substances to the product.
- Only the gel is sterile, but not the outside of the syringe.
- BELOTERO® Intense Lidocaine must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.
- Discard the syringe and the remaining product after use.
- Do not re-sterilize and do not reuse due to the associated risks including infection.
- The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

8. Incompatibilities

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that BELOTERO® Intense Lidocaine does not come into contact with such substances.
- There is no known interaction with other local or loco-regional anesthetics.

9. Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain, sometimes accompanied by itching in the treated area. These reactions may last for a week.
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate and lidocaine hydrochloride.
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

10. Assembly of needle to syringe

For optimal use of BELOTERO® Intense Lidocaine, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

- (1) **Firmly hold** the glass cylinder of the **syringe and the Luer-lock adaptor** between the thumb and the index fingers.
- (2) Grasp the protective cap with the other hand and unscrew it.
- (3) **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.
- (4) Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

12. Disposal



Used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with medical procedures adopted in medical institutions and local, state or federal regulations.

13. Requirements for environmental protection


The device should be disposed of in accordance with any applicable local/national regulations. With proper disposal, there is no harm to the environment.

14. References

Updated documentation may be available from ANTEIS SA, Switzerland.

Backstop in the right position during injection

Caution

Consult instructions for use

Do not use if package is damaged

Single use product. Do not re-use.

Open the blister by pulling the Tyvek lid following the arrow

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.

Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Temperature limit of storage: 2 °C – 25 °C

Date of manufacture

Use-by date

.. Batch code

CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number.

Manufacturer

Manufacturer of BELOTERO® Intense Lidocaine:

ANTEIS SA
18 Chemin des Aulx
CH-1228 Plan-les-Ouates
Geneva, Switzerland

Distributed by:

Authorized Representative in RF

LLC «Merz Pharma»: 125040, Moscow, str. 3-ya Yamskogo polya, 2, bld. 26, 2nd floor,
premises VII, room 29; tel. +7 (495) 653 8555; www.merz.ru

Date of the instructions for use: 2015-04-13

BELOTERO® LIPS SHAPE

MERZ AESTHETICS

INSTRUCTIONS FOR USE FOR BELOTERO® LIPS SHAPE

Belotero line viscoelastic implant based on cross-linked sodium hyaluronate for correction of depressions of tissue, in the following versions:

1) Belotero SOFT Lidocaine, consisting of:

- implantable gel, in a 1.0 ml syringe - 1 pc.;
- sterile needle 30G ½ (0,3 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

2) Belotero BALANCE Lidocaine:

Belotero BALANCE Lidocaine;

Belotero LIPS Contour,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) and 30G ½ (0,3 x 13) - 1 pc. each;
- label - 2 pcs.;
- Instructions for use – 1 pc.

3) Belotero INTENSE Lidocaine:

Belotero INTENSE Lidocaine;

Belotero LIPS Shape,

consisting of:

- implantable gel, in a 1.0 ml or 0.6 ml syringe - 1 pc.;
- sterile needle 27G ½ (0,4 x 13) - 2 pcs.;
- label - 2 pcs.;
- Instructions for use – 1 pc.

Intended Use

Belotero line in the following versions: Belotero SOFT Lidocaine, Belotero BALANCE Lidocaine, Belotero LIPS Contour, Belotero INTENSE Lidocaine, Belotero LIPS Shape is an injectable biodegradable implant indicated for filling of depressions due to facial atrophic scars, filling of facial wrinkles and folds as well as for lip enhancement. The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

For use by medical practitioner.

For successful treatment, it is essential that the practitioner has received a special training on the injection techniques for soft tissue augmentation and for increasing the volume of skin tissues. A good knowledge of the anatomy and physiology of the site to be treated is required.

1. Description

Belotero® Lips Shape is a sterile, non-pyrogenic, viscoelastic, colourless, transparent cross-linked sodium hyaluronate gel of non-animal origin in a physiological phosphate buffer. Belotero® Lips Shape contains 0.3% of lidocaine hydrochloride.

2. Presentation

Belotero® Lips Shape is presented in a single use pre-filled sterile glass syringe sterilized by moist heat. Each box contains one instruction leaflet, one syringe, two traceability labels and two sterile CE-marked needles for single use only. The dimensions of needles are stated on the external box.

3. Composition

Cross-linked sodium hyaluronate:	25.5 mg/ml
Lidocaine hydrochloride:	3.0 mg/ml
Phosphate buffer pH 7 q.s.:	gel volume

The volume of the gel in each syringe is as stated on the external box.

4. Indications

- Belotero® Lips Shape is an injectable resorbable implant indicated to correct facial atrophic scars and to enhance soft tissue volume (e.g. lip volume, severe oral commissures).
- The presence of lidocaine aims to reduce local pain associated with the injection of the gel and to improve patient comfort.

5. Posology and administration method

- Belotero® Lips Shape is designed to be injected into the deep dermis and the mucous membrane of the lips by a legally approved practitioner.
- Belotero® Lips Shape can be used for all skin types.
- For successful treatment it is essential that the practitioner has received a specific training on the injection techniques for soft tissue augmentation. A good knowledge of the anatomy and physiology of the site to be treated is required.
- The treatment must be carried out under appropriate aseptic conditions.
- Belotero® Lips Shape must be injected into a healthy, non-inflamed and previously rigorously disinfected skin.
- It is recommended to use one of the supplied needles.

- To ensure optimal use of Belotero® Lips Shape, it is recommended to assemble the needle according to the diagrams below.
- The use of the enclosed 27G½ needle is recommended, as a smaller needle diameter would require a greater force to inject the implant.
- Belotero® Lips Shape can be used in combination with other Belotero® products during the same session. Instructions for use of each product should be followed.
- If the needle becomes obstructed and the injection pressure becomes too high, stop the injection and change the needle.
- Belotero® Lips Shape should be injected slowly. The quantity of gel to be injected depends on the area to be treated and the correction to be achieved. Do not over-correct.
- Gently massage the treated area after the injection to distribute the product uniformly.

6. Contra-indications

Belotero® Lips Shape is contra-indicated:

- In case of known hypersensitivity to one of the product's components, especially to sodium hyaluronate, lidocaine hydrochloride and amide-type local anesthetics,
- In pregnant or breast-feeding women,
- In young patients under 18 years old,
- In patients presenting a general infection.
- Do not inject Belotero® Lips Shape into blood vessels.
- Do not inject Belotero® Lips Shape into areas presenting cutaneous problems of an inflammatory or infectious type (acne, herpes...).
- Do not inject Belotero® Lips Shape into an area previously treated with a permanent dermal filler.
- Do not inject Belotero® Lips Shape in the glabellar region.

7. Precautions for use

- Before treatment the patient must be informed about the device, its contra-indications and possible side effects.
- In the absence of available clinical data on tolerance and efficacy of the injection of Belotero® Lips Shape in patients with antecedents or with an active autoimmune disease or in patients presenting a history of severe multiple allergies or anaphylactic shock, the practitioner must decide whether to inject Belotero® Lips Shape on a case-by-case basis depending on the nature of the disease as well as the associated treatment. It is recommended to propose a prior double test to these patients and to not inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection.
- It is recommended not to inject Belotero® Lips Shape in patients with a history of streptococcal diseases or in patients pre-disposed to hypertrophic scars or cheloids.

- Belotero® Lips Shape can be injected in the periorbital region only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may be associated with an increased frequency and severity of side effects.
- Belotero® Lips Shape can be injected in the nose only by trained experienced practitioners with a deep knowledge of the anatomy. Injection of dermal fillers into this area may lead to local vascular complications such as ischemia or necrosis.
- Belotero® Lips Shape must not be used in association with other aesthetic medicine techniques such as peeling, dermabrasion, or any type of laser treatment before complete healing of the last treatment. In any case, even if the healing occurs earlier, Belotero® Lips Shape must not be used earlier than 2 weeks after the last treatment. No clinical data is available on the combined use of Belotero® Lips Shape with the above-mentioned treatments.
- No clinical data is available on the injection of Belotero® Lips Shape into an area already treated with another filling product.
- Check the integrity of the inner packaging and the expiry date for both the syringe and the needle prior to use. Do not use these products if the expiry date has lapsed or if the packaging has been opened or damaged.
- Patients using antithrombotic substances such as aspirin or non-steroidal anti-inflammatory medications may have increased reactions of hematomas, nodules or bleeding at the injection site.
- In cases of patients suffering from epilepsy, impaired cardiac conditions, severely impaired hepatic function or severe renal dysfunction or porphyria, the practitioner must decide whether to inject Belotero® Lips Shape on a case-by-case basis depending of the nature of the disease as well as the associated treatment.
- Practitioners and athletes should consider that lidocaine may produce positive results to anti-doping tests.
- It should be noted that the presence of lidocaine may cause local redness or hypersensitivity.
- For normal healthy adults, it is recommended that the maximum total dose of lidocaine HCl (without epinephrine) does not exceed 300 mg (or 4.5 mg/kg) per session. Overdosage of lidocaine HCl usually results in sign of the central nervous system or cardiovascular toxicity.
- When using concurrently (topical administration...), the total administered dose of lidocaine should be considered. The concomitant use of other local anesthetic agents or agents structurally related to amide-type local anesthetics should also be considered since the systemic toxic effects may be additive.
- Care should be taken for patients with congenital methemoglobinemia, with glucose-6-phosphate dehydrogenase deficiencies and patients who are receiving concomitant treatment with methemoglobin-inducing agents.
- Do not transfer Belotero® Lips Shape into another container and do not add other substances to the product.
- Only the gel is sterile, but not the outside of the syringe.

Instructions for use –Belotero® Lips Shape - Russia

- Belotero® Lips Shape must not be used with an automated injection system not recommended by MERZ/ANTEIS. If an automatic system is used, it is recommended that the practitioner has previously read its instructions for use and is trained on the use of the system.
- Discard the syringe and the remaining product after use.
- Do not re-sterilize and do not reuse due to the associated risks including infection.
- The patient must avoid applying makeup for at least 12 hours after treatment as well as avoid saunas, Turkish baths and prolonged exposure to the sun or UV rays for 2 weeks after the treatment. Patients should also avoid putting pressure on and/or handling the treated area.

8. Incompatibilities

- Sodium hyaluronate precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride). It is therefore recommended that Belotero® Lips Shape does not come into contact with such substances.
- There is no known interaction with other local or loco-regional anesthetics.

9. Side effects

Patients must be informed by the practitioner about possible side effects before treatment. A slight bleeding may occur during the injection, and it disappears spontaneously as soon as the injection is finished. In occasional cases one or more of the following may occur either immediately or as a delayed reaction (list not exhaustive):

- Reactions usually associated with injections such as redness, erythema, oedema or pain, sometimes accompanied by itching in the treated area. These reactions may last for a week,
- Hematomas in the treated area,
- Swelling in the treated area,
- Indurations or nodules in the treated area,
- Coloration or discoloration in the treated area,
- Allergy to one of the product's components, especially to sodium hyaluronate and lidocaine hydrochloride,
- Cases of necrosis, abscesses and granulomas after sodium hyaluronate injections have been reported in the literature. These rare potential risks must nevertheless be considered.

Patients should be instructed to report any side effects which last for more than one week to his/her practitioner. The practitioner may then prescribe the patient appropriate treatment.

10. Assembly of needle to syringe

For optimal use of Belotero® Lips Shape, it is important that the needle is properly connected to the syringe. See diagrams 1, 2, 3 and 4.

(1) **Firmly hold** the glass cylinder of the syringe and the Luer-lock adaptor between the thumb and the index fingers.

(2) Grasp the protective cap with the other hand and unscrew it.

(3) **Push & Twist** the needle on the syringe **until a resistance is felt**. Do not over-tight. Over-tightening of the needle may lead to the Luer-lock moving and dislodging from the syringe.

(4) Keep holding the Luer-lock and remove the sheath from the needle.

11. Storage

Store between 2°C and 25°C. Protect from light and freezing. Avoid mechanical shocks.

12. Disposal


Used and partially used syringes and needles for injection can be biologically dangerous and therefore should be considered and disposed of in accordance with medical procedures adopted in medical institutions and local, state or federal regulations.

13. Requirements for environmental protection


The device should be disposed of in accordance with any applicable local/national regulations. With proper disposal, there is no harm to the environment.

14. References

Updated documentation may be available from ANTEIS SA, Switzerland.

Backstop in the right position during injection

Caution

Consult instructions for use

Do not use if package is damaged

Single use product. Do not re-use

Open the blister by pulling the Tyvek lid following the arrow

Sterile. Sterilized by moist heat. Only the gel is sterile, but not the outside of the syringe.

Sterile. Sterilized by irradiation. Only the needle itself is sterile, but not the outside of the needle packaging.

Temperature limit of storage: 2 °C – 25 °C

Date of manufacture

Use-by date

Batch code

CE mark in accordance with Directive 93/42/EEC relating to medical devices. This mark is followed by the notified body number.

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Manufacturer of Belotero® Lips Shape:

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