
BE IT KNOWN that I, JONATHAN RICHARD KEMP, of Newcastle upon Tyne, England, NOTARY PUBLIC duly admitted and sworn, authorised to practise throughout England and Wales, CERTIFY that the attached Instructions for Use was signed before me in Newcastle upon Tyne today by Steven Elstob, who was identified to me and is a Director and the Finance Director of **LEICA BIOSYSTEMS NEWCASTLE LIMITED** (the "Company"). I also CERTIFY that the Company exists as a United Kingdom company with company number 02163063.

SIGNED and SEALED with my seal of office in Newcastle upon Tyne, England, on this nineteenth day of December in the year two thousand and twenty-three.

Jonathan Richard Kemp
Notary Public, England and Wales
Protocol No. 2023/663.B

For Russian Federation only

Instructions for Use

BOND Primary monoclonal Antibody Anaplastic Lymphoma Kinase (5A4) for the qualitative identification by immunohistochemical method for in vitro diagnostics for BOND-III and Bond-maX immunohistochemistry instruments

Designer:

Leica Biosystems Newcastle Limited
Balliol Business Park West, Benton Lane,
Newcastle upon Tyne, NE12 8EW,
United Kingdom
Tel: +44 191 215 0567
Fax: + 44 191 215 1152
e-mail: msds@leicabiosystems.com

Name and address of Production site:

Leica Biosystems Newcastle Limited
Balliol Business Park West, Benton Lane,
Newcastle upon Tyne, NE12 8EW,
United Kingdom

Legal Manufacturer:

Leica Biosystems Newcastle Limited
Balliol Business Park West, Benton Lane,
Newcastle upon Tyne, NE12 8EW,
United Kingdom
Tel: +44 191 215 0567
Fax: + 44 191 215 1152
e-mail: msds@leicabiosystems.com

Authorized Representative in the Russian Federation

BioLine LLC
197022 Saint-Petersburg, municipal district Aptekarsky ostrov, Professora Popova street 23 lit. E,
structure 3H/5H, room 208, Russia


Tel: +7 (812) 320 49 49
Fax: +7 (812) 320 49 40
e-mail: main@bioline.ru

Approved by:

Full name: Steven Elstob

Position: Director and Finance Director, Advanced Staining Reagents

Company Name: Leica Biosystems Newcastle Limited

Signature: 

Date: 19th December 2023

Leica
BIOSYSTEMS

Leica Biosystems Newcastle Ltd
Balliol Business Park
Benton Lane
Newcastle upon Tyne
NE12 8EW

NAME

BOND Primary monoclonal Antibody Anaplastic Lymphoma Kinase (5A4) for the qualitative identification by immunohistochemical method for in vitro diagnostics for BOND-III and Bond-maX immunohistochemistry instruments, composed of:

- BOND Primary Antibody Anaplastic Lymphoma Kinase (5A4), 7.0 ml, container - 1 pc.
- Instructions for Use (included) - 1 pc.

(Hereafter referred to as the BOND Primary Antibody Anaplastic Lymphoma Kinase (5A4), Anaplastic Lymphoma Kinase (5A4), the product, IVD medical device)

INTENDED USE

Anaplastic Lymphoma Kinase (5A4) monoclonal antibody for BOND-III and Bond-maX is intended to be used for the qualitative identification by light microscopy of endogenous ALK and ALK fusion proteins in formalin-fixed, paraffin embedded tissue by immunohistochemical staining using the BOND-III and Bond-maX.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

INTRODUCTION

Immunohistochemical (IHC) staining techniques allow for the visualization of antigens via the sequential application of a specific antibody to the antigen (primary antibody), a secondary antibody to the primary antibody and an enzyme complex with a chromogenic substrate with interposed washing steps. The enzymatic activation of the chromogen results in a visible reaction product at the antigen site. The specimen may then be counterstained and coverslipped. Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular antigen.

PRINCIPLE OF PROCEDURE

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Anaplastic Lymphoma Kinase (5A4) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of endogenous ALK and ALK fusion proteins is achieved by first allowing the binding of Anaplastic Lymphoma Kinase (5A4) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

SCOPE OF APPLICATION

Clinical laboratory medicine.

For in vitro diagnostics.

For professional use.

POTENTIAL USER

This product is intended for professional use in laboratories and clinicodiagnostic departments of medical institutions by medical specialists or other persons as intended, only for *in vitro* diagnostics.

Professional level of potential users

The medical product is intended only for professional use by pathologists, laboratory technicians, medical laboratory assistants or other specialists in the field of clinical laboratory medicine.

INDICATIONS FOR USE

Human ALK protein in formalin-fixed, paraffin embedded tissue

COUNTER-INDICATIONS

Not applicable.

There are no contra-indications to use the medical product for in vitro diagnostics.

EXPECTED AND FORECASTED SIDE EFFECTS ASSOCIATED WITH THE INTENDED USE OF THE MEDICAL PRODUCT.

Not applicable.

There are no side effects when using the medical product for in vitro diagnostics.

REAGENTS PROVIDED

The product is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin 950 as a preservative.

Total volume = 7 mL.

Clone

5A4

Immunogen

Prokaryotic recombinant protein corresponding to a region which spans the tyrosine kinase catalytic domain and part of the C-terminus of the NPM-ALK transcript (419-520 amino acids).

Specificity

Human ALK protein and C-terminal ALK fusion proteins containing the immunized target sequence.

Ig Class

IgG1

Total Protein Concentration

Approx. 10 mg/mL.

Antibody Concentration

Greater than or equal to 2.6 mg/L as determined by ELISA.

Dilution and Mixing

Anaplastic Lymphoma Kinase (5A4) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

MATERIALS OF HUMAN AND (OR) ANIMAL ORIGIN IN THE MEDICAL PRODUCT

There are no materials of human origin (cell culture products) in the medical product.

There are materials of animal origin (cell culture products) in the medical product:

- mouse anti-human monoclonal antibodies

Antibody concentration: ≥ 2.6 mg/L as determined by ELISA.

Reagents of animal origin (e.g., BSA and animal sera) are provided with research reports and certificates of origin confirming that they are virus-free and collected from healthy animals.

Animals used to produce monoclonal and polyclonal antibodies are of British origin and raised under highly regulated live-stock and veterinary regimes.

DELIVERY SET

BOND Primary monoclonal Antibody Anaplastic Lymphoma Kinase (5A4) for the qualitative identification by immunohistochemical method for in vitro diagnostics for BOND-III and Bond-maX immunohistochemistry instruments, composed of:

- BOND Primary Antibody Anaplastic Lymphoma Kinase (5A4), 7.0 ml, container - 1 pc.

- Instructions for Use (included) - 1 pc.

MATERIALS REQUIRED BUT NOT PROVIDED

Product name	Registration certificate No.
BOND-III IHC/ISH stainer with accessories	P3H 2016/3560 dd. 07.04.2021
Bond-maX IHC/ISH stainer with accessories	ΦC3 2009/04611 dd. 07.04.2021
Bond Polymer Refine Detection	ΦC3 2011/10032 dd. 09.11.2018
In vitro IHC reagents: Bond Dewax Solution Bond Wash Solution 10X Concentrate, 1 L Bond Epitope Retrieval 1 Bond Epitope Retrieval 2	ΦC3 2011/11352 dd. 25.07.2022
Dyes and ancillary reagents for specimen preparation and IVD staining of biospecimens for medical purposes	P3H 2013/1170 dd. 28.04.2017
Tissue processor ASP6025 S	P3H 2017/5542 dd. 05.10.2022
Rotary Microtome HistoCore BIOCUT	P3H 2020/12056 dd. 18.10.2022

Products that are not medical, but intended for use in combination with the declared medical product: BOND Cleaning system

USE, STORAGE AND TRANSPORTATION CONDITIONS, STABILITY

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Anaplastic Lymphoma Kinase (5A4) are: turbidity of the solution, odor development, and presence of precipitate. Return to 2–8 °C immediately after use.

Shelf life: 12 months from the date of production when stored at 2-8°C.

After opening, the product is ready to be used for the whole stated shelf life when stored at 2–8 °C.

Transportation can be carried out in any way subject to storage conditions. The products should be transported in disposable insulated containers with cold elements, or in reusable insulated containers with automatically maintained temperature or in refrigerated trucks with temperature indicators. The products transported with incorrect temperature regime are not subject to use.

SAMPLE REQUIREMENTS

Type of sample analyzed: sections (on slides) of formalin-fixed and paraffin-embedded (FFPE) human tissues.

To collect and transport the samples, it is recommended to use a fixative - 10% neutral-buffered formalin (for example, Dyes and ancillary reagents for specimen preparation and IVD staining of biospecimens for medical purposes, cl. 14. 10% NBF, Leica Biosystems Richmond, Inc., USA. RC No. P3H 2013/1170 dd. 28.04.2017) for formalin-fixed and paraffin-embedded patient tissues.

SPECIAL NOTES AND PRECAUTIONS

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth. Avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

INSTRUCTIONS FOR USE OF THE IVD MEDICAL DEVICE

PROCEDURE

Refer to BOND-III and Bond-maX Operator Manuals for more information on procedure.

Anaplastic Lymphoma Kinase (5A4) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection.

The recommended staining protocol for Anaplastic Lymphoma Kinase (5A4) is IHC Protocol F. Heat induced epitope retrieval is recommended using the product "IVD IHC reagents, cl. 25 BOND Epitope Retrieval Solution 2" for 30 minutes.

Results Expected

Normal Tissues

Clone 5A4 detected Anaplastic Lymphoma Kinase protein in some neuroendocrine cells of the pituitary gland and small intestine, alveolar macrophages in lung and a proportion of prostatic stromal cells. All other tissues tested were negative. (Total number of normal cases evaluated = 130).

Tumor Tissues

PA0831 stained 13/35 anaplastic large cell lymphomas and 2/5 large cell carcinoma of the lung. No staining was observed in a variety of other tumors, including lung tumors (0/68), breast tumors (0/5), thyroid tumors (0/5), tumors of the colon (0/4), liver tumors (0/4), metastatic tumors (0/4), tumors of the esophagus (0/3), stomach tumors (0/3), head and neck tumors (0/3), ovarian tumors (0/3), tumors of the adrenal gland (0/2), bladder tumors (0/2), brain tumors (0/2), tumors of the small intestine (0/2), rectal tumors (0/2), kidney tumors (0/2), tumors metastatic to the lymph node (0/2), salivary gland tumors (0/2), testicular tumors (0/2), cervical tumors (0/2) endometrial tumors (0/2), an osteosarcoma (0/1), a pancreatic adenocarcinoma (0/1), an adenocarcinoma of the prostate (0/1), a squamous cell carcinoma of the skin (0/1) and a melanoma (0/1). (Total number of abnormal tissues evaluated = 169).

The product is recommended to aid in the evaluation and expression of endogenous ALK and ALK fusion proteins, in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

PRODUCT SPECIFIC LIMITATIONS

Anaplastic Lymphoma Kinase (5A4) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

FUNCTIONAL SPECIFICATIONS

Analytical characteristics of the IVD medical product

Analytical sensitivity

Not applicable

Analytical specificity and interference

Not applicable.

Repeatability and reproducibility

Estimation of variability between sites, analyzers, users, batches

Repeatability – 100%

Reproducibility – 100%

Linearity

Not applicable.

Assay measurement range

Not applicable.

Threshold value

Not applicable.

Accuracy

Not applicable.

Diagnostic characteristics of the IVD medical product

Diagnostic sensitivity

Anaplastic large cell lymphoma and non-small cell lung cancer. Note: Expected prevalence of ALK and ALK-fusion proteins in anaplastic large cell lymphoma and non-small cell lung cancer is ~60% and ~3%, respectively.

An observed rate of 37% (13/35) of tested anaplastic large cell lymphoma cases stained positive, which reflects the positivity prevalence of the cases evaluated rather than true diagnostic sensitivity because the number of cases harboring ALK-translocation is unknown. All 35 tested anaplastic large cell lymphoma cases were verified with the Kreatech break-apart FISH probe (KI-10747 – ALK (2p23) Break); all 13 cases (100%) positive by ALK IHC (PA0831) were positive by FISH, and all 22 cases negative by ALK IHC (PA0831) were negative by FISH.

In non-small lung cancer cases, 2/72 (2/5 large cell carcinoma of the lung and 0/67 other non-small cell lung cancer cases) stained positive, giving an overall positivity prevalence of 3% (2/72).

Diagnostic sensitivity – 37%

Diagnostic specificity

Other neoplasms (excluding anaplastic large cell lymphoma and non-small cell lung cancer)

PA0831 was assessed in a range of clinical samples from abnormal (neoplastic) tissue specimens. The results obtained were as follows:

10

PA0831 stained 13/35 anaplastic large cell lymphomas and 2/5 large cell carcinoma of the lung. No staining was observed in a variety of other tumors, including lung tumors (0/68), breast tumors (0/5), thyroid tumors (0/5), tumors of the colon (0/4), liver tumors (0/4), metastatic tumors (0/4), tumors of the esophagus (0/3), stomach tumors (0/3), head and neck tumors (0/3), ovarian tumors (0/3), tumors of the adrenal gland (0/2), bladder tumors (0/2), brain tumors (0/2), tumors of the small intestine (0/2), rectal tumors (0/2), kidney tumors (0/2), tumors metastatic to the lymph node (0/2), salivary gland tumors (0/2), testicular tumors (0/2), cervical tumors (0/2) endometrial tumors (0/2), an osteosarcoma (0/1), a pancreatic adenocarcinoma (0/1), an adenocarcinoma of the prostate (0/1), a squamous cell carcinoma of the skin (0/1) and a melanoma (0/1). (Total number of abnormal tissues evaluated = 169).

An observed rate of 100% (62/62) of other tumor cases (excluding anaplastic large cell lymphoma, large cell carcinoma of the lung, and 67 lung tumors which were non-small cell lung cancer) stained negative, with ALK protein staining limited to only the anaplastic large cell lymphoma and large cell carcinoma of the lung cases.

Diagnostic specificity – 100%

DISPOSAL AND ENVIRONMENTAL PROTECTION

When using and storing, the BOND Primary monoclonal Antibody Anaplastic Lymphoma Kinase (5A4) for the qualitative identification by immunohistochemical method for in vitro diagnostics for BOND-III and Bond-maX immunohistochemistry instruments does not affect the patient, medical staff and environment if rules described in Instructions for Use are followed.

Waste disposal

Materials used in this analysis, including reagents, specimens and products used, should be disposed of in accordance with federal, state or local regulations.

Reagents used and unused, as well as any other contaminated disposable materials should be disposed of in accordance with rules for disposal of infected or potentially infected products.

The laboratory staff shall be responsible for handling and disposing of waste and liquid waste according to their type and hazard level, in accordance with all applicable rules and regulations. Disposal of the products shall be carried out by organizations, that have the appropriate license, at specially equipped sites, landfills and premises in accordance with the requirements stipulated by existing Federal laws and in compliance with mandatory environmental protection requirements.

Disposal shall be carried out on a general basis and in accordance with the requirements of SanPiN 2.1.3684-21.

Used products, unused products with expired shelf life, products with damaged packaging should be disposed of as medical waste of class B in accordance with the classification according to SanPiN 2.1.3684-21.

Methods of waste disposal

Recommendations:

Chemicals should be disposed of in accordance with current national regulations.

Packaging contaminated:

Recommendations:

Waste should be disposed of in accordance with local regulations. Packaging that cannot be cleaned must be disposed of in the same way as the product.

Detergents recommended:

Water, if necessary, in combination with non-oxidizing detergents. Do not use bleaching agents.

MANUFACTURER WARRANTY

The warranty is valid for the entire period until the expiration date (12 months), subject to compliance with the recommendations for transportation, storage and use specified on the external labeling of the medical device and in the Instructions for Use.

MANUFACTURER (PRODUCER)

Leica Biosystems Newcastle Limited

Balliol Business Park West, Benton Lane, Newcastle upon Tyne, NE12 8EW, United Kingdom

Tel: +44 191 215 0567

Fax: + 44 191 215 1152

E-mail: msds@leicabiosystems.com

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION:

Limited Liability Company "BioLine" (BioLine LLC)

Russia, 197022 Saint Petersburg, Municipal district Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, structure 3H/5H, room 208

Tel: +7 (812) 320 49 49

Fax: +7 (812) 320 49 40

E-mail: main@bioline.ru

CLAIM PROCEDURE

To receive technical advice and support, please contact an authorized representative of the manufacturer in the Russian Federation:

Limited Liability Company "BioLine" (BioLine LLC)

Russia, 197022 Saint Petersburg, Municipal district Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, structure 3H/5H, room 208

Tel: +7 (812) 320 49 49

Fax: +7 (812) 320 49 40

e-mail: main@bioline.ru

In case of side effects not specified in the Instructions for Use or operator manual for the medical product, undesirable reactions during its use, features of the interaction of medical products with each other, facts and circumstances that pose a threat to the life and health of citizens and medical specialists during the use and operation of medical products, it is necessary to send a message containing the specified information to the Federal Service for Surveillance in Healthcare in accordance with current legislation.

BIBLIOGRAPHY

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Hisaoka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumors: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721
5. Savic S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. Journal of Thoracic Oncology. 2013; 8(8):1004-11
6. Savic S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. Lung Cancer. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. Virchows Archiv: an international journal of pathology. 2016; 469(5):533-540
8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. Journal of Molecular Endocrinology. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. Cancer Biology and Therapy. 2015; 16(12):1691-701

CALIBRATION PRINCIPLE

Not applicable.

LIST AND DESCRIPTION OF THE IVD MEDICAL PRODUCT MATERIALS THAT COME IN CONTACT WITH THE PATIENT'S BODY (HUMAN BODY)

Not applicable. Components of this medical product do not come into contact with the patient's body (human body) either directly or indirectly.

INFORMATION ON MEDICINES FOR MEDICAL PURPOSE, PHARMACEUTICAL SUBSTANCES FOR MEDICAL PURPOSE USED FOR THE FINAL PRODUCT, INCLUDING PRINCIPLES AND METHODS OF EVALUATION

Not applicable. There are no medicines for medical purpose, pharmaceutical substances for medical purpose in this medical product.

METHOD OF FINAL STERILIZATION

Not applicable. This medical product is not sterile and not subject to be sterilized.

15

INFORMATION ON PREVIOUSLY REGISTERED VERSIONS OF THE MEDICAL PRODUCT

Other versions have not previously been registered in Russia.

RISKS OF USING THIS MEDICAL PRODUCT

Risk management is carried out in accordance with the requirements of harmonized standard EN ISO 14971, "Medical devices – Application of risk management to medical devices".

After conducting a risk analysis and obtaining risk reduction indicators in accordance with the Quality Management System, corrective and preventive action plans were reviewed. Residual risks are considered to be acceptable.

SYMBOLS ON THE MARKING

List of symbols used on the marking:

Symbol	Meaning
	Catalogue number
	Medical device for <i>in vitro</i> diagnostics.
	Manufacturer
	Temperature range
	Best before (year-month-day)
	Lot number
	Refer to the Instructions for Use
	EU conformity mark
	Caution!
	Recycling symbol

LIST OF NATIONAL AND INTERNATIONAL REGULATIONS/STANDARDS WITH WHICH THE IVD MEDICAL PRODUCT COMPLIES

No.	DOCUMENT TITLE
1	EN ISO 13485 - Medical devices. Quality Management Systems. Requirements for regulatory purposes

2	EN ISO 13612 - Performance evaluation of in vitro diagnostic medical devices
3	EN ISO 23640 - In vitro medical devices. Evaluation of stability of in vitro diagnostic reagents
4	EN ISO 18113-1 - In vitro medical devices. Information supplied by the manufacturer (labelling). Part 1. Terms, definitions and general requirements
5	EN ISO 18113-2 - In vitro medical devices. Information supplied by the manufacturer (labelling). Part 2. In vitro diagnostic reagents for professional use
6	EN ISO 14971 - Medical devices. Application of risk management to medical devices
7	EN ISO 9001 - Quality management systems. Requirements
8	EN ISO 15223-1 - Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements
10	ISO 14644-5 - Cleanrooms and associated controlled environments. Part 5. Operations
11	GOST R 51088 - In vitro diagnostic medical devices. Reagents, kits, the test - systems, control materials, culture media. Requirements to devices, and to supporting
12	GOST R 51352 - In vitro diagnostic medical devices. Test methods

INFORMATION ON INITIAL ISSUE OR LATEST REVISION OF THE OPERATIONAL DOCUMENTATION

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