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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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
THIS IS TO CERTIFY THAT: the seal of Guilin Hbm Health Protections, Inc. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明 网网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

恒保防護
HBM Protections

Guilin HBM Health Protections, Inc.

No. 1-2, Shuijing East Road, Economic and Technological
Development Area, 541805 Guilin, Guangxi, China

Tel.: 0773-8569067

Email: guilinhbm@hbmchina.com

APPROVED BY

Quality Director

Mr. Pu Lei

INSTRUCTION FOR USE
of medical device

Surgical disposable sterile natural latex gloves

Revision: 00 (April 2024)

1. NAME

Surgical disposable sterile natural latex gloves

consisting of:

I. Design versions:

1. Surgical disposable sterile natural latex gloves textured powdered (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
2. Surgical disposable sterile natural latex gloves textured powder-free with inner gel-based polymer coating, with humidifier (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
3. Surgical disposable sterile natural latex gloves textured powder-free with inner gel-based polymer coating, antimicrobial (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
4. Surgical disposable sterile natural latex gloves textured powder-free with inner gel-based polymer coating, elongated cuff 350 mm (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
5. Surgical disposable sterile natural latex gloves textured powder-free specialized with inner gel-based polymer coating, for gynecology elongated cuff 400 mm (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
6. Surgical disposable sterile natural latex gloves textured powder-free specialized with inner gel-based polymer coating, for gynecology elongated cuff 480 mm (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
7. Surgical disposable sterile natural latex gloves textured powder-free specialized with inner gel-based polymer coating, for microsurgery (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: brown, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
8. Surgical disposable sterile natural latex gloves textured powder-free specialized with inner gel-based polymer coating, for orthopedics (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: brown, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
9. Surgical disposable sterile natural latex gloves textured powder-free specialized with inner gel-based polymer coating, with puncture indication (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: inner glove – green, outer glove – white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
10. Surgical disposable sterile natural latex gloves textured powder-free specialized with inner gel-based polymer coating, when working in a puncture indication system (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: green, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.
11. Surgical disposable sterile natural latex gloves textured powder-free with gel-based inner polymer coating (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.

12. Surgical disposable sterile natural latex gloves textured powder-free with inner gel-based polymer coating, with low protein content (beaded cuff/cuff without beading/cuff without beading with adhesive strip), color: white, size: 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0 – 10/25/50/100/200 pairs/pack.

II. Instructions for use – 1 pc (on the consumer package).

Manufacturer and producer of the medical device

Manufacturer	Guilin HBM Health Protections, Inc. No. 1-2, Shuijing East Road, Economic and Technological Development Area, 541805 Guilin, Guangxi, China
Address of the manufacturing facility	No. 1-2, Shuijing East Road, Economic and Technological Development Area, 541805 Guilin, Guangxi, China

Authorized representative of the manufacturer

SKANMASTERS Company Limited

SKANMASTERS Co. Ltd.

Russia, 198261, Saint-Petersburg, ext. ter. c. municipal district Ulyanka, Stoykosti st., 38,
bldg. 1, litera A, office 5-N, room 2.

Tel.: +7 (911) 251-54-76, e-mail: masalex7@yandex.ru

The scopes of delivery of the medical device:	- 10 pairs in a consumer package; - 25 pairs in a consumer package; - 50 pairs in a consumer package; - 100 pairs in a consumer package; - 200 pairs in a consumer package.
--	---

2. PURPOSE

The gloves are intended for single use in protecting medical personnel's hands and creating a barrier against bacterial and/or viral contamination between the health worker and the patient during surgical interventions and other sanitary actions that require compliance with sterility conditions.

The gloves can be used by medical personnel involved in surgical intervention, as well as performing any other sanitary actions that require compliance with aseptic requirements.

Indications for use of the medical device is the necessity for surgical interventions and other sanitary actions that require compliance with sterility conditions. Gloves represent the personal protective equipment for medical personnel's hands. The use of gloves in healthcare organizations is a mandatory requirement.

3. FUNCTIONAL SPECIFICATION

The operational principle of the medical device is to create a physical barrier between the skin surfaces of the medical personnel's hands and the surgical field or the other surfaces of the patient, that require compliance with sterility conditions. Thus, the medical device prevents bacterial and/or viral contamination between the medical personnel and the patient and increases the effectiveness of aseptic conditions during surgical interventions, protecting the health worker from possible infection during surgery and other sanitary actions.

Medical device design versions allow covering almost all areas of general and specialized surgery, including microsurgery, orthopedics, traumatology, gynecology and operations with

an increased risk of infection.

The microtextured surface of all design versions of the medical device contributes to the reliable grip of medical instruments and thus increases the probability of a successful surgical outcome.

4. RISK OF USING, CONTRAINDICATIONS AND SIDE EFFECTS

Potential risk of using of the gloves – Class 2a.

Contraindications of the medical device use is a presence of allergic reaction to powdered gloves and/or to natural latex gloves in a patient or a health worker.

Healthcare providers may observe the following adverse effects of using gloves: non-allergic contact dermatitis, allergic contact dermatitis, and immediate allergic reactions accompanied by hives, pruritus and systemic allergic reactions. Allergic reactions when using gloves are most often associated with a latex allergy, individual sensitization of a medical professional to the chemical components of latex material. The use of powder to facilitate donning of gloves also increases the risk of allergies in health workers and may also cause a number of postoperative complications in patients.

The manifestation of dermatitis and other skin irritations is possible as a result of improper hands skincare and harsh conditions of hand skin isolation from natural air ventilation.

If skin irritation occurs, stop using gloves and consult a doctor.

Shelf life of surgical gloves is 3 years.

5. TECHNICAL SPECIFICATION

Form Surgical gloves have curved fingers, with division into the right and left hands (anatomical), the cuff can be beaded or without beading, or without beading with adhesive strip.

Surface The gloves have a textured surface over the entire area of fingers and palm.

Dimensions Width, length and thickness of gloves, depending on their design versions and size, must comply with the values given in table 1-5.

The thickness of the bead for all design versions of the medical device is not more than 2.5 mm.

Table 1. Dimensions of the medical device for design versions Nos. 1-3, 7, 9-12.

Size	Length <i>l</i> , mm, min	Width <i>W</i> , mm	Thickness, mm, min		
			finger	palm	cuff
6.0	260	77 ± 5	0.13	0.13	0.10
6.5		83 ± 5			
7.0		89 ± 5			
7.5	270	95 ± 5			
8.0		102 ± 6			
8.5		108 ± 6			
9.0	280	114 ± 6			

Table 2. Dimensions of the medical device for design version No. 4.

Size	Length <i>l</i> , mm, min	Width <i>W</i> , mm	Thickness, mm, min		
			finger	palm	cuff
6.0	350	77 ± 5	0.13	0.13	0.10
6.5		83 ± 5			

7.0		89 ± 5			
7.5		95 ± 5			
8.0		102 ± 6			
8.5		108 ± 6			
9.0		114 ± 6			

Table 3. Dimensions of the medical device for design version No. 5.

Size	Length <i>l</i> , mm, min	Width <i>W</i> , mm	Thickness, mm, min		
			finger	palm	cuff
6.0	400	77 ± 5	0.13	0.13	0.10
6.5		83 ± 5			
7.0		89 ± 5			
7.5		95 ± 5			
8.0		102 ± 6			
8.5		108 ± 6			
9.0		114 ± 6			

Table 4. Dimensions of the medical device for design version No. 6.

Size	Length <i>l</i> , mm, min	Width <i>W</i> , mm	Thickness, mm, min		
			finger	palm	cuff
6.0	480	77 ± 5	0.13	0.13	0.10
6.5		83 ± 5			
7.0		89 ± 5			
7.5		95 ± 5			
8.0		102 ± 6			
8.5		108 ± 6			
9.0		114 ± 6			

Table 5. Dimensions of the medical device for design version No. 8.

Size	Length <i>l</i> , mm, min	Width <i>W</i> , mm	Thickness, mm, min		
			finger	palm	cuff
6.0	260	77 ± 5	0.20	0.20	0.15
6.5		83 ± 5			
7.0		89 ± 5			
7.5	270	95 ± 5			
8.0		102 ± 6			
8.5	280	108 ± 6			
9.0		114 ± 6			

Physico-mechanical
characteristics

Physico-mechanical characteristics of the gloves must comply with
the values given in table 6.

Table 6.

Characteristics	Before accelerated aging	After accelerated aging
Tensile strength, MPa, min	24	18
Elongation at break, %, min	750	560
Stress at 500% elongation, MPa, max	5.5	—
Force at break, N, min	12.5	9.5
Force at 300% elongation, N, max	2.0	—

Inspection levels and acceptable quality level

For the purposes of quality control, gloves are selected and inspected according to the requirements of ISO 2859-1 standard. Inspection levels and acceptable quality levels AQL must comply with the values given in table 7.

Table 7.

Parameter	Defect related to parameter	Inspection level	AQL
Watertightness	Presence of holes	G1	1.5
Dimensions	Mismatch of length, width and thickness	S-2	4.0
Physico-mechanical parameters	Mismatch of parameters before and after accelerated aging	S-2	4.0
Residual powder concentration	Exceeds recommended maximum value	N=5	—
Water-extractable protein content	Exceeds recommended maximum value	N=3	—
Antigenic protein content	Exceeds recommended maximum value	N=1	—

Watertightness	Gloves are watertight. The acceptable quality level (AQL) is ≤ 1.5 .
Powder content	In powder-free gloves the residual powder content does not exceed 2 mg per one glove. In powdered gloves, the powder content does not exceed 10 mg/dm ² .
Water-extractable protein content	The concentration of water-extractable proteins is not more than 200 µg/dm ² .
Antigenic protein content	The concentration of antigenic proteins is not more than 10 µg/dm ² .
Resistance to viral penetration	The medical device is resistant to viral penetration.

6. APPLICATION PROCEDURE

Before using the gloves health workers carry out hygienic treatment of their hands having performed surgeon's hand treatment technology.

It is necessary to strictly follow the rules for donning gloves so as not to violate the sterility of their work surface.

Open the individual outer package of sterile gloves (away from the sterile table), remove ("shake out") from it the inner wallet with gloves on a sterile surface. Then open the inner wallet with sterile forceps or aseptically treated hands. Initially, the glove is donned on the dominant (right) hand. Grab the edge of the right glove cuff turned inside out with the thumb and forefinger of the left hand and put it on the right hand without touching the outer sterile surface of the glove. The cuff remains inverted. Next, place fingers of the right hand, on which the glove is already worn, under the cuff lapel of the left glove and don it on the left hand, without touching the skin and the wrong side of the cuff. Spread the cuff of the left glove on the gown in circular motions. The cuff of the glove should overlap the cuff of the gown sleeve by more than 5-10 cm. Only then spread the inverted cuff edge of the right glove in the same circular movements. Pull on the cuff with your four fingers and finally adjust the gloves, ensuring their comfortable use.

Visually the process of surgical gloves donning is shown in Fig.1.


Figure 1. Surgical gloves donning algorithm

There are features of work in surgical gloves.

Hands with sterile gloves should not be lowered below the sterile surfaces of the instrument and operating table. It is not recommended to treat gloves with alcohol-containing and other antiseptic agents before the operation and during surgery, as they increase the porosity/permeability of the glove material and create conditions for overcoming the protective barrier by microorganisms. To maintain an optimal level of protection, the gloves must be changed every 60-120 min, depending on the nature of the surgical intervention. Before donning new gloves, it is necessary to treat the hands with an alcohol-containing antiseptic.

Gloves should be replaced in the following situations:

- after damage (puncture);
- after defect detection;
- after an accidental electric shock from an electrosurgical instrument;
- if any liquid gets (absorbs) under the glove;
- upon appearance of stickiness;
- upon sensation of "glove juice" appearance;
- in the transition from a "dirty" surgical stage to a "clean".

If one glove is damaged, both gloves must be changed. When changing gloves after removing them from the hands, the hands should be treated with a skin antiseptic. Don a new pair of gloves on dry hands.

Remove used gloves carefully, avoiding splashes from the surfaces of the glove, which could cause microbial contamination of the hands and the environment.

To remove the glove, it is necessary to take the left glove cuff from the outside by fingers of the right hand in the glove. Remove the glove from the left hand turning it inside out and leave it holding the lapel in the right hand. Then place the fingers of the free left hand in the glove cuff on the right hand and pull it off similarly to the first one. Having turned the glove inside out, use it as a bag for another glove. Immerse both gloves in a container with a disinfectant solution for disinfection or in a disposable container for temporary waste storage for subsequent centralized disinfection/neutralization.

Visually the process of surgical gloves removing is shown in Fig.2.


Figure 2. Surgical gloves removal algorithm

After finishing the work and gloves removal, it is necessary to carry out hygienic hand treatment. Disposable gloves after use must be decontaminated/neutralized as a medical waste of the appropriate class. To prevent dryness and adverse effects on the hand skin, it is recommended for health workers to always use professional nourishing hand cream.

7. WARNINGS AND PRECAUTIONS

- The device is intended for single use only.
- For individual use only.
- Every time it is necessary to wash hands before donning or changing gloves, as well as after removing them.
- It is necessary to choose gloves of the right size and depending on the nature of the performed manipulations – they should be easily donned and should not hamper movements or squeeze the hands.
- Touching is permissible only on the limited area of the glove corresponding to the wrist (upper edge of the cuff).
- Sterile. Do not use if individual package is damaged.
- Sterilized by radiation method.
- If the integrity of the individual package is broken, the gloves must be disposed of as General Wastes (Class A Wastes according to the regulatory guidelines of the Russian Federation SanPiN 2.1.3684-21).
- Do not wet and protect from moisture.
- Avoid exposure to direct sunlight or bright artificial light that has a high UV content.
- Store in a place protected from sunlight. Storage rooms should not contain any equipment capable of producing ozone, such as mercury lamps or high voltage equipment, causing electric sparks or inaudible electric discharges. Combustible gases and organic vapors should also be excluded from the premises.
- The product is made of natural rubber latex, which may cause an allergic reaction.
- Powder may cause allergic reactions, as well as inflammation and the formation of adhesions in the patient's surgical wound.
- The gloves do not guarantee 100% hand protection against punctures, cuts, and damage caused by exposure to chemicals.
- If skin irritation occurs, stop using the gloves and consult a doctor.

8. APPLICATION CONDITIONS

Gloves must be used as per the instruction on the medical device. The use of the device does not require special qualification or special education.

Medical gloves are used in medical or healthcare facilities indoors, in a dry place, away from direct sunlight, at a temperature of 15°C to 25°C and relative humidity from 65% up to 85%.

For single use only. For individual use only. Do not use in case of size mismatch.

9. TRANSPORTATION, STORAGE AND DISPOSAL CONDITIONS OF THE MEDICAL DEVICE

The transportation of the device may be carried out by all types of covered vehicles in the original package as per the transport regulations, which are applied to this type of transport. During transportation the storage conditions must not be violated.

The devices must be stored in a dry place, protected from direct sunlight or bright artificial light having a high ultraviolet content, at a temperature of 0°C to 25°C and relative air humidity – from 65 % to 85%.

The devices should be handled with caution when they are stored at temperature below 15°C, as during careless handling, they can be stiff and sensitive to deformation. After storage at a low temperature, gloves should be warmed up to 30°C before use throughout all device volume.

Storage rooms should not contain ozone-generating equipment, such as mercury lamps or high voltage electrical equipment, generating excessive sparks or inaudible electric discharges.

The devices must be stored away from sources of ignition. Combustible gases and organic vapors should be excluded from the storage facilities.

An unused medical device that has become useless, including due to the expiration date, belongs to the category of General waste (Class A waste according to the regulatory guidance of the Russian Federation SanPiN 2.1.3684-21).

The used medical device belongs to the category of Infectious waste (Class B waste according to the regulatory guidance of the Russian Federation SanPiN 2.1.3684-21).

Collection, temporary storage, and waste disposal should be carried out following the medical waste management scheme adopted in the given organization carrying out medical and/or pharmaceutical activities.

10. MANUFACTURER'S WARRANTIES

The manufacturer guarantees the safety of all parameters and characteristics of the surgical gloves for 3 years from the manufacturing date upon the condition of package integrity, subject to the conditions of transportation, storage and application set forth in the instruction.

11. LIST OF NATIONAL STANDARDS, APPLIED BY MANUFACTURER OF THE MEDICAL DEVICE

EN 455-1:2020	Medical gloves for single use. Part 1. Requirements and testing for freedom from holes
EN 455-2:2015	Medical gloves for single use. Part 2. Requirements and testing for physical properties
EN 455-3:2015	Medical gloves for single use. Part 3. Requirements and testing for biological evaluation
EN 455-4:2009	Medical gloves for single use. Part 4. Requirements and testing for shelf life determination

EN 556-1:2002	Sterilization of medical devices. Requirements for medical devices to be designated "sterile". Part 1. Requirements for terminally sterilized medical devices
ISO 2859-1:1999	Sampling procedures for inspection by attributes – Part 1: Sampling schemes indexed by acceptance quality limit (AQL) for lot-by-lot inspection
ISO 11137-1:2006	Sterilization of health-care products – Radiation – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
ISO 11137-2:2013	Sterilization of health-care products – Radiation – Part 2: Establishing the sterilization dose
ISO 11137-3:2017	Sterilization of health-care products – Radiation – Part 3: Guidance on dosimetric aspects of development, validation and routine control
ISO 10993-1:2018	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process
ISO 10993-10:2010	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization
ISO 10993-12:2012	Medical devices. Biological evaluation of medical devices – Part 12: Sample preparation and reference materials
ISO 11607-1:2019	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems
ISO 11607-2:2019	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
ISO 11737-1:2018	Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms on products
ISO 11737-2:2019	Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
ISO 13485:2016	Medical devices – Quality management systems – Requirements for regulatory purposes
ISO 14971:2007	Medical devices – Application of risk management to medical devices
ISO 15223-1:2016	Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements
ASTM D412	Standard Test Methods for Vulcanized Rubber and Thermoplastic Elastomers – Tension
ASTM D573	Standard Test Method for Rubber – Deterioration in an Air Oven
ASTM D3767	Standard Practice for Rubber – Measurement of Dimensions
ASTM D5151	Standard Test Method for Detection of Holes in Medical Gloves
ASTM D5712	Standard Test Method for Analysis of Aqueous Extractable Protein in Latex, Natural Rubber, and Elastomeric Products Using the Modified Lowry Method
ASTM D6124	Standard Test Method for Residual Powder on Medical Gloves
ASTM D6499	Standard Test Method for Immunological Measurement of Antigenic Protein in Hevea Natural Rubber (HNR) and its Products
ASTM F1671	Standard test method for resistance of materials used in protective clothing to penetration by blood-borne pathogens using Phi-X174 bacteriophage penetration as a test system
CE	European Conformity. Requirements

Instruction for use of the medical device is released for the first time (Rev: 00 (April 2024)).

In case of an undesirable event showing the signs of an adverse event (incident), it is necessary to notify thereof the manufacturer or an authorized representative of the manufacturer.

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